


General Disclaimer

One or more of the Following Statements may affect this Document

- This document has been reproduced from the best copy furnished by the organizational source. It is being released in the interest of making available as much information as possible.
- This document may contain data, which exceeds the sheet parameters. It was furnished in this condition by the organizational source and is the best copy available.
- This document may contain tone-on-tone or color graphs, charts and/or pictures, which have been reproduced in black and white.
- This document is paginated as submitted by the original source.
- Portions of this document are not fully legible due to the historical nature of some of the material. However, it is the best reproduction available from the original submission.

NASA CR 167737



SPECIAL REPORT ON THE
DATA COLLECTION PROGRAMS FOR THE
GROUND BASED NITROGEN WASHOUT EXPERIMENT
VOLUME TWO: DETAILED PROGRAM DESCRIPTIONS, LISTINGS,
EXAMPLES AND HARDWARE SPECIFICATIONS

Prepared for the NASA Johnson Space Center
Life Sciences Medical Directorate
by Penelope A. Bueker

July 30, 1982

Contract NAS9-14880
Project 0185-40

(NASA-CR-167737) SPECIAL REPORT ON THE DATA
COLLECTION PROGRAMS FOR THE GROUND BASED
NITROGEN WASHOUT EXPERIMENT. VOLUME 2:
DETAILED PROGRAM DESCRIPTIONS, LISTINGS,
EXAMPLES AND HARDWARE (Technology, Inc.,

N83-10760

HCAID

Unclas
G3/52 35790



**Technology
Incorporated**
Life Sciences Division



P.O. Box 58827
Houston, Texas 77058

SPECIAL REPORT
ON THE DATA COLLECTION PROGRAMS
FOR THE GROUND BASED NITROGEN WASHOUT EXPERIMENT

VOLUME TWO: DETAILED PROGRAM DESCRIPTIONS, LISTINGS, EXAMPLES
AND HARDWARE SPECIFICATIONS

The Nitrogen Washout System measures nitrogen elimination on a breath basis from the body tissues of a subject breathing pure oxygen. The system serves as a prototype for a Space Shuttle Life Sciences experiment and in the Environmental Physiology Laboratory. Typically, a subject washes out body nitrogen for three hours while breathing oxygen from a mask enclosed in a positive pressure-oxygen tent. A nitrogen washout requires one test operator and the test subject.

A DEC LSI-11/02 computer is used to (1) control and calibrate the mass spectrometer and Skylab spirometer, (2) gather and store experimental data (3) and provide limited real time analysis and more extensive post-experiment analysis. Five programs are used to gather and store the experimental data and perform all the real time control and analysis.

The Personal Data Input Program collects the data on a subject's physical characteristics and body measurements as entered by the test operator. The body composition and body nitrogen content are calculated by the program for one or more processes: (1) Water Immersion, (2) Skinfold Process I, (3) Skinfold Process II and (4) Nutrition Journal Formula. Body composition and body nitrogen content are also calculated from the subject's known per cent body fat, if available. The program produces a subject characteristics profile on the printing terminal and stores the same information in a data file.

The Decompression Data Program collects data from a decompression experience experiment, then computes the final tissue nitrogen and the decompression ratio at the time of final decompression. The data is written to a data file suitable for analysis. It types a short report on final tissue nitrogen partial pressure and decompression ratio at time of final decompression for specified tissue saturation half times.

The Nitrogen Washout Program controls the gas valving system for the mass spectrometer, calibrates the mass spectrometer and spirometer, and records all the subject and test data entered by the test operator. It samples the spirometer and the gas channels from the mass spectrometer, calculates inspired/expired air gas concentrations and breath volumes, then writes the data to a floppy disk for storage. It also produces a typed report of the nitrogen washout parameters along with one minute summaries as the washout proceeds.

The Nitrogen Data Program reads the information from the nitrogen washout and produces a typed report of one of four available reports to a printing terminal. The Update Program takes a previously collected, unformatted data file generated by the old nitrogen washout program and rewrites it to a formatted data file for use with N2DA TA report writing program or for transfer to the VAX.

DATE: 01/10/1973
OF POOR QUALITY

APPROVAL SHEET FOR THE
DATA COLLECTION PROGRAMS FOR THE
GROUND BASED NITROGEN WASHOUT EXPERIMENTS

VOLUME TWO: DETAILED PROGRAM DESCRIPTIONS, LISTINGS, EXAMPLES
AND HARDWARE SPECIFICATIONS

Approved by:

Edward C. Moseley

Edward C. Moseley, Ph.D.

Medical Sciences Division

NASA/Johnson Space Center

Joseph T. Baker

Joseph T. Baker

Supervisor, Laboratory Research

Support Section

Technology Incorporated

Life Sciences Division

Houston, Texas

Harry F. Walbrecher

Harry F. Walbrecher

Project Manager

Technology Incorporated

Life Sciences Division

Houston, Texas

Prepared by:

STATEMENTS
OF POOR QUALITY




Craig E. Litton

Task Manager, Biomathematics Services

Technology Incorporated

Life Sciences Division

Houston, Texas



Penelope Bueker

Junior Programmer

Biomathematics Services

Technology Incorporated

Life Sciences Division

Houston, Texas

PERSONNEL

Craig E. Litton - Project Design and Development

Penelope Bueker - Programming

Barry Levitan - Engineering, Programming

Linda Weaver - Graphics

Scott Thompson - Technical Assistance

Lita Holt - Typing

TABLE OF CONTENTS

	ABSTRACT	i
	APPROVAL SHEET	ii
	PERSONNEL	iv
	TABLE OF CONTENTS	v
	LIST OF FIGURES	vii
1.0	PROGRAM DESIGN AND CONCEPTS	1
1.1	PERSONAL DATA INPUT PROGRAM	1
1.1.1	PURPOSE	1
1.1.2	INPUT FILES AND DATA	1
1.1.3	OUTPUT FILES AND REPORTS	1
1.1.4	CALCULATIONS	1
1.2	DECOMPRESSION DATA PROGRAM	4
1.2.1	PURPOSE	4
1.2.2	INPUT FILES AND DATA	4
1.2.3	OUTPUT FILES AND REPORTS	4
1.2.4	CALCULATIONS	4
1.3	NITROGEN WASHOUT PROGRAM	5
1.3.1	PURPOSE	5
1.3.2	INPUT FILES AND DATA	5
1.3.3	OUTPUT FILES AND REPORTS	5
1.3.4	CALCULATIONS	6
1.4	NITROGEN DATA PROGRAM	8
1.4.1	PURPOSE	8
1.4.2	INPUT FILES AND DATA	8
1.4.3	OUTPUT FILES AND REPORTS	8
1.4.4	CALCULATIONS	9
1.5	UPDATE PROGRAM	11
1.5.1	PURPOSE	11
1.5.2	INPUT FILES AND DATA	11
1.5.3	OUTPUT FILES AND REPORTS	11
2.0	INPUT DATA AND FORMATS	12
2.1	PERSONAL DATA INPUT PROGRAM	12
2.2	DECOMPRESSION DATA PROGRAM	14
2.3	NITROGEN WASHOUT PROGRAM	16
2.4	NITROGEN DATA PROGRAM	17
2.5	UPDATE PROGRAM	17
3.0	PROGRAM OUTPUT, REPORTS AND DATA	18
3.1	PERSONAL DATA INPUT PROGRAM	18
3.2	DECOMPRESSION DATA PROGRAM	23
3.3	NITROGEN WASHOUT PROGRAM	28
3.4	NITROGEN DATA PROGRAM	30
3.5	UPDATE PROGRAM	32
4.0	PROGRAM FLOWCHARTS	34
4.1	PERSONAL DATA INPUT PROGRAM	34
4.2	DECOMPRESSION DATA PROGRAM	41
4.3	NITROGEN WASHOUT PROGRAM	49
4.4	NITROGEN DATA PROGRAM	64
4.5	UPDATE PROGRAM	69

5.0	PROGRAM LISTINGS	70
5.1	PERSONAL DATA INPUT PROGRAM	70
5.2	DECOMPRESSION DATA PROGRAM	90
5.3	NITROGEN WASHOUT PROGRAM	104
5.4	NITROGEN DATA PROGRAM	125
5.5	UPDATE PROGRAM	137
6.0	SAMPLE RUNS WITH INPUT AND OUTPUT PAGES	142
6.1	PERSONAL DATA INPUT PROGRAM	142
6.2	DECOMPRESSION DATA PROGRAM	148
6.3	NITROGEN WASHOUT PROGRAM	151
6.4	NITROGEN DATA PROGRAM	157
6.5	UPDATE PROGRAM	164
7.0	HARDWARE OPERATION	165
7.1	OVERVIEW	165
7.2	MASS SPECTROMETER	167
7.3	CALIBRATION VALUES	168
7.4	OPERATING THE N2WASH PROGRAM	170
8.0	ENGINEERING DIAGRAMS AND TABLES	172
9.0	APPENDIX: RT-11 COMMANDS	193
9.1	SPECIAL FUNCTION KEYS	193
9.2	USING THE SYSTEM UTILITIES FOR FILE MAINTENANCE.	194
9.3	GLOSSARY	200

LIST OF FIGURES

NITROGEN WASHOUT SYSTEM	173
BLOCK DIAGRAM OF N2WASHOUT HARDWARE	174
BOARD CONFIGURATION OF LSI-1102 MICROCOMPUTER	175
N2WASHOUT ELECTRONICS CARD CONNECTOR DESIGNATIONS	176
A/D CONVERTER	179
DRV/11 PARALLEL BOARD	180
FACE MASK AND GAS PLUMBING	182
POWER SUPPLIES FOR N2WASHOUT SYSTEM ELECTRONICS	183
CALIBRATION GAS SELECT CIRCUIT	184
MASS SPECTROMETER INLET SELECT CIRCUIT	185
MODIFICATIONS TO MODEL 1100 MASS SPECTROMETER	186
SPIROMETER LED INDICATORS	187
SPIROMETER CALIBRATION SWITCHES	188
SPIROMETER DUMP VALVE ELECTRONICS	189
SPIROMETER ELECTRONICS CHASSIS WIRING	190
SPIROMETER PRIM/FAIL SAFE BOARD	191
SPIROMETER POWER SUPPLY AND VOLUME READOUT	192

1.0 PROGRAM DESIGN AND CONCEPTS

1.1 PERSONAL DATA INPUT PROGRAM

1.1.1 PURPOSE

The PDINPT program is designed to collect descriptive information on the subject's physical characteristics and data from one or more of three skin measurement processes for determining body composition. This data is then used to compute nitrogen content information by the water immersion process, by each of the two skinfold processes used, by a nutrition journal formula and by a known per-cent body fat from the personal characteristics data.

1.1.2 INPUT FILES AND DATA

The data is entered by the test operator in response to computer prompts. Refer to section 2.1 INPUT DATA AND FORMATS.

1.1.3 OUTPUT FILES AND REPORTS

The program stores the data in a file on a floppy disk in device DY1: and produces a typed report on the printing terminal. Refer to section 3.1 PROGRAM OUTPUT, REPORTS AND DATA.

1.1.4 CALCULATIONS

1. Body Composition by Water Immersion

- a. Average under water weight from $(1 \leq N \leq 10)$ (kg)
- b. Water Weight = Average under water weight - under water chair weight (kg)

$$c. \text{ Body volume} = \frac{\text{Nude Weight} - \text{Water Weight}}{\text{Water Density}} - \text{Residual Vol (L)}$$

$$d. \text{ Body Density} = \frac{\text{Nude Weight}}{\text{Body Volume}}$$

$$e. \% \text{ Body Fat} = \left(\frac{4.95}{\text{Density}} - 4.50 \right) \times 100$$

2. Body Composition by Skinfold Process I

- a. Biceps skinfold thickness average from three trials (mm)
- b. Triceps skinfold thickness average from three trials (mm)
- c. Subscapular skinfold thickness average from three trials (mm)
- d. Suprailiac skinfold thickness average from three trials (mm)
- e. $S = \text{Log of the sum of Biceps average, Triceps average, Subscapular average, and Suprailiac average}$
- f. Body Density
 - i. If Male = $1.1610 - (.0632 \times S)$
 - ii. If Female = $1.1581 - (.0720 \times S)$
- g. $\% \text{ Body Fat} = \left(\frac{4.95}{\text{Density}} - 4.50 \right) \times 100$

3. Body Composition by Skinfold Process II

- a. Midclavicular skinfold thickness average from three trials (mm)
- b. Umbilicus skinfold thickness average from three trials (mm)
- c. Anterior mid thigh skinfold thickness average from three trials (mm)
- d. $S = \text{skinfold thickness sum of Midclavicular average, Umbilicus average, and Anterior mid thigh average (mm)}$
- e. Body Density
 - i. If Male = $1.10938 - (8.267 \times 10^{-4}) \times S + (1.6 \times 10^{-6}) \times (S)^2 - (2.574 \times 10^{-4}) \times \text{Age}$
 - ii. If Female = $1.105339 - (1.1964 \times 10^{-3}) \times S + (3.8 \times 10^{-6}) \times (S)^2 - (1.069 \times 10^{-4}) \times \text{Age}$
- f. $\% \text{ Body Fat} = \left(\frac{4.95}{\text{Density}} - 4.50 \right) \times 100$

4. Formula estimation of Body Characteristics from American Journal of Clinical Nutrition

a. Volume = 1.015 x Nude Weight - 4.937 (L)

b. Density = $\frac{\text{Nude Weight}}{\text{Volume}}$

c. % Body Fat = $52.67 - \frac{2454.38}{\text{Weight}}$

d. % Body Fat Free Mass = 100 - % Body Fat

e. % Total Body Water = $34.93 + \frac{1811.33}{\text{Weight}}$

5. Estimated Total Nitrogen Content of Body Before Nitrogen Washout

a. Fat Weight = % Body Fat x Nude Weight (kg)

b. Fat Free Mass = Nude Weight - Fat Weight (kg)

c. Essential Lean Body Water Mass = Fat Free Mass x .72 (kg)

d. Volume of Total Fat Tissue = $\frac{\text{Fat Weight (kg)}}{.90 \text{ kg/Liter at } 37 \text{ deg C}}$ (L)

e. Volume of Essential Lean Body Water =

$$\frac{\text{Essential Lean Body Water Mass (kg)}}{.9933 \text{ kg/Liter at } 37 \text{ deg C}} \text{ (L)}$$

f. Quantity of N2 in Total Fat Tissue = Fat Volume x .067 (L)

g. Quantity of N2 in Essential Lean Body Water =

$$\text{Lean Body Water Volume} \times .013 \text{ (L)}$$

h. Total Body N2 = Fat N2 + Lean Body Water N2 (L)

i. Nitrogen/kg Body = $\frac{\text{Total Body N2}}{\text{Nude Weight}}$ (L)

1.2 DECOMPRESSION DATA PROGRAM

1.2.1 PURPOSE

The DECDA T program records the data from a decompression experience in the chamber and types a report on the final tissue partial pressure of nitrogen in psi for each of three types of body tissue and the decompression ratio for each type of tissue at the time of final decompression. Tissues of half saturation time 180 minutes, 240 minutes and 360 minutes are used.

1.2.2 INPUT FILES AND DATA

The test operator enters the data on the subject's identification, the nitrogen partial pressure and duration of washout for each phase of a nitrogen washout, the experimental conditions during the chamber run, and the subject's symptoms during decompression in the hypobaric chamber. Refer to section 2.2 INPUT DATA AND FORMATS

1.2.3 OUTPUT FILES AND REPORTS

The computer calculates the final tissue partial pressure of nitrogen in PSI for each of three types of body tissue and the decompression ratio for each type of tissue at the time of final decompression. The program produces a summary report on the printing terminal, containing subject identification, the partial pressure of nitrogen in each of the three types of tissue and the decompression ratio for each of the three types of tissue. It also produces a data file on a floppy disk which contains all the data entered by the test operator and the computed information in the summary on partial pressure of nitrogen in the three types of tissue and the decompression ratio for each of the three types of tissue. Refer to section 3.2 PROGRAM OUTPUT, REPORTS AND DATA

1.2.4 CALCULATIONS

1. Computation of the Final Tissue Partial Pressure in psi of Nitrogen after an Exposure of t Minutes

$$FP = P_o + [(P_a - P_o) (1 - e^{-kt})]$$

FP: Final Tissue Nitrogen Pressure for a specific tissue type at the end of a washout phase

P_o: The original or beginning N₂ partial pressure (equivalent to ambient N₂ pressure in phase 1 of washout)

Pa: The N2 partial pressure in the breathing medium during the washout phase.

e: Base of natural logarithms

0.693

k: $\frac{0.693}{t_{1/2}}$

t 1/2: Tissue half saturation time in minutes

0.693: Logarithm to the base e of 2

2. Equation for the Computation of the Decompression Ratio

$$\text{Decompression Ratio} = \frac{\text{Final Tissue N2 Pressure (psi)}}{\text{Ambient (or Suit) Pressure (psi)}}$$

1.3 NITROGEN WASHOUT PROGRAM

1.3.1 PURPOSE

The N2WASH program is the main program for running the nitrogen washout experiment. A read only data file named DATA.CON.VAL is located on the system floppy and contains the initial calibration values used by the N2WASH program. The program prints the data constants and allows the test operator to make changes if necessary. It collects data on the subject's identification and the experimental conditions. The computer calibrates the spectrometer and spirometer and collects the test data as the experiment runs.

1.3.2 INPUT FILES AND DATA

The program prompts the test operator to enter data on the subject's identification and the experimental conditions. It collects the test data as the experiment runs. Refer to section 2.3 INPUT DATA AND FORMATS.

1.3.3 OUTPUT FILES AND REPORTS

During the washout, test data and one minute summaries are produced on the printing terminal. The test data is also being stored on a floppy disk along with the data that was entered by the test operator. Dummy zeroes mark the end of data input into the storage file. If a printing terminal is used, the hard copy provides a preliminary report on the nitrogen washout. A final report is obtained by using the N2DATA Program. Refer to section 3.3 PROGRAM OUTPUT, REPORTS AND DATA.

1.3.4 CALCULATIONS

1. Calculate Slopes and Intercepts for Hi and Lo range N2 and for O2 and CO2.

CALCNT is the digital count equivalent to the voltage recorded.

CALPCT is the concentration of a gas in a medium.

GASSLP is the slope calculated for a gas.

GSNTCP is the intercept calculated for a gas.

Nitrogen-Hi Range

$$\text{GASSLP} = \frac{\text{CALPCT N2 in Air} - \text{CALPCT N2 in Gas Mix}}{\text{CALCNT N2 in Air} - \text{CALCNT N2 in Gas Mix}}$$

$$\text{GSNTCP} = \text{CALPCT N2 in Gas Mix} - (\text{GASSLP Hi N2} * \text{CALCNT N2 in Gas Mix})$$

Nitrogen-Lo Range

$$\text{GASSLP} = \frac{\text{CALPCT N2 in Gas Mix} - \text{CALPCT N2 in 100\% O2}}{\text{CALCNT N2 in Gas Mix} - \text{CALCNT N2 in 100\% O2}}$$

$$\text{GSNTCP} = \text{CALPCT N2 in 100\% O2} - (\text{GASSLP Lo N2} * \text{CALCNT N2 in 100\% O2})$$

Oxygen

$$\text{GASSLP} = \frac{\text{CALPCT O2 in 100\% O2} - \text{CALPCT O2 in Gas Mix}}{\text{CALCNT O2 in 100\% O2} - \text{CALCNT O2 in Gas Mix}}$$

$$\text{GSNTCP} = \text{CALPCT O2 in Gas Mix} - (\text{GASSLP O2} * \text{CALCNT O2 in Gas Mix})$$

Carbon Dioxide

$$\text{GASSLP} = \frac{\text{CALPCT CO2 in Gas Mix} - \text{CALPCT CO2 in 100\% O2}}{\text{CALCNT CO2 in Gas Mix} - \text{CALCNT CO2 in 100\% O2}}$$

$$\text{GSNTCP} = \text{CALPCT CO2 in 100\% O2} - (\text{GASSLP CO2} * \text{CALCNT CO2 in 100\% O2})$$

- Calculate Slope and Intercept for Calibrating the Spirometer. One liter was measured at one volt and one volt equals 410.5 counts.

$$\text{Hi Cal Cnt} = \text{Lo Cal Cnt} + 410.5$$

$$\text{Slope} = \frac{\text{Hi Volume} - \text{Lo Volume}}{\text{Hi Cal CNT} - \text{Lo Cal Cnt}}$$

$$\text{Slope} = \frac{1 \text{ liter} - 0 \text{ liters}}{410.5} = \frac{1}{410.5}$$

$$\text{Intercept} = \text{Lo Volume} - (\text{Slope} \times \text{Lo Cal Cnt})$$

- Calculating the Average Digital Computer Count

The computer samples the gas channels 25/sec and reads the digital count equating to the analog voltage. At the end of a breath (determined by no increase in voltage in the spirometer) the counts are summed and divided by the number of samples occurring in the breath. The average digital count is computed for the spirometer and is used to compute the volume of a breath. An average digital count is also computed for each gas channel and is used to determine the concentration of a gas.

- Volume of Breath in Liters

$y = mx + b$, where x is the average digital computer count corresponding to average voltage in spirometer during breath.

- Gas Concentrations in Expired Breath

$y = mx + b$, where x is the average digital computer count corresponding to average voltage in gas channel during breath.

- $\text{N}_2 \text{ Volume of Breath} = \text{N}_2 \text{ Concentration} * \text{Volume of Breath}$
- $\text{N}_2 \text{ Volume for the Minute} = \text{Sum of N}_2 \text{ volumes for the minute}$
- $\text{Total N}_2 \text{ Volume} = \text{Sum of N}_2 \text{ volumes}$

1.4 NITROGEN DATA PROGRAM

1.4.1 PURPOSE

The N2DATA Program uses a data file from the N2WASH program to type a formatted report of the Nitrogen Washout Experiment.

1.4.2 INPUT FILES AND DATA

The input data for N2DATA is stored on a floppy disk in a data file created by the N2WASH program. Refer to section 2.4 INPUT DATA AND FORMATS

1.4.3 OUTPUT FILES AND REPORTS

There are four types of reports produced by the N2DATA program.

(1) A short report types the data on subject identification and experimental conditions with a summary of the average breath volume, average nitrogen volume per breath, average oxygen concentration per breath, average nitrogen concentration per breath, average carbon dioxide concentration per breath, total nitrogen expired and total number of breaths.

(2) A full report types the data on subject identification and experimental conditions with a breath by breath account of breath number, elapsed time in seconds, breath volume, nitrogen volume per breath, oxygen concentration per breath, nitrogen concentration per breath and carbon dioxide concentration per breath. It also types a summary of the average breath volume, average nitrogen volume per breath, average oxygen concentration per breath, average nitrogen concentration per breath, average carbon dioxide concentration per breath, total nitrogen expired and total number of breaths.

(3) A One Minute Interval report types the data on subject identification and experimental conditions with elapsed time in minutes and one minute readings of CO2 ml/min, O2 ml/min, N2 ml/min, N2 ml/kg of body-mass/min, N2 Out (Torr)/minute, cumulative N2 washed out, cumulative % of N2 estimate washed out, and the metabolic rate for the minute as $CO_2 * 1362.4$. It also types a summary of the average breath volume, average nitrogen volume per breath, average oxygen concentration per breath, average nitrogen concentration per breath, average carbon dioxide concentration per breath, total nitrogen expired and total number of breaths.

(4) A Five Minute Interval report types the data on subject identification and experimental conditions with five minute readings of CO2 ml/five-minute interval, O2 ml/five-minute interval, N2 ml/five-minute interval, N2 ml/kg of body-mass/five minute interval, N2 Out (Torr)/minute, cumulative N2 washed out, cumulative % of N2 estimate washed out, and the metabolic rate for five minutes as CO2 * 1362.4. It also types a summary of the average breath volume, average nitrogen volume per breath, average oxygen concentration per breath, average nitrogen concentration per breath, average carbon dioxide concentration per breath, total nitrogen expired and total number of breaths. Refer to section 3.4 PROGRAM OUTPUT, REPORTS AND DATA

1.4.4 CALCULATIONS

1. CO2 and O2 Volume/Breath in Liters x Breath Volume

CO2 Volume O2 Concentration x Breath Volume

O2 Volume = O2 Concentration * Breath Volume

2. The N2, O2, and CO2 Volume in ml for Interval

N2 = Sum of N2 x 1000

O2 = Sum of O2 x 1000

CO2 = Sum of CO2 x 1000

3. The 1 Minute Average N2 Concentration

$$= \frac{\text{Sum of N2 Concentration Readings}}{\text{Number of Readings}}$$

4. The 1 Minute N2 ml/kg of Body Weight

$$= \frac{\text{N2 Volume in ml}}{\text{Body Weight in kg}}$$

5. The Cumulative N2 Eliminated by 1 Minute Intervals

= Sum of N2 Eliminated in Previous Intervals

6. The % of Estimated Body N2 Content Eliminated by 1 Minute Intervals

$$= \left(\frac{\text{Cumulative N2 / 1000}}{\text{Estimated N2 Content of Body}} \right) \times 100$$

7. The Metabolic Rate (BTU/hr) during the 1 Minute Interval

$$= \text{Liters of CO}_2 \text{ Eliminated per Minute} \times 1362.4$$

8. The 5 Minute Average N₂ Concentration

$$= \frac{\text{Sum of N}_2 \text{ Concentration Readings}}{\text{Number of Readings}}$$

9. The 5 Minute N₂ ml/kg of Body Weight

$$= \frac{\text{N}_2 \text{ Volume in ml}}{\text{Body Weight in kg}}$$

10. The Cumulative N₂ Eliminated by 5 Minute Intervals

$$= \text{Sum of N}_2 \text{ Eliminated in Previous Intervals}$$

11. The % of Estimated Body N₂ Content Eliminated by 5 Minute Intervals

$$= \left(\frac{\text{Cumulative N}_2 \text{ (ml)} / 1000}{\text{Estimated N}_2 \text{ Content of Body}} \right) \times 100$$

12. The Metabolic Rate (BTU/hr) during the 5 Minute Interval

$$= \text{Average Liters of CO}_2 \text{ Eliminated per Minute} \times 1362.4$$

13. Average of Breath Volumes

14. Average of N₂ Volume of Breath

15. Average N₂ Concentration per Breath

16. Average O₂ Concentration per Breath

17. Average CO₂ Concentration per Breath

1.5 UPDA TE PROGRAM

1.5.1 PURPOSE

The UPDA TE program was designed for limited use. It reads the unformatted data stored on a floppy disk file from an earlier version nitrogen washout program and writes it to another file so that it is formatted. It also records blanks or zeroes for the header information data which was not included in the early nitrogen washouts. The updated data floppy disk file format is then identical to the ones being produced by the current version of the nitrogen washout program.

1.5.2 INPUT FILES AND DA TA

The input for UPDA TE is the unformatted data stored on a floppy disk file from an earlier version of the nitrogen washout program. Refer to 2.5 INPUT DA TA FORMATS.

1.5.3 OUTPUT FILES AND REPORTS

UPDA TE creates a formatted file on a floppy disk which is identical to the format of the files being produced by the current version of the nitrogen washout program. Refer to section 3.5 PROGRAM OUTPUT REPORTS AND DA TA.

2.0 INPUT DATA AND FORMATS

2.1 PERSONAL DATA INPUT PROGRAM

The following table lists the variables which the test operator will be asked to enter.

INPUT VARIABLE	FORMAT	UNITS	RANGE
FILENAME	(FILENAME.EXT, e.g., SMITH.PDI)		
SUBJECT ID CODE	(XXX-XX-XXXX)		
DATE OF DATA COLLECTION	(MM/DD/YY)		
DATE OF BIRTH	(MM/DD/YY)		
AGE	(XX)	yrs	0-99
SEX	(A)		(M/F)
HEIGHT	(XXX.X)	cm	0.0-250.0
NUDE WEIGHT	(XXX.X)	kg	0.0-250.0
KNOWN % BODY FAT	(XX.X)		0.0-99.9

WATER IMMERSION PROCESS

WATER TEMP	(XX.XX)	deg C	0.0-99.99
WATER DENSITY	(X.XXX)	kg/L	0.0-9.999
RESIDUAL VOLUME	(XX.XX)	L	0.0-99.99
UNDER WATER CHAIR WEIGHT	(XX.XX)	kg	0.0-99.99
NUMBER OF TRIALS	(XX)		1-10

(Let N represent the number of the trial and let N increase from 1 to the number of trials.)

UNDER WATER WEIGHT N	(XX.XX)	kg	0.0-99.99
MEASUREMENT DATE	(MM/DD/YY)		

SKINFOLD PROCESS I

BODY SIDE USED (A) (R/L)

NUMBER OF TRIALS (X) 1-3

(Let N represent the number of the trial and let N increase from 1 to the number of trials.)

BICEPS N (XX.XX) mm 0.0-99.99

TRICEPS N (XX.XX) mm 0.0-99.99

SUBSCAPULAR N (XXX.X) mm 0.0-250.0

SUPRAILIAC N (XXX.X) mm 0.0-250.0

DATE OF MEASUREMENTS (MM/DD/YY)

SKINFOLD PROCESS II

BODY SIDE USED (A) (R/L)

NUMBER OF TRIALS (X) 1-3

(Let N represent the number of the trial and let N increase from 1 to the number of trials.)

MIDCLAVICULAR N (XXX.X) mm 0.0-250.0

UMBILICUS N (XXX.X) mm 0.0-250.0

ANTERIOR MID-THIGH N (XXX.X) mm 0.0-250.0

DATE OF MEASUREMENTS (MM/DD/YY)

2.2 DECOMPRESSION DATA PROGRAM

The following table lists the variables which the test operator will be asked to enter.

INPUT VARIABLES	FORMAT	UNITS	RANGE
FILENAME	(FILENAME .EXT, e.g., SMITH.DEC)		
SUBJECT ID CODE	(XXX-XX-XXXX)		
TEST NUMBER	(XX)		(1-99)
TEST DATE	(MM-DD-YY)		
INITIAL PARTIAL PRESSURE OF NITROGEN	(XX.XX)	PSI	(4.0-20.99)
NUMBER OF WASHOUT PHASES	(XX)		(1-10)
(Let N represent the number of the washout phase and let N increase from 1 to the number of washout phases so that these next two items are entered for each phase)			
PARTIAL PRESSURE OF NITROGEN IN WASHOUT PHASE N	(XX.XX)	PSI	(0.0-12.99)
DURATION OF WASHOUT PHASE N	(XXXX)	MIN	(0-1440)
FINAL PRESSURE	(XX.XX)	PSIA	(0.0-10.99)
TEMPERATURE	(XX)	deg C	(0-40)
TIME OF FINAL DECOMPRESSION	(HH:MM)	hrs:min	0-24:0-59
TIME OF RECOMPRESSION	(HH:MM)	hrs:min	0-24:0-59
EXERCISE LEVEL	(L/M/H)		
BUBBLES MEASURED?	(Y/N)		
BUBBLES PRESENT?	(Y/N)		
TIME OF ONSET OF BUBBLES	(HH:MM)	hrs:min	0-24:0-59
TIME OF REMISSION OF BUBBLES	(HH:MM)	hrs:min	0-24:0-59

PRESSURE AT REMISSION OF BUBBLES	(XX.X)		
BENDS PRESENT?	(Y/N)		
LOCATION OF BENDS	(Text)		
TIME OF ONSET OF BENDS	(HH:MM)	hrs:min	0-24:0-59
TIME OF REMISSION OF BENDS	(HH:MM)	hrs:min	0-24:0-59
PRESSURE AT REMISSION OF BENDS	(XX.X)		

2.3 NITROGEN WASHOUT PROGRAM

The following table lists the variables which the test operator will be asked to enter.

INPUT VARIABLES	FORMAT	UNITS	RANGE
FILENAME	(FILENAME.EXT, e.g., SMITH.N2W)		
DATE	(MM/DD/YY)		
NAME	(Up TO 30 CHAR)		
SUBJECT IDCODE	(xxx-xx-xxxx)		
SUBJECT AGE	(xx)	yrs	(0-99)
EXTERNAL TEMP	(xx.xx)	deg C	(0.0-99.99)
EXTERNAL GRAVITY	(1/0)		
SUBJECT IMMERSED	(Y/N)		
BODY POSITION CODE	(xx)		(0-99)
AMBIENT PRESSURE	(xxx.xx)	mm Hg	(0.0-999.99)
DRUG CODE	(xx)		(0-99)
EXERCISE CODE	(xx)		(0-99)
BREATHING MIXTURE CODE	(xx)		(0-99)
BODY HYDRATION CODE	(xx)		(0-99)
AMBIENT GAS CODE	(XX)		(0-99)
STARTING TIME	(HH:MM)	hrs:min	0-24:0-59
SPECIAL CONDITION CODE	(XX)		
COMMENTS	(UP TO 4 LINES OF TEXT)		

2.4 NITROGEN DATA PROGRAM

The input data for N2DATA is the formatted data stored on a floppy disk from a nitrogen washout. The test operator supplies the name of the data file in which the data is stored and answers yes or no to the question, "Do you want a full report?" If the answer is "no" a short report is produced.

2.5 UPDATE PROGRAM

The input data for the UPDATE program is the unformatted data stored on a floppy disk from an early nitrogen washout. The test operator enters the name of the data file where the unformatted data is stored and the name of the file where he would like the newly formatted data to be stored. The program puts zeros or blanks in the new data file for the new variables created in the modified nitrogen washout program.

3.0 PROGRAM OUTPUT, REPORTS AND DATA

3.1 PERSONAL DATA INPUT PROGRAM

The following table lists all the variables which are output in the hard copy report and which are stored on a data floppy for later use.

OUTPUT VARIABLE	FORMAT	UNITS
DESCRIPTIVE OUTPUT		
SUBJECT ID CODE	(XXX-XX-XXXX)	
FILENAME	(FILENAME.EXT, e.g., SMITH.PDI)	
% BODY FAT	(XX.X)	
NUDE WEIGHT	(XXX.X)	kg
DATE OF BIRTH	(DD/MON/YY)	
BODY COMPOSITION OUTPUT BY WATER IMMERSION		
NUDE WEIGHT	(XXX.XX)	kg
WATER TEMP	(XX.XX)	deg C
RESIDUAL VOLUME	(XX.XX)	L
WATER DENSITY	(X.XXX)	kg/L
UNDER WATER CHAIR WEIGHT	(XX.XX)	kg
AVERAGE UNDER WATER WEIGHT	(XX.XX)	kg
(Let N represent the number of the trial and let N increase from 1 to the number of trials.)		
TRIAL N UNDER WATER WEIGHT	(XX.XX)	kg
WATER WEIGHT	(XX.XX)	kg
BODY VOLUME	(XX.XX)	L
BODY DENSITY	(XX.XX)	kg/L
% BODY FAT	(XX.XX)	

BODY COMPOSITION OUTPUT BY SKINFOLD PROCESS I

BODY SIDE USED (RIGHT/LEFT)

DATE (DD/MON/YY)

(Let N represent the number of the trial and let N increase from 1 to the number of trials.)

TRIAL N BICEPS	(XX.XX)	mm
TRIAL N TRICEPS	(XX.XX)	mm
TRIAL N SUBSCAPULAR	(XX.XX)	mm
TRIAL N SUPRAILIAIC	(XX.XX)	mm
AVERAGE BICEPS	(XX.XX)	mm
AVERAGE TRICEPS	(XX.XX)	mm
AVERAGE SUBSCAPULAR	(XX.XX)	mm
AVERAGE SUPRAILIAIC	(XX.XX)	mm
BODY DENSITY	(X.XXX)	kg/L
% BODY FAT	(XX.XX)	

BODY COMPOSITION OUTPUT BY SKINFOLD PROCESS II

BODY SIDE USED (RIGHT/LEFT)

DATE (DD/MON/YY)

(Let N represent the number of the trial and let N increase from 1 to the number of trials.)

TRIAL N MIDCLAVICULAR	(XX.XX)	mm
TRIAL N UMBILICUS	(XX.XX)	mm
TRIAL N ANTERIOR MID-THIGH	(XX.XX)	mm
AVERAGE MIDCLAVICULAR	(XX.XX)	mm
AVERAGE UMBILICUS	(XX.XX)	mm
AVERAGE ANTERIOR MID-THIGH	(XX.XX)	mm

BODY DENSITY	(X.XXX)	kg/L
% BODY FAT	(XX.XX)	

BODY COMPOSITION OUTPUT BY NUTRITION JOURNAL FORMULA

BODY VOLUME	(XXX.XX)	L
BODY DENSITY	(XX.XXX)	kg/L
% BODY FAT	(XX.XX)	
% FAT FREE MASS	(XX.XX)	
% TOTAL BODY WATER	(XX.XX)	

NITROGEN CONTENT OUTPUT FOR EACH OF THE 5 METHODS

(1) DESCRIPTIVE DATA - KNOWN % BODY FAT

% BODY FAT	(XX.XX)	
FAT WEIGHT	(XX.XX)	kg
FAT FREE MASS	(XX.XX)	kg
LEAN BODY WATER MASS	(XX.XX)	kg
FAT VOLUME	(XX.XX)	L
LEAN BODY WATER VOLUME	(XX.XX)	L
NITROGEN IN FAT	(X.XXX)	L
NITROGEN IN LEAN BODY WATER	(X.XXX)	L
TOTAL BODY NITROGEN	(X.XXX)	L
NITROGEN/KG BODY WEIGHT	(X.XXXX)	L

(2) WATER IMMERSION PROCESS

% BODY FAT	(XX.XX)	
FAT WEIGHT	(XX.XX)	kg
FAT FREE MASS	(XX.XX)	kg
LEAN BODY WATER MASS	(XX.XX)	kg
FAT VOLUME	(XX.XX)	L
LEAN BODY WATER VOLUME	(XX.XX)	L
NITROGEN IN FAT	(X.XXX)	L
NITROGEN IN LEAN BODY WATER	(X.XXX)	L
TOTAL BODY NITROGEN	(X.XXX)	L
NITROGEN/KG BODY WEIGHT	(X.XXXX)	L

(3) SKINFOLD PROCESS I

% BODY FAT	(XX.XX)	
FAT WEIGHT	(XX.XX)	kg
FAT FREE MASS	(XX.XX)	kg
LEAN BODY WATER MASS	(XX.XX)	kg
FAT VOLUME	(XX.XX)	L
LEAN BODY WATER VOLUME	(XX.XX)	L
NITROGEN IN FAT	(X.XXX)	L
NITROGEN IN LEAN BODY WATER	(X.XXX)	L
TOTAL BODY NITROGEN	(X.XXX)	L
NITROGEN/KG BODY WEIGHT	(X.XXXX)	L

(4) SKINFOLD PROCESS II

% BODY FAT	(XX.XX)	
FAT WEIGHT	(XX.XX)	kg
FAT FREE MASS	(XX.XX)	kg
LEAN BODY WATER MASS	(XX.XX)	kg
FAT VOLUME	(XX.XX)	L
LEAN BODY WATER VOLUME	(XX.XX)	L
NITROGEN IN FAT	(X.XXX)	L
NITROGEN IN LEAN BODY WATER	(X.XXX)	L
TOTAL BODY NITROGEN	(X.XXX)	L
NITROGEN/KG BODY WEIGHT	(X.XXXX)	L

(5) NUTRITION JOURNAL FORMULA

% BODY FAT	(XX.XX)	
FAT WEIGHT	(XX.XX)	kg
FAT FREE MASS	(XX.XX)	kg
LEAN BODY WATER MASS	(XX.XX)	kg
FAT VOLUME	(XX.XX)	L
LEAN BODY WATER VOLUME	(XX.XX)	L
NITROGEN IN FAT	(X.XXX)	L
NITROGEN IN LEAN BODY WATER	(X.XXX)	L
TOTAL BODY NITROGEN	(X.XXX)	L
NITROGEN/KG BODY WEIGHT	(X.XXXX)	L

3.2 DECOMPRESSION DATA PROGRAM

The DECDAT Output consists of the hard copy from the printing terminal of all the input variables as they are entered by the test operator and the summary generated by the computer of the final tissue nitrogen partial pressure and decompression ratio at time of final decompression for each of the three types of tissues. Tissues of half saturation time 180 minutes, 240 minutes and 360 minutes will be used. The following table lists the variables calculated by the computer which make up the summary.

OUTPUT VARIABLE	FORMAT
PARTIAL PRESSURE OF NITROGEN IN 180 TISSUE	(XX.X)
PARTIAL PRESSURE OF NITROGEN IN 240 TISSUE	(XX.X)
PARTIAL PRESSURE OF NITROGEN IN 360 TISSUE	(XX.X)
DECOMPRESSION RATIO FOR 180 TISSUE	(XX.X)
DECOMPRESSION RATIO FOR 240 TISSUE	(XX.X)

DECOMPRESSION RATIO
FOR 360 TISSUE (XX.X) DECDAT also writes a formatted data file containing the variables on subject identification, washout data, experimental conditions and the data on the chamber run.

RECORD 1

ITEM	FORMAT	STARTING COLUMN
IDCODE	11A1	1
TEST NUMBER	I2	12
TEST DATE	8A1	14
NUMBER OF PHASES	I2	22
INITIAL PARTIAL PRESSURE OF N2	F5.2	24
FINAL PRESSURE	F5.2	29
TEMPERATURE	I2	34
TIME OF FINAL DECOMPRESSION	5A1	36

TIME OF RECOMPRESSION	5A1	41
EXERCISE LEVEL	A1	46
BUBBLES MEASURED?	A1	47
BUBBLES PRESENT?	A1	48

RECORD 2

ITEM	FORMAT	STARTING COLUMN
------	--------	-----------------

(Let N represent the number of the washout phase and let N increase from 1 to the number of phases so that partial pressure is entered for each phase.)

PARTIAL PRESSURE OF N2 DURING PHASE N	F5.2	1
---------------------------------------	------	---

RECORD 3

ITEM	FORMAT	STARTING COLUMN
------	--------	-----------------

(Let N represent the number of the washout phase and let N increase from 1 to the number of phases so that duration of phase is entered for each phase.)

DURATION OF PHASE N	I4	1
---------------------	----	---

RECORD 4

ITEM	FORMAT	STARTING COLUMN
------	--------	-----------------

BUBBLES GRADE 1 ONSET TIME	5A1	1
----------------------------	-----	---

BUBBLES GRADE 2 ONSET TIME	5A1	6
----------------------------	-----	---

BUBBLES GRADE 3 ONSET TIME	5A1	11
----------------------------	-----	----

BUBBLES GRADE 4 ONSET TIME	5A1	16
----------------------------	-----	----

RECORD 5

ITEM	FORMAT	STARTING COLUMN
------	--------	-----------------

BUBBLES GRADE 1 REMISSION TIME	5A1	1
--------------------------------	-----	---

BUBBLES GRADE 2 REMISSION TIME	5A1	6
--------------------------------	-----	---

BUBBLES GRADE 3
REMISSION TIME 5A1 11

BUBBLES GRADE 4
REMISSION TIME 5A1 16

RECORD 6

ITEM FORMAT STARTING COLUMN

BUBBLES GRADE 1
REMISSION PRESSURE F4.1 1

BUBBLES GRADE 2
REMISSION PRESSURE F4.1 5

BUBBLES GRADE 3
REMISSION PRESSURE F4.1 9

BUBBLES GRADE 4
REMISSION PRESSURE F4.1 13

RECORD 7

ITEM FORMAT STARTING COLUMN

BENDS PRESENT? A1 1

BENDS LOCATION? 15A1 2

RECORD 8

ITEM FORMAT STARTING COLUMN

BENDS GRADE 1
ONSET TIME 5A1 1

BENDS GRADE 2
ONSET TIME 5A1 6

BENDS GRADE 3
ONSET TIME 5A1 11

BENDS GRADE 4
ONSET TIME 5A1 16

RECORD 9

ITEM FORMAT STARTING COLUMN

BENDS GRADE 1
REMISSION TIME 5A1 1

BENDS GRADE 2 REMISSION TIME	5A1	6
BENDS GRADE 3 REMISSION TIME	5A1	11
BENDS GRADE 4 REMISSION TIME	5A1	16

RECORD 10

ITEM	FORMAT	STARTING COLUMN
BENDS GRADE 1 REMISSION PRESSURE	F4.1	1
BENDS GRADE 2 REMISSION PRESSURE	F4.1	5
BENDS GRADE 3 REMISSION PRESSURE	F4.1	9
BENDS GRADE 4 REMISSION PRESSURE	F4.1	13

RECORD 11

ITEM	FORMAT	STARTING COLUMN
PARTIAL PRESSURE OF NITROGEN IN 180 TISSUE	F5.2	1
PARTIAL PRESSURE OF NITROGEN IN 240 TISSUE	F5.2	6
PARTIAL PRESSURE OF NITROGEN IN 360 TISSUE	F5.2	11
DECOMPRESSION RATIO FOR 180 TISSUE	F4.2	16
DECOMPRESSION RATIO FOR 240 TISSUE	F4.2	20
DECOMPRESSION RATIO FOR 360 TISSUE	F4.2	24

3.3 NITROGEN WASHOUT PROGRAM

THE NITROGEN WASHOUT program output consists of a data file and the data which was output to the printing terminal during the run.

The data file contains all of the input variables on subject identification and experimental conditions plus the values collected in the run. For each breath there is :

Record 1

ITEM	FORMAT	STARTING COLUMN
NAME	30A1	1
DATE	8A1	31
IDCODE	11A1	39
AGE	I2	50
EXTERNAL TEMPERATURE	F7.2	52
EXTERNAL GRAVITY	I1	59
IMMERSION?	A1	60
BODY POSITION CODE	I2	61
AMBIENT PRESSURE	F7.2	63
DRUGS CODE	I2	70
EXERCISE CODE	I2	72
BREATHING MIXTURE CODE	I2	74
BODY HYDRATION CODE	I2	76
AMBIENT GAS CODE	I2	78
STARTING TIME	5A1	80
SPECIAL CONDITION CODE	I2	85
COMMENTS?	A1	87

RECORD 2

ITEM	FORMAT	STARTING COLUMN
COMMENT	80A1	1

RECORD 3

ITEM	FORMAT	STARTING COLUMN
COMMENT	80A1	1

RECORD 4

ITEM	FORMAT	STARTING COLUMN
COMMENT	80A1	1

RECORD 5

ITEM	FORMAT	STARTING COLUMN
COMMENT	80A1	1

RECORD 6 thru N, where N is the number of breaths

ITEM	FORMAT	STARTING COLUMN
BREATH NUMBER	I4	1
ELAPSED TIME	F10.1	5
BREATH VOLUME	F10.5	15
NITROGEN VOLUME IN BREATH	F10.5	25
CONCENTRATION OF O2 IN BREATH	F10.5	35
CONCENTRATION OF N2 IN BREATH	F10.5	45
CONCENTRATION OF CO2 IN BREATH	F10.5	55

The data file end is marked by dummy zero values for the variables in the last record.

The output to the printing terminal consists of the input variables on subject identification and experimental conditions and the following parameters for each breath:

1. Breath number
2. Elapsed time in seconds
3. Breath volume in liters
4. Nitrogen volume per breath in liters
5. Concentration of nitrogen in breath

At the end of every minute there is a summary of elapsed minutes, nitrogen eliminated for the minute and total nitrogen eliminated.

3.4 NITROGEN DATA PROGRAM

The N2DATA Program produces a hard copy report on the printing terminal of the data from a nitrogen washout.

The output for a short report is:

1. Subject identification
2. Experimental conditions
3. Average breath volume
4. Average nitrogen volume per breath
5. Average oxygen concentration per breath
6. Average nitrogen concentration per breath
7. Average carbon dioxide concentration per breath
8. Total nitrogen expired
9. Total number of breaths

A full report requires about 20 minutes. The output for a full report is the subject's identification, experimental conditions and the breath by breath account of:

1. Breath number
2. Elapsed time
3. Breath volume
4. Nitrogen volume in breath
5. Oxygen concentration in breath
6. Nitrogen concentration in breath
7. Carbon dioxide concentration in breath
8. Average breath volume
9. Average nitrogen volume per breath
10. Average oxygen concentration per breath
11. Average nitrogen concentration per breath
12. AVERAGE carbon dioxide concentration per breath
13. Total nitrogen eliminated
14. Total number of breaths

The output for a one minute report is the subject's identification, experimental conditions and the one minute interval parameters for:

1. Elapsed time in minutes
2. Nitrogen volume eliminated during the minute
3. Oxygen volume eliminated during the minute
4. Carbon dioxide volume eliminated during the minute
5. Nitrogen volume per kilogram of body weight during the interval
6. Average nitrogen concentration during the interval
7. Cumulative nitrogen volume eliminated
8. Per-cent of estimated nitrogen content eliminated
9. Metabolic rate during the interval
10. Average breath volume
11. Average nitrogen volume per breath
12. Average oxygen concentration per breath
13. Average nitrogen concentration per breath
14. Average carbon dioxide concentration per breath
15. Total nitrogen eliminated

16. Total number of breaths

The output for a five minute report is the subject's identification, experimental conditions and the one minute interval parameters for:

1. Elapsed time in minutes
2. Nitrogen volume eliminated during the five minute interval
3. Oxygen volume eliminated during the five minute interval
4. Carbon dioxide volume eliminated during the five minute interval
5. Nitrogen volume per kilogram of body weight during the interval
6. Average nitrogen concentration during the interval
7. Cumulative nitrogen volume eliminated
8. Per-cent of estimated nitrogen content eliminated
9. Metabolic rate during the interval
10. Average breath volume
11. Average nitrogen volume per breath
12. Average oxygen concentration per breath
13. Average nitrogen concentration per breath
14. Average carbon dioxide concentration per breath
15. Total nitrogen eliminated
16. Total number of breaths

3.5 UPDATE PROGRAM

The output of the UPDATE program is a formatted data file identical to the data files generated by the modified nitrogen washout program and suitable for use with the N2DATA report writing program or for transfer to the VAX. The updated data file is stored on a floppy disk. The data file contains zeros and blanks for all of the new input variables on subject identification and experimental conditions.

Record 1

ITEM	FORMAT	STARTING COLUMN
NAME	30A1	1
DATE	8A1	31
IDCODE	11A1	39
AGE	I2	50
EXTERNAL TEMPERATURE	F7.2	52
EXTERNAL GRAVITY	I1	59
IMMERSION?	A1	60
BODY POSITION CODE	I2	61
AMBIENT PRESSURE	F7.2	63
DRUGS CODE	I2	70
EXERCISE CODE	I2	72
BREATHING MIXTURE CODE	I2	74
BODY HYDRATION CODE	I2	76
AMBIENT GAS CODE	I2	78
STARTING TIME	5A1	80
SPECIAL CONDITION CODE	I2	85
COMMENTS?	A1	87

RECORD 2

ITEM	FORMAT	STARTING COLUMN
COMMENT	80A1	1

RECORD 3

ITEM	FORMAT	STARTING COLUMN
COMMENT	80A1	1

RECORD 4

ITEM	FORMAT	STARTING COLUMN
COMMENT	80A1	1

RECORD 5

ITEM	FORMAT	STARTING COLUMN
COMMENT	80A1	1

The data file contains a record for each breath of the following variables collected in the run.

RECORD 6


ITEM	FORMAT	STARTING COLUMN
BREATH NUMBER	I4	1
ELAPSED TIME	F10.1	5
BREATH VOLUME	F10.5	15
NITROGEN VOLUME IN BREATH	F10.5	25
CONCENTRATION OF O2 IN BREATH	F10.5	35
CONCENTRATION OF N2 IN BREATH	F10.5	45
CONCENTRATION OF CO2 IN BREATH	F10.5	55

The data file end is marked by dummy zero values for the variables in record 6.


ORIGINAL PAGE IS
OF POOR QUALITY

4.0 PROGRAM FLOWCHARTS


4.1 PERSONAL DATA INPUT PROGRAM




ORIGINAL PART 19
OF POOR QUALITY




ORIGINAL PAGE IS
OF POOR QUALITY




ORIGINAL
OF POOR QUALITY




ORIGINAL FIGURE 13
OF POOR QUALITY




ORIGINAL PAGE IS
OF POOR QUALITY




ORIGINAL PAGE IS
OF POOR QUALITY




4.2 DECOMPRESSION DATA PROGRAM




ORIGINAL PAGE IS
OF POOR QUALITY




ORIGINAL PAGE IS
OF POOR QUALITY




ORIGINAL PAGE IS
OF POOR QUALITY







ORIGINAL BY THE IS
OF POOR QUALITY




CLASSIFICATION
OF FLOW CHARTS




UNION OF
OF POOR QUALITY



4.3 NITROGEN WASHOUT PROGRAM






OPERATOR CAN CHANGE SELECTED
DATA VALUES

ORIGINAL PERCENTAGES
OF POOR QUALITY


CAL GAS PERCENTAGES

INTERRUPT RATE i.e. 'SAMPLE'
EXECUTION RATE


SPIROMETER NOISE, MINIMUM
BREATH DIFFERENCE



ORIGINAL VERSIONS
OF POOR QUALITY



ORIGIN
OF POINT




SELECT MASS SPEC PORT &
OPEN A VALVE FOR VENTING

SELECT MASS SPEC N2 RANGE
START CAL. SEQUENCE

OPEN A VALVE

FLUSH FOR 10 SEC'S

ORIGINAL DOCUMENT
OF POOR QUALITY



ZERO INTEGRATOR

READ A D CHANNEL

ADD TO INTEGRATOR

AVERAGE VALUE TO GET
ANSWER

INCREMENT INDEX

REPEAT IF NECESSARY

SELECT MASS SPEC N2 RANGE
START CAL SEQUENCE

OPEN A VALVE

FLUSH FOR 10 SEC'S

ORIGINAL PAGE NO
OF POOR QUALITY

MASS
SPEC
2

GASSUM(1)=0.0

ZERO INTEGRATOR

CALL IPOKE (IADCSR,
IADCGC(1))
INT = IPEEK IADBUF
GASSUM (1) =
GASSUM (1) +
FLOAT (INT)

READ A-D CHANNEL
ADD TO INTEGRATOR

CALCNT (1,J) =
GASSUM (1)/100

AVERAGE VALUE TO GET
ANSWER

INCREMENT INDEX

REPEAT IF NECESSARY

J = 3

(1) = [CALPCT(1,1) - CALPCT(1,3)]/[CALCNT(1,1) - CALCNT(1,3)]
(2) = " (2,2) " (2,3) " (2,2) " (2,3)
(3) = " (3,3) " (3,2) " (3,3) " (3,2)
(4) = " (1,3) " (1,2) " (1,4) " (1,2)

GSNTCP(1) = CALPCT(1,3) - [GASSLP(2) * CALCNT(1,3)]
(2) " (2,3) " (2) " (2,3)
(3) " (3,2) " (3) " (3,2)
(4) " (1,2) " (4) " (1,2)

CALCULATE SLOPES
AND INTERCEPTS

TYPE CALCNT(1,1),CALCNT(1,3)
GASSLP(1) GSNTCP(1)
TYPE CALCNT(1,3),CALCNT(1,2)
GASSLP(4) GSNTCP(4)
TYPE CALCNT(2,2),CALCNT(2,3)
GASSLP(2) GSNTCP(2)
TYPE CALCNT(3,3),CALCNT(3,2)
GASSLP(3) GSNTCP(3)

PRINT RESULTS


CALL IPOKE
(IDRVOT,0)
PRINT 'END OF PHASE'

SHUT VALUES, CHANGE MASS
SPEC PORT

92

GET NEW PHASE ENTRY

ORIGINAL FILE
OF POOR QUALITY



READ SPIROMETER 0

ZERO ACCUMULATOR

READ A-D
ADD TO ACCUMULATOR


AVERAGE RESULTS

VALUE FOR HIGH CAL. VOLUME

GET SLOPE AND INTERCEPT

PRINT RESULTS

GET NEW PHASE ENTRY



START OF PROGRAM
OF POOR QUALITY

B

CALL INTSET ("104,4,ID,SAMPLE)
ICLCNT = 10000/ICLRAT
CALL IPOKE ("172542, ICKNT)

ATTACH SUBROUTINE
SAMPLE TO INTERRUPT

IDRIFT = 2.0	ICON = 0	BRTST = FLOZER + BRGAT
N2TOT = 0.0	IFLO = 0	VOLUME = 0.0
1MIN = 1	IBTH = 0	VOL1 = BRTST
	N2VLM = 0.0	VOL2 = 0.0
	LWAT = .f.	1CHECK = 0
	LATCH = .t.	

CALCULATE BUFFER COUNT
FOR PROGRAM CLOCK

TINT = 1.0/FLOAT (ICDRAT)
DELAY = FLOAT (ICLRAT)
 *1.1
IDL M = IFIX (DELAY)

INITIALIZE SOME COUNTER'S
AND LATCHES

TSEC = SECONDS (0.0)
TOCK = TSEC

CALCULATE INTEPRATING
FACTOR


CALCULATE M.S. TRANSPORT
LAG (=1.1 SEC'S) COUNTERS

CALL IPOKE ("172540, "0113)

INITIALIZE TIME TAG'S
AND EXPERIMENT CLOCK

510

START CLOCK GENERATING
INTERRUPTS



CHECK MASS SPEC. N2 RANGE

LOOK FOR KEYBOARD INPUT

"S car "
set


CHECK NEED TO PURGE
SPIROMETER

PURGE SPIROMETER (OPEN DUMP)


CLOSE DUMP VALVE

IS DATA READY

ORIGINAL IS
OF POOR QUALITY



ORIGINAL PAGE IS
OF POOR QUALITY



INTERRUPT SERVICE ROUTINE

SAMPLE N2 DRIFT OF
INSPIRED O2


CHECK FOR MAIN
PROGRAM READY

READ MASS SPEC
N2 RANGE


GET VOLUME SAMPLE
AND TEST FOR END OF
BREATH

TURN ON INDICATOR LIGHT

INCREMENT VOLUME INDEX .
CALCULATE VOLUME



ORIGINAL PAGE IS
OF POOR QUALITY



INCREMENT TRANSPORT DELAY
COUNTER

IS DELAY COUNTER TRANSPORT
LAG TIME


INCREMENT GAS CONCENTRATION
INDEX

READ MASS SPE AND CALCULATE GAS
CONCENTRATIONS FOR N2, O2, AND
CO2.


ARE INDEX'S EQUAL

ALL DATA IS GATHERED: SET
FLAG'S TO GIVE MAIN PROGRAM
ANALYSIS CONTROL.


4.4 NITROGEN DATA PROGRAM




ORIGINAL PAGE IS
OF POOR QUALITY




ORIGINAL COPY
OF FOUR (4)...



ORIGINAL PAGE IS
OF POOR QUALITY




ORDER OF
OF PRIORITY



4.5 UPDATE PROGRAM

ORIGINAL PAGE IS
OF POOR QUALITY



ORIGINAL DOCUMENT
OF POOR QUALITY

```

C
0001 CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
      PROGRAM PDINPT
C
      TABLE OF VARIABLES
C
C     VARIABLE                               USE
C     AGE                                     SUBJECT'S AGE
C     AVGBIC                                 AVERAGE OF BICEPS MEASUREMENTS
C     AVGMCL                                 AVERAGE OF MIDCLAVICULAR MEASUREMENTS
C     AVGMTH                                 AVERAGE OF MID THIGH MEASUREMENTS
C     AVGSUB                                 AVERAGE OF SUBSCAPULAR MEASUREMENTS
C     AVGSUP                                 AVERAGE OF SUPRAILIAC MEASUREMENTS
C     AVGTRI                                 AVERAGE OF TRICEPS MEASUREMENTS
C     AVGUMB                                 AVERAGE OF UMBILICUS MEASUREMENTS
C     AVSUM1                                 SUM OF AVERAGES OF BICEPS,
C                                           TRICEPS, SUBSCAPULAR AND SUPRAILIAC
C     AVSUM2                                 SUM OF AVERAGES OF
C                                           MIDCLAVICULAR, UMBILICUS AND MIDTHIGH
C     AVUWW7                                 AVERAGE OF UNDERWATER WEIGHT TRIALS
C     BDP1                                    BODY DENSITY IN SKINFOLD PROCESS I
C     BDP2                                    BODY DENSITY IN SKINFOLD PROCESS II
C     BDWI                                    BODY DENSITY IN WATER IMMERSION PROCESS
C     BFKNOW                                  % BODY FAT IF KNOWN
C     BFPC T                                  TEMPORARY FOR % BODY FAT
C     BSID1                                   BODY SIDE USED IN SKINFOLD PROCESS I
C     BSID2                                   BODY SIDE USED IN SKINFOLD PROCESS II
C     BVOLWI                                  BODY VOLUME IN WATER IMMERSION PROCESS
C     DEN                                     BODY DENSITY IN NUTRITION JOURNAL FORMULA
C     HEIGHT                                  SUBJECT'S HEIGHT
C     J                                       INDEX VARIABLE
C     K                                       INDEX VARIABLE
C     KSP1                                    DO SKINFOLD PROCESS I?
C     KSP2                                    DO SKINFOLD PROCESS II?
C     KWI                                     DO WATER IMMERSION PROCESS?
C     K1                                       INTERMEDIATE VALUE
C     K2                                       INTERMEDIATE VALUE
C     M                                       INDEX VARIABLE
C     MEAS1                                    NUMBER OF MEASUREMENTS MADE IN PROCESS I
C     MEAS2                                    NUMBER OF MEASUREMENTS MADE IN PROCESS II
C     N                                       INDEX VARIABLE
C     NC                                      NUMBER OF CHARACTERS
C     NF                                      NUMBER OF CHARACTERS
C     NUDEWT                                  SUBJECT'S WEIGHT
C     PCBFP1                                  % BODY FAT IN SKINFOLD PROCESS I
C     PCBFP2                                  % BODY FAT IN SKINFOLD PROCESS II
C     PCBFWI                                  % BODY FAT IN WATER IMMERSION PROCESS
C     PCF                                       % BODY FAT IN NUTRITION JOURNAL FORMULA
C     PCFFM                                   % FAT FREE MASS IN NUTRITION JOURNAL FORMULA
C     PCTBW                                  % TOTAL BODY WATER , NUTRITION JOURNAL FORMULA
C     RESVOL                                  RESIDUAL VOLUME OF AIR IN LUNGS
C     SEX                                     SUBJECT'S SEX
C     S1                                       LOG OF AVSUM1
C     S2                                       SQUARE OF AVSUM2

```

C	TEMP	TEMPORARY FOR SEX
C	TEMP1	TEMPORARY FOR BSID1
C	TEMP2	TEMPORARY FOR BSID2
C	TRIAL	NUMBER OF UNDER WATER WEIGHT TRIALS
C	UWCHWT	UNDER WATER CHAIR WEIGHT
C	VOL	BODY VOLUME IN NUTRITION JOURNAL FORMULA
C	WATERD	WATER DENSITY
C	WATERT	WATER TEMPERATURE
C	WATRW	WATER WEIGHT

C		
C		
C		
C	ARRAY	USE
C	BICEP(3)	BICEPS SKINFOLD MEASUREMENT
C	BIRTHD(6)	DATE OF BIRTH
C	FATWT(4)	FAT WEIGHT
C	FFMASS(4)	FAT FREE MASS
C	FILENM(11)	SUBJECT'S IDENTIFICATION CODE
C	LBWM(4)	LEAN BODY WATER MASS
C	MEASDT(6)	DATE OF WATER IMMERSION PROCESS
C	MIDCLA(3)	MIDCLAVICULAR SKINFOLD MEASUREMENT
C	MTHIGH(3)	MID THIGH SKINFOLD MEASUREMENT
C	NI TDEN(4)	NITROGEN/KG OF BODY WEIGHT
C	NI TFAT(4)	NITROGEN IN TOTAL FAT TISSUE
C	NI TLBW(4)	NITROGEN IN LEAN BODY WATER
C	PRIDT(6)	DATE OF SKINFOLD PROCESS I
C	PR2DT(6)	DATE OF SKINFOLD PROCESS II
C	SUBSCA(3)	SUBSCAPULAR SKINFOLD MEASUREMENT
C	SUPRAI(3)	SUPRAILIAC SKINFOLD MEASUREMENT
C	TBNIT(4)	TOTAL BODY NITROGEN
C	TRICEP(3)	TRICEPS SKINFOLD MEASUREMENT
C	UMBILI(3)	UMBILICUS SKINFOLD MEASUREMENT
C	UWWT(10)	UNDER WATER WEIGHT
C	VOLFAT(4)	VOLUME OF FAT
C	VOLLBW(4)	VOLUME OF LEAN BODY WATER

NITROGEN WASHOUT PERSONAL DATA INPUT PROGRAM

0002 LOGICAL*1 IDCOD(11), BIRTHD(8), SEX, DATE(8),
 1 KWI, MEASDT(8), KSP1, PRIDT(8),
 2 BSID2, KSP2, PR2DT(8), FILENM(31),
 3 BSID1

0003 INTEGER*2 TRIAL, MEAS1, MEAS2, YORN, AGE

0004 REAL*4 HEIGHT, BFKNOW, NUDEWT,
 1 WATERT, WATERD, RESVOL, MONTH,
 2 UWCHWT, UWWT(10), BICEP(3),
 3 TRICEP(3), SUBSCA(3), SUPRAI(3),
 4 MIDCLA(3), UMBILI(3), MTHIGH(3),
 5 FATWT(5), FFMASS(5), LBWM(5),
 6 VOLFAT(5), VOLLBW(5), NI TFAT(5),
 7 NI TLBW(5), TBNIT(5), NI TDEN(5),BFPC T

ORIGINAL
OF P...

```
C
0005 REAL*8    TEMP, TEMP1, TEMP2
C
C
C    PRINT START MESSAGE
C
0006 WRITE(7,901)
C
C    GET DATA FILENAME
C
0007 101 WRITE(7,902)
0008 READ(5,800, END=101) NF, (FILENM(J),J=1,NF)
0009 IF (NF.LE.0) GO TO 101
0011 FILENM (NF+1)=0
0012 OPEN(UNIT=1, NAME=FILENM, TYPE='NEW', ACCESS='SEQUENTIAL',
1 FORM='FORMATTED', ERR=80, DISP='KEEP',
2 CARRIAGECONTROL='FORTRAN')
0013 REWIND 1
0014 GO TO 90
0015 80 WRITE(7,936) (FILENM(J), J=1,NF)
0016 GO TO 101
C
C    GET DESCRIPTIVE INFORMATION
C
0017 90 AGE=0
0018 SEX=' '
0019 HEIGHT=0.
0020 NUDEWT=0.
0021 BFKNOW=0.
0022 DO 91 J=1,11
0023 91 IDCODE(J)=' '
0024 DO 93 J=1,8
0025 93 BIRTHD(J)=' '
0026 DO 94 J=1,8
0027 94 DATE(J)=' '
0028 102 WRITE(7,903)
0029 READ(5,801, END=102) NC ,(IDCODE(J),J=1,NC)
0030 IF(NC.NE.11) GO TO 102
0032 110 WRITE(7,924)
0033 READ(5,802, END=110) NC ,(DATE(J),J=1,NC)
0034 IF(NC.NE.8) GO TO 110
0036 103 WRITE(7,904)
0037 READ(5,802, END=103) NC ,(BIRTHD(J),J=1,NC)
0038 IF (NC.NE.8) GO TO 103
0040 104 WRITE(7,935)
0041 READ(5,*, END=104) AGE
0042 105 WRITE(7,905)
0043 READ(5,803, END=105) SEX
0044 IF (SEX.NE.'M'.AND.SEX.NE.'F') GO TO 105
0046 IF(SEX.EQ.'M') TEMP='MALE '
0048 IF(SEX.EQ.'F') TEMP='FEMALE'
0050 106 WRITE(7,906)
0051 READ(5,*, END=106) HEIGHT
0052 202 WRITE(7,910)
0053 READ(5,*, END=202) NUDEWT
```

OF POOR QUALITY

```
0054 107 WRITE(7,907)
0055      IF (YORN(0)) 200, 107, 108
0056 108 WRITE(7,908)
0057      READ(5,*, END=108) BFKNOW
0058 200 CONTINUE
0059      WA TERT=0.
0060      WA TERD=0.
0061      RESVOL=0.
0062      UWCHWT=0.
0063      TRIAL=0
0064      DO 230 N=1,10
0065 230 UWWT(N)=0.
0066      DO 210 J=1,8
0067 210 MEASD T(J)= ' '

C
C      WATER IMMERSION
C

0068 201 WRITE(7,909)
0069      KWI=YORN(0)
0070      IF (KWI) 300, 201, 203
0071 203 WRITE(7,911)
0072      READ(5,*, END=203) WA TERT
0073 204 WRITE(7,912)
0074      READ(5,*, END=204) WA TERD
0075 205 WRITE(7,913)
0076      READ(5,*, END=205) RESVOL
0077 206 WRITE(7,914)
0078      READ(5,*, END=206) UWCHWT
0079 207 WRITE(7,915)
0080      READ(5,*, END=207) TRIAL
0081      IF (TRIAL.LT.1.OR. TRIAL.GT.10) GO TO 207
0083      DO 220 N=1, TRIAL
0084 208 WRITE(7,916)N
0085      READ(5,*, END=208) UWWT(N)
0086 220 CONTINUE
0087 221 WRITE(7,917)
0088      READ(5,802, END=221) NC ,(MEASD T(J),J=1,NC)
0089      IF(NC.NE.8)GO TO 221
0091 300 CONTINUE
0092      MEAS1=0
0093      BSIU1=' '
0094      DO 330 N=1,3
0095      BIC EP(N)=0.
0096      TRIC EP(N)=0.
0097      SUBSCA(N)=0.
0098 330 SUPRAI(N)=0.
0099      DO 340 J=1,8
0100 340 PRID T(J)= ' '

C
C      SKINFOLD PROC ESS I
C

0101 301 WRITE(7,918)
0102      KSP1=YORN(0)
0103      IF (KSP1) 400, 301, 302
0104 302 WRITE(7,919)
```

CHANGE LOG
OF FOUR QUARTERS

```

0105 READ(5,803, END=302) BSID1
0106 IF(BSID1.NE.'R'.AND.BSID1.NE.'L') GO TO 302
0108 IF(BSID1.EQ.'R') TEMP1='RIGHT'
0110 IF(BSID1.EQ.'L') TEMP1='LEFT'
0112 303 WRITE(7,920)
0113 READ(5,*, END=303) MEAS1
0114 IF (MEAS1.LT.1.OR.MEAS1.GT.3) GO TO 303
0116 DO 305 N=1,MEAS1
0117 304 WRITE(7,921)N
0118 READ(5,*, END=304) BIC EP(N)
0119 305 CONTINUE
0120 DO 307 N=1,MEAS1
0121 306 WRITE(7,923)N
0122 READ(5,*, END=306) TRIC EP(N)
0123 307 CONTINUE
0124 DO 309 N=1,MEAS1
0125 308 WRITE(7,925)N
0126 READ(5,*, END=308) SUBSCA(N)
0127 309 CONTINUE
0128 DO 320 N=1,MEAS1
0129 310 WRITE(7,927)N
0130 READ(5,*, END=310) SUPRAI(N)
0131 320 CONTINUE
0132 321 WRITE(7,917)
0133 READ(5,802, END=321) NC,(PRIDT(J),J=1,NC)
0134 IF(NC.NE.8) GO TO 321
0136 400 CONTINUE
0137 MEAS2=0
0138 BSID2=' '
0139 DO 430 N=1,3
0140 MIDCLA(N)=0.
0141 UMBILI(N)=0.
0142 430 MTHIGH(N)=0.
0143 DO 440 J=1,8
0144 440 PR2DT(J)=' '

```

C
C
C

SKINFOLD PROCESS II

```

0145 401 WRITE(7,928)
(146 KSP2=YORN(0)
J147 IF (KSP2) 500, 401, 402
0148 402 WRITE(7,919)
0149 READ(5,803, END=402) BSID2
0150 IF(BSID2.NE.'R'.AND.BSID2.NE.'L') GO TO 402
0152 IF(BSID2.EQ.'R') TEMP2='RIGHT'
0154 IF(BSID2.EQ.'L') TEMP2='LEFT'
0156 403 WRITE(7,920)
0157 READ(5,*, END=403)MEAS2
0158 IF (MEAS2.LT.1.OR.MEAS2.GT.3) GO TO 403
0160 DO 405 N=1,MEAS2
0161 404 WRITE(7,930)N
0162 READ(5,*, END=404) MIDCLA(N)
0163 405 CONTINUE
0164 DO 407 N=1,MEAS2
0165 406 WRITE(7,932)N

```

OF PAPER QUALITY

```

0166 READ(5,*, END=406) UMBILI(N)
0167 407 CONTINUE
0168 DO 409 N=1,MEAS2
0169 408 WRITE(7,934)N
0170 READ(5,*, END=408) MTHIGH(N)
0171 409 CONTINUE
0172 410 WRITE(7,917)
0173 READ(5,802, END=410) NC,(PR2DT(J),J=1,NC)
0174 IF(NC.NE.8) GO TO 410
0176 500 CONTINUE

```

```

C
C COMPUTATION FORMULAE
C
C
C BODY COMPOSITION BY WATER IMMERSION
C

```

```

0177 AVUWWT=0.
0178 WA TRWT=0.
0179 BVOLWI=0.
0180 BDWI=0.
0181 PCBFWI=0.
0182 IF (KWI) 600,600,501
0183 501 DO 510 N=1, TRIAL
0184 510 AVUWWT=AVUWWT + UWWT(N)
0185 AVUWWT=AVUWWT/ TRIAL
0186 WA TRWT=AVUWWT - UWCHWT
0187 BVOLWI=( NUDEWT - WA TRWT)/WA TERD) -RESVOL
0188 BDWI=NUDEWT/BVOLWI
0189 PCBFWI=((4.95/BDWI) - 4.50)*100

```

```

C
C BODY COMPOSITION BY SKINFOLD PROCESS 1
C

```

```

0190 600 AVGBIC=0.
0191 AVGTRI=0.
0192 AVGSUB=0.
0193 AVGSUP=0.
0194 AVSUM1=0.
0195 S1=0.
0196 PCBFP1=0.
0197 BDP1=0.
0198 IF (KSP1) 700,700,601
0199 601 DO 603 N=1,MEAS1
0200 603 AVGBIC=AVGBIC + BICEP(N)
0201 AVGBIC=AVGBIC/MEAS1
0202 DO 605 N=1,MEAS1
0203 605 AVGTRI=AVGTRI + TRICEP(N)
0204 AVGTRI=AVGTRI/MEAS1
0205 DO 607 N=1,MEAS1
0206 607 AVGSUB=AVGSUB + SUBSCA(N)
0207 AVGSUB=AVGSUB/MEAS1
0208 DO 609 N=1,MEAS1
0209 609 AVGSUP=AVGSUP + SUPRAI(N)
0210 AVGSUP=AVGSUP/MEAS1
0211 AVSUM1=AVGBIC+AVGTRI+AVGSUB+AVGSUP
0212 S1=ALOG10(AVSUM1)

```

OF POOR QUALITY

0213 IF (SEX.EQ.'M') BDP1=1.1610-.0632*S1
 0215 IF (SEX.EQ.'F') BDP1=1.1581-.0720*S1
 0217 PCBFP1=((4.95/BDP1)-4.50)*100

C
 C BODY COMPOSITION BY SKINFOLD PROCESS 2

C
 0218 700 AVGMCL=0.
 0219 AVGUMB=0.
 0220 AVGMTH=0.
 0221 AVSUM2=0.
 0222 S2=0.
 0223 BDP2=0.
 0224 PCBFP2=0.
 0225 IF (KSP2) 1000,1000,701
 0226 701 DO 703 N=1,MEAS2
 0227 703 AVGMCL=AVGMCL+MIDCLA(N)
 0228 AVGMCL=AVGMCL/MEAS2
 0229 DO 705 N=1,MEAS2
 0230 705 AVGUMB=AVGUMB+UMBILI(N)
 0231 AVGUMB=AVGUMB/MEAS2
 0232 DO 707 N=1,MEAS2
 0233 707 AVGMTH=AVGMTH+MTHIGH(N)
 0234 AVGMTH=AVGMTH/MEAS2
 0235 AVSUM2=AVGMCL+AVGUMB+AVGMTH
 0236 S2=AVSUM2*AVSUM2
 0237 IF (SEX.EQ.'M') BDP2=1.10938-((8.267E-4)*AVSUM2)+((1.6E-6)*S2)-
 1 ((2.574E-4)*AGE)
 0239 IF (SEX.EQ.'F') BDP2=1.105339-((1.1964E-3)*AVSUM2)+((3.8E-6)*S2)-
 1 ((1.069E-4)*AGE)
 0241 PCBFP2=((4.95/BDP2)-4.50)*100

C
 C FORMULA ESTIMATION OF BODY CHARACTERISTICS FROM
 C 1 AMERICAN JOURNAL OF CLINICAL NUTRITION

C
 0242 1000 VOL=0.
 0243 DEN=0.
 0244 PCF=0.
 0245 PCFFM=0.
 0246 PCTBW=0.
 0247 VOL=(1.015*NUDEWT)-4.937
 0248 DEN=NUDEWT/VOL
 0249 PCF=52.67-(2454.38/NUDEWT)
 0250 PCFFM=100-PCF
 0251 PCTBW=34.93+(1811.33/NUDEWT)

C
 C NITROGEN CONTENT CALCULATIONS

C
 0252 1101 DO 1900 N=1,5
 0253 FATWT(N)=0.
 0254 FFMAS(N)=0.
 0255 LBWM(N)=0.
 0256 VOLFAT(N)=0.
 0257 VOLLBW(N)=0.
 0258 NITFAT(N)=0.
 0259 NITLBW(N)=0.

ORIGINAL PAGE IS
OF POOR QUALITY

```

0260      TBNIT(N)=0.
0261      NITDEN(N)=0.
0262      BFPC T=0.
0263      IF(N.EQ.1.AND.BFKNOW.EQ.0.) GO TO 1900
0265      IF(N.EQ.2.AND.PCBFWI.EQ.0.) GO TO 1900
0267      IF(N.EQ.3.AND.PCBFP1.EQ.0.) GO TO 1900
0269      IF(N.EQ.4.AND.PCBFP2.EQ.0.) GO TO 1900
0271      IF(N.EQ.5.AND.PCF.EQ.0.) GO TO 1900
0273      IF(N.EQ.1) BFPC T=BFKNOW
0275      IF(N.EQ.2) BFPC T=PCBFWI
0277      IF(N.EQ.3) BFPC T=PCBFP1
0279      IF(N.EQ.4) BFPC T=PCBFP2
0281      IF(N.EQ.5) BFPC T=PCF
0283      FAWT(N)=BFPC T*NUDEWT/100.
0284      FFMAS(N)=NUDEWT-FAWT(N)
0285      LBWM(N)=FFMAS(N)*.72
0286      VOLFA T(N)=FAWT(N)/.90
0287      VOLLBW(N)=LBWM(N)/.9933
0288      NITFA T(N)=VOLFA T(N)*.067
0289      NITLBW(N)=VOLLBW(N)*.013
0290      TBNIT(N)=NITFA T(N)+NITLBW(N)
0291      NITDEN(N)=TBNIT(N)/NUDEWT
0292      1900 CONTINUE
0293      CALL LINE(0)

C
C      REPORT
C
C
C      PERSONAL CHARACTERISTICS PROFILE
C
0294      1901 WRITE(7,937)
0295      READ(5,804)
0296      WRITE(7,940)
0297      CALL LINE(5)
0298      WRITE(7,941) IXCDE,(FILENM(J),J=1,NF)
0299      CALL LINE(2)
0300      CALL SWAPDT(DATE,MONTH)
0301      WRITE(7,942) DATE(4),DATE(5),DATE(3),MONTH,
1 DATE(6),DATE(7),DATE(8),TEMP,AGE
0302      CALL LINE(2)
0303      IF(BFKNOW.NE.0) GO TO 1920
0305      WRITE(7,943)HEIGHT
0306      CALL LINE(2)
0307      GO TO 1925
0308      1920 WRITE(7,944) HEIGHT,BFKNOW
0309      CALL LINE(2)
0310      CALL SWAPDT(BIRTHD,MONTH)
0311      1925 WRITE(7,938) NUDEWT,BIRTHD(4),BIRTHD(5),BIRTHD(3),
1 MONTH,BIRTHD(6),BIRTHD(7),BIRTHD(8)
0312      CALL LINE(3)

C
C      WATER IMMERSION DATA
C
0313      1930 IF (KWI) 1931, 1931, 1932
0314      1931 WRITE(7,971)

```

ORIGINAL PRINTING
OF POOR QUALITY

```

0315     CALL LINE(5)
0316     GO TO 1945
0317 1932 CALL SWAPDT(MEASDT,MONTH)
0318     WRITE(7,945) NUDEWT, WATERT, RESVOL, WATERT,
1     UWCHWT, MEASDT(4), MEASDT(5), MEASDT(3), MONTH,
2     MEASDT(6), MEASDT(7), MEASDT(8)
0319     CALL LINE(12)
0320     WRITE(7,946) AVUWWT
0321     CALL LINE(4)
0322     N=(TRIAL+1)/2
0323     DO 1940 J=1,N
0324     K1=((J-1)*2)+1
0325     K2=K1+1
0326     IF(K2.LE. TRIAL) WRITE(7,947) K1,UWWT(K1),K2,UWWT(K2)
0328     IF(K2.GT. TRIAL) WRITE(7,947) K1,UWWT(K1)
0330     CALL LINE(2)
0331 1940 CONTINUE
0332     WRITE(7,948) WATRWT, BVOLWI, BDWI, PCBFWI
0333     CALL LINE(4)
0334 1945 CONTINUE
C
C     SKINFOLD PROCESS I DATA
C
0335     IF (KSP1) 1946, 1946, 1947
0336 1946 CALL ESTIM(4)
0337     WRITE(7,972)
0338     CALL LINE(4)
0339     GO TO 1955
0340 1947 CALL ESTIM(14+(2*MEAS1))
0341     CALL SWAPDT(PR1DT,MONTH)
0342     WRITE(7,950) TEMP1, PR1DT(4), PR1DT(5), PR1DT(3), MONTH,
1     PR1DT(6), PR1DT(7), PR1DT(8)
0343     CALL LINE(10)
0344     DO 1950 J=1, MEAS1
0345     WRITE(7,951) J, BIC EP(J), TRIC EP(J), SUBSCA(J),
1     SUPRAI(J)
0346 1950 CALL LINE(2)
0347     WRITE(7,952) AVGBIC, AVGTRI, AVGSUB, AVGSUP
0348     CALL LINE(2)
0349     WRITE(7,953) BDP1, PCBFP1
0350     CALL LINE(2)
C
C     SKINFOLD PROCESS II DATA
C
0351 1955 IF (KSP2) 1956, 1956, 1957
0352 1956 CALL ESTIM(4)
0353     WRITE(7,973)
0354     CALL LINE(4)
0355     GO TO 1965
0356 1957 CALL ESTIM(14+(MEAS2*2))
0357     CALL SWAPDT(PR2DT,MONTH)
0358     WRITE(7,954) TEMP2, PR2DT(4), PR2DT(5), PR2DT(3), MONTH,
1     PR2DT(6), PR2DT(7), PR2DT(8)
0359     CALL LINE(10)
0360     DO 1960 J=1, MEAS2

```

ORIGINAL PAGE IS
OF POOR QUALITY

```
0361 WRITE(7,955) J, MIDCLA(J), UMBILI(J), MTHIGH(J)
0362 1960 CALL LINE(2)
0363 WRITE(7,956) AVGMCL, AVGUMB, AVGMTH
0364 CALL LINE(2)
0365 WRITE(7,953) BDP2, PCBFP2
0366 CALL LINE(2)
C
C NUTRITION JOURNAL FORMULA DATA
C
0367 1965 CALL ESTIM(10)
0368 WRITE(7,960) VOL, DEN, PCF, PCFFM, PC TBW
0369 CALL LINE(10)
C
C NITROGEN CONTENT DATA
C
0370 1970 CALL ESTIM(32)
0371 WRITE(7,961)
0372 CALL LINE(7)
0373 WRITE(7,957) BFKNOW, PCBFWI, PCBFP1, PCBFP2, PCF
0374 CALL LINE(2)
0375 WRITE(7,962) (FATWT(K), K=1,5)
0376 CALL LINE(2)
0377 WRITE(7,963) (FFMASS(K), K=1,5)
0378 CALL LINE(2)
0379 WRITE(7,964) (LBWM(K), K=1,5)
0380 CALL LINE(3)
0381 WRITE(7,965) (VOLFAT(K), K=1,5)
0382 CALL LINE(2)
0383 WRITE(7,966) (VOLLBW(K), K=1,5)
0384 CALL LINE(3)
0385 WRITE(7,967) (NITFAT(K), K=1,5)
0386 CALL LINE(2)
0387 WRITE(7,968) (NITLBW(K), K=1,5)
0388 CALL LINE(3)
0389 WRITE(7,969) (TBNIT(K), K=1,5)
0390 CALL LINE(3)
0391 WRITE(7,970) (NITDEN(K), K=1,5)
0392 CALL LINE(3)
C
C STORE DATA ON DATA FILE
C
0393 WRITE(1,981) IDCODE, DATE, AGE, SEX, BIRTHD,
1 HEIGHT, BFKNOW, NUDEWT
0394 WRITE(1,985) NUDEWT, RESVOL, BVOLWI, BDWI, PCBFWI
0395 WRITE(1,983) BDP1, PCBFP1
0396 WRITE(1,983) BDP2, PCBFP2
0397 WRITE(1,982) VOL, DEN, PCF, PCFFM, PC TBW
0398 DO 2000 M=1,5
0399 WRITE(1,984) FATWT(M), FFMASS(M), LBWM(M), VOLFAT(M),
1 VOLLBW(M), NITFAT(M), NITLBW(M), TBNIT(M), NITDEN(M)
0400 2000 CONTINUE
0401 CLOSE(UNIT=1,DISPOSE='KEEP')
0402 CALL EXIT
0403 901 FORMAT(///10X,'NITROGEN WASHOUT PERSONAL DATA INPUT PROGRAM'///)
0404 902 FORMAT(10X,'FILENAME ON WHICH TO STORE DATA ',
```

ORIGINAL FILED
OF POOR QUALITY

1 /10X,'DEVICE:FILENAME.TYPE = ',,\$)
0405 903 FORMAT(/1X,'IDENTIFICATION NUMBER OF ',
1 'THE SUBJECT(XXX-XX-XXXX): ',,\$)
0406 924 FORMAT(/1X,'DATE (MM-DD-YY) : ',,\$)
0407 904 FORMAT(/1X,'DATE OF BIRTH (MM-DD-YY): ',,\$)
0408 905 FORMAT(/1X,'SEX (M/F): ',,\$)
0409 906 FORMAT(/1X,'HEIGHT IN CM: ',,\$)
0410 907 FORMAT(/1X,'IS % BODY FAT EXPLICITLY KNOWN (Y/N)? ',,\$)
0411 908 FORMAT(/1X,'ENTER % BODY FAT: ',,\$)
0412 909 FORMAT(/1X,'IS BODY COMPOSITION TO BE MEASURED BY WATER ',
1 'IMMERSION (Y/N)? ',,\$)
0413 910 FORMAT(/1X,'NUDE WEIGHT IN KG: ',,\$)
0414 911 FORMAT(/3X,'WATER TEMPERATURE IN CENTIGRADE: ',,\$)
0415 912 FORMAT(/3X,'WATER DENSITY IN KG/L: ',,\$)
0416 913 FORMAT(/3X,'RESIDUAL VOLUME IN L: ',,\$)
0417 914 FORMAT(/3X,'WEIGHT OF UNDERWATER CHAIR IN KG: ',,\$)
0418 915 FORMAT(/3X,'ENTER THE NUMBER OF TRIALS FOR UNDERWATER ',
1 'WEIGHT: ',,\$)
0419 916 FORMAT(/5X,'TRIAL ',I2,' UNDERWATER WEIGHT IN KG: ',,\$)
0420 917 FORMAT(/3X,'DATE OF MEASUREMENTS (MM-DD-YY): ',,\$)
0421 918 FORMAT(/1X,'IS BODY COMPOSITION TO BE MEASURED BY ',
1 'SKINFOLD PROCESS I (Y/N)? ',,\$)
0422 919 FORMAT(/3X,'BODY SIDE USED (R/L)? ',,\$)
0423 920 FORMAT(/3X,'ENTER THE NUMBER OF MEASUREMENTS RECORDED: ',,\$)
0424 921 FORMAT(/5X,'MEASUREMENT ',I1,' BICEPS IN MM: ',,\$)
0425 923 FORMAT(/5X,'MEASUREMENT ',I1,' TRICEPS IN MM: ',,\$)
0426 925 FORMAT(/5X,'MEASUREMENT ',I1,' SUBSCAPULAR IN MM: ',,\$)
0427 927 FORMAT(/5X,'MEASUREMENT ',I1,' SUPRAILIAC IN MM: ',,\$)
0428 928 FORMAT(/1X,'IS BODY COMPOSITION TO BE MEASURED BY ',
1 'SKINFOLD PROCESS II (Y/N)? ',,\$)
0429 930 FORMAT(/5X,'MEASUREMENT ',I1,' MIDCLAVICULAR IN MM: ',,\$)
0430 932 FORMAT(/5X,'MEASUREMENT ',I1,' UMBILICUS IN MM: ',,\$)
0431 934 FORMAT(/5X,'MEASUREMENT ',I1,' ANTERIOR MID-THIGH IN MM: ',,\$)
0432 999 FORMAT(1X)
0433 935 FORMAT(/1X,'AGE IN YEARS: ',,\$)
0434 936 FORMAT(/1X,'ERROR IN ENTERING FILENAME, FILENAME= ', 30A1)
0435 937 FORMAT(/1X,'ALIGN PAPER TO TOP OF PAGE AND HIT CARRIAGE ',
1 'RETURN'/)
0436 940 FORMAT(/20X,'PERSONAL CHARACTERISTICS PROFILE'/)
0437 941 FORMAT(/2X,'SUBJECT',T11,I1A1,T36,'DATA FILE',T52,30A1)
0438 942 FORMAT(/2X,'DATE',T11,3A1,A3,3A1,T36,A8,T52,'AGE',4X,I2)
0439 944 FORMAT(/2X,'HEIGHT',T11,F6.1,' CM',T36,'% BODY FAT',T51,F5.1,' %')
0440 943 FORMAT(/2X,'HEIGHT',T11,F6.1,' CM',T36,'% BODY FAT',
1 T49,'UNKNOWN')
0441 938 FORMAT(/2X,'NUDE',
1 /2X,'WEIGHT',T10,F7.1,' KG',T36,'DATE OF BIRTH',T52,
2 3A1,A3,3A1)
0442 945 FORMAT(/T23,'BODY COMPOSITION EVALUATIONS',
1 ///1X,22(' '), 'WATER IMMERSION PROCESS ',22(' '),
2 //2X,'NUDE WEIGHT',T26,F7.2,' KG',T40,'WATER TEMPERATURE',T57,
3 F6.2,' C',
4 //2X,'RESIDUAL VOLUME',T26,F7.2,' L',T40,'WATER DENSITY',T57,
5 F6.2,' KG/L',
6 //2X,'UNDERWATER CHAIR WEIGHT',T26,F7.2,' KG',T40,'DATE',T58,
7 3A1,A3,3A1)

ORIGINAL LISTING
OF PCGR QUALITY

0443 946 FORMAT(/10X,'UNDERWATER WEIGHT (KG)',T41,
1 'AVERAGE UNDERWATER WEIGHT =' ,F7.2,' KG' ,
2 //10X,' TRIAL' ,T21,' WEIGHT (KG)' ,T41,' TRIAL' ,T51,' WEIGHT (KG)')

0444 947 FORMAT(/T13,I1,T19,F7.2,T43,I2,T49,F7.2)

0445 948 FORMAT(/2X,' WATER WEIGHT' ,T16,F7.2,' KG' ,T40,' BODY VOLUME' ,
1 T53,F7.2,' L' ,
2 //2X,' BODY DENSITY' ,T16,F8.3,' KG/L' ,T40,' % BODY FAT' ,
3 T53,F7.2,' %')

0446 950 FORMAT(///1X,22('-')),' SKINFOLD PROCESS I ' ,22('-') ,
1 //2X,' BODY SIDE USED' ,2X,A5,T48,' DA TE' ,2X ,
2 3A1,A3,3A1
3 //2X,' TRIAL' ,T21,' MEASUREMENTS IN MM' ,
4 //T11,' BIC EPS' ,T25,' TRIC EPS' ,T38,' SUBSCAPULAR' ,T54 ,
5 'SUPRAILIAC')

0447 951 FORMAT(/T5,I1,T11,F5.2,T26,F5.2,T41,F5.2,T56,F5.2)

0448 952 FORMAT(/2X,' MEAN' ,T11,F5.2,T26,F5.2,T41,F5.2,T56,F5.2)

0449 953 FORMAT(/2X,' BODY DENSITY' ,T17,F7.3,' KG/L' ,T48 ,
1 '% BODY FAT' ,T60,F7.2,' %')

0450 954 FORMAT(///1X,22('-')),' SKINFOLD PROCESS II ' ,22('-') ,
1 //2X,' BODY SIDE USED' ,2X,A5,T48,' DA TE' ,2X,3A1,A3,3A1 ,
2 //2X,' TRIAL' ,T21,' MEASUREMENTS IN MM' ,
3 //T7,' MIDCLAVICULAR' ,T24,' UMBILICUS' ,T37 ,
4 'ANTERIOR MID-THIGH')

0451 955 FORMAT(/T5,I1,T11,F5.2,T26,F5.2,T41,F5.2)

0452 956 FORMAT(/2X,' MEAN' ,T11,F5.2,T26,F5.2,T41,F5.2)

0453 957 FORMAT(/2X,' % BODY FAT' ,T26,F7.2,T38,F7.2 ,
1 T49,F7.2,T60,F7.2,T72,F7.2)

0454 960 FORMAT(///1X,22('-')),' NUTRITION JOURNAL FORMULA ' ,
1 22('-') ,
2 //2X,' BODY VOLUME' ,T23,F7.2,' L' ,T44,' BODY DENSITY' ,
3 T60,F8.3,' KG/L' ,
4 //2X,' % BODY FAT' ,T23,F7.2,' %' ,T44,' % FAT FREE MASS' ,
5 T60,F7.2,' %' ,
6 //2X,' % TOTAL BODY WATER' ,T23,F7.2,' %')

0455 961 FORMAT(///1X,22('-')),' NITROGEN CONTENT DATA ' ,32('-') ,
1 //T28,' KNOWN' ,T39,' WATER' ,T48,' SKINFOLD' ,T59,' SKINFOLD' ,T71 ,
2 'NUTRITION' ,
3 /T25,' % BODY FAT' ,T37,' IMMERSION' ,T48,' PROCESS I' ,
4 T59,' PROCESS II' ,T72,' JOURNAL')

0456 962 FORMAT(/2X,' FAT WEIGHT (KG)' ,T26,F7.2,T38,F7.2,T49,F7.2,T60,F7.2 ,
1 T72,F7.2)

0457 963 FORMAT(/2X,' FAT FREE MASS (KG)' ,T26,F7.2,T38,F7.2,T49,F7.2,T60 ,
1 F7.2,T72,F7.2)

0458 964 FORMAT(/2X,' LEAN BODY' ,
1 /2X,' WATER MASS (KG)' ,T26,F7.2,T38,F7.2,T49,F7.2,T60 ,
2 F7.2,T72,F7.2)

0459 965 FORMAT(/2X,' VOLUME OF FAT (L)' ,T26,F7.2,T38,F7.2,T49,F7.2,T60 ,
1 F7.2,T72,F7.2)

0460 966 FORMAT(/2X,' VOLUME LEAN' ,
1 /2X,' BODY WATER (L)' ,T26,F7.2,T38,F7.2,T49,F7.2,T60 ,
2 F7.2,T72,F7.2)

0461 967 FORMAT(/2X,' NITROGEN IN FAT (L)' ,T26,F8.3,T38,F8.3,T49,F8.3,T60 ,
1 F8.3,T72,F8.3)

0462 968 FORMAT(/2X,' NITROGEN IN' ,
1 /2X,' LEAN BODY WATER (L)' ,T26,F8.3,T38,F8.3,T49,F8.3,T60 ,

ORIG: 11/13/63
OF POU: 11/13/63

0463 2 F8.3,T72,F8.3)
969 FORMAT(/2X,'TOTAL BODY',
1 /2X,'NITROGEN (L)',T26,F8.3,T38,F8.3,T49,F8.3,T60,
2 F8.3,T72,F8.3)
0464 970 FORMAT(/2X,'NITROGEN/KG',
1 /2X,'BODY WEIGHT (L)',T26,F9.4,T38,F9.4,T49,F9.4,T60,F9.4,
2 T72,F9.4)
0465 971 FORMAT(//T23,'BODY COMPOSITION EVALUATIONS',
1 //1X,22('-'),'WATER IMMERSION PROCESS NOT DONE',22('-'))
0466 972 FORMAT(///1X,22('-'),'SKINFOLD PROCESS I NOT DONE',22('-'))
0467 973 FORMAT(///1X,22('-'),'SKINFOLD PROCESS II NOT DONE',22('-'))
0468 981 FORMAT(11A1,8A1,I2,A1,8A1,F6.1,F5.1,F7.2)
0469 982 FORMAT(F7.2,F7.3,F7.2,F7.2,F7.2)
0470 983 FORMAT(F7.3,F7.2)
0471 984 FORMAT(F7.2,F7.2,F7.2,F7.2,F7.2,F7.3,F7.3,F7.3,F7.4)
0472 985 FORMAT(F7.2,F7.2,F7.2,F7.3,F7.2)
0473 800 FORMAT(Q,30A1)
0474 801 FORMAT(Q,11A1)
0475 802 FORMAT(Q,8A1)
0476 803 FORMAT(A1)
0477 804 FORMAT(1X;
0478 END

ORIGINAL PRINT IS
OF POOR QUALITY

C
C

FORTRAN IV Storage Map for Program Unit PDINPT
Local Variables, .PSECT \$DATA, Size = 001326 (363. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
AGE	I*2	001040	AVGBIC	R*4	001172	AVGMCL	R*4	001232
AVGMTH	R*4	001242	AVGSUB	R*4	001202	AVGSUP	R*4	001206
AVGTRI	R*4	001176	AVGUMB	R*4	001236	AVSUM1	R*4	001212
AVSUM2	R*4	001246	AVUWW1	R*4	001146	BDP1	R*4	001226
BDP2	R*4	001256	BDWI	R*4	001162	BFKNOW	R*4	001046
BFPC T	R*4	001102	BSID1	L*1	001031	BSID2	L*1	001027
BVOLWI	R*4	001156	DEN	R*4	001272	HEIGHT	R*4	001042
J	I*2	001140	K	I*2	001316	KSP1	L*1	001026
KSP2	L*1	001030	KWI	L*1	001025	K1	I*2	001312
K2	I*2	001314	M	I*2	001320	MEAS1	I*2	001034
MEAS2	I*2	001036	MONTH	R*4	001072	N	I*2	001144
NC	I*2	001142	NF	I*2	001136	NUDEWT	R*4	001052
PCBFP1	R*4	001222	PCBFP2	R*4	001262	PCBFWI	R*4	001166
PCF	R*4	001276	PCFFM	R*4	001302	PCTBW	R*4	001306
RESVOL	R*4	001066	SEX	L*1	001024	S1	R*4	001216
S2	R*4	001252	TEMP	R*8	001106	TEMP1	R*8	001116
TEMP2	R*8	001126	TRIAL	I*2	001032	UWCHWT	R*4	001076
VOL	R*4	001266	WATERD	R*4	001062	WATERT	R*4	001056
WATRWT	R*4	001152						

Local and COMMON Arrays:

Name	Type	Section	Offset	-----Size-----	Dimensions
BIC EP	R*4	\$DATA	000172	000014 (6.)	(3)
BIRTHD	L*1	\$DATA	000013	000010 (4.)	(8)
DATE	L*1	\$DATA	000023	000010 (4.)	(8)
FATWT	R*4	\$DATA	000316	000024 (10.)	(5)
FFMASS	R*4	\$DATA	000342	000024 (10.)	(5)
FILENAME	L*1	\$DATA	000063	000037 (16.)	(31)
IDCODE	L*1	\$DATA	000000	000013 (6.)	(11)

ORIGINAL DOCUMENT
OF PROGRAM UNIT

LBWM	R*4	\$DA TA	000366	000024 (10.) (5)
MEASDT	L*1	\$DA TA	000033	000010 (4.) (8)
MIDCLA	R*4	\$DA TA	000252	000014 (6.) (3)
MTHIGH	R*4	\$DA TA	000302	000014 (6.) (3)
NITDEN	R*4	\$DA TA	000556	000024 (10.) (5)
NITFAT	R*4	\$DA TA	000462	000024 (10.) (5)
NITLBW	R*4	\$DA TA	000506	000024 (10.) (5)
PR1DT	L*1	\$DA TA	000043	000010 (4.) (8)
PR2DT	L*1	\$DA TA	000053	000010 (4.) (8)
SUBSCA	R*4	\$DA TA	000222	000014 (6.) (3)
SUPRAI	R*4	\$DA TA	000236	000014 (6.) (3)
TBNIT	R*4	\$DA TA	000532	000024 (10.) (5)
TRIC EP	R*4	\$DA TA	000206	000014 (6.) (3)
UMBILI	R*4	\$DA TA	000266	000014 (6.) (3)
UWWT	R*4	\$DA TA	000122	000050 (20.) (10)
VOLFAT	R*4	\$DA TA	000412	000024 (10.) (5)
VOLLBW	R*4	\$DA TA	000436	000024 (10.) (5)

FORTRAN IV Storage Map for Program Unit PDINPT
Subroutines, Functions, Statement and Processor-Defined Functions:

Name	Type	Name	Type	Name	Type	Name	Type	Name	Type
ALOG10	R*4	ESTIM	R*4	EXIT	R*4	LINE	I*2	SWAPDT	R*4
YORN	I*2								

ORIGINAL
OF POOR QUALITY

```

0001     INTEGER FUNCTION YORN*2(K)
0002     LOGICAL*1 X
0003     801 FORMAT(A1)
0004     100 READ(5,801, END=200) X
0005     YORN=0
0006     IF (X.EQ.1HN.OR.X.EQ.1HN) YORN=-1
0008     IF (X.EQ.1HY.OR.X.EQ.1HY) YORN=1
0010     IF (YORN.EQ.0) GO TO 200
0012     RETURN
0013     200 WRITE(7,901)
0014     901 FORMAT(' PLEASE ANSWER Y OR N')
0015     GO TO 100
0016     END

```

C

C

FORTRAN IV Storage Map for Program Unit YORN
Local Variables, .PSECT \$DATA, Size = 000005 (3. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
K	I*2 @	000000	X	L*1	000004	YORN	I*2	000002 Eqv

CONFIDENTIAL
OF THE U.S. GOVERNMENT

```
0001      SUBROUTINE LINE(N)  
0002      COMMON/LC,COUNT/LC  
0003      IF (N.EQ.0) LC=0  
0005      IF (N.NE.0) LC=LC+N  
0007      RETURN  
0008      END
```

C
C

FORTRAN IV Storage Map for Program Unit LINE
Local Variables, .PSECT \$DATA, Size = 000002 (1. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
N	I*2	@ 000000						

COMMON Block /LC,COUNT/, Size = 000002 (1. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
LC	I*2	000000						

ORIGINAL
OF PROGRAM

```

0001      SUBROUTINE ESTIM(N)
0002      COMMON/LCOUNT/LC
0003      IF (N+LC.LT.51) RETURN
0005      ND=51-LC
0006      DO 100 I=1,ND
0007      100 WRITE(7,901)
0008      901 FORMAT(1X)
0009      LC=0
0010      RETURN
0011      END

```

FORTRAN IV Storage Map for Program Unit ESTIM
Local Variables, .PSECT \$DATA, Size = 000006 (3. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
I	I*2	000004	N	I*2 @	000000	ND	I*2	000002

COMMON Block /LCOUNT/, Size = 000002 (1. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
LC	I*2	000000						

ORIGINAL PAGE IS
OF POOR QUALITY

```

0001 SUBROUTINE SWAPDT(X,MONTH)
0002 LOGICAL*1 X(8)
0003 REAL*4 MONTH
0004 IF(X(1).EQ.'0'.AND.X(2).EQ.'1') MONTH='JAN'
0006 IF(X(1).EQ.'0'.AND.X(2).EQ.'2') MONTH='FEB'
0008 IF(X(1).EQ.'0'.AND.X(2).EQ.'3') MONTH='MAR'
0010 IF(X(1).EQ.'0'.AND.X(2).EQ.'4') MONTH='APR'
0012 IF(X(1).EQ.'0'.AND.X(2).EQ.'5') MONTH='MAY'
0014 IF(X(1).EQ.'0'.AND.X(2).EQ.'6') MONTH='JUN'
0016 IF(X(1).EQ.'0'.AND.X(2).EQ.'7') MONTH='JUL'
0018 IF(X(1).EQ.'0'.AND.X(2).EQ.'8') MONTH='AUG'
0020 IF(X(1).EQ.'0'.AND.X(2).EQ.'9') MONTH='SEP'
0022 IF(X(1).EQ.'1'.AND.X(2).EQ.'0') MONTH='JAN'
0024 IF(X(1).EQ.'1'.AND.X(2).EQ.'1') MONTH='FEB'
0026 IF(X(1).EQ.'1'.AND.X(2).EQ.'2') MONTH='MAR'
0028 IF(X(1).EQ.'1'.AND.X(2).EQ.'3') MONTH='APR'
0030 IF(X(1).EQ.'1'.AND.X(2).EQ.'4') MONTH='MAY'
0032 IF(X(1).EQ.'1'.AND.X(2).EQ.'5') MONTH='JUN'
0034 IF(X(1).EQ.'1'.AND.X(2).EQ.'6') MONTH='JUL'
0036 IF(X(1).EQ.'1'.AND.X(2).EQ.'7') MONTH='AUG'
0038 IF(X(1).EQ.'1'.AND.X(2).EQ.'8') MONTH='SEP'
0040 IF(X(1).EQ.'1'.AND.X(2).EQ.'9') MONTH='OCT'
0042 IF(X(1).EQ.'1'.AND.X(2).EQ.'0') MONTH='NOV'
0044 IF(X(1).EQ.'1'.AND.X(2).EQ.'1') MONTH='DEC'
0046 RETURN
0047 END

```

FORTRAN IV Storage Map for Program Unit SWAPDT
Local Variables, .PSECT \$DATA, Size = 000004 (2. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
MONTH	R*4	@ 000002						

Local and COMMON Arrays:

Name	Type	Section	Offset	-----Size-----	Dimensions
X	L*1	@ \$DATA	000000	000010 (4.)	(8)

C-2

5.2 DECOMPRESSION DATA PROGRAM

```

CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
C
C *****
C *****
C *
C *
C *              TECHNOLOGY INCORPORATED
C *              LIFE SCIENCES DIVISION
C *    DEPARTMENT OF BIOMATHEMATICS SERVICES
C *
C *
C *****
C *
C *
C *    PROGRAM NAME:.....DECDA T
C *    DESIGNER/ANALYST:.....CRAIG E. LITTON
C *    PROGRAMMER:.....PENNY BUEKER
C *    DATE:.....31 JANUARY 1982
C *
C *-----*
C *
C *
C *    COMPUTER SYSTEM:.....LSI-11
C *    OPERATING SYSTEM:.....RT-11v4
C *
C *-----*
C *
C *
C *    COMPILING SEQUENCE:
C *
C *      LSI:  FORTRAN DECDA T
C *
C *-----*
C *
C *
C *    LINKING SEQUENCE:
C *
C *      LSI:  LINK DECDA T
C *
C *-----*
C *
C *
C *    EXECUTION SEQUENCE:  RUN DECDA T
C *
C *****
C *****
C

```


ORIGINAL LOG
OF POOR QUALITY

```

C
0002 LOGICAL*1 IDCODE(11), TESTDT(8), TIMFD(5), TIMR(5), EXERLV,
      1 BBG10T(5), BBG20T(5), BBG30T(5), BUBMS, BUBPRS,
      2 BBG1RT(5), BBG2RT(5), BBG3RT(5), BBG4RT(5),
      3 BNG10T(5), BNG20T(5), BNG30T(5), BNG40T(5),
      4 BNG1RT(5), BNG2RT(5), BNG3RT(5), BNG4RT(5),
      5 BNDPRS, BNDLC(15), BBG40T(5), FILENM(31)

C
0003 INTEGER*2 TESTNO, TEMP, DURPH(10), YORN, NPH

C
0004 REAL*4 IPN2, PN2PH(10), FP, MONTH,
      1 BBG1RP, BBG2RP, BBG3RP, BBG4RP,
      2 BNG1RP, BNG2RP, BNG3RP, BNG4RP,
      5 PPT180, PPT240, PPT360,
      5 DR180, DR240, DR360

C
C
C PRINT START MESSAGE
C
C
0005 100 WRITE(7,900)

C
C GET DATA FILE NAME
C
C
0006 50 WRITE(7,901)
0007 READ(5,805, END=50) NF, (FILENM(J), J=1, NF)
0008 IF (NF.LE.0) GO TO 50
0010 FILENM (NF+1)=0
0011 OPEN(UNIT=1, NAME=FILENM, TYPE='NEW', ACCESS='SEQUENTIAL',
      1 FORM='FORMATTED', ERR=80, DISP='KEEP',
      2 CARRIAGECONTROL='FORTRAN')
0012 REWIND 1
0013 GO TO 102
0014 80 WRITE(7,905) (FILENM(J), J=1, NF)
0015 GO TO 50

C
C TEST IDENTIFICATION
C
0016 102 WRITE(7,902)
0017 READ(5,800, END=102) NC, (IDCODE(J), J=1, NC)
0018 IF(NC.NE.11) GO TO 102
0020 103 WRITE(7,903)
0021 READ(5,*, END=103) TESTNO
0022 104 WRITE(7,904)
0023 READ(5,801, END=104) NC, (TESTDT(J), J=1, NC)
0024 IF (NC.NE.8) GO TO 104

C
C N2 WASHOUT
C
0026 DO 112 I=1,10
0027 PN2PH(I)=0.
0028 112 DURPH(I)=0
0029 105 WRITE(7,911)
0030 READ(5,*, END=105) NPH
0031 106 WRITE(7,906)

```

ONE
OF POOR QUALITY

```
0032      READ(5,*,END=106) IPN2
0033      DO 111 I=1,NPH
0034 108 WRITE(7,907)
0035      READ(5,*,END=108) PN2PH(I)
0036 110 WRITE(7,908)
0037 111 READ(5,*,END=110) DURPH(I)
```

C
C
C

CHAMBER RUN

```
0038 116 WRITE(7,912)
0039      READ(5,*,END=116) FP
0040 118 WRITE(7,913)
0041      READ(5,*,END=118) TEMP
0042 120 WRITE(7,914)
0043      READ(5,802,END=120) NC, (TIMFD(J),J=1,NC)
0044      IF(NC.NE.5) GO TO 120
0046 122 WRITE(7,915)
0047      READ(5,802,END=122) NC, (TIMR(J),J=1,NC)
0048      IF(NC.NE.5) GO TO 122
0050 124 WRITE(7,916)
0051      READ(5,803,END=124) EXERLV
0052      IF(EXERLV.NE.'L'.AND.EXERLV.NE.'M'.AND.EXERLV.NE.'H')GO TO 124
```

C
C
C

SYMPTOMS

```
0054      BUBMS=0
0055      BUBPRS=0
0056      DO 125 J=1,5
0057      BUBG10T(J)= ' '
0058      BUBG20T(J)= ' '
0059      BUBG30T(J)= ' '
0060      BUBG40T(J)= ' '
0061      BUBG1RT(J)= ' '
0062      BUBG2RT(J)= ' '
0063      BUBG3RT(J)= ' '
0064 125 BUBG4RT(J)= ' '
0065      BUBG1RP=0.
0066      BUBG2RP=0.
0067      BUBG3RP=0.
0068      BUBG4RP=0.
0069 126 WRITE(7,918)
0070      BUBMS=YORN(0)
0071      IF (BUBMS) 160, 126, 128
0072 128 WRITE(7,919)
0073      BUBPRS=YORN(0)
0074      IF (BUBPRS) 160, 128, 130
0075 130 WRITE(7,920)
0076 132 WRITE(7,921)
0077      READ(5,802,END=132) NC, (BUBG10T(J),J=1,NC)
0078      IF(NC.LE.0) GO TO 132
0080 134 WRITE(7,922)
0081      READ(5,802,END=134) NC, (BUBG20T(J),J=1,NC)
0082      IF(NC.LE.0) GO TO 134
0084 136 WRITE(7,923)
0085      READ(5,802,END=136) NC, (BUBG30T(J),J=1,NC)
```

```

0086      IF(NC.LE.0) GO TO 136
0088 138 WRITE(7,924)
0089      READ(5,802,END=138) NC , (BBG40T(J),J=1,NC)
0090      IF(NC.LE.0) GO TO 138
0092 140 WRITE(7,920)
0093 142 WRITE(7,925)
0094      READ(5,802,END=142) NC , (BBG1RT(J),J=1,NC)
0095      IF(NC.LE.0) GO TO 142
0097 144 WRITE(7,926)
0098      READ(5,802,END=144) NC , (BBG2RT(J),J=1,NC)
0099      IF(NC.LE.0) GO TO 144
0101 146 WRITE(7,927)
0102      READ(5,802,END=146) NC , (BBG3RT(J),J=1,NC)
0103      IF(NC.LE.0) GO TO 146
0105 148 WRITE(7,928)
0106      READ(5,802,END=148) NC , (BBG4RT(J),J=1,NC)
0107      IF(NC.LE.0) GO TO 148
0109 150 WRITE(7,920)
0110 152 WRITE(7,929)
0111      READ(5,*,END=152) BBG1RP
0112 154 WRITE(7,930)
0113      READ(5,*,END=154) BBG2RP
0114 156 WRITE(7,931)
0115      READ(5,*,END=156) BBG3RP
0116 158 WRITE(7,932)
0117      READ(5,*,END=158) BBG4RP
0118 160 BNDPRS=0
0119      DO 162 J=1,15
0120 162 BNDLC(J)=' '
0121      DO 165 J=1,5
0122      BNG10T(J)=' '
0123      BNG20T(J)=' '
0124      BNG30T(J)=' '
0125      BNG40T(J)=' '
0126      BNG1RT(J)=' '
0127      BNG2RT(J)=' '
0128      BNG3RT(J)=' '
0129 165 BNG4RT(J)=' '
0130      BNG1RP=0.
0131      BNG2RP=0.
0132      BNG3RP=0.
0133      BNG4RP=0.
0134 161 WRITE(7,933)
0135      BNDPRS=YORN(0)
0136      IF (BNDPRS) 194, 161, 163
0137 163 WRITE(7,935)
0138      READ(5,804,END=163) NC ,(BNDLC(J),J=1,NC)
0139      IF(NC.LE.0) GO TO 163
0141 164 WRITE(7,934)
0142 166 WRITE(7,921)
0143      READ(5,802,END=166) NC , (BNG10T(J),J=1,NC)
0144      IF(NC.LE.0) GO TO 166
0146 168 WRITE(7,922)
0147      READ(5,802,END=168) NC , (BNG20T(J),J=1,NC)
0148      IF(NC.LE.0) GO TO 168

```

ORIGINAL POINT IS
OF POOR QUALITY

```

0150 170 WRITE(7,923)
0151      READ(5,802,END=170) NC , (BNG30T(J),J=1,NC)
0152      IF(NC.LE.0) GO TO 170
0154 172 WRITE(7,924)
0155      READ(5,802,END=172) NC , (BNG40T(J),J=1,NC)
0156      IF(NC.LE.0) GO TO 172
0158 174 WRITE(7,934)
0159 175 WRITE(7,925)
0160      READ(5,802,END=176) NC , (BNG1RT(J),J=1,NC)
0161      IF(NC.LE.0) GO TO 176
0163 178 WRITE(7,926)
0164      READ(5,802,END=178) NC , (BNG2RT(J),J=1,NC)
0165      IF(NC.LE.0) GO TO 178
0167 180 WRITE(7,927)
0168      READ(5,802,END=180) NC , (BNG3RT(J),J=1,NC)
0169      IF(NC.LE.0) GO TO 180
0171 182 WRITE(7,928)
0172      READ(5,802,END=182) NC , (BNG4RT(J),J=1,NC)
0173      IF(NC.LE.0) GO TO 182
0175 184 WRITE(7,934)
0176 186 WRITE(7,929)
0177      READ(5,*,END=186) BNG1RP
0178 188 WRITE(7,930)
0179      READ(5,*,END=188) BNG2RP
0180 190 WRITE(7,931)
0181      READ(5,*,END=190) BNG3RP
0182 192 WRITE(7,932)
0183      READ(5,*,END=192) BNG4RP
0184 194 CONTINUE

C
C      CALCULATIONS
C

0185      CALL DECRAT (IPN2, NPH,
1      PN2PH, DURPH, FP,
2      PPT180, PPT240, PPT360, DR180,
3      DR240, DR360)

C
C      REPORT
C
C      DECOMPRESSION DATA OUTPUT
C

0186      WRITE(7,937)
0187      READ(5,806)
0188      CALL SWAPDT( TESTDT, MONTH)
0189 200 WRITE(7,940) IDCODE, TESTDT(4), TESTDT(5), TESTDT(3),
1      MONTH, TESTDT(6), TESTDT(7), TESTDT(8)
0190 202 WRITE(7,941)PPT180, PPT240, PPT360
0191 204 WRITE(7,942)DR180, DR240, DR360

C
C      STORE DATA ON FILE
C

0192      WRITE(1,981) IDCODE, TESTNO, TESTDT, NPH, IPN2,
1      FP, TEMP, TIMFD, TIMR, EXERLV, BUBMS, BUBPRS
0193      WRITE(1,986) (PN2PH(J), J=1,NPH)
0194      WRITE(1,987) (DURPH(J), J=1,NPH)

```

ORIGINAL PAGE IS
OF POOR QUALITY

```
0195 WRITE(1,982) BBG10T, BBG20T, BBG30T, BBG40T
0196 WRITE(1,982) BBG1RT, BBG2RT, BBG3RT, BBG4RT
0197 WRITE(1,983) BBG1RP, BBG2RP, BBG3RP, BBG4RP
0198 WRITE(1,984) BNDPRS, BNDLC
0199 WRITE(1,982) BNG10T, BNG20T, BNG30T, BNG40T
0200 WRITE(1,982) BNG1RT, BNG2RT, BNG3RT, BNG4RT
0201 WRITE(1,983) BNG1RP, BNG2RP, BNG3RP, BNG4RP
0202 WRITE(1,985) PPT180, PPT240, PPT360, DR180, DR240, DR360
0203 CLOSE(UNIT=1,DISP='SAVE')
0204 CALL EXIT
0205 900 FORMAT(///10X,'DECOMPRESSION EXPERIENCE DATA PROGRAM'/)
0206 901 FORMAT(10X,'FILENAME ON WHICH TO STORE DATA ',
1 /10X,'DEVICE:FILENAME.TYPE= ',,$)
0207 902 FORMAT(/1X,'SUBJECT IDENTIFICATION NUMBER (XXX-XX-XXXX): ',,$)
0208 903 FORMAT(/1X,'TEST NUMBER (XX): ',,$)
0209 904 FORMAT(/1X,'TEST DATE (MM-DD-YY): ',,$)
0210 905 FORMAT(/1X,'ERROR IN ENTERING FILENAME, FILENAME= ', 30A1)
0211 906 FORMAT(/1X,'INITIAL PN2 (PSI) (XX.XX): ',,$)
0212 907 FORMAT(/1X,'PN2 DURING WASHOUT, PHASE I (PSI) (XX.XX): ',,$)
0213 908 FORMAT(/1X,'DURATION OF WASHOUT, PHASE I (MIN) (XXXX): ',,$)
0214 911 FORMAT(/1X,'HOW MANY PHASES (XX): ',,$)
0215 912 FORMAT(/1X,'FINAL PRESSURE (PSIA) (XX.XX): ',,$)
0216 913 FORMAT(/1X,'TEMPERATURE IN CENTIGRADE (XX): ',,$)
0217 914 FORMAT(/1X,'TIME OF FINAL DECOMPRESSION (HH:MM): ',,$)
0218 915 FORMAT(/1X,'TIME OF RECOMPRESSION (HH:MM): ',,$)
0219 916 FORMAT(/1X,'EXERCISE LEVEL ',
1 /16X,'LOW 400- 800 BTU/HR',
2 /16X,'MED 800-1600 BTU/HR',
3 /16X,'HIGH 1600 BTU/HR AND ABOVE',
4 /46X,'ENTER (L/M/H): ',,$)
0220 918 FORMAT(/1X,'WERE BUBBLES MEASURED (Y/N)? ',,$)
0221 919 FORMAT(/1X,'WERE BUBBLES PRESENT (Y/N)? ',,$)
0222 920 FORMAT(/6X,'GRADE OF BUBBLES : ENTER 0 FOR NOT OBSERVED ')
0223 921 FORMAT(/1X,'TIME OF ONSET GRADE 1 (HH:MM): ',,$)
0224 922 FORMAT(/1X,'TIME OF ONSET GRADE 2 (HH:MM): ',,$)
0225 923 FORMAT(/1X,'TIME OF ONSET GRADE 3 (HH:MM): ',,$)
0226 924 FORMAT(/1X,'TIME OF ONSET GRADE 4 (HH:MM): ',,$)
0227 925 FORMAT(/1X,'TIME OF REMISSION GRADE 1 (HH:MM): ',,$)
0228 926 FORMAT(/1X,'TIME OF REMISSION GRADE 2 (HH:MM): ',,$)
0229 927 FORMAT(/1X,'TIME OF REMISSION GRADE 3 (HH:MM): ',,$)
0230 928 FORMAT(/1X,'TIME OF REMISSION GRADE 4 (HH:MM): ',,$)
0231 929 FORMAT(/1X,'PRESSURE AT REMISSION GRADE 1 (PSIA) (XX.XX): ',,$)
0232 930 FORMAT(/1X,'PRESSURE AT REMISSION GRADE 2 (PSIA) (XX.XX): ',,$)
0233 931 FORMAT(/1X,'PRESSURE AT REMISSION GRADE 3 (PSIA) (XX.XX): ',,$)
0234 932 FORMAT(/1X,'PRESSURE AT REMISSION GRADE 4 (PSIA) (XX.XX): ',,$)
0235 934 FORMAT(/6X,'GRADE OF BENDS : ENTER 0 FOR NOT OBSERVED ')
0236 933 FORMAT(/1X,'WERE BENDS PRESENT (Y/N)? ',,$)
0237 935 FORMAT(/1X,'BENDS LOCATION (TEXT): ',,$)
0238 937 FORMAT(/1X,'ALIGN PAPER TO TOP OF PAGE ',
1 'AND HIT CARRIAGE RETURN')
0239 940 FORMAT(//5X,'DECOMPRESSION DATA OUTPUT',
1 //1X,'SUBJECT #',T12,11A1,T26,'DATE',T32,
2 3A1,A3,3A1)
0240 941 FORMAT(/1X,'NITROGEN PRESSURE AT TIME OF FINAL DECOMPRESSION',
1 //10X,'FOR 180 MIN TISSUE : ',F6.2,' PSI',
```

ORIGINAL PAGE IS
OF POOR QUALITY

```
      2 //10X,'FOR 240 MIN TISSUE : ',F6.2,' PSI',  
      3 //10X,'FOR 360 MIN TISSUE : ',F6.2,' PSI')  
0241 942 FORMAT(/1X,'DECOMPRESSION RATIO AT TIME OF FINAL DECOMPRESSION',  
      1 //10X,'FOR 180 MIN TISSUE : ',F6.2,  
      2 //10X,'FOR 240 MIN TISSUE : ',F6.2,  
      3 //10X,'FOR 360 MIN TISSUE : ',F6.2)  
0242 981 FORMAT(11A1,I2,8A1,I2,2F5.2,I2,2(5A1),3(A1))  
0243 982 FORMAT(4(5A1))  
0244 983 FORMAT(4(F4.1))  
0245 984 FORMAT(A1,15A1)  
0246 985 FORMAT(3(F5.2),3(F4.2))  
0247 986 FORMAT(10F5.2)  
0248 987 FORMAT(10I4)  
0249 800 FORMAT(Q,11A1)  
0250 801 FORMAT(Q,8A1)  
0251 802 FORMAT(Q,5A1)  
0252 803 FORMAT(A1)  
0253 804 FORMAT(Q,70A1)  
0254 805 FORMAT(Q,30A1)  
0255 806 FORMAT(1X)  
0256      END
```

C
C

ORIGINAL COPY IS
OF POOR QUALITY

Local Variables, .PSECT \$DATA, Size = 000612 (197. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
BBG1RP	R*4	000512	BBG2RP	R*4	000516	BBG3RP	R*4	000522
BBG4RP	R*4	000526	BNDPRS	L*1	000467	BNG1RP	R*4	000532
BNG2RP	R*4	000536	BNG3RP	R*4	000542	BNG4RP	R*4	000546
BUBMS	L*1	000465	BUBPRS	L*1	000466	DR180	R*4	000566
DR240	R*4	000572	DR360	R*4	000576	EXERLV	L*1	000464
FP	R*4	000502	I	I*2	000610	IPN2	R*4	000476
J	I*2	000604	MONTH	R*4	000506	NC	I*2	000606
NF	I*2	000602	NPH	I*2	000474	PPT180	R*4	000552
PPT240	R*4	000556	PPT360	R*4	000562	TEMP	I*2	000472
TESTNO	I*2	000470						

Local and COMMON Arrays:

Name	Type	Section	Offset	-----Size-----	Dimensions
BBG10T	L*1	\$DA TA	000035	000005 (3.)	(5)
BBG1RT	L*1	\$DA TA	000054	000005 (3.)	(5)
BBG20T	L*1	\$DA TA	000042	000005 (3.)	(5)
BBG2RT	L*1	\$DA TA	000061	000005 (3.)	(5)
BBG30T	L*1	\$DA TA	000047	000005 (3.)	(5)
BBG3RT	L*1	\$DA TA	000066	000005 (3.)	(5)
BBG40T	L*1	\$DA TA	000167	000005 (3.)	(5)
BBG4RT	L*1	\$DA TA	000073	000005 (3.)	(5)
BNDLC	L*1	\$DA TA	000150	000017 (8.)	(15)
BNG10T	L*1	\$DA TA	000100	000005 (3.)	(5)
BNG1RT	L*1	\$DA TA	000124	000005 (3.)	(5)
BNG20T	L*1	\$DA TA	000105	000005 (3.)	(5)
BNG2RT	L*1	\$DA TA	000131	000005 (3.)	(5)
BNG30T	L*1	\$DA TA	000112	000005 (3.)	(5)
BNG3RT	L*1	\$DA TA	000136	000005 (3.)	(5)
BNG40T	L*1	\$DA TA	000117	000005 (3.)	(5)
BNG4RT	L*1	\$DA TA	000143	000005 (3.)	(5)
DURPH	I*2	\$DA TA	000234	000024 (10.)	(10)
FILENM	L*1	\$DA TA	000174	000037 (16.)	(31)
IDCODE	L*1	\$DA TA	000000	000013 (6.)	(11)
PNZPH	R*4	\$DA TA	000260	000050 (20.)	(10)
TESTD	L*1	\$DA TA	000013	000010 (4.)	(8)
TIMFD	L*1	\$DA TA	000023	000005 (3.)	(5)
TIMR	L*1	\$DA TA	000030	000005 (3.)	(5)

Subroutines, Functions, Statement and Processor-Defined Functions:

Name	Type	Name	Type	Name	Type	Name	Type	Name	Type

DEC RAT R*4 EXIT R*4 SWAPDT R*4 YORN I*2

OF POOR QUALITY

```

0001      SUBROUTINE DECRAT (IPN2, NPH, PN2, DUR,
1      FP, PPT180, PPT240, PPT360, DR180, DR240, DR360)
      TABLE OF VARIABLES
C
C
C      VARIABLE              USE
C      DR180                DECOMPRESSION RATIO AT 180 MIN
C      DR240                DECOMPRESSION RATIO AT 240 MIN
C      DR360                DECOMPRESSION RATIO AT 360 MIN
C      DUR                  DURATION OF WASHOUT PHASE
C      FP                  FINAL PRESSURE
C      IPN2                 INITIAL PRESSURE OF N2
C      KT180                VALUE OF EXPONENT IN CALCULATION
C      KT240                VALUE OF EXPONENT IN CALCULATION
C      KT360                VALUE OF EXPONENT IN CALCULATION
C      NPH                  NUMBER OF PHASES
C      PN2                  PRESSURE OF N2 DURING WASHOUT PHASE
C      PP                   TEMPORAY FOR NEW INITIAL PRESSURE
C      PPT180               PRESSURE OF N2 AT 180 MIN
C      PPT240               PRESSURE OF N2 AT 240 MIN
C      PPT360               PRESSURE OF N2 AT 360 MIN
C
0002      INTEGER*2 DUR(10), NPH
0003      REAL*4 IPN2, PN2(10),
1      FP, PPT180, PPT240, PPT360, DR180, DR240, DR360,
2      KT180, KT240, KT360, PP
C
0004      PP=IPN2
0005      DO 1100 I=1,10
0006      1001 KT180=(DUR(I)/180.)*0.693
0007      PPT180=PP+((PN2(I)-PP)*(1.-EXP(-KT180)))
0008      IF(I.EQ.NPH) GO TO 1200
0010      PP=PPT180
0011      1100 CONTINUE
0012      1200 DR180=PPT180/FP
C
0013      PP=IPN2
0014      DO 2100 I=1,10
0015      2000 KT240=(DUR(I)/240.)*0.693
0016      PPT240=PP+((PN2(I)-PP)*(1.-EXP(-KT240)))
0017      IF(I.EQ.NPH) GO TO 2200
0019      PP=PPT240
0020      2100 CONTINUE
0021      2200 DR240=PPT240/FP
C
0022      PP=IPN2
0023      DO 3100 I=1,10
0024      3000 KT360=(DUR(I)/360.)*0.693
0025      PPT360=PP+((PN2(I)-PP)*(1.-EXP(-KT360)))
0026      IF(I.EQ.NPH) GO TO 3200
0028      PP=PPT360
0029      3100 CONTINUE
0030      3200 DR360=PPT360/FP
C
0031      RETURN
0032      END

```

FORTRAN IV Storage Map for Program Unit DECRA T
Local Variables, .PSECT \$DATA, Size = 000100 (32. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
DR180	R*4	@ 000020	DR240	R*4	@ 000022	DR360	R*4	@ 000024
FP	R*4	@ 000010	I	I*2	000062	IPN2	R*4	@ 000000
KT180	R*4	000042	KT240	R*4	000046	KT360	R*4	000052
NPH	I*2	@ 000002	PP	R*4	000056	PPT180	R*4	@ 000012
PPT240	R*4	@ 000014	PPT360	R*4	@ 000016			

Local and COMMON Arrays:

Name	Type	Section	Offset	-----Size-----	Dimensions
DUR	I*2	@ \$DATA	000006	000024 (10.)	(10)
PN2	R*4	@ \$DATA	000004	000050 (20.)	(10)

Subroutines, Functions, Statement and Processor-Defined Functions:

Name	Type	Name	Type	Name	Type	Name	Type	Name	Type
EXP	R*4								

ORIGINAL PAGE IS
OF POOR QUALITY

```
0001      INTEGER FUNCTION YORN*2(K)
0002      LOGICAL*1 X
0003      801 FORMAT(A1)
0004      100 READ(5,801,FND=200) X
0005      YORN=C
0006      IF (X.EQ.1HN.OR.X.EQ.1HN) YORN=-1
0008      IF (X.EQ.1HY.OR.X.EQ.1HY) YORN=1
0010      IF (YORN.EQ.0) GO TO 200
0012      RETURN
0013      200 WRITE(7,901)
0014      901 FORMAT(' PLEASE ANSWER Y OR N')
0015      GO TO 100
0016      END
```

C

C

Local Variables, .PSECT \$DATA, Size = 000005 (3. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
K	I*2	000000	X	L*1	000004	YORN	I*2	000002 Eqv

```

0001 SUBROUTINE SWAPD T(X,MONTH)
0002 LOGICAL*1 X(8)
0003 REAL*4 MONTH
0004 IF(X(1).EQ.'0'.AND.X(2).EQ.'1') MONTH='JAN'
0006 IF(X(1).EQ.'0'.AND.X(2).EQ.'2') MONTH='FEB'
0008 IF(X(1).EQ.'0'.AND.X(2).EQ.'3') MONTH='MAR'
0010 IF(X(1).EQ.'0'.AND.X(2).EQ.'4') MONTH='APR'
0012 IF(X(1).EQ.'0'.AND.X(2).EQ.'5') MONTH='MAY'
0014 IF(X(1).EQ.'0'.AND.X(2).EQ.'6') MONTH='JUN'
0016 IF(X(1).EQ.'0'.AND.X(2).EQ.'7') MONTH='JUL'
0018 IF(X(1).EQ.'0'.AND.X(2).EQ.'8') MONTH='AUG'
0020 IF(X(1).EQ.'0'.AND.X(2).EQ.'9') MONTH='SEP'
0022 IF(X(1).EQ.'1'.AND.X(2).EQ.'1') MONTH='JAN'
0024 IF(X(1).EQ.'1'.AND.X(2).EQ.'2') MONTH='FEB'
0026 IF(X(1).EQ.'1'.AND.X(2).EQ.'3') MONTH='MAR'
0028 IF(X(1).EQ.'1'.AND.X(2).EQ.'4') MONTH='APR'
0030 IF(X(1).EQ.'1'.AND.X(2).EQ.'5') MONTH='MAY'
0032 IF(X(1).EQ.'1'.AND.X(2).EQ.'6') MONTH='JUN'
0034 IF(X(1).EQ.'1'.AND.X(2).EQ.'7') MONTH='JUL'
0036 IF(X(1).EQ.'1'.AND.X(2).EQ.'8') MONTH='AUG'
0038 IF(X(1).EQ.'1'.AND.X(2).EQ.'9') MONTH='SEP'
0040 IF(X(1).EQ.'1'.AND.X(2).EQ.'0') MONTH='OCT'
0042 IF(X(1).EQ.'1'.AND.X(2).EQ.'1') MONTH='NOV'
0044 IF(X(1).EQ.'1'.AND.X(2).EQ.'2') MONTH='DEC'
0046 RETURN
0047 END

```

Local Variables, .PSECT \$DATA, Size = 000004 (2. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
MONTH	R*4	@ 000002						

Local and COMMON Arrays:

Name	Type	Section	Offset	-----Size-----	Dimensions
X	L*1	@ \$DATA	000000	000010 (4.)	(8)

ORIGINAL MANUSCRIPT
OF POOR QUALITY

TABLE OF VARIABLES

C		
C		
C		
C		
C		
C	AGE	SUBJECT'S AGE
C	AMGAS	AMBIENT GAS COMPOSITION
C	AMPRES	AMBIENT PRESSURE (MM HG)
C	BODHYD	BODY HYDRATION
C	BODPOS	BODY POSITION
C	BRGAT	MINIMUM NOISE LEVEL OF SPIROMETER
C	BRN2V	TOTAL NITROGEN VOLUME EXHALED PER BREATH
C	BRTMIX	BREATHING MIXTURE
C	BRTST	MINIMUM VOLUME TO DETERMINE A BREATH
C	BRVOL	TOTAL BREATH VOLUME
C	CALCNT	NOT USED
C	CALPC T	FROM DATA ON VAL CAL GAS CONCENTRATIONS
C	CALZER	VALUE EQUIVALENT TO CALCNT(1,2) FOR ADJUSTING N2 DRIFT BETWEEN BREATHS
C	COMMTS	ARRAY TO HOLD OPTIONAL TEXT OF COMMENTS
C	COSUM	ACCUMULATOR FOR CARBON DIOXIDE END TIDAL CONCENTRATION CALCULATION
C	COTID	AVERAGE CARBON DIOXIDE CONCENTRATION OF BREATH
C	DELAY	FLOATING CALCULATION OF MASS SPEC TRANSPORT LAG
C	DMY	DUMMY ARGUMENT FOR IFFINR SUBROUTINE
C	DRUGS	CODE FOR DRUGS
C	DT	ELAPSED TIME FOR CAL GAS VALUE, PLUMBING AND FLUSHING
C	ELPST	EXPERIMENTAL ELAPSED TIME
C	EXERC I	EXERCISE PROTOCOL
C	EXGRAV	EXTERNAL GRAVITY
C	EXTEMP	EXTERNAL TEMPERATURE IN CENTIGRADE
C	FLB	SPIROMETER CAL CURVE INTERCEPT
C	FLM	SPIROMETER CAL CURVE SLOPE
C	FLCAL	HIGH CALIBRATION VALUE OF SPIROMETER (1 LITER)
C	FLOSUM	ACCUMULATOR FOR SPIROMETER ZERO CALIBRATION
C	FLOZER	LOW CAL VALUE OF SPIROMETER (0 LITERS)
C	GASSLP	MASS SPEC CAL CURVE SLOPE
C	GASSUM	ACCUMULATOR FOR CAL GAS SAMPLE READINGS
C	GSNTCP	MASS SPEC CAL CURVE INTERCEPT
C	I	GAS TYPE INDEX
C	IADBUF	ADDRESS OF A/D OUTPUT BUFFER
C	IADC G	A/D CHANNEL AND GAIN CODES
C	IADC SR	ADDRESS OF A/D CONTROL STATUS REGISTER
C	IANS	OPERATOR RESPONSE CHARACTER (1-6)
C	IBTH	BREATH NUMBER I.D.
C	ICHECK	FLAG TO SET NEW VOLUME FOR BREATH OVER DETERMINATION
C	ICLCNT	PROGRAM CLOCK PRE-SET COUNT BUFFER
C	ICLRAT	INTERRUPT RATE OF SAMPLE PROGRAM (TIMES SPER SECOND)
C	ICON	NUMBER OF CONCENTRATION SAMPLES MADE
C	ID	PARAMETER OF RT-11 SUBROUTINE INTSET

ORIGINAL IN 10
OF PGCN QUALITY

C	IDA T	ALPHA-NUMERIC ARRAY FOR DATE TAG
C	IDCODE	SUBJECT IDENTIFICATION #
C	IDLM	MASS SPEC TRANSPORT LAG 1.1 SEC LAG + F
C	IDMY	INTEGER ZERO END OF DATA TAG (SAMPLE RATE)
C	IDRIFT	VALUE OF NITROGEN IN INSPIRED OXYGEN
C	IDRVIN	PARALLEL BOARD INPUT REGISTER ADDRESS
C	IDRVOT	PARALLEL BOARD OUTPUT REGISTER ADDRESS
C	IFIL	ALPHA-NUMERIC ARRAY FOR DATA FILE NAMING
C	IMIN	EXPERIMENT MINUTE COUNT
C	IMMERS	IMMERSION? (Y/N)
C	INAM	ALPHA-NUMERIC ARRAY FOR NAME TAG
C	INO	OPERATOR RESPONSE NO (FOR INSTRUCTIONS N, CARRIAGE RETURN)
C	INT	INTGER RESULT OF A/D CONVERSION OR PARALLEL BOARD READING
C	INTSET	PROCESSOR DEFINED SUBROUTINE
C	IPHS	PHASE NUMBER ENTERED BY OPERATOR
C	ISTP	STOP RESPONSE ENTERED BY OPERATOR (S, CARRIAGE RETURN)
C	IVAL	OPERATOR INPUT OF DATA CONSTANTS (INTEGERS)
C	IVLVM	BIT STRUCTURE TO CONTROL GAS VALUES
C	IYS	OPERATOR RESPONSE YES (Y, CARRIAGE RETURN)
C	J	CAL GAS SOURCE INDEX
C	KC	NUMBER OF CHARACTERS IN COMMENT LINE
C	KCMMS	COMMENTS? (Y/N)
C	LATCH	LOGICAL FLAG FOR NITROGEN RANGE SWITCHING
C	LBTH	LOGICAL LATCH FOR BREATH DATA GATHERING DONE
C	LWAT	LOGICAL LATCH OF DATA ANALYSIS DONE F -> DATA DONE T -> ANALYSIS IN PROGRESS
C	M	DO LOOP DUMMY COUNTER
C	NTSUM	ACCUMULATOR FOR NITROGEN END TIDAL CONCENTRATION CALCULATION
C	N2TID	AVERAGE NITROGEN CONCENTRATION OF BREATH
C	N2TOT	TOTAL VOLUME OF NITROGEN ELIMINATED
C	N2VLM	MINUTE ACCUMULATOR VOLUME OF NITROGEN ELIMINATED
C	OTSUM	ACCUMULATOR FOR OXYGEN END TIDAL CONCENTRATION CALCULATION
C	O2TID	AVERAGE OXYGEN CONCENTRATION OF BREATH
C	PCCO	CONCENTRATION READING OF CARBON DIOXIDE
C	PCN2	CONCENTRATION READING OF NITROGEN
C	PCO2	CONCENTRATION READING OF OXYGEN
C	RDMY	FLOATING POINT ZERO END OF DATA TAG
C	SPCOND	OTHER SPECIAL CONDITIONS
C	STRTIM	STARTING TIME OF DAY
C	T	TIME FOR DELAY OF CAL GAS PLUMBING, VALVING, PURGING
C	TEST	CHANGE IN BREATH VOLUME TO DETERMINE END OF BREATH
C	TICK	EXPERIMENTAL TEST 1-MINUTE COUNTER
C	TINT	TIME INTERVAL (DT) INTEGRATOR (VOL=DV/DT)
C	TUCK	EXPERIMENTAL TEST START TIME (FOR ELAPSED TIME)
C	TSEC	CPU MONITOR TIME
C	VALUE	OPERATOR INPUT OF DATA CONSTANTS (FLOATING POINT)
C	VOLUME	BREATH VOLUME = BRVOL

ORIGINAL PAGE IS
OF POOR QUALITY

```

C      VOL1      SAMPLE 1 OF VOLUME TO COMPARE TO 5TH SAMPLE
C      TO DETERMINE IF BREATH IS COMPLETE
C      VOL2      SAMPLE 2 OF VOLUME TO COMPARE TO 5TH SAMPLE
C      TO DETERMINE IF BREATH IS COMPLETE
C
C      TYPE
0002     IMPLICIT LOGICAL(L),REAL(N)
0003     LOGICAL*1 IANS,IYS,INO,ISTP,IDCODE(11),IDAT(8),
1     IMMERS,STRTIM(5),KCMMS,COMMS(4,80),
2     INAM(30),IFIL(31)
0004     REAL*4 EXTEMP,AMPRES
0005     INTEGER*2 AMGAS,BODHYD,BODPOS,BRTMIX,
1     DRUGS,EXERCI,EXGRAV,SFCOND,KC,AGE
C
C      SIZE
0006     DIMENSION CALCNT(3,4),CALPC T(3,3)
0007     DIMENSION GASSUM(3),GASSLP(4),GSNTP(4)
0008     DIMENSION IVLVMD(3),IADC GC(4)
0009     DIMENSION PC TN2(300),PC TCO(300),PC T02(300)
C
C      COMMON
0010     COMMON/FLOCOM/FLOZER,FLM,FLB,IFLO,BRTST
0011     COMMON/ADCON/IADSR,IADBUF,IADC GC,VOLUME,VOL1,VOL2
0012     COMMON/GAS/GASSLP,GSNTP,LBTH,BRGAT,TEST
0013     COMMON/PARPRS/ICON,PC TCO,PC TN2,PC T02,CALZER
0014     COMMON/DRV11/IDRVIN,IDRVOT,IDLML,LWAT,IDRIFT
C
0015     EXTERNAL SAMPLE
C
C      DATA
0016     DATA IDRVT/"167772/,IDRVIN/"167774/
0017     DATA IVLVMD/"030,"420,"220/
0018     DATA IYS/1HY/,INO/1HN/,ISTP/1HS/
0019     DATA IADSR/"176770/,IADBUF/"176772/
0020     DATA IADC GC/"0431,"1031,"1431,"2021/
C
C      READ DATA CONSTANTS
0021     OPEN (UNIT=2,TYPE='OLD',NAME='SY:DATCON.VAL',READONLY)
0022     READ (2,*) ((CALPC T(I,J),J=1,3),I=1,3),ICLRAT,
*     BRGAT
0023     CLOSE (UNIT=2,DISPOSE='SAVE')
C
C      TYPE HEADER
0024     TYPE 10
0025     ACCEPT 11,IANS
0026     IF(IANS.EQ.INO) GO TO 91
C
C      GET PHASE & BRANCH
0028     90 CONTINUE
0029     TYPE 12
0030     GO TO 92
0031     91 TYPE 15

```

ORIGINAL PAGE IS
OF POOR QUALITY

0032 92 READ(5,13) IPHS
0033 93 IF(IPHS.LT.1.OR.IPHS.GT.6)GO TO 91
0035 94 GO TO (100,200,300,400,500,600) IPHS

C
C
C
C

PHASE 1--DATA CONSTANTS

0036 100 CONTINUE
0037 TYPE 14

C
C

CAL-GAS %'S

0038 DO 110 J=1,3
0039 DO 110 I=1,3
0040 TYPE 16,I,J,CALPC T(I,J)
0041 ACCEPT 11,IANS
0042 IF(IANS.NE.IYS) GO TO 110
0044 TYPE 18
0045 ACCEPT 19,VALUE
0046 CALPC T(I,J)=VALUE
0047 110 CONTINUE

C
C

CLOCK RATE

0048 TYPE 161,ICLRAT
0049 ACCEPT 11,IANS
0050 IF(IANS.NE.IYS) GO TO 115
0052 TYPE 18
0053 ACCEPT 191,IVAL
0054 ICLRAT=IVAL
0055 115 CONTINUE

C

BREATH GATE

0056 TYPE 162,BRGAT
0057 ACCEPT 11,IANS
0058 IF(IANS.NE.IYS) GO TO 120
0060 TYPE 18
0061 ACCEPT 19,VALUE
0062 BRGAT=VALUE
0063 120 CONTINUE

C
C

END OF PHASE 1--GET NEW PHASE

0064 TYPE 20,IPHS
0065 GO TO 92

C
C
C

PHASE 2-- HEADER INFO

0066 200 CONTINUE

C
C

INITIALIZE HEADER VARIABLES

0067 DO 201 M=1,8
0068 201 IDAT(M)=' '
0069 DO 202 M=1,30
0070 202 INAM(M)=' '
0071 DO 203 M=1,11
0072 203 ICODE(M)=' '
0073 EXTEMP=0.0
0074 EXGRAV=0

```

0075     IMMERS=' '
0076     BODPOS=0
0077     AMPRES=0.0
0078     DRUGS=0
0079     EXERC I=0
0080     BRTMIX=0
0081     BODHYD=0
0082     AMGAS=0
0083     AGE=0
0084     DO 204 M=1,5
0085     204 STRTIM(M)=' '
0086     SPCOND=0
0087     KCMMTS=' '
0088     DO 216 J=1,4
0089     DO 216 I=1,80
0090     216 CUMMTS(J,I)=' '

```

C
C

```

ENTER THE EXPERIMENTAL CONDITIONS
0091     TYPE 22
0092     ACCEPT 23,( IDA T(M) ,M=1,8)
0093     TYPE 24
0094     ACCEPT 25, IC ,( INAM(M) ,M=1, IC )
0095     TYPE 26
0096     ACCEPT 27,( IDC ODE(M) ,M=1,11)
0097     TYPE 28
0098     ACCEPT 29, AGE
0099     TYPE 900
0100     ACCEPT 800, EXTEMP
0101     TYPE 901
0102     ACCEPT 801, EXGRAV
0103     TYPE 902
0104     ACCEPT 802, IMMERS
0105     TYPE 903
0106     ACCEPT 803, BODPOS
0107     TYPE 904
0108     ACCEPT 800, AMPRES
0109     TYPE 905
0110     ACCEPT 803, DRUGS
0111     TYPE 906
0112     ACCEPT 803, EXERC I
0113     TYPE 907
0114     ACCEPT 803, BRTMIX
0115     TYPE 908
0116     ACCEPT 803, BODHYD
0117     TYPE 909
0118     ACCEPT 803, AMGAS
0119     TYPE 910
0120     ACCEPT 810, ( S TRTIM(M) ,M=1,5)
0121     TYPE 911
0122     ACCEPT 803, SPCOND
0123     TYPE 912
0124     ACCEPT 802, KCMMTS
0125     IF (KCMMTS.EQ.'N') GO TO 225
0127     TYPE 913
0128     ACCEPT 813

```

ORIGINAL COPY
OF POOR QUALITY

```
0129      DO 220 J=1,4
0130      TYPE 31
0131      READ(5,814) KC, (COMMITTS(J,I), I=1,KC)
0132      220 CONTINUE
C
C          END OF PHASE 2--GET NEW PHASE
0133      225 TYPE 20,IPHS
0134      GO TO 92
C
C
C          PHASE 3--MASS SPEC TROMETER CALIBRATION
0135      300 CONTINUE
0136      CALL IPOKE(IDRVOT,IVLVMD(2))
0137      TYPE 30
C
0138      PAUSE 'SELECT N-2 0-10%. RETURN'
0139      J=2
C          OPEN VALVE J--WAIT 15 SECONDS
0140      305 CALL IPOKE(IDRVOT,IVLVMD(J))
0141      T=SECNDS(0.0)
0142      308 DT=SECNDS(T)
0143      IF(DT.LT.15.0) GO TO 308
C
0145      DO 315,I=1,3
0146      GASSUM(I)=0.0
C
0147      DO 310,M=1,100
0148      CALL IPOKE(IADC SR,IADC GC(I))
0149      INT=IPEEK(IADBUF)
0150      310 GASSUM(I)=GASSUM(I)+FLOA T(INT)
C
0151      CALC NT(I,J)=GASSUM(I)/100.
0152      315 CONTINUE
0153      J=J+1
0154      IF(J.EQ.3) GO TO 305
0156      CALC NT(1,4)=CALC NT(1,3)
C
0157      PAUSE 'SELEC T N-2 0-100%. RETURN'
0158      J=1
C          OPEN VALVE J--WAIT 15 SECONDS
0159      325 CALL IPOKE(IDRVOT,IVLVMD(J))
0160      T=SECNDS(0.0)
0161      328 DT=SECNDS(T)
0162      IF(DT.LT.15.) GO TO 328
0164      GASSUM(1)=0.0
0165      DO 330 M=1,100
0166      CALL IPOKE(IADC SR,IADC GC(1))
0167      INT=IPEEK(IADBUF)
0168      330 GASSUM(1)=GASSUM(1)+FLOA T(INT)
0169      CALC NT(1,J)=GASSUM(1)/100.
0170      J=J+2
0171      IF(J.EQ.3) GO TO 325
C
C
C          GET SLOPES & INTERCEPTS
```

ORIGINAL L
OF FOUR

```

0173      GASSLP(1)=(CALPC T(1,1)-CALPC T(1,3))
* /(CALCNT(1,1)-CALCNT(1,3))
0174      GSNTCP(1)=CALPC T(1,3)-(GASSLP(1)*CALCNT(1,3))
C
0175      GASSLP(2)=(CALPC T(2,2)-CALPC T(2,3))
* /(CALCNT(2,2)-CALCNT(2,3))
0176      GSNTCP(2)=CALPC T(2,3)-(GASSLP(2)*CALCNT(2,3))
C
0177      GASSLP(3)=(CALPC T(3,3)-CALPC T(3,2))
* /(CALCNT(3,3)-CALCNT(3,2))
0178      GSNTCP(3)=CALPC T(3,2)-(GASSLP(3)*CALCNT(3,2))
C
0179      GASSLP(4)=(CALPC T(1,3)-CALPC T(1,2))
* /(CALCNT(1,4)-CALCNT(1,2))
0180      GSNTCP(4)=CALPC T(1,2)-(GASSLP(4)*CALCNT(1,2))
C
0181      TYPE 321
0182      TYPE 32,CALCNT(1,1),CALCNT(1,3),GASSLP(1),GSNTCP(1)
0183      TYPE 322
0184      TYPE 32,CALCNT(1,4),CALCNT(1,2),GASSLP(4),GSNTCP(4)
0185      TYPE 323
0186      TYPE 32,CALCNT(2,2),CALCNT(2,3),GASSLP(2),GSNTCP(2)
0187      TYPE 324
0188      TYPE 32,CALCNT(3,3),CALCNT(3,2),GASSLP(3),GSNTCP(3)
0189      CALZER=CALCNT(1,2)
C
C      END OF PHASE 3
0190      CALL IPOKE(IDRVOT,"0")
0191      TYPE 20,IPHS
0192      GO TO 92
C
C
C      PHASE 4--SPIROMETER CALIBRATION
0193      400 CONTINUE
0194      TYPE 40
C
0195      PAUSE'ZERO SPIROMETER.RETURN THANK YOU'
C      SAMPLE "ZERO FLOW"
0196      FLOSUM=0.0
0197      DO 415,M=1,1000
0198      CALL IPOKE(IADC SR,IADC GC(4))
0199      INT=IPEEK(IADBUF)
0200      415 FLOSUM=FLOSUM+FLOAT(INT)
0201      FLOZER=FLOSUM/1000.
C      SIMULATE 1 LITER VOLUME
0202      FLOCAL=FLOZER+410.5
C      GET FLOW SLOPES & INTERCEPTS
0203      FLM=1.0/(FLOCAL-FLOZER)
0204      FLB=-FLOZER
0205      TYPE 421
0206      TYPE 32,FLOCAL,FLOZER,FLM,FLB
C
C      END OF PHASE 4
0207      TYPE 20,IPHS
0208      GO TO 92

```

ORIGINAL FILED
OF POOR QUALITY

```
C
C
C
0209 PHASE 5--TEST DATA
0210 499 PAUSE 'INSERT DATA FLOPPY WITH AT LEAST 400 FREE BLOCKS, RETURN'
      500 CONTINUE
C
C
0211 OPEN DISK FILE & WRITE HEADER
0212 501 WRITE(7,920)
0213 READ(5,820,END=501) MC, (IFIL(M),M=1,MC)
0214 IF (MC.LE.0) GO TO 501
0215 IFIL(MC+1)=0
0216 OPEN(UNIT=2,TYPE='NEW',NAME=IFIL,FORM='FORMATTED',ERR=499,
      1 INITIALSIZE=400)
0217 WRITE (2,914) (INAM(I),I=1,30),(IDAT(I),I=1,8),
      1 (IDCODE(I),I=1,11), AGE,
      2 EXTEMP,EXGRAV,IMMERS,BODPOS,AMPRES,DRUGS,
      3 EXERC I,BRTMIX,BODHYD,AMGAS,(STRTIM(M),M=1,5),SPCOND,KCMNTS
0218 DO 502 I=1,4
0219 WRITE(2,915)(COMMTS(I,J),J=1,80)
0220 502 CONTINUE
0221 TYPE 505
0222 ACCEPT 11,IANS
0223 IF(IANS.EQ.ISTP) GO TO 560
C
C
0225 PREPARE INTERRUPT ROUTINE 'SAMPLE'
0226 CALL INTSET("104,4,ID,SAMPLE)
0227 ICLCNT=10000/ICLRAT
      CALL IPOKE("172542,ICLCNT)
C
C
0228 INITIALIZE SOME TERMS
0229 ICON=0
0230 IFLO=0
0231 IBTH=0
0232 N2VLM=0.0
0233 BRTST=FLOZER+BRGAT
0234 VOLU=0.0
0235 VOL1=BRTST
0236 VOL2=0.0
0237 ICHECK=0
0238 IDRIFT=20
0239 LWAT=.FALSE.
0240 LATCH=.TRUE.
0241 N2TOT=0.0
      IMIN=1
C
C
0242 SET TIME AND DELAY TERMS DUE TO CATHETER LAG
0243 TINT=1.0/FLOAT(ICLRAT)
0244 DELAY=FLOAT(ICLRAT)*.36
      IDLM=FIX(DELAY)
C
0245 TSEC=SECONDS(0.0)
0246 TOCK=TSEC
C
C
0247 START SAMPLE ROUTINE
      CALL IPOKE("172540,"0113)
```

```

C
C      START TO MONITOR SAMPLE ROUTINE
C
0248 510 CONTINUE
0249     INT=IPEEK("167774)
0250     IF(INT.EQ."177777) LATCH=.FALSE.
0252     IANS=ITINR(DMY)
C
0253     IF(IANS.EQ.ISTP) GO TO 560
C
C      CHECK NEED TO PURGE
0255     IF((VOL1.LE.BRTST).OR.(ABS(TEST).GT..5)) GOTO 520
0257     DO 515,J=1,3000
0258     CALL IPOKE("176760,"2050)
0259     515 CONTINUE
0260     CALL IPOKE("176760,"0010)
C      CHECK FOR DATA READY
0261 520 IF(.NOT.LWAT) GO TO 530
C
C      BREATH DATA IS READY FOR PROCESSING
0263     IBTH=IBTH+1
0264     ELPST=SECNDS(TSEC)
C
C
0265     DO 540,J=2,ICON
C
0266     NTSUM=NTSUM+PC T12(J)
0267     OTSUM=OTSUM+PC T02(J)
0268     COSUM=COSUM+PC T0(J)
C
0269 540 CONTINUE
C
0270     N2TID=NTSUM/ICON
0271     O2TID=OTSUM/ICON
0272     COTID=COSUM/ICON
0273     BRVOL=VOLUME
0274     BRN2V=BRVOL*N2TID
0275     N2VLM=N2VLM+BRN2V
0276     N2TOT=N2TOT+BRN2V
C
0277 542 CONTINUE
C
0278     WRITE(2,916) IBTH,ELPST,BRVOL,BRN2V,O2TID,N2TID,COTID
0279     TYPE 55,IBTH,ELPST,BRVOL,BRN2V,N2TID
0280 55 FORMAT(5X,I4,F10.1,3F10.5)
C      RESET COUNTERS & INTEGRATORS
0281     ICON=0
0282     IFLO=0
0283     ICHECK=0
0284     VOL2=0.0
0285     BRVOL=0.0
0286     BRN2V=0.0
0287     NTSUM=0.0
0288     OTSUM=0.0
0289     COSUM=0.0

```

ORIGINAL PAGE IS
OF POOR QUALITY

```
0290     IDRIFT=20
0291     LWAT=.FALSE.

C
C       CHECK FOR 1-MIN OPERATIONS
0292 530 TICK=SECNDS(TOCK)
0293     IF(TICK.LE.60.) GO TO 510

C
0295     TOCK=TOCK+60.
0296     TYPE 550,IMIN,N2VLM,N2TOT
0297     IMIN=IMIN+1
0298     N2VLM=0.0
0299     IF((N2TID.LT.0.070).AND.(LATCH)) TYPE 545
0301     GO TO 510

C
C       END OF PHASE 5
0302 560 CONTINUE
0303     CALL IPOKE("172540,"000)
0304     CALL IPOKE(IDRVOT,"000)
0305     IDMY=0
0306     RDMY=0.0
0307     WRITE (2,916) IDMY,RDMY,RDMY,RDMY,RDMY,RDMY,RDMY
0308     CLOSE(UNIT=2,DISP='SAVE')
0309     TYPE 20,IPHS
0310     GO TO 92

C
C
C
C
C       PHASE 6-- TERMINATION
0311 600 CONTINUE
0312     STOP 'BYE'

C
C
C
C
0313     9 FORMAT(9F7.2,I5,F5.0,F7.2)
C
0314     10 FORMAT(/,17X,'N 2  W A S H O U T',//1X,
* 'INSTRUCTIONS (Y/N) ? ',$,)
C
0315     11 FORMAT(A1)
C
0316     12 FORMAT(X,//,X,' TESTLOAD N2WASH HAS 6 PHASES',
* //,5X,'1  DATA CONSTANTS--ALLOWS USER',
* /,' MODIFICATION OF CAL-GAS %'S.  OTHERWISE',
* /,' PROGRAM USES DEFAULT VALUES.',
* //,5X,'2  INITIALIZATION--ENTER HEADER INFORMATION',
* //,5X,'3  CALIBRATE MASS SPECTROMETER',
* //,5X,'4  CALIBRATE SPIROMETER',
* //,5X,'5  TEST DATA--SAMPLES TEST DATA CONTINUOUSLY',
* /,' UNTIL USER TERMINATES',
* //,5X,'6  TERMINATE--RETURNS CONTROL TO SYSTEM MONITOR',
* ///,10X,'ENTER PHASE # & RETURN ',$,)
C
0317     13 FORMAT(I1)
C
0318     14 FORMAT(X,//,10X,'PHASE 1  DATA CONSTANTS',
```

```

* /,5X,' TO LEAVE VALUE UNCHANGED, ENTER',
* /,5X,' H & RETURN. TO CHANGE, ENTER Y',
* /,5X,' & RETURN. GET "NEW VALUE" PROMPT,' ,
* /,5X,' ENTER NEW VALUE & RETURN',/)

C
0319 15 FORMAT(5X,'ENTER PHASE #',,$)
0320 16 FORMAT(5X,'CALPC T(',I1,',',I1,')=',F6.2,$)
C
0321 161 FORMAT(5X,'SAMPLE RATE=',I4,'/SEC',,$)
C
0322 162 FORMAT(5X,'BREATH GATE=',F5.0,$)
C
0323 18 FORMAT(3X,'NEW VALUE =',,$)
C
0324 19 FORMAT(F6.2)
C
0325 191 FORMAT(I4)
C
0326 20 FORMAT(X,/,5X,'END OF PHASE ',I3,/1X,
* 'ENTER NEW PHASE # ',,$)
C
0327 22 FORMAT(X,/,10X,'PHASE 2  INITIALIZATION',/,
* 2X,'DATE (MM-DD-YY): ',,$)
C
0328 23 FORMAT(8A1)
C
0329 24 FORMAT(/2X,'SUBJECT NAME : ',,$)
C
0330 25 FORMAT(Q,30A1)
0331 26 FORMAT(/2X,'SUBJECT IDENTIFICATION # (XXX-XX-XXXX): ',,$)
0332 27 FORMAT(11A1)
0333 28 FORMAT(/2X,'SUBJECT AGE (YEARS AS OF THIS DATE): ',,$)
0334 29 FORMAT(I2)
C
0335 30 FORMAT(X,/,10X,'PHASE 3  MASS SPEC CALIBRATION',
* /,5X,'VERIFY CAL GASES ARE OPEN',
* /,5X,'SET MASS SPEC FUNCTION TO "AUTO"')
0336 31 FORMAT(1X,'?',,$)
C
0337 32 FORMAT(5X,'C TS-1',6X,'C TS-2',5X,'SLOPE',
* 7X,'INTCPT',/,2(3X,F8.2),2X,2(3X,F9.3))
C
0338 321 FORMAT(/,10X,'NITROGEN, HI-RANGE')
0339 322 FORMAT(/,10X,'NITROGEN, LO-RANGE')
0340 323 FORMAT(/,10X,'OXYGEN')
0341 324 FORMAT(/,10X,'CARBON DIOXIDE')
C
0342 40 FORMAT(X,/,10X,'PHASE 4  SPIROMETER CALIBRATION')
C
0343 421 FORMAT(/,10X,'SPIROMETER')
C
0344 505 FORMAT(5X,'SELECT N-2 0-100% RANGE.',
* /,5X,'RETURN TO START.',
* /,5X,'TO STOP AT ANY TIME, ENTER S & RETURN.',,$)
C

```

OF POOR QUALITY

0345 545 FORMAT(5X, 'SELECT N-2 0-10% RANGE')

C

0346 550 FORMAT(2X, I3, 2X, 'MIN N2(L)', F8.5, 3X, 'N2TOT', 2X, F8.5)

0347 900 FORMAT(/2X, 'EXTERNAL TEMPERATURE IN CENTIGRADE (XX.X): ', \$)

0348 901 FORMAT(/2X, 'EXTERNAL GRAVITY (1,0): ', \$)

0349 902 FORMAT(/2X, 'IS SUBJECT IMMersed (Y/N): ', \$)

0350 903 FORMAT(/2X, 'ENTER CODE FOR BODY POSITION (0 FOR NO DATA ',
* 'OR 1-99): ', \$)

0351 904 FORMAT(/2X, 'AMBIENT PRESSURE IN MM HG (XXX.X): ', \$)

0352 905 FORMAT(/2X, 'ENTER CODE FOR DRUGS (0 FOR NO DATA ',
* 'OR 1-99): ', \$)

0353 906 FORMAT(/2X, 'ENTER CODE FOR EXERCISE PROTOCOL (0 FOR NO DATA ',
* 'OR 1-99): ', \$)

0354 907 FORMAT(/2X, 'ENTER CODE FOR BREATHING MIXTURE (0 FOR NO DATA ',
* 'OR 1-99): ', \$)

0355 908 FORMAT(/2X, 'ENTER CODE FOR BODY HYDRATION (0 FOR NO DATA ',
* 'OR 1-99): ', \$)

0356 909 FORMAT(/2X, 'ENTER CODE FOR AMBIENT GAS COMPOSITION (0 FOR NO ',
* 'DATA OR 1-99): ', \$)

0357 910 FORMAT(/2X, 'STARTING TIME OF DAY (HH:MM): ', \$)

0358 911 FORMAT(/2X, 'ENTER CODE FOR OTHER SPECIAL CONDITIONS (0 FOR ',
* 'NO DATA OR 1-99): ', \$)

0359 912 FORMAT(/2X, 'ARE THERE ANY COMMENTS TO RECORD(Y/N)? ', \$)

0360 913 FORMAT(/2X, 'TO ENTER COMMENTS HIT CARRIAGE RETURN THEN ',
* 'ENTER UP TO 4 LINES OF COMMENT ', \$)

0361 914 FORMAT(30A1, 8A1, 11A1, I2, F7.2, I3, A1, I2, F7.2, 5(I2), 5A1, I2, A1)

0362 915 FORMAT(80A1)

0363 916 FORMAT(I4, F10.1, 5(F10.5))

0364 920 FORMAT(10X, 'FILENAME ON WHICH TO STORE DATA ',
1 /10X, 'DEVICE:FILENAME. TYPE = ', \$)

0365 800 FORMAT(F7.2)

0366 801 FORMAT(I1)

0367 802 FORMAT(A1)

0368 803 FORMAT(I2)

0369 810 FORMAT(5A1)

0370 814 FORMAT(Q, 80A1)

0371 813 FORMAT(1X)

0372 820 FORMAT(Q, 30A1)

0373 END

Local Variables, .PSECT \$DATA, Size = 001326 (363. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
AGE	I*2	001106	AMGAS	I*2	001064	AMPRES	R*4	001060
BODHYD	I*2	001066	BODPOS	I*2	001070	BRN2V	R*4	001270
BRTMIX	I*2	001072	BRVOL	R*4	001264	COSUM	R*4	001244
COTID	R*4	001260	DELAY	R*4	001210	DMY	R*4	001224
DRUGS	I*2	001074	DT	R*4	001136	ELPST	R*4	001230
EXERCI	I*2	001076	EXGRAV	I*2	001100	EXTEMP	R*4	001054
FLOCAL	R*4	001150	FLOSUM	R*4	001144	I	I*2	001110
IANS	L*1	001050	IBTH	I*2	001162	IC	I*2	001130
ICHECK	I*2	001170	ICLCNT	I*2	001160	ICLRAT	I*2	001114
ID	I*2	001156	IDMY	I*2	001300	IMIN	I*2	001202
IMMERS	L*1	001051	INO	L*1	000775	INT	I*2	001142
IPHS	I*2	001116	ISTP	L*1	000776	IVAL	I*2	001124
IYS	L*1	000774	J	I*2	001112	KC	I*2	001104
KCMMS	L*1	001052	LATCH	L*4	001172	M	I*2	001126
MC	I*2	001154	NTSUM	R*4	001234	N2TID	R*4	001250
N2TOT	R*4	001176	N2VLM	R*4	001164	OTSUM	R*4	001240
O2TID	R*4	001254	RDMY	R*4	001302	SPCOND	I*2	001102
T	R*4	001132	TICK	R*4	001274	TINT	R*4	001204
TOCK	R*4	001220	TSEC	R*4	001214	VALUE	R*4	001120

COMMON Block /FLOCOM/, Size = 000022 (9. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
FLOZER	R*4	000000	FLM	R*4	000004	FLB	R*4	000010
IFLO	I*2	000014	BRTST	R*4	000016			

COMMON Block /ADCON/, Size = 000030 (12. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset

COMMON BLOCKS
OF POOR QUALITY

IADC SR I*2 000000 IADBUF I*2 000002 IADC GC I*2 000004

VOLUME R*4 000014 VOL1 R*4 000020 VOL2 R*4 000024

COMMON Block /GAS /, Size = 000054 (22. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
GASSLP	R*4	000000	GSNTCP	R*4	000020	LB TH	L*4	000040
BRGAT	R*4	000044	TEST	R*4	000050			

COMMON Block /PARPRS/, Size = 007026 (1803. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
ICON	I*2	000000	PC TC0	R*4	000002	PC TN2	R*4	002262
PC T02	R*4	004542	CALZER	R*4	007022			

COMMON Block /DRV11 /, Size = 000014 (6. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
IDRVIN	I*2	000000	IDRVOT	I*2	000002	IDLM	I*2	000004
LWAT	L*4	000006	IDRIFT	I*2	000012			

CONTAINS
OF POOR QUALITY

Local and COMMON Arrays:

Name	Type	Section	Offset	-----Size-----	Dimensions
CALCNT	R*4 Vec	\$DA TA	000626	000060 (24.)	(3,4)
CALPC T	R*4 Vec	\$DA TA	000706	000044 (18.)	(3,3)
COMMTS	L*1	\$DA TA	000030	000500 (160.)	(4,80)
GASSLP	R*4	GAS	000000	000020 (8.)	(4)
GASSUM	R*4	\$DA TA	000752	000014 (6.)	(3)
GSN TCP	R*4	GAS	000020	000020 (8.)	(4)
IADC GC	I*2	ADC ON	000004	000010 (4.)	(4)
IDAT	L*1	\$DA TA	000013	000010 (4.)	(8)
IDCODE	L*1	\$DA TA	000000	000013 (6.)	(11)
IFIL	L*1	\$DA TA	000566	000037 (16.)	(31)
INAM	L*1	\$DA TA	000530	000036 (15.)	(30)
IVLVMD	I*2	\$DA TA	000766	000006 (3.)	(3)
PC TCO	R*4	PARPRS	000002	002260 (600.)	(300)
PC TN2	R*4	PARPRS	002262	002260 (600.)	(300)
PC TO2	R*4	PARPRS	004542	002260 (600.)	(300)
S TRTIM	L*1	\$DA TA	000023	000005 (3.)	(5)

Subroutines, Functions, Statement and Processor-Defined Functions:

Name	Type	Name	Type	Name	Type	Name	Type	Name	Type
ABS	R*4	FLOAT	R*4	IFIX	I*2	INTSET	I*2	IPEEK	I*2
IPOKE	I*2	ITTINR	I*2	SAMPLE	R*4	SECNDS	R*4		

```

0001      SUBROUTINE SAMPLE
          C
          C
          C              TABLE OF VARIABLES
          C
          C      IBAR      ACCUMULATOR FOR NITROGEN DRIFT IN INSPIRED AIR
          C      IBARRY    AVERAGE COUNT FOR NITROGEN IN INSPIRED AIR
          C      ICHECK    FLAG TO SET NEW VOLUME FOR BREATH OVER DETERMINATION
          C      IDLY      MASS SPEC TRANSPORT LAG COUNTER
          C      IFLLT     BIT CONFIGURATION FOR BREATH INDICATOR LAMPS
          C      K         DO LOOP DUMMY COUNTER
          C      KFLOW     ACCUMULATOR FOR VOLUME BREATH END DETERMINATION
          C      NRNG      INDEX FOR NITROGEN HI/LO SCALE USE
          C      TYPE
0002      IMPLICIT LOGICAL(L)
          C
          C      SIZE
0003      DIMENSION PC TCO(300),PC TN2(300),PC TU2(300)
0004      DIMENSION IADC GC(4),GASSLP(4),GSNTCP(4)
          C
          C      COMMON
0005      COMMON/FLOCOM/FLOZER,FLM,FLB,IFLO,BRTST
0006      COMMON/ADC ON/IADC SR,IADBUF,IADC GC,VOLUME,VOL1,VOL2
0007      COMMON/GAS/GASSLP,GSNTCP,LB TH,BRGAT,TEST
0008      COMMON/PARPRS/IC ON,PC TCO,PC TN2,PC TU2,CALZER
0009      COMMON/DRV11/IDRVIN,IDRVOT,IDLM,LWAT,IDRIFT
          C      SAMPLE N2 DRIFT BETWEEN BREATHS
0010      IF(LB TH) GO TO 45
0012      IBAR=0
0013      DO 44,J=1,5
0014      CALL IPOKE(IADC SR,IADC GC(1))
0015      44 IBAR=IBAR+IPEEK(IADBUF)
0016      IBARRY=IBAR/5
0017      IF(IDRIFT.GT.IBARRY)IDRIFT=IBARRY
          C
0019      45 IF(LWAT)GO TO 500
          C
          C
          C      CHECK N-2 RANGE SWITCH
0021      INT=IPEEK("167774)
0022      IF(INT.EQ."177777) GO TO 60
0024      IF(INT.EQ."177677) GO TO 50
0026      TYPE 55
0027      55 FORMAT(5X,'READ ERROR')
0028      GO TO 500
          C
          C      RANGE = 0-100%
0029      50 IFLLT="100
0030      NRNG=1
0031      GO TO 70
          C
          C      RANGE = 0-10%
0032      60 IFLLT="040
0033      NRNG=4

```

```

0034      GO TO 70
0035      70 CONTINUE

C
C          SET 'BREATH OVER' TEST EVERY 5th SAMPLE
0036      ICHECK=ICHECK+1
0037      IF( ICHECK.LT.10)GOTO 75
0039      VOL2=VOL1
0040      ICHECK=0
0041      75 KFLOW=0
0042      DO 80,K=1,15
0043      CALL IPOKE(IADC SR, IADC C(4))
0044      80 KFLOW=KFLOW+IPEEK(IADBUF)
0045      VOL1=FLOAT(KFLOW)/15.0
C          CHECK FOR BREATH
0046      TEST=VOL1-VOL2
0047      IF((VOL1.GE.BRTST).AND.(TEST.GT..5))GOTO 100

C
0049      CALL IPOKE(IDRVOT,"0")
0050      IF(LBTH) GO TO 110
0052      GO TO 500

C
0053      100 CONTINUE
0054      CALL IPOKE(IDRVOT,IFLLT)
0055      LBTH=.TRUE.
0056      IFLO=IFLO+1
0057      VOLUME=(VOL1+FLB)*FLM

C
0058      110 IDLY=IDLY+1
C
C          CHECK DELAY
0059      IF(IDLY.LT.IDLM) GO TO 500

C
C          GET GAS CONCENTRATIONS
0061      ICON=ICON+1
0062      INT=0
0063      DO 111 J=1,5
0064      CALL IPOKE(IADC SR, IADC C(1))
0065      111 INT=IPEEK(IADBUF)+INT
0066      INT=INT/5
0067      PC TN2(ICON)=(GASSLP(NRNG)*(INT+(CALZER-IDRIFT))
* +GSNTCP(NRNG))/100.

C
0068      INT=0
0069      DO 112 J=1,5
0070      CALL IPOKE(IADC SR, IADC C(2))
0071      112 INT=IPEEK(IADBUF)+INT
0072      INT=INT/5
0073      PC T02(ICON)=(GASSLP(2)*INT
* +GSNTCP(2))/100.0

C
0074      INT=0
0075      DO 113 J=1,5
0076      CALL IPOKE(IADC SR, IADC C(3))
0077      113 INT=IPEEK(IADBUF)+INT
0078      INT=INT/5

```

ORIGINAL PAGE IS
OF POOR QUALITY

0079 PC TCO(IC ON)=(GASSLP(3)*INT
* +GSNTCP(3))/100.

C
C
C

CHECK END OF DELAY

0080 IF(IC ON.LT.IFLO)GO TO 500
0082 LB TH=.FALSE.
0083 LWA T=.TRUE.
0084 IDLY=0

C
C

0085 500 CONTINUE
0086 RETURN
0087 END

Local Variables, .PSECT \$DATA, Size = 000024 (10. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
IBAR	I*2	000000	IBARRY	I*2	000004	ICHECK	I*2	000014
IDLY	I*2	000022	IFLLT	I*2	000010	INT	I*2	000006
J	I*2	000002	K	I*2	000020	KFLOW	I*2	000016
NRNG	I*2	000012						

COMMON Block /FLOCOM/, Size = 000022 (9. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
FLOZER	R*4	000000	FLM	R*4	000004	FLB	R*4	000010
IFLO	I*2	000014	BRTST	R*4	000016			

COMMON Block /ADCON /, Size = 000030 (12. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
IADC SR	I*2	000000	IADBUF	I*2	000002	IADC GC	I*2	000004
VOLUME	R*4	000014	VOL1	R*4	000020	VOL2	R*4	000024

COMMON Block /GAS /, Size = 000054 (22. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
GASSLP	R*4	000000	GSNTP	R*4	000020	LBTH	L*4	000040
BRGAT	R*4	000044	TEST	R*4	000050			

COMMON Block /PARPRS/, Size = 007026 (1803. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
ICON	I*2	000000	PC TCO	R*4	000002	PC TN2	R*4	002262
PC T02	R*4	004542	CALZER	R*4	007022			

COMMON Block /DRV11 /, Size = 000014 (6. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
IDRVIN	I*2	000000	IDRVOT	I*2	000002	IDLM	I*2	000004
LWAT	L*4	000006	IDRIFT	I*2	000012			

CRITICAL QUALITY
OF POOR QUALITY

Local and COMMON Arrays:

Name	Type	Section	Offset	-----Size-----	Dimensions
GASSLP	R*4	GAS	000000	000020 (8.)	(4)
GSNTCP	R*4	GAS	000020	000020 (8.)	(4)
IADC CC	I*2	ADC ON	000004	000010 (4.)	(4)
PC TC0	R*4	PARPRS	000002	002260 (600.)	(300)
PC TN2	R*4	PARPRS	002262	002260 (600.)	(300)
PC T02	R*4	PARPRS	004542	002260 (600.)	(300)

Subroutines, Functions, Statement and Processor-Defined Functions:

Name	Type	Name	Type	Name	Type	Name	Type	Name	Type
FLOAT	R*4	IPEEK	I*2	IPOKE	I*2				

5.4 NITROGEN DATA PROGRAM

0001

PROGRAM N2DATA

```
C *****
C *****
C *
C *
C *           TECHNOLOGY INCORPORATED
C *           LIFE SCIENCES DIVISION
C *
C *
C *
C *****
C *
C *
C *   PROGRAM NAME:.....N2DATA
C *   PROGRAMMERS,VERSION 1.0 A:.....BARRY LEVITAN AND
C *                                  KEN KOESTER
C *   PROGRAM MODIFIED,VERSION 1.0 B.....BARRY LEVITAN
C *   PROGRAM MODIFIED,VERSION 2.0.....PENNY BUEKER
C *   DATE:VERSION 1.0 A.....AUGUST 1980
C *   DATE:VERSION 1.0 B.....AUGUST 1981
C *   DATE:VERSION 2.0.....DECEMBER 1981
C *
C *
C *-----*
C *
C *   COMPUTER SYSTEM:.....LSI-11
C *   OPERATING SYSTEM:.....RT-11v4
C *
C *
C *-----*
C *
C *   COMPILING SEQUENCE:
C *
C *       LSI:  FORTRAN N2DATA
C *
C *
C *-----*
C *
C *   LINKING SEQUENCE:
C *
C *       LSI:  LINK N2DATA
C *
C *
C *-----*
C *
C *   EXECUTION SEQUENCE:  RUN N2DATA
C *****
```

N2DA TA PROGRAM WRITES A REPORT (FULL, 1 MINUTE, 5 MINUTE OR SHORT)
ON THE NITROGEN WASHOUT EXPERIMENT

VARIABLE	USE
AGE	SUBJECT'S AGE
AMGAS	CODE FOR AMBIENT GAS MIXTURE
AMPRES	AMBIENT PRESSURE IN MM HG
ANSWER	TEMPORARY FOR VALUE OF IMMERS
BMR	METABOLIC RATE
BODHYD	CODE FOR BODY HYDRATION
BODPOS	CODE FOR BODY POSITION
BRAVG	AVERAGE BREATH VOLUME
BRN2V	AVERAGE BREATH NITROGEN VOLUME
BRSM	INTERVAL SUM OF BREATH VOLUMES
BRSUM	SUM OF BREATH VOLUMES
BRTMIX	CODE FOR BREATHING MIXTURE
BRVOL	BREATH VOLUME
BTH	RUNNING TOTAL NUMBER OF BREATHS
COTDAV	AVERAGE CO2 CONCENTRATION
COTDSM	SUM OF CO2 CONCENTRATIONS
COTID	CO2 CONCENTRATION OF BREATH
COTSM	SUM OF CO2 CONCENTRATION OF INTERVAL
COJHT	BREATH COUNTER
COVOL	INTERVAL VOLUME OF CO2 IN LITERS
CO2	INTERVAL VOLUME OF CO2 IN ml
CUMTIM	ELAPSED TIME BEYOND INTERVAL
DRUGS	CODE FOR DRUGS
ELPST	ELAPSED TIME
EXERCI	CODE FOR EXERCISE PROTOCOL
EXGRAV	EXTERNAL GRAVITY (1 OR 0)
EXTEMP	EXTERNAL TEMPERATURE DEGREES CENTIGRADE
I	INDEX COUNTER
IAN5	OPERATOR SUPPLIED ANSWER
IBTH	BREATH NUMBER
IMMERS	IMMERSION (YES/NO)
IMIN	ELAPSED MINUTES
IYS	YES
J	INDEX COUNTER
KCM	TEMPORARY FOR VALUE OF KCMMTS
KCMMTS	COMMENTS (YES/NO)
M	INDEX COUNTER
MONTH	MONTH
N2	INTERVAL VOLUME OF N2 IN ML
N2AVG	AVERAGE N2 VOLUME
N2AVT	INTERVAL AVERAGE OF N2 CONCENTRATION
N2CONT	ESTIMATED N2 CONTENT OF BODY
N2CUM	CUMMULATIVE N2 ELIMINATED
N2SM	INTERVAL VOLUME OF N2 IN LITERS
N2SUM	SUM OF N2 VOLUME
N2TID	N2 CONCENTRATION OF BREATH
N2TSM	SUM OF N2 CONCENTRATIONS OF INTERVAL
N2VOL	INTERVAL VOLUME OF N2 IN LITERS

ORIGINAL QUALITY
OF POOR QUALITY

C	N2VWT	ELIMINATED N2 ml/kg BODY WEIGHT
C	NTDAVG	AVERAGE OF AVERAGE N2 CONCENTRATION
C		PER BREATH
C	NTDSM	SUM OF AVERAGE N2 CONCENTRATION
C		PER BREATH
C	O2	INTERVAL VOLUME OF O2 IN ml
C	O2TID	O2 CONCENTRATION OF BREATH
C	O2TSM	SUM OF O2 CONCENTRATIONS OF INTERVAL
C	O2VOL	INTERVAL VOLUME OF O2 IN LITERS
C	OTDAVG	AVERAGE O2 CONCENTRATION
C	OTDSM	SUM OF O2 CONCENTRATIONS
C	PC TN2	% OF ESTIMATED N2 CONTENT ELIMINATED
C	PULN2	PULMONARY N2
C	SPCOND	CODE FOR OTHER SPECIAL CONDITIONS
C	TIMCNT	ELAPSED TIME IN INTERVAL
C	WEIGHT	SUBJECT'S BODY WEIGHT IN kg
C		
C	ARRAY	USE
C	COMMTS	OPTIONAL SECTION FOR COMMENTS
C	IDA T	DATE
C	IDCODE	SUBJECT IDENTIFICATION NUMBER
C	IFIL	FILE NAME
C	INAM	NAME OF SUBJECT
C	STRTIM	STARTING TIME OF DAY
C		
C		
0002		IMPLICIT REAL(N)
0003		LOGICAL*1 IANS,IYS,COMMTS(4,80),IMMERS,KCMMTS, 1 STRTIM(5),IDA T(8),IDCODE(11), 2 INAM(30),IFIL(31)
0004		INTEGER*2 AMGAS,BODHYD,BODPOS,BRTMIX,DRUGS,EXERC I, 1 EXGRAV,SPCOND,AGE,BTH,IMIN
0005		REAL*4 ANSWER,KCM,MONTH,BRSM,BRVOL,COAVT,COTDAV,COTDSM, 1 COTSM,COVOL,CUMTIM,N2AVT,N2CONT,N2CUM,N2SM,N2TSM, 2 N2VOL,N2VWT,O2AVT,O2TSM,O2VOL,OTDAVG,OTDSM,PC TN2, 3 TIMCNT,WEIGHT,COSM,O2SM,BMR,PULN2
C		
0006		DATA IYS/1HY/
C		
0007		5 WRITE(7,910)
0008		READ(5,810) MC,(IFIL(M),M=1,MC)
0009		IFIL(MC+1)=0
C		
0010		OPEN(UNIT=2,TYPE='OLD',NAME=IFIL,FORM='FORMATTED')
C		
0011		READ(2,814) (INAM(I),I=1,30),(IDA T(I),I=1,8), 1 (IDCODE(J),J=1,11),AGE, 2 EXTEMP,EXGRAV,IMMERS,BODPOS,AMPRES,DRUGS, 3 EXERC I,BRTMIX,BODHYD,AMGAS,STRTIM,SPCOND,KCMMTS
0012		IF(IMMERS.EQ.'N') ANSWER='NO'
0014		IF(IMMERS.EQ.'Y') ANSWER='YES'
0016		IF(IMMERS.EQ.' ') ANSWER='NO'
0018		IF(KCMMTS.EQ.'N') KCM='NO'
0020		IF(KCMMTS.EQ.'Y') KCM='YES'
0022		IF(KCMMTS.EQ.' ') KCM='NO'

```

0024      DO 10 I=1,4
0025      10 READ(2,815) (COMMTS(I,J),J=1,80)
0026      11 WRITE(7,939)
0027      READ(5,877) IANS
0028      IF (.NOT.((IANS.EQ.'F').OR.(IANS.EQ.'1').OR.(IANS.EQ.'5')
1 .OR.(IANS.EQ.'S'))) GO TO 11
0030      FAUSE ' ADVANCE CARRIAGE IF DESIRED. RETURN'
C
C      N2DATA REPORT
C      WRITE HEADER INFORMATION
0031      WRITE(7,920)
0032      CALL SWAPDT(IDAT,MONTH)
0033      WRITE(7,921)(INAM(I),I=1,30),AGE, IDAT(4), IDAT(5),
1 IDAT(3),MONTH, IDAT(6), IDAT(7), IDAT(8)
0034      WRITE(7,922)(IDCODE(J),J=1,11),(STRTIM(J),J=1,5)
0035      WRITE(7,923)EXTEMP,EXGRAV
0036      WRITE(7,924)ANSWER,BODPOS
0037      WRITE(7,925)AMPRES,DRUGS
0038      WRITE(7,926)EXERC I,BRTMIX
0039      WRITE(7,927)BODHYD,AMGAS
0040      WRITE(7,928)SPCOND,KCM
0041      IF(KCMMS.EQ.'N') GO TO 100
0043      WRITE(7,929)
0044      DO 20 I=1,4
0045      20 WRITE(7,930)(COMMTS(I,J),J=1,80)
C
C      BREATH BY BREATH PARAMETERS FOR FULL REPORT
C
0046      100 IF(IANS.NE.'F')GO TO 200
0048      WRITE(7,938)
C      INITIALIZE VARIABLES
0049      BRSUM=0.0
0050      N2SUM=0.0
0051      NTDSM=0.0
0052      OTDSM=0.0
0053      COTDSM=0.0
C      READ A RECORD AND WRITE A RECORD
0054      DO 150,J=1,10000
0055      READ(2,816) IBTH,ELPST,BRVOL,BRN2V,02TID,N2TID,COTID
0056      IF(IBTH.EQ.0) GO TO 500
0058      BTH=IBTH
0059      WRITE(7,935)IBTH,ELPST,BRVOL,BRN2V,02TID,N2TID,COTID
C      SUM THE PARAMETERS
0060      BRSUM=BRSUM+BRVOL
0061      N2SUM=N2SUM+BRN2V
0062      NTDSM=NTDSM+N2TID
0063      OTDSM=OTDSM+02TID
0064      COTDSM=COTDSM+COTID
0065      150 CONTINUE
C
C      1 MINUTE INTERVAL PARAMETERS FOR 1 MINUTE INTERVAL REPORT
C
0066      200 IF (IANS.NE.'1') GO TO 300
0068      WRITE(7,941)
0069      READ(5,*)WEIGHT

```

ORIGINAL DOCUMENT
OF POOR QUALITY

```

0070      WRITE(7,942)
0071      READ(5,*)N2CONT
0072      WRITE(7,944)
0073      READ(5,*)PULN2
0074      WRITE(7,940)
      C    INITIALIZE VARIABLES
0075      IMIN=0
0076      COUNT=0
0077      N2CUM=0.0
0078      CUMTIM=0.0
0079      BRSM=0.0
0080      N2SM=0.0
0081      O2SM=0.0
0082      COSM=0.0
0083      N2TSM=0.0
0084      O2TSM=0.0
0085      COTSM=0.0
0086      BRSUM=0.0
0087      N2SUM=0.0
0088      NTDSM=0.0
0089      OTDSM=0.0
0090      COTDSM=0.0
      C    READ A 1 MINUTE INTERVAL
0091      DO 220,J=1,10000
0092  210  READ(2,816)IBTH,ELPST,BRVOL,BRN2V,O2TID,N2TID,COTID
0093      IF(IBTH.EQ.0) GO TO 215
0095      TMCNT=ELPST-CUMTIM
0096      IF(TMCNT.GT.60) GO TO 213
0098  211  BTH=IBTH
      C    SUM THE BREATH AND N2 VOLUMES IN LITERS
0099      BRSM=BRSM+BRVOL
0100      N2SM=N2SM+BRN2V
      C    COMPUTE THE O2 AND CO2 VOLUME IN LITERS
0101      O2VOL=O2TID*BRVOL
0102      COVOL=COTID*BRVOL
      C    SUM THE O2 AND CO2 VOLUME
0103      O2SM=O2SM+O2VOL
0104      COSM=COSM+COVOL
      C    SUM THE O2, CO2 AND N2 CONCENTRATIONS
0105      O2TSM=O2TSM+O2TID
0106      COTSM=COTSM+COTID
0107      N2TSM=N2TSM+N2TID
0108      COUNT=COUNT+1
0109      GO TO 210
0110  213  IMIN=IMIN+1
0111      CUMTIM=CUMTIM+60
      C    CHANGE THE VOLUMES TO ml
0112      N2=N2SM*1000
0113      O2=O2SM*1000
0114      CO2=COSM*1000
      C    COMPUTE THE 1 MINUTE AVERAGE N2 CONCENTRATION
0115      N2AVT=N2TSM/COUNT
      C    COMPUTE THE 1 MINUTE N2 ml/kg BODY WEIGHT
0116      N2VWT=N2/WEIGHT
      C    COMPUTE THE CUMMULATIVE N2 ELIMINATED

```

COMPUTATION OF PULN2

```

0117      N2CUM=N2CUM+N2
C        COMPUTE THE % OF ESTIMATED BODY N2 CONTENT ELIMINATED
0118      PC TN2=((N2CUM/1000)-PULN2)/N2CONT)*100
C        COMPUTE THE METABOLIC RATE
0119      BMR=COSM*1362.4
C        TRANSFER THE SUMS TO SUMMARY VARIABLES
0120      215 BRSUM=BRSUM+BRSM
0121      N2SUM=N2SUM+N2SM
0122      OTDSM=OTDSM+O2TSM
0123      NTDSM=NTDSM+N2TSM
0124      COTDSM=COTDSM+COTSM
0125      IF(IBTH.EQ.0)GO TO 500
C        WRITE THE 1 MINUTE INTERVAL AVERAGES
0127      WRITE(7,943)IMIN,N2,O2,CO2,N2VWT,N2AVT,N2CUM,PC TN2,BMR
C        INITIALIZE VARIABLES
0128      COUNT=0
0129      BRSM=0.0
0130      N2SM=0.0
0131      N2TSM=0.0
0132      O2TSM=0.0
0133      COTSM=0.0
0134      COSM=0.0
0135      O2SM=0.0
0136      GO TO 211
0137      220 CONTINUE
C
C        5 MINUTE INTERVAL PARAMETERS FOR 5 MINUTE INTERVAL REPORT
C
0138      300 IF(IANS.NE.'5') GO TO 400
0140      WRITE(7,941)
0141      READ(5,*)WEIGHT
0142      WRITE(7,942)
0143      READ(5,*)N2CONT
0144      WRITE(7,944)
0145      READ(5,*)PULN2
0146      WRITE(7,940)
C        INITIALIZE VARIABLES
0147      IMIN=0
0148      COUNT=0
0149      N2CUM=0.0
0150      CUMTIM=0.0
0151      BRSM=0.0
0152      N2SM=0.0
0153      COSM=0.0
0154      O2SM=0.0
0155      N2TSM=0.0
0156      O2TSM=0.0
0157      COTSM=0.0
0158      BRSUM=0.0
0159      N2SUM=0.0
0160      NTDSM=0.0
0161      OTDSM=0.0
0162      COTDSM=0.0
C        READ A 5 MINUTE INTERVAL AND SUM THE PARAMETERS
0163      DO 320,J=1,10000

```

ORIGINAL D. 113
OF POOR QUALITY

```

0164 310 READ(2,816)IB TH,ELPST,BRVOL,BRN2V,O2TID,N2TID,COTID
0165     IF(IB TH.EQ.0) GO TO 315
0167     TMCNT=ELPST-CUMTIM
0168     IF(TMCNT.GT.300)GO TO 313
0170 311 B TH=IB TH
      C   SUM THE BREATH AND N2 VOLUMES IN LITERS
0171     BRSM=BRSM+BRVOL
0172     N2SM=N2SM+BRN2V
      C   COMPUTE THE O2 AND CO2 VOLUME IN LITERS
0173     O2VOL=O2TID*BRVOL
0174     COVOL=COTID*BRVOL
      C   SUM THE O2 AND CO2 VOLUME
0175     O2SM=O2SM+O2VOL
0176     COSM=COSM+COVOL
      C   SUM THE O2, CO2 AND N2 CONCENTRATIONS
0177     O2TSM=O2TSM+O2TID
0178     CO TSM=CO TSM+COTID
0179     N2 TSM=N2 TSM+N2 TID
0180     COUNT=COUNT+1
0181     GO TO 310
0182 313 IMIN=IMIN+5
0183     CUMTIM=CUMTIM+300
      C   CHANGE THE VOLUMES TO ml
0184     N2=N2SM*1000
0185     O2=O2SM*1000
0186     CO2=COSM*1000
      C   COMPUTE THE 5 MINUTE AVERAGE N2 CONCENTRATION
0187     N2AVT=N2TSM/COUNT
      C   COMPUTE THE 5 MINUTE N2 ml/kg BODY WEIGHT
0188     N2VWT=N2/WEIGHT
      C   COMPUTE THE CUMMULATIVE N2 ELIMINATED
0189     N2CUM=N2CUM+N2
      C   COMPUTE THE % OF ESTIMATED BODY N2 CONTENT ELIMINATED
0190     PC TN2=((N2CUM/1000)-PULN2)/N2CONT)*100
      C   COMPUTE THE METABOLIC RATE
0191     BMR=(COSM/5)*1362.4
      C   TRANSFER THE SUMS TO SUMMARY VARIABLES
0192 315 BRSUM=BRSUM+BRSM
0193     N2SUM=N2SUM+N2SM
0194     O2SUM=O2SUM+O2TSM
0195     N2TSM=N2TSM+N2TSM
0196     CO TSM=CO TSM+COTSM
0197     IF(IB TH.EQ.0)GO TO 500
      C   WRITE THE 5 MINUTE INTERVAL AVERAGES
0199     WRITE(7,943)IMIN,N2,O2,CO2,N2VWT,N2AVT,N2CUM,PC TN2,BMR
      C   INITIALIZE VARIABLES
0200     COUNT=0
0201     BRSM=0.0
0202     N2SM=0.0
0203     N2TSM=0.0
0204     O2TSM=0.0
0205     CO TSM=0.0
0206     COSM=0.0
0207     O2SM=0.0
0208     GO TO 311

```

ORIGINAL DOCUMENT
OF POOR QUALITY

```

0209      320 CONTINUE
          C
          C      SHORT REPORT
          C
0210      400 IF(IANS.NE.'S') GO TO 11
          C      INITIALIZE VARIABLES
0212      BRSUM=0.0
0213      N2SUM=0.0
0214      NTDSM=0.0
0215      OTDSM=0.0
0216      COTDSM=0.0
          C      READ RECORDS AND SUM PARAMETERS
0217      DO 440,J=1,10000
0218      READ(2,816)IBTH,ELPST,BRVOL,BRN2V,O2TID,N2TID,COTID
0219      IF(IBTH.EQ.0) GO TO 500
0221      BTH=IBTH
0222      BRSUM=BRSUM+BRVOL
0223      N2SUM=N2SUM+BRN2V
0224      OTDSM=OTDSM+O2TID
0225      NTDSM=NTDSM+N2TID
0226      COTDSM=COTDSM+COTID
0227      440 CONTINUE
          C      CHECK FOR END OF DATA
0228      500 IF(ELPST.EQ.0.0) GO TO 600
0230      WRITE(7,955)
0231      GO TO 670
          C      COMPUTE AVEARAGES FOR SUMMARY
0232      600 BRAVG=BRSUM/BTH
0233      N2AVG=N2SUM/BTH
0234      OTDAVG=OTDSM/BTH
0235      NTDAVG=NTDSM/BTH
0236      COTDAV=COTDSM/BTH
0237      WRITE(7,961)
0238      WRITE(7,962)
0239      WRITE(7,963)BRAVG,N2AVG,OTDAVG,NTDAVG,COTDAV
0240      WRITE(7,964)N2SUM
0241      WRITE(7,965)BTH
          C
0242      670 CLOSE(UNIT=2,DISP='SAVE')
          C
0243      IANS=' '
0244      700 WRITE(7,972)
          C      CHECK TO SEE IF THEY WANT TO READ ANOTHER FILE
0245      ACCEPT 877,IANS
0246      IF(IANS.EQ.IYS) GO TO 5
0248      910 FORMAT(//10X,'NITROGEN WASHOUT DATA ANALYSIS',
1 /,5X,'ENTER FILENAME ON WHICH DATA IS STORED ',
2 /5X,'(DEVICE:FILENAME.TYPE): ',S)
0249      920 FORMAT(/25X,'N2DATA REPORT'/)
0250      921 FORMAT(/2X,30A1,T41,'AGE: ',I2,T52,'DATE: ',2X,3A1,A3,3A1)
0251      922 FORMAT(/2X,'SUBJECT ID NUMBER: ',T24,11A1,T41,'STARTING ',
1 'TIME: ',4X,5A1)
0252      923 FORMAT(/2X,'EXTERNAL TEMPERATURE: ',T28,F7.2,T41,'EXTERNAL ',
1 'GRAVITY: ',T68,I1)
0253      924 FORMAT(/2X,'IMMERSION? ',T32,A3,T41,'BODY POSITION CODE: ',

```

ORIGINAL PRINTING
OF POOR QUALITY

0254 1 T67,I2)
 925 FORMAT(/2X,'AMBIENT PRESSURE: ',T28,F7.2,T41,'DRUG CODE: ',
 1 T67,I2)
 0255 926 FORMAT(/2X,'EXERCISE CODE: ',T33,I2,T41,'BREATH MIXTURE CODE: ',
 1 T67,I2)
 0256 927 FORMAT(/2X,'BODY HYDRATION CODE: ',T33,I2,T41,'AMBIENT GAS CODE: ',
 1 T67,I2)
 0257 928 FORMAT(/2X,'SPECIAL CONDITION CODE: ',T33,I2,T41,'COMMENTS? ',
 1 T66,A3)
 0258 929 FORMAT(/2X,'COMMENTS: ')
 0259 930 FORMAT(2X,80A1)
 0260 935 FORMAT(5X,I4,F10.1,5(F10.5))
 0261 938 FORMAT(/2X,'BREATH # ',2X,'ELAPSED ',2X,'BREATH ',
 1 2X,'N2 VOL ',2X,'O2 CONC ',2X,'N2 CONC ',
 2 2X,'CO2 CONC ',
 3 /13X,' TIME ',2X,' VOLUME',2X,'/BREATH ',2X,
 4 '/BREATH ',2X,'/BREATH ',2X,'/BREATH ',
 5 /13X,' (SEC) ',2X,' (L) ',2X,' (L) ')
 0262 939 FORMAT(/5X,'A FULL REPORT (F) GIVES BREATH BY BREATH PARAMETERS',
 1 /5X,'A 1 MINUTE REPORT (1) GIVES PARAMETERS FOR 1 MINUTE ',
 2 'INTERVALS',
 3 /5X,'A 5 MINUTE REPORT (5) GIVES PARAMETERS FOR 5 MINUTE ',
 4 'INTERVALS',
 5 /5X,'A SHORT REPORT (S) GIVES THE OVERALL AVERAGES',
 6 //5X,'DO YOU WANT A FULL REPORT (F), A 1 MINUTE REPORT (1),',
 7 /5X,'A 5 MINUTE REPORT (5),OR A SHORT REPORT (S) ?',
 8 /5X,'TYPE (F/1/5/S): ',S)
 0263 940 FORMAT(/3X,' TIME',2X,'N2 VOL/',2X,'O2 VOL/',1X,'CO2 VOL/',
 1 1X,'N2/WEIGHT',1X,'N2 CONC/',1X,'CUMMUL.',2X,'% OF EST.',2X,
 2 'BMR',/2X,'(min)',1X,'INTERVAL',1X,'INTERVAL',1X,'INTERVAL',
 3 2X,'(ml/kg)',2X,'INTERVAL',3X,'N2',6X,'N2 CONT.','(BTU/hr)',
 4 /10X,'(ml)',5X,'(ml)',5X,'(ml)',23X,'(ml)',6X,'(%)')
 0264 941 FORMAT(/2X,39HENTER THE SUBJECT'S BODY WEIGHT IN kg: ,S)
 0265 942 FORMAT(/2X,52HENTER THE SUBJECT'S ESTIMATED N2 CONTENT IN LITERS:
 1 ,S)
 0266 944 FORMAT(/2X,42HENTER THE SUBJECT'S ESTIMATED PULMONARY N2,
 1 20H CONTENT IN LITERS: ,S)
 0267 943 FORMAT(2X,I3,1X,F9.2,1X,F9.2,F9.2,F9.5,F9.5,F9.2,F9.2,F9.2)
 0268 955 FORMAT(10X,'BAD END DATA')
 0269 961 FORMAT(/25X,'AVERAGES: ')
 0270 962 FORMAT(/4X,' BREATH ',4X,'N2 VOLUME ',4X,
 1 'O2 CONC ',4X,'N2 CONC ',4X,'CO2 CONC ',
 2 /4X,' VOLUME(L)',4X,'/BREATH(L)',4X,'/BREATH',
 3 4X,'/BREATH',4X,'/BREATH')
 0271 963 FORMAT(/5X,F7.5,6X,F7.5,5X,F8.5,3X,F8.5,3X,F8.5)
 0272 964 FORMAT(/2X,' TOTAL N2 VOLUME IN LITERS= ',F8.5)
 0273 965 FORMAT(/2X,' TOTAL # OF BREATHS= ',I6)
 0274 972 FORMAT(/5X,'END OF DATA..READ ANOTHER FILE (Y/N): ',S)
 0275 810 FORMAT(Q,30A1)
 0276 814 FORMAT(30A1,8A1,11A1,I2,F7.2,I1,A1,I2,F7.2,5(I2),5A1,I2,A1)
 0277 815 FORMAT(80A1)
 0278 816 FORMAT(I4,F10.1,5(F10.5))
 0279 877 FORMAT(1A1)
 0280 END

C

UNIVERSITY
OF POOR

Local Variables, .PSECT \$DATA, Size = 001220 (328. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
AGE	I*2	000670	AMGAS	I*2	000650	AMPRES	R*4	001102
ANSWER	R*4	000676	BMR	R*4	001056	BODHYD	I*2	000652
BODPOS	I*2	000654	BRAVG	R*4	001174	BRN2V	R*4	001130
BRSM	R*4	000712	BRSUM	R*4	001106	BRMIX	I*2	000656
BRVOL	R*4	000716	BTH	I*2	000672	COAVT	R*4	000722
COSM	R*4	001046	COTDAV	R*4	000726	COTDSM	R*4	000732
COTID	R*4	001144	COTSM	R*4	000736	COUNT	R*4	001150
COVOL	R*4	000742	CO2	R*4	001170	CUMTIM	R*4	000746
DRUGS	I*2	000660	ELPST	R*4	001124	EXERC I	I*2	000662
EXGRAV	I*2	000664	EXTEMP	R*4	001076	I	I*2	001072
IANS	L*1	000644	IBTH	I*2	001122	IMIN	I*2	000674
IMMERS	L*1	000645	IYS	L*1	000625	J	I*2	001074
KCM	R*4	000702	KCMMTS	L*1	000646	M	I*2	001070
MC	I*2	001066	MONTH	R*4	000706	NTDAVG	R*4	001204
NTDSM	R*4	001116	N2	R*4	001160	N2AVG	R*4	001200
N2AVT	R*4	000752	N2CONT	R*4	000756	N2CUM	R*4	000762
N2SM	R*4	000766	N2SUM	R*4	001112	N2TID	R*4	001140
N2TSM	R*4	000772	N2VOL	R*4	000776	N2VWT	R*4	001002
OTDAVG	R*4	001022	OTDSM	R*4	001026	O2	R*4	001164
O2AVT	R*4	001006	O2SM	R*4	001052	O2TID	R*4	001134
O2TSM	R*4	001012	O2VOL	R*4	001016	PC TN2	R*4	001032
PULN2	R*4	001062	SPCOND	I*2	000666	TIMCNT	R*4	001036
TMCNT	R*4	001154	WEIGHT	R*4	001042			

ORIGINAL FILED IN
OF POOR QUALITY

Local and COMMON Arrays:

Name	Type	Section	Offset	-----Size-----	Dimensions
COMMTS	L*1	\$DA TA	000000	000500 (160.)	(4,80)
IDAT	L*1	\$DA TA	000505	000010 (4.)	(8)
IDCODE	L*1	\$DA TA	000515	000013 (6.)	(11)
IFIL	L*1	\$DA TA	000566	000037 (16.)	(31)
INAM	L*1	\$DA TA	000530	000036 (15.)	(30)
STRTIM	L*1	\$DA TA	000500	000005 (3.)	(5)

Subroutines, Functions, Statement and Processor-Defined Functions:

Name	Type	Name	Type	Name	Type	Name	Type	Name	Type
SWAPDT	R*4								

```

0001 SUBROUTINE SWAPD T(X,MONTH)
0002 LOGICAL*1 X(8)
0003 REAL*4 MONTH
0004 IF(X(1).EQ.'0'.AND.X(2).EQ.'1') MONTH='JAN'
0006 IF(X(1).EQ.'0'.AND.X(2).EQ.'2') MONTH='FEB'
0008 IF(X(1).EQ.'0'.AND.X(2).EQ.'3') MONTH='MAR'
0010 IF(X(1).EQ.'0'.AND.X(2).EQ.'4') MONTH='APR'
0012 IF(X(1).EQ.'0'.AND.X(2).EQ.'5') MONTH='MAY'
0014 IF(X(1).EQ.'0'.AND.X(2).EQ.'6') MONTH='JUN'
0016 IF(X(1).EQ.'0'.AND.X(2).EQ.'7') MONTH='JUL'
0018 IF(X(1).EQ.'0'.AND.X(2).EQ.'8') MONTH='AUG'
0020 IF(X(1).EQ.'0'.AND.X(2).EQ.'9') MONTH='SEP'
0022 IF(X(1).EQ.'1'.AND.X(2).EQ.'1') MONTH='JAN'
0024 IF(X(1).EQ.'1'.AND.X(2).EQ.'2') MONTH='FEB'
0026 IF(X(1).EQ.'1'.AND.X(2).EQ.'3') MONTH='MAR'
0028 IF(X(1).EQ.'1'.AND.X(2).EQ.'4') MONTH='APR'
0030 IF(X(1).EQ.'1'.AND.X(2).EQ.'5') MONTH='MAY'
0032 IF(X(1).EQ.'1'.AND.X(2).EQ.'6') MONTH='JUN'
0034 IF(X(1).EQ.'1'.AND.X(2).EQ.'7') MONTH='JUL'
0036 IF(X(1).EQ.'1'.AND.X(2).EQ.'8') MONTH='AUG'
0038 IF(X(1).EQ.'1'.AND.X(2).EQ.'9') MONTH='SEP'
0040 IF(X(1).EQ.'1'.AND.X(2).EQ.'0') MONTH='OCT'
0042 IF(X(1).EQ.'1'.AND.X(2).EQ.'1') MONTH='NOV'
0044 IF(X(1).EQ.'1'.AND.X(2).EQ.'2') MONTH='DEC'
0046 RETURN
0047 END

```

Local Variables, .PSECT \$DATA, Size = 000004 (2. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
------	------	--------	------	------	--------	------	------	--------

MONTH	R*4	@ 000002						
-------	-----	----------	--	--	--	--	--	--

Local and COMMON Arrays:

Name	Type	Section	Offset	-----Size-----	Dimensions
X	L*1	@ \$DATA	000000	000010 (4.)	(8)

ORIGINAL PART 1
OF FOUR PARTS

PROGRAM UPDATE TAKES PREVIOUSLY COLLECTED, UNFORMATTED
N2 DATA AND WRITES A FORMATTED FILE FOR A REPORT BY
N2DATA PROGRAM

VARIABLES	USE
AGE	SUBJECT'S AGE
AMGAS	CODE FOR AMBIENT GAS COMPOSITION
AMPRES	AMBIENT PRESSURE IN MM HG
BODHYD	CODE FOR BODY HYDRATION
BODPOS	CODE FOR BODY POSITION
BRN2V	N2 VOLUME OF BREATH
BRTMIX	CODE FOR BREATHING MIXTURE
BRVOL	VOLUME OF BREATH
COTID	CO2 CONCENTRATION OF BREATH
DRUGS	CODE FOR DRUGS
ELPST	ELAPSED TIME
EXERCI	CODE FOR EXERCISE PROTOCOL
EXGRAV	EXTERNAL GRAVITY (1, 0)
EXTEMP	EXTERNAL TEMPERATURE DEGREES CENTIGRADE
I	INDEX COUNTER
IBTH	BREATH NUMBER
IDMY	DUMMY INTEGER ZERO END OF DATA TAG
IMMERS	IMMERSION (YES/NO)
J	INDEX COUNTER
KCMMS	COMMENTS (YES/NO)
N2TID	N2 CONCENTRATION OF BREATH
O2TID	O2 CONCENTRATION OF BREATH
RDY	DUMMY REAL ZERO END OF DATA TAG
SPCOND	CODE FOR OTHER SPECIAL CONDITIONS

ARRAY	USE
COMMS	SECTION FOR OPTIONAL COMMENTS
IDA T	OLD INTEGER ARRAY FOR DATE
IDCODE	SUBJECT IDENTIFICATION NUMBER
IFIL	FILE NAME
INAM	NAME OF SUBJECT
LDA T	NEW LOGICAL ARRAY FOR DATE
STRTIM	STARTING TIME OF DAY

```

0002 LOGICAL*1 IDC CODE(11),IMMERS,STRTIM(5),KCMMS,
1 COMMS(4,80),IFIL(20),IFIL2(20),LNAM(30)
0003 INTEGER*2 INAM(15),IDA T(3),EXGRAV,BODPOS,
1 DRUGS,EXERCI,BRTMIX,BODHYD,AMGAS,SPCOND,IBTH,AGE
0004 REAL*4 EXTEMP,AMPRES,ELPST,BRVOL,BRN2V,O2TID,
1 N2TID,COTID

```

```

0005 EQUIVALENCE (INAM,LNAM)

```

```

INITIALIZE NEW VARIABLES

```

```

0006 DO 20 J=1,11
0007 20 IDC CODE(J)=' '

```

ORIGINAL PAGE IS
OF POOR QUALITY

```

0008     IMMERS=' '
0009     DO 30 J=1,5
0010     30 STRTIM(J)=' '
0011     KCMMTS=' '
0012     DO 40 J=1,30
0013     40 LNAM(J)=' '
0014     DO 80 I=1,4
0015     DO 80 J=1,80
0016     80 COMMTS(I,J)=' '
0017     AGE=0
0018     EXGRAV=0
0019     BODPOS=0
0020     DRUGS=0
0021     EXERC I=0
0022     BR TMIX=0
0023     BODHYD=0
0024     AMGAS=0
0025     SPCOND=0
0026     EXTEMP=0.0
0027     AMPRES=0.0

C
C
0028     WRITE(7,917)
0029     READ(5,810) MC, (IFIL(J),J=1,MC)
0030     IFIL(MC+1)=0
0031     WRITE(7,918)
0032     READ(5,810) MC, (IFIL2(J),J=1,MC)
0033     IFIL2(MC+1)=0
0034     OPEN(UNIT=1,TYPE='OLD',NAME=IFIL,FORM='UNFORMATTED')
0035     OPEN(UNIT=2,TYPE='NEW',NAME=IFIL2,FORM='FORMATTED')

C
C
0036     READ(1)INAM,IDA T
0037     WRITE(2,914)(LNAM(I),I=1,30),IDA T(1),'-',
1     IDA T(2),'-',IDA T(3),(IDCODE(I),I=1,11),
2     AGE,EXTEMP,EXGRAV,IMMERS,BODPOS,AMPRES,DRUGS,EXERC I,
3     BR TMIX,BODHYD,AMGAS,(STRTIM(I),I=1,5),SPCOND,KCMMTS
0038     DO 90 I=1,4
0039     90 WRITE(2,915)(COMMTS(I,J),J=1,80)
0040     100 READ(1,END=150,ERR=150)IB TH,ELPST,BRVOL,BRN2V,O2TID,
1     N2TID,COTID
0041     IF(IB TH.EQ.0)GO TO 200
0043     WRITE(2,916)IB TH,ELPST,BRVOL,BRN2V,O2TID,N2TID,COTID
0044     GO TO 100
0045     150 IB TH=0
0046     ELPST=0.
0047     BRVOL=0.
0048     BRN2V=0.
0049     O2TID=0.
0050     N2TID=0.
0051     COTID=0.
0052     200 WRITE(2,916)IB TH,ELPST,BRVOL,BRN2V,O2TID,N2TID,COTID

C
C
0053     300 CLOSE(UNIT=1,DISP='SAVE')

```

```

0054      CLOSE(UNIT=2,DISP='SAVE')
      C
      C
0055      914 FORMAT(30A1,2(I2,A1),I2,11A1,I2,F7.2,I1,A1,I2,F7.2,5(I2),
          1 5A1,I2,A1)
0056      915 FORMAT(80A1)
0057      916 FORMAT(I4,F10.1,5(F10.5))
0058      917 FORMAT(/5X,'ENTER FILENAME ON WHICH DATA IS STORED ',
          1 /10X,'(DEVICE:FILENAME.TYPE): ','$)
0059      918 FORMAT(/5X,'ENTER FILENAME ON WHICH DATA IS TO BE PLACED ',
          1 /10X,'(DEVICE:FILENAME.TYPE): ','$)
0060      919 FORMAT(I4,F10.1,4(F10.5))
0061      810 FORMAT(Q,20A1)
0062      END

```

COPIES OF THIS FILED
 OF POOR QUALITY

FORTRAN IV Storage Map for Program Unit UPDA TE
 Local Variables, .PSECT \$DATA, Size = 000754 (246. words)

Name	Type	Offset	Name	Type	Offset	Name	Type	Offset
AGE	I*2	000700	AMGAS	I*2	000672	AMPRES	R*4	000706
BODHYD	I*2	000670	BODPOS	I*2	000660	BRN2V	R*4	000722
BRTMIX	I*2	000666	BRVOL	R*4	000716	COTID	R*4	000736
DRUGS	I*2	000662	ELPST	R*4	000712	EXERCI	I*2	000664
EXGRAV	I*2	000656	EXTMP	R*4	000702	I	I*2	000744
IBTH	I*2	000676	IMMERS	L*1	000654	J	I*2	000742
KCMMTS	L*1	000655	MC	I*2	000746	N2TID	R*4	000732
O2TID	R*4	000726	SPCOND	I*2	000674			

Local and COMMON Arrays:

Name	Type	Section	Offset	-----Size-----	Dimensions
COMMTS	L*1	\$DATA	000020	000500 (160.)	(4,80)
IDAT	I*2	\$DATA	000626	000006 (3.)	(3)
IDCODE	L*1	\$DATA	000000	000013 (6.)	(11)
IFIL	L*1	\$DATA	000520	000024 (10.)	(20)
IFIL2	L*1	\$DATA	000544	000024 (10.)	(20)
INAM	I*2	\$DATA	000570	000036 (15.)	(15)
LNAM	L*1	\$DATA	000570	000036 (15.)	(30)
STRTIM	L*1	\$DATA	000013	000005 (3.)	(5)

ORIGINAL FILE
 OF POOR QUALITY

6.0 SAMPLE RUNS WITH INPUT AND OUTPUT

6.1 PERSONAL DATA INPUT PROGRAM

ORIGINAL PAGE IS
OF POOR QUALITY.

.RUN DY1:PDIHPT

NITROGEN WASHOUT PERSONAL DATA INPUT PROGRAM

FILENAME ON WHICH TO STORE DATA
DEVICE FILENAME.TYPE = SY:TEST.PDI

IDENTIFICATION NUMBER OF THE SUBJECT (XXX-XX-XXXX): [REDACTED]

DATE (MM-DD-YY) : 07-08-82

DATE OF BIRTH (MM-DD-YY): [REDACTED]

AGE IN YEARS: 36

SEX (M/F): M

HEIGHT IN CM: 178

NUDE WEIGHT IN KG: 78

IS % BODY FAT EXPLICITLY KNOWN (Y/N)? YU

ENTER % BODY FAT: 15

IS BODY COMPOSITION TO BE MEASURED BY WATER IMMERSION (Y/N)? Y

WATER TEMPERATURE IN CENTIGRADE: 34

WATER DENSITY IN KG/L: .988

RESIDUAL VOLUME IN L: 1.3

WEIGHT OF UNDERWATER CHAIR IN KG: 3.1

ENTER THE NUMBER OF TRIALS FOR UNDERWATER WEIGHT: 10

TRIAL 1 UNDERWATER WEIGHT IN KG: 4

TRIAL 2 UNDERWATER WEIGHT IN KG: 4

TRIAL 3 UNDERWATER WEIGHT IN KG: 5

TRIAL 4 UNDERWATER WEIGHT IN KG: 5

TRIAL 5 UNDERWATER WEIGHT IN KG: 4

TRIAL 6 UNDERWATER WEIGHT IN KG: 5

TRIAL 7 UNDERWATER WEIGHT IN KG: 5

TRIAL 8 UNDERWATER WEIGHT IN KG: 6

TRIAL 9 UNDERWATER WEIGHT IN KG: 5

TRIAL 10 UNDERWATER WEIGHT IN KG: 4

DATE OF MEASUREMENTS (MM-DD-YY): 12-22-81

IS BODY COMPOSITION TO BE MEASURED BY SKINFOLD PROCESS I (Y/N)? Y

BODY SIDE USED (R/L)? R

ENTER THE NUMBER OF MEASUREMENTS RECORDED: 3

MEASUREMENT 1 BICEPS IN MM: 15

MEASUREMENT 2 BICEPS IN MM: 15

MEASUREMENT 3 BICEPS IN MM: 14

MEASUREMENT 1 TRICEPS IN MM: 9

MEASUREMENT 2 TRICEPS IN MM: 9

MEASUREMENT 3 TRICEPS IN MM: 10

MEASUREMENT 1 SUBSCAPULAR IN MM: 5

MEASUREMENT 2 SUBSCAPULAR IN MM: 5

MEASUREMENT 3 SUBSCAPULAR IN MM: 4

MEASUREMENT 1 SUPRAILIAC IN MM: 18

MEASUREMENT 2 SUPRAILIAC IN MM: 19

MEASUREMENT 3 SUPRAILIAC IN MM: 18

DATE OF MEASUREMENTS (MM-DD-YY): 12-18-82

ORIGINAL PAGE IS
OF POOR QUALITY

IS BODY COMPOSITION TO BE MEASURED BY SKINFOLD PROCESS II (Y/N)? Y

BODY SIDE USED (R/L)? R

ENTER THE NUMBER OF MEASUREMENTS RECORDED: 3

MEASUREMENT 1 MIDCLAVICULAR IN MM: 5

MEASUREMENT 2 MIDCLAVICULAR IN MM: 5

MEASUREMENT 3 MIDCLAVICULAR IN MM: 4

MEASUREMENT 1 UMBILICUS IN MM: 11

MEASUREMENT 2 UMBILICUS IN MM: 10

MEASUREMENT 3 UMBILICUS IN MM: 10

MEASUREMENT 1 ANTERIOR MID-THIGH IN MM: 11

MEASUREMENT 2 ANTERIOR MID-THIGH IN MM: 11

MEASUREMENT 3 ANTERIOR MID-THIGH IN MM: 12

DATE OF MEASUREMENTS (MM-DD-YY): 12-25-82

ORIGINAL PAGE IS
OF POOR QUALITY

ALIGN PAPER TO TOP OF PAGE AND HIT CARRIAGE RETURN

PERSONAL CHARACTERISTICS PROFILE

SUBJECT	██████████	DATA FILE	SY:TEST.PDI
DATE	08-JUL-82	MALE	AGE 36
HEIGHT	178.0 CM	% BODY FAT	15.0 %
NUDE WEIGHT	78.0 KG	DATE OF BIRTH	██████████

BODY COMPOSITION EVALUATIONS

----- WATER IMMERSION PROCESS -----

NUDE WEIGHT	78.00 KG	WATER TEMPERATURE	34.00 C
RESIDUAL VOLUME	1.30 L	WATER DENSITY	0.99 KG/L
UNDERWATER CHAIR WEIGHT	3.10 KG	DATE	22-DEC-81

UNDERWATER WEIGHT (KG) AVERAGE UNDERWATER WEIGHT = 4.70 KG

TRIAL	WEIGHT (KG)	TRIAL	WEIGHT (KG)
1	4.00	2	4.00
3	5.00	4	5.00
5	4.00	6	5.00
7	5.00	8	6.00
9	5.00	10	4.00

WATER WEIGHT	1.60 KG	BODY VOLUME	76.03 L
BODY DENSITY	1.026 KG/L	% BODY FAT	32.48 %

ORIGINAL PAGE IS
OF POOR QUALITY

----- SKINFOLD PROCESS I -----

BODY SIDE USED RIGHT

DATE 18-DEC-82

TRIAL MEASUREMENTS IN MM

	BICEPS	TRICEPS	SUBSCAPULAR	SUPRAILIAC
1	15.00	9.00	5.00	18.00
2	15.00	9.00	5.00	19.00
3	14.00	10.00	4.00	18.00
MEAN	14.67	9.33	4.67	18.33

BODY DENSITY 1.055 KG/L

% BODY FAT 19.05 %

----- SKINFOLD PROCESS II -----

BODY SIDE USED RIGHT

DATE 25-DEC-82

TRIAL MEASUREMENTS IN MM

	MIDCLAVICULAR	UMBILICUS	ANTERIOR MID-THIGH
1	5.00	11.00	11.00
2	5.00	10.00	11.00
3	4.00	10.00	12.00
MEAN	4.67	10.33	11.33

BODY DENSITY 1.079 KG/L

% BODY FAT 8.57 %

----- NUTRITION JOURNAL FORMULA -----

BODY VOLUME 74.23 L

BODY DENSITY 1.051 KG/L

% BODY FAT 21.20 %

% FAT FREE MASS 78.80 %

% TOTAL BODY WATER 58.15 %

----- NITROGEN CONTENT DATA -----

	KNOWN % BODY FAT	WATER IMMERSION	SKINFOLD PROCESS I	SKINFOLD PROCESS II	NUTRITION JOURNAL
% BODY FAT	15.00	32.48	19.05	8.57	21.20
FAT WEIGHT (KG)	11.70	25.34	14.86	6.68	16.54
FAT FREE MASS (KG)	66.30	52.66	63.14	71.32	61.46
LEAN BODY WATER MASS (KG)	47.74	37.92	45.46	51.35	44.25
VOLUME OF FAT (L)	13.00	28.15	16.51	7.42	18.38
VOLUME LEAN BODY WATER (L)	48.06	38.17	45.77	51.70	44.55
NITROGEN IN FAT (L)	0.871	1.886	1.106	0.497	1.231
NITROGEN IN LEAN BODY WATER (L)	0.625	0.496	0.595	0.672	0.579
TOTAL BODY NITROGEN (L)	1.496	2.383	1.701	1.169	1.810
NITROGEN/KG BODY WEIGHT (L)	0.0192	0.0305	0.0218	0.0150	0.0232

ORIGINAL PAGE IS
OF POOR QUALITY

6.2 DECOMPRESSION DATA PROGRAM

.RUN BY:DECDAT

DECOMPRESSION EXPERIENCE DATA PROGRAM

FILENAME ON WHICH TO STORE DATA
DEVICE:FILENAME,TYPE= SY:TEST,DEC

SUBJECT IDENTIFICATION NUMBER (XXX-XX-XXXX): ██████████

TEST NUMBER (XX): 01

TEST DATE (MM-DD-YY): 07-08-82

HOW MANY PHASES (XX): 2

INITIAL PN2 (PSI) (XX,XX): 11.2

PN2 DURING WASHOUT, PHASE I (PSI) (XX,XX): 9

DURATION OF WASHOUT, PHASE I (MIN) (XXXX): 80

PN2 DURING WASHOUT, PHASE I (PSI) (XX,XX): 8.3

DURATION OF WASHOUT, PHASE I (MIN) (XXXX): 80

FINAL PRESSURE (PSIA) (XX,XX): 3.7

TEMPERATURE IN CENTIGRADE (XX): 32

TIME OF FINAL DECOMPRESSION (HH:MM): 12:00

TIME OF RECOMPRESSION (HH:MM): 14:40

EXERCISE LEVEL

LOW 400- 800 BTU/HR
MED 800-1600 BTU/HR
HIGH 1600 BTU/HR AND ABOVE

ENTER (L/M/H): L

WERE BUBBLES MEASURED (Y/N)? Y

WERE BUBBLES PRESENT (Y/N)? Y

GRADE OF BUBBLES : ENTER 0 FOR NOT OBSERVED

TIME OF ONSET GRADE 1 (HH:MM): 13:00

TIME OF ONSET GRADE 2 (HH:MM): 0

TIME OF ONSET GRADE 3 (HH:MM): 0

TIME OF ONSET GRADE 4 (HH:MM): 0

GRADE OF BUBBLES : ENTER 0 FOR NOT OBSERVED

TIME OF REMISSION GRADE 1 (HH:MM): 15:00

TIME OF REMISSION GRADE 2 (HH:MM): 0

TIME OF REMISSION GRADE 3 (HH:MM): 0

TIME OF REMISSION GRADE 4 (HH:MM): 0

GRADE OF BUBBLES : ENTER 0 FOR NOT OBSERVED

PRESSURE AT REMISSION GRADE 1 (PSIA) (XX.XX): 10

PRESSURE AT REMISSION GRADE 2 (PSIA) (XX.XX): 0

PRESSURE AT REMISSION GRADE 3 (PSIA) (XX.XX): 0

PRESSURE AT REMISSION GRADE 4 (PSIA) (XX.XX): 0

WERE BENDS PRESENT (Y/N)? Y

BENDS LOCATION (TEXT): LEFT ELBOW

GRADE OF BENDS : ENTER 0 FOR NOT OBSERVED

TIME OF ONSET GRADE 1 (HH:MM): 14:00

TIME OF ONSET GRADE 2 (HH:MM): 0

TIME OF ONSET GRADE 3 (HH:MM): 0

TIME OF ONSET GRADE 4 (HH:MM): 0

GRADE OF BENDS : ENTER 0 FOR NOT OBSERVED

TIME OF REMISSION GRADE 1 (HH:MM): 14:50

TIME OF REMISSION GRADE 2 (HH:MM): 0

TIME OF REMISSION GRADE 3 (HH:MM): 0

TIME OF REMISSION GRADE 4 (HH:MM): 0

GRADE OF BENDS : ENTER 0 FOR NOT OBSERVED

PRESSURE AT REMISSION GRADE 1 (PSIA) (XX.XX): 10

PRESSURE AT REMISSION GRADE 2 (PSIA) (XX.XX): 0

PRESSURE AT REMISSION GRADE 3 (PSIA) (XX.XX): 0

PRESSURE AT REMISSION GRADE 4 (PSIA) (XX.XX): 0

ALIGN PAPER TO TOP OF PAGE AND HIT CARRIAGE RETURN

ORIGINAL PAGE IS
OF POOR QUALITY

DECOMPRESSION DATA OUTPUT

SUBJECT # [REDACTED] DATE 08-JUL-82

NITROGEN PRESSURE AT TIME OF FINAL DECOMPRESSION

FOR 180 MIN TISSUE : 10.00 PSI

FOR 240 MIN TISSUE : 10.24 PSI

FOR 360 MIN TISSUE : 10.52 PSI

DECOMPRESSION RATIO AT TIME OF FINAL DECOMPRESSION

FOR 180 MIN TISSUE : 1.49

FOR 240 MIN TISSUE : 1.53

FOR 360 MIN TISSUE : 1.57

6.3 NITROGEN WASHOUT PROGRAM

.RUN N2WASH

N 2 W A S H O U T

INSTRUCTIONS ?

TESTLOAD N2WASH HAS 6 PHASES

- 1 DATA CONSTANTS--ALLOWS USER
MODIFICATION OF CAL-GAS %'S. OTHERWISE
PROGRAM USES DEFAULT VALUES.
- 2 INITIALIZATION--ENTER HEADER INFORMATION
- 3 CALIBRATE MASS SPECTROMETER
- 4 CALIBRATE FLOW METER
- 5 TEST DATA--SAMPLES TEST DATA CONTINUOUSLY
UNTIL USER TERMINATES
- 6 TERMINATE--RETURNS CONTROL TO SYSTEM MONITOR

ENTER # & RETURN

PHASE 1 DATA CONSTANTS
TO LEAVE VALUE UNCHANGED, ENTER
N & RETURN. TO CHANGE, ENTER Y
& RETURN. GET "NEW VALUE" PROMPT,
ENTER NEW VALUE & RETURN

CALPCT(1,1)= 79.00
CALPCT(2,1)= 21.00
CALPCT(3,1)= 0.05
CALPCT(1,2)= 0.00
CALPCT(2,2)= 99.99
CALPCT(3,2)= 0.00
CALPCT(1,3)= 9.31Y
NEW VALUE =8.010
CALPCT(2,3)= 82.69Y
NEW VALUE =86.99
CALPCT(3,3)= 7.99Y
NEW VALUE =5.004
SAMPLE RATE= 25/SECY
NEW VALUE =15
BREATH GATE= 15.

END OF PHASE 1
ENTER NEW PHASE # 2

ORIGINAL PAGE IS
OF POOR QUALITY

PHASE 2 INITIALIZATION

DATE (MM-DD-YY): 03-09-82

SUBJECT NAME : JOHN SMITH

SUBJECT IDENTIFICATION # (XXX-XX-XXXX): [REDACTED]

SUBJECT AGE (YEARS AS OF THIS DATE): 22

EXTERNAL TEMPERATURE IN CENTIGRADE (XX.X): 35

EXTERNAL GRAVITY (1,0): 1

IS SUBJECT IMMersed (Y/N): N

ENTER CODE FOR BODY POSITION (0 FOR NO DATA OR 1-99): 0

AMBIENT PRESSURE IN MM HG (XXX.X): 760

ENTER CODE FOR DRUGS (0 FOR NO DATA OR 1-99): 0

ENTER CODE FOR EXERCISE PROTOCOL (0 FOR NO DATA OR 1-99): 0

ENTER CODE FOR BREATHING MIXTURE (0 FOR NO DATA OR 1-99): 0

ENTER CODE FOR BODY HYDRATION (0 FOR NO DATA OR 1-99): 0

ENTER CODE FOR AMBIENT GAS COMPOSITION (0 FOR NO DATA OR 1-99): 0

STARTING TIME OF DAY (HH:MM): 08:00

ENTER CODE FOR OTHER SPECIAL CONDITIONS (0 FOR NO DATA OR 1-99): 0

ARE THERE ANY COMMENTS TO RECORD(Y/N)? Y

TO ENTER COMMENTS HIT CARRIAGE RETURN THEN ENTER UP TO 4 LINES OF COMMENT

- ?THIS IS THE FIRST LINE OF A TEST
- ?THIS IS THE SECOND LINE OF A TEST
- ?THIS IS ANOTHER TEST LINE
- ?THIS IS THE LAST LINE FOR ENTERING TEXT.

END OF PHASE 2

ENTER NEW PHASE # 3

ORIGINAL PAGE IS
OF POOR QUALITY

PHASE 3 MASS SPEC CALIBRATION
VERIFY CAL GASES ARE OPEN
SET MASS SPEC FUNCTION TO "AUTO"
PAUSE -- SELECT N-2 0-10%. RETURN

PAUSE -- SELECT N-2 0-100%. RETURN

NITROGEN, HI-RANGE			
CTS-1	CTS-2	SLOPE	INTCPT
1601.49	171.30	0.050	-0.493

NITROGEN, LO-RANGE			
CTS-1	CTS-2	SLOPE	INTCPT
1651.81	-0.27	0.005	0.001

OXYGEN			
CTS-1	CTS-2	SLOPE	INTCPT
2047.00	1813.27	0.056	-13.864

CARBON DIOXIDE			
CTS-1	CTS-2	SLOPE	INTCPT
1036.52	-0.64	0.005	0.003

END OF PHASE 3

ORIGINAL PAGE IS
OF POOR QUALITY

ENTER NEW PHASE # 4

PHASE 4 FLOW METER CALIBRATION
PAUSE -- ZERO SPIROMETER.RETURN THANK YOU

SPIROMETER		SLOPE	INTCPT
CTS-1	CTS-2		
460.93	50.43	0.002	-50.426

END OF PHASE 4

ENTER NEW PHASE # 5

SELECT N-2 0-100% RANGE.
RETURN TO START.

TO STOP AT ANY TIME, ENTER S & RETURN.

1	9.0	2.51719	0.32105	0.12754
2	17.0	2.77444	0.32822	0.11830
3	29.0	3.02292	0.20593	0.06812
4	40.0	2.95325	0.13357	0.04523
5	50.0	2.46571	0.06869	0.02786

1 MIN N2(L) 1.05747 N2TOT 1.05747

SELECT N-2 0-10% RANGE

6	63.0	2.08422	0.03345	0.01605
7	72.0	2.46863	0.05014	0.02031
8	85.0	2.54334	0.03602	0.01416
9	97.0	2.02381	0.02119	0.01047
10	108.0	2.05873	0.01837	0.00892
11	120.0	2.19856	0.01782	0.00811

2 MIN N2(L) 0.17699 N2TOT 1.23446

12	128.0	1.82779	0.01440	0.00788
13	137.0	1.39287	0.01011	0.00726
14	146.0	1.56421	0.00942	0.00602
15	155.0	1.17964	0.00696	0.00590
16	157.0	1.17736	0.00011	0.00009
17	165.0	1.65856	0.00764	0.00461
18	173.0	1.78037	0.00832	0.00467

3 MIN N2(L) 0.05695 N2TOT 1.29141

19	183.0	1.56372	0.00624	0.00399
20	193.0	1.42941	0.00547	0.00383
21	201.0	1.08268	0.00432	0.00399
22	204.0	1.08333	0.00005	0.00005
23	211.0	1.45686	0.00474	0.00324
24	220.0	1.31297	0.00423	0.00322
25	228.0	1.27058	0.00403	0.00317
26	237.0	1.32840	0.00382	0.00287

4 MIN N2(L) 0.03292 N2TOT 1.32433

27	248.0	1.32580	0.00391	0.00295
28	256.0	1.22592	0.00348	0.00284
29	264.0	1.21098	0.00328	0.00271
30	273.0	1.19766	0.00338	0.00282
31	282.0	1.27091	0.00365	0.00288
32	290.0	1.34025	0.00370	0.00276
33	299.0	1.55187	0.00416	0.00268

5 MIN N2(L) 0.02556 N2TOT 1.34989

34	308.0	1.44159	0.00384	0.00266
35	316.0	1.20903	0.00307	0.00254
36	325.0	1.21488	0.00306	0.00252
37	334.0	1.27188	0.00322	0.00253
38	343.0	1.14439	0.00287	0.00251
39	350.0	0.99904	0.00230	0.00230
40	360.0	1.10623	0.00208	0.00189

6 MIN N2(L) 0.02044 N2TOT 1.37033

41	367.0	1.12036	0.00225	0.00201
42	373.0	1.24167	0.00274	0.00221
43	383.0	1.14131	0.00191	0.00168
44	390.0	1.20010	0.00261	0.00217
45	396.0	0.16559	0.00022	0.00132
46	398.0	0.16851	0.00001	0.00004
47	406.0	1.79222	0.00388	0.00216
48	414.0	1.18597	0.00273	0.00231

ORIGINAL PAGE IS
OF POOR QUALITY

ORIGINAL PAGE --
OF POOR QUALITY

ORIGINAL PAGE IS
OF POOR QUALITY

7	MIN	N2(L)	0.01635	N2TOT	1.38668
49		421.0	1.11711	0.00249	0.00223
50		430.0	1.19636	0.00246	0.00205
51		438.0	1.30323	0.00290	0.00222
52		446.0	1.13595	0.00210	0.00185
53		454.0	1.18483	0.00228	0.00193
54		465.0	1.25126	0.00208	0.00167
55		472.0	1.34025	0.00270	0.00202
56		480.0	1.20887	0.00229	0.00189
8	MIN	N2(L)	0.01930	N2TOT	1.40597
57		489.0	1.19718	0.00230	0.00192
58		496.0	1.15430	0.00216	0.00187
59		505.0	1.13238	0.00209	0.00185
60		513.0	1.24928	0.00244	0.00192
61		515.0	1.26457	0.00010	0.00008
62		522.0	1.33863	0.00241	0.00180
63		531.0	1.33928	0.00233	0.00174
64		540.0	1.27822	0.00192	0.00150
9	MIN	N2(L)	0.01576	N2TOT	1.42173
65		547.0	1.12426	0.00186	0.00166
66		555.0	1.04370	0.00175	0.00168
67		562.0	0.97712	0.00165	0.00169
		571.0	1.04906	0.00182	0.00173
		578.0	1.02259	0.00172	0.00168
			1.19571	0.00207	0.00173
			07034	0.00180	0.00168
				N2TOT	1.43440
1251				00209	0.00172
1252		10611.			0.00172
177	MIN	N2(L)	0.00174		1.96
1253		10622.0	1.44962		
1254		10632.0	1.21001	v.	
1255		10640.0	1.20221	0.00021	
1256		10653.0	1.55836	0.00032	0.
1257		10664.0	1.48544	0.00030	0.00020
1258		10673.0	1.13205	0.00028	0.00020
178	MIN	N2(L)	0.00174	N2TOT	2.17932
1259		10683.0	0.82219	0.00015	0.00018
1260		10685.0	0.00000	0.00000	0.00000
1261		10696.0	1.76007	0.00043	0.00025
1262		10709.0	1.52344	0.00031	0.00021
1263		10721.0	1.42649	0.00029	0.00020
1264		10731.0	1.27074	0.00026	0.00020
1265		10733.0	0.08211	0.00000	0.00003
179	MIN	N2(L)	0.00145	N2TOT	2.18077
1266		10742.0	0.87919	0.00016	0.00018
1267		10752.0	1.79904	0.00051	0.00028
1268		10763.0	1.34236	0.00030	0.00023
1269		10774.0	2.00254	0.00047	0.00024
1270		10785.0	1.68634	0.00034	0.00020
1271		10797.0	1.75503	0.00036	0.00021
S180	MIN	N2(L)	0.00215	N2TOT	2.18292

END OF PHASE 5
ENTER NEW PHASE # 6

STOP -- BYE

ORIGINAL PAGE IS
OF POOR QUALITY

6.4 NITROGEN DATA PROGRAM

RUN BY: N2DATA

NITROGEN WASHOUT DATA ANALYSIS
ENTER FILENAME ON WHICH DATA IS STORED
(DEVICE:FILENAME.TYPE): DY1:IERWW1.N2W

A FULL REPORT (F) GIVES BREATH BY BREATH PARAMETERS
A 1 MINUTE REPORT (1) GIVES PARAMETERS FOR 1 MINUTE INTERVALS
A 5 MINUTE REPORT (5) GIVES PARAMETERS FOR 5 MINUTE INTERVALS
A SHORT REPORT (S) GIVES THE OVERALL AVERAGES

DO YOU WANT A FULL REPORT (F), A 1 MINUTE REPORT (1),
A 5 MINUTE REPORT (5), OR A SHORT REPORT (S) ?
TYPE (F/1/5/S): F

PAUSE --- ADVANCE CARRIAGE IF DESIRED. RETURN

N2DATA REPORT

IERWW1		AGE: 0	DATE: 1-DEC-81
SUBJECT ID NUMBER:		STARTING TIME:	
EXTERNAL TEMPERATURE:	0.00	EXTERNAL GRAVITY:	0
IMMERSTION?	NO	BODY POSITION CODE:	0
AMBIENT PRESSURE:	0.00	DRUG CODE:	0
EXERCISE CODE:	0	BREATH MIXTURE CODE:	0
BODY HYDRATION CODE:	0	AMBIENT GAS CODE:	0
SPECIAL CONDITION CODE:	0	COMMENTS?	NO
COMMENTS:			

ORIGINAL PAGE IS
OF POOR QUALITY

BREATH #	ELAPSED TIME (SEC)	BREATH VOLUME (L)	N2 VOL /BREATH (L)	O2 CONC /BREATH	N2 CONC /BREATH	CO2 CONC /BREATH
1	9.0	2.51719	0.32105	0.81261	0.12754	0.03284
2	17.0	2.77444	0.32822	0.81574	0.11830	0.03991
3	29.0	3.02292	0.20593	0.88528	0.06812	0.05382
4	40.0	2.95325	0.13357	0.90448	0.04523	0.03693
5	50.0	2.46571	0.06869	0.92569	0.02786	0.03502
6	63.0	2.08422	0.03345	0.94523	0.01605	0.02927
7	72.0	2.46863	0.05014	0.93448	0.02031	0.03974
8	85.0	2.54334	0.03602	0.94588	0.01416	0.03561
9	97.0	2.02381	0.02119	0.95183	0.01047	0.03204
10	108.0	2.05873	0.01837	0.95308	0.00892	0.03265
11	120.0	2.19856	0.01782	0.94887	0.00811	0.03480
12	128.0	1.82779	0.01440	0.94779	0.00788	0.03720
13	137.0	1.39287	0.01011	0.94576	0.00726	0.03748
14	146.0	1.56421	0.00942	0.95191	0.00602	0.03308
15	155.0	1.17964	0.00696	0.94392	0.00590	0.03781
16	157.0	1.17736	0.00011	0.89991	0.00009	0.00064
17	165.0	1.65856	0.00764	0.95518	0.00461	0.03209
18	173.0	1.78037	0.00832	0.95117	0.00467	0.03601
19	183.0	1.56372	0.00624	0.95605	0.00399	0.03333
	193.0	1.42941	0.00547	0.95516	0.00383	0.03344
		0.08268	0.00432	0.94535	0.00399	0.03817
			0.00005	0.88880	0.00005	0.00093
				0.95675	0.00326	0.03187
				0.95177	0.00322	0.03221
1258	10683.0				0.00317	0.03512
1259	10683.0	0.00000				0.03327
1260	10685.0	0.00000				
1261	10696.0	1.76007	0.00045			
1262	10709.0	1.52344	0.00031	0.96111		
1263	10721.0	1.42649	0.00029	0.96334	0.00000	
1264	10731.0	1.27074	0.00026	0.95785	0.00020	
1265	10733.0	0.08211	0.00000	0.79992	0.00003	0.00120
1266	10742.0	0.87919	0.00016	0.96392	0.00018	0.02311
1267	10752.0	1.79904	0.00051	0.96155	0.00028	0.03352
1268	10763.0	1.34236	0.00030	0.96183	0.00023	0.03021
1269	10774.0	2.00254	0.00047	0.96173	0.00024	0.03223
1270	10785.0	1.68634	0.00034	0.95896	0.00020	0.03037
1271	10797.0	1.75503	0.00036	0.96345	0.00021	0.02912

AVERAGES:

BREATH VOLUME (L)	N2 VOLUME /BREATH (L)	O2 CONC /BREATH	N2 CONC /BREATH	CO2 CONC /BREATH
1.27147	0.00172	0.94703	0.00095	0.03036

TOTAL N2 VOLUME IN LITERS= 2.18273

TOTAL # OF BREATHS= 1271

END OF DATA..READ ANOTHER FILE (Y/N): Y

ORIGINAL PAGE IS
OF POOR QUALITY

NITROGEN WASHOUT DATA ANALYSIS
ENTER FILENAME ON WHICH DATA IS STORED
(DEVICE:FILENAME.TYPE): DY1:IERWW1.N2W

A FULL REPORT (F) GIVES BREATH BY BREATH PARAMETERS
A 1 MINUTE REPORT (1) GIVES PARAMETERS FOR 1 MINUTE INTERVALS
A 5 MINUTE REPORT (5) GIVES PARAMETERS FOR 5 MINUTE INTERVALS
A SHORT REPORT (S) GIVES THE OVERALL AVERAGES

DO YOU WANT A FULL REPORT (F), A 1 MINUTE REPORT (1),
A 5 MINUTE REPORT (5), OR A SHORT REPORT (S) ?
TYPE (F/1/5/S): 1

PAUSE -- ADVANCE CARRIAGE IF DESIRED. RETURN

N2DATA REPORT

IERWW1 AGE: 0 DATE: 1-DEC-81
SUBJECT ID NUMBER: STARTING TIME:
EXTERNAL TEMPERATURE: 0.00 EXTERNAL GRAVITY: 0
IMMERSION? NO BODY POSITION CODE: 0
AMBIENT PRESSURE: 0.00 DRUG CODE: 0
EXERCISE CODE: 0 BREATH MIXTURE CODE: 0
BODY HYDRATION CODE: 0 AMBIENT GAS CODE: 0
SPECIAL CONDITION CODE: 0 COMMENTS? NO
COMMENTS:

ENTER THE SUBJECT'S BODY WEIGHT IN kg: 77.2

ENTER THE SUBJECT'S ESTIMATED N2 CONTENT IN LITERS: 2.3

TIME (min)	N2 VOL/ INTERVAL (ml)	O2 VOL/ INTERVAL (ml)	CO2 VOL/ INTERVAL (ml)	N2/WEIGHT (ml/kg)	N2 CONC/ INTERVAL	CUMMUL. N2 (ml)	% OF EST. N2 CONT. (%)	BMR (BTU/hr)
1	1057.46	11938.49	491.04	13.69767	0.07741	1057.46	45.98	668.99
2	176.99	12657.25	458.25	2.29262	0.01300	1234.45	53.67	624.32
3	56.96	9989.33	334.63	0.73782	0.00520	1291.41	56.15	455.90
4	32.90	9966.83	319.79	0.42617	0.00305	1324.31	57.58	435.88
5	25.56	8694.38	307.06	0.33109	0.00281	1349.87	58.69	418.34
6	20.44	7980.53	279.52	0.26477	0.00242	1370.31	59.58	380.82
7	16.35	7618.70	249.28	0.21179	0.00174	1386.66	60.29	339.63
8	19.30	9279.59	321.17	0.25000	0.00198	1405.96	61.13	437.56
9	15.75	9319.81	284.61	0.20402	0.00159	1421.71	61.81	387.76
10	12.67	7109.67	247.16	0.16412	0.00169	1434.38	62.36	336.73
11	14.64	8037.31	284.84	0.18964	0.00173	1449.02	63.00	388.07
12	14.31	8019.60	277.13	0.18536	0.00170	1463.33	63.62	377.57

13	12.08	7249.82	247.84	0.15648	0.00158	1475.41	64.15	337.66
14	12.62	8022.86	285.11	0.16347	0.00149	1488.03	64.70	388.43
15	10.45	6435.98	228.16	0.13536	0.00154	1498.48	65.15	310.84
16	12.50	8785.69	308.66	0.16192	0.00136	1510.98	65.69	420.51
17	11.20	8726.81	299.41	0.14508	0.00123	1522.18	66.18	408.20
18	10.47	8156.81	287.66	0.13562	0.00122	1532.65	66.64	391.91
19	10.93	8746.92	297.07	0.14158	0.00119	1543.58	67.11	404.73
20	9.29	7559.20	267.16	0.12034	0.00117	1552.87	67.52	363.98
21	10.61	10228.19	344.54	0.13744	0.00098	1563.48	67.98	469.41
22	8.56	8421.31	290.79	0.11088	0.00097	1572.04	68.35	396.17
23	7.93	8158.42	273.69	0.10272	0.00092	1579.97	68.69	372.87
24	7.88	9452.10	281.93	0.10207	0.00079	1587.85	69.04	384.10
25	8.98	9396.11	320.56	0.11632	0.00091	1596.83	69.43	436.73
26	6.77	6709.58	233.98	0.08769	0.00096	1603.60	69.72	318.78
27	7.45	8046.30	258.78	0.09650	0.00088	1611.05	70.05	366.19
28	7.81	9113.00	305.07	0.10117	0.00082	1618.86	70.39	415.63
29	6.32	7008.47	225.62	0.08187	0.00086	1625.18	70.66	307.39
30	9.64	11420.86	362.66	0.12487	0.00081	1634.82	71.08	494.09
31	11.23	14616.33	462.58	0.14547	0.00073	1646.05	71.57	630.22
32	10.67	11642.66	421.87	0.13821	0.00088	1656.72	72.03	574.76
33	8.68	9911.68	320.73	0.11244	0.00081	1665.40	72.41	436.96
34	6.80	7733.59	268.99	0.08808	0.00085	1672.20	72.70	366.47
35	6.37	7373.51	263.86	0.08251	0.00082	1678.57	72.98	359.48
36	6.24	6916.55	247.06	0.08083	0.00086	1684.81	73.25	336.59
37	7.34	7512.04	244.87	0.08212	0.00081	1691.15	73.53	333.62
38	7.27	7077.26	216.42	0.06813	0.00079	1696.41	73.76	294.85
161			260.60	0.07850	0.00077	1702.47	74.02	355.04
162	1.82			0.07176	0.00064	1708.01	74.26	368.85
163	2.19	8122.80			0.00071	1713.33	74.49	327.49
164	2.29	8234.90	267.97			1719.88	74.78	393.75
165	1.99	7510.80	237.47	0.02570			75.04	384.78
166	1.97	7009.65	220.03	0.02552	0.00000			359.74
167	2.14	8772.87	270.56	0.02772	0.00023	2157.00		360.00
168	1.89	8361.51	242.34	0.02448	0.00021	2159.76		
169	2.27	9930.67	311.13	0.02940	0.00022	2162.03	94.00	
170	2.04	9507.63	265.19	0.02642	0.00020	2164.07	94.09	361.29
171	1.81	7592.87	229.31	0.02345	0.00023	2165.88	94.17	312.42
172	2.31	9223.87	280.34	0.02992	0.00024	2168.19	94.27	381.93
173	2.36	10104.62	304.94	0.03057	0.00022	2170.55	94.37	415.45
174	1.83	8429.64	243.98	0.02370	0.00021	2172.38	94.45	332.40
175	1.90	9398.94	291.70	0.02461	0.00019	2174.28	94.53	397.41
176	1.72	9246.69	291.06	0.02228	0.00018	2176.00	94.61	382.91
177	1.41	6099.01	179.05	0.01826	0.00018	2177.41	94.67	243.93
178	1.74	7718.02	241.92	0.02254	0.00022	2179.15	94.75	329.60
179	1.44	6590.36	209.12	0.01865	0.00015	2180.59	94.81	284.90

AVERAGES:

BREATH VOLUME(L)	N2 VOLUME /BREATH(L)	O2 CONC /BREATH	N2 CONC /BREATH	CO2 CONC /BREATH
1.27147	0.00172	0.94703	0.00095	0.03036

TOTAL N2 VOLUME IN LITERS= 2.18273

TOTAL # OF BREATHS= 1271

END OF DATA..READ ANOTHER FILE (Y/N): Y

ORIGINAL PAGE IS
OF POOR QUALITY

ORIGINAL PAGE IS
OF POOR QUALITY

NITROGEN WASHOUT DATA ANALYSIS
ENTER FILENAME ON WHICH DATA IS STORED
(DEVICE:FILENAME.TYPE): DY1:IERWW1.N2W

A FULL REPORT (F) GIVES BREATH BY BREATH PARAMETERS
A 1 MINUTE REPORT (1) GIVES PARAMETERS FOR 1 MINUTE INTERVALS
A 5 MINUTE REPORT (5) GIVES PARAMETERS FOR 5 MINUTE INTERVALS
A SHORT REPORT (S) GIVES THE OVERALL AVERAGES

DO YOU WANT A FULL REPORT (F), A 1 MINUTE REPORT (1),
A 5 MINUTE REPORT (5), OR A SHORT REPORT (S) ?
TYPE (F/1/5/S): 5

PAUSE -- ADVANCE CARRIAGE IF DESIRED. RETURN

N2DATA REPORT

IERWW1		AGE: 0	DATE: 1-DEC-81
SUBJECT ID NUMBER:		STARTING TIME:	
EXTERNAL TEMPERATURE:	0.00	EXTERNAL GRAVITY:	0
IMMERSION?	NO	BODY POSITION CODE:	0
AMBIENT PRESSURE:	0.00	DRUG CODE:	0
EXERCISE CODE:	0	BREATH MIXTURE CODE:	0
BODY HYDRATION CODE:	0	AMBIENT GAS CODE:	0
SPECIAL CONDITION CODE:	0	COMMENTS?	NO
COMMENTS:			

ORIGINAL PAGE IS
OF POOR QUALITY

ENTER THE SUBJECT'S BODY WEIGHT IN kg: 77.2

ENTER THE SUBJECT'S ESTIMATED N2 CONTENT IN LITERS: 2.3

TIME (min)	N2 VOL/ INTERVAL (ml)	O2 VOL/ INTERVAL (ml)	CO2 VOL/ INTERVAL (ml)	N2/WEIGHT (ml/kg)	N2 CONC/ INTERVAL	CUMMUL. N2 (ml)	% OF EST. N2 CONT. (%)	BMR (BTU/hr)
5	1349.87	53246.28	1910.77	17.48537	0.01653	1349.87	58.69	520.65
10	84.51	41308.30	1381.75	1.09469	0.00187	1434.38	62.36	376.50
15	64.10	37765.59	1323.08	0.83031	0.00161	1498.48	65.15	360.51
20	54.39	41975.43	1460.17	0.70453	0.00123	1552.87	67.52	397.87
25	43.96	45656.13	1511.51	0.56943	0.00091	1596.83	69.43	411.86
30	37.99	42298.21	1396.12	0.49210	0.00087	1634.82	71.08	380.41
35	43.75	51277.77	1738.03	0.56671	0.00081	1678.57	72.98	473.58
40	29.44	36199.19	1239.68	0.38135	0.00077	1708.01	74.26	337.79
45	41.90	41562.05	1353.97	0.54275	0.00100	1749.91	76.08	368.93
50	30.13	41909.43	1349.61	0.39028	0.00065	1780.04	77.39	367.74
55	26.74	46125.07	1452.00	0.34637	0.00054	1806.78	78.56	395.64
60	23.51	43617.86	1393.54	0.30453	0.00051	1830.29	79.58	379.71
65	24.94	40958.88	1405.17	0.32306	0.00058	1855.23	80.66	382.88
70	19.91	37431.40	1269.87	0.25790	0.00051	1875.14	81.53	346.01
75	18.88	39210.10	1230.53	0.24456	0.00046	1894.02	82.35	335.30
80	18.62	40979.50	1244.67	0.24119	0.00043	1912.64	83.16	339.15
85	20.39	45240.64	1426.25	0.26412	0.00043	1933.03	84.04	388.63
90	21.10	54872.51	1620.45	0.27332	0.00038	1954.13	84.96	441.54
95	17.45	41760.23	1213.42	0.22604	0.00040	1971.58	85.72	330.63
100	16.98	42485.95	1355.14	0.21995	0.00037	1988.56	86.46	369.25
105	14.79	40231.99	1347.70	0.19158	0.00035	2003.35	87.10	367.22
110	15.91	43093.18	1414.17	0.20609	0.00035	2019.26	87.79	385.33
115	15.45	41953.68	1310.54	0.20013	0.00036	2034.71	88.47	357.09
120	14.87	45553.75	1400.65	0.19262	0.00030	2049.58	89.11	381.65
125	13.82	47164.87	1477.40	0.17902	0.00028	2063.40	89.71	402.56
130	12.19	39932.64	1288.64	0.15790	0.00029	2075.59	90.24	351.13
135	11.67	41307.94	1264.48	0.15117	0.00026	2087.26	90.75	344.54
140	11.05	39403.11	1274.23	0.14313	0.00027	2098.31	91.23	347.20
145	11.13	40149.90	1334.43	0.14417	0.00026	2109.44	91.71	363.61
150	10.94	44346.65	1327.57	0.14171	0.00022	2120.38	92.19	361.74
155	12.65	48831.21	1485.62	0.16386	0.00025	2133.03	92.74	404.80
160	10.70	40833.64	1302.66	0.13860	0.00024	2143.73	93.21	354.95
165	10.03	36201.57	1179.06	0.12992	0.00026	2153.76	93.64	321.27
170	10.31	43582.33	1309.25	0.13355	0.00022	2164.07	94.09	356.74
175	10.21	44749.94	1350.27	0.13225	0.00022	2174.28	94.53	367.92

AVERAGES:

BREATH VOLUME(L)	N2 VOLUME /BREATH(L)	O2 CONC /BREATH	N2 CONC /BREATH	CO2 CONC /BREATH
1.27147	0.00172	0.94703	0.00095	0.03036

TOTAL N2 VOLUME IN LITERS= 2.18273

TOTAL # OF BREATHS= 1271

END OF DATA..READ ANOTHER FILE (Y/N): Y

NITROGEN WASHOUT DATA ANALYSIS
ENTER FILENAME ON WHICH DATA IS STORED
(DEVICE:FILENAME.TYPE): DY1:IERWW1.N2W

A FULL REPORT (F) GIVES BREATH BY BREATH PARAMETERS
A 1 MINUTE REPORT (1) GIVES PARAMETERS FOR 1 MINUTE INTERVALS
A 5 MINUTE REPORT (5) GIVES PARAMETERS FOR 5 MINUTE INTERVALS
A SHORT REPORT (S) GIVES THE OVERALL AVERAGES

DO YOU WANT A FULL REPORT (F), A 1 MINUTE REPORT (1),
A 5 MINUTE REPORT (5), OR A SHORT REPORT (S) ?
TYPE (F/1/5/S): S

PAUSE -- ADVANCE CARRIAGE IF DESIRED. RETURN

N2DATA REPORT

IERWW1	AGE: 0	DATE: 1-DEC-81
SUBJECT ID NUMBER:	STARTING TIME:	
EXTERNAL TEMPERATURE: 0.00	EXTERNAL GRAVITY: 0	
IMMERSION? NO	BODY POSITION CODE: 0	
AMBIENT PRESSURE: 0.00	DRUG CODE: 0	
EXERCISE CODE: 0	BREATH MIXTURE CODE: 0	
BODY HYDRATION CODE: 0	AMBIENT GAS CODE: 0	
SPECIAL CONDITION CODE: 0	COMMENTS? NO	
COMMENTS:		

AVERAGES:

BREATH VOLUME(L)	N2 VOLUME /BREATH(L)	O2 CONC /BREATH	N2 CONC /BREATH	CO2 CONC /BREATH
1.27147	0.00172	0.94703	0.00095	0.03036
TOTAL N2 VOLUME IN LITERS= 2.18273				
TOTAL # OF BREATHS= 1271				

END OF DATA..READ ANOTHER FILE (Y/N): N

STOP --

ORIGINAL PAGE IS
OF POOR QUALITY

6.5 UPDA TE PROGRAM

R UPDATE

ENTER FILENAME ON WHICH DATA IS STORED
(DEVICE:FILENAME,TYPE): DY1:ENTNW1.N2W

ENTER FILENAME ON WHICH DATA IS TO BE PLACED
(DEVICE:FILENAME,TYPE): SY:TEST.UPD

STOP --

7.0 HARDWARE OPERATION

7.1 OVERVIEW

The volume of nitrogen eliminated from a subject is measured on a breath by breath basis by carefully measuring the volume and the gas composition of each expired breath. A Perkin-Elmer Medical Gas Analyzer 1100 Mass Spectrometer is used to measure the concentrations of nitrogen, oxygen, and carbon dioxide (in percent of total gas concentration) in expired breaths and the nitrogen concentration of the inspired air. The volume of nitrogen eliminated is equal to the breath volume times the concentration of nitrogen in the expired air less the volume of nitrogen inspired (this should be zero when breathing pure oxygen). The mass spectrometer has two nitrogen range settings, 0-100% and 0-10% to allow for greater resolution of nitrogen concentration at low levels. The output signal for each channel is 0-10 volts, including the two nitrogen ranges.

Breath volume is determined by a Skylab exhalation spirometer. The spirometer is a dry roller seal type which outputs an increasing voltage (approximately 0-7 volts) as the measured breath volume increases. The maximum measurable volume is five liters. After each breath, a five volt signal from the computer activates a solenoid which dumps the spirometer. The dump solenoid is gas driven and is connected to a high pressure (160 PSI) nitrogen source. A return spring empties the spirometer chamber, forcing over four liters of air out in less than 0.25 seconds. The breath is determined to be over when five consecutive samples of the spirometer input have not increased. The accuracy of the spirometer has been measured and is within 2% of the actual volume (1.000 liter reads 1.000 ± 0.0201).

Calibration of the mass spectrometer is performed by analyzing calibration gases of known purity and concentrations and determining slope and intercept values for each of the three gases analyzed (N_2 , CO_2 , and O_2). Real time analysis solves the linear equation y (% concentration) = slope \times mass spectrometer output + b (intercept). Three N_2 concentrations are used to provide high and low range calibration curves for each N_2 range.

The LSI-11/02 measures the gas concentrations and breath volume by means of an ADAC 1030 analog to digital converter. The ADAC 1030 has 8 channels for differential analog input signals. Sampling of the analog input determines a digital output to the computer of an integer from 0-2047. The gain and channel selection is under computer software control. It is this count which is converted to the appropriate measurement (% or liters) by means of the slope and intercept determined by the calibration procedure.

The spirometer is calibrated by emptying the spirometer and having the computer sample the zero volume output. A one liter volume is simulated by the computer and determines the high calibration volume. Slope and intercept values are then calculated from the zero and one liter values.

Breath volume is continuously updated during a breath; a breath being defined as a spirometer analog output being greater than a minimum noise gate over the zero volume output. Consecutive volume samples are compared to each other to determine if the breath volume is still increasing. If

there is no additional volume measured between five samplings, the breath is considered to be over. The computer then loads one ADAC 1030 digital to analog channel with a 5 volt signal to the spirometer controller that dumps the spirometer. The spirometer purges to the ambient air.

The parallel interface board (ADAC 1632TTL) output controls the three 12 volt calibration gas solenoids and valves, and the LED breath indicators. The input buffer reads the mass spectrometer nitrogen range select switch and inlet valve select position.

The MDB KW11P Programmable Clock board regulates the interrupt timings for the subroutine SAMPLE. The two MDB DLV11 serial interface boards allow the use of either the printing terminal or a CRT terminal.

7.2 MASS SPECTROMETER

The Perkin-Elmer MGA 1100 Mass Spectrometer has been found to have a severe baseline drift (approximately $\pm 0.5\%$) due to a nitrogen exposure problem. When subjected to more than several minutes of gas samples with a high (ambient) N_2 concentration, several hours may be required to read a true 0.00% nitrogen when pure oxygen is being sampled. It is necessary to flush pure oxygen into the mass spectrometer for at least two days prior to and up to the time of a test. With this flushing, a pure oxygen source will read 0.00% N_2 . The 30 seconds of nitrogen during the calibration does not affect the baseline drift for more than the first fifteen minutes of the washout, during which time the drift is a small (less than 2%) fraction of the expired breath nitrogen concentration. If the mass spectrometer is not adequately flushed out, pure oxygen will show up to 0.5 - 0.6% nitrogen, which is up to twenty times the value that may actually exist after a three hour washout. The software subtracts out low level nitrogen baseline drifts measured in the inspired air, but cannot cancel the effect of large drifts that overwhelm the desired measurement. Following lengthy oxygen flushing and minimizing the mass spectrometer exposure to nitrogen has given repeatable three hour washouts on a subject that are within 50 ml. of nitrogen of each other.

Typically, an easy method to flush the mass spectrometer is to insert the face mask catheter (inlet #1) into the oxygen tent O_2 line set at a low flow. The inlet select switch must be switched to inlet #2 during calibration and returned to inlet #1 for the test and no longer than one or two minutes of ambient nitrogen exposure will occur.

7.3 CALIBRATION VALUES

A read only data file named DA TCON.VAL is stored on the system floppy (DY0:), and contains all the initial calibration values used by the main program N2WASH.SAV. The three gas mixture calibration gas tank should be 9% N2, 86% O2, and 5% CO2 (all + 1 %). Certification must be precise to at least 20 parts per million (0.002%). The first three values in DA TCON.VAL represent the concentrations of N2 in ambient air, in 100% O2 and in the three gas mixture. The next three values in DA TCON.VAL represent the concentrations of O2 in ambient air, in 100% O2 and in the three gas mixture. The next three values in DA TCON.VAL represent the concentrations of CO2 in ambient air, in 100% O2 and in the three gas mixture. These nine values, all separated by commas, are followed by the sample rate and the breath gate.

Each time a new tank of three gas calibration mixture gas is received in the Environmental Laboratory it will be necessary to change the gas percentages in DA TCON.VAL.

In order to change the values listed in DA TCON.VAL (originally 79.0,0.0,9.313,21.0,99.99,82.69,0.05,0.00,7.99,25,15.0),

Type

R LEDITV

The computer will type:

```
LINE-ORIENTED TEXT EDITOR
VERSION 2.0/A
FILE (DEVICE:FILENAME.TYPE) =
```

Type

DYU:DA TCON.VAL

The computer will type:

NEW FILE (Y/N)?

Type

N (for no)

The computer will type:

```
FILE:DA TCON.VAL           CONTAINS 1 LINES
BEGIN TEXT EDITING
```

Type

L (to list the line)

The computer will type the line. The N2, O2, and CO2 in the three gas mixture are the 3rd, 6th, and 9th values, respectively.

To replace the current value with a new value type:

RS:/(current value)/,/(new value)/

To check that the new value has been entered correctly type:

L (to list the line)

Repeat the previous 2 steps for each value which needs to be replaced. When you are satisfied that the line is correct with the new values, type:

END

This will make the changes permanent. Typing stop will abort the editing process and leave DATA.CON.VAL in its original form.

7.4 OPERATING THE N2WASH PROGRAM

N2WASH collects all the header information prior to the washout, prompting the test operator to enter all the information. Next, it controls and calibrates all of the subsystems necessary to the measurement of the washout, determining slopes and intercepts from calibration values. Finally, the program reads all the data during the experiment, performs some real time analysis and loads the data to the floppy disk and prints out the measured parameters on a breath by breath basis. One minute summary reports are also typed out.

This report assumes that the test operator is familiar with the general operation of the DEC LSI-RT11 operating system. If this is not so, refer to the APPENDIX: RT11 COMMANDS.

Before starting the program, the following must be available and, where appropriate, connected correctly (see Interconnections):

- Perkin-Elmer Mass Spectrometer MGA 1100
- Skylab Exhalation Spirometer and control electronics
- LSI-11/02 Computer console w/floppy disk drive and printing terminal
- Breathing O2 gas
- Calibration gases (100% O2, and three gas mixture)
- Mask assembly and hoses to spirometer
- Floppy disk containing N2WASH.SAV on RT-11 version 4.0 Operating System (System floppy)
- An initialized, double density, single sided floppy disk containing at least 400 free blocks for 3 hours worth of data (Data Floppy)
- N2 gas with approximately 160 psi (depends on the actuating pressure listed on the air diverter 3-way valve) to drive spirometer dump solenoid

1. Turn power on the following, (verify with on lights)

- Computer console
- LSI-11/02 computer (Line Time Clock, Run/Halt switches in up position)
- Disk drive (write protect off)
- Spirometer controller
- Mass spectrometer (on switch depressed)
- Printing terminal

2. Open all calibration gas valves (low flow, approximately 5 psi)

3. Instructions

- a. Insert the system floppy into the left hand slot DY1: with the label on the top. The exposed portion is inserted first. Push the floppy in all the way and close the door.

- b. Insert the data floppy in the right hand slot DY1: in the same manner and close the door. (Note: DY1: is vital to program startup).
- c. The disk drive will engage (with a series of audible clicks and the door LED flashing) and the terminal will respond by printing:

```
RT-11 SJ V04.00A  
.D 56 = 1012  
SET TT:SCOPE  
ASS DY1: DK  
TYPE SY:REMIND.TXT  
ENTER DATE AND TIME
```

Enter the current date and time by typing:

DATE dd-mmm-yy (where dd=day 0-31, mmm=month three letters of the month JAN-DEC and yy=year 81,82..., separated by hyphens).

Press Return

Type TIME hh:mm (where hh=hours 00-24 and mm=minutes 00-59, separated by a colon).

Press Return

Verify the time by typing TIME and return, and the computer will respond by printing the updated time.


To switch the printing terminal to a CRT terminal, one must type

R TTY and press return.

The prompt character will then appear on the CRT (if turned on and connected). There will be no hard copy of the test generated unless R TTY is again typed in to return control to the printing terminal. This must be done while under monitor control, not in the N2WASH program.

8.0 ENGINEERING DIAGRAMS AND TABLES


NITROGEN WASHOUT SYSTEM	173
BLOCK DIAGRAM OF N2WASHOUT HARDWARE	174
BOARD CONFIGURATION OF LSI-1102 MICROCOMPUTER	175
N2WASHOUT ELECTRONICS CARD CONNECTOR DESIGNATIONS	176
A/D CONVERTER	179
DRV/11 PARALLEL BOARD	180
FACE MASK AND GAS PLUMBING	182
POWER SUPPLIES FOR N2WASHOUT SYSTEM ELECTRONICS	183
CALIBRATION GAS SELECT CIRCUIT	184
MASS SPECTROMETER INLET SELECT CIRCUIT	185
MODIFICATIONS TO MODEL 1100 MASS SPECTROMETER	186
SPIROMETER LED INDICATORS	187
SPIROMETER CALIBRATION SWITCHES	188
SPIROMETER DUMP VALVE ELECTRONICS	189
SPIROMETER ELECTRONICS CHASSIS WIRING	190
SPIROMETER PRIM/FAIL SAFE BOARD	191
SPIROMETER POWER SUPPLY AND VOLUME READOUT	192



ORIGINAL PAGE IS
OF POOR QUALITY

QTY	CODE	PART OR	MANUFACTURE	
REQD	IDENT	IDENTIFYING NO.	OR DESCRIPTION	
PARTS LIST				
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE		CONTRACT NO.	Technology Incorporated LIFE SCIENCES DIVISION HOUSTON, TEXAS	
FUNCTIONS	REVISIONS	APPROVALS	NITROGEN WASHOUT SYSTEM	
1	AS SHOWN	DATE	SHEET NO. 1 OF 2	
2			DRAWN / INCH / PARTS / 1/27/67	
3		CHECKED	THRIB7-OHOI-J300	
MATERIAL			SCALE N/A	
FINISH			SHEET NO. 2 OF 2	
NEXT ASST	URGE OR			
APPLICATION				

ORIGINAL PAGE IS
OF POOR QUALITY



REVISIONS		DATE	APPROVED
FORM	LIB		
QTY REQD		CODE IDENT	PART OR IDENTIFYING NO
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES FRACTIONS DECIMALS ANGLES		CONTRACT NO	MANUFACTURE OR DESCRIPTION
1/2	1/16	1/16	1/16
MATERIAL		DATE	DATE
INSTR		APPROVALS	DATE
BEST ABST		GREEN, LINDA	10-14-71
USED (M)		CHECKED	
APPLICATION		SIZ	DOC IDENT NO
		C	THRU67-1E03-ED000
		SCALE	N/A
			SHEET 3 OF 13

technology incorporated
LIFE SCIENCES DIVISION HOUSTON, TEXAS

BLOCK DIAGRAM OF N₂
WASHOUT HARDWARE

BOARD CONFIGURATION OF THE LSI-11 MICROCOMPUTER

<u>Slot Location</u>	<u>Manufacturer</u>	<u>Device #</u>	<u>Description</u>
1A	DEC	M7270	LSI-11/2 microprocessor
1B	DEC	M8044	32K memory
2A	MDB	DLV11	Serial interface (1)
2B	Data Systems	DSD480	Floppy disk interface
3A & B	MDB	KW11P	Programmable clock
4A & B	ADAC	1030	A/D, D/A Programmable Gain
5A	ADAC	1632 TTL	Parallel interface
6A	MDB	DLV11	Serial interface (2)
6B-7B	EMPTY		
8A	MDB	SMU	Line clock/system monitoring unit
	MDB	M2689A	Power Supply 5/12/12
	Data Systems	480	Dual Floppy Disc Drive
	ITT	43	Teletype
	Hazeltine	2000	CRT

ORIGINAL PAGE IS
OF POOR QUALITY

N₂ WASHOUT ELECTRONICS CARD CONNECTOR DESIGNATIONS

<u>Top of Plug</u>		2	1		
		4	3		
		6	5		
		8	7		
		10	9		
		12	11		
		14	13	<u>Brn #2</u>	S2 No.
Clear Ind. Sig.	<u>Cable Blk</u>	16	15	<u>Red #2</u>	S1 No.
		18	17	<u>Rib.Cable</u>	DRV11 J2 Pin 29 (HH) In 5 Flow Meter Cal Done
S1 N.C. Cal Done	<u>Yel #2</u>	20	19	<u>Or #2</u>	S2 N.C. Hi Cal Done
		22	21	<u>Rib.Cable</u>	DRV11 J1 Pin 10 (L) Out 4 Inlet Select
Flow Meter Hi Cal Done DRV11 J1 Pin 31 (KK) In 4	<u>Rib Cable</u>	24	23	<u>Cable Yel</u>	Select Inlet #2
		26	25	<u>Cable Red</u>	24 VDC Grnd. from Mass Spec
Select Inlet #1	<u>Cable Wh</u>	28	27	<u>Rib.Cable</u>	DRV11 J1 Pin 12 (N) Out 5 Do Hi Cal L.E.D.
		30	29	<u>Sol. Wire</u>	Cal #3 83% O ₂ Solenoid
Do Cal L.E.D. DRV11 J1 Pin 14 (R) Out 6	<u>Rib Cable</u>	32	31	<u>Gry</u>	Red L.E.D. L3 83% O ₂
		34	33		
Capacitor C5 3K FD	<u>Brn</u>	36	35	<u>Blk</u>	Grnd
		38	37	<u>Blk</u>	
		40	39		
		42	41	<u>Lt. Grn</u>	Transformer
Transformer	<u>Lt. Grn.</u>	44	43		
		46	45	<u>Red</u>	VR #1 Pin 2
+5V D.C. VR #1 Pin 1	<u>Grn #2</u>	48	47	<u>Red</u>	3K FD Capacitor C5 +5 V.D.C.
		50	49	<u>Red</u>	VR #1 Pin 2

ORIGINAL PAGE IS
OF POOR QUALITY

		52	51	<u>Or</u>	VR #2 Pin 2
		54	53	<u>Or</u>	VR #2 Pin 2 +12 V.D.C.
Capacitor 3K FD C6 +12 V.D.C.	<u>Red</u>	56	55	<u>Or</u>	VR #2 Pin 1
		10	9		
		12	11		
		14	13	<u>Brn #2</u>	S2 No.
Clear Ind. Sig.	<u>Cable Blk</u>	16	15	<u>Red #2</u>	S1 No.
		18	17	<u>Rib.Cable</u>	DRV11 J2 Pin 29 (HH) In 5 Flow Meter Cal Done
S1 N.C. Cal Done	<u>Yel #2</u>	20	19	<u>Or #2</u>	S2 N.C. Hi Cal Done
		22	21	<u>Rib.Cable</u>	DRV11 J1 Pin 10 (L) Out 4 Inlet Select
Flow Meter Hi Cal Done DRV11 J1 Pin 31 (KK) In 4	<u>Rib Cable</u>	24	23	<u>Cable Yel</u>	Select Inlet #2
		26	25	<u>Cable Red</u>	24 VDC Grnd. from Mass Spec
Select Inlet #1	<u>Cable Wh</u>	28	27	<u>Rib.Cable</u>	DRV11 J1 Pin 12 (N) Out 5 Do Hi Cal L.E.D.
		30	29	<u>Sol. Wire</u>	Cal #3 83% O ₂ Solenoid
Do Cal L.E.D. DRV11 J1 Pin 14 (R) Out 6	<u>Rib Cable</u>	32	31	<u>Gry</u>	Red L.E.D. L3 83% O ₂
		34	33		
Capacitor C5 3K FD	<u>Brn</u>	36	35	<u>Blk</u>	Grnd
		38	37	<u>Blk</u>	
		40	39		
		42	41	<u>Lt. Grn</u>	Transformer
Transformer	<u>Lt. Grn.</u>	44	43		
		46	45	<u>Red</u>	VR #1 Pin 2
+5V D.C. VR #1 Pin 1	<u>Grn #2</u>	48	47	<u>Red</u>	3K FD Capacitor C5 +5 V.D.C.
		50	49	<u>Red</u>	VR #1 Pin 2
		52	51	<u>Or</u>	VR #2 Pin 2
		54	53	<u>Or</u>	VR #2 Pin 2 +12 V.D.C.
Capacitor 3K FD C6 +12 V.D.C.	<u>Red</u>	56	55	<u>Or</u>	VR #2 Pin 1

ORIGINAL PAGE IS
OF POOR QUALITY

		58	57		
		60	59	<u>Rib.Cable</u>	DRV11 J1 Pin 16 (T) Out 7 Cal #3 On-Off
Amber L.E.D. L4 Do Cal	<u>Yel</u>	62	61	<u>Grn</u>	Green L.E.D. L5 Do Hi Cal
		64	63	<u>Sol.Wire</u>	Cal #2 100% N ₂ Solenoid
		66	65	<u>Vio</u>	Red I.E.D. L2 100% N ₂
Pin 21 (Y) DRV11 J1 Comp. Grnd	<u>Rib Cable</u>	68	67	<u>Rib.Cable</u>	DRV11 J1 Pin 17 (U) Out 3 Cal # 2 On-Off
		70	69	<u>Sol.Wire</u>	Cal #1 100% O ₂ Solenoid
Red L.E.D. L1 100% O ₂	<u>Blu</u>	72	71	<u>Rib.Cable</u>	DRV11 J1 Pin 19 (W) Out 8 Cal #1 On-Off

A/D CONVERTER

Command Status Register	"176770 ₈	
Output Buffer	"176772 ₈	
Bit 0	1 on	
(enable)	0 off	
Bits 3-4	00 x 10	
gain	01 x 5	
	10 x 2	
	11 x 1	
Bits 8-10	001 ch 2	
channel	010 3	
selection	011 4	
	100 5	
	101 6	
	110 7	
	111 8	
IADCGC (1)	"0431	Nitrogen channel enable
IADCGC (2)	"1031	Oxygen channel enable
IADCGC (3)	"1431	Carbon Dioxide channel enable
IADCGC (4)	"2021	Spirometer channel enable

DRV11 PARALLEL BOARD

Output buffer 167772

Input buffer 167774

Output control

bit 0

1

2

bit 5 4 3 2 1 0

3 cal gas 1

*

4 mass spectrometer
port 1 select

*

*

*

5 #1 lamp

6 #2 lamp

*

7 cal gas 3

*

8 cal gas 2

IVLVMD (1) = "030

IVLVMD

IVLVMD (2) = "420

IVLVMD (3) = "220

Input data

bit 0

1

2

3

4 #1 switch

switches are not used

5 #2 switch

6 = 0 100% N₂ range

7 = 1 10% N₂ range

8

LEGEND

1. 3-Way Rudolf Valve
 - 1A. One way valve from regulator (4) to face mask inhalation port (open on inspiration)
 - 1B. One way valve from face mask exhalation port (3) to spirometer (open on exhalation)

2. 3-Way Rudolf Valve
 - 2A. One way valve from 1 to spirometer (open on exhalation)
 - 2B. Sample port for mass spectrometer catheter #1
 - 2C. Adapter to demand regulator (4)

3. Sierra Firefighters Face Mask 651-02-4
 - 3A. Inhalation port
 - 3B. Exhalation port

4. 900-002-074-01 Regulator
 - 4A. High pressure O₂ supply hose to regulator
 - 4B. High pressure O₂ supply hose to O₂ source

5. Mass spectrometer catheter #1 from mass spectrometer inlet #1 to sample port (2B)


6. Oxygen tent O₂ supply hose (from O₂ source to inside O₂ positive pressure tent)

7. Air Diverter 3-Way Valve
 - 7A. Opening to face mask
 - 7B. Opening to ambient air (note: connect only between A+B or A+C)
 - 7C. Opening to spirometer intake port



8. Nitrogen gas solenoid drive connector (3/8" swagelock stainless fitting)


9. Spirometer dump port

ORIGINAL PART IS
OF POOR QUALITY




ORIGINAL PAGE IS
OF POOR QUALITY




CONTRACT NO.		 Technology Incorporated LIFE SCIENCES DIVISION HOUSTON, TEXAS
APPROVALS	DATE	
DRAWN		POWER SUPPLIES FOR N ₂ WASHOUT SYS. ELECTRONICS
CHECKED		
SIZE	CODE IDENT NO.	DRAWING NO.
A		TH8187-5E07-E000
SCALE	N/A	SHEET 7 OF 13

ORIGINAL PAGE IS
OF POOR QUALITY.



DESIGN NO.	REV.	DATE	BY	APP'D.
technology Incorporated LIFE SCIENCES DIVISION HOUSTON TEXAS				
THERMAL GAS SELECT CIRCUIT			K. V. WASHGUT	
TITLE: THERMAL GAS SELECT CIRCUIT PROJECT NO.: DRAWING NO.: REV. NO.: DATE:			SHEET NO.: 1 OF: 13	



NOTES:
 →→→
 3
 EDGE
 CONN. PIN NO.

CAUTION:
 USE ADAC
 PARALLEL
 BOARD ONLY

CONTRACT NO.	
APPROVALS	DATE
DRAWN BY LINDA WEAVER	9-11-81
CHECKED	


Technology Incorporated
 LIFE SCIENCES DIVISION HOUSTON, TEXAS


**MASS SPEC INLET CIRCUIT
 N₂ WASHOUT SYSTEM**

SIZE A	CODE IDENT NO.	DRAWING NO. TH8187-5E09-E000
SCALE N/A	SHEET 4 OF 13	



ORIGINAL PAGE IS
OF POOR QUALITY

MODIFICATIONS TO
MODEL 1100 MASS SPEC
FOR
N₂ WASHOUT SYSTEM




CONTRACT NO.		 Technology Incorporated LIFE SCIENCES DIVISION HOUSTON, TEXAS MODIFICATIONS TO MODEL 1100 MASS SPEC FOR N ₂ WASHOUT SYSTEM SIZE: A CODE IDENT NO.: DRAWING NO.: TH8187-5E08-E000 SCALE: N/A SHEET 8 OF 13
APPROVALS	DATE	
DRAWN		
CHECKED		



C-3




NOTES:
 BIT 6 USED FOR 0-100% OF J2
 0-10%
 SWIN 0-100%
 POSIT.
 BIT = 0 OR GND.

CONTRACT NO.		 Technology Incorporated LIFE SCIENCES DIVISION HOUSTON, TEXAS
APPROVALS	DATE	
DRAWN ZINDA WEAVER CHECKED	10/5/81	SPIROMETER LED INDICATORS N ₂ WASHOUT SYSTEM
SIZE	CODE IDENT NO.	DRAWING NO.
A		TH8187-5E05-E000
SCALE	N/A	SHEET 5 OF 13

ORIGINAL PAGE IS
OF POOR QUALITY




NOTES:



EDGE CONN.
PIN NO.


CAUTION:

USE ADAC PARALLEL
BOARD ONLY

CONTRACT NO.		 Technology Incorporated LIFE SCIENCES DIVISION HOUSTON, TEXAS
APPROVALS	DATE	
DRAWN	LINDA WEAVER	9-22-81
CHECKED		
		SPIROMETER CAL SWITCHES (PUSHBUTTON)
		N ₂ WASHOUT SYSTEM
SIZE	CODE IDENT NO.	DRAWING NO.
A		TH8187-5E04-E000
SCALE	N/A	SHEET 4 OF 13


ORIGINAL PAGE IS
OF POOR QUALITY

- NOTES:
1. Q1 THROUGH Q8 ARE 2N305A
 2. Q9 & Q10 ARE 2N3055 OR 2N3530P
 3. D1 & D2 ARE IN 92B .30V
 4. D3 IS IN 302020B .50V 1W
 5. O.I.1 & O.I.2 ARE .4N33
 6. S1 SHOWN IN PRIME POSITION
 7. S2 SHOWN IN NC POSITION
 8. TP READINGS MADE WITH +5V APPLIED FIRST, THEN ON FOR THE EXPIRATION CMD, S1=53, AND S2=54



TEST POINT READINGS

- TP1 (5V) +5V ON ON
 TP2 (GND) 0V ON ON
 TP4 +5V ON ON
 TP5 +5V ON ON
 TP6 +5V ON ON
 TP7 +5V ON ON
 TP8 +5V ON ON
 TP9 +5V ON ON
 TP10 +5V ON ON
 TP11 +5V ON ON



REVISIONS		DATE		AUTHOR	
DRAWING NO.		PARTS LIST		PROJECT NO.	
PROJECT NO.		CIRCUIT NO.		DESIGNER	
DRAWN BY		CHECKED BY		DATE	
TITLE		APPLICATION		REVISIONS	

technology Incorporated
LIFE SCIENCES DIVISION HOUSTON TEXAS


DUMP VALVE ELECTRONICS

THEIR PART NO. TH-1017-5E13-0003
DATE 12-13-78

ORIGINAL PAGE IS
OF POOR QUALITY

POWER SUPPLY FILTER — DRIVER BOARD

EDGE CONNECTOR



CONTRACT NO	
APPROVALS	DATE
DRAWN LINDA WEAVER	9-22-81
CHECKED	

Technology Incorporated
LIFE SCIENCES DIVISION HOUSTON, TEXAS


SPIROMETER ELECTRONICS
CHASSIS WIRING (N₂)

SIZE A	CODE IDENT NO.	DRAWING NO. TH8187-5E06-E000
SCALE N/A	SHEET 6	OF 13

4 3 2 1

REV. NO.	DESCRIPTION	DATE	APPROVED

ORIGINAL PART IS OF POOR QUALITY



- + 28 V
- NC/NO SW 52A (NC)
- NC/NO SW 52A (NO)
- BASE Q4
- COLLECTOR Q4 SW 51B
- NO/NC SW 52A WIPER
- PRIM/FAIL SW 51B WIPER
- + 10 V
- CONTROL SIGNAL
- PRIM/FAIL SW 51A (CONT SIG)
- PRIM/FAIL SW 51A (CONT SIG)
- EMITTER Q4
- PRIM/FAIL SW 51C (PRIM GRND)
- PRIM/FAIL SW 51E (FAIL GRND)
- PRIM/FAIL SW 51C WIPER
- EMITTER Q10
- NC/NO SW 52B (NO)
- NC/NO SW 52B (NC)
- Q10 BASE
- NO/NC SW 52B WIPER
- COLLECTOR Q10 i B/W 51B
- PRIM/FAIL SW 51A WIPER

QTY	CODE	PART OR	NOMENCLATURE
REQD	IDENTIFYING	NO	OR DESCRIPTION
PARTS LIST			
UNLESS OTHERWISE SPECIFIED			
DIMENSIONS ARE IN INCHES			
TOLERANCES ARE			
FRACTIONS	DECIMALS	ANGLES	
XX =		=	
XXX =		=	
MATERIAL			
FINISH			
NEXT ASSY USED ON			
APPLICATION			
CONTRACT NO			DATE
DRAWN / INCH / MEAS / CR			CHECKED
SIZE CODE BEAT NO			DRAWING NO
C			TH8187-5G12-E-000
SCALE			SHEET G. OF D

Technology Incorporated
LIFE SCIENCES DIVISION HOUSTON, TEXAS


DUMP VALVE ELECTRONICS
PRIM/FAIL SAFE (N₂)

3 2 1

4

4	3	2	1	REVISIONS	DATE	APPROVED
				ZONE	LTB	DESCRIPTION

ORIGINAL PAGE IS
OF POOR QUALITY



- INSPIRATION VALVE OUT
- INSP SPIRO POT WIPER PIN F
- BASE, Q3
- COLLECTOR, Q3
- +10 V OUT
- 200Ω, 5 W RESISTOR
- +15 V
- +28V GND
- 15 V
- 10K POT WIPER
- 15V GND
- EMITTER, Q3
- +28V OUT
- +15V TO INSP. TO TERM STRIP (POT HIGH)
- 15 V GND TO TERM STRIP (POT LOW)
- EXPIRATION VALVE OUT
- EXP SPIRO POT WIPER PIN F

QTY		CODE IDENT	PART OR IDENTIFYING NO	MANUFACTURE OPTION	PARTS LIS
UNLESS OTHERWISE SPECIFIED TOLERANCES ARE IN INCHES		CONTRACT NO			
FRACTIONS DECIMALS ANGLES		DATE			
SERIAL		DRAWING REVISION			
MATERIAL		CHECKED			
FINISH		SIZE CODE IDENT NO (DRAWING NO)			
NEXT ASSY USED ON		SCALE N/A			
APPLICATION		SHEET N. OF 13			

Technology Incorporated
LIFE SCIENCES DIVISION HOUSTON, TEXAS

DUMP VALVE ELECTRONICS
POWER SUPPLY AND
VOLUME READ OUT (N₂)

THE 167-5611-E-000

9.0 APPENDIX: RT-11 COMMANDS

9.1 SPECIAL FUNCTION KEYS

The special functions of certain terminal keys you need for communication with the keyboard monitor are explained below. A more complete list of special function keys and commands is provided in section 3 of the RT-11 System User's Guide.

Enter CTRL commands by holding the CTRL key down while typing the appropriate letter.

CTRL/C terminates program execution and returns control to the keyboard monitor. CTRL/C echoes C on the terminal. You must type two CTRL/C's to terminate execution unless the program to be terminated is waiting for terminal input or is using the TT handler for input. In these cases, one CTRL/C terminates execution.

CTRL/O causes RT-11 to suppress terminal output while continuing program execution. CTRL/O echoes O on the terminal. RT-11 enables terminal output when one of the following occurs:

1. You type a second CTRL/O.
2. You return the control to the monitor by typing CTRL/C when program terminates.

Note that when you are using CTRL/O the system can print an extraneous character after the monitor echoes the CTRL/O and a carriage return/line feed.

CTRL/Q resumes printing characters on the terminal from the point previously stopped because of a CTRL/S.

CTRL/S temporarily suspends output to the terminal until you type a CTRL/Q. CTRL/S does not echo. You can alternate between CTRL/S and CTRL/Q to temporarily stop and then resume the display of output on a CRT terminal. Type CTRL/S when the screen is almost full, then CTRL/Q to continue when you are ready.

DELETE OR RUBOUT deletes the last character from the current line and echoes a backlash plus the character deleted. Each RUBOUT succeeding DELETE deletes and echoes another character. The system prints an enclosing backlash when you type a key other than DELETE. This erasure is performed from right to left up to the beginning of the current line.

If you are using a video display terminal and you have issued the SET TT SCOPE command, DELETE erases the characters with a backspace, space, backspace sequence. Your corrections appear on the screen; RUBOUT does not enclose them with backlash characters.

9.2 USING THE SYSTEM UTILITIES FOR FILE MAINTENANCE

Keyboard commands let you communicate with the RT-11 monitor to allocate system resources, start programs, and use various services.

The following keyboard commands are the ones most frequently used with our LSI-11 computers.

They should be used regularly to do such things as obtaining listings of the directory information for the data floppies and making backup copies of files. For more detailed information about using these or other commands, refer to the DEC RT-11 System User's Guide.

Some of the more commonly used RT-11 commands for performing file maintenance tasks are listed below. They may be used in any order, whenever the system has typed a "." to let you know it is ready.

All information that you enter is underlined here. Anything which is not underlined is typed by the computer.

A. Initializing a Floppy Disk

This is done either:

1. When you want to use a new floppy disk which has never been used before, or
2. When a floppy has old data on it which is of no value, and you want to remove the old information so that the floppy can be used again.

To initialize the floppy, put it in the right hand disk slot and type:

DIR DY1:

The system will then type: DIR-F-Illegal Directory (if the floppy is new and has never been initialized before), or it will type a list of the files already stored in the floppy. Make sure that there are none that you might want to save.

Next type:

INIT/BAD/VOL DY1:

When the system asks "Are you sure?", check to make sure that you have the correct floppy in DY1 and that you have not erroneously typed DY0 instead of DY1. The type Y (for Yes) followed by a carriage return. If you are not sure, type N instead, and repeat the INIT command.

The /BAD option with the INIT command will cause the floppy to be checked for any bad blocks that might be present on it. This will take a minute or so. When finished, the system will type ?DUP-I-No bad blocks detected DY1: if the floppy has no bad spots, it will then ask for a volume (or floppy) identification

code, followed by the owner name. You may type in anything you want for these two entries, but each must be 12 characters or less in length, including spaces, hyphens, etc. Typically, you would use a floppy ID number or code for the Volume ID, and the laboratory name (or your own name) for the owner. For example:

Volume ID? Disk 52A

Owner? G Washington

The INIT command causes the directory on the floppy in DY1 to be zeroed out. This effectively erases the information on it. Although it is sometimes possible to retrieve information from a floppy after it has had its directory zeroed (if nothing on it has been overwritten), it is rather difficult and generally not practical.

B. Getting a Floppy Directory (Contents) Listing

To get a complete, formatted summary of what is stored on a floppy disk, type the following:

DIR/VOL/FULL/BLOCKS DY1: (Carriage Return)

This causes a listing of the file names, their creation dates, and other information concerning such things as their physical locations on the floppy in DY1 to be printed on the terminal. The date printed for each file is the date the file was created or last modified, whichever is later.

To get a simplified directory listing, type:

DIR DY1: (Return).

A listing of each subject data floppy, using the DIR command, should be obtained regularly if data files are being created on it. This should preferably be done at the end of each day in which tests were run. The latest listing for each floppy should be saved in order to help keep track of what information is on each floppy, and also because it can sometimes be useful in retrieving files which have accidentally been deleted or overwritten. Each such listing should be labeled to correspond with the particular floppy it summarizes, and then taped on the special protective envelope for that floppy. Alternatively, the listings can be kept in a notebook.

If you need to see if a particular file is on a floppy without getting a complete listing, put the floppy in DY1 and type:

DIR DY1:name.suf (Return)

Here name is the filename that you used when creating the file, and "suf" is the 3-letter suffix (or extension or file type) for that file. If you are not sure what the suffix is, just type DIR.DY1: name with no period following the name. If the file

exists on the floppy, the system will print a line containing the directory information for that file, preceded by the date if it has been entered into the system, and followed by the total number of files listed (not the total number on the floppy) and the number of free (unused) blocks on it. The number of free blocks can be used to estimate how many more subject data files can be stored on the floppy. If the requested filename does not exist on the floppy, then only the date, "0 files", and number of free blocks will be printed. Therefore, this is also an easy way of determining how much unused space is on the floppy. Note that the unused space may not all be in one contiguous area, however, the /FULL option can be used to tell you where the unused blocks are distributed.

C. Deleting Files

To delete a single file which has erroneous or useless information on it, put the floppy in DY1 (the right hand disk slot), and type:

```
DELETE DY1:name.suf followed by a carriage return,
```

where "name" is the filename that you want to delete. This can be repeated for several files, if desired. Be sure that the filename is entered correctly, because a different file could be inadvertently deleted if the incorrect filename is typed in. The system will type the file name followed by "?".

Respond with Y (Return) if you are sure that you want to delete a file.

D. Renaming a File

If a test has been run and it is later discovered that an incorrect filename was used throughout, then it can be renamed, without changing any of the data stored in it, by inserting the floppy into DY1 and typing:

```
RENAME DY1:oldname.suf DY1:newname.suf
```

where "newname" is the new filename under which you want the data to be stored, and "oldname" is the old filename that was previously used.

For example:

```
RENAME DY1:T00900.DAT DY1:T01000.DAT
```

will cause the filename for file T00900.DAT to be changed to T01000.DAT and the old filename will no longer exist.

E. Copying Data Files (Making Disk Backups)

Data files may be copied from one floppy to another very easily. This enables one or more backup copies to be made of all important data.

To copy all files from one floppy disk to another, do the following:

1. Initialize the blank floppy you are going to copy onto as described above in "A. Initializing a Floppy Disk". Note that this will destroy any information stored on it, so be sure you have used a floppy which does not have any data that you want to save.

This step can be skipped if you are only going to add or replace files on a floppy which already has files on it.

2. Several options are now available depending on what you want to do.

- (a) To copy all files from DY0 to DY1, type:

```
COPY/WAIT/SYS DY0:*. * DY1:*. *
```

(You do not need to include /SYS unless you are copying system files from a system floppy.)

The system will then let you take out the system floppy (if it is not the one you are copying from) and replace it with whatever floppy has the files you wish to copy.

(The /WAIT option is not necessary if you are copying from your system floppy.) This process can take from several seconds to several minutes, depending on how many files are being copied. When the system has finished, it will tell you to replace (mount) the original system floppy back in DY0: (if you used the /WAIT option), and then let you type Y when you are ready. The system will type a period when it is ready for a new command.

- (b) To copy only those files with the current date (which were created or last modified today), include /NEWFILES in the list of options with the copy command.

```
COPY/NEWFILES/WAIT DY0:*. * DY1:*. *
```

The /NEWFILES causes the system to compare the dates of each file in the disk directory with the current date entered previously with the DATE command. Therefore, the current date must have been entered correctly when the computer was powered up in order for this to work properly.

- (c) To selectively copy certain files, include /Q (for QUERY) in the list of options. For example:

COPY/SYS/Q DY0:*. * DY1:*. *

The /Q causes the system to individually list each filename on DY0 to see if it is to be copied. To bypass copying a particular file, simply hit the Carriage Return Key after the filename is printed. For files that you do want to copy, type Y(Return) (Y for yes) after their names. For example: if DY0 has files named A.DAT, B.DAT, and C.DAT on it and you want to copy only A.DAT and C.DAT do the following (your responses are underlined as usual):

Type:

COPY/Q DY0:*. * DY1:*. *

A .DAT? Y(Return) (yes, copy it)

B .DAT? (Return) (no, do not copy it)

C .DAT? Y(Return) (yes, copy it)

(d) To copy one or a few files from DY0 to DY1, type:

COPY/SYS/WAIT DY0:name.suf DY1:*. *

Where "name" is the filename. This can be repeated for several files, if desired.

3. A backup system floppy can be created in a similar manner. First, put a blank floppy in DY1 and perform step (1) above. Next, type:

COPY/SYS DY0:*. * DY1:*. * (A list of files copied will follow.)

COPY/BOOT DY1:RT11SJ DY1:

DY1 now contains a floppy which has all of the files that the original system floppy had on it. Note that the copying caused by the first line above may take a minute or so to perform, so wait until it is finished and the system responds with another "." before typing the second line.

4. These procedures can be used regularly to maintain backup copies of all data files. The recommended procedure is:
 - (a) First, perform step (1) above.
 - (b) Second, perform step (2) (a) to copy all the files onto a new backup floppy.
 - (c) Third, perform step (2) (b) at the end of each day to copy new files onto the backup floppy. As an alternative, the new files from one or more days may be copied using the procedures in steps (2) (c) or (2) (d).

It is probably best to keep two backup floppies, in addition to the original, for important files, since it takes only a few minutes to perform the copying operation.

If you try to copy a data file from one floppy onto another that already contains a file with the same name, then the file on the floppy being copied into will be deleted and the copy operation will be performed as requested. Normally, this is not really a problem, since files with the same name should generally be identical. However, if an erroneous filename was entered when creating a file, for example, then it is possible that this file might replace one on the backup floppy unintentionally. This could happen if the file was copied before correcting the filename using the Rename command (step D above). Therefore, be particularly careful when entering filenames or copying with possibly erroneous filenames, since valuable data might be deleted. Frequently the /NOREPLACE option with the COPY command will help to avoid this potential problem.

F. Getting Help with System Commands

If you do not remember which command or option to use for a particular function, or the syntax for a particular command, use the HELP command. If you type only HELP followed by a carriage return, information on the help command itself will be typed. For a list of all the possible commands (some are not always available, depending on which system disk you are presently using), type HELP (Return). For more information on one of these commands, including possible options, type HELP (space) command, where "command" is the name of the command (or the first few letters of it). For example, for information on the DIRECTORY (or DIR) command, type:

HELP DIR followed by a carriage return.

If the information typed is not clear or you need more details, refer to the RT-11 System User's Guide. The first half of that manual contains an alphabetized list of all commands, with complete information and examples on most options.

9.3 GLOSSARY

Breath gate - the minimum noise level of the spirometer. Also referred to as minimum noise gate.

DATAUN.VAL - a data file on the system disk which contains the default values for calibrating the spirometer and mass spectrometer.

DECDAT - a program on the system disk which accepts information from the test operator on a decompression experience experiment, types a report on the printing terminal, then stores the data in a data file.

DY1: - the disk drive device on the right where a floppy disk can be inserted, e.g., for a data file.

External gravity - the gravity outside of the experimental environment, one if the experiment is conducted at ground level and zero if conducted in space.

N2DATA - a program on the system disk which reads the nitrogen washout data from a data file and types a report on the printing terminal.

N2WASH - the main program on the system disk which calibrates the spirometer and spectrometer, records the data during a nitrogen washout experiment and stores it in a data file.

PDINPT - a program on the system disk which accepts information from the test operator on the physical characteristics of a subject, types a report on the printing terminal and stores the data in a data file.

UPDATE - a program on a system disk which reads an unformatted nitrogen washout data file, formats the data, then stores it in a data file.