

INITIAL BLOOD STORAGE EXPERIMENT

.

NASA 9-17222

10 21 - 5 - 4 - 3

(NARA-CR-185513-Pt-2) INITIAL BLOOP STORAGE N91-14701 EXPERIMENT (Fenter for cloop Research) 301 p CSCL 06C Unclus

63/51 0235683

The Center for Blood Research

Contractor

Subcontractors Arthur D. Little, Inc. The Children's Hospital Lahey Clinic Medical Center University of Massachusetts Medical Center

NAS 9-17222 D.M. SURGENOR, Ph.D.

INITIAL BLOOD STORAGE EXPERIMENT

TABLE OF CONTENTS

- 1. SPACE SHUTTLE BLOOD STORAGE EXPERIMENT ADDENDUM
- 2. SCIENTIFIC EXPERIMENTAL PROTOCOL
- 3. SCHEDULES FOR BLOOD COLLECTION, COMPONENT PREPARATIONS AND DISTRIBUTION
- 4. PRE-LAUNCH CLOSEOUT PROCEDURES (REVISION 1.0)
- 5. OPERATING, MAINTENANCE AND HANDLING PROCEDURES
- 6. DD 1149 FORM WITH EQUIPMENT AND SUPPLY LIST ATTACHED
- 7. SUMMARY OF MAIN INITIAL BLOOD STORAGE EXPERIMENT ACTIVITY
- 8. POWER DISSIPATION OF THE INITIAL BLOOD STORAGE EXPERIMENT
- 9. DATA ANALYSIS

The Center for Blood Research

SFACE SHUTTLE BLOOD STORAGE EXPERIMENT

ADDENDUM

The Center for Blood Research Children's Hospital Lahey Clinic U. of Mass. Medical Center Arthur D. Little, Inc.

7.0 ADDENDUM

7.1 Introduction

This addendum provides further justification for the use of microgravity for the Blood Storage Experiments. The addendum is comprised of five portions: 1) significance, 2) working hypothesis, 3) storage lesions at earth's gravity, 4) biomaterials and 5) summary.

7.2. Significance

The possibility of conducting experiments with the formed elements of the blood under conditions of microgravity opens up important opportunities to improve our understanding of basic formed element physiology as well as contribute to improved preservation of the formed elements for use in transfusion. In attempting to identify the possible effects of reduced gravitational attraction on formed elements during preservation, it is essential to recognize that 1) the formed elements - red cells, white cells and platelets - are complex living systems which depend upon multiple interactions with their extracellular environment, of which gravitational force is but one component; and 2) the surivial within the circulation and the functional capacity of stored formed elements after transfusion are markedly diminished when compared to cells from fresh blood.

Multiple factors affect blood cells during storage such as sedimentation, diffusion of nutrients and metabolites, gas transport, mechanical trauma due to mixing and leaching of chemical compounds from the container material. Storage under standard blood bank conditions invariably exposes blood to Earth's gravity, and leads to injuries, i.e. storage lesions (see below) to the stored formed elements.

Exposure of the blood at microgravity in the space shuttle is expected to result in a greatly reduced perturbation of cell functions and permit the separation of sedimentation, mixing and container lesions which are associated with terrestrial blood storage.

Information gained from storage experiments under microgravity as proposed will add a new and important dimension to our understanding of blood cell membrane functions and the storage defects, and serve as a guideline for design of the improved methodology for blood storage on Earth. Further, the knowledge gained can serve as a guide not only for the polymer industry to develop alternative blood bags but also for the pharmaceutical industry to study drug effects on cell membranes.

An Experimental Hypothesis

The objective of this proposal is a comparative evaluation of the physiological, biochemical and physical changes of the membrane of the erythrocyte, platelet, and leukocyte during storage under two specific conditions: standard blood bank conditions and microgravity, utilizing three plastic formulations.

Given the complexity of living biological systems, the wide range of possible measurements which could be made and the limited set of measurements which have been chosen in the first set of experiments, the following are some of the results we expect from our studies:

Studies comparing the state of human blood following prolonged preservation at a) Earth's gravity, and b) the microgravity within the shuttle, should shed important new light on the fundamental cell physiology of each type of blood cells, and, if we are fortunate in our choice of measurements, should provide significant new clues to improved preservation of these important materials on earth.

For cells preserved at microgravity, the adverse effects of sedimentation and mixing should be less apparent against the background of container lesions and aging lesions.

Red cells preserved at microgravity should exhibit improved pH control and more nearly physiological levels of important metabolites such as ATP and 2,3 DPG.

Platelets and leukocytes which are more sensitive than red cells to harmful effects from cell-cell interactions, should show superior indices of cell viability and potential to function after storage at microgravity in comparison to controls kept on Earth.

Adverse effects of synthetic polymers which originate from cell-surface interactions should be reduced during microgravity preservation.

ORICHAL RACE IS

OF PORY QUALITY

7.3

7.4 Storage Lesions at Earth's Gravity

In this portion, the discussion is divided into four sections: 1) storage lesions overview, 2) red cell storage on Earth, 3) platelet storage on Earth, and 4) leukocyte storage on Earth.

7.4.1 Storage Lesion Overview

In discussing the science of blood formed element preservation, it has proved useful to classify observed results into several categories of effects, termed lesions. It must be noted, however, that this does not necessarily mean that we can identify a specific area of physical damage or injury with each lesion.

7.4.1.1 Sedimentation Lesions

The sedimentation lesion results from the dense packing of the formed elements, with resulting intensity of injurious cell-cell and cell-container surface interactions. Under these circumstances, metabolic injury can occur as cells are deprived of their substrates or as cell wastes concentrate in the immediate vicinity rather than being dissipated and removed as in the physiological state.

Within the blood stream, erythrocytes, leukocytes and platelets, like all cells, are metabolically active. They are supplied with energy in the form of small organic molecules (substrates), which are contained in the fluid medium, the plasma, which surrounds them, and they depend upon the plasma to accept the waste products of their metabolism. If these waste products, particularly carbon dioxide, are not removed, they can quickly alter the extracellular environment and injure the cell.

The blood cells utilize energy to carry out their physiological functions, and to maintain their own integrity. Thus, the red cell requires a constant supply of energy derived primarily from the metabolism of glucose to maintain its unique biconcave discoid shape. If the metabolic system fails, the cells invariably die.

The rationale of preservation of the red blood cell, the formed element about which we know the most, takes advantage of what is known about red cell metabolism. A preservative diluent solution provides glucose, the main energy source, and other chemical intermediates needed

-3-

to foster a viable metabolic system. Another constituent of the diluent lowers the pH of the red cell suspension. And finally, and most important, the temperature is reduced from 37.5°C to between 1 and 6°C. Under these conditions, the metabolic system is put into a quasi-hibernative state, in which it survives remarkably well. Glucose utilization proceeds at a very slow rate, as does carbon dioxide production; and no oxygen is required.

It follows that an experiment at microgravity will shed light on cellular metabolism and other parameters of cellular viability which are most sensitive to the harmful effects of gravitational force.

7.4.1.2 Aging Lesions

Aging lesions arise from the simple fact of biologic aging which all biological cells are subject to. Red Cells, for example, normally live for 120 days in the circulating blood. Thus, if preserved outside the body, red cells may continue to age. In the blood bank, at $1-6^{\circ}C$, this aging contributes to the time limit for preservation. The specific cause of biologic aging is not known but is thought to include membrane changes of cells.

The ability of blood cells (red cells, platelets and leukocytes) to maintain their size and shape against the stress of gravity requires the expenditure of energy. Reduction or elimination of this stress in microgravity will lead to marked energy conservation. As a result, the characteristics of the aging processes of blood cells can be better defined.

7.4.1.3 Container Lesions

Container lesions are defined as those effects brought about by exposure of the blood or its formed elements to foreign surfaces such as glass tubes, metal needles and plastic bags. These include but are not limited to cell-foreign surface effects and cell-foreign surface constituent effects. A whole field of endeavor is concerned with cell-foreign surface effects; these being particularly important in respect to implanted prosthetic devices, such as artificial heart valves, prosthetic aortas, etc. In these cases, biocompatibility is probably determined by physical forces at the interface between the solid and liquid phases of the blood forces which may adsorb and rupture formed elements, or which may damage the complex composition of lipids and proteins in the cell membrane. An example of a cell-foreign surface constituent effect is the leaching of the organic chemical plasticizer, DEHP, from the plastic blood bag onto the red cell surface.

7.4.1.4 Mixing Lesions

Mixing lesions represent the fourth type of storage lesions which are encountered in blood preservation. These result from the use of mixing procedures to improve the preservation of separated formed elements. In a simplistic view, mixing is usually used to prevent the deleterious effects of sedimentation.

Mixing ostensibly reduces the deleterious effects of "crowding" of cells, results in the accessibility of essential nutrients to the cell, and eliminates the undesirable effects on the cell of accumulated waste products.

But mixing unfortunately can introduce undesired effects as well. Mixing subjects cells to energetic forces of turbulence, cell-container surface and cell-cell interactions. Not surprisingly, the mixing can produce variable effects which are often brought to light when mixing is carried out in different polymer bags and modes of movement (see Table I in section 7.4.3.1.).

Microgravity which maintains the cells in suspension, provides a unique way of avoiding both the sedimentation and the mixing lesions.

- 7.4.2 Red Cell Storage on Earth
- 7.4.2.1 Sedimentation Effects on Red Cell Storage

Storage of whole blood under standard blood bank conditions leads to sedimentation of the cells which produces storage lesions because metabolites accumulate around the cells, pH decreases, and substrates are depleted. These lesions develop as a result of the limitation of diffusion which in turn is determined by the thickness of sedimented cell layers.

Recently, Beutler et al. (1) compared 24-hour posttransfusion survival and hemolysis of red cell concentrates collected into CPD-A2 (anticoagulant with increased glucose) and stored for 42 days either standing vertically or lying horizontally. It was reasoned that the horizontal position would facilitate diffusion by increasing the plasma to cell surface ratio. They found that hemolysis was significantly greater and survival significantly lower in the concentrates stored in the vertical position than those stored in the horizontal position.

The significant differences in hemolysis may be due to the inability of substrates to diffuse from the overlying plasma into the deeper layers of the sedimented cells in the vertical position during the storage period. However, the percent of hemolysis in the units stored in the horizontal sedimented position was still unacceptably high (>12).

7.4.2.2 Mixing Injury of Red Cells During Storage

Sedimentation can be prevented on Earth by continuous mixing. Dern et al. (2) studied the effect of intermittent mixing (5 times weekly), single mixing (once midway through storage) and non-mixing on the storage of CPD whole blood for 28 days. They demonstrated that intermittent mixing gave the best in vivo and in vitro results after storage, followed by single mixing and then non-mixing. Units with intermittent mixing showed an increase in ATP level and in post-infusion survival but a decrease in plasma hemoglobin when compared to the single- or non-mixing units.

However, injury caused by agitation (or mixing) to red cells was described as early as 1884 by Meltzer et al. (3). These authors observed that agitation of whole blood in glass tubes produced the release of hemoglobin from red cells and demonstrated that this was dependent on shear stress and on interactions between red cells and container wall surfaces.

Although the mechanism of the mixing lesion is poorly defined, it has been shown that mechanical mixing can result in loss of membrane lipid (4), reduction in the surface to volume ratio (5), increase in cation permeability (6) and decrease in deformability which correlates with poor post-infusion survival (7).

7.4.3.1 Platelet Storage on Earth: Mixing Injury

Storage of platelet concentrates under standard blood bank conditions requires continuous mixing. Work from several laboratories has shown that mixing of platelets during storage prevented extensive deterioration of platelet functions as measured by activation, yield and survival (8-11). for example, Murphy et al. (8) studied the effect of mixing on the survival of 48 hour stored platelet concentrates. The mean survival $(T_{1/2})$ was 3.6 days for the mixed, 2.6 days for the non-mixed concentrates, and 4.2 days for fresh platelets. Morphology of platelets and the pH of the concentrate were also better maintained in the mixed units (9).

Several commercial mixing devices are currently available for storage of platelet concentrates. There are two general types, flatbed (horizontal) and rotary. Results from storage of platelet concentrates using two different five-day bags and four different mixing devices are summarized in Table I. As illustrated in Table I, the mode of agitation is extremely important for both in vivo and in vitro function of platelets. At the present time, the mechanism as to why only minor variations in the mode of mixing should have such significant effect on platelet storage has not been defined. It should be emphasized, however, that even with continuous mixing, the recovery, survival and function of stored platelets are decreased when compared to fresh platelets.

7.4.3.2 Contact Activation of Platelets

It is well established that platelet-platelet contact as well as platelet-container surface contact can cause alteration and activation of these cells. In other words, contact among platelets themselves and between platelets and foreign surfaces can induce platelet injury. Continuous mixing during storage facilitates contact activation of platelets. This may be one of the mechanisms which leads to decreased recovery, survival and function of stored platelets.

7.4.4 Leukocyte Storage on Earth

To date, it is necessary to transfuse leukocytes as soon as possible after collection (less than 6 hours) because these cells rapidly lose their integrity and function during storage. There are many unsolved problems with regard to leukocyte storage (18-21). For example, leukocyte collection under standard blood bank conditions invariably contains a large number of platelets and red cells. Interactions among these cells during sedimentation result in a fall in pH lack of diffusion of nutrients and clumping of leukocytes by platelet activation.

Little is known about the effect of mixing on leukocytes during storage. McCullough et al. (18) reported that mixing during storage, although it had no effect on the bacterial killing function of leukocytes, caused a Table I. Effects of Methods of Agitation on in Vivo and in Vitro Characteristics of Platelets Stored for 5 Days at 22°C

					Polymers			
Method of		PL 732	ja			crxp		
Agitation	Function	Survival	Recovery	Ref.	Function	Survival	Recovery	Ref.
HORIZONTAL								
Flathed	++C	* * *	+ + +	12,13	+++	+ + +	+ + +	14-i6
ROTARY								
Eliptical Turbler	+1	÷1	÷1 .	13	*	:	•	14,15
rumuter Ferrig Wheel	• • • •	+ + + +I	+ + + +1	13 12,13	+ + + +	+ + +1 + +1	+ + +	14-16 17
^a polyolefic plas c +++: acceptable	tic; ^D PVC ; ++: fair	plasticiz ; +: equiv	ed with ocal; ::	TOTM; unaccept	able			

ORIGINAL PAGE IS OF POOR QUALITY 1

, ,

÷

reduction in chemotaxis (migration of leukocytes to a noxious stimulus such as bacteriae). The chemotactic response of leukocytes, an energy-dependent process, deteriorates first and to the greatest extent during storage when compared to other functions (19-21).

We have previously proposed that storage of blood at microgravity leads to marked energy conservation (see section 7.4.1.2.). Therefore, microgravity will better maintain energy-dependent functions of leukocytes, e.g. chemotaxis. Leukoctye function is maintained through a narrow pH range (7.0 to 7.5) (18). Storage at microgravity which prevents sedimentation would reduce the steep pH gradients to which leukocytes are subjected under standard blood bank conditions may in fact improve the viability of leukocytes.

7.5.1 Biomaterials and Blood Storage on Earth

A phthalate plasticizer is currently employed in the manufacture of blood bags. The evidence of the harmful consequence of the phthalate plasticizer has been reviewed in our original proposal (section 2.1.4, Biomaterials). Accordingly, industry has been engaged in extensive studies to develop alternative polymers.

In comparison to the standard PVC-DEHP blood bags, all alternative polymers developed to date have been deleterious to the red cell in that: 1) there is increased hemolysis; 2) the red cell potassium does not leak out to the same degree as the hemoglobin; 3) there is increased osmotic fragility; 4) the red cell survival at 21 days in CPD anticoagulant is reduced and 5) the red cell survival in all new polymer formulations utilizing CPD adenine anticoagulant is less than 70% (FDA Regulations require greater than 70% survival) after 35 days storage. These factors point to a selective membrane alteration. All attempts to define the mechanism of this alteration on earth have failed.

7.5.2 Biomaterials and Microgravity Blood Storage

The interaction of biomaterials and blood cells during storage is one of the multiple variables that can affect blood storage on earth (see above). Microgravity as provided by the space shuttle would eliminate the effect of two most important variables, i.e. sedimentation and the mechanical mixing injury to stored blood cells (see previous section), and is the only means by which the separation of the sedimentation, container and mechanical mixing lesions can be achieved. A comparison of the membrane changes of blood cells stored in different

OF POOR QUALITY

polymer formulations at microgravity will give new insight into the membrane characteristics of blood cells upon storage.

In this project, we will utilize three FDA approved blood bags, one, the currently used PVC-DEHP bag, and two alternative polymer formulations, the PVC-TOTM and polyolefin bags. PVC-DEHP and PVC-TOTM contain the plastizers phthalate and trimellitate, respectively, whereas, polyolefin bags do not contain a plasticizer. The studies summarized below have been carried out in the Blood Bank Research Laboratory at Children's Hopspital (under the leadership of Drs. S. Kevy and M. Jacobson) as part of a cooperative program with industry for the in vitro and in vivo evaluation of alternative polymer formulations. Our results demonstrate the existence of a membrane abnormality due to the interaction of the polymer with the membrane of both the red cell and platelet as well as the affinity of the plasticizer for lipid.

7.5.3 Leaching of Plasticizer

The leaching of DEHP from "VC during blood storage is dependent on diffusion and is directly proportional to the lipoprotein concentration. Two aspects of diffusion are involved: 1) migration through the plastic to the internal surface of the container, and 2) diffusion of lipoprotein molecules toward and away from the container wall. Although only a fraction of the leached DEHP is associated with the red cell, this plasticizer has a membrane stabilizing effect (see below). In a previous study, we compared the degree of leaching of TOTM and DEHP and their effect on human fibroblast tissue culture. The leaching of TOTM by plasma is onehundredth that of DEHP and is independent (unlike DEHP) of lipoprotein concentration. Furthermore, results from tissue culture studies showed that TOTM did not inhibit the growth of human fibroblasts, whereas DEHP was a potent inhibitor (22).

7.5.4 Lipid Association of Plasticizer

During storage there is an active exchange between the plasma and erythrocyte membrane lipids (23). It is thought that the neutral lipids other than cholesterol are not an integral part of the red cell membrane but represent small amounts of plasma lipoprotein material adsorbed to the surface of the membrane (24). Utilizing polyvinyl chloride plasticized with ¹⁴C-carbonyl labeled DEHP we were able to determine that the degree of leaching was limited by the lipoprotein concentration. Ultracentrifugation studies revealed that 93% of the DEHP found in plasma was in the lipoprotein fraction. As shown in Table II, agarose gel chromatography demonstrated that the low density component bound more plasticizer than either the very low or high density components (25). There was little bound to albumin. The plasticizer is probably solubilized and transported in plasma in a manner similar to that of the triglycerides. ٢

Table II. Association of ¹⁴C-DEHP with Lipoprotein Components of Plasma

	2 0	f Total	Recovery	
DEHP Addition	VLDL	LDL	HDL	Vi*
23.3 mg/ml	19	47	27	7

*VLDL: very low density lipoprotein; LDL: low density lipoprotein; HDL: high density lipoprotein; Vi: non-lipoprotein fraction

Membrane Stabilization Effect of DEHP

Red blood cells when exposed to DEHP are more resistent to hemolysis as compared to those cells which are stored in non-DEHP containers (26). The mechanism of the "protective" effect of DEHP on hemolysis is unknown. However, Seeman (27) demonstrated that there is a 50% inhibition of hemolysis when the membrane area of the intact erythrocyte expands by 2 to 37. Possible mechanisms by which compounds could expand the membrane area are: 1) occupation of space within the membrane, 2) disordering and expansion of the lipid regions of the membrane, 3) induction of conformational changes in membrane proteins, and 4) adsorption to the membrane, thus altering the amount of water in contact with the membrane. Furthermore, Roth and Seeman (28) demonstrated that lipophilic anesthetics stabilize the red cell membrane against hypotonic hemolysis while drugs with low lipid solubility had no effect. It is conceivable, therefore, that DEHP, due to its lipophilic nature may stabilize the red cell membrane in a like manner.

7.5.6 Red Cell Studies

Storage of blood in alternative polymer containers (PVC-TOTM, the polyolefins, and etc.) resulted in increased hemolysis and osmotic fragility and decreased post-infusion survival as compared to PVC-DEHP bags. In

7.5.5

this section additional data from studies of hemolysis and osmotic fragility are included whereas the survival data have been reviewed in our original proposal.

7.5.6.1 Hemolysis Studies

The rate of hemolysis in stored blood is generally accepted as the primary initial screening procedure to indicate biocompatibility with the erythrocyte. In 1954, Gibson and Thorn (29) compared blood stored in glass and plasticized polyvinyl chloride containers, and demonstrated that the plasma hemoglobin level and the susceptibility to hemolysis in hypetonic 0.67 saline was evident at one week and markedly increased during the third week of storage in the "inert" glass container as compared to polyvinyl chloride formulations, whereas all other biochemical parameters measured were identical.

Paired studies were carried out in our laboratory to compare plasma hemoglobin, osmotic fragility, red cell electrolytes and glycolytic intermediates in blood stored in PVC-DEHP bags and alternative polymers. The average plasma hemoglobin level at one week was markedly different for PVC-DEHP and polyolefin but similar for the polyolefin and a polyolefin copolymer. By 3 weeks the plasma hemoglobin level for the polyolefin bag was 64.2 mg/dl with its PVC-DEHP control at 34.1 mg/dl; the polyolefin co-polymer had a level of 44.0 mg/dl with its PVC control at 27.5 mg/dl (means of 5 paired studied). Similar results were obtained with 5 other polyolefin derivatives. All exhibited significant differences by one week.

Studies were also performed comparing PVC-DEHP with an identical PVC formulation using a TOTM plasticizer. The average plasma hemoglobin for the PVC-DEHP controls was 33.4 mg/dl and the PVC-TOTM was 66.5 mg/dl. The rate of hemolysis of the PVC-DEHP and PVC-TOTM packs were quite different. This was evident within the first week of storage.

7.5.6.2 Osmotic Fragility Studies

Results from paired studies on osmotic fragility are shown in Figs. 1 and 2. We demonstrated that significant losses in stability as measured by hypotonic hemolysis occurred in red cells stored in containers fabricated with non-plasticized polymers and TOTM as compared to the PVC-DEHP containers (Figs. 1 and 2). Fragility curves were generated for each sample but only the percent hemolysis in 0.6% NaCl is shown. By the end of four weeks the percent hemolysis is 11 for PVC-DEHP and 27.5 and 29.5, respectively for two of the experimental polymers (Fig. 1). After 3 weeks of storage the PVC-DEHP units had 12.5% hemolysis and after four weeks 18.5% hemolysis whereas the PVC-TOTM units had 23.5 and 34.8 percent hemolysis at the corresponding storage times (Fig. 2). It should be emphasized that as early as one week after storage, significant differences were noted.

In the same paired studies, red cell electrolytes and glycolytic intermediate levels were identical for all polymers tested. All the aforementioned polymers have excellent gas transfer and two of them have been approved for 5 day platelet storage.

7.5.7 Platelet Studies

Therapeutic effectiveness of platelets is judged by three parameters: 1) Yield or Recovery: This is the term applied to the percent of the platelets that circulate following transfusion. 2) Survival: The duration in days that the platelets circulate. 3) Correction of the patient's bleeding time.

Previous studies from our laboratory has demonstrated that varying formulations of PVC bags plasticized with DENP exert a profound effect on platelet yield and survival during storage (Table III). As shown in Table III, polymers #2 and #3 were suitable for storage up to 72 hours with normal survival whereas polymers #1 and #4 had markedly decreased survival within 48 hours. Both polymers 2 and 3 are currently used in blood banks for platelet storage up to 72 hours. It should be emphasized, however, that platelet recovery in the two latter polymers was only 40% of that of fresh platelets.

Plastic	Storage		Recovery	Survival
Pack	Temp (*C)	Hrs.	(Z)	Days
#1	22	0	51.2 + 1.5	8.6
#2	22	0	61.3 + 2.2	8.5
#3	22	0	58.6 ∓ 2.7	8.3
#4	22	0	42.5 🛨 3.1	7.9
#1	22	48	43.5 + 2.9	2.9 ⁸
#2	22	48	52.5 + 2.2	7.5
#3	22	48	54.8 + 2.1	6.9
#4	22	48	41.1 + 2.6	1.8 ^a
#2	22	72	35.0 + 1.6	7.2
#3	22	72	35.2 <u>+</u> 1.4	7.1

Table III. Effect of Four PVC-DEHP Formulationson Platelet Recovery and Survival
(n=10)

^aThe 72-hour recovery and survival were not determined due to poor survival at 48 hours.

7.6

Summary

Clearly, from the above discussion storage of whole blood, platelet concentrates and leukocytes under earth's gravity, either in the sedimented state or with continuous mixing is not optimal. Storage of blood at microgravity in the space shuttle is expected to result in a greatly reduced perturbation of cell function. The continuous suspension of blood during storage in the shuttle at a microgravity environment will more closely mimic the relationship of the cells to the plasma in the circulatory system than does ground storage of blood. Continuous suspension of blood will 1) eliminate the accumulation of metabolites immediately surrounding blood cells, 2) increase the availability of substrates, 3) increase gas transport, 4) better maintain pH, and 5) reduce or eliminate the chances of contact activation for platelets.

REFERENCES

- 1. Beutler, E. and C. West: J. Lab. Clin. Med. 102:53, 1983.
- 2. Dern, R.J., Wiorkowski, J.J. and T. Matsade: J. Lab. Clin. Med. 75:37, 1970.
- 3. Meltzer, S.J. and W.H. Welch: J. Physiol. (London) 5:255, 1884.
- 4. Haradin, A.R., Weed, R.I. and C.F. Reed: Transfusion 9:229, 1969.
- 5. Lux, S.E. and K.M. John: Cell Shape and Surface Architecture (Alan R. Liss, ed.), New York, 1977.
- Bernstein, E.F., Indeglia, R.A., Shea, M.A. and R.L. Varco: Monograph No. 16, Am. Heart Assoc., Cardiovasc. Surg. 1:226, 1967.
- 7. Card, R.T., Mohandas, N. and P.L. Mollison: Br. J. Haematol. 53:237, 1983.
- 8. Murphy, S., Sayer, S.N. and F.H. Gardner: Blood 35:549, 1970.
- 9. Runicki, J.J., Tuccelli, M., Becker, G.A. and R.H. Aster: Transfusion 15:414, 1975.
- 10. Slichter, S.J. and L.A. Harker: Br. J. Haematol. 34:403, 1976.
- 11. Holme, S., Vaidja, S. and S. Murphy: Blood 52:425, 1978.
- 12. Murphy S. and S. Holme: Transfusion 20:624, 1980.
- 13. Murphy, S., Kahn, R.A., Holme, S., et al.: Blood 21:637, 1981.
- 14. Murphy, S. and T. Simon: Transfusion 21:637, 1981.
- Moroff, G., Friedman, A. and L. Robkin-Rline: Vox Sang. 42:33, 1982.
- 16. Murphy, S., Holme, S., Nelson, E. and R. Carmen: Transfusion 24:31, 1984.
- 17. Champion, A.B.: Effect of mode of agitation on in vitro platelet function during 7-day storage in the CLX system in platelet symposium, Cutter Biological, 1982.
- 18. McCullough, J., Weibblen, B.J. and D. Fine: Transfusion 23:20, 1983.
- 19. English, D. and B.R. Andersen: J. Immunol. Methods 5:249, 1974.
- 2C. McCullough, J.: Transfusion 20:129, 1980.
- 21. McCullough, J., Weibblen, B.J., Petersen, P.K. and P.G. Quie: Blood 52:301, 1978.

ACT IN CALL AND LODIE

- 22. Jacobson, M.S. and S.V. Kevy: Transfusion 20:443, 1980.
- Beel, F.P. and C.J. Schwartz: Biochem. Biophys. Acta 231:553, 1972.

ł.

- 24. Nelson, C.J.: Biochem. Biophys. Acta 144:221, 1967.
- Kevy, S.V., Button, L.N. and M.S. Jacobson: National Heart Lung Institute, NIH contract PHS-2968. Comprehensive Report National Technical Information Service Ref. No. PB-247,168 857, 1978.
- 26. Stern, I.J. and R.A. Carmen: Abstract International Soc. Hemat. and Trans. Montreal, 1980.
- 27. Seeman, P.: Biochem. Pharmac. 15:1632, 1966.
- 28. Roth, S. and P. Seeman: Nature, London 231:284, 1971.
- 29. Gibson, J.G. and G.W. Thorn: Report 11, NIH, 1954.



-17-

Initial Blood Storage Experiment (IBSE)

Scientific Experimental Protocol

Draft 5/85

:

Table of Co	ntents	
1.0	Introduction	
2.0	Blood Collection	
3.0	Component Preparation	
4.0	Records	
5.0	Pooling	
6.0	Aliquoting of Pooled Blood	
7.0	Samples for Baseline Studies	
3.0	Package and Shipment of Blood	X
9.0	Storage of Pooled and Aliquoted Bags	
10.0	Loading of Dewars	
11.0	Summary of Pre-Flight Timeline	
12.0	ADL FLight Preparations	
13.0	Postflight Aliquoting of Samples for Assays	
14.0	Detailed Experimental Procedures	
15.0	Acquisition of Data	
16.0	Data Analysis	
17.0	Safety Procedures	

-

Initial Blood Storage Experiment

Experimental Protocol

1.0 INTRODUCTION

The Initial Blood Storage Experiment (IBSE) will involve pre and post flight biomedical operations. Blood will be collected from a group of blood donors, and sterile suspensions of blood cells will be separated into sealed plastic blood bags. The bags will be placed in temperature controlled environmental stainless steel dewars. The dewars will be housed in three modified middeck lockers. Two lockers comprise the flight payload; the third dewar and a functional IBSE unit will be held on the ground in the laboratory facility.

This experimental protocol details the procedures that will be performed at the onsite facility, Central Florida Blood Bank, and base laboratories. All procedures, where applicable, are in accordance with the American Association of Blood Banks Technical Standards, Eleventh Edition and the Good Manufacturing Practices for Blood and Blood Components of the Code of Federal Regulations.

2.0 BLOOD COLLECTION

- 2.1 Blood will be drawn from a group of 36 volunteer donors under standard operating procedures and aegis of the Central Florida Blood Bank (CFBB), which will be responsible for collection and preparation of the components used in the experiment. A member of the IBSE team will be present.
 - 2.1.2 The donors will be preselected by CFBB according to the following guidelines:

a.) having successfully donated blood previously for transfusion

b.) having previously been shown to be negative for HBsAg

c.) having previously been shown to be negative for irregular antibodies against red cell antigens

d.) being of a single ABO and Rh type or at the discretion of the CFBB, two blood groups 12 units of one ABO and Rh for use as whole blood and 24 units of a second ABO and Rh for use as platelets and buffy coats.

e.) having not ingested aspirin the day of blood collection or 5 days previous.

f.) non-transfused males.

g.) having previously been shown to be negative for anti-HTLV III antibodies, ie. suspected AIDS virus antibodies.

- 2.1.3 FDA approved standard CPD anticoagulant blood bags will be used.
- 2.1.4 CFBB criteria will govern the acceptance or rejection of all donors at the time of phlebotomy; as well as the disposition of any units found to be unsatisfactory for any reason.
- 2.1.5 Standard CFBB laboratory testing will be conducted by CFBB after the collection, from pilot samples obtained at the time of phlebotomy, i.e. ABO, Rh, antibody screening, syphilis testing, hepatitis testing and anti-HTLV III testing. All testing will be complete prior to pooling. All donor and donation records will be retained by CFBB.
- 2.1.6 Units collected as whole blood and not intended for component preparation will be stored at 1-6 C until the time of pooling and aliquoting.
- 2.1.7 Units intended for production of components will be stored at room temperature (20-24 C) for no more than six hours.

3.0 COMPONENT PREPARATION

- 3.1 All component preparation will be completed by the CFBB according to their standard operating procedures with the exception of the final plasma volume of platelet concentrates. These concentrates will be resuspended to volume of 65 ml.
- 3.2 The red cells and plasma from the units used for platelet and buffy coat preparation will be stored at 1-6 C for 21 days.

4.0 RECORDS

4.1 Pooling Records:

A Pooling Data Sheet is prepared prior to pooling each blood component. The accuracy of all data must be verified by a second person, who initials the sheet, before pooling. The following information must be documented:

- Component (whole blood, platelets, or leukocytes)
- Collection date
- Pooling date and time
- Temperature of pooling area (room temp, 5°C, etc.)
- Lot number and description of pooling bag
- Whole blood number of each component
- ABO and Rh type of each component
- Identity of person pooling
- Identity of person verifying records.

The time pooling is completed will be documented.

Initial Blood Storage Experiment POOLING DATA SHEET

Component: Collection D Pooling Bag Whole Blood	Whole BI Platelet Leukocyt ate: Descripti Lot Numb	ood s es on: er: Technician	Temp. of Technician	Pooling Date: Time Begun: Time Completed: f Pooling Area:
Number	ABO/Rh	Pooling	Reviewing	Remarks
<u></u>				
<u></u>				
· <u> </u>				

Initial Blood Storage Experiment PRE-FLIGHT ALIQUOTING DATA SHEET

Whole Blood 250 ml

Date:_____

Polymer Code: PVC-DEHP = D PVC-TOTM = T Polyolefin = PO Time Begun:_____

Time Completed:_____

Temperature of Aliquoting Area:_____

.

	Bag ID	Baluman	Gross	Net Weicht(c)	Tech	Roview	Remarks
<u>Uraer</u>	_Code#	Forymer	Heighildy	_ METQHI (QZ_		INCY ICH	
	1	P					
_2	8	Ţ					
	15	P0					
_4		D					
5	9						
6	16	PO					
_7	3	D					
_8	10	T					
9	17	Po					
10	4	D					
11	11	T					
12	18	PO					Spare Bag
13	5	D					
14	12	T			,		Spare Bag
15	19	P0					
16	66	D					
<u>17</u>	13						
18	20	P0					
19	7	D					Spare Bag
20	14						·
21	21	P0					

Initial Blood Storage Experiment PRE-FLIGHT ALIQUOTING DATA SHEET

Leukocytes 75 ml

Date:_____

Polymer Code: PVC-DEHP = D PVC-TOTM = T Polyolefin = P0 Time Begun:_____

Time Completed:_____

Temperature of Aliquoting Area:_____

Order	Bag ID Code#	Polymer	Gross Weight(g)	Net Weiaht(a)	Tech.	Review	Remarks
1	22	D					
 °	20						
				-			
	36	<u> </u>					
_4	23	D				•	
_5		ļ					Spare Bag
_6	37	P0					
7	24	D					
8		т					
9	38	PO					
10	25	D					· ·
11	32	T					
12	39	P0					
13	26	D					
14	33	Т					
15	40	PO					
<u>16</u>	27	D					Spare Bag
17	34	L_T					
18	41	PO					
19	28	D					
20	35	т					
21	42	PO					Spare Bag

Initial Blood Storage Experiment PRE-FLIGHT ALIQUOTING DATA SHEET

Platelets 65 ml

Date:_____

Polymer Code: PVC-DEHP = D PVC-TOTM = T Polyolefin = PO Time Begun:_____

Time Completed:_____

Temperature of Aliquoting Area:_____

Order	Bag ID Codei	Polymer	Gross Weight(g)	Net Weight(g)	Tech.	Review	Remarks
<u>vi.ver</u>	43	0					
	<u> </u>				· · · · · · · · · · · · · · · · · · ·		
		<u> </u>					Spare Bag
<u> </u>		P0					
_4	44						
_6	58	<u> </u>					
_7	45	D					
_8	52	<u> </u>					Spare Bag
_9		P0					
10	46						
11	53	Ţ					
12	60	PO					Spare Bag
<u>13</u>	47						
14	54	T					
15	61	P0					
16	48	D					
17	55	Ţ					
18	62	P0					
<u> 19</u>	49	D					Spare Bag
20	56	<u></u> т					
21.	63	P0					
22	98	Т					
07	00	т					

4.2 Records of Pre-flight Aliquoting of Pools:

A Pre-Flight Aliquoting Data Sheet is prepared prior to aliquoting each pool into the coded experiment blood bags. The order in which the bags are to filled is pre-deterimined on the data sheet. The aliquoting technician will verify the labeled bags as being of the correct polymer and having the correct code. After the aliquoting, a second person will review the records and initial to verify accuracy.

5.0 POOLING

All pooling and aliquoting will be performed using aseptic techniques. All bench area will be washed down with a 10% bleach solution prior to use and after any blood spills.

5.1.0 Whole Blood Pool

- 12 units of whole blood will be pooled in groups of 6 units to make two pools. The two pools will then be mixed together to make one pool. Whole blood pooling may be carried out at room temperature, but after aliquoting the bags should be placed immediately at 1-6°C.
- 5.1.1 Remove the segments from each of the units of whole blood. Place the segments from each unit in separate plastic bags labeled with the corresponding whole blood number and store between 1-6°C for 21 days.
- 5.1.2 Using the Dupont Sterile Connection Device (SCD), connect the integral tubing from a unit of whole blood to the tubing of a cytoagglomerator bag. Allow the whole blood to empty into the cytoagglomerator bag. Once the bag is empty, clamp the tubing with a hemostat and heat seal the integral tubing in three places close to the port of the whole blood bag.
- 5.1.3 Disconnect tubing at the center seal. Repeat with 5 more units of whole blood.
- 5.1.4 Repeat the above procedure with the remaining 6 units of whole blood and a second cytoagglomerator bag.
- 5.1.5 Once the second bag is filled, the two pooling bags must be mixed together. Using the SCD, connect the two cytoagglomerator bags.
- 5.1.6 Hang one bag above the other. Allow the blood in the hanging bag to flow into the second bag, mixing gently until the hanging bag is empty. This should take about 10 minutes.
- 5.1.7 Reverse positions of the two bags and repeat until the the two bags are mixed a total of ten times.
- 5.1.8 Store whole blood pool at 1-6°C until ready to aliquot.

5.2 Platelets

All platelet concentrates must rest for 1 hour and 15 minutes and be rotated for a minimum of 1 hour and 15 minutes prior to pooling. All platelet pooling will be performed at room temperature.

- 5.2.1 Pool 24 units of platelet concentrates into a 2 liter plasma transfer pack.
- 5.2.2 Using the SCD, connect a 2 liter plasma transfer pack with a unit of platelet concentrate. Hang the platelet unit above the transfer pack and allow the platelets to gravity feed into the transfer pack. When the platelet bag is empty, clamp the integral tubing with a hemostat. Heat seal the integral tubing 3 times close to the port of the platelet bag. Disconnect at the center seal.
- 5.2.3 Repeat until all 24 platelet units have been pooled.
- 5.2.4 Mix by inversion.
- 5.2.5 Store pooled platelets at room temperature until ready to aliquot.
- 5.3 Leukocyte Pool
 - Preparation of leukocytes will be performed under a laminar flow hood.
 - 5.3.1 24 units of buffy coats will be sedimented with dextran to separate the leukocytes. Suspensions of leukocyte will be washed in plasma transfer packs. See white cell protocol, Section 14.3.
 - 5.3.2 Connect a 2 liter plasma transfer bag to a unit of leukocytes, using the SCD. Hang the leukocytes above the plasma transfer pack and allow the leukocytes to gravity feed into the transfer pack. When the leukocyte bag is empty, clamp the integral tubing with a hemostat. Heat seal the integral tubing 3 times close to the port of the leukocyte bag. Disconnect the 2 bags at the center seal.
 - 5.3.3 Repeat until all leukocytes have been pooled.
 - 5.3.4 Mix by inversion.
 - 5.3.5 Store pooled leukocytes at 1-6°C until ready to pool.

5.4 Plasma

- 5.4.1 Select 8 units of plasma from units used for platelet preparations for pooling into a 2 liter plasma transfer pack. The pooled plasma will be used for tests on platelets.
- 5.4.2 Using the SCD, connect a unit of plasma to a 2 liter plasma transfer pack. Hang the plasma above the transfer pack and allow the plasma to gravity feed into the transfer pack. When the plasma bag is empty, hemostat the integral tubing. Heat seal the tubing

close to the port of the plasma bag. Disconnect the 2 bags at the center seal.

5.4.3 Repeat until 8 bags of plasma have been pooled.

- 5.4.4 Mix by inversion.
- 5.4.5 Store pooled plasma at 1-6°C until needed.

6.0 ALIQUOTING

Each product pool will be aliquoted into 3 different, FDA OoB approved plastic formulations: PVC-DEHP, PVC-TOTM, and Polyolefin. Each whole blood and leukocyte pool will be aliquoted into 21 samples in 300 ml bags. The platelet pool will be aliquoted in 23 samples in 300 ml bags. There will be 7 bags of each plastic formulation for each whole blood and leukocyte pool. For the platelet pool, there will be 7 bags PVC-DEHP and 7 bags polyolefin and 9 bags PVC-TOTM. The total number of bags will be sixty-five. The volumes vary according to each product.

6.1 Prelabeling of bags

Each bag will be given a unique number which will identify the blood component, the type of plastic, the dewar and locker used for storage.

	PVC-DEHP	PVC-TOTM	Polyolefin
Whole Blood	1-7	8-14	15-21
Leukocytes	22-28	2 9– 35	36-42
Platelets	43-49	50-56,98,99	57-63

Spare bags are: Whole Blood: 7,12,18

Leukocytes: 27,30,42

Platelets: 49,52,57

	6.1.	2 C	ode	for lo	cati	on and	stor	age	e cor	nditio	n			
			ORBI	TER							EA	RTH		
	WB	LEUK	WB	LEUK	WB	LEUK			WB	LEUK	WB	LEUK	WB	LEUK
A	6	22	10	29	20	37		D	4	40	13	32	15	28
B	8	41	17	26	5	35		E	11	24	19	38	2	33
с	2 1	34	3	39	14	23		F	16	31	1	25	9	36
P	56	48						P*	60	50				
L A	62	53						A T	47	98				
T B	43	58						E L	51	63				
E	54	99						E	59	46				
S	44	61						Ŝ	45	55				
FF ()	NONT				B (B	ACK SOTTOM)		F) ('	RONT TOP)				(BACK (BOTTOM)

Letter codes are dewars. A and B are cold dewars stored in one orbiter locker. C is a cold dewar and PLATELETS is a warm dewar stored in the orbiter locker. D and E are cold dewars stored in one earth locker. F is a cold dewar and PLATELETS* is a warm dewar stored in the other earth locker.

6.2 Whole Blood

6.2.1 Whole blood will be aliquoted in 250 ml volumes in the following prelabeled 300 ml bags and stored at 1-6°C until loaded in the dewars:

1.	PVC DEHP	7	bags
2.	Polyolefin	7	bags
3.	PVC TOTM	7	bags

6.3 Platelets

6.3.1 Platelets will be aliquoted in 65 ml volumes in the following prelabeled 300 ml bags and stored at room temperature (20-24°C) until loaded in dewar:

1.	PVC DEHP	7 bags	
2.	Polyolefin	7 bags	
3.	PVC TOTM	9 bags	

- 6.4 Leukocytes
 - 6.4.1 Leukocytes will be aliquoted in 75 ml volumes in the following prelabeled 300 ml bags and stored at 1-6°C until loaded in the dewars:
 - 1. PVC DEHP7 bags2. PVC TOTM7 bags3. Polyolefin7 bags
- 6.5 Procedure for aliquoting
 - 6.5.1 Whole blood
 - 6.5.1.2 Mix the whole blood pool well prior to aliquoting and again between each bag.
 - 6.5.1.3 Place a 300 ml bag to be filled on scale and tare for the weight of the bag.
 - 6.5.1.4 Using the SCD, connect one cytoagglomerator bag containing whole blood pool with the 300 ml bag for aliquoting.
 - 6.5.1.5 Allow 250 g of whole blood to flow into the above prelabeled bag.
 - 6.5.1.6 When bag is filled with the appropriate volume, clamp with a hemostat. Heat seal tubing 3 times close to the port of the aliquoted bag. Disconnect the 2 bags at the center seal.
 - 6.5.1.7 Repeat with each prelabeled bag for whole blood until 11 bags are filled or cytoagglomerator bag is empty.
 - 6.5.1.8 Using the sterile docking device connect second cytoagglomerator bag to a bag for aliquoting. Repeat steps 6.5.1.2-6.5.1.7.
 - 6.5.1.9 The blood remaining in the cytoagglomerator bag will be aliquoted into plastic tubes for use in baseline studies. See 7.0.
 - 6.5.1.10 Store whole blood aliquots at 1-6°C until ready to load in the dewar.

6.5.2 Platelets

- 6.5.2.1 Mix platelet pool well before aliquoting and between each bag during aliquoting.
- 6.5.2.2 Place bag for aliquot on scale and tare for the weight of the bag.
- 6.5.2.3 Using the SCD, connect the platelet pool with a prelabeled bag for aliquoting.

- 6.5.2.4 Allow 65 g of platelet suspension to flow into aliquoting bag. Clamp with hemostats when the appropriate volume has been aliquoted. Heat seal the integral tubing 3 times, close to the port of the aliquot bag. Disconnect bags at the center seal.
- 6.5.2.5 Repeat until all 23 bags are filled.
- 6.5.2.6 The platelet suspension remaining in the pooled bag will be aliquoted into tubes for baseline studies. See 7.0.
- 6.5.2.7 Aliquoted platelets will be stored at room temperature (20-24°C) until loaded in the dewar.

6.5.3 Leukocytes

- 6.5.3.1 Mix leukocyte pool well before aliquoting and between each bag during aliquoting.
- 6.5.3.2 Place bag for aliquot on scale and tare for the weight of the bag.
- 6.5.3.3 Using the SCD, connect the leukocyte pool with a prelabeled bag for aliquoting.
- 6.5.3.4 Allow 75 g of leukocytes to flow into aliquoting bag. Clamp with hemostats when the appropriate volume has been aliquoted. Heat seal the integral tubing 3 times, close to the port of the aliquot bag. Disconnect bags at the center seal.
- 6.5.3.5 Repeat until all 21 bags are filled.
- 6.5.3.6 The leukocyte suspension remaining in the pooled bag will be aliquoted into tubes for baseline studies. See 7.0.
- 6.5.3.7 Aliquoted leukocytes will be stored at 1-6 until loaded in the dewar.

7.0 SAMPLES FOR BASELINE STUDIES

Time zero samples are those taken at the time of aliquoting. Their purpose is to verify the condition of the products pre flight. All assays outlined in the protocol will be performed on time zero samples. Each investigator will receive samples according to the description below.

7.1 Whole blood

Whole blood pool must be mixed well prior to aliquoting.

7.1.1 Using a hemostat, clamp the integral tubing of the cytoagglomerator bag. Cut the end of the tubing. Place the tubing in a prelabeled

centrifuge tube. Release the hemostat and allow the require volume (listed in 7.1.2) to flow into the tube. When the tube is filled with the require volume reclamp the hemostat and repeat until all whole blood tubes are filled.

7.1.2 Whole blood will be aliquoted into the following tubes labeled whole blood and place on ice or to maintain a temperature of $1-6^{\circ}C$.

# of Tubes	Tube Volume	Volume Needed	Investigator
3 1 2	50 ml 50 ml 15 ml	50 ml 30 ml 10 ml	M. Jacobson (A.B.[1]) I. Szymanski M. Jacobson D. Ausprunk
2 2	10 ml 15 ml	10 ml 5 ml	M. Jacobson (A.S.) M. Jacobson (S.K.) W. Curby

7.2 Platelets

Platelets must be well mixed prior to any aliquoting.

- 7.2.1 Using a hemostat, clamp the integral tubing of the platelet bag. Cut the end of the tubing. Place the tubing in a prelabeled centrifuge tube. Release the hemostat and allow the require volume (listed in 7.2.2) to flow into the tube. When the tube is filled with the require volume reclamp the hemostat and repeat until all platelet tubes are filled.
- 7.2.2 Platelets will be aliquoted into the following tubes labeled platelets and left at room temperature:

# of	Tube	Volume	Investigator
Tubes	Volume	Needed	
1 1 1 1 1	3 ml syringe 50 ml 50 ml 15 ml 15 ml 15 ml	1.5 ml 30 ml 15 ml 5 ml 4 ml 3 ml	M. Jacobson F. Chao F. Chao I. Szymanski D. Ausprunk M. Jacobson W. Curby

7.3 Leukocytes

Leukocyte pool must be well mixed prior to aliquoting.

7.3.1 Using a hemostat, clamp the integral tubing of the leukocyte bag. Cut the end of the tubing. Place the tubing in a prelabeled centrifuge tube. Release the hemostat and allow the require volume
(listed in 7.3.2) to flow into the tube. When the tube is filled with the require volume reclamp the hemostat and repeat until all whole blood tubes are filled.

7.3.2 Leukocytes will be aliquoted into the following tubes labeled leukocytes and place on ice or to maintained at a temperature between 1-6°C:

# of Tubes	Tube Volume	Volume Needed	Investigator
2	50 ml	30 ml	F. Lionetti
2	15 ml	5 ml	D. Ausprunk
1	2 ml	l ml	W. Curby

- 8.0 PACKAGING AND SHIPMENT OF BLOOD
- 8.1 Whole blood, red blood cells, leukocyte preparations and plasma will be packed and shipped on wet ice to maintain at a temperature of 1-6°C.
- 8.2 Platelets will be packed in IBSE shipping racks and shipped at 20-24°C.
- 8.3 A member of the IBSE team will observe the packing of blood components prior to shipment. This person will accompany the blood in the helicopter to KSC. A second person will accompany the ground stored blood to Boston.

9.0 STORAGE OF ALIQUOTED POOLED SAMPLES

Of the sixty-three aliquoted bags, fifty-six will be used in the experiment. One extra bag of each plastic and each component will be available in the event a replacement is required. There are 2 dewars to each locker and 4 lockers total. 2 lockers will remain in the laboratory and 2 lockers will be placed on the orbiter. The following table is a the description of how the products and the specific bag types will be divided.

Earth stored aliquots Orbiter stored aliquots

Whole blood:

3 PVC-DEHP 3 PVC-TOTM 3 Polyolefin	3 PVC-DEHP 3 PVC-TOTM 3 Polyolefin
Leukocytes:	
3 PVC-DEHP	3 PVC-DEHP
3 PVC-TOTM	3 PVC-TOTM
3 Polvolefin	3 Polyolefin

Platelets:

3 PUC-DEHP	3 PVC-DEHP
	4 PVC-TOTM
4 FVC-1010 9 Deluciofin	3 Polyolefin
2 LOTAOTELTU	

A total of 56 bags, 18 whole blood, 18 leukocytes and 20 platelets.

9.1 Whole blood and leukocytes

Whole blood and leukocytes will be stored in a total of 6 cold dewars. Three of the cold dewars will be left in the laboratory under controlled conditions and three will be placed on the orbiter. Each cold dewar will contain 3 whole blood aliquots and three leukocyte aliquots in the configuration described in 6.1.2.

9.2 Platelets

Platelets will be stored in a total of 2 warm dewars (22°C). One warm dewar will be left in the laboratory under controlled conditions and one will be placed on the orbiter. Each dewar will contain 9 platelet aliquots and one "dummy" bag in the configuration described in 6.1.2.

The "dummy" bag will contain platelets from the pool and aliquoted into a PVC-TOTM bag.

10.0 LOADING OF DEWARS

- 10.1 This activity requires 3 IBSE staff members (2 CBR and 1 ADL staff) to manipulate the hardware. These three individuals will work in a resricted area, 0 + C Building, room 2249 or 2251, until all dewars are loaded and ADL closeout of the lockers is complete. All clerical entries will be verified and initialed by a second person.
- 10.2 Two Curby cranes for suspending IBSE hardware over the lockers will be in place for the loading procedure.
- 10.3 Each dewar will be prelabeled by ADL (A-D for Orbiter, E-H for Earthstored dewar). The bag types, blood components, and temperature of each dewar by label are defined in section 6.2.
- 10.4 The microgravity-stored dewars will be loaded at KSC and the ground stored will be loaded simultaneously in Boston. The procedure is identical for microgravity and Earth-stored dewars.
- 10.5 Selection of blood components:
 - 10.5.1 Prior to loading the dewars, filled, coded bags will be separated into 6 identical groups for whole blood and leukocytes, each containing 3 whole blood units (one in each type of plastic bag), 3 leukocyte units (one in each type of plastic bag), and 2 identical groups for platelet units.

- 10.5.2 Before loading commences, the code number of each bag to go in the dewar will be recorded on the appropriate worksheet.
- 10.5.3 A plastic protective drape will be placed over the lower locker assembly to prevent spillage into the hardware in the unlikely event of a broken blood bag during the loading procedure.
- 10.6 Loading will begin with platelet bags since they are stored at room temperature.
 - 10.6.1 Platelet bags will be secured in vertical ____m holders as described in section ___, arranged by bag type as designated in 6.1.1. The final position of each coded bag will be recorded on the appropriate worksheet.
 - 10.6.2 The platelet holders will be joined by rubber bands maintaining the relative positions defined in section 6.1, and the entire unit slipped into the 22 + 1°C dewar compartment.
- 10.7 Cold-stored components intended for a single dewar (3 whole blood and 3 leukocytes) will be removed from refrigerated storage and brought to the dewar loading site. Each bag will be removed from its protective plastic overwrap and wiped dry with towelling to remove condensation prior to loading.
 - 10.7.1 Whole blood and leukocytes will be loaded alternately in the arrangement described in section 6.1.2, beginning at the bottom of the 5 + 1°C dewar.
 - 10.7.2 Each bag must be folded at the top, label side up, maintaining maximum surface area against the air flow grid surface. The top of the bag is folded back onto the label and is positioned at the top of the individual compartment.
 - 10.7.3 The folded blood bag is placed horizontally in a parallel folded Teflon sheet which facilitates placement of the blood bag into the hardware without sliding the bag across any rough surfaces. Then, while the blood bag is held in place from the opposite side, the teflon sheet "envelope" is withdrawn. Corners of the plastic bag may not protrude out of the prescribed area in order for the hardware to lower properly into the locker.
 - 10.7.4 To fully contain any protruding bag corners and to ensure ease of lowering, each loaded bag will be encircled in the hardware by Kapton tape (Connecticut Hard Rubber). Whole blood will be encircled by 1" tape and leukocytes will be encircled by 1/2" tape. Precut strips of tape will be attached to a vertical metal bar in the hardware and circled in the direction of the 3 adjacent bars and finally around the open space allowed for loading.
- 10.8 The 2-dewar unit is ready for ADL closeout. Repeat with Earth-stored lockers and dewars.

11.0 SUMMARY OF PRE-FLIGHT TIMELINE

The launch is scheduled for 7 a.m., 1985. The IBSE lockers must be turned over to NASA 12 hours prior to that, 7 p.m., ______ 1985. In order to minimize pre-flight exposure to 1 x g, the donor blood will be collected early on the morning before and all subsequent preparations are defined within that 12-13 hour period.

- 11.1 Day Prior to Launch
 - 6:00 9:00 a.m. Collection of blood from CFBB donors. 6:00 - 8:00 24 donors for platelet and leukocyte production. 8:00 - 9:00 12 donors for whole blood.
 - 8:00 10:30 Platelet concentrates rest and rotate.
 - 8:00 11:00 Leukocytes prepared from buffy coats at CFBB.
 - 9:00 2:00 p.m. Whole blood units pooled and divided into coded IBSE bags, at CFBB, 1 person.
 - 10:30 2:00 Platelets and leukocytes pooled and divided into coded IBSE bags, at CFBB, 2 people.
 - 2:00 3:30 Prepared Blood components and 1 IBSE staff transported to KSC by helicopter. Other IBSE personnel needed for loading the dewars at KSC must travel by car.
 - 3:30 7:00 Load orbiter and ground dewars. ADL closeout of orbiter and ground lockers.
 - ** 7:00 p.m. Turn over orbiter dewars to NASA.

12.0 ADL FLIGHT PREPARATIONS

The experiment requires eight nearly identical environmental chambers. During the preflight preparatory phase, suspensions of human blood formed elements will be prepared, introduced and sealed into plastic bags (containers), which will then be mounted in fixed positions within the environmental chambers as described above. These chambers will then be closed, the associated environmental control systems (temperature control and air supply) will be activated and the entire system will be enclosed in modified NASA storage lockers for flight. The environmental chambers must then operate continuously from the time of insertion and closure of the blood preparation into the chambers until the locker assembly has been removed from the orbiter and returned to IBSE personnel.

12.1 The design of the chambers and their associated environmental control systems must assure that the following environmental specifications will be met from the point of preflight loading and closure to the termina-

tion of flight and the lockers have been turned over to the IBSE personnel (see 12.1 Table I).

12.1 TABLE I

Chambers needed	Whole Blood	Platelets	Leukocytes
Capacity	250ml	65ml	75ml
Thickness of bags with blood	3.0cm	0.5cm	1.5cm
Temperature Range	5° <u>+</u> 1°C	22° <u>+</u> 1°C	5° <u>+</u> 1°C
Air	Cabin air	Cabin air	Cabin air
Orientation of the bag during take off	Horizontal	Vertical	Horizontal

The air flow through the chambers will be a minimum of 85 cc/min.

- 12.2 All blood preparations will be placed in 300 ml volume plastic bags. Three different formulations of plastics will be used for the bags. The dimension of the filled blood bag will vary with the component place in it due to differing volumes being inserted into the bags. The configuration of the chambers must ensure that the bags do not come in contact with one another, are held in position to avoid any shifting, and are spaced in uniform distances. The chambers must provide a constant flow of air to maintain the exchange of oxygen and carbon dioxide across the wall of the blood bags containing the formed elements of the blood.
- 12.3 The power and air supply must have capability of maintaining the restrictions stated in the Table I for a total of 10 days. The temperature must remain within stated ranges. Should there be an interruption of power, the temperature within the chambers containing red cells and leukocytes should not rise above TBD for more than TBD at any one time. The temperature of the chamber containing platelets should not rise at above TBD for more than TBD at any one time.
- 12.4 In addition to the design of the environmental chambers, ADL is responsible for the following:
- 12.5 Proper placement of the insert holding the blood bags in the dewars, activation of the environmental system and closure of the lockers.
- 12.6 Packaging and delivery of the lockers to KSC personnel.
- 12.7 Assuring that the proper environmental conditions are maintained pre and post flight.

12.8 Placement of the lockers on and off the orbiter

The lockers will be loaded and delivered to KSC personnel a minimum of 12 hours prior to launch. During this period, specified environmental conditions must be maintained within the lockers.

13.0 POSTFLIGHT ALIQUOTING OF SAMPLES FOR TESTING

13.1 Labeling of plastic tubes for aliquoted samples

All blood will be aliquoted into prelabeled prearranged plastic centrifuge tubes, test tubes or syringes. All labeling of tubes will be completed prior to the landing of the orbiter. All labeling and aliquoting will be performed in a secluded area by authorized personnel.

13.2 Whole Blood

For each of the 18 whole blood units label with corresponding code numbers the following tubes:

# of Tubes	Tube Volume	Volume Needed	Investigator
3.	50 ml 50 ml	50 ml 30 ml	M. Jacobson (A.B.[1]) I. Szymanski
2	15 ml	10 ml	M. Jacobson D. Ausprunk
2	10 ml	10 m1	M. Jacobson (A.S.)
1	15 ml	5 ml	M. Jacobson (S.K.)
1	15 ml	5 ml	W. Curby

13.3 Leukocytes

For each of the 18 leukocyte units label with corresponding code numbers the following tubes:

# of	Tube	Volume	Investigator
Tubes	Volume	Needed	
2	50 ml	30 ml	F. Lionetti
1	15 ml	5 ml	D. Ausprunk
1	15 ml	5 ml	W. Curby

13.4 Platelets

For each of the 20 platelet units label with corresponding code numbers the following tubes:

# of Tubes	Tube Volume	Volume Needed	Investigator
1	3 ml syringe	1.5 ml	M. Jacobson
1	50 ml	30 ml	F. Chao
1	50 ml	15 ml	F. Chao
1	15 ml	5 ml	I. Szymanski
1	10 == 15 ml	4 ml	D. Ausprunk
1	15 ml	3 ml	M. Jacobson
1	15 ml	l ml	W. Curby

13.5 Removal of bags from the dewars

TBD

- 13.6 ADL will unload the lockers (separate procedure TBD).
- 13.7 CBR personnel will unload the dewars (separate procedure TBD).

The platelet dewar will be unloaded and platelets place on a rotator for approximately 30 minutes prior to aliquoting.

13.8 Aliquoting

- 13.8.1 For each component a specified volume will be aliquoted into pre labeled tubes according to the charts in 13.2-13.4
- 13.8.2 Using a coupler with a spike on each end, enter a bag. Be sure the roller clamp is closed. Hang the bag above the rack of prelabeled tubes. Open the clamp and allow the required volume of sample to flow into premarked tubes.
- 13.9 Once the earth stored dewars are unloaded and samples aliquoted the same procedure will be repeated for the orbiter stored dewars.
- 13.10 Aliquoted samples will be stored at the following temperatures until received by the investigator:

Whole blood and leukocytes 1-6°C Platelets 20-24°C

13.11 Only when all aliquoting operations are complete, will each investigator will receive a rack of coded samples. All assays will be performed as described in 14.0 Detailed Experimental Procedures

14.0 DETAILED EXPERIMENTAL PROCEDURES

PRECEDING PAGE BLANK NOT FILMED

```
IBSE Protocol
Page 21
```

Page 21 14.1 Erythrocyte Studies (Technical Proposal 3.2.2)

Sample Aliquoting and Designated Responsibility for the Erythrocyte Study

Personnel:

D. Ausprunk (D.A.); A. Byrnes (A.B.); W. Curby (W.C.); M. Jacobson (M.J.); S. Kevy (S.K.); A. Scanlan (A.S.); I. Szymanski (I.S.); L. Wolfe (L.W.).

	Whole Blood	Responsible	Performance
TESTS	Volume (ml)	Person	
Physical Measurements			
Mean Critical Hemol. vol	5	W.C.	on site
Osmotic Fragility	5	S.K., A.S.	on site
Morphology and Shape Change (E.M.)	10	D.A.	on site and Boston
pH, P02, PC02		A.S.	on site
<u>Membrane Properties</u> Intact Cells - Lipids	10	M.J.	on site and Boston
Ghosts	10	A.B., L.W.	on site
Spectrin Extraction	50	Α.Β.	on site and Boston
Spectrin-Actin-4.1 Bindi	ing	Α.Β.	Boston
Changes in Membrane Prot	teins	Α.Β.	Boston
Membrane Glycoproteins	50	M.J.	on site and Boston
<u>Metabolites</u>			
ATP, DPG, Glucose	(red t 10 tube	op) A.S.	on site and Boston
Electrolytes	(.	1	
WB K+, Na+, plasma K+, I	(autoc Na+ 10 tube) A.S.	on site and Boston
plasma Hgb, WB, Hgb		A.S.	on site and Boston
<u>Plasma</u> Lipids	50	M.J. M.J.	on site on site and Boston
plasticizer quant.		M.J.	on site and Boston
coagulation factors		M.J.	on site and Boston
membrane dust		M.J., A.B.	
Immunology	30	I.S.	on site and Boston

ON-SITE STUDIES

18 Samples

A. MEMBRANE PROPERTIES

1. Intact Cells

a) 10 cc of each whole blood sample will be washed by centrifugation 3X with cold buffered saline and the buffy coat removed by passing through a cellulose column. The cells will be resuspended in 3 cc of PBS.

b) 0.5 cc will be saved for electronic Hematology parameters by Coulter analysis.

(1) Lipid Analysis: 1 ml cells will be extracted with chloroform-isopropanol and back extracted with 0.05 M KCl. The extracts will be stored frozen (-20°C) and will be transported back to the laboratory at CHMC in Boston for completion of analysis (see below for methodology).

2. Red Cell Ghosts

a) 50 cc of whole blood will be prepared into ghosts by a modification of the method of Dodge (see below). The ghosts will be counted electronically by Coulter analysis.

(1) Protein Structure will be analyzed by two-dimensional SDS gel electrophoresis. Ghosts will be solubilized in 0.5% SDS plus 0.5 mM EDTA and with and without 0.02 M dithiothreitol, incubated 1 hr (37°C), frozen and brought back to Boston (-40°C). The ghosts will be analyzed in Boston by 5% polyacrylamide gel electrophoresis, as described by Fairbanks et al. The resulting gels will be stained for protein with Coomassie Blue and for glycoproteins with PAS (see below).

(2) <u>Spectrin Extraction</u>: Spectrin dimer will be isolated from 10 cc of ghosts by low ionic strength extraction at 37°C and purified by gel chromatography at 4-6°C (see below). Samples will be stored at 4-6°C and transported back to CHMC laboratory in Boston. See protocol below for Spectrin-Actin Protein-4.1 assay.

(3) <u>Glycoprotein Labeling</u>: See section under platelets.

B. RED CELL METABOLISM

1. 10 ml of whole blood will be used to determine the metabolic state of the red cell.

a) 2,3 Diphosphoglyceric Acid: 1 ml of whole blood will be immediately deproteinized by pipetting into 3.0 ml Trichloroacetic Acid in a 15 ml polypropylene centrifuge tube.

(1) Vortex immediately and sit in ice for 5 min. Centrifuge at 2500 x g for 10 min and take off clear supernatant. Freeze (-20°C) and transport back to CHMC laboratory in Boston for completion. See protocol.

b) ATP: Adenosine-5-triphosphate

(1) Pipet into a 15 ml centrifuge tube 4.0 ml of 0.6 N perchloric acid (ice cold) and 1.0 ml whole blood. Mix well (Vortex).

(2) Allow to stand in ice for 10 min then centrifuge for 10 min at 2500 x g at 4° C.

(3) Take off supernatant and freeze at -40°C for transport back to CHMC lab in Boston for completion. See protocol.

c) Glucose

(1) Pipet 0.4 ml whole blood into 15 ml centrifuge tube into 3.6 ml 3% trichloroacetic acid. Centrifuge for 10 min at 2500 x g. Take off supernatant and freeze for transportation back to CHMC lab in Boston for completion.

d) Plasma Hemogloblin

(1) 3 ml of whole blood will be placed into a plastic 7 ml centrifuge tube and centrifuged in a Sorvall at 4°C for 30' at 5000 rpm. The supernatant is taken off and respun. The supernatant will be stored at 4-6°C for transport back to the CHMC lab in Boston. See protocol.

e) The remaining whole blood will be transported at 4-6°C back to the CHMC lab for K+, HgB, Na+ determinations.

2. <u>Plasma</u>: 25 ml of whole blood in a 50 ml plastic centrifuge tube will be centrifuged at 4200 x g for 5 min. The plasma will be separated.

a) 3 ml aliquot will be frozen at -20°C for factor analysis and transported back to the CHMC laboratory.

b) 3 ml aliquot will be extracted for lipid analysis in chloroform-methanol and back extracted with 0.2 M KCl. Extracts will be stored frozen (-20°C) for return to the CHMC laboratory. See protocol for analysis.

c) The remaining aliquot will be frozen at -20°C and transported back to the CHMC lab for plasticizer quantitation. See protocol for analysis.

C. Physical Properties

1. <u>Osmotic Fragility</u>: 1 ml of whole blood will be used for this test. A unopette osmotic fragility kit will be used for these determinations. Fragilities will be determined with and without incubation for 1 hour at 22°C in buffered saline. See protocl for analysis.

2. <u>pH</u>: 0.5 ml of whole blood will be placed into a 1 ml plastic vial and stored frozen at -20°C until pH is read.

PREPARATION OF RBC GHOSTS

Modification of method of Dodge, J.T., Mitchell, C., and Hanahan, D.J.: Arch. Biochem. Biophys. <u>100</u>:119-130 (1963).

<u>Media</u>

5mM Isotonic phosphate buffered saline, pH 8.0. 0.1 mm EGTA

Procedure

- Approximately 4 cc anticoagualted blood washed 3x with cold buffered saline. Resuspended in ~3 cc of B.S. Send ~.5 cc or less to lab for Hgb and RBC. WBC must be removed when cells hemolyzed.
- 2. 2.5 cc of resuspended RBC's diluted to 35 cc in a graduated cylinder with cold 30 mosm buffer and inverted several times to ensure complete mixing and hemolysis of the red cells. (Ghosts should be kept at ~4°C at all times during the procedure). Approximately 9 cc apportioned to each of 4 centrifuge tubes. Graduate rinsed once with ~5 cc 30 mosm buffer, which is added to tubes.
- 3. Centrifuge tubes balanced, centrifuged at 9000 x g for 10 min in cold.
- 4. Hemolysate aspirated and button of ghosts resuspended with stream of cold 30 mosm buffer from a wash bottle.
- Ghosts in each centrifuge tube are resuspended with ~10 ml of 30 mosm buffer. Tubes are inverted several times to ensure complete suspension.
- 6. Centrifugation, aspiration of supernatant and reusspension is repeated, using cold 30 mosm buffer, until supernatant is completely clear and ghost button is white. At this point all ghosts are combined in one centrifuge tube, resuspended in 10 cc, 30 mosm buffer and collected a final time (usually a total of 4-5 centrifugations will be required).
- 7. If a small button of debris and/or unhemolyzed BRC's remains below the ghost button, this is left behind when ghosts are combined in the last step (6 above).
- 8. Ghosts are resuspended at the desired volume (usually 5 ml) with 30 mosm media. They may be washed and stoerd in other media as desired. Recovery = usually 90-95% Yield = 5 ml of ghosts with a concentration of ~3x10° ghosts/ml.

Ghosts may be counted in Coulter counter or in RBC counting chamber (phase microscopy) after standing 30 min.

SDS GEL ELECTROPHORESIS

REFERENCE: Fairbanks, G., Steck, T.L., and Wallach, D.F.H. Electrophoretic analysis of the major polypeptides of the human erythrocyte membraned. Biochemistry <u>10</u>:2606-2617 (1971).

as modified by -- Steck, T.L. Cross-linking of the major proteins of the isolated erythrocyte membrane. J. Mol. Biol. <u>66</u>:295-305 (1972).

as modified by -- Laemmlie, U.K. Cleavage of structural problems during the assemblyof the head of bacteriophage T4. Nature 227:680 (1970).

STOCK SOLUTIONS:

1.	Conc A+B:	Acrylamide (20 g) (USE GLOVES) Bis (0.75 g)
		Water to 100 ml. Freeze in 150 ml lots (-20°C.

2. <u>Gel Buffer, pH 7.4 (x10)</u>: Tris 24.25 g Final conc = 0.4 M Na0H 4.0 g " " = 0.2 M as NaAc Na2EDTA 3.7 g " = 0.02 M HAc to pH 7.4 Water to 500 ml. Freeze in 10 ml lots.

- 3. Solubilizing Solution (x5): Tris 1.21 g Final conc = 0.1 M NazEDTA 0.37 g " = 0.01 M Sucrose 50.0 g " = 50% Titrate to pH 7.4 with 1M HCl Water to 100 ml. Freeze in 5 ml lots.
- 4. 20% SDS: SDS 100 g Water to 500 ml. Freeze in 100 ml lots.
- 5. <u>1.5% Ammonium Persulfate</u>: Amm. Per. 0.15 g Water to 10 ml. (USE GLASS TO PIPET) (make up fresh)
- 6. <u>0.5 TEMED</u>: TEMED 0.05 ml Water to 10 ml. Freeze in 5 ml lots. . (USE GLASS TO PIPET) (make up fresh)

7. <u>Coomassie Blue Stain</u>: MeUH 454 ml H2U 454 ml Glacial Acetic Acid 92 ml Coomassie Brilliant Blue 1.25 g Filter with Whatman No. 1 filter paper.

Fixative:	MeOH	500 ml	
	H2 Û	450	ml
	HAC	50	сc
	<u>Fixative</u> :	<u>Fixative</u> : MeOH H2O HAc	<u>Fixative</u> : Me0H 500 ml H20 450 HAC 50

- 9. Destain: HAC 10%
- B. Making Gel
 - 1. Set up plates (use medium spacers and 20 well comb).
 - 2. Make the following solution:

Stock solution	Amount	
A + B	12.0 cc	
X10	6.0 cc	
20% SDS	.6 cc	
Water	32.4 CC	
1.5% AP	6.0 cc	
.05% TEMED	3.0 cc	

3. Pour solution between plates, place comb, and wait for polymerization.

4. Remove comb and place gel into chamber.

5. Fill chamber top and bottom with electrophoresis buffer.

X10	200	ml
10% SDS	20	ml
Water to	2	1

- 6. Place 30 ug of protein in each well a. Solubilizing Solution X5 Sol. soln. = .8 ml 20 % SDS = .4 ml Dithiothreitol = ~ 25 mg
 - Pyronin Y = dash
 - b. To solubilize protein add: 33 ul Sol. soln. to 100 ul of sample Boil 3 minutes.

7. Run gel at \sim 75 volts until tracking dye reaches the bottom of the gels (takes about 2 hrs).

8. Put gel in fixative overnight shaking.

9. Stain with Coomassie Blue, shaking, 1 hr.

10. Destain with 10% HAc, shaking.

ATP: ADENOSINE -5'-TRIPHOSPHATE ASSAY

REF: Method Bucher, Th., Biochim. Biophys. Acta (Amst) 1, 129, (1947 a Denneman, H.A., ges. exp. Med. 134, 335) (1961).

REAGENTS:

Solutions I:

 0.151 gms Glycerat-3-phosphate (BM#127116) dissolved in 100 ml of Triethanolamine buffer (Sigma #665-5).

Solution stable for four months at 4°C.

Solution II:

 0.018 gms -Nicotinamid-adenin-dinucleotide NADH (BM#127345) in 10 ml of distilled water.

Solution stable for four weeks at 4°C.

- a) Glycerin-3Phosphat-dehydrogenase/Triosephosphat Isomerase (BM#127787).
 - b) 3-Phosphoglycerat-Kinase/Glycerinaldehyd-3-phosphate-Dehydrogenase (BM#108456).
- 0.6N Perchloric Acid: 86.1 ml of 70% Perchloric Acid diluted to 1000ml with distilled water in a volumetric flask.
- 5. ATP standard (Sigma #A-3127) 100 um/ml .0623 gms/ml of distilled water Aliquots stored in freezer for daily use. Dilute 25 ul of standard in 2.5 ml distilled water just prior to use.

SPECIMEN:

Since ATP is rapidly decomposed, specimen should be deproteinized as soon as possible.

DEPROTEINIZATION:

1. Pipet into a 15 ml centrifuge tube:

4.0ml 0.6N Perchloric Acid (ice cold) 1.0 ml blood and mix well

 Allow to stand at room temperature for 10 min. then centrifuge for 10 min at approximately 3000 rpm (Sorvall). Use supernatant for assay. Store clear supernatant at -80°C.

SPECTROPHOTOMETRIC MEASUREMENTS:

Wavelength: 340 nm Glass cuvette: 1 cm light path Temperature: 20-25°C

PROCEDURE:

Pipette into cuvette:

2.0 ml Solution I 0.2 ml Solution II 0.2 ml supernatant fluid, unknowns, diluted standards, etc.

Mix by inverting with parafilm cover. Read optical density E1. Ad 10 ul of BM suspension #127787 and 10 ul of BM suspension #108456.

Mix by inverting with parafilm cover, wait until reaction stops (approx. 10 min) and read optical density E2.

 $E_1 - E_2 = E_1$

:

CALCULATIONS:

239 x E 340 nm = ATP/100 ml <u>mg ATP/100 ml</u> -- gms% of hgb = uMoles ATP/gm hgb. 507 (MW of ATP)

2,3 - DIPHOSPHOCYLCERIC ACID ASSAY

REF: Sigma Technical Bulletin No. 35-UV

REAGENTS:

- A. -<u>DPHN, Preweighted Vial, Stock No. 340-101</u> Each vial contains 1.0 mg -DPNH Stable for over one year at room temperature while kept dark and dessicated.
- B. <u>Triethanolamine Buffer Solution, Stock No. 665-5</u> Triethanolamine Buffer containing Magnesium ions and EDTA. Stable 6 months or longer when stored at 0-5°C.
- C. <u>Adenosine 5' Triphosphate (ATP)</u> Vial containing 500 mg of Product No. A-5394 Sigma Grade ATP, Disodium Salt Store below 0°C. Stabel one year or longer. c'. Reconstitute vial with 5 ml Triethanolamine Buffer, Stock No. 655-5 (reagent B). Stable several months when stored frozen.
- D. <u>GAPD/PGK Enzymes, Stock No. 366-2</u> Glyceraldhyde Phosphate Dehydrogenase and Phosphoglyceric Phosphokinase suspension in Ammonium Sulfate. Mix by inverting before each use (DO NOT SHAKE). Stable at least 6 months when stored at 0-5°C.
- E. <u>Phosphoglycerate Mutase (PGM), Stock No. 665–3</u> Phosphoglycerate Mutase Enzyme suspension in Ammonium Sulfate. Mix by inverting before each use (DU NOT SHAKE). Stable at least 6 months when stored at 0–5°C.
- F. <u>Phosphoglycolic Acid, Stock No. 665-2</u> 50 mg preweighted vial. Stable when stored below 0°C. f'. Reconstitute vial with 5.0 ml water. Stable for several months when stored frozen.
- G. <u>Trichloroacetic Acid 82, Stock No. 665-8</u> Approximately 82 Trichloroacetic Acid (w/v). Stable for at least 6 months when stored at 0-5°C.
- H. <u>DPG Standard (BM#105821) 1 um/m1</u> 0.0083 gms/10 ml of distilled water 0.0083 gms/100 ml of distilled water Aliquots stored in freezer for daily use.

2,3 DIPHOSPHOCYLCERIC ACID ASSAY

SPECIMEN:

2,3-DPG appears to be stable in whole blood for at least 2 hours when stored in an ice bath. Trichloroacetic filtrates are stable up to 2 weeks when stored at $0-5^{\circ}C$.

DEPROTEINIZATION:

Pipett 1.0 ml of freshly drawn blood (heparin, ACD or CPD) into 3.0 ml of cold 8% Trichloroacetic Acid, Stock No. 655-8.

Shake vigorously for several seconds. Keep the mixture cold for an additional 5 minutes to assure complete protein precipitation.

Centrifuge 5 to 10 minutes at approximately 3000 rpm (Sorvall) to obtain a clear supernatant.

SPECTROPHOTOMETRIC MEASUREMENTS:

Wavelength: 340 mu Glass cuvette: 1 cm light path Temperature: 20-25°C

PROCEDURE:

- Into a 1 mg DPNG vial, Stock No. 340-101 pipette:
 8.0 ml Triethanolamine Buffer Solution (Stock No. 665-5).
 Cap and invert several times to dissolve the DPNH.
 Stable up to one week when stored at 0-5°C. (D0 NOT FREEZE).
- Into two cuvettes marked BLANK and TEST, pipette the following: (NOTE: ONLY ONE BLANK IS NEEDED FOR EACH SERIES OF TESTS)

Mixture from Step 1 ATP Solution (Reagent C) Water Protein Free Filtrate (From Step B under Preparation of Protein Free Filtrate)	BLANK 2.5 ml 0.1 ml 0.25 ml	TEST 2.5 ml 0.1 ml 0.25 ml
Mix by inversion	and add:	
GAPD/PGK Enzymes (Reagent D) PGM Enzyme (Reagent E)	0.02 ml 0.02 ml	0.02 ml 0.02 ml
Mix by inversion		

Wait approximately 5 minutes

- Read and record OD340 mu of both the BLANK and TEST. Cuvette using water as reference.
- 4. To each cuvette add: 0.1 ml phosphoglycolic acid (Reagent F). Mix by inversion. Let stand for 25-30 minutes (to allow reaction to go to completion).
- 5. Read and record OD340 mu of BLANK and TEST using water as reference.

CALCULATIONS:

- a) OD BLANK = OD of BLANK (step 3) OD of BLANK (step 5).
- b) OD TEST = OD of TEST (step 3) OD of TEST (step 5).
- c) Calculate OD by subtracting a) from b);

OD = (OD TEST) - (OD BLANK)

d) uMoles of 2,3-DPG = $(0D) \times (3.0)$ = (0D) $\times (7.7)$ per ml of blood (6.22) $\times (0.0625)$

The factors are explained as follows:

a) 3.0 represents the volume of liquid in the cuvette.

b) 6.22 is the millimolar extinction coefficient for -DPNH at 340 mu.
 c) 0.0625 represents the volume of original samplein ml in the reaction mixture.

NOTE

You can also use above answer to calculate 2,3-DPG levels on the basis of Packed Cells, or Hemoglobin as follows:

1)	uMoles of 2,3-DPG per ml Packed Cells :	(uMoles of 2,3-DPG per ml Whole Blood X (100)
		Hematocrit in percent
	Ex: If 2,3-DPG level blood and Hemato	is 2.2 uMoles per ml Whole crit is 45%, then:
	uMoles of 2,3-DPG = per ml of Cells	(2.2) (100) = 4.9

(uMoles of 2,3-DPG per ml Whole Blood X (100)

2) uMoles of 2,3-DPG per grams Hemoglobin =

grams Hemoglobin per 100 ml

Ex: If 2,3-DPG level is 2.2 uMoles per ml Whole blood and Hemoglobin level is 15 grams, then:

uMoles of 2,3-DPG = (2.2) (100) per grams Hemoglobin = 14.7

(15)

Tetramethylbenzidine (TMB) Method for Measuring Supernatant Hemoglobin

Reagents:

1. TMB Reagent, 1% (Aldrich Chemical Co.)

Dissolve 1 gram of TMB in 90 ml of glacial acetic acid and dilute up to 100 ml in a volumetric flask with distilled water. The reagent's reaction with hemoglobin will remain stable for at least <u>5 weeks</u> despite some darkening of color if stored at 4°C.

2. Hydrogen Peroxide, 1% (30%, Fisher #H-325)

Dilute 3.3 ml of 30% hydrogen peroxide to 100 ml in a volumetric flask with distilled water. Prepare daily.

3. Acid Diluent Solution, 10% Glacial Acetic Acid

A 10% by volume acid diluent solution is prepared by adding 100 ml of glacial acetic acid to 900 ml distilled water. Stability of reagent is at least one month at R.T.

4. Platelet-Poor Hemologlobin "Free" Plasma (PPHFP), Diluent for standards and samples

Prepare from fresh citrated whole blood clarified by centrifugation and frozen in aliquots at -20°C. Prior to use, thaw and centrifuge. In order to compensate for the hemoglobin in this plasma diluent, a plasma blank is included with every run. The absorbance reading of this blank is then subtracted from the standards and the diluted samples.

Preparation of Hemoglobin Standards:

- Packed red cells (hematocrit 70-80%) prepared from a unit of citrated whole blood should be lysed with distilled water (3 ml H₂0 per ml red cells). Stroma and other debris should be removed by centrifugation and the hemaglobin concentration of the lysate should be determined by the cyanmethemoglobin method (Protocol #9). Aliquots should be frozen at -20°C.
- Prepare a 200 mgs % standard from lysed cells:

Example:

W.B. Hgb of lysed cells = 11.1 gms%

$$\frac{11,000}{200} = x$$

or .181 10 ml distilled water 1.81 100 ml distilled water = 200 mg % solution

 Prepare a working standard of 100 mg% with equal amounts of PPHFP and the 200 mg% standard solution.

Conc. of Standard	Amt. of 100 mg% Soln.
5 mgz	1 ul
25 mg%	5 ul
50 mg ž	10 ul
75 mgZ	15 ul
100 mg%	20 ul

Plot 0.D. versus concentration.

Preparation of Samples for Assay:

Supernatant or plasma samples will be prepared by filling a sufficient number of 12 mm polypropylene tubes to provide enough plasma or supernatant for assay. The amount of blood required will, of course, be dependent upon the hematocrit of the blood.

Centrifuge the tubes in the Sorval at 5000 RPM for 20' at a temperature of 10° C.

Pipet the supernatant/plasma into another 12 mm polypropylene tube using a plastic pipette. Clear any entrained red cells from the sample by repeating the centrifugation, and pipetting the supernatant/plasma into a separate glass tube.

Spectrophotometric Measurements:

Instrument	:	Unicam spectrophotometer
Wavelength	:	515 nm
Glass cuvette	:	1 cm light path
Temperature	:	20-25°C

Assay Procedure:

- For each unknown, blank, plasma bank, and standards, measure 1 ml of 1% TMB into a test tube and add 20 ul of plasma or standard to each tube.
- Add 1 ml of 1% H2U2 solution to each tube. Mix and let stand at room temperature for 20 minutes.
- 3. Following incubation add 10 ml of 10% CH3COOH, mix well and wait 10 minutes before measuring absorbancy at 515 nm.
- 4. Calculate concentration of hemoglobin from standard curve.

pH Determination

Instrument:

Radiometer pH meter 4 Micro Electrode Unit Circulation Thermostat

Operation:

- 1. Be sure the electrode water jacket is full of circulating water and the thermometer is set for the desired temperature.
- Remove the glass electrode off its holder and press the suction cap against the plastic tip, with your thumb and suck distilled water through the capillary electrode until a few drops have passed the suction cap.
- 3. Fill the polyethylene tip and the capillary electrode with the buffer solution (7.381 pH).
- 4. Stop the suction and lift the polyethylene capillary out of the buffer solution.
- 5. Put the glass electrode back in its holder with the polyethylene tip immersed in the KCl pool of the calomet electrode (this is the measuring position).
- 6. Set the pH meter at READING.
- Rotate the BUFFER ADJUSTMENT knob of the pH meter until the meter reads the pH of the buffer solution.
- 8. Set the pH meter to STAND BY.
- Flush the electrode with distilled water, by sucking it through the capillary.
- 10. Suck air through the electrode for 3-5 seconds to dry its capillary.
- 11. Place the glass electrode in its resting position but with the polyethylene tip outside the beaker.
- 12. Mix the blood and insert the polyethylene tip into the tube or segment.
- 13. Press the suction cap against the plastic tip and suck blood into the capillary electrode until one drop is visible in the suction cap. Air bubbles must not be present in the polyethylene tip and in the glass capillary electrode.

- 14. Place the glass electrode in its measuring position, and read the pH when the meter has come to rest. Take down the result (the actual pH).
- 15. Immerse the polyethylene tip in a buffer solution for rinsing purposes. (Buffer solution is used for the rinsing because it prevents poisoning of the capillary electrode. The pH of the rinsing buffer should be between 7.3 and 7.4). Now create the suction and the rinsing buffer will clear the capillary electrode of the blood.
- 16. Rinse the glass capillary electrode with distilled water, refill with buffer of pH 7.38 and place the electrode in its resting position.

DETERMINATION OF RED CELL HEMOGLOBIN CONTENT

Instrument: Coleman Jr. Spectrophotometer or Unicam Spectrophotometer

The standard curve is set up by increasing dilutions of the Hycel Cyanmethemoglobin Standard. For the 5.0 ml volume employed in the unknown samples the undiluted standard corresponds to 20.0 Gm/dl hemoglobin. Dilutions are made to correspond to 15 Gm/dl, 10 Gm/dl, 5 Gm/dl and zero. Dilutions must be made with the Hycel Cyanmethemoglobin Reagent -- never with water.

Place five test tubes in a rack. Mark the tubes 20, 15, 10, 5 and 8. Dilutions must correspond to the following table:

Gm/dl Hemoglobin	20	15	10	5	Blank
Volume of Standard	6.0 ml	4.5 ml	3.0 ml	1.5 ml	none
Volume of Reagent	none	1.5 ml	3.0 ml	4.5 ml	6.0 ml

Transfer the dilutions to well matched cuvettes. Set the instrument to the proper wave length, 540 mu or filter. Adjust the instrument so the Blank tube has zero Optical Density or 100% Transmission. Take the readings for the standard and plot Optical Density on straight graph paper and % Transmission on semi-log graph paper.

Determining Unknown Sample:

Place 5.0 ml of Reagent in 19 x 150 mm cuvette test tube. Add exactly 0.02 ml of blood and mix contents. Read against Reagent Blank. Transfer reading to standard curve and obtain hemoglobin concentration in Gm/dl.

REF: Sigma Technical Bulletin No. 635

THE COLORIMETRIC DETERMINATION OF GLUCOSE

REAGENTS:

- A. <u>O. Toluidine Reagent, Stock No. 635-6</u> Store in dark at room temperature. Caution: Do not pipet by mouth. Avoid contact with skin. Do not inhale fumes.
- B. <u>Trichloroacetic Acid Solution, 3% (W/V), Stock No. 635-3</u> Store in refrigerator at 0-5°C.
- C. <u>Glucose Standard Solution, Stock No. 635-100</u> Standardized at 1.0 mg/ml (100 mg/100 ml or 5.5 m mol/1) with Benzoic Acid added as preservative Store in refrigerator at 0-5°C.
- D. <u>Hyland Q-Pak--Chemistry Control Serum 1</u> List No. 045-030 Store between 2° and 8°C before and after reconstitution.

SPECIMEN:

Glucose can be determined on whole blood, plasma, or serum. At room temperature, glucose in blood undergoes glycolysis at a rate of approximately 5% per hour.

DEPROTEINIZATION:

Into a centrifuge tube, pipet: 0.4 ml sample 3.6 ml of 3% Trichloroacetic Acid, Stock No. 635-3 Mix well by shaking. Allow to stand approximately 5 minutes to precipitate proteins. Centrifuge (5-10 minutes) until clear supernatant is obtained in the RC-3 at 3000 rpm. Store clear supernatant at 0-5°C, if determination is to be within a few hours. Otherwise freeze clear supernatant.

INSTRUMENT:

Practically any photoelectric colorimeter that transmits light in the range of 620-650 nm can be used.

Instrument: Coleman Jr. Wavelength: 635 nm Cuvette: 19 x 150 mm

> NOTE: Instrument readings should be made in terms of Absorbance (A) for use in calculations and for preparation of calibration curves.

PROCEDURE WITH DEPROTEINIZATION

For whole blood and markedly icteric or hemolyzed samples.

CAUTION: Avoid inhaling fumes resulting from the O-Toluidine Reagent. It is suggested that the procedure be carried out in a <u>well ventilated</u> area or hood. DO NOT pipet this reagent by mouth.

1. Label test tubes or cuvets BLANK, STANDARD, CONTROL, TEST 1, TEST 2, etc.

TO BLANK ADD:	TO STANDARD ADD:	TO CONTROL & TESTS ADD:
0.1 ml water 0.9 ml 3% TCA Solution Stock No. 635-3	0.1 ml Glucose Standard Solution Stock No. 635-100 0.9 ml 3% TCA Solution	1.0 ml supernatant
	JLUCK NU. DJJ J	

- To each tube add: CAUTION: D0 NOT PIPET BY MOUTH
 5.0 ml 0-Toluidine Reagent, Stock No. 635-6 Mix by lateral shaking.
- 3. Place all tubes in a vigorously boiling water bath for exactly 10 minutes. NOTE: If results are calculated from a standard that is assayed with the Test, the boiling period can be 10 ± 1 minute.
- Quickly remove all tubes and cool to room temperature by placing in tap water for approximately 3 minutes.
- 5. Transfer contents of tubes to cuvets and read Absorbance of STANDARD and Tests at 635 ± 15 nm, using BLANK as reference. Complete readings within 30 minutes.

CALCULATIONS:

Use of Standard: Glucose (mg/100 ml) = <u>*Test</u> *Standard × 1000*

*Represents the concentration of the Glucose Standard (mg/100 ml).

NOTE: If the reading of Test indicates a Glucose concentration greater than 250 mg/100 ml, dilute the test with an equal volume of 0-Toluidine Reagent, Stock No. 635-6. Read Absorbance of the diluted Test and multiply result by 2.

FLAME PHOTOMETRY

SODIUM AND POTASSIUM DETERMINATIONS

REF: Instrumentation Laboratory Manual 143

EQUIPMENT:

IL Model 143 Flame Photometer Compressed Air - Supply with a minimum pressure of 25#/in² site gauge showing clean oil and with a flow of moisture-free air of 0.5 ft³/min. Fuel: Propane, (instrument grade) Dade Dilutor Model 200

REAGENTS: (Instrumentation Laboratory, Inc.)

#35000 Stock Lithium Concentrate 1500 mEq Li/liter dilute 20 ml stock to 2,000 ml with distilled water #33203 150 mEq K/l #35050 50 mEq Na/liter/100 mEq K/liter #35100 100 mEq Na/liter/5 mEq K/liter #97517-50 140 mEq/Na/liter/2 mEq K/liter #97518-50 120 mEq/Na/liter/2 mEq K/liter #97519-50 160 mEq/Na/liter/8 mEq K/liter

OPERATION:

 Turn the front panel toggle switch to the ON position. This allows ignition to take place.

 Open the propane Shut-off valve located at the top of the fuel cylinder, one-half counter clockwise turn.

b. Upen laboratory air manifold fully.

- Turn ON/OFF switch to its upper position. The Flame On indicator light should come on and the faint sparkling sound cease within a matter of a few seconds.
- 3. Place a lithium blank solution or distilled water on the sample stand, raise the stand and aspirate for at least five minute before calibrating.

SAMPLE AND STANDARD PREPARATION:

 Lithium Diluent - It is imperative that the lithium concentration in all calibrating standards and unknown samples be kept constant. Dilute (20 ml) #3500 Stock Lithium Concentrate 1500 mEq Li/liter to 2,000 ml with distilled water.

2. Calibrating Standards - Once the lithium diluent is prepared, dilute working standard can be made. Dilute all calibrating standards and also all biological unknown TWO HUNDRED TIMES with the lithium diluent. This is done via the Dade Dilutor with 50 ul of sample and 9.95 ml of diluent. NOTE: It is recommended that approximately one ounce disposable plastic cups be used to hold standards and samples for analysis.

ANALYSIS OF DILUTED SAMPLES:

- As recommended aspirate lithium diluent through the atomizer for several minutes before calibration. Doing this allows the diluent to thoroughly circulate around and wet the walls of the atomizing chamber, thus preventing the collection of droplets on the atomizer walls and instability in the readout.
- 2. Change the lithium solution and aspirate a fresh lithium diluent.
- 3. Set the needle on the Lithium Response Meter at the center of the reference triagle with the Lithium Set Control.
- Set both the Sodium and Potassium Digital Concentration Displays to 0000 with their respective Zero Controls.
- 5. Aspirate an appropriate working standard 140 Na/5 K and switch the Potassium Range Selector Switch to the required position (20 for Serum).
- 6. Set both the Sodium and Potassium Digital Concentration Displays to the corresponding values of this standard with the respective Balance Controls. For example, the 140 Na/5 K working standard would be calibrated at 1400 and 0500, respectively.
- 7. Replace the standard with an unknown sample correctly diluted with lithium diluent and aspirate. Read and record the concentration of the unknown from the Sodium and Potassium Digital Concentration Display. NOTE: Initially, it is advisable to recheck the calibration of the instrument after every five to ten unknowns. After ten to fifteen minutes of use, the calibration, typically, need be check less frequently.

INSTRUMENT SHUT DOWN:

- 1. After completing your determinations, aspirate lithium diluent through the unit for several minutes. This serves to thoroughly wash the atomizer and burner system.
- Turn the propane Shut-Off Valve at the top of the fuel cylinder, clockwise (closed). Do not overtighten this valve, finger tight is adequte.
- 3. Close the air manifold fully.
- Place the ON/OFF Switch in its lower position.

MAKING GHOSTS IN RC 2B OR 5B SS34 ROTOR (All done at 4°C)

1. Solutions for making ghosts:

- a. <u>0.15 M NaPÚ4 Stock</u> titrate .015 M Na2HPÚ4 7H2Ú with .15 M NaH2PÚ4 H2Ú until pH = 7.4
- b. <u>PBS</u> 50 ml 0.15 M NaPÚ4 stock 8.77 g NaCl

.

- c. 5P8 33.3 ml of 0.15 M P04 can also make it 0.5 mM EGTA = 2.5 cc of 200 mM stock 1 l = total volume pH to 8 with Na0H
- 2. Spin down whole blood (go up to 12 K rpm, then cut out rotor).
- 3. Remove plasma.
- 4. Suspend and mix RBC's with PBS, spin down RBC (go up to 12 K rmp, then cut out rotor), aspirate off supernatant.
- 5. Repeat #4 2 more times.
- 6. Suspend and mix cells with 5P8, spin down ghosts at 20 K rpm for 10-15 min. Aspirate supernatant and any white cells at the bottom of the tube.
- Repeat #5 until ghosts are pure white (4-8 times depending on the amount of blood in the tube).

MAKING SPECTRIN EXTRACT

Unless otherwise indicated, done at 4°C in SS34 rotor in RC2B or RC5B centrifuge.

1. Make ghosts.

- 2. Suspend and mix ghosts in ice cold 0.1 mM PO4 (20 cc of 5P8/1).
- 3. Spin at 20 K rpm = 50 kg for 20 min, aspirate off supernatant.
- 4. Repeat #2 and 3.
- 5. Incubate ghosts, shaking at 37°C for .5 hr.
- Place ghosts into 60 TI rotor and spin for 1 hr at 39 K rpm = 110 kg in Sorvall or Beckman ultracentrifuge.
- 7. Aspirate and save supernatant = Spectrin extract.
- 8. Read $0D^{289}$ against H₂0 = Spectrin concentrate.

MAKING SPECTRIN DIMER (All done at 4°C)

1. Make spectrin from 50 cc of whole blood.

- 2. Concentrate spectrin to about 5 mg/ml.
- Make extract 10% sucrose.
- 4. Load extract on 12 x 1000 mm column (column gel = Biorad A-15).

A-15 Buffer 150 mm NaCl 10 mm Tris .1 mm EDTA .02% Na Azide

- 5. Run at 3 cc/hr, collect 20 min fractions.
- 6. Chromatograph.

.

SPECTRIN ACTIN 4.1 BINDING ASSAY

A. Solutions

Actin Buffer = Buffer A	
2 mM Tris-HCl pH 8	2 ml of 1 M stock
0.2 mM ATP	~0.111 g (variable M.W.)
0.5 mM DTT	77 mg
0.2 mM CaClz	0.2 ml of 1 M stock
•••	total vol = 1 liter

High Salt Buffer =	HS Buffer
1 M NaCl	5.84 g
59 M Tris	0.7139 g
(pH 7.4)	
2.7 mM DTT	0.0416 g
0.9 mM ATP	0.53 <u>6 g</u>
•••	total vol = 100 ml

[store at -20°C]

4.:	1 Bi	<u>iffer</u>	
10	۵M	PÚ4	
.5	mΜ	DTT	
			[pH = 7.4]

X10 Sucrose Buffer 1.5 M NaCl 40 mM PO4 2 mM MgCl2 5 mM DTT 100 mM Tris 1.0 mM ATP 50 mM KCl

```
[store at -20°C]
```

B. Polymerize Actin

1. Take Actin frozen at -20°C in glycerol and dialyze X3 in Buffer A.

2. Make Actin 50 mM KCl, 2 mM MgCl2.

3. Let sit 2 hrs at room temperature.

<u>S</u>

C. Mixes a la Wolfe

S-A-4.1 Mix F Actin 240 \ of 4 4.1 600 \ of 4 BSA 20 \ of 5 HS Buffer 160 \ 1020 \ 1020 \ 1000 \	OD280 = 2.0 in Buffer A OD280 = 0.2 in 4.1 Buffer 35% Sterile
--	---

<u>SA. Mix</u> F Actin A 1 Buffer	240 \	of	ÚD280	=	2.0	in	Buffer	A
BSA	20 \	of	35% S	tei	rile			
HS Buffer	$\frac{160}{1020}$							

<u>S Mix</u> F Actin Buffer 4.1 Buffer		
BSA	20 \ of 35% Ster	ile
HS Buffer	<u>160 \</u>	
	1020 \	

D. The Assay

<u>SA 4.1</u>

SA

1.	50 \ of SA 4.1 mix	50 \ of SA mix	50 \ of S mix
h	2-40 up 1257 spectrin	2-40 µg ¹²⁵ I spectrin	2-40 ug ¹²⁵ I spectrin

2.	dimer up to 30 \ in A-15 Buffer	dimer up to 30 \ in A-15 Buffer	dimer up to 30 \ in A-15 Buffer
3.	A-15 Buffer	A-15 Buffer	A-15 Buffer
	total vol = 80 \	total vol = 80 \	total vol = 80 \

- 4. Incubate at 23°C 60-90 min.
- Layer 60 \ of each sample above 5% sucrose in sucrose buffer leaving an air space in polyethylene tubes.

6. Centrifuge 2 hrs at 20 K rpm in SS34.

7. Freeze tubes and clip 7 mm from bottom. Count supernatant and pellet.

PROTEIN 4.1 PREPARATION BY TWEEN 20

All done at 4°C in SS34 rotor in RC-2B or RC-5B unless otherwise indicated.

- A. Make ghosts.
- B. Band 6 extraction (Glyceraldehyde-3-phosphate dehydrogenase)
 - 1. Stir membranes in: 0.155 M NaCl = 9 g/l 0.5 mM EGTA = 2.5 ml/l of 200 mM stock 5 mM Na phosphate, pH 8 for 30 min on ice, i.e. add salt and EGTA to ghosts.
 - 2. Centrifuge at 20 K rpm for 10 min.
 - 3. Wash with PBS and recentrifuge.
 - 4. Wash with glycine buffer (part C) and recentrifuge.

C. Spectrin and actin (Band 5) extraction

- 1. Suspend membranes in 6 vols of: (for unit bring vol to 1 1)
 5 mM glycine = 1.13 g/3 l
 0.5 mM EGTA = 7.5 ml of 200 mm stock/3 l
 5 mM beta-mercaptoethanol = 1.07 ml/3 l
 3 mM azide = .585 g/3 l
 pH to 9.5
- 2. Stir on ice 6-12 hrs and spin at 20 K rpm for 20 min.
- 3. Repeat extraction 2 more times.
- 4. Suspend membranes in 6 vols of: (for unit bring to 1 1) 0.6 mM Na Azide = 0.78 g/2 1 0.5 mM EGTA = 5 mls of 200 mm stock/2 1 pH to 3
- 5. Stir on ice 6-12 hrs and spin at 20 K rpm for 20 min.
- 6. Repeat Azide extraction once more and recentrifuge.

D. Tween 20 4.1 Extraction

1.	Dilute membrane with an	equal volume of:
	0.2 M glycine	7.5 g
	2% Tween 20	10 cc
	2 mM Na tetrathionate	.306 g
	1 mM EGTA	<u>2.5 ml</u> of 200 mm stock
	pH 9.8	500 ml

2. Stir on ice 10-15 hrs.

. _ . . _ ____

3. Spin at 39 K rpm for 1 hr in 60 Ti rotor in Sorval or Beckman ultracentrifuge.

- E. Purification at 4.1
 - 1. Layer sample on DE52 column (.15 ml DE52/mg protein) equilibrated in: 0.05 M glycine 3.75 g 0.5 mM EGTA 2.5 cc of 200 mm stock 0.5 mM dithiothreitol <u>.08 g</u> pH 9.8 1 liter

2. Wash with equilibration buffer until OD280 falls to 0.05 or lower using white-white tubing (36 cc/hr).

- 3. Apply gradient
- 4. Collect using black-black tubing (19cc/hr).
- 5. Collect 10 min fractions.
- 6. Run on SDS-PAGE gels to check purity.

METHOD FOR RBC LIPID EXTRACTION

DO EVERYTHING IN DUPLICATE

- Wash anticoagulated RBC three times in saline, remove WBC with cellulose acetate columns. Final WBC should be less than 400/mm³.
- 2. Extract approximately 0.2 cc of packed RBC and determine RBC count on Coulter.
- 3. Use disposable 16 x 150 mm glass tubes to do the extract.

.

- Add 0.5 cc distilled water to hemolyse samples. Vortex. After lysis is completed, place samples in ice bath and let cool for 5-10 minutes.
- 5. Add 5.5 cc of <u>cold</u> isopropanol to lysate and vortex. The isopropanol must be cold or the lipids will not be completely extracted from the RBC. After approximately 15 minutes in ice, let samples remain at room temperature with frequent vortexing for approximately 1 hour. The RBC membranes should look fluffy and will settle by this time.
- 6. Place the samples back in the cold ice bucket and add 3.5 cc <u>COLD</u> chloroform. Again vortex, remove from the cold in 15 minutes, and keep at room temperature for 1 hour. Frequent vortexing is necessary.
- 7. After 1 hour, centrifuge at ~1,5000 x g, filter supernatant with glass wool pipets into 16 x 150 mm tube, and wash pellet.
- 8. Remove inorganic phosphorus by adding 2 cc .5M KCl/10 cc extract. Vortex and spin down. Remove aqueous layer (top). Repeat 2 more times. If doing just cholesterol assay skip this step.
- 9. Dry down extract with nitrogen in 40-45°C bath.
- 10. Bring up to known volume (9.25 cc) with chloroform/methanol 2/1 and aliquot for assay.
- NOTE: All glassware must be phosphorus free for phosphorus assay. No rubber stoppers should be used for cholesterol assay. They are not recommended for phosphorus either.

May stop at 2 points:

- 1) after addition of isopropanol
- 2) after addition of chloroform and store at -20°C overnight.
METHOD FOR PHOSPHORUS DETERMINATION

 Reagents: Perchloric Acid 69-72% as purchased from J.T. Baker. Ammonium Molybdate 2.5%. Dissolve 2.5 g per 100 ml water.* Ascorbic Acid 10% (Purchased from Fisher). Dissolve 10.0 g per 100 ml water.

*Keep at 4°C in dark bottle, make 500 cc, storage okay for 1 month.

- 2. Tubes used: 16 x 150 mm disposable borosilicate glass culture tubes.
- 3. Specimen:

Plasma: Make lipid extract from 0.2 ml plasma or serum. Use 1.0 ml extract for phosphorous determination.

RBC: After washing in isotonic saline 3X and removing WBC's, use 1 cc at 20% Hct, spin and discard supernatant or use 0.3 cc red cells from 80% Hct to make lipid extract. Use 1.0 of the extract for phosphorous determination.

4. Procedure:

1. Evaporate all aliquoted extracts to dryness with N2 in 40-45°C H20 bath.

2. Add 0.5 ml perchloric acid to each.

3. Heat in 160°C heating block for 30 min; if samples are not clear, heat longer.

- 4. Let cool ~5 min.
- 5. Add 3.3 ml H20 to each.
- 6. Add 0.5 ml ammonium molybdate and vortex.
- 7. Add 0.5 ml ascorbic acid and vortex.
- 8. Boil in 100°C water bath for 10 min; put marbles on top of tubes.
- 9. After cooling ~10 min read at 797 nm.

5. Standard Curve:

Make .8650 g Na₂HPO₄ . 7H₂O stock solution = 100 ug P/cc. Make working stocks of 1 ug/ml, 2 ug/ml, 3 ug/ml from above stock. For daily standard curve, use 1 ug, 2 ug, 3 ug by drying 1 ml from each of the above working stocks in oven. These may be dried in a batch and stored at room temp indefinitely.

Alt: Pre-ashed phosphorus STD may be purchased.

NOTE: All glassware must be new or acid washed.

TOTAL CHOLESTEROL

Method of Zlatkis

REAGENTS

1. Standard Cholesterol solution (1 mg/ml).

a. Dissolve 100 mg pure, dry, ash-free cholesterol in 100 ml glacial acetic acid.

b. Alternative: Purchase Sigma's Chol. Standard Solution in Glacial Acetic Acid (1 mg/ml).

c. Daily set up fresh standard curve with standards of 50 ug, 75 ug, 100 ug.

- FeCl₃ · 6H₂O Solution: 10 g/100 ml HAc (#2). Make fresh daily.
- 3. Color Reagent (#3): Slowly dilute 1 ml of #2 to 100 ml with conc. H2S04.

PROCEDURE

- Use either 0.2 ml plasma or 0.2 cc packed RBC. Perform isop/chlor lipid extract as indicated.
 Use 3.0 ml of the extract for cholesterol determination. Use 16 x 150 mm borosilicate glass tubes. Evaporate extract to dryness with N2 in H20 bath at 40-45°C.
- 2. Add 3 ml HAc to sample tubes. Vortex 30 sec.
- Add 2.0 ml Color Reagent (#3), slowly down the side of the tube. This will form two layers. Vortex 1 min to mix <u>very</u> well.
- 4. Let tubes come to RT (30 min).
- 5. Read at 560 nm in same order as added Color Reagent.

FEST 5830



BECTON-DICKINSON

FRAGILITY DETERMINATION FOR MANUAL METHODS **RBC OSMOTIC**

INFORMATION AND PROCEDURE

PRODUCT INFORMATION

USE

tions equivalent to 0.85%, 0.65%, 0.60%, 0.55%, 0.50% 0.45%, 0.40%, 0.35%, 0.30%, and 0.00% sodium chloride. These 10 saline concentrations are contained in UNOPETTE Reservoirs correquindingly labeled 85. reagent system for the determination of erythrocyte UNOPETTE Test 5830 is a stable in vitra diagnostic ournotic fragility, using 10 buffered satine concentra 55. 60, 55, 50, 45, 40, 35, 30, and 00

SUMMARY

buffering system to maintain pH stability.³ Thunerosal The method is based on an adaptation of the original ery throcyte osumtic. In such procedure using various Dacie," and inviditied by the addition of a phosphate concentrations of sodium chlorale, as described by s added to inhibit bucieruit growth

PRINCIPLE

Whole blood is wided to the buffrind subne diluents which have concentrations equivalent to 0.00% to 0.85% sodium chloride.

lent to isotonic (0.85°°) safine will retain their shape and Erythrocytes suspended in the buffered diluent equiva become spheroid and more fragile, and eventually lyse, will not hemolyze. Erythrocytes suspended in the hypotonic saline dilucints will take up water, swill, releasing humoglobin into the solution.

Characteristically abriormal curves will be generated by The percent hemolysis that occurs within a 20 minute incubation period in each concentration is determined strengths of the soline diluents. In normal subjects, an almost symmetrical signoid shaped curve is obtained. different erythrocyte dyscrasias. (See PROCEDURE, formula and plotted on a graph against the osmotic by absorbance measurements utilizing appropriate Step 9.1

Results obtained with UNOPETTE Test 5830 compare favorably with those obtained using Dace's method with sodium chiloride alone in the diluent, ^{3,4}

REAGENTS

 UNOPETTE Reversions containing 3.98 mLot dilurnt misture

0 020 gm 0 100 gm U 140 gm us to 1 liter Dibase antivotions sortium phosphate Dibase anti-trous sodium phosphate Alonobasic sodium photohale Memohasic sortium phosphate Reservoir Labelist 00 Reservoir Labelad 30 Distributivities Thunerosat

0 407 am 0 057 am 0 100 am[•] 2 700 am us to 1 hier 0 4 / 4 gm Ditasic antivitious sodium photohate Atomotosic sodium phosphate Hereron Calieled 35 Software Interaction United voter Thimeros.d

0 162 gm 0 100 gm 7 650 gm qs to 1 liter ...1.152 gm

> Sortium chilorate **Uistillet water** hunerosal

Dibuse antigenes sedium phosphare Monobuse sedium phoquhute

Rewrynie Latieliei 84

IN UNOPETTE Capitlary Pyneties ... 20 µl capacity

N

0.100 gm 0.100 gm 3.150 ym qs to 1 liter Sodium chloriste Distribution water Thimerosal

Herervoir Lideful 40

Dubates with stores at a second store at a second store sto
Alimit and a satium phosphale U.542 gm
www.combec.scy.houn phosphale
f humerosul
Surfum chlorade 3 600 gm
Untilled water and an and a set of the
Heservoir Lubeluri 45.
Othistic antivities softium phosphate 0.610
Monthuse subjects about its 2000
0 100 gm
Utstituted water
Reservour Listuelert 50:
Dahasic antiyeleous sodiumi phosphale 0.678 am
filonobasic sodium phosphale
I hinterosal 0.100 gm
Sotium chloride
Distilled water
Reservoir Laheled 55:
Dibasic anhydrous sodium phosphate 0.745 am
Monobasic sodium ohomhale
Sodum chlorida
Desition (1990 gm
Unstituted water and a second se
Reservoir Lubeleri 60
Didiate anhydrous sodium phosphate 0.813 pm
f-lonobasic sodium phosphate
Thimerosal 0.100 am
Sodium chloride
Distilled water
Reservoir Labeleri 65:
Ditusic anhydrous sodium phosphate 0 881 am
Monobasic sodium phosphate 0.124 am
Thimerosal 0.100 am
Sotium chlorwie 5 850 an
Distilled water

ORIGINAL PAGE IS OF POOR QUALITY TI

F

7

[7]

Π

17

1

,

•



:



- pressure will thaw blood reservoir Then remove Release pressure on Inder from pipette papette ser orely an openning Regarive 4 YOU HOADEN
- but not out at, overllow Squerze reservoir grintly forcing diluent up into, Iwo or three times to pressure each time to rinse capillary bore. churcher, releasing return misture to into diluent

Ξ

- CAUTION See WAHN ING Section Place unlex fuger over invert several times to ILLIA IDS BI
- Approximiting the second second thoroughly mix bloot with diluent. đ

3

Repeat PROCEDURE 2 (a.g) for all two labored reservoirs using a new paterte for each didution

3. INCUBATE

Incubate diduted sumples at room temperature for twenty (20) minutes.

4. TRANSFER CONTENTS

- **Transfer thor**oughly inixed contents of each reser voir to appropriately labeled test tube is follows Convert to thropper assembly by withdiswarut
 - pipette from reservoir und revuit securely in reverse position.

•

 Place capitlary tup into approve code, based test tube and squeeze reservor to e-tail entre Contents

5. CENTRIFUGE

Centrifuge all samples at 2000 gran for their full mmules.

6. TRANSFER CONTENTS

tabeled curvettes which accommodute 4.0 ml of Hund Be careful not to decant centrifuged cell button supernate from test tubes into correspondently. Keeping surples in proper surprises of property

7. SET WAVELÈNGTH

nin or use a litter photometer with appropriate filter Set the wavelength on a spectrophy output to fold

B. READ ABSORBANCE

ċ. ð. 3 2

- Wipe outside of cuvettes liefond physical in ŗ
 - instrument. •
- Check zero stability of instrumment using coverte fulled with distribution water to act a Blank ف
 - Aleasure absorbance of all trim (10). Unknown concentrations against Water Blank set at 2000 absorbance, and record ú

9. CALCULATE

(O.D.)) for the desired concentration in the follow Determine percent housing sis in our transmission. substituting absorbance value [equit_defends] ing formula

- X IIII 0.0.9 % 00. 00000 • Hemolysis 0 D 0.
- osmotic strength of sectoric chipade 1 e 0.65°+, 0.60°+, 0.55°+, 0.50°-0.45°-0.40 O.D. Absorbance of solution of an aven 0.35 -, 0.30

#

-

17

7

1

"1

_ 1

*

--

-,

EXAMPLE:

. С.,

51M) (1	4 Menudy	
	- 2-	
	HH.	
	్	<u>ب</u>
	000	000
.,# 0	0.248	0 305
¥.	فالمستحداث ويتو	

÷

10. CONSTRUCT ERYTHROCYTE FRAGILITY GRAPH

Protyner entilled of your state structure approximation of the development of the developments on the orderate. Enst point on the ordinate must be O 85 -

Typical Erythinscyte Fragility Graph with Normal Limits Defined and Two (2) Dyscrasias Illustrated¹



ORIGINAL PAGE 1\$ OF POOR QUALITY

:

LIMITATIONS OF PROCEDURE

and resugnered both sumples so that the ratio of cells to solve is the same for both 1. Then proceed with It specimen of tow hemoglobin content is to be treated. sumple in isotomic saline. Centrifuge, decant subernativ wigh erythrocytes in the specimen and in the control analysis.

EXPECTED NORMAL VALUES

		1111			1.6.06	100.05	5-45	2	•	=	=	-
Complete Strength of	Sodium Chloride Salution	0.00		. 1)E O	0 35'.	040	0.45%	0.50 0	0 55'+	0.60 ,	0.65	0 85
	Latiolat		8	20	15	40	45	50	55	60	65	85

EXPECTED PERFORMANCE

. .

Precision of Values of Erythrocyte Fragility Tration ŝ E 3

ing UNOPETTE To 1 5630.1	trength of the strands for ute Solution in ange	1001	1411 7 6.	2010		81 6-	-	=	-	-	-
vi (n. 10) u	Sodium	0 00:	0 30 -	0 35 0	0 10	0.45	050	0.55.0	0 CO °	0 65 🖕	
Dold Luni	beled	8	90	35	40	45	50	55	09	65	;

TECHNICAL NOTES

 It is advisable to currantematible and a control adve-regime an unknown somple for explore syneric. fragility

BECTON-DICKINSON

В О

•••••

1.00

11,101,114-1

Ξ

٢Į

}

Π

[]

ŧ

 \mathcal{F} . A testic sumple is generally developed for this proa web it show and show and owned to be a second change with time

.

- within the work area, reservoirs should be on ubited Respondence II. One to choose in temperature 3.051 in a room temperature water buth for consitent results.² for diation percet doubt be conducted at contra-_
- blood is modulish in erro under sterile combinats af 4. To disconcertain subile abiomatices or explore eyre osmotic feapfity, defidiomative or hep-neived 37. C. for 24 hours. Consult Reference 7 for
- $b_{\rm c}$. If we determine the constraint $0.35_{\rm c}$, i.e. we do not a number of the determination of each lie. munited and a first strength ter homense

. A DESCRIPTION OF A DESC 111 11 11 11 11 1 100 C 100 C 100 C . • ÷ A more principle sector. 1 1 1 1 PH 1 1 1 1 1 1 1 REFERENCES . ••• 1 I 1.1 -----1.1 .

14.2 Platelet Studies

Sample Aliquoting and Designated Responsibility for the Platelet Study.

Persons Involved:

D. Ausprunk (D.A.); W. Curby (W.C.); M. Jacobson (M.J.); I. Szymanski (I.S.); A. England (A.E.); Lisa Rice (L.R.); B. Kim (B.K.); D. Kenney(D.K.); F. Chao (F.C.).

Tests	Responsible Person	Blood Volume(ml	Performance) Site	Reference Section in Original Proposal
pH viscosity	M.J.	1	on site	4.2.2.1
platelet count volume size shape	W.C.	1	on site	4.2.2.1 4.2.2.2
morphology & shape change	D.A.	7	on site	4.2.2.2 4.2.2.3
platelet aggregati & ATP_release	on B.K. & A.	E. 4	on site	4.2.2.4 4.2.2.7
Thromboxane generation	L.R.	2	on site & in Boston	4.2.2.5
uptake & release of 5 HT	F.C.	2	on site & in Boston	4.2.2.6
&-thromboglobulin release	A.E.	2	on site & in Boston	4.2.2.7
PF-3 measurement	F.C. & L.	R. 3	on site	4.2.2.8
y Høpotonic Stress	В.К.	2	on site	4.2.2.9 .
membrane glycoproteins	D.K. & F.	C. 15	on site & in Boston	4.2.2.10
cytoskeleton	D.K. & F.	.C. 15	on site & in Bos	ton 4.2.2.11
Immunology	I.S.	3	on site	

X

Assays of Platelet Aggregations and Release Reaction

A. Working Space

50" x 25" bench space 4 plugs of the power line with 115 V AC, 60 Hz

- B. Equipment and Accessories:
 - 1. Lumi-Aggregometer (Model 1020, Payton Scientific, Inc., Buffalo, NY)
 - Two Dual Channel Strip Chart Recorders (Model 500, Payton Sci., Inc., Buffalo, NY).
 - 3. Chart Paper (#0100-025, Payton Sci., Inc., Buffalo, NY); one roll
 - 4. Cuvettes, disposable (CU312, Payton Sci., Inc., Buffalo, NY); two boxes
 - 5. Stirring Bar, (SB312, Payton Sci., Buffalo, NY); 100 bars
 - 6. Pipettes, 1 ml, 200 ul, 10 ul capacities; two each
 - 7. Pipette tips, blue and yellow

C. Reagents:

- 1. ADP stock solution, 10^{-3} M in 0.87% NaCl, 10 ml in total
- 2. Collagen stock solution, 20 mg %, 5 ml in total
- 3. ATP standard solution, 10^{-3} M, 3 ml in total
- 4. Luciferase-Luciferin Reagent (Part #395, Chrono-Log Co.); four vials
- D. Assay Procedures:

 Preparation of samples: Platelet concentrate of 3 ml aliquot sampled from each storage unit will be adjusted to platelet count 300,000 per ul by diluting with pooled citrate plasma. For this, 60 ml of pooled plasma will be obtained of which 10 ml will be used for control samples (before storage; to) and 50 ml for stored samples. The plasma for stored samples will be kept at 4°C until used.

- 2. Simultaneous assays for aggregation and release reaction:
 - Platelet samples of 0.4 ml plus the luciferase-luciferin reagent of 0.05 ml will be prewarmed to 37° C in the incubator of lumiaggregometer. The recording chart will be adjusted to a full scale with the sample and plasma as 0° and 100 $^{\circ}$ light transmittance, respectively. The assay will be started by adding aggregating agent of 0.05 ml into stirring sample. Aggregating pattern and the intensity of luciferin light produced by the released ATP will be recorded in a dual channel recorder. At the end of reaction (usually 4-5 mins), 2 ul of ATP standard solution containing 2 nmol will be injected and light intensity will be recorded and used as an internal standard for calculation. Two different concentrations of both aggregating agents, ADP (10 uM and 100 uM) and collagen (6.25 ug and 12.5 ug per assay) will be tested for each sample. Maximal degress of aggregation will be recorded and the value of release reaction will be expressed in umol ATP released from the aggregating platelets, e.g. per 10^{11} cells.

PRECEDING PAGE BLANK NOT FILMED

IBSE Protocol Page 54

UPTAKE AND RELEASE OF SERVIONIN

Equipment:

Microfuge	pH meter
Vortex	Scintillation Counter
Water Bath 37° C	Cell counters
Pipettes	Heater

Materials:

Distilled H20	pooled plasma; platelet concentrates
NaCl	TES
EDTA	HCl
Horm collagen	14C - 5HT
Protozol	NaOH
plastic tubes for incubation	Scintillation fluid
plastic vials for counting	microfuge tube, polycarbonate outer tubes

Plasma preparation:

Dilute platelet concentrates with pooled plasma to 300,000 cells/ul.

Solutions:

1.2% paraformaldehyde solution (PFA)

- Dissolve 6 gm of PFA in 100 ml of distilled H20; heat to 60-70°C; add 5N Na0H until solution becomes clear.
- Dilute (1 part) 6% formalin solution with 5 parts of TES-saline-EDTA (0.01 M TES, 0.14 M NaCl, 5 mM EDTA; pH 7.4).
- 3. Adjust pH to 7.4 if necessary by adding HCl.

Collagen solution (Horm Collagen)

Dilute collagen stock (100 ul) with 400 ul buffer (conc. 200 ug/ml).

Procedure:

- Mix PRP (2 ml, 300,000 cells/ul) with ¹⁴C-5HT (final conc. 5 uM) and incubate at 37°C.
- At 1, 3, 5, 10 and 30 min. after incubation, remove 150 ul of the reaction mixture and put in microfuge tubes (hard tube); centrifuge at 12000 rpm for 2 min. in a microfuge.
- Remove 50 ul (x2) of supernatant plasma and add 0.2 ml of protozol. Incubate overnight. Add scintillation fluid and count.

- 4. Remove 50 ul of <u>PRP</u> in duplicates and mix with 0.2 ml of protozol to determine total radioactivity.
- To the remainder of PRP, add 10 ul of diluted collagen (final conc. 2 ug/ml), mixed by vortex.
- 6. At 1, 3, 5, 10 and 15 min after addition of collagen, 100 ul of the reactive mixture will be removed and added to 400 ul of 1.2% paraformaldehyde solution.
- 7. Centrifuge at 12000 rmp for 2 min. in a microfuge, remove 50 ul of supernatant, mix with 0.2 ml of protozol overnight and 5 ml of Liquiscint solution (National Diagnostics; Liquiscint will be acidified by mixing 10 ul of glacial acetic acid with 1.5 ml of Liquiscint prior to use). The samples will be counted for radioactivity in a liquid scintillation counter using an established counting program.

•

SAMPLE PREPARATIONS FOR β -TG AND TXB MEASUREMENT BY RIA

- Gently add 2.5 ml of platelet concentrates to the blood collection tubes (which contain anticoagulent and antiplatelet agents) provided with the radioimmunoassay kit. The tubes are pre-marked at 2.5 ml.
- 2. Stopper the tubes and mix the content by gentle inversion two or three times.

.

- Cool the sample by placing the blood collecting tubes in cooling bath for 15 min.
- 4. Centrifuge the sample at 1500-2000 g at 2-4°C for 30 min.
- 5. Remove the top 0.5 ml of plasma using a pipette with a disposable plastic tip.
- 6. Transfer the sample to a separate labeled specimen tube.
- 7. Store the sample at -20°C and transport the sample back to Boston for assay.

Storage conditions:	24 hrs at room temperature
	7 days at 2-4°C
	4 weeks at −20°C

8. See attached package inserts.

INITIAL BLOOD STORAGE EXPERIMENT - PLATELETS PF3-1

PF-3 Measurements

Equipment:

Microfuge		Wire loop
Conculation assay wat	er bath	Ice buckets
Stop watches		Cell counter

Materials:

Platelet concentrates	Distilled water
glass clotting tubes	CaCl2, TES, NaCl
Russell's vipor venom	Pipette
Microfuge tube (hard tubes)	Ice
PPP prepared from platelet concentra	ates

Solutions:

- 1. Russell's vipor venom, working stock 10 ug/ml in TES-saline.
- 2. CaCl2. 0.05 M CaCl2 in H20.

Platelet and plasma preparations:

- Centrifuge 2 ml of platelet concentrate at 12,000 rmp at room temperature for 2 min to obtain PPP.
- 2. Save 0.5 ml of PPP for PF-3 measurement.
- Dilute platelet concentrate with PPP to obtain PRP at cell count of 300,000/ul.
- 4. During preparation of PPP and PRP, the tubes will be capped.

Procedures:

- Prewarm CaCl2 solution at 37° C. Store PRP and PPP at room temperature. Keep RVV on ice.
- Add 0.1 ml of PRP or PPP into a glass clotting tube.
 Add 0.1 ml of CaCl2, incubate at 37° C for 1 min.
- 3. Add 0.1 ml of RVV. Start the stop watch.
- Mix the reaction mixture constantly with a wire loop until a film clot forms over the loop; stop the watch; record the time.

Assay of Platelet Response to Hypotonic Stress (PRHS)

Working Space: Α.

> 65" x 25" bench space 2 plugs of the power with 115 V AC, 60 Hz

- B. Equipment and Accessories
 - 1. Spectrophotometer (Beckman DU, Model #2400)
 - Strip Chart Recorder (Omni-Scribe, Model #B-5000)
 - 3. Chart Paper (Cat. #13-939-35, Fisher Sci. Co.); one roll
 - 4. Cuvettes, semi-micro uv cell, 1 ml capacity; one set of four
 - 5. Pipettes; two 1 ml and one 0.2 ml capacity
 - Pipette tips: blue and yellow 6.
 - Other facilities for platelet countings (Coulter count, sample diluting 7. vials and Isotonic diluent, etc.)
- С. Reagents: NONE
- Assay Procedures: D.
 - 1. Preparation of samples: Platelet concentrate of 2 ml aliquot sampled from each storage unit will be adjusted to platelet count 400,000 per ul by diluting with pooled donor plasma. The pooled plasma of 50 ml will be needed, 5 ml for the to sample and 45 ml for the stored samples.
 - 2. Assay:

Platelet sample of 0.5 ml will be admixed rapidly with an equal volume of distilled water in the cuvette and the change of optical density of the platelets will be recorded in a spectrophotometer at the light wave length of 420 nM. A sudden drop in optical density occurs by exposing the platelets to a hypotonic condition which is followed by a gradual return to the orginal density. The latter phase, an increase in optical density units during the first 2 minutes, will be the value of the PRHS test.

Introduction of ¹²⁵I into Tyrosines on Platelet Membrane Proteins and Glycoproteins by Lactoperoxidase Catalyzed Radioiodination

Literature: Modified from Phillips, D.R., <u>Biochem</u>. <u>11</u>:4582-4588, 1972; Phillips, D.R. & Agin, P.P., <u>J. Biol</u>. <u>Chem</u>. <u>252</u>:2121-26, 1977; Holihan, J.R. & White, G.C., <u>Blood</u> <u>57</u>:174-181, 1981.

REAGENTS:

- H202 -Fisher Cat. #H-325; 30% stock (8.8 M) diluted 1:1000 into cold (4°C) PBS immediately before use. Stock is kept refrigerated and discarded approximately 1 month after opening.
- Lactoperoxidase Cal. Biochem Cat. #427488; dissolve in PBS pH 7.4 at 200 ug or .2 U/ml store frozen at -20°C in 600 ul aliquots. Molecular weight 80,000.
- 3. Na¹²⁵I Amersham Cat. #IMS-30; purchased in 1 mCi quantities takes 1-2 days for delivery after order; store at room temp. in radioactive area.
- 4. 10 mM KI prepared as a stock, kept frozen in aliquots.
- 5. Solution of TES/Tyrode's with 5 mM EDTA + glucose and 200 uM KI.
- 6. Solution of TES/Tyrode's with 5 mM EDTA + glucose and 100 uM KI.
- 7. Solution of TES/Tyrode's with 5 mM EDTA, no glucose, no KI.
- Solution of 1% Nonidet P40 in TES/Tyrode's with 5 mM EDTA (no glucose), 5 mM Iodoacetamide, .2 mM PMSF.

EQUIPMENT:

Hood equipped for radioiodinationVortexIntermediate Speed CentrifugepH MeterFacilities for disposal of 125I wasteCell Counter

PROCEDURE FOR 5 ml LABELING-REACTION MIXTURE:

- 1. Count washed platelets and adjust concentration to 0.5 x $10^{9}/ml$ (range 3-7 x $10^{9}/ml$) in 5 ml.
- 2. Bring solution to 10 uM KI (5 ul of 10 mM/5 ml).
- 2.5 ul of NaI¹²⁵ (or 250 uCi).
- 4. Add 1.25 nMoles LPO (500 ul of frozen stock).
- 5. H202 diluted in ice-cold PBS just before use.

- 6. Add H202 at 30 sec. intervals in 5,6 ul aliquots. Platelet suspension is stirred constantly on a magnetic stirrer with Fisher microbar (Cat. No. 14511-69) magnet.
- 7. Dilute to 35 ml with TES/Tyrodes, 5 mM EDTA, + 200 uM KI.
- 8. Pellet in IEC centrifuge at 2,000 rpm for 10 min at room temp.
- 9. Resuspend pellet in 15 ml TES/Tyrodes with 5 mM EDTA + 100 uM KI.
- 10. Pellet in IEC at 2,000 rpm for 10 min at room temp.
- Resuspend in TES/Tyrodes/EDTA (no glucose) at a final concentration of 1.0 x 10⁹/ml (cells must be counted).
- 12. Extract platelets for 5 min at room temp. with 0.5% NP40 by adding an equal volume of TES/Tyrodes + 5 mM EGTA containing 5 mM iodoacetamide and 1% NP40.
- 13. Freeze at -20°C for electrophoretic analysis.

Analysis of 125I-labeled Membrane Proteins (Boston)

- A. One dimensional SDS Polyacrylamide Gel Electrophoresis (SDS-PAGE) will allow detection of gross changes (i.e., cleavage, loss from membrane) in major platelet surface glycoproteins, particularly gpIIb and gpIIIa. Samples will be analyzed under both reducing and non-reducing conditions to detect changes in S-S bonding and SH-groups (and effect of microgravity) of membrane proteins that occur during storage. SDS-PAGE will be carried out using the discontinuous Laemmli System (Laemmli, Nature 227:680-85, 1970).
 - Solubilize extracts (prepared in #12 above containing 1 mg platelet protein/ml) with 1/10 volume 10% SDS, boil for 2 min at 100°C, mix with Laemmli Sample Buffer containing 1.8% SDS, 125 M Tris-HCl pH 6.8, 20% glycerol, bromophenol blue, and for analysis of proteins under reducing conditions, 5% beta-mercaptoethanol.
 - Apply aliquots of solubilized samples containing 50-120 ug of protein to precast SDS-slab gels consisting of a 4% stacking gel and 7.5% resolving gel.
 - 3. Electrophorese at 24 mAmps/slab until bromophenol blue tracking dye reaches the leading edge of the gel. Run off for 45 min additional minutes at 24 mAmps to improve resolution of components in the 100-200,000 molecular weight range.
 - *4. Stain gels in 0.025% Coomassie Brilliant Blue for 2 hrs.
 - 5. Destain in 2-3 changes of 30% methanol/10% acetic acid overnight.
 - 6. Photograph wet gel using an orange filter to record profile of total platelet polypeptides.

- 7. Dry down gel on vaccum slab gel dryer.
- Detect radiolabeled proteins by autoradiography using X-OMAT (Kodak) film and intensifying screens. 3-5 days exposure at -80°C will give bands of sufficient intensity.
- Quantitation of radiolabeled proteins will be performed by densitometric scanning of autoradiograms using a Hoeffer gel scanner.

*Duplicate gels will be stained by the Alcian Blue method which selectively stains glycoproteins (Wardi and Mechow, Anal. Biochem 49:607-09, 1972).

1. Fix gels in 12.5% Trichloroacetic acid for 30 min.

- Rinse gels 5X with distilled water then oxidize glycoproteins with 1% periodic acid.
- Rinse 2X, with changes 30 min each distilled water.
- 4. Soak gels in 0.5% potassium metabisulfite for 30 min.
- 5. Rinse gels in distilled water.
- 6. Stain overnight in 1% Alcian Blue in 3% acetic acid.
- 7. Destain with 7% acetic acid.
- 8. Photograph wet gels.
- 9. Dry down and subject to autoradiography as described above.

B. Two-Dimensional Electrophoretic Analysis of ¹²⁵I-labeled surface proteins involving separation of proteins first on the basis of charge by Isoelectric Focusing then on the basis of molecular weight by SDS-PAGE. The method used is that of O'Farrell (J. Biol. Chem 250:4007, 1975) as modified by Linck et al. (Cell Motility 1:127, 1982). The analyses will permit detection of subtle changes in charge distribution of membrane proteins.

- NP-40 extracts of ¹²⁵I-labeled platelets are prepared for isoelectric focusing by adding solid urea (ultrapure), Triton-X-100, dithiothreitol to final concentrations of 9.5 M, 3.2% and 100 mM respectively.
- 2. Isoelectric focusing gels, cast in 200 ul glass capillary pipettes, consisting of 3.5% polyacrylamide, 3.2% Triton-X-100, 4% ampholines (pH range 5-8), 0.4% ampholines (pH range 4-6), 0.6% ampholines (pH range 3.5-10), are prefocused at 200 V for 15 min, 300 V and 400 V each for 30 min. Electrode solutions are 0.02 M Na0H at the cathode and 0.01 M H3P04 at the anode.
- 3. Samples containing 30-60 ug total platelet proteins are applied to the basic end of the isoelectric focusing gels and focused for 16 hrs at 400 V and 2 hrs at 1000 volts (8400 volt hours).

- 4. The isoelectric points of relevant polypeptides are determined from duplicate IEF gels electrophoresed in parallel without protein. The blank gradient gel is cut into 5 mm length and eluted with 0.5 ml deionized water; the pH of the eluted ampholines is measured with a pH meter. The pH values of eluted ampholines are plotted with respect to gel length. The isoelectric points are then correlated with the stained bands on another duplicate gel after correction for expansion of the stained/destained gel.
- 5. For SDS-PAGE, the IEF gel containing fractionated proteins is equilibrated for 20 min at room temperature in Laemmli sample buffer containing 5% betamercaptoethanol and placed across a precast slab gel consisting of a 4% stacking gel and a 7.5% resolving gel.
- 6. Electrophoresis of the 2nd dimension SDS-PAGE is as described above for one dimensional SDS-PAGE.
- Gels are stained with Coomassie Blue, destained and labeled proteins detected by autoradiography as described above.

3H-Labeling of Sialic Acid Residues of Platelet Membrane Glycoproteins

It is necessary to examine possible storage-induced changes in platelet surface membrane proteins (and the effects of microgravity on these changes) by both the ¹²⁵I and ³H methods since different membrane proteins become radiolabeled by these techniques. For example, platelet membrane glycoprotein Ib, the vWF receptor, which is lost from the platelet membrane during storage, labels well with ³H-NaBH4 technique, but does not incorporate significant ¹²⁵I or stain with Coomassie Blue.

Literature:

Steiner et al., Thrombo. Res. 29: 43, 1983.

Reagents:

Modified Tyrode's Buffer containing 5 mM EDTA buffered with Hepes at pH 6.8 and 7.6 EDTA: .25 M stock [³H] NaBH4 (5 mCi in 50 ul cold 0.01 M Na0H) HEPES: 1 M stock iodoacetamide - .5 M frozen stock NP-40: 10% stock sol. (w/w) sodium metaperiodate - prepared fresh as 4 mM Na0H: .1 N stock solution in Hepes/Tyrodes EDTA at pH 6.8

Equipment:

Intermediate Speed Centrifuge	Ice Machine pH Meter
Fume hood Barrel for disposal of ³ H contaminated waste	Vortex

Procedure:

- Wash platelets twice by differential centrifugation in a modified Tyrodes' buffer containing 5 mM EDTA and buffered at pH 6.8 with 10 mM HEPES.
- Resuspend washed platelets in 1-2 ml of the same buffer.
- Adjust platelet count to 0.5 x 10⁹/ml. Equilibrate at 4°C in the dark for 20 mins.
- 4. Add freshly prepared 4 mM sodium metaperiodate to the cooled platelet suspension to a final concentration of 2 mM and allow the reaction mixture to stand in the dark for 10 mins.
- Dilute platelet suspension to 15 ml with HEPES/Tyrode's/EDTA at pH 6.8 and pellet.
- 6. Wash platelets twice with HEPES/Tyrode's/EDTA at pH 7.6 and resuspend in 1-2 ml of the same buffer; adjust platelet count to 0.5 x 10⁹/ml.

- Add [³H] NaBH4 (final concentration 500 uCi/0.5 x 10⁹ platelets) to the 7. washed platelet suspension and incubate with periodic agitation, for 30 mins at room temperature.
- Wash radio-labeled platelets twice with HEPES/Tyrode's/EDTA pH 6.8 and 8. resuspend in the same buffer.
- Extract platelets with 0.5% NP-40 in the presence of 5 mM iodacetamide for 5 9. mins at room temperature.

Store extracts of radio-labeled platelets frozen at -20°C. 10.

Analysis of ³H-labeled platelet membrane proteins will be performed in Boston using:

1) SDS-PAGE and

2) Two dimensional isoelectric focusing and SDS-PAGE, essentially as described in the previous section for ¹²⁵I-labeled membrane proteins. In the case of the ³H-labeled samples, radiolabeled glycoproteins will be detected by fluorography using diphenyloxazole as described by Bonner and Laskey (Eur. J. Biochem. 46:83, 1974).

Cytoskeleton Studies

These experiments will analyze both qualitative and quantitative changes in the composition and assembly state of platelet cytoskeletal components.

Approach:

That percentage of the total platelet 1) actin and 2) tubulin present in the cytoskeletal fraction of platelet lysates after extraction with Triton-X-100 will be quantitated. The cytoskeletal fraction which is detergent insoluble will be harvested by selective centrifugation. Actin filament content will be determined by quantitative densitometry of Coomassie blue stained SDS-gels of electrophoretically fractionated cytoskeletons. Microtubule content will be determined as "sedimentable tubulin" also by densitometry of Coomassie blue stained SDS-gels of electrophoretically fractionated cytoskeletons but from platelets that were pretreated with Taxol to stabilize assembled microtubules prior to extraction with Triton-X-100.

Literature:

Fox et al., J. Cell Biol. 98: 1985, 1984. Kenney and Linck, submitted for publication, 1985.

Equipment:

Low speed centrifuge Ultracentrifuge and/or Air fuge (Beckman Instruments) Cell counter 37°C waterbath 60°C waterbath boiling waterbath

Platelet Preparation for Assays Described Below:

- Wash platelets twice in TES/Tyrode's buffer, pH 7.4 containing 10 ng/ml prostaglandin E1, and 5.5 mM glucose maintained at 22-35°C.
- Resuspend platelets at a final concentration of 0.5 x 10°/ml in TES/Tyrode's pH 7.2 with 5.5 mM glucose.
- Equilibrate for 1-1.5 hours at 37°C.

Assembled Microtubules:

Reagents:

- MT-Cytoskeleton Extraction Buffer: TES/Tyrode's pH 6.8; 12 mM EGTA; Leupeptin 10 ug/ml; 1% Triton-X-100 at 22°C.
- 2. Taxol: 20 mg/ml in dimethylsulfoxide stored at -20°C.
- 3. 27% Sucrose: in TES/Tyrode's pH 6.8 with 6 mM EGTA.
- 4. Solubilization Buffer: 1% SDS in 0.125 M Tris-HCl pH 6.8; 6 mM EGTA; 10 m MDTT; 0.5% Na azide.

Procedure:

- Treat 10 ml aliquot of washed, 37°C equilibrated platelets with Taxol (4 ug/10⁸ platelets) for 2 mins at 22 + 2°C.
- Extract platelets by mixing with 10 ml MT-cytoskeleton Extraction Buffer for 5 mins at 22 ± 2°C.
- Layer 20 ml extract on top of 27% Sucrose.
- 4. Pellet assembled microtubules at 100,000 x g for 30 min at 25°C.
- 5. Discard supernatant and sucrose layer.
- 6. Solubilize pellet (prepared from 5.0 x 10° platelets) in 400 ul Solubilization Buffer, heat at 60°C for 30 min and then boil for 2 min. Store samples frozen at -20°C for analysis in Boston.
- 7. To determine total platelet tubulin (microtubule protein), extract a parallel 10 ml aliquot of each platelet suspension, with 10 ml extraction buffer. Add Taxol to a final concentration of 50 ug/ml of extract (assuming that each ml contains protein from .25 x 10° platelets), and incubate for 30 mins at 37°C; Pellet and solubilize microtubules as described above in #4 and #6.

Actin Filament:

Reagents:

- A-Cytoskeleton Extraction Buffer: 2% Triton-X-100; 10 mMM EGTA; .1 M Tris-HCl, pH 7.4 at 22°C.
- Solubilization Buffer: (see Microtubules).
- 3. 4X Solubilization Buffer: 0.48 M Tris-HCl, pH 6.8; 9.6% SDS; 40 mM MDTT.

Procedure:

- Treat 5 ml washed, 37°C equilibrated platelets with 5 ml A-Cytoskeleton Extraction Buffer at 22 <u>+</u> 2°C.
- 2. Pellet insoluble cytoskeletons at 150,000 x g at 0°C for 45 mins.
- Wash sedimented material with 1:1 mixture of TES/Tyrode's pH 7.4 and A-Cytoskeleton Extraction Buffer.
- 4. Solubilize actin cytoskeletons from 2 x 10° platelets in 200 ul solubilization buffer, boil for 2 min, then freeze at -20°C for analysis in Boston.
- 5. Extract 1.5 ml aliquot of each platelet suspension directly with 500 ul 4X Solublilization Buffer for determination of the total actin in platelet extract. Boil the extract for 2 min, then freeze at -20°C for analysis in Boston.

Analysis of Cytoskeleton Polypeptides

Analysis of cytoskeleton polypeptides will be carried out on polypeptides which are electrophoretically fractionated using the SDS system of Laemmli (Nature 227:680, 1970) with 4% stacking gel and 8% resolving gels.

- 1. Cytoskeleton pellets (and platelet extracts used for determination of total platelet tubulin and actin) are solubilized with 1% SDS, brought to a final concentration of 0.0625 M Tris HCl, pH 6.8, and the disulfide bonds reduced with 5% beta-mercaptoethanol.
- 2. Electrophoresis is carried out at 24 mAmps per gel for 4-5 hrs.
- Gels are stained with 0.025% Coomassie Brilliant Blue and destained with 30% methanol/10% acetic acid and photographed.
- 4. Molecular weights of cytoskeletal polypeptides are determined based on migration of the standard polypeptides: spectrin heterodimer 240,000 and 220,000; beta galactosidase, 130,000; phosphorylase a, 94,000; Bovine serum albumin, 68,000; glutamate dehydrogenase, 56,000; creatine kinase, 40,000 and carbonic anhydrase, 29,000.
- 5. Major polypeptides of the platelet cytoskeleton to be examined for qualitative and quantitative changes are: Actin binding protein, 260,000; Myosin, 200,000; alpha actinin, 90,000; Tubulin alpha, 58,000; Tubulin beta, 53,000 and actin, 43,000.
- 6. Amounts of individual polypeptides in cytoskeletons and supernatant fractions will be quantitated by densitometry of Coomassie Blue Stained gels using a Helena Quick Scan Densitometer with a 570 nm filter. Mean values will be determined from three scans of each sample.
- 7. Percent of the total cytoskeletal protein which is assembled in organized cytoskeletal structures will be determined by quantitating that particular polypeptide in the insoluble cytoskeletal fraction relative to the total

quantity of that platelet protein. It is anticipated that during storage 1) the content of cytoskeleton associated microtubule protein will decrease and 2) the content of cytoskeleton associated actin will increase.

Viscosity

Viscosity determinations will be made on a modified cone plate Wells-Brookfield microviscometer. The modifications facilitate measuring shear stress at low shear rates. The cone and plate contain radial etchings 60 microns deep to reduce the plasma layer effect at the containing surfaces. All determinations are made at 37°C. Viscosity will be measured at shear rates from 45-2.2 sec -1. Each determination requires 0.5 cc sample and is done in duplicate. The instrument is standardized with samples of known viscosity obtained from Wells-Brookfield.

THE COLORIMETRIC DETERMINATION OF GLUCOSE

REF: Sigma Technical Bulletin No. 635

REAGENTS:

A. <u>0. Toluidine Reagent, Stock No. 635-6</u>
 Store in dark at room temperature.
 Caution: Do not pipet by mouth. Avoid contact with skin.
 Do not inhale fumes.

.

- B. <u>Trichloroacetic Acid Solution, 3% (W/V), Stock No. 635-3</u> Store in refrigerator at 0-5°C.
- C. <u>Glucose Standard Solution, Stock No. 635-100</u> Standardized at 1.0 mg/ml (100 mg/100 ml or 5.5 m mol/1) with Benzoic Acid added as preservative Store in refrigerator at 0-5°C.
- D. <u>Hyland Q-Pak--Chemistry Control Serum 1</u> List No. 045-030 Store between 2° and 8°C before and after reconstitution.

DEPROTEINIZATION OF PLATELET CONCENTRATE:

Into a centrifuge tube, pipet: 0.4 ml sample 3.6 ml of 3% Trichloroacetic Acid, Stock No. 635-3 Mix well by shaking. Allow to stand approximately 5 minutes to precipitate proteins. Centrifuge (5-10 minutes) until clear supernatant is obtained in the RC-3 at 3000 rpm. Store clear supernatant at 0-5°C, if determination is to be within a few hours. Otherwise freeze clear supernatant.

INSTRUMENT:

Practically any photoelectric colorimeter that transmits light in the range of 620-650 nm can be used.

Instrument: Coleman Jr. Wavelength: 635 nm Cuvette: 19 x 150 mm

> NOTE: Instrument readings should be made in terms of Absorbance (A) for use in calculations and for preparation of calibration curves.

ASSAY PROCEDURE:

1. Label test tubes or cuvets BLANK, STANDARD, CONTROL, TEST 1, TEST 2, etc.

TO BLANK ADD:TO STANDARD ADD:0.1 ml water0.1 ml Glucose Standard0.9 ml 3Z TCA SolutionSolutionStock No. 635-3Stock No. 635-1000.9 ml 3Z TCA Solution0.9 ml 3Z TCA SolutionStock No. 635-3Stock No. 635-3	TO CONTROL & TESTS ADD: 1.0 ml supernatant
--	---

- To each tube add: CAUTION: DO NOT PIPET BY MOUTH
 5.0 ml O-Toluidine Reagent, Stock No. 635-6 Mix by lateral shaking.
- 3. Place all tubes in a vigorously boiling water bath for exactly 10 minutes. NOTE: If results are calculated from a standard that is assayed with the Test, the boiling period can be 10 ± 1 minute.
- Quickly remove all tubes and cool to room temperature by placing in tap water for approximately 3 minutes.
- Transfer contents of tubes to cuvets and read Absorbance of STANDARD and Tests at 635 <u>+</u> 15 nm, using BLANK as reference. Complete readings within 30 minutes.

CALCULATIONS:

Use	of	Sta	ndard:					
Gl	uco	se	(mg/100	ml)	=	<u>ATest</u> AStandard	×	1000*

*Represents the concentration of the Glucose Standard (mg/100 ml).

NUTE: If the reading of Test indicates a Glucose concentration greater than 250 mg/100 ml, dilute the test with an equal volume of U-Toluidine Reagent, Stock No. 635-6. Read Absorbance of the diluted Test and multiply result by 2.

DETERMINATION OF LACTIC ACID REF: Sigma Technical Bulletin No. 826-UV

REAGENTS: Α.

- Lactic dehydrogenase (Sigma stock #826-6).
- 2. Glycine Buffer (Sigma stock #826-3).
- NAD preweighed vial (Sigma stock #260-110).
- 4. Lactic acid standard solution (Sigma stock #826-10).

EQUIPMENT: Β.

1. Spectrophotometer

DEPROTEINIZATION OF PLATELET SAMPLE: C. Trichloroacetic acid extract of sample prepared for glucose assay will be used.

ASSAY PROCEDURE: D.

1. Pipet into each of the NAD vials

2 ml Glycine Buffer 4 ml water 0.1 ml LDH

The number of vials has to be determined based on each vial per 2 assays.

2. Combine the reaction mixture prepared above in a flask.

3. Label an appropriate number of tubes and pipet into each tube; 2.8 ml.

4. To blank add 0.2 ml of 3% TCA. To sample tube add 0.2 ml of respetive extract, and mix gently.

Incubate at 37°C for 30 min or at 25°C for 45 min. 5.

6. Read absorbance of sample tubes at 340 nm vs. blank as reference. At this point the reading should be stable (less than 0.001 per min increment). If it is changing more than 0.002 per min, an additional incubation for 15 min should be done.

 $A_{340} \times 7.23 \times \frac{10}{3} = mmol L A/L$ 7. Calculation: $A_{340} \times 65.1 \times \frac{10}{3} = L A mg/100 ml$

or

B-Thromboglobulin (B-TG) RIA Kit

Code IM.88

 β -Thromboglobulin (B-TG) RIA Kit

Intended Use

the Amersham β . Thromboglobulin (β -TG) RIA Kit provides a quantitative method for the direct measurement of β -TG in human <u>platelet poor plasma over</u> the range 10-225ng /J-TG/ml

Summary and Explanation of the Test

The platelet-specific protein. B-Thromboglobulin is released into the circulation when blood platelets undergo the release reaction.

blood vessel starts a sequence of events including platelet release. aggregation and formation of a platelet "plug" which stems blood loss and is followed by The platelet release reaction is essential in the primary role of blood platelets in controlling hemostasis. Under normal circumstances, injury to the wall of a the formation of a stable fibrin clot. Under certain conditions, these processes can occur within the circulation without injury to the vessel wall.

The appearance of a platelet-specific material in the circulation produced by the platelet release reaction could provide a means of monitoring this reaction.

by Moore. Pepper and Cash.³ The protein has a molecular weight of 36.000 is no function has yet been ascribed to β -TG. β -Thromboglobulin is released when platelets undergo their release reaction⁴ and is present in only minute smounts in other tissues.⁵ β -<u>Ihromboolobulin can therefore</u> he considered ss <u>a platelet specific protein. The release of which is a marker for the platelet</u> quently. It was isolated, characterized, and named eta-Thromboglobulin (eta-TG) thought to comprise six identical sub-units and is probably located in the alphagranules of platelets. It is the most abundant platelet-specific protein, though The presence of a platelet-specific β -globulin had been suggested 1.2 Subserelease reaction

Principles of the Procedure

C-2

12.1 β -1G for a limited number of binding sites on a β -1G specific antibody. The amount of ¹³³1-labeled β -1G bound by the antibody will be inversely progation and removal of the supernatant, the precipitated radioactivity is measured in a gamma counter. By measuring the proportion of 1311 β -TG bound in the Portional to the concentration of unlabeled B-TG present. The antibody-bound 171 pt-TG is separated by precipitation with ammonium sulfate. After contribupresence of a series of eta-TG standards, the concentration of eta-TG in unknown The radioimmimoassay method depends on competition between A-IG and samples can be interpolated from a standard dose response curve.

> 2636 South Clearbrook Drive Arlington Heights, Illinois 60005 1 (800) 323-0668/(312) 364 7523 Amersham Corporation Printed in U.S.A. Amersham is a trademark.

1 (800) 268 5061/(416) 842-2720 Oakville, Ontario L6H 2H3 505 troquois Shore Road

In Canada

Amersham

2789608

Amersham Corporation 2636 South Clearbrook Drive Artington Heights, Ittinois 60005 1-(800) 323 0668/(312) 364-7523

In Canada

Mmersham

505 troquors Shore Road Oakville. Ontario 1.6H 2R3 1 (800) 268 5061((416) 842 2720

November, 1983

ORIGINAL PAGE IS OF POOR QUALITY

2. Recovery of added $[^{3}H]$ -TXB ₂ has been consistently >90% and it is considered unnecessary to employ a $[^{3}H]$ -TXB ₂ recovery marker. 3. It is good practice to assay, in each run, an aliquot of extracted distilled water to ensure that non-specific interfering materials have not been introduced from solvents, etc. Any detectable method blank must be corrected for as described for plasma (Section II C). 4. It is, of course, the user's responsibility to validate any sample preparation procedure used for uring other than the one recommended here.	In the event the user employs a plasma sample preparation procedure other than that recommended in Section V B and Appendix I, he should determine: (1) the recovery of (³ H)-TXB ₂ through his procedure, and (2) the effect of his solvents and sample handling techniques with regard to the introduction of non-specific interfering materials. Suggested methods for accomplishing these studies are as follows: A. [³ H)-TXB ₃ Recovery	 Add -2000cpm [³H]-TXB₂ per ml of plasma sample. Process samples by your selected method with final reconstitution prior to assay being made with assay buffer provided in the kit. After allowing at least 30 minutes for complete dissolution, the reconstituted extract should be aliquoted for the RIA and for counting of the [³H]-TXB₂ recovered. Note: Contributions of recovery marker to mass of TXB₂ measured in RIA must be subtracted if the user chooses to add [³H]-TXB₂ to samples.) B. <u>Matrix Correction</u> 	 Extract, reconstitute, and assay a prostaglandin-free human plasma sample by the same method as planned for your unknowns. Prostaglandin-free human plasma can be made in the following manner: Collect a sufficient volume of normal human plasma following the procedures and precautions given in Part V. Strip the sample overright with 1000mm activated charcoal*. Centrifuge to remove the bulk of the charcoal, then filter down to 0.22µm to remove fine particulates. If this results in 6958 B/Bo, it will undoubtedly be necessary to correct for this interference by one of the two methods given in #3 of the procedural notes in Appendix 1. 	*Norit A, the trade name of a decolorizing carbon which is commercially available from Amend Drug and Chemical Co., Irvington, NJ 07111, is recommended.
 C. Procedural Notes 1. Solvents should be prepared fresh on the day of use and stored in glass containers. 2. Recoveries of [³H]-TXB₂ have routinely been >90% using this method and no [³H]-recovery marker is necessary. 3. Although this method results in considerably cleaner eluates than conventional methods. a small method blank (>90% B/Bo) may be obtained. This must be accounted for by preparing for each assay, by the same procedure as described above, an extract of distilled 	 Fun the extract along with other samples in assay. Extrapolate the apparent TXB, value (pg added) using log-logit graph paper and subtract that from each sample value (pg added) before subsequent final calculations to pg/ml. b. Alternatively, add an aliquot of the distilled water extract to the entire standard curve. Keep the reaction volume constant throughout the assay by adding an equal aliquot of assay buffer to each sample tube. Using this method, any method blank will automatically be corrected for when reading sample values of this standard curve. Assay binding (Bo) will be somewhat decrased because of the large to obtained. Fallure to correct for any differ the sample volume. 	 B/BO), which may vary slightly from day to day, may result in the over-estimation of samples having very low TXB₂ levels. 4. Concentration of samples 2-4 times should be adequate for most samples tested. This can be accomplished by using 2-4ml of plasma for the procedure and reconstituting the final extract to 1.0ml. 5. Care must be taken to ensure complete reconstitution of the dried extract by adequate sourcexing, since the final volume is less than the original solvent volume of the eluate. 6. This procedure has been validated for single use of BOND-ELUT columns only. 	 III. <u>URINE</u> A. <u>Preparation and Extraction</u> 1. Acidify an aliquot of the urine to pH 3-3.5 with 2N HCl or 2N Citric Acid. 2. Extract with two volumes of ethyl acetate, vortexing 30 seconds. 3. Centrifuge 5 minutes at 1000 x g to separate phases. 4. Carefully remove 1ml of the organic (upper) phase and elute through a BOND-ELUT 5I column prepped with 5.0ml Benzene: Ethyl Acetate (80:20). 5. Continue with 5I column extraction as described for plasma (Steps 2e-2h). 	The concentration of TXB ₂ in urine is generally high enough to make sample concentration unnecessary. Extraction of 2ml of urine by the procedure above with final reconstitution to 1.0ml results in a final 1:2 sample dilution for assay and should fall within a readable portion of the standard curve for the majority of samples assayed.

ORIGINAL PAGE IS OF POOR QUALITY

-16-

-11-

-15-	 PLASMA PLASMA	with assay burger, mix and allow to support the transmission of the extract through a 0.22µm Miller filter using a polypropylene or siliconized glass syringe into a suitable tube or vial.
-14-	 Sors, H., Pradella, P. and Dray, F., Prostaglandin, 16:277 (1938). Dray, F., et, al. Accos. In Prost. and Thromb. Res. (18) (1990). Wandou, S. Friso, (Ed.). Academic Verk, pp. 137 (1931). Sors, J., Friso, (Ed.). Academic Verk, pp. 137 (1931). Sors, J., Friso, (Ed.). Academic Verk, pp. 137 (1931). Sors, M. Caran, D. Toharsty, S. and Sorkin, R., Clin. Chem. J. 21 (1971). Merris, H.G., Sharman, N.A. and Shapperdon, F.T., Prostaglanding, 21 (1971) (1981). Prosti, and Thromb, R.J., Sustanon, N.A. and Shapperdon, F.T., Prostaglanding, 21 (1971). Strinka, V. and Ramel, N.A. and Shapperdon, F.T., Prostaglanding, 21 (1971). Prosti, and Thromb. Res., 19 (1980). Prosti, and Thromb. Res., 19 (1980). Strinka, V. and Dray, F., J. and Oles, J.A. J. 1900. Prosti, and Thromb. Res., 19 (1980). Resta, L.J., Sweatman, B., and Oles, J.A. J. 1900. Prosti, and Thromb Res., I. and Oles, J.A. J. 1900. Prosti, and Thromb Res., F. Prostaglandin, 22 (1991). Prosti, and Thromb Res., J.A. J. 2000. Resta, J.L.J., Sweatman, B. and Oles, J.A. J. 1900. Resta, J.L.J., Sweatman, B. and Oles, J.A. J. 1900. Resta, J.L.J., Sweatman, B. and Oles, J.A. J. 1900. Resta, J.L.J., Sweatman, B. and Oles, J.A. J. 1900. Resta, J.L.J., Sweatman, B. and Oles, J.A. J. 1900. Resta, J.L.J., Sweatman, B. and Oles, J.A. J. 1900. Resta, J.L.J., Sweatman, B. and Oles, J.A. J. 1900. Resta, J.L.J., South P. Prostaglandin, P. A. Kanda, P. J. 2010. Resta, S.L.J., Shara, P. Kandari, A. Kanda, K. Kand	ORIGINAL PAGE IS OF POOR QUALITY

.

-

.

 D. Conparative Sample Values: 1261 vs. 341 Kits Sample values were compared by assaying an extracted plasma sample at various dilutions in both the 1261 kit (NEK-025) and the 341 kit (NEK-025) and the 1341 kit (NEK-025). The sample extract was concentrated four times in order to cover the standard curve range of both kits. Extract Dilution 1281 - Kit 341 - Kit 341 - Kit 	straight off curve 955 112 999 999 113 999 999 114 999 999 115 919 919 115 919 919 1112 919 919 1112 910 Resentiuity of the system was found to be approximately 1112 Ramwell, P. W., Biol. Reprod. 152.00 (197) 919 (1918) 1112 Ramwell, P. W., Biol. Reprod. 152.00 (197) 919 (1918) 1112 Ramwell, P. W., Biol. Reprod. 152.00 (197) 919 (1918) 111 Ramwell, P. W., Biol. Reprod. 152.00 (197) 919 (1918) 111 Ramvelsson K. et al., Ann. Reprod. 152.00 (197) 919 (1918) 112 Honerada, S. and Vana, J. R. Reprod. 152.00 (197) 919 (1918) 113 Honerada, S. et al., Ann. Reprod. 190 (1918) 91000 114 Honerada, S. et al., Ann. Reprod. Reprod. 2019 (1918) 91000 115 Honerada, S. et al., Ann. Reprod. Reprod. 2019 (1918) 91000 116 Honerada, S. et al., Ann. Reprod. Reprod. 2019 (1918) 91000 113 Honerada, S. et al., Ann. Repro. 10000 910000 <	
 B. The thromboxane B₂ standards are prepared in the assay phosphate buffer. The effect of other sample matrices upon the assay system must be determined by the investigator. C. Any exogenous radioactivity will lead to erroneous results. Proper controls should be employed to determine if samples are contaminated. IX. <u>PERFORMANCE CHARACTERISTICS</u> A. <u>Reproducibility of the Assay</u> Reproducibility was determined by running multiple duplicate analyses of 	PRECEDING PAGE BLANNING MORE DI VINIG AND COMPACT FLAMED Sample A Strand branch and bra	

-12-

-13-

RECEDING PAGE BLANK NOT FILMED

ORIGINAL PAGE IS OF POOR QUALITY

I

ا 1

mj i

-1

[]] ÷.

ī ī

7

1

-.



÷

thromboxane B_2 added, sample values must then be corrected for method blank, aliquots, dilution and recovery to determine the original concentration semi-logarithmic graph paper, plot & B/Bo for each standard versus The values If all tubes have been counted for the same period of time, use the total Determine the normalized percent bound (\$ B/Bo) for each standard and See Figure 3 for a typical standard curve using the standard protocol.) the corresponding amounts of thromboxane B_2 added in picograms (pg) Determine the pg thromboxane B_2 in each sample by interpolation from <u>_</u> of the the standard curve. Since the standard curve is expressed as pg determined from a standard curve. The following method After counting has been completed, the concentration of thromboxane \mathbf{B}_2 the samples is determined from a standard curve. The following method accumulated counts; otherwise, correct all raw counts to counts per Average the counts for each set of duplicates. Calculate the average NET counts for all standards and samples by Any samples with concentrations which are above the range with assay buffer and re-assayed. subtracting from each the average blank counts (tubes 3-4) obtained are then multiplied by the appropriate dilution factor **&** B/Bo = Net cpm of standard or sample X 100 Net cpm of $\frac{100}{100}$ standard (See Table 2 for sample calculations.) VII. PROCEDURE FOR CALCULATING UNKNOWNS curve may be diluted sample as follows: sample. minute (CPM). in the suggested. Using standard NOTE: ວ່ ຕໍ ÷ ŝ ف . 20-30 minutes 3 To minimize S ŝ

÷

÷

solution into each tube and mix Pipet 100ul of

- antiserum into all tubes beginning with tube 5 and seconds. vortex thoroughly for 2-5 Pipet 100ul <u>ы о</u>
 - Incubate overnight (16-24 hours) at 2-8°C Ξ
- At the end of the overnight incubation, place all tubes in an ice hath ž
- Pipet 1mi of cold precipitating reagent into all tubes beginning with tube 3. Vortex each tube thoroughly for 2-5 seconds. Allow tubes to incubate at 2-8°C for 20-30 minutes. З.
 - Centrifuge the tubes in a refrigerated centrifuge at 1000-2000 x g <u>15</u>.
 - Decant the supernatants of all tubes beginning with tube 3 and for 30 minutes. ē.
- allow to drain for ~ one minute on absorbant paper. Blot the tubes to remove residual liquid at the rim.
 - Count all tubes in a gamma counter. At normal efficiencies, a one minute count time is sufficient. 17.
 - Calculate results as described in Section VIII. 18.

Table 1 - Thromboxane B, Assay Protocol Schematic (All volumes are in microliters)

i

.

1

	I UDE NO.	BULLEL	Standard	Samples	I Lacer	Antibooy
Total Counts	1-2	;	:		001	ţ
Blank	9-F	200	:	1	81	ţ
0 Standard	2-6 5	100	:	:	100	100
Standards	7-20	;	100	1	100	100
Samples	21,22 et		-	00	001	9 0

Add 1mi of cold Precipitating and centrifuge at refrigerated temperatures for 30 minutes at Decant all tubes except total count tubes, blot and Mix, incubate Incubate overnight (16-24 hours) at 2-8°C. Reagent to all tubes except total counts. 1000-2000 × g. at 2-8°C. count.

Precautions ò

- polypropylene or siliconized glass pipets, pipet tips and tubes when Pipetting must be done reproducibly and accurately. Thromboxane assay interference from this "sticking" problem, use only B, has the tendency to adhere to many surfaces. transferring diluted materials or incubating. -
- for any method blank that may arise as a result of the extraction or An aliquot of an appropriate control should be assayed to account responsibility to check and correct for non-specific matrix and chromatographic steps (see Appendix 1). It is the user's solvent effects. e,
 - centrifugation time or speed may result in incomplete precipitation Inadequate incubation time with the precipitating reagent or bound counts. ٦ ÷
- inadequate centrifugation speed or prolonged inversion after decanting may cause the peliets to become dislodged from the bottom of the i het ÷
- Since the assay is sensitive to sub-picogram quantities of TXB_2 , extreme care must be taken to avoid contamination of laboratory equipment (such as pipettors) with TXB₂. ທ່

PRECEDING PAGE BLACK NOT FILMED

ORIGINAL PAGE IS OF POOR QUALITY

η

1

"]

Ţ

11

1.8

-1-	 B. Preparation of Thromboard B, Working Standard concentrate is diluted with sample varies but for the reconstituted TXB, standard concentrate is diluted with sample varies but the sample varie of pote to folgo added per (1)m() is shown below. Other dilution schemes below). A suggested for dilution scheme to cover a standard during regime to be used but the samp buffer provided with the sample varies of the sample va	
	d by d by strutulate lon to strutulate lon to struture te his atlats atlats atlats dispenser	
	1, pp. 14 samples p an of TXB LUT C ₁₆ c ins followe reing particulation of added ssay. P assay. P as	Ē
	Appendix ication of purification move inter rent/colum action and covery inter action and coveries inter action acti	
	i the purif the purif on ND-ELU on ND-ELU on ND-ELU on ND-ELU on the solver i from solver	
-9-	and the set of the set	•
	three major but of procession of concer- smainto of purification of recon- apparentij mear result of recon- mear result near result free the ti- ed. Our says the resp inter- pitation i increate ass test ti- increate ass test ti- tors, suc-	-
	e recommer d involves Extraction from pla from pla from pla from pla from pla fruction (rutation (rutation (rification (rification)	~
	Image: State of the state	
	Street to the second s	

-

—

ORIGINAL PAGE IS OF POOR QUALITY

E. <u>Precipitating Reagent</u> The solution contains 16% polyethylene glycol (PEG 6000) and 0.05% sodium azide in 50mM phosphate buffer, pH 6.8. Stored at 4°C, it is	stable for at least two months. This reagent is quite viscous and use of a positive displacement device facilitates dispensing. V. <u>SAMPLE HANDLING</u> A. <u>Collection and Storage</u> it is recommended that all samples be processed immediately after collection	Blood samples should be collected in pre-chilled siliconized glass or polypropyremeters tubes coaled with a solution of 4.5mM EDTA combining a prostaglandin synthetase inhibitor such as indomethacin or aspirin (24) indomethacin has been reported to be very effective at concentrations of up to 10µg/ml and we have been able to demonstrate the lack of interference in the assay at these concentrations. The plasma fraction should be tricks in the assay at these concentrations. The plasma fraction and informated from the whole blood as soon as possible after collection and frozen at -70°C if it is not assayed on same day. Tissue samples should be deep frozen on dry ice or in liquid nitrogen samples should be stored at -70°C or lower. Tissue samples should be processed in the presence of prostaglandin synthetase inhibitors (25), such as indomethacin, at the concentration, cited above.	B. Preparation for Assay In many plasma samples, the low level of analyte will necessitate extraction and subsequent concentration prior to assay. Normal extraction and subsequent concentration prior to assay. Normal There is considerable ambiguity in the literature concerning the extraction of prostaglandins. Useful solvent extraction procedures have been reported (25,28), but considerable care must be taken to produce a sample clean enough for assay. At NEN, solfd-phase extraction using modifications of procedures reported by Skrinska and Lucas (29) has been found preferable to liquid-liquid extraction methods. Commercially available chemically bonded extraction have been used by others (30,13). These have yielded consistently have been used by others (40,14).	*BOND-ELUT C _{1a} and SI extraction columns are available from Analytichem International, Harbor City, CA. SEP-PAK C _{1a} and Si extraction columns are available from Waters Associates, Milford, MA.
INSTRUCTIONS RELATING TO THE HANDLING, USE, STORAGE, AND DISPOSAL OF THIS RADIOACTIVE MATERIAL.	This radioactive material may be received, acquired, possessor, and used only by research laborators for <u>in vitro</u> laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use, and transfer are subject to the regulations and a general license of the U.S. Nuclear Regulatory Commission or of a State with which the Commission has entered into an agreement for the exercise of regulatory authority. 1. All radioactive materials should be stored in specifically designated 2. All work with these materials should be carried out only in authorized	 areas. No pipetting should be done by mouth. There should be no smoking or eating within the work area. There should be no smoking or eating within the work area. Any spilled material should be wiped up quickly and thoroughly and the contaminated substances transferred to a suitable receptacle. The surfaces involved should be washed thoroughly with an appropriate decontaminate the vial. This radioactive material can be discarded into the sanitary sewerage system, provided the discharge concentration of ¹²⁹ does not exceed 4 x 10 µc/im. (See 10 CFR 20.303 for maximum daily, monthly and yearly quantity discharge limits.) Prior to disposal of the empty, uncontaminated kit and tracer containers in unrestricted areas, remove or deface the radioactive material should be waterial can be disconted with an average solution of the samilary severage system. 	 B. TXB, Antibody B. TXB, Antibody The rabbit anti-TXB, is supplied ready to use. Stored at 4°C, the solution is stable for at least two months. C. <u>TXB, Standard Concentrate</u> To reconstitute, add exactly 1.0ml of distilled water. The reconstituted solution contains 100ng/ml of thromboxane B, in buffer. Immediately before use in the assay, dilute an aliquot of the stock solution to prepare standards. A suggested procedure for preparing diluted standards. Pipels and/or pipet tips used to transfer diluted standards. Pipels and/or pipet tips used to transfer diluted standard at 4°C. Under these 	D. <u>Assay Buffer</u> Assay Buffer Assay Buffer The solution consists of 0.94 NaCl, 0.01M EDTA, 0.34 bovine y-globulin, 0.0054 Triton-X-100 and 0.054 sodium azide in 50mM phosphate buffer, pH 6.8. When stored at 4°C, it is stable for at least two months.

+

÷.

ORIGINAL PAGE IS OF POOR QUALITY

Ļ

m 1

77

7 ŝ -----1 ١ ----Π 1 ----1 ł

In the New England Nuclear Thromboxane B ₂ [¹²⁶] RIA Kit, separation of the antibody-antigen complexes from free antigen is achieved by precipitation of the antibody-bound tracer with polyethylene glycol in the presence of carrier immunoglobulin. After centrifugation, the supernatant, containing the unbound antigen, is decanted and the pellet, containing the antibody-antigen complex, is counted in a gamma counter. Results obtained for the studends are used to construct a standard (dose-response) curve from which the unknowns are read by interpolation.	 AGEND ANTIGEN • SPECIFIC ANTIGEN • (In standard solutions Ages) UNLABELED ANTIGEN • MATTRODY COMPLEX • Ages UNLABELED ANTIGEN • MATTRODY COMPLEX • Ages Ages Ages	
Figure 1 Schematic of Arachidonic Acid Metabolic Pathway LINOLEIC ACID ARACHIDONIC ACID	Prostaglardin Synthetase PGG, THROMBOXANE B, Thromboxane Synthetase Front Front Prostacyclin Synthetase Fromboxane Synthetase From Prostagname out PROSTACLANDIN F. Fromboxane B, TTXB, PCE, PCE, PCE, PCE, PCE, PCE, PCE, PCE	

÷

÷

PRECEDING PAGE BEANK NOT FILMED

ORIGINAL PAGE IS OF POOR QUALITY 1

·• ,

₹* 3

~ •]

-1

. .

1

I INTRODUCTION	The schematic for the "in vivo" metabolic pathways of the major prostaglandins, shown in Figure 1, represents only the most general outline of a very complex series of blochemical interconversions. Linoleic acid is an essential fatty acid which is metabolized to archidonic acid, the starting compound of the cascade. Arachidonic acid is stored in cell walls esterified in phospholipids (1). Upon demand, arachidonic acid may be released from the cell wall by phospholipase A_2 (2).	Prostaglandin synthetase contains both cyclo-oxygenase and peroxidase activities to convert arachidonic acid to prostaglandin endoperoxide. Cyclo-oxygenase enzymatically metabolizes arachidonic acid to prostaglandin G ₂ (PGG ₂), a cyclic endoperoxide, while peroxidase reduces PGG ₂ to another cyclic endoperoxide, prostaglandin H ₂ (PGH ₂) (3,4). Both PGG ₂ and PGH ₂ have a hair-life of about 5 minutes at 37°C in acueous buffer at pH 7.4 (5). Prostaglandin endoperoxide PGH ₂ is considered a pivotal compound, since it is metabolized by three different, reaction pathways; whown to be present in introbolized by three different, reaction pathways; whown to be present in the Thromboxane A ₂ is rapidly hydrolyzed to thromboxane B ₂ (TXB ₂). In the prestacyclin synthetiase, which has been demonstrated in the microsomal prostacyclin is unstable and converts to 6-heto-prostagjandin F ₁ 3. Classical prostaglandins, PGE ₂ , PGE ₃ , can be formed from the prostaglandin (6-keto-PGF ₁).	Because of the opposing effects of prostacyclin and thromboxane A ₂ synthesis, there is the potential for the delicate control of hemostasis and arterial thrombosis in <u>vivo</u> . ⁴ That is, thromboxane A ₂ promotes plateiet aggregation and is a vasoconstrictor, while prostacyclin prevents plateiet aggregation is a vasoconstrictor. By quantitatively measuring both stable conversion products - TXB ₂ and 6-keto-PGF ₁ - researchers can study normal as well as pathological mechanisms. For example, pathological states can be generated in the experimental model by mainpulating the activity of selected anzyme(s) (8-11). Pharmacological agents can then be assessed for their effectiveness in preventing or alleviating various pathological cudute thromboembolis. (8), atheromatous plaques (12) diabetes (13) and a variety of congenital and acquired bleeding plaques (12) diabetes (13) and a variety of congenital and acquired bleeding	defects (14-16). There presents such as the lung and kidney and possibly certain cells in tissue culture can produce both thromboxane A_2 and a_1 prostacyclin (1,4,8). In addition, there appears to be a species dependent predisposition for producing larger quantities of one of the two compounds (6,17). In order to fully understand the biological roles of thromboxane A_2 and prostacyclin, it is first necessary to be able to monitor their production. For thromboxane A_2 and brostacyclin, it is first necessary to be able to monitor their production. For thromboxane A_2 and brostacyclin, their respective stable conversion products (TXB ₂ and 6-keto-PGF _{1d}) are the parameters of choice.	
					-
	- - -				L
• *		PRECEDING PAGE BLANK NOT	f filmed	original pa of poor qu	GE IS ALITY




OPICINAL FACE IS OF FOUR QUALITY Ī

ī

13

1

•

÷

Bibliography 1) Sokal, G., Etude morphologique des plaquettes sanguines et de la metamor- chose viscueuse au moven d'antiserums fluorescents antifibrinogene et anti-	plaquettes, Acta Haemat., 28, 313-325 (1962). 2) Davey, M.G. and Luscher, E.F., Release reactions of human platelets induced	by thrombin and other agents, <i>Brochim, Brophys.</i> Acta, 100, 490-909 (1990). 3) Moore, S., Pepper, D.S. and Cash, J.D., The isolation and characterisation of a matelet-specific B-dobulin (B-thrombodiobulin) and the detection of anti-	urokinase and antiplasmin released from thrombin aggregated washed human platelets, <i>Biochim, Biophys.</i> Acta, 379, 360-369 (1975).	4) Ludlam, C.A., Moore, S., Bolton, A.E., Pepper, D.S. and Cash, J.D., The release of a human platelet-specific protein measured by a radioimmunoassay. <i>Thrombosis Res.</i> 6, 543-548 (1975).	5) Ludiam, C.A. and Anderton, J.L., Platelet & thromboglobulin in The signifi-	cance of platelet function tests in the evaluation of hemostatic and thrombotic tendencies, Day, J.J., Holmsen, H. and Zucker, M. (eds.), in press.	b) Luciam, C.A. and Cash, J.D., Stoures on the menantic of Principle Section 19 (1976). from human platelets in vitro, British Journal of Haematology, 33, 239-247 (1976).	 Additional netrements regraming out-intromogramming (β-TG) - a new approach to the diagnosis of thrombosis, in <i>Platelets and Thrombosis</i> Symposium, Milan (1976). Ludlam, C.A., Bolton, A.E., Moore, S. and Cash, J.D., New rapid method for diagnosis of deep versions thrombosis, <i>Lancet</i>, 2, 259-260 (1975). Denham, M.J., Fisher, M., James, G. and Hassan, M., Plasma concentration of β-Thromboglobulin in venous and arterial thrombosis, <i>Lancet</i>, 1, 1154 (1977). 			
Effect	Lowers #-TG level by 3.5ng/ml	Lowers /J-TG level by 2ng/ml	Lowers /}-TG level by 1-2.5ng/mt	Lowers /}-TG level by 14-15ng/ml	- 	.hed Irom zero with					
ai Dose Tested	100µg/ml	10mg/ml	1010/m	50IU/ml		can be distinguis					
Max. Physiologic. Dose	53 µg/ml	1.2mg/ml	0-30IU/m			vel of <i>J</i> J-TG that .) is 5ng/ml.					-
Compound	Pontopon ^e – Roche	Aspirin	Heparin		Sensitivity	The fowest detectable le 95% confidence (± 2 S.D	Accuracy				
•							. 1		and the second	ORIGINAL P	
				I	PR	ECE	DIN	PAGE BLANK NOT FIL	MCD	OF FOOR Q	

- - •

.....

•

i

Sample Calculation	c			interest should be verified in the u	ser's laboratory. It sh	ould be r	noted that	elevat-
Tube No.	Standard	Counts/60 sec.	ng/J-TG/ml	ed p. 1G levels may result from im	proper sample cone			Þ
•	10ng/ml	12082		Performance Characteristics				
2	:	11862		Reproducibility				
ю	2 lng/ml	10339		The reproducibility of the /3-TG I	IA Kit has been eva	luated ar	id the res	ults ob-
4	:	10209		tained in a hospital laboratory are	Snown in Labia 4.			
Q	54ng/mł	1651		Table 2. Reproducibility of the /	I-TG RIA KIL	Ċ		
ø	:	1732				3		
7	101ng/ml	5880				<	-	c
æ	:	5829		Number of assays		66	54	39
ŋ	218ng/ml	3739		Man no //. TG/ml		9.6	43.6	108.0
10	:	3900				-	A R	16.5
11	Unknown 1	10543	20	Standard Deviation (each result the mean of duplic	ates)		 7	
12	:	10356		Confirment of variation &		6.6	10.7	15.3
13	Unknown 2	8302	44				- Hoom C	
14	:	8103		4 operators, 14 kit batches a	jed trom 2 - 14 week			
15	Unknown 3	6323	611	: : :				
16	•	5392		Specificity The following compounds were	tested and were for	nd not to	o interfere	with the
Details of Calibrat	tion			assay in concentrations up to th	ose calculated to be	the max	inum poi	sible cir-
The standards sup	pplied with the β -T	G RIA Kit have been a	accurately calibrated					
against primary β .	-TG standards at Ai	mersham Internationa A visition of 421 G/ml	i pic, and the exact		Max. Physiolog	ical		•
value of each stam	dard is printed on In	e viai in rig p- i colini.		Compound	Dose		Dose Te	sted
Limitations				Pethidine	0.5mg/ml		10mg	/ml
Care should be tak	ken to ensure that th	he patient samples do	not contain radioac-	Morbhine	13µg/ml		50040	/ml
tivity sufficient to	cause interference /	with the assay loce m		Sentrate - B Wellcome	1mg/ml		0.4mg	/mlt
				Detrait	17ուց/ուլ		20ing	/ml
EXPECTED VAI	LUES				. 1		0.4m	
It is postidated the	at the measurement	of 13-TG levels may b	je useful in assessing	Platelet-factor 4				
platelet release, w	which is a characteri	istic of numerous disc	orders. However, cur- erefore, this kit is in-	Ethylene-diamine-tetra-acetic	: 3mg/mf		20m	lm/g
tended for research	the use only and the	e results obtained sh	ould not be used for	Theophylline	I		2.56m	0/ml
clinical purposes.				Warlarin	. 0.1mg/ml		100m	۹/ml
For informative p	purposes, the mear	n values for //-TG le	vels were within the 31 annarently normal	troi soluble m higher c	oncentrations			

. .

2767

15

. .

OF POCK QUALITY

individuals obtained at 4 centers. Ninety-five percent of the values were below 52ng β -TG/ml. This does not imply that all normal persons will have values in research applications, the values to be found in controls and populations of range 24-28ng β -TG/ml on 124 β -TG determinations on 31 apparently normal this range or that persons with platelet disorders will fall outside this range For For info

PRECEDING PAGE BLANK NOT FILMED

and a second second second second

OF P008 QUALT

ORIGINAL PACE

The interference levels recorded by Heperin. Pontopon^{Θ} and sspirin are in excess of those normally observed in patients, therefore, unlikely to cause problems in measurement of patient samples

2

i

ţ

97.75 E

2

يە 2 مەر بەر مەر

والمحاصر والمسولاتين والمتعدية , . . and the second second m -_

ł

i ÷

:

η

7

n

η

11

Ī

۹

F

-.

5. Decantation

The supernatants should be decanted within 10 minutes of the end of centrifugation. After decantation, the tubes must be dramed for at least five minutes it is important to ensure that the tubes are not placed upright once they have been inverted until draining is complete.

6. Radioactivity Counting

The following points should be noted in the radioactivity counting:

- Sufficient counts should be accumulated to reduce the statistical variation to within acceptable limits for example, the standard deviation on 10,000 counts is 1%.
- b. The use of glass containers for secondary containment of the assay tubes in the gamma counter should be avoided, since the attenuation of the low energy ¹³³1 gamma radiation by glass is relatively high. Not only is the count rate markedly reduced but serious errors may arise due to variations in the wall thickness of the glass container. Where secondary containment is required, molded plastic tubes or vials are preferred.
 - c. Possible variations in the court rates due to the positioning of the assart tube in the gamma counter should also be checked.

Quality Control

The user may apply the following checks on assay performance

- The time required to accumulate 10,000 counts in the 10ng/ml standard (tubes 1 and 2) should be not more than 1-4 minutes depending on the efficiency of the counter and the age of the kit.
- The shape of the standard curve should be similar to that of the example shown in Figure 1.
 - 3. The ratio: Mean counts in tubes 1, 2 should be ≥ 1.7 Mean counts in tubes 7, 8

(Background counts should be subtracted in calculating this ratio.)

Assay Modification

In accordance with good laboratory procedure. <u>It least one control sample of</u> h<u>uman platelet poor plasma should be included in each assay run</u>.

Blanks

Corrections for samples containing radioactivity can be made by assaying the sample and a blank tube for the sample

.=

The original sample is assayed according to the usual protocol glanks are assaved according to the sample procedure substituting 200 at of waiter to the antiserum solution. In addition, a blank must be run on a standard (preterably the 10ng/ml standard) to use this correction. If the courts for the standard blank and the blank for the patient sample are within a small range (5%), then the acconcivity in the sample does not interfere with the assay because the radioactive material is not precipitated by ammonium suffate and is decanted in the supernatant If the counts for the patient sample blank are higher than that for the standard, the radioactivity present can interfere.

A corrected standard curve should be drawn plotting the standard counts - the

1

2017

į.

والمحتمد فالمعالية فالمعاد يتفكك كالمسابق وكالشا

:

f

the second second

i.

[]

1

π±

_]]

-

standard blank counts versus β - IG concentration and the sample β -IG concentration determined using the sample counts minus the sample blank counts.

na stan sen sen an La La Cala sen sen distan atau

Results

- 1. Using the linear graph paper provided, plot a curve of 13 l counts for the five standards against the β -TG concentration printed on the standard vials. Draw a smooth curve through the means of the duplicate points to give a standard curve. Grossly aberrant counts should be rejected. A typical curve is given in Figure 1.
- Using the mean of duplicate counts for the unknowns, read off their *B*-TG concentrations from the standard curve.



ORIGINAL PAGE IS OF POOR COMMITY

- the assay tubes (See Note 3). Inumediately vortex mix each tube thor-oughly until the contents are completely homogeneous. It is important to ensure efficient mixing at this stage
 - Centrifuge the tubes for 10-15 minutes at 1000-15000 (See Note 4). Remove the tubes from the centrifuge, place them in decantation racks. 1 ∼ œ
- DOUT Off and discard the supernatant liquids. The tubes should be gentry inverted in one continuous movement so that the precipitate will remain undisturbed at the bottom of the tube. Avoid tapping or shaking the tubes during decantation (See Note 5)
 - Keeping the tubes inverted, place them so that they drain onto paper issues for 5 minutes. ຫ່
- After draiming, gently touch the rims of the inverted tubes on paper lissues to remove any liquids remaining at the rims õ
 - Count the precipitates in a gamma counter (See Note 6) 11.

NOTES

1. Pipetting

As with other tests of this type, precision pipetting is essential. The use of a micropipette with disposable tip is recommended as they provide the er's operating instructions for use. The pipette tip should be moistened with the reagent before pipetting a series of aliquots. Clean tips should be used for each new sample or reagent. The pipettes used for the assay should have a precision of less than 1% standard deviation. It is recommended that the reproducibility of pipetting is regularly checked according to directions desired precision and convenience it is essential to follow the manufacturprovided by the manufacturer

- Incubation Times and Temperatures N
- The temperature range recommended for the assay is 15.30° C. The incubation time should b<u>e 60 ± 10 minutes</u>
- Separation Procedure e

The addition of the ammonium sulfate solution is time independent up to 8 minutes. After the addition of the ammonium sulfate solution, the assay neous. After mixing, there is no significant effect on the assay values if the lubes should be vortex mixed until the contents are completely homoge. tubes are left for up to 10 minutes before centrifugation

- Centrilugation
- It is recommended that centrifugation is performed at room temperature [15-30°C) However, the temperature of the centrifuge can be reduced to 2-4°C without affecting assay results.

For information, the g force (Relative Centrifugal Force or RCF) for a horizon tal head centrifuge can be calculated from the following formula:

RCF (g/= 11.2 × 10-6 × B × N²

where R = radius in centimeters from the center of the centrifuge shalt to the outside tip of the centrifuge tube when in the horizontal position

N = revolutions per minute

			nseaM	rivitosoiber s	ty in treetpi	.29161							
		Centrituge	snim 21-01	91-000L 18 -	nscan. Decan	reternaque t	sbiupil II						
			•	lz xim xərrov	νοιοπδμιλ'								
esettuč munommA noizutoč													
			Incuba	n te nuori file	noom temper	ature.							
				хэлоV	.xim								
P-TG ^{1 2 2 1} Anti-9-TG Serum	500												
		500						**************************************					
Unknown		—		—		09	09	09					
Disbraid	09	09	09	09	09								
Tube Numbers	1'5	3'4	9'9	8.7	01 '6	21,11	13, 14	15, 16 etc.					
	01~	02∿	09~	001~	\$ZZ2	L	2	3 etc.					
			Unicornes										
.r eldet		AIR DT-4	VezzA JiM A	· 1000101	•								

risticiation of all semular llA

-15-55

Second street

Π

7

7

[]

1

	The following sample collection protocol has been developed to allow prepara- tion of platelet poor plasma samples as simply and conveniently as possible. The protocol is a modification of that recommended by Ludiam and Cash ⁴ . The system employs a mixture of anticosguinat and antiplatelet reagents. Which, in	ent if the count rate for 50µJ patient sample is greater than 10% of the counts in the highest standard (after assay). Activity may be significant if the patient has received ¹²⁵ 1 thinnogen for throm- hosis detection. Preferably, samples for tJ-1G measurement should be obtained	RAGE IS
	complication with reduced temperature, immines p^{-1} is reasoning the processing <i>(in vitro)</i> and separates platelets and other formed elements of blood from the plasma by centrifugation techniques. In this way, the measured β -IG value reflects the <i>in vivo</i> circulating level of β -IG at the time of sampling.	prior to or 24 hours or more after ¹³⁵ I fibrinogen administration Procedure	POOR BOOR
	The tubes containing the antiplatelet and anticoagulant egents are provided in a separate container which has been designed to act as an ice bath for cooling samples.	Materials Provided Standard JJ-1 G in Buller — 5 vials	0 <u>F</u>
	Sample Collection and Preparation 1. Tap the sampling tube gently on a hard surface to ensure that the liquid con-	μ-TG ¹²⁵ 1 — 1 vial Anti-β-TG Serum — 1 vial	
	tents are at the bottom. 2. Place a mixture of CRUSHED ICE AND WATER in the base of the tube box and replace the rack containing the sampling tubes. Sufficient water is required	Ammonium Sulfate Solution — 1 vial Blood Sampling Tubes — 24	·
	to ensure good thermal contact and rapid cooling of the blood samples. 3. Using standard venipuncture techniques, withdraw a 2.5 ml blood sample using a 200 x 1" needle and svringe (preferably polystyrene). EVACUATED	Tube Rack/Cooling Bath Container Data Sheet	
	SAMPLE COLLECTION TUBES MUST NOT BE USED FOR COLLECTION OF BLOOD SAMPLES. A good veripuncture with minimal venous occlusion must be performed. When blood for other tests is drawn at the same time.	Materials Required For Sample Collection and Separation: Tube marker (waterproof)	
<u>.</u>	 4. IMMEDIATELY after collection, the needle should be removed from the syring e and 2.5ml of blood should be transferred gently to a pre-cooled labeled sampling tube (supplied with the kit). The tube is marked at 2.5ml 	Refrigerated centrifuge capable of 2000g — norizontal nead Crushed ice 5 or 10ml syringes — preferably polystyrene Needles — 20 gauge (1" length preferred)	ುಕ್ಕಾರ್ ಎಂ ಕ್ಷಮ್ಮ ಕೇರೆ ಎಂಕಿಕ್ಸಾ,
د ال الحمد و هذا الله مورد الحمد و الح	for convenience. 5. Stopper the tube and mix the contents by gentle inversion two or three times and IMMEDIATELY cool the sample by placing it in the rack in the cool- ing bath. THE BLOOD SAMPLE MUST BE MIXED WITH THE ANTICOAGULANT AND	Sample storage tubes For Assay: Precision pipettes – 50, 200 and 500 µl sizes Assay tubes – 12 x 75mm polystyrene tubes, round bottomed fother types of	<u>4</u>
18777 18 - ™ 1	ANTIPLATELET REAGENT MIXTURE AND <i>COOLED</i> AS RAPIDLY AS POSSIBLE AFTER COLLECTION. 6. Allow sufficient time for the sample to cool. Ensure that the ice remains in the cooling container throughout this time. 7. Centrifuge the sample at 1500-2000g and 2.4°C for 30 minules. Centrifuge the sample out within three hours after collection. 8. After centrifugation. carefully remove the top 0.5ml of plasme using a pipette with a disposable plastic tip. Transfer this sample to a separate	tubes are not recommended) Vortex mixer Centrifuge capable of 1000 <i>g</i> Absorbent towels Decantation racks Gamma scintillation counter Graduated 10ml glass pipette	
	labeled specimen tube. The sample for assay can be stored at room temperature for up to 24 hours or at 2-4°C for up to seven days. Storage for longer periods should be at -20°C for up to 4 weeks. If the sample is stored frozen, it should not be repeatedly frozen and thawed. Thaved samples should be mixed well prior to use. Interfering Substances Samples should be checked for radioactivity prior to assay since significant radioactivity may cause interference with the assay depending on the behavior of the radioactive material in the assay system. Significant activity may be pres-	Tube racks Distuiled water Distuiled water Details of Procedure 1 Number the assay tubes as shown in the protocol (Table 1). 2 Pipette duplicate 50µd aliquots of the reconstituted standards and unknowns into the protocol (See Note 1). 3 Pipette 200µd aliquots of the 0.1G sectum solution into each assay tube Pipette 200µd aliquots of the 0.1G sectum solution into each assay tube tube and vortex mix the tubes.	an an an an ann an an an an an an an an
			·

ORIGINAL PAGE IS OF POOR QUALITY +

_

5
2
ē
8
2

two standard curves and a total of 14 unknowns in duplicate on two occasions or any combination of standards and unknowns to a total of 48 tubes. A stan-Each β -TG RIA Kit contains sufficient reagents for the construction of a 5 point standard curve and the assay of 19 unknowns in duplicate on one occasion, or dard curve must be constructed on both occasions If a large number of samples are to he assayed on a single occasion, reagents from two or three kits of the same lot may be pooled before use. Only a single standard curve is required in this case

Description

- ¹³⁵ β -TG (human, freeze-dried) contains up to 2 μ Ci, 74kBq⁻¹³⁵ in 10ml solution after reconstitution.
- Anti-B-TG (rabbit, freeze-dried) contains antiserum sufficient to bind at least 40% of 0.5ng β -TG in 10ml of solution after reconstitution.
 - Ammonium Suifate Solution contains >25ml of a 3.3M solution of ammonium sulfate (sulphate)
- tains 0 5ml of a solution with *B*-TG concentrations of approximately 10, 20. eta-TG (human) in Bulter (freeze-dried) — after reconstitution, each vial con-Blood Sampling Tubes — contains anticoagulant and antiplatelet agents. 50, 100 and 225ng/ml. The exact values are stated on the vial labels.

Precautions

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES. Caution: Radioactive Material

Not for Internal or External Use in Humans or Animals.

radioactive material which should be handled with appropriate precautions in or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and a general license of the U.S. Nuclear Regulatory Commission, or of a state with which the Commission has entered into an agreement for the exercise of regulatory authority, or license by the Canadian Atomic Energy Control Board. This kit contains aboratory tests not involving internal or external administration of the material. This radioactive material may be received, acquired, possessed and used only by physicians, clinical laboratories, or hospitals and only for in vitro clinical or use and disposal.

instructions Relating to the Handling, Use and Storage of Radioactive Materials

Introduction Into Foods. Beverages, Cosmetics, Drugs, or Medicinals, or Into Products Manufactured for Commercial Distribution is Prohibited – Exempt Quantities Should Not Be Combined.

- Radioactive material should be stored in specially designated areas not normally accessible to unauthorized personnel. _
- Radioactive material should be used only by responsible persons in authorized areas. Care should be exercised to prevent ingestion or contact with the skin or clothing In the event that contact is made with radioactive material, The Si unit of radioactivity is the bacqueriel (Bq) equal to one distribution per second 1μ Cr = 37 the contaminated areas should be thoroughly washed with detergent.

No smoking, drinking, or eating should be allowed in areas where radioactive Pipetting of radioactive solutions must not be done by mouth

e

- materials are used
 - Hands should be washed after using radioactive materials
- Work should be carried out on a surface covered with absorbent materials. ف ی
- lammated materials disposed as radioactive waste. Contaminated surfaces Any spills of radioactive material should be cleaned immediately and all conshould be washed with a detergent

Reconstitution

Standards

the rubber stoppers and place them inverted on a clean surface Gently add 500 ± 10 µl of freshly distilled water onto the inside surface of each of the vials ules invert each vial gently to remove any particles of desiccute from the cap from the stoppers. Remove the tear-off aluminnum closures, carefully remove and replace the stoppers. Allow to dissolve at room temperature for two nun-Tap the vials of freeze-dried builter standards to distodge any large particles and to obtain a homogeneous solution. Avoid vigorous agitation and foaming

Antiserum and 1351

0.2ml distilled water. Ensure that homogeneous solutions are obtained by +I Reconstitute the treeze-dried antiserum and 1251 JJ-TG reagents with 10 gentle inversion Avoid vigorous agitation and loaming

Storage

The β -TG RIA Kit should be stored at 2-4°C.

The reagents may be stored at 2-4°C for up to one week after reconstriution The ammonium suitate solution can be stored at 2-4°C or at room temperature

2.5

Betore use on a second occasion, reagents should be allowed to come to room temperature and be mixed to ensure homogeneity

Availability

From stock.

Expiration

The expiration date is stated on the package and will normally be 3-6 weeks from the date of shipment.

Do not use the kit after the expiration date shown on the label

Preparation and Storage of Samples

Sample preparation is the critical stage in this assay. THE FOLLOWING PROTOCOL MUST BE STRICTLY ADHERED TO Improper sample preparation may result in elevated μ -TG levels. The user should be familiar with this protocol prior to collecting any samples for assay.

platelets during the blood clotting process, resulting in serum β -TG levels in the cause plateter refease can occur during the collection and processing of plasma samples, resulting in β.TG levels of 7-10μg/ml Platelet poor plasma processed Serum samples are not acceptable for this test because β -IG is released by $(0.25 \mu g/m l$) in addition, samples of normal plasma are unsuitable beso as to minimize platelet release during handling must be used for this assay raivge

Ĺ

and the second states and with the second second

;

I

•)

F1

7

1

Π

1

1

e

PRECEDING PAGE BLANK NOT FILMED

4/17/85

4.3 LEUKOCYTES

OUTLINE OF ISOLATION AND TESTING

13 Bags of Sedimented Buffy Coat (150.0 ml each, 9 in space, 9 on ground) 6 each in PVC-DEHP, PVC-TOTM, and Polyolefine Plastic



B & T Lymphocytes H.C. J.C.

LEUKUCYTES

Abstract

Buffy coat white cells will be isolated from whole blood and sedimented with Dextran. Leukocyte concentrates will be centrifuged, and the supernatants containing Dextran and platelets discarded. Nine bags of leukocytes will be reconstituted to 150 ml with buffered plasma and sent into space. An additional nine bags will be maintained on the ground for an equivalent storage period. All bags will be coded and randomized post-flight.

Polymorphonuclear cells without further purification will be examined by electron and fluorescence microscopy, and electronically for volume distribution. Phagocytic and respiratory burst activity also will be assayed.

Mononuclear cells will be separated from leukocyte concentrates with Ficoll-Hypaque. They will be examined without further purification by electron and fluorescence microscopy, electronically for size distributions and by flow cytometry for differential counts and T cell stimulation responses.

A. POLYMORPHONUCLEAR CELLS

1. Preparation of Leukocyte Concentrates.

Using the standard procedures for obtaining platelets, 75.0 ml of buffy coats will be harvested into a 300 ml transfer pack. Via a sterile connector 15.0 ml of sterile 2% Dextran will be added from a stock previously filtered through 0.44 u membrane filter. The bag contents are mixed briefly manually and the bag suspended upside down. Sedimentation of red cells is allowed to proceed for 45 minutes at room temperature. The red cells are carefully drained via the original sterile connector and discarded. The bag containing the leukocyte concentrate is centrifuged in a RC-3 Sorval centrifuge in a blood bank head HG at 4°C and 300 x g for 8-10 minutes. The supernatant fluid containing Dextran and platelets is expressed and discarded. The white cell pellet is resuspended to 150 ml in the bag with PBS buffered plasma (1:1).

Microfluorescence of Cytoplasm and Nuclei of Granulocytes.

The reaction of granulocytes with fluorescein diacetate (FDA) and ethidium bromide (EB) in Hank's balanced salt solution (HBSS) without calcium and magnesium is performed. About 1 x 10⁶ granulocytes are mixed in proportions of 0.05 to 0.25 ml with 0.50 ml of a mixture of fluorescein diacetate and ethidium bromide. Wet mounts are made at room temperature and the cells are viewed within one minute with an Olympus Vanox transmission microscope. A green exciter filter (G553) is used for identification of the granulocytes and fluorescence is viewed after switching to a UV exciter filter (Schott BG-12) and a blue barrier filter (Schott 06-530). Esterase activity in the cytoplasm of granulocytes is measured as the percent of cells showing green fluorescence of the fluorescein liberated from fluorescein diacetate. The percentage of cells with red fluorescence nuclei due to uptake of ethidium bromide is also recorded. Two hundred cells are counted.

The percentage of viable cells is obtained from the fraction of green cells to the total of green plus red fluorescent cells.

3. Determination of Size Distributions

Cell suspensions containing 100,000 to 150,000 granulocytes per ml are sized in a Coulter Model ZH with a C-1000 Channelyzer and an X-Y recorder. Prior to testing the granulocyte samples, the Coulter aperture is flushed with the buffered solution that is used to dilute the granulocytes. The sample volume is 0.1 ml. All solutions are rendered particle-free by filtration through both the 0.45 and the 0.20 micron filters in a Coulter filtration system or with 0.22 micron Falcon filters. A 7084 aperture is used. This is set at an amplification of 1, an aperature current of 4, an exclusion of 4, a lower threshold of 5, and an upper threshold of 99. The settings for the Channelyzer are: count control on External, count range 1 K, base channel threshold 1, and window width 100. Calibration is with 2.0 micron and 9.8 micron polystyrene particles. Granulocyte counts are made by integration of the number of cells within the granulocyte distribution. Granulocyte counts in whole blood and other samples are also made with a Coulter Counter Model F blood cell counter after lysis of red blood cells with Zap-Isoton. Red blood cell contamination of isolated granulocytes is counted as the difference in counts between the nontreated and Zap-Isoton treated aliquots at a 500 fold dilution.

4. Phagocytosis Assay.

a. Preparation of 1251 labeled Staph aureus.

Staphyloccus aureus (IgGsorb) is reconstituted with 10 ml cold distilled $H_2\bar{0}$ yielding a 10% (v/v) cell suspension in phosphate buffered saline (PBS) (150 mM NaCl, 40 mM phosphate pH 7.2 containing 0.05% azide). A 0.5 ml aliquot of the stock suspension is washed three times with 3 ml of PBS pH 7.0 at 2200 rpm for 10

minutes at 4°C and resuspended to 0.3 ml. Three hundred ul of washed IgGsorb is then added to a 10 x 100 glass tube which contains 50 ug of Iodogen evaporated on the vessel bottom and the original tube is washed with 100 ul buffered saline (pH 7.0) and combined. Five hundred uCi of ¹²⁵I in 100 ul of buffered saline is then added, gently rotated to initiate contact with the Iodogen pellet and allowed to incubate with occasional rotation on ice bucket under constant rocking. One half milliliter of 150 mM KI is added to stop the reaction. The labeled cell suspension is washed four times with 3 ml volumes of buffered saline (pH 7.0) and the cpm/ul of *Staph aureas* and the ¹²⁵I incorporated is determined.

b. Phagocytosis of Iodinated Staphyloccus aureus by Human Neutrophils.

In each sample to be tested 100 ul of a 20-fold diluted ¹²⁸I-Staph aureus suspension is added to 12 x 75 polypropylene tubes containing (in duplicate) 0.5 ml of fresh normal autologous, heterologous, or heat-inactivated serum. Controls will contain buffered saline plus Mg⁺⁺ and Ca⁺⁺. The suspensions are incubated at 37° C for 30 minutes with constant shaking to keep cells suspended. The opsonized particles are then washed 2X in ice cold buffered saline at 2200 rpm for 10 min at 4° C. After the final centrifugation, the cell pellets are resuspended to 0.4 ml in KRP pH 7.4 plus 5 mM glucose. Neutrophils (5 x 10⁶) are added to each tube and the final volume adjusted to 0.5 ml. The suspensions are mixed continuously for 20 minutes at 37° C. Ice cold buffered saline is added to stop ingestion and the cells are washed to remove free uningested particles (150 x g for 10 minutes at 4° C). The cells are then counted in a gamma counter. The phagocytes indices are determined as counts of radioactivity ingested per 10⁶ cells (cpm in normal sera cpm in PBS = cpm ingested by *Staph aureus*).

5. <u>Glucose Oxidation by Human Leukocytes</u>. <u>C1402 from 14C-1 labeled</u> glycose.

Human leukocytes from control and 5 day old test samples in a concentration of 1-1.5 x 10⁷ cells/ml will be suspended in 3 mls of Kreb's Ringer Phosphate pH 7.4 containing 1 mM glucose and 0.5 uCi of [1-1⁴C]glucose (approx. 50 mCi/mM). Then 1-2 ug of PMA or f-MLP (10⁶ M), both soluble stimuli for neutrophils will be added. The suspension will be incubated in a 25 ml side-arm erlenmeyer flask stoppered by a cap fitted with a cup containing 0.2 ml of hyamine hydroxide. After a 30 min incubation at 37°C in a shaking water bath, the mixture will be acidified by injection of 1 ml of 5 N H2SO4. The cell suspension will be equilibrated for an additional 30 min and then the cup removed and the contents assayed for radioactivity by liquid scintillation counting using a Beta Counter.

The amount of 14CO2 released from [1-14C]glucose will be expressed as the percent of (control) unstimulated leukocytes. Fresh preparations of purified human neutrophils and leukocyte rich plasma will also be included as controls for each test group.

6. <u>Electron Microscopy - Diana Ausprunk</u>.

Granulocytes and mononuclear cells will be studied by transmission and scanning electron microscopy following fixation in 2.5% glutaraldehyde buffered to pH 7.4 with sodium cacodylate. The samples will be transported to Boston at 4°C following post fixation in 1% osmium tetroxide in sodium cacodylate buffer.

a. <u>Scanning</u>: An aliquot of fixed leukocytes will be passed through a nucleopore filter. The filters with attached cells will be dehydrated in ethanol. Following deposition of gold-palladium the cells will be examined and photographed in a JSM-35 scanning microscope. This will provide a three-dimensional view of the shape and surface structure of the cells.

- b. <u>Transmission</u>: An aliquot of fixed cells will be dehydrated in increasing concentration of ethanol. The dehydrated cells are infiltrated with plastic resin, cut on an ultra microtome, placed on grids and stained. The stained sections are studied and photographed in a JEM 100 B electron microscope. This will provide a high resolution view of the structure of the cell cytoplasm.
- c. <u>Whole Mount Preparation</u>: A suspension of cells is spread on a Formvarcoated electron microscope grid and processed by critical point drying. The entire complement of organelles and cytoskeleton within each cell can be observed simultaneously. Any defects in cell spreading or cytoplasmic organization can be readily detected.

7. Cell shape analysis of granulocytes and mononuclear cells - W.M. Curby.

- a. <u>Electronic Analysis</u>: The shape of stored white cells will be analyzed using the Curby Biodetector (CBD). The ability of each cell type to maintain size and shape against known osmotic gradients over fixed time intervals will be monitored.
- b. <u>Photomicrographs</u> of wet slide preparations will be made of cells from each formed element population before and after the application of gradients.

Stop action photomicrographs of each cell population will be done using Xenon strobe, point source light illumination. The cells will be photographed in hanging drop slide preparations to obviate distortion associated with the adherence of a blood cell to a glass microscope slide.

B. MONONUCLEAR CELLS*

1. Isolation of Mononuclear Cells

Ten ml of leukocyte concentrate are overlaid onto 3 ml of Lymphoprep in a 15 ml conical tube and centrifuged at 1500 rpm (400 x g) for 30 min at 20°C. The top 5 ml is removed and discarded. The mononuclear cells are removed from the interface between Lymphoprep and saline and placed in a 15 ml test tube. The remaining fluid above the pellet is discarded. All cells are kept on ice (2°C).

2. Tests on Mononuclear Cells

Monocytes and lymphocytes in Lymphoprep isolates will be tested for size and membrane integrity by microfluorescence as described in A (Polymorphonuclear Cells). Lymphocyte differential counts will be made with an EPICS cell sorter using monoclonal antibodies to specific cell types. Helper-suppressor T cell ratios will be determined as well as the T cell stimulation response.

*Note: Protocols for mononuclear cell studies from collaborators, Carter, Taylor and Yunis are in discussion. WHITE CELL EXPERIMENT

MATERIALS AND EQUIPMENT

Processing Buffy Coat to isolate PMS and MNES Quantity per unit processed

MATERIALS:

2 Sampling site coupler (Fenwal 4C2405) 50 ml syringes (BD 5663) 1 box 18 g needles 1.5" x 18 g (100/box) (BD 5196) 2 Styrofoam rack for 15 ml and 50 ml T.T. DEXTRAN HÍ M.W. Sigma D-52S1 200,000-300,000 M.W. NaCl 10 ml pipettes 25/unit 2 pipette bulbs 15 and 50 ml polyproplene test tubes/conival tubes Lymphoprep 1.077 gr/ml 1 case (100 ml bottles) PBS pH 7.4 + glucose 3 liters glass microscope slides 75mm x 25 1 box glass microscope slips 22mm x 22m 1.5 thickness 1 box powder Fluorescein Diacetate powder Ethidium Bromide Pooled Plasma 3.0 liters 12 x 75 mm test tubes polyproplene with caps 2 cases Coulter accuvettes for cell counting&ZAP-oglobin II (Azide free) 2 boxes Coulter chart paper for C-1000 xy plotter Glass Pasteur pipette 5.75" disposable glass Plasma expressor and 3 hemostats & 1 pair scissors Isoton, Coulter, 1 cube

EQUIPHENT

Sorvall RC 3B refrigerated centrifuge with head (6 place) IEC CRU 5000 centrifuge with head #253 Refrigerator for blood storage 4-6°C range & reagents Fluorescence microscope w/Schott 0G-530 barrier and Schott BG-12 filters Water bath, shaking, 37°C Coulter ZH counter, C1000 channelyzer and X-Y plotter Vortex mixer Magnetic stirrer and magnetic stir bars Balance - triplebeam Vacuum and vacuum trap set up with small air pump Double beam spectrophotmeter, Cuvettes J isotope counter (gamma)

Leukocyte-81

SPACE, PLUMBING AND ELECTRICAL REQUIREMENTS FOR WHITE CELLS' STORAGE EXPERIMENT

Fabian Lionetti F. Wm. Luscinskas James Cohn

For Microscopy:

Water bath, 37°C, 20" x 15" x 6" approximately Olympic Vanox fluorescence microscope Illuminator

Bench space 2' x 6' x 2.5' = 14 sq' 2 Electrical outlets 110 V 1 sink, stainless, 20" x 10" x 10"

For Cell Counting and Volume Studies:

Coulter Model ZH Coulter Channelyzer Coulter X-Y Recorder Dade Diluter Vortex Mixer Coulter Filtration System

Bench Space 2' x 10' x 3' = 20 sq' 6 Electrical outlets 110V Coulter S+4 Diff. Counter 10 sq' LEUKOCYTES SUMMARY OF BENCH SPACE, PLUMBING AND ELECTRICAL REQUIREMENTS

Bench Space: 70 sq'

40 sq'Work 30 sq'For Equipment

- Desks: (2) 2' x 5' = 20 sq'
- Sinks: (2) 2' x 2' x 1.5'

Outlets: (16) 110 V

Special power requirements for Coulter S+4 Diff. Counter: An independent, protected circuit three wire outlet, 105-125 Vac. single phase input power, 20 A. 1500 W power consumption.

- 14.4 Electron Microscopy Protocol
- Part I. To be completed at on-site facility
 - A. Fixation of RBC
 - 1) 10 ml whole blood
 - 2) ghosts
 - B. Fixation of Platelets
 - 1) 5 ml platelet concentrate
 - 2) ADP-stimulated platelets
 - 3) serotonin-stimulated platelets
 - C. Fixation of Leukocytes
 - D. Preparation of Platelet Whole Mounts
 - E. Preparation of Leukocyte Whole Mounts

Part II. To be completed at Children's Hospital, Boston

- A. Transmission Electron Microscopy of RBC, Platelets and Leukocytes
 - 1) postfixation, dehydration, embedding
 - 2) sectioning and microscopy
- B. Scanning Electron Microscopy of RBC, Platelets and Leukocytes
- C. Freeze-fracture of RBC
- D. TEM & SEM of Whole Mounts of Platelets and Leukocytes

EQUIPMENT AND SUPPLIES FOR ELECTRON MICROSCOPY

Equipment:

Tabletop centrifuge

Top loading balance

Refrigerator

CO2 incubator

Phase contract microscope

Tank of 5% CO2: 95% air mixture

Supplies:

Pipettes (1, 5 & 10 ml)
Pasteur pipettes (6", plastic)
Conical test tubes (plastic w/caps) - 100 ml + 50 ml
Parafilm (roll)
Scissors
Timer (1 hr)
Racks for test tubes
Gold grids (formvar coated)
Petri dishes (35 mm, sterile, tissue culture plastic)
Forceps
Hemacytometer
Pocket calculator
Graduated cylinders (50 + 100 ml)
Reagent bottles

<u>Reagents</u>:

Distilled H₂O

Gluteraldehyde (10%, nitrogen sealed vials)

Sodium cacodylate

CaCl2

Sucrose

Dulbecco's PBS

Eagle's medium

30 müsm buffer

DETAILED EM PROTOCOL

Part I. To be completed at on-site facility

A. Fixation of RBC (9 samples of 5 ml each)

100 ml of 1% gluteraldehyde in Eagle's medium (pH 7.0) will be added to each 5 ml sample of RBC. After 30 min at room temperature, RBCs are centrifuged at room temperature and 100 rpm for 10 min. After removing the supernatant, the cells are washed in Eagle's medium, recentrifuged and resuspended in 100 ml of Eagle's medium containing 6% sucrose. The cells are stored at 4°C and transported to CHMC. Red cell ghosts (obtained from Dr. Jacobson) will be similarly prepared except that 30 m 0sm buffer replaces Eagle's medium.

B. Fixation of Platelets (9 samples of 2 ml each)

40 ml of 2.5% gluteraldehyde in 0.1 M sodium cacodylate buffer (pH 7.4) will be added to each 2 ml sample of platelets. After 30 min at room temperature, platelets are centrifuged at room temperature and 5000 g for 7 min. After removing the supernatant, the platelets are washed in 0.1 M sodium cacodylate buffer, recentrifuged and resuspended in 40 ml of cacodylate buffer containing 6% sucrose. The cells are stored at 4°C and transported to CHMC. Platelets stimulated by ADP and serotonin will be similarly prepared.

C. Fixation of Leukocytes (9 samples 3 ml each)

60 ml of 2.5% gluteraldehyde in 0.1 M sodium cacodylate buffer (pH 7.4) will be added to each 3 ml of leukocytes. After 30 min at room temperature, cells are centrifuged at room temperature, 100 rpm for 10 min. After removing the supernatant, the cells are washed in cacodylate buffer, recentrifuged and resuspended in 60 ml of cacodylate buffer containing 6% sucrose. The cells are stored at 4°C and transported to CHMC.

D. Preparation of Platelet Whole Mounts (9 samples of 1 ml each)

Gold electron microscopy grids will be prepared at CHMC in advance of arriving at Cape Kennedy by cleaning them in dilute nitric acid, coating the grids with a thin film of 0.7%. Formvar plastic and evaporating a thin layer of carbon onto the Formvar film. Platelets are washed free of plasma proteins and resuspended at a concentration of 50,000/ml in 0.9% NaCl. One drop of platelet suspension from each sample is placed onto each of four grids laying on a piece of Parafilm contained in a 35 mm Petri culture dish. The dishes are covered and placed in a himidified CO2 incubator (5% CO2: 95% air) at 37°C for 10-20 min. Each grid is removed from the incubator, rinsed in Dulbecco's phosphate buffered saline and placed in a vial containing 1 ml of 2.5% gluteraldehyde in 0.1 M sodium cacodylate buffer (pH 7.4). After fixation for 30 min at room temperature, the grids are rinsed in cacodylate buffer and then stored at 4°C in 1 ml of cacodylate buffer containing 6% sucrose.

E. Preparation of Leukocyte Whole Mounts (9 samples of 1 ml each)

Gold electron microscopy grids are prepared as described for platelet whole mounts. Leukocytes are diluted to a concentration of 200,000 cells/ml in 0.9% saline. A drop of cell suspension is placed onto each of four grids contained in a Petri dish as described for platelets. The covered dishes are placed in the 37°C incubator for 2 hr and then fixed and processed using identical methods described for the platelet whole mounts.

Part II. To be completed at Children's Hospital, Boston

A. Transmission Electron Microscopy of RBC, Platelets and Leukocytes (27 samples)

One ml of each sample of fixed RBCs, RBC ghosts, platelets and leukocytes will be centrifuged into a pellet and postfixed in buffered 1 ? osmium tetroxide (in Eagle's medium for RBC, 30 m0sm buffer for ghosts, 0.1 M sodium cacodylate buffer for platelets and leukocytes) for 60 min at 4°C. After washing 3X in distilled water, cells are dehydrated in increasing concentrations of ethanol (70%, 80%, 95%, 100%) for 5 min in each alcohol. The dehydrated cells are placed in 2 changes of propyline oxide (20 min each) and embedded in plastic resin (Polybed). Sections 1 mm thick are stained in methylene blue-azure II examined by light microscopy and photographed using Kodak Plus-X film. Thin sections are cut with diamond knives on a Reichert ultramicrotome, doubly stained with uranyl acetate and lead citrate and examined and photographed in a JEM 100 B electron microscope.

B. Scanning Electron Microscopy of RBC, Ghosts, Platelets and Leukocytes (27 samples)

0.25 ml of each sample of fixed RBCs, ghosts, platelet and leukocytes will be passed through a Nucleapore filter (0.45 um pore size) by pressure filtration. Filters with attached cells are dehydrated in ethanol and critical point dried from liquid CO₂. A thin layer gold-polledium is deposited on the cells in a sputter coater. After attaching the filters to specimen stubs, the cells will be examined and photographed in a JEUL JSM-35 scorning electron microscope.

C. Transmission and Scanning Electron Microscopy of Whole Mounts of Platelets and Leukocytes

Whole mounts of fixed platelets and leukocytes are postfixed for 5 min at 4°C 0.5% osmium tetroxide buffered to pH 7.4 with 0.1 M sodium cacodylate. The cells are dehydrated in ethanol and critical point dried from liquid CO2. A thin layer of carbon is evaporated over the cells. The pre-

parations will be studied, by transmission electron microscopy and photographed in a JEM 100 B electron microscope. The microscope is equipped with a goneiometer stage so that pairs of stereophotographs will be obtained by tilting the stage *6°-7° from the horizontal. After transmission electron microscopy, the same specimen will be coated with gold-palledium in a sputter coater and examined in the scanning electron microscope.

D. Freeze Fracture of RBCs and Ghosts

One ml aliquots of fixed RBCs and ghosts are centrifuged into a pellet, washed with Dulbecco's PBS and the pellet resuspended in 20% glycerol in PBS for 15 min. Glycerinated cells will be fractured in a Bolzers apparatus at -115°C and the fracture face will be replicated by evaporating platinum and carbon onto its surface. Replicas are cleaned by floating over methanol and chlorine bleach to digest away the cells. The replicas are transferred to copper grids and examined in the transmission electron microscope.

Time Frame for Completion of Electron Microscopy

Part I. Work to be completed on-site

1.	Fixation of all cells	2 hours
2.	Centrifugation and washing of cells	3 hours
3.	Resuspension of cells in storage buffer	2 hours
4.	Packing of specimens for transport to Boston	3 hours
	TOTAL	10 hours

Part II. Work to be completed at Children's Hospital, Boston

The Initial Blood Storage Experiment

NASA-IBSE

Protocols to be Performed at the Launching Site

Immunology

Introduction

The following protocol describes the tests to be done by UHMC staff at the Kennedy Space Center in the context of the NASA-IBSE (Initial Blood Storage Experiment). The protocol is divided into two sections, <u>I Materials</u> and <u>II Methods</u>. Section I will describe the <u>essential laboratory space</u> <u>required including</u> the <u>major equipment</u>. Some of the equipment already exists at KSC laboratories and is available for our use. Other equipment is provided by the manufacturer for the period of the experiment. The remainder is to be either brought in by UMMC staff or ordered to be sent to KSC directly. These items are indicated by *. In the Materials section there will be a list of <u>the laboratory supplies and the reagents and buffers</u> required for processing the RBC and platelet specimens for the immunological assays. Special reagents needed for individual tests will be given in the Methods section, under the heading of each special test. The II section will describe the following procedures to be conducted at KSC:

1.0 <u>Glycerolization</u> and <u>Freezing</u> of <u>RBC</u>

2.0 Processing Platelet Concentrate Samples

3.0 <u>Quantitative Antiglobulin Consumption</u> (QAC) for Platelet-Associated C3c and C3d

4.0 <u>QAC</u> for <u>Platelet-Associated</u> IgG

I. MATERIALS

LABORATORY SPACE AND MAJOR EQUIPMENT:

A -80°C freezer must be available in the laboratory at the site.

A 4°C cold room with sufficient space and accessible electrical outlet to accommodate a 0.5 amp, 110 V powered test tube agitator* (12'x9'x24').

A sink with 2 faucets for the attachment of two aspirator pumps. An RC-5 centrifuge (Sorval) with SS-34 rotor.

Bench space (approx. 15'x2.5') to accommodate the following:

2 IEC centr-7 centrifuges equipped with #210 rotors and carriers for 13 mm, 16 mm, and 30 mm diameter test tubes.

*2 vortex mixers.

An electronic cell counter (capable of counting platelets) and associated reagents and supplies.

Crushed ice supply.

Dry ice supply.

LABORATORY SUPPLIES:

*Adjustable FinnpipettesTH and tips.

*Test tube agitator.

*20 and 5 ml graduated syringes.

*4" and 250' roll of parafilm.

*250 ml beakers.

*2000 ml and 500 ml graduated cylinders.

*Plastic and glass transfer pipets.

*Insulated containers for sample transport.

*Electronic 3-channel timer.

*Test Tubes.
*12 x 75 mm glass tubes.
*16 x IDo mm siliconized glass tubes (vacutainer red tops).
*30 x 100 mm plastic centrifuge tubes.
*2 ml and 13 ml capped plastic freezing tubes.
*Assorted reagent containers.
*Tape, markers, etc. for labeling.
*Racks for test tubes.
*500 ml volumetric flasks.

REAGENTS, BUFFERS:

100 1 reagent Hz0
*Dry chemicals, pre-weighed, for making the following:
0.9 g/dl NaCl with 10 mM phosphate, 0.1 g/dl NaN3, pH 7.2 (PBS).
PBS containing 0.15 g/dl Na2EDTA (PBSA-EDTA).

GlycerolyteTH 57 solution.

II. METHODS

1.0 PROCESSING WHOLE BLOOD SAMPLES FOR SHIPMENT TO UMMC

Reference: I.O. Szymanski, J.M. Harper, P.R. Odgren, C.R. Valeria. Freezing red blood cells prepared for quality control of antiglobulin sera. Transfusion 21:498-501, 1981.

- 1.1 Sample required: 30 ml whole blood
- 1.2 Procedure

1.2.1 The sample(s) is divided into two 50 ml centrifuge tubes and washed 3 times with PBS and the RBC packed after final wash.
1.2.2 After the final wash add to the packed cells a volume of glycerolyte equal to 0.4 times the volume of the packed cells. The

glycerolyte is added slowly and with constant gentle mixing on vortex mixer.

1.2.3 The RBC and glycerolyte are allowed to equilibrate at room temperature for 10 min.

1.2.4 Slowly and with constant gentle mixing on vortex mixer, a volume of glycerolyte equal to 1.6 times the original volume of the packed RBC is added.

1.2.5 The contents of both tubes are then pooled, mixed and aliquoted as follows:

<u>Two</u> aliquots of at least 1.5 ml each (for detection of IgG and complement by direct agglutination tests)

Four aliquots of at least 6 ml each (for quantitation of C3 and IgG on RBC membrane)

1.2.6 10 min. after the final addition of glycerolyte, the samples are frozen at -80°C.

1.2.7 Samples are stored in -80°C freezer until they are transported to UMMC, at which time they are packed in an insulated container with dry ice, and shipped.

Man hours for the processing of whole blood sample: one person can handle 18 samples in approximately 6 hours.

2.0 PROCESSING PLATELET CONCENTRATE SAMPLES

2.1 Sample required: A 10 ml sample is obtained on the day of launch, and 18 5 ml samples of platelet concentrate are obtained on the day of return of the satellite. A larger sample is required initially since platelet-bound complement is low in that sample.

2.2 Platelets are pelleted by centrifugation at 2000 RCF for 7 minutes, the supernatant plasma decanted, and they are then washed 6 times with 5 ml PBSA-EDTA.

2.3 After the final wash, platelets are suspended with a volume of PBSA-EDTA sufficient to yield a concentration of approximately 1 x 10^6 /microliter in day 0 sample and approximately 5 x 10^5 /microliter in day 5 samples. The exact platelet count in the washed platelet suspension is determined using the cell counter.

3.0 QUANTITATIVE ANTIGLOBULIN COMSUMPTION TEST (QAC) FOR PLATELET-ASSOCIATED C3c AND C3d

Reference: I.O. Szymanski, R.E. Swanton, and P.R. Odgren. Quantitation of the third component of complement on stored red cells. *Transfusion* 24:194-197, 1984.

3.1 Additional reagents required:

3.1.1 <u>Anti-C3 (C3c-specificf)</u>: Goat anti-human C3, IgG fraction, obtained from Atlantic Antibodies (Scarborough, ME). For this assay, it is diluted 1:100,000 in PBSA, and has been given lot #C3(G).

3.1.2 <u>Anti-C3d</u>: Goat anti-human C3d, obtained from Netherlands Red Cross (Amsterdam). Heteroagglutinins have been removed from this serum by absorption with bromelin-treated RBC. For this assay, it is diluted 1:5,000 in PBSA and has been given lot #17. 3.1.3 <u>Zymosan Standards</u> (ZyC3) consist of zymosan A particles which have been coated with C3, washed, and diluted in PBSA. This ZyC3 has been analyzed for C3 quantity as described before and sets of aliquots have been prepared for both C3c and C3d assays. Each set consists of ten 1.0 ml aliquots, each containing a dif-

ferent concentration of ZyC3. For C3d standardization, the C3 content in aliquots range from 1.81 x 10^{10} molecules/ml through 7.24 x 10^{11} mol/ml. For C3c standardization, the C3 range is from 8.06 x 10^9 mol/ml through 1.6 x 10^{11} mol/ml.

3.2 Procedure

3.2.1 Sample dilution:

3.2.1.1 Day O suspension is tested both undiluted and diluted 1/2 in PBSA-EDTA.

3.2.1.2 Day 5 suspensions are tested undiluted and diluted 1/2 and 1/5 in PBSA-EDTA.

3.2.2 Sample pipetting: 1 ml of each platelet suspension and required dilution(s) of ZyC3 is dispensed into a 2 ml tube in triplicate. One set is pipetted for C3c, another one for C3d. 3.2.3 OZ neutralization controls: 1 ml PBSA is dispensed into 20 2 ml tubes.

3.2.4 Antibody pipetting: To each tube of the C3 set (i.e., including dilutions of various platelets, three sets of ZyC3 standards, and 10 PBSA controls) exactly 1 ml of anti-C3c will be added; to each tube of the C3d set exactly 1 ml of anti-C3d will be added.

3.2.5 The tubes are all tightly capped, placed on test tube agitator, and allowed to incubate overnight (18-24 hrs.) at 4°C.

3.2.6 The incubated samples are centrifuged for 30 minutes at 12,000 RCF and the supernates transferred to appropriately labeled 2 ml freezing tubes.

3.2.7 The supernates are frozen at -80°C and stored until shipped to UMMC in insulated containers with dry ice.

4.0 WACT FOR PLATELET-ASSOCIATED IGG

<u>Reference</u>: N. Levitan, R. Tano, I. Szymanski. An autoanalyzer test for the quantitation of platelet-associated IgG. Submitted for publication.

4.1 Additional reagents required:

*4.1.1 <u>Anti-IgG (Fc specific)</u>, obtained from Atlantic Antibodies (Scarborough, ME) was absorbed with bromelin-treated RBC to remove heteroagglutinins and given BBR Lot #2. For this assay it is diluted 1:1,000,000 with PBSA.

*4.1.2 <u>IgG Standards</u> consist of Beckman ICS Calibrator I serum (Beckman Instruments) which has been calibrated against a World Health Organization IgG Standard. It is diluted with PBSA (in two-fold fashion) to a minimum of 8 dilutions, each 1.0 ml in volume. The concentration of IgG in these 8 standards is to range from 31.8 ng/ml to 0.24 ng/ml. Four "0"% neutralization targets are prepared also by dispensing 1.0 ml of PBSA into each of four tubes.

*4.1.3 <u>IgG Control</u> consists of Beckman ICS Calibrator serum which is diluted in PBSA (in two-fold fashion) to 4 dilutions, each 1.0 ml in volume. The concentration of IgG in these controls is to range from 14.2 to 1.77 ng/ml.

4.1.4 <u>Platelet Dilution</u>: Platelet suspension (from 5.3 above) is diluted to contain approximately 200,000 platelets/ml in PBSA-EDTA.

4.2 Procedure

4.2.1 <u>Anti-IgG neutralization by the IgG Standards</u>: To each of the 1.0 ml standards (from 7.1.2 above), add 1.0 ml of the 1:1,000,000 dilution of anti-IgG (from 7.1.1 above).

4.2.2 <u>Anti-IgG neutralization by Platelets</u>: Set up the consumption of anti-IgG by testing four concentrations of the platelet suspension (7.1.4 from above) as follows:

Dispense 1.0 ml of PBSA into each of four 12 x 75 mm tubes. To one of each of the tubes add 10, 20, 30 and 40 ul of the platelet suspension. Next add 1.0 ml of the diluted anti-IgG reagent (4.1.1 from above) and mix thoroughly. (<u>Since</u> <u>these are set up in duplicate, the final number of tubes per</u> <u>test is eight</u>.) (It has been found that with such small volumes of the platelet suspensions it is not necessary to take into account the platelet volume when making these dilutions.)

4.2.3 <u>Anti-IgG neutralization by the IgG controls</u>: To each of the 1.0 ml controls (from 4.1.3 above), add 1.0 ml of the diluted anti-IgG reagent (4.1.1 from above).

4.2.4 After the above have been set up, all tubes are stoppered and incubated at 4°C for at least 18 hr with constant, gentle agitation on the tube shaker (2.1). It is important that the unknowns be incubated simultaneously with the standards. After incubation, the tubes are centrifuged at 2000 x g for 10 min and the supernates transferred to the 2.0 ml plastic freezing tubes, capped tightly, and placed in the -80°C freezer. These are eventually transported in the frozen state to UMMC where the assay will be completed.

PRECEDING PAGE PLANK NOT FILMED

Immunology-9

TEST TO BE PERFORMED FOR NASA AT THE UNIVERSITY OF MASSACHUSETTS MEDICAL CENTER

INTRODUCTION:

The purpose of these studies will be to measure red blood cell (RBC) and platelet associated immunoglobulin and complement on samples of human blood which have been stored in a micro-gravity environment and on the Earth's surface. This is done with the intent of establishing if deposition of these plasma proteins increases (or decreases) upon the storage of human blood in space. A sensitive, automated <u>direct agglutination assay</u> will be used to determine relative amounts of immunoglobulin (Ig)G and the complement fragments, C3a, C3c and C3d on intact RBC. Precise quantitation of IgG, C3c and C3d on intact RBC and platelets will be made with the <u>quantitative antiglobulin consumption test (QACT)</u>. In addition, the QACT for IgG will be performed on isolated RBC stroma and lysates. Both of these assays use the Autoanalyzer (Technicon Corp., Tarrytown, NY) and are described in detail below.

1.0 <u>The automated direct agglutination assay</u>: This assay detects the presence of specific protein(s) on the surface of RBC and is based on the phenomenon of agglutination of said RBC by specific antibody. The results, expressed in percent agglutination, are proportional to the amount of protein of the RBC. Therefore, it is possible to express relative concentrations of cell-bound protein among different RBC populations. This assay will be performed as described in detail elsewhere (1). Briefly, washed RBC, suspended to a 10% hematocrit in Ficoll and 0.9% NaCl are introduced continuously into an Autoanalyzer manifold together with 1.0% polyvinylpyrrolidone K-90 (PVP) in 0.9% NaCl while appropriately diluted antibodies (to each of the four proteins) are sampled discontinuously. RBC bearing the specific protein aglutinate in the presence of corresponding antibody. The agglutinated RBC are
removed by gravity and the remaining unagglutinated RBC hemolyzed with Triton X-100 and water. The ODsso of the lysate is recorded and the percent agglutination which reflects directly the amount of protein on the RBC membrane is then calculated.

1.1 Materials.

1.1.1 Phosphate buffered saline (PBS): 0.9% NaCl, 0.1% NaN3, buffered to pH 7.2 w/0.01 M sodium phosphate.

1.1.2 Manifold cleaning solution: 2.5% urea in 0.05% NaOH in H2O.

1.1.3 0.5% Ficoll (Pharmacia) in PBS.

1.1.4 0.5% bovine serum albumin (BSA) in PBS.

1.1.5 1.0% polyvinylpyrrolidone K-90 (PVP) in PBS containing 2 drops of Tween-20 (Technicon) per 50 ml of solution.

1.1.6 Antisera to:

IgG (Fc specific) (Atlantic Antibodies, Scarborough, ME) (diluted 1:5000) C3a (Cappel) (diluted 1:200) C3c (Atlantic Antibodies) (diluted 1:2000) C3d (Neatherlands Red Cross) (diluted 1:200)

1.1.7 Lysing solution: 1.0% Triton X-100 (Sigma Chemical Co.,
St. Louis, M0) in H20 with 0.1% NaN3.
1.1.8 The glassware, instruments and facilities for performing the assay are available at the UMMC Blood Bank Research Laboratories.

1.2 Samples: Samples of glycerolyzed RBC will be thawed and deglycerolyzed (Appendix). The recovered RBC will be washed in PBS and suspended to a 10.0% hematocrit in 0.5% Ficoll solution.

1.3 Procedure:

1.3.1 An Autoanalyzer (Technicon, Tarrytown, NY) is outfitted with a sampler, tubing, manifold, reservoirs, pump, and colorimeter as shown in Figure 1.

1.3.2 The reservoirs for PBS and Triton are filled.

1.3.3 The platen on the proportionating pump is engaged and the sampling probes placed in the cleaning solution.

1.3.4 The pump, colorimeter and recorder are turned on.

1.3.5 After 20-30 min the probes are placed in PBS containing 2 drops of Tween-20 per liter of solution for 20-30 additional min.

1.3.6 The electronic scale is adjusted to read between 00 and 100T(T).

1.3.7 The antibody probe is placed in the BSA solution. The Analyzer is allowed to run for 10-20 min. The appropriate sampling probes are placed in the RBC suspension (which is agitated continuously with a magnetic stirrer and always maintained in an ice bath) and PVP solution. The antibody probe is engaged to the sampler and sampling is begun. The reaction products will reach the photo cell of the colorimeter in approximately 25 min.

1.3.8 After the fluids have reached the colorimeter, the baseline is adjusted with the "baseline" knob to 80-90%T.

1.3.9 After the baseline is established, sampling of the antibody solu-

1.3.10 The change in optical density from the baseline for each sample is determined from the chart paper and the % agglutination calculated from the following formula:

> %T of baseline-%T of the sample X 100 = %T of the baseline % agglutination

The percent agglutination is proportional to the number of molecules of antigen on the surface of the RBC and provides a measure of the relative quantities of membrane bound protein among different cell populations.

2.0 <u>The quantitative antiglobulin consumption test (QACT)</u>: This assay quantitates the molecules of a specific protein (IgG, C3c ro C3d) that accumulate on or within RBC and platelets. In this test varying amounts of test sample or known standards are incubated with a known amount of specific antibody. Thereafter, the antibody remaining in the supernate is measured by direct agglutination of indicator RBC coated with the specific protein. The percent neutralization is calculated and dose response curves constructed. The total number of molecules of protein per cell is then determined by using the dose response curves, Avogadro's number and the known molecular weight of the protein.

- 2.1 Materials:
 - 2.1.1 All of the material described in section 1.1.
 - 2.1.2 RBC products to be tested:

2.1.2.1 Intact RBC: RBC samples received will be deglycerolyzed (Appendix), washed three times with PBS, resuspended to about 50% hematrocrit and the total cell count determined electronically with a Coulter Counter (Coulter Electronics, Hialeah, FL).
2.1.2.2 Lysates for internal IgG determinations: For each 1.0 ml of the 50% RBC suspension, 0.5 ml of a 0.5% solution of digitonin in PBS is added while gently mixing on a Vortex mixter. 1.5 ml of PBS is added and after 2 min the mixture is centrifuged at 2000 x g for 7 min and about 1.0 ml of the supernatant hemolysate is transferred to dialysis tubing (12,000 to 14,000 molecular weight)

sieve) and dialysed against 250 volumes of PBS for 1 hr. (The pellet is saved for the preparation of the stroma as outlined below). The hemoglobin concentration of the lysate is determined on a Coulter Counter and kept at 4°C until tested for IgG. 2.1.2.3 RBC stroma: The pellet from step 2.1.2.1 is washed with PBS five times or until there is no visible hemoglobin and resus-

pended in an equal volume of PBS. A 100 ul aliquot is mixed with 2.0 ml of 1.0 x sodium dodecylsulfate in HzO for 10 min. The OD of the solution at 260 and 230m nm is obtained and the total stromal protein determined from the following formula:

(21(1.55(0Dzso)-0.77(0Dzoo) = mg protein/ml stromal suspension

(<8.1 mg/ml = 10¹⁰ RBCm stroma>>

2.1.3 Indicator, IgG coated RBC: Approximately 250 ml of packed group O Rho positive RBC are combined with 60 ml of anti-Rho (Ortho Diagnostics, Raritan, NJ). After 30 min at 37°C the RBC are glycerolyzed (2), aliquoted and sotred at -80°C.

2.1.4 Indicator, C3-coated RBC: These indicator RBC are prepared as previously described (3). These are also glycerolyzed and stored at -80°C.

2.1.5 IgG soluble standard: Beckman ICS Calibrator serum (Beckman Instruments, Irvine, CA) with known concentration of IgG is diluted with PBS to a concentration of approximately 30 ng/ml. This stock solution is aliquoted and stored at -80°C. Working standards (to range in concentration from 30 to 0.025 ng/ml are prepared by serially diluting the stock solution.

2.2 Procedure for the QACT:

2.2.1 The first step of this assy consists of combining different amounts of the antigens (intact RBC, lysates, soluble IgG standards and ZyC3) with a known amount of antiserum according to the following table:

Antigen	Volumes of the and Anti-IgG 4 ml (1:2M)	tisera Anti-C3c O.5 ml (1:50K)	Anti-C3d 0.5 ml (1:4.5K)
	20*	_	-
SUZ KBC	40	-	-
15% RBC	-	588	588
7.5% RBC	-	541	541
Stroma	200 100	-	-
Lysate	20 40	-	-
IgG standards	20	-	-
ZyC3 standards	-	5000	500

*Volume of antigen in ul.

These mixtures are incubated ovenight at 4°C with gentle agitation whereafter they are centrifuged and the supernatants tested for remaining antibody activity against the appropriate indicator RBC by the direct agglutination test described in 1.0.

NOTE ON PLATELET WACT: The incubation of the platelets with the antibodies will have been done on NASA site and the frozen supernatants containing the remaining antibody will be tested at UMMC.

2.2.2 Calculation of results: The percent agglutination is calculated and a dose-response curve depicting the percent neutralization of antibody by known standard is constructed. Molecules per cell are determined by using the dose response curves, Avogadro's number and the molecular weight of the given antigen (IgG, C3c, or C3d).

Appendix:

Deglycerolization of frozen RBC:

1. Frozen RBC are thawed at 37°C.

- 2. 12% NaCl is added* in a volume equal to 25% of the glycerolyzed sample and allowed to equilibrate at room temperature for 5-10 min.
- 3. 1.6% NaCl is added in a volume equal to the glycerolyzed sample and allowed to equilibrate for at least 2 min.

4. The sample is centrifuged and the supernatant decanted.

- 5. 1.6% NaCl is added as in step 3.
- 6. Step 4 is repeated.
- 7. PBS is added in a volume equal to the original glycerolyzed sample volume and allowed to equilibrate for at least 2 min.
- 8. Step 4 is repeated.
- 9. Step 7 is repeated.
- 10. Step 4 is repeated.
- 11. The RBC are resuspended in the desired diluent.

*All solutions are added slowly and with constant stirring with Vortex stirrer.

References cited:

- Szymanski, I.O., Huff, S.R., Delsignore, R.: An autoanalyzer test to determine immunoglobulin class and IgG subclass of blood group antibodies. Transfusion <u>22</u>:90 (1982).
- <u>Technical Manual of the American Association of Blood Banks</u>. J.B.
 Lippincott, Philadelphia, pp. 60-61 (1981).
- 3. Chaplin, H. Jr., Freedman, J., Massey, A., Monroe, M.C.: Characterization of red blood cells strongly coated <u>in vitro</u> by C3 via the alternative pathway. Transfusion <u>20</u>:254 (1980).
- 4. Szymanski, I.O., Odgren, P.R.: Measurements of fragments of the third component of human complement on erythrocytes by a new immunochemical method. Vox Sang. <u>46</u>:9 (1984).

Figure 1.

INITIAL BLOOD STORAGE EXPERIMENT

Biophysics Protocol

All data taking to be completed at on-site facility.

- A. Total Cellular Count
 - 1. From each red blood cell sample, a 10^{-2} ml aliquot will be added to a 10.0 ml of 0.9% w/v particle free saline and mixed by capping and inversion ten times. From this stock solution, 1.0 ml will be taken and added to 9 ml of the 0.9% particle free saline solution. This final concentration will be read directly on the Curby Biodetecter (CBD) using a short tunnel 5 x 10^{-2} mm diameter aperture. Cell counts and cell shape analysis will be made on 5 x 10^{-2} ml volumes of the final concentration fluid.
 - 2. White blood cell counts will be made on a final solution concentration containing 10^{-2} ml of each white blood cell sample in 10.0 ml of 0.9% particle free saline. Cell counts and cell shape analysis will be made on 5 x 10^{-2} ml volumes of the final concentration fluid.
 - 3. Platelet counts concentration will be arrived at by first making a stock solution containing 10⁻² ml of the platelet preparation in 10 ml of 0.9% particle free saline. 2 ml of the stock solution is then added to 3 ml of the 0.9% saline to make the final counting solution. All pipetting will be done using automatic microliter adjustable pipets and certified pipet tips.
 - 4. For all cell counts, the CBD will be set to count one peak per generated pulse. The size will be determined by placing the add one count in the memory storage channel equal to the maximum pulse height within the pulse. All pulse height measurements will be made when a negative pressure equivalent of 7.0 inches of mercury exists across the aperture tunnel of the 5 x 10⁻² mm aperture.
 - B. Cell Shape Analysis
 - Pulse shape analysis will be done using 1 inch of mercury pressure differential across the aperture. Aliquots of each sample will be drawn through the aperture, and in sequence, reading first the peak height and then the pulse height generated by each cell at 4, 8, 16 and 32 microseconds after the initial rise of the voltage pulse. All samples will be normalized by allowing the cell frequency distribution to be accumulated in the computer memory until the highest peak in the storage channel vs. frequency distribution reaches 100 counts full scale.
 - Distributions will be plotted by the CBD plotter and filed for analysis.

C. Red Blood Cell Action to Applied Unit Stress

- Pulse height data will be collected as in Section A1. Immediately after the initial reading 1 ml of particle free, pyrogen free distilled water will be added to the sample and peak readings will be repeated once every 2 minutes until a peak shift of more than 10 channels is observed between subsequent readings.
- The cell counts and size frequency data will be recorded and printed out by the CBD.

D. Photomicrographs

- Photomicrographs of each sample will be taken. Tungsten source light and Xenon source light will be used.
- 2. Two photomicrographs of each sample of red blood cells will be made using a Nikon Model S-Ke Microscope, an anachromatic dark contrast optical combination, viz., a 100X DLL achromatic objective and a matching 100 annular diaphragm in the focal plane of the condenser. Dark-low-low phase contrast images will be recorded on 35 mm, 1000 ASA speed print film. The field coverage film frame indexing will be done using a 10⁻² mm etched spacing stage micrometer. Red blood cells will also be photographed using a 100X plan apochromat lens and direct lighting under similar tungsten light source conditions.
- 3. A slide for each sample will be prepared by dropping exactly 10⁻² ml of sample onto a cleaned microscope slide and tipping a 1 cm x 1 cm cover slip onto the sample so that the center of the cover slip lands on the center of the drop. A drop of immersion oil having a refractive index of 1.515 will be centered on top of the cover slip. No oil will be used to fill in below the microscope slide and the microscope condenser.
- 4. Samples of platelets and white blood cells will be applied to make up microscope slide preparations in the manner described in D3. DLL phase contract pictures will be made from these samples.
- 5. Red blood cell samples will also be prepared for hanging drop stopaction photomicrographs by placing a 10⁻² ml of a one in ten diluted drop of each sample on a 1 cm x 1 cm cover slip and carefully inverting the cover slip and placing it over the center of a depression slide. The cover slip will be waxed in place.
- 6. Stop action (100X objective) photomicrographs will be made using E. Leitz Microscope Equipment and a Leica "Microibso" unit with a 1/3 lens to reform the image on the camera film plane. The light source will be a Xenon Strobe Unit which is centered along the optical axis of the microscope 15 centimeters below the objective focal plane of the microscope. Five centimeters above the strobe light bulb a diaphragm having a 4 mm opening is placed centered on

the optical axis parallel to the film plane. (This diaphragm effectively forms the point for formation of a new wave front.) No condenser is needed. A columnated tungsten light source and 45° removable mirror is placed above the diaphragm when a focusing light is needed. The microscope is focused onto the center of the hanging drop of sample, the mirror is removed and a flash picture is made of the cells drifting in the drop.

7. Color print film (35 mm) will be ASA matched to the optical system. As before a 10⁻² mm spacing stage micrometer will be used for fram identification and size indexing on the films.

- 15.0 ACQUISITION OF DATA
- 16.0 ANALYSIS OF DATA
- 17.0 SAFETY PROCEDURES
- 17.1 Dr. Fabian Lionetti from The Center for Blood Research, is the Safety Officier.

.

INITIAL BLOOD STORAGE EXPERIMENT

,

.....

ł

Schedules for Blood Collection, Component Preparations and Distribution

1.0 The purpose of this document is to define the timelines and responsibility for IBSE personnel to collect, prepare, pool and distribute whole blood and its components in preparation for loading into flight hardware and shuttle launch. The activity described in this document is detailed in Annex 8 of the Payload Integration Plan.

011140

2.0 Manpower Requirements

M S Jacobson		LHML
n Van Pelt		CBR
		CHMC
S. Kevy		CBR
F. Lionetti		CBD
T. Curran	•	
C. Smith		CEBB
A. Scanlon		CHMC

3.0 Schedule

. •

The activity of blood collection, component preparation and sample distribution is scheduled in relation to the shuttle launch. The schedule is developed on the basis of a launch on December 18, 1985 at 7:00 a.m. The activity begins at twenty-five hours prior to launch (T-25, 6:00 a.m., 12/17/85). The schedule of the activity is as follows:

See following page for table

4.0 Activity Record

The activity related to blood collection, platelet and leukocyte preparation will be entered on the attached record sheets.

			8 - + - · · · · · · · ·	Personnei Responsible
Code		Time	Activity	
3.1	T-25	(6:00 a.m. 12/17)	Begin blood collection	May Jacobson
3.2	T-24.5	(6:30 a.m. 12/17)	Begin platelet and leukocyte preparation	H
3.3	T-22.5	(8:30 a.m. 12/17)	Begin whole blood pooling	u
3.4	T-22	(9:00 a.m. 12/17)	Completion of blood collection	
3.5	T-20	(11:00 a.m. 12/17)	Completion of whole blood pool	11
3.6	T-20	(11:00 a.m. 12/17)	Completion of individual leukocyte units	64
3.7	T-20	(11:00 a.m. 12/17)	Begin distribution of whole blood	**
3.8	T-20	(11:00 a.m. 12/17)	Begin leukocyte pooling & distribution	11
3.9	T-19.5	(11:30 a.m. 12/17)	Completion of individual platelet units	
3.10	T-19.5	(11:30 a.m. 12/17)	Begin platelet pooling & distribution	u
3.11	T-19	(12:00 p.m. 12/17)	Completion of distribution of whole blood, leukocytes, platelets into designated ba	" gs
3.12	T-19	(12:00 p.m 12/17)	Begin packing blood bags into 8 boxes for transportation to KSC	n
3.13	T-18.5	. (12:30 p.m. 12/17)	Completion of packing blood bags	61
2 14	T-18	(1:00 p.m. 12/17)	Departure of blood from CFBB	11
3 15	T-17	(2:00 p.m. 12/17)	Arrival of blood at KSC	**
3.12	,-,, T_17	(2:00 p.m. 12/17)	Begin loading of dewars	**
5.10	1 - 71	, · · ·		Closeout

•

3.17 Continuation of pre-launch activity as described in "Prelaunch Closeout Procedures" (ADL Reference C-53282-45).

٠

ACTIVITY RECORD - WHOLE BLOOD COLLECTION

. .

UNITS COLLECTION	TIME STARTED	TIME COMPLETED	OPERATOR
1ST DONOR	•		:
2ND DONOR			
3RD DONOR			
4TH DONOR			
5TH DONOR .			
6TH DONOR			
7TH DONOR			
8TH DONOR			
9TH DONOR			
10TH DONOR			
11TH DONOR			
12TH DONOR			
13TH DONOR			
14TH DONOR			
15TH DONOR			
16TH DONOR			
17TH DONOR			
18TH DONOR ·			
19TH DONOR			
20TH DONOR			
21ST DONOR			
22ND DONOR			
23RD DONOR			
24TH DONOR			
25TH DONOR			

ACTIVITY RECORD - WHOLE BLOOD COLLECTION

UNITS COLLECTION	TIME STARTED	TIME COMPLETED	OPERATOR
26TH DONOR			
27TH DONOR			
28TH DONOR			
29TH DONOR		•	
30TH DONOR			
31ST DONOR			
32ND DONOR			
33RD DONOR			
34TH DONOR			
35TH DONOR			
36TH DONOR			
37TH DONOR			
38TH DONOR			
39TH DONOR			
40TH DONOR			

WHOLE BLOOD POOLING AND DISTRIBUTION

	TIME	TIME	OPERATOR
ACTIVITY	SIARIED	0014 0019	
POOLING OF FIRST 6 UNITS OF WHOLE			
BLOOD (FIRST POOL)			
MIXING 1ST			
MIXING 2ND			
MIXING 3RD			
MIXING 4TH			
MIXING 5TH			
MIXING 6TH			
POOLING OF SECOND 6 UNITS OF WHOLE			

ż

BLOOD (SECOND POOL)

MIXING	1ST
MIXING	2ND
MIXING	3RD
MIXING	4TH
MIXING	5th
MIXING	бтң

POOLING OF FIRST AND SECOND WHOLE BLOOD POOL

MIXING 1ST
MIXING 2ND
MIXING 3RD
MIXING 4TH
MIXING 5TH
MIXING 6TH

•

PLATELET PREPARATION

	17 ·	TIME STARTED	TIME COMPLETED	STATUS OF PREPARATION	RESPONBILE PERSONNEL
ACTIVIT PREPARA INDIVID	TION OF UAL UNITS				
UNIT	PRP Preparation				
	Packed platelets				
	Resting				
	Rotation				
UNIT	PRP Preparation				
	Packed platelets				
	Resting				
	Rotation	-			
UNIT	PRP Preparation				
	Packed platelets				
	Resting				
	Rotation				
UNIT	PRP Preparation	•			
	Packed platelets				
	Resting				
	Rotation				
UNIT	PRP Preparation				
	Packed platelets				
	Resting				
	Rotation				
UNIT	PRP Preparation				
	Packed platelets			•	
	Resting				
	Rotation ·				

+ •

LEUKOCYTE PREPARATION

.

	TIME STARTED	TIME COMPLETED	STATUS OF PREPARATION	RESPONBILE PERSONNEL
ACTIVITY DECRAPATION OF				
INDIVIDUAL UNITS				
UNIT 1				
UNIT 2				
UNIT 3				
UNIT 4				
UNIT 5				
UNIT 6				
UNIT 7				
UNIT 8				
UNIT 9				
UNIT 10				
UNIT 11				
UNIT 12				
UNIT 13				
UNIT 14				
UNIT 15				
UNIT 16				
UNIT 17				
UNIT 18				
UNIT 19				
UNIT 20				
UNIT 21	•			
UNIT 22				
UNIT 23				
UNIT 24				

· · · · ·

INITIAL BLOOD STORAGE EXPERIMENT PRE-LAUNCH CLOSEOUT PROCEDURES

(Revision 1.0)

submitted to

The Center for Blood Research Boston, Massachusetts

submitted by

Arthur D. Little, Inc. Cambridge, Massachusetts

December 31, 1985

ADL Reference C-53282-45

A Arthur D Little, Inc.

.

TABLE OF CONTENTS

- ----

-

1.0	Purpose	-
2.0	Scope	Ţ
3.0	Manpower Requirements	1
4.0	Schedule	2

ii

Arthur D. Little, Inc.

LIST OF FIGURES

		Page
1	IBSE Organization	Ĵ
2	Closeout Schedule	ú.
3	Loading/Closeout Sequence	5

iii

Arthur D. Little, Inc.

.

1.0 PURPOSE

The purpose of this document is to define the procedures for Arthur D. Little, Inc. (ADL), personnel to "close-out" the flight and ground-based lockers containing the Initial Blood Storage Experiment (IBSE) in preparation for the shuttle launch.

2.0 SCOPE

The scope of this document is to describe the hardware-related activities and the loading procedures as they interface with the hardware. The scientific aspects of the loading are described in the IBSE Scientific Protocol and in Annex 8 of the Payload Integration Plan. Performance checkout of the IBSE hardware is described in the <u>Opera-</u> ting. <u>Maintenance and Handling Procedures</u>.

3.0 MANPOWER REQUIREMENTS

The following IBSE personnel will be involved in the IBSE closeout:

D. W. Almgren, Arthur D. Little, Inc. - Engineering Manager

E. J. Boudreau, Arthur D. Little, Inc. - Technician

W. A. Curby, Lahey Clinic Medical Center - Science/Engineering Interface Manager

J. Young, Jr., Lahey Clinic Medical Center - Technician

D. Elwood, Lahey Clinic Medical Center - Technician

M. S. Jacobson, Children's Hospital - Scientist in Charge of Blood Bag Loading

F. Chao or D. Van Pelt, Children's Hospital - Observer

In addition to loading and closeout personnel, additional personnel will be available to:

- o observe procedures and schedule,
- authorize proposed changes to the procedures
 (e.g., use of spare unit or parts, change of loading sequence, etc.),
- o record any procedure changes.
- o interface with NASA,

o interface with Central Florida Blood Bank.

ORIGINAL PAGE IS OF POOR QUALITY

Arthur D. Little, Inc.

1

- o interface with working press,
- o interface with point of contact of closeout team (W. Curby).

The personnel for these activities include:

E. Michel, NASA JSC

- W. Paton, NASA KSC
- D. Surgenor, Center for Blood Research
- S. Kevy, Children's Hospital
- P. Glaser, Arthur D. Little, Inc.

Figure 1 illustrates the IBSE organization during loading and closeout.

4.0 SCHEDULE

The activities related to IBSE loading and closeout occur in relation to the shuttle launch. This schedule assumes a launch of January 6. 1986, at 7 A.M. The schedule of events 3 days prior to launch are shown on the timeline in Figure 2 and are described below. Figure 3 shows the loading sequence of units.

- 4.1 Three Days Prior to Launch (1/3/86)
 - 4.1.1 Post access list to closeout room.
 - 4.1.2 D. Almgren and E. Boudreau check inventory of IBSE Hardware and Equipment in Hangar L.
 - 4.1.3 D. Almgren and E. Boudreau check DC power supplies in Hangar L.
 - 4.1.4 D. Almgren and E. Boudreau inspect each of four IBSE module, (2) P/N 4215-1200 and (2) P/N 4215-1400, in its locker.
 - 4.1.4.1 D. Almgren and E. Boudreau install Kapton air flow passages in three platelet dewars.
 - 4.1.4.2 Perform airflow check (see paragraph 6.1 of <u>Operating, Maintenance and Handling Procedures</u>).

4.1.4.3 Set elapsed time meters and record in log book.

ORIGE STATES IS



FIGURE 1 IBSE ORGANIZATION

24 HOURS TO LAUNCH	18 HOURS TO LAUNCH	17 HOURS TO LAUNCH	15 HOURS TO LAUNCH	12 HOURS TO LAUNCH		IAUNCH 1/6/86 7 A.M.	Transfer to Orbiter 7 P.M. 1/5/86	Completion of loading 4 P.M. 1/5/86	L-Final Preparation for loading 2 P.M. 1/5/86	-Final Preparation for loading begins 1 P.M. 1/5/86	Preliminary Preparation for loading 7 A.M. 1/5/86	I DOPOUT SCHEDULE
			IECIKS	reters	6 1/5/86	ch to launch						FIGURE 2 0
		AIR FLOW	PERFORMANCE CH	ELAPSED TIME M	1/3/86 1/4/86	Day 3 Day 2 to launch to launc						

Arthur D. Little, Inc.



* Ground Unit Loading Sequence will be determined by need for flight substitution.

FIGURE 3

5

- 4.1.4.4 Closeout module without Kapton tape on perimeter of air passages.
- 4.1.4.5 Install module in locker.
- 4.1.4.6 Verify spare fuse in foam.
- 4.1.4.7 Perform thermal performance test (see paragraph 6.2 of Operating, Maintenance and Handling Procedures).
- 4.1.5 Repeat airflow and thermal tests with IBSE P/N 4215-1200 spare C/D module (S/N 0003)
- 4.1.6 D. Almgren and E. Boudreau store modules in Hanger L.
 - 4.1.6.1 Disconnect DC power supplies from IBSE modules by unplugging 25 foot cable from pigtail to DC power supply.
 - 4.1.6.2 Cover middeck lockers containing IBSE modules (located in Curby Crane) and associated 25-foot power cable with polyethylene bags.
 - 4.1.6.3 Cover spare IBSE module in cardboard box with polyethylene bag.
- 4.1.7 Power each locker for approximately 1 minute/day until launch preparations. Record data in logbooks.
- 4.2 Twenty-Four Hours Prior to Launch (1/5/86, 07:00 A.M.)
 - 4.2.1 E. Boudreau unplugs IBSE module A/B, P/N 4215-1400 (S/N 0002) from power source.
 - 4.2.2 E. Boudreau removes two electrical subpanel screws.
 - 4.2.3 E. Boudreau opens locker door to 180° position and supports it with Lucas door support.
 - 4.2.4 E. Boudreau removes top foam block.
 - 4.2.5 E. Boudreau removes belly band foam.
 - 4.2.6 E. Boudreau removes screws (approximately 16) from control box/dewar flange and verifies that the control box/bag support structure can be lifted from dewar.
 - 4.2.7 E. Boudreau repeats 4.2.1 to 4.2.6 for IBSE module P/N 4215-1200 (S/N 0004), dewars E/F.

- 4.2.8 E. Boudreau repeats 4.2.1 to 4.2.6 for IBSE module P/N 4215-1200 (S/N 0001), dewars J/K.
- 4.2.9 E. Boudreau repeats 4.3.1 to 4.3.7 for IBSE module P/N 4215-1200 (S/N 0005), dewars G/H.
- 4.3 Eighteen hours prior to launch (1/5/86 1:00 P.M.)
 - 4.3.1 D. Almgren and D. Elwood lift the blood storage rachs out of the locker and suspend them for P/N 4215-1-00 (S/N 0002), dewars A/B.
 - 4.3.2 D. Almgren removes the two thermal shields and stores them inside the dewars.
 - 4.3.3 D. Almgren covers the open dewars with polyethylene.
 - 4.3.4 Process repeated for P/N 4215-1200 (S/N 0004), dewars E/F, one thermal shield.
 - 4.3.5 Process repeated for P/N 4215-1200 (S/N 0001), dewars J/K, one thermal shield.
 - 4.3.6 Process repeated for P/N 4215-1400 (S/N 0005), dewars G/H, two thermal shields.
- 4.4 Seventeen hours prior to launch (1/5/86 2:00 P.M.)
 - 4.4.1 CBR personnel load blood bags into racks of P/N 4215-1400 (S/N 0002), dewars A/B, according to protocol.
 - 4.4.2 D. Almgren and D. Elwood install two thermal shields on P/N 4215-1400 (S/N 0002), dewars A/B, and apply Kapton tape to exterior of shields.
 - 4.4.3 D. Almgren and D. Elwood close out Kapton box around platelet grids.
 - 4.4.4 D. Almgren and D. Elwood lower bag support structure into dewars.
 - 4.4.5 E. Boudreau and J. Young closeout P/N 4215-1400 (S/N 0002), dewars A/B.
 - 4.4.5.1 E. Boudreau installs screws holding control box to dewars.
 - 4.4.5.2 J. Young raises the IBSE modules.
 - 4.4.5.3 E. Boudreau applies Kapton tape to perimeter of air passage.

- 4.4.5.4 J. Young lowers the IBSE modules.
- 4.4.5.5 J. Young removes eyebolt.
- 4.4.5.6 E. Boudreau puts in belly band foam.
- 4.4.5.7 E. Boudreau installs foam block while feeding subpanel through cutout in door.
- 4.4.5.8 E. Boudreau closes locker door.
- 4.4.5.9 E. Boudreau puts back subpanel screws.
- 4.4.5.10 E. Boudreau connects 25 foot cable to electric subpanel with laboratory dc power supply. J. Young notes time in log book.
- 4.4.5.11 E. Boudreau verifies current draw and temperature. J. Young records in log book.
- 4.4.5.12 Turnover of loaded fight module P/N 4215-1400 (S/N 0002) to NASA.
- 4.4.6 Repeat 4.4.3 to 4.4.5.12 for P/N 4215-1200 (S/N 0004), dewars E/F.
- 4.4.7 Closeout process is repeated for ground based units: P/N 4215-1200 (S/N 0001), dewars J/K, and P/N 4215-1400 (S/N 0005), dewars G/H.
- 4.5 Twelve hours prior to launch (1/5/86, 6:30 P.M.)
 - 4.5.1 NASA personnel remove P/N 4215-1400 (S/N 0002) and P/N 4215-1200 (S/N 0004) for cleaning of exterior, final weighing, transportation and installation in Orbiter.
 - 4.5.2 E. Boudreau and J. Young record temperature of units: o as they set out to Orbiter.

o as they arrive at Orbiter,

Activity record forms for each pre-launch time step (24 hours. 18 hours, 17 hours, and 12 hours) are attached.

24 HOUR ACTIVITY RECORD

			TIME	i
	FUNCTION	PERSON(S)	STARTED/COMPLETED	NOTES
1.	Unplug P/N 4215-1400 (S/N 0002) from power source (A/B)	Boudreau		
2.	Remove two electrical subpanel screws	Boudreau		
3.	Open locker to 180° pos.+ support w/ Lucas door support	Boudreau		
4.	Remove top foam block	Boudreau		
5.	Remove belly band foam	Boudreau		
6.	Remove (16) screws from control box and dewar flange	Boudreau		
7.	Verify control box/bag support structure can be lifted from dewar	Boudreau		
8.	Unplug P/N 4215-1200 (S/N 0004) from power source (E/F)	Boudreau		
<u>9</u> .	Remove two electrical subpanel screws	Boudreau		
10.	Open locker to 180° pos.+ support w/ Lucas door support	Boudreau		
11.	Remove top foam block	Boudreau		
12.	Remove belly band foam	Boudreau		
13	Remove (16) screws from control box and dewar flange	Boudreau		
 →	Verify control box/bag support structure can be lifted from dewar	Boudreau		

Arthur D. Little, Inc.

24 HOUR ACTIVITY RECORD

		TIME					
	FUNCTION	PERSON(S)	STARTED/COMPLETED	NOTES			
15.	Unplug P/N 4215-1200 (S/N 0001) from power source (J/K)	Boudreau					
16.	Remove two electrical subpanel screws	Boudreau					
 - ·	Open locker to 180° pos.+ support w/ Lucas door support	Boudreau					
18.	Remove top foam block	Boudreau					
19.	Remove belly band foam	Boudreau					
20.	Remove (16) hot plate screws from control box and dewar flange	Boudreau					
21.	Verify control box/bag support structure can be lifted from dewar	Boudreau					
22.	Unplug P/N 4215-1400 (S/N 0005) from power source (G/H)	Boudreau					
23.	Remove two electrical subpanel screws	Boudreau					
24.	Open locker to 180° pos.+ support w/ Lucas door support	Boudreau					
25.	Remove top foam block	Boudreau					
26.	Remove belly band foam	Boudreau					
27.	Remove (16) screws from control box and dewar flange	Boudreau					
28.	Verify control box/bag support structure can be lifted from dewar	Boudreau					
_							

RECORD COMPLETED BY:

18 HOUR ACTIVITY RECORD

		TIME					
	FUNCTION	PERSON(S)	STARTED/COMPLETED	NOTES			
1.	Lift blood storage racks and suspend P/N 4215-1400 (S/N 0002)	Almgren/Elwood					
2.	Remove thermal shields (2) and store in dewar	Almgren					
3.	Cover dewars with polyethylene	Elwood					
4.	Lift blood storage rack and suspend P/N 4215-1200 (S/N 0004) (E/F)	Almgren/Elwood					
5.	Remove thermal shield (1) and store in dewar	Almgren					
6.	Cover dewars with polyethylene	Elwood					
7.	Lift blood storage rack and suspend P/N 4215-1200 (S/N 0001) (J/K)	Almgren/Elwood					
8.	Remove thermal shield (1) and store in dewar	Almgren					
9.	Cover dewars with polyethylene	Elwood					
10	Lift blood storage rack and suspend P/N 4215-1400 (S/N 0005) (G/H)	Almgren/Elwood					
.....

			TIME	
	FUNCTION	PERSON(S)	STARTED/COMPLETED	NOTES
11.	Remove thermal shields (2) and store in dewars	Almgren		
12.	Cover dewars with polyethylene	Elwood		
	1		· · · · · · · · · · · · · · · · · · ·	

RECORD COMPLETED BY:_____

.

			TIME		
	FUNCTION	PERSON(S)	STARTED/COMPLETED	NOTES	
1.	Load blood bags into P/N 4215-1400 (S/N 0002) (A/B)	Jacobson/Van Pelt			
2.	Install two thermal shields	Almgren/Elwood			
3.	Apply Kapton tape to shields' exterior	Almgren/Elwood			
4.	Lower bag support structure into dewars	Almgren/Elwood			
5.	Install screws holding control box to dewars	Boudreau			
6.	Raise the IBSE modules	Young			
7.	Apply Kapton tape to air passage perimeter	Boudreau			
8.	Lower IBSE modules	Young			
<u>9</u> .	Remove eyebolts	Young			
10.	Put in bellyband foam	Boudreau			
11.	Install foam block while feeding sub- panel through cut- out in door	Boudreau			
12.	Close locker door	Boudreau			
13.	Put back subpanel screws	Boudreau			

,

			TIME	
	FUNCTION	PERSON(S)	STARTED/COMPLETED	NOTES
15.	Connect 25-foot cable to electric subpanel w/laboratory dc power supply	Boudreau		
16.	Note time of No. 15 in log book	Young		
17.	Verify current draw and temperature	Boudreau		
18.	Turnover loaded flight unit P/N 4215-1400 (S/N 0002) to NASA	Curby		
19.	Load blood bags into P/N 4215-1200 (S/N 0004) (E/F)	Jacobson/Van Pelt		
20.	Install thermal shield (one)	Almgren/Elwood		
21.	Apply Kapton tape to shield exterior	Almgren/Elwood		
22.	Close out Kapton box around platelet grid	Almgren/Elwood		
23.	Lower bag support structure into dewars	Almgren/Elwood		
24.	Install screws holding control box to dewars	Boudreau		
25.	Raise the IBSE Modules	Young		
26.	Apply Kapton tape to air passage perimeter	Boudreau		
27	Lower IBSE modules	Young		ĺ
28	. Remove eyebolts	Young		
29	Put in bellyband foam	Boudreau		

Arthur D. Little, Inc.

•

			TIME	
	FUNCTION	PERSON(S)	STARTED/COMPLETED	NOTES
30.	Install foam black while feeding sub- panel through cut- out in door	Boudreau		
31.	Close locker door	Boudreau		
32.	Put back subpanel screws	Boudreau		
33.	Connect 25-foot cable to electric subpanel w/labora- tory dc power supply	Boudreau		
34.	Note time of No. 33 in log book	Young		
35.	Verify current draw and temperature	Boudreau		
36.	Turnover loaded flight unit P/N 4215-1200 (S/N 0004) to NASA	Curby		
37.	Load bag into P/N 4215-1200 (SN 0001) (J/K)	Jacobson/Van Pel	5	
38.	Install thermal shield (one)	Almgren/Elwood		
39.	Apply Kapton tape to shield exterior	Almgren/Elwood		
40.	Close out Kapton box around platelet grid	Almgren/Elwood		
41.	Lower bag support structure into dewars	Almgren/Elwood		
42	Install screws holding control box to dewars	Boudreau		
43	. Raise the IBSE Modules	Young		

		TIME				
	FUNCTION	PERSON(S)	STARTED/COMPLETED	NOTES		
44.	Apply Kapton tape to air passage perimeter	Boudreau				
45.	Lower IBSE modules	Young				
46.	Remove eyebolts	Young		!		
; ∔7,	Put in bellyband foam	Boudreau		!		
48.	Install foam block while feeding sub- panel through cut- out in door	Boudreau				
49.	Close locker door	Boudreau				
50.	Put back subpanel screws	Boudreau				
51.	Connect 25-foot cable to electric subpanel w/labora- tory dc power supply	Boudreau				
52.	Note time of No. 51 in log book	Young				
53.	Verify current draw and amps	Boudreau				
54.	Load blood bags into P/N 4215-1400 (SN 0005) (G/H)	Jacobson/Van Pelt				
55.	Install two thermal shields	Almgren/Elwood				
56.	Apply Kapton tape to shields' exterior	Almgren/Elwood				
57	Lower bag support structure into dewars	Almgren/Elwood				
58	Install screws holding control box to dewars	Boudreau				

.

			TIME			
	FUNCTION	PERSON(S)	STARTED/COMPLETED	NOTES		
59.	Raise the IBSE modules	Young				
60.	Apply Kapton tape to air passage perimeter	Boudreau				
<u>51</u> .	Lower IBSE modules	Young	,	1		
52.	Remove eyebolts	Young				
63.	Put in bellyband foam	Boudreau				
64.	Install foam block while feeding sub- panel through cut- out in door	Boudreau				
65.	Close locker door	Boudreau				
66.	Put back subpanel screws	Boudreau				
67.	Connect 25-foot cable to electric subpanel w/laboratory dc power supply	Boudreau				
68.	Note time of No. 67 in log book	Young				
69.	Verify current draw and temperature	Boudreau				

RECORD COMPLETED BY:_____

Arthur D. Little, Inc.

			TIME			
	FUNCTION	PERSON(S)	STARTED/COMPLETED	NOTES		
1.	Convert modules P/N 4215-1400 (S/N 0002) and P/N 4215-1200 (S/N 0004) from laboratory dc power supply to battery power source	NASA				
2.	Units P/N 4215-1400 (S/N 0002) and P/N 4215-1200 (S/N 0004) are exterior cleaned, weighed, transported installed in Orbiter	NASA				
3.	Temperature of P/N 4215-1400 (S/N 0002) and P/N 4215-1200 (S/N 0004) logged as they set out to, and arrive at, the Orbiter	Boudreau				

RECORD COMPLETED BY:_____

•

Initial Blood Storage Experiment

OPERATING, MAINTENANCE AND HANDLING PROCEDURES

submitted as part of the

Acceptance Data Package Data Item 017 Contract No. NAS9-17222

submitted to

The Center for Blood Research Boston, Massachusetts 02115

submitted by

Arthur D. Little, Inc. Cambridge, Massachusetts 02140

October 30, 1985

ADL Reference C-53282

Arthur D. Little, Inc.

Table of Contents

٠

1.0	Purpose	1
2.0	Scope	1
3.0	Special Tools Required	1
4.0	Maintenance Schedule	1
5.0	Condition of Environment	1
6.0	Maintenance Procedure	1
7.0	Fault Analyses	8
8.0	Operational Instructions	10

٩.

.

1.0 Purpose

The purpose of this procedure is to provide guidance and instruction to personnel who will be maintaining and operating the IBSE modules.

2.0 Scope

This procedure covers maintenance and operation of the five, flight IBSE modules.

3.0 Special Tools Required

No special tools are required.

4.0 Maintenance Schedule

Maintenance of the IBSE modules is not required on a regularly scheduled basis. It is recommended that a dewar air flow test and a thermal performance test of a module be conducted at least 60 days prior to any reflight of the module.

5.0 Condition of Environment

The thermoelectric coolers within the IBSE module reject their heat to the cooling air flowing through the module with more power dissipation required at higher air temperatures. Satisfactory operation of the modules has been measured with inlet air temperature up to 80° F. It is recommended, however, that all testing of the IBSE module be performed at ambient air temperatures between 65 and 75°F to ensure a reduced power dissipation of the module during testing.

6.0 Maintenance Procedure

Maintenance procedures for the IBSE hardware is based on two basic concerns: the dewar air flow rate and the thermal performance of the thermoelectric devices. In the following description, reference should be made to Figures 6-1 through 6-4 which are the exploded views of the IBSE hardware.

6.1 Dewar Air Flow Maintenance Procedure

The 40-50 ml/min specified air flow rate through the dewars is provided by a small, plastic, gear pump powered by a 12 v dc, brush-type motor. This pump is located in the center compartment of the control box (see assembly drawing 4215, sheet 5). The air flow rate being provided to each dewar is



FIGURE 6-1 OVERALL VIEW OF THE COLD/WARM MODULE

Arthur D. Little, Inc.



FIGURE 6-2 DETAILED EXPLODED VIEW OF THE COLD/WARM MODULE

Arthur D. Little, Inc.

3



FIGURE 6-3 OVERALL VIEW OF THE COLD/COLD MODULE

A Arthur D. Little, Inc.



FIGURE 6-4 A DETAILED EXPLODED VIEW OF THE COLD/COLD MODULE

Arthur D. Little, Inc.

ORIGINAL PAGE IS OF POOR QUALITY

best measured at the exit of the air tube which terminates inside each dewar at its closed end. The following procedure is required to measure the air flow rate:

- 6.1.1 Disassemble control box from flanges of two dewars by removing eight, allen head screws around the perimeter of each of the two dewars.
- 6.1.2 Lift control box assembly and bag support structures from dewar assembly and lay on its side on a piece of soft material, e.g., a silicone rubber mat.
- 6.1.3 Remove the open end of the short horizontal section of TFE air tubing at the bottom of each dewar from inside the thermal shield or bag support structure.
- 6.1.4 Attach a calibrated (air) orifice or similar device to the exit of each of the two exposed ends of the air tubing. An air filter identical to the one in the cold plate can be used.
- 6.1.5 Install an inclined water manometer, or a similar pressure differential measuring device, to read the pressure drop across the orifice or similar device. The air flow pressure drop characteristics of the orifice should be approximately 0.1 inch of water at 100 ml/min air flow.
- 6.1.6 Apply 28 v dc power to the module via the connector on the electrical subpanel and measure the flow rate of air exiting each air tube.
- 6.1.7 If adjustments are required to the air flow rate, turn off the power to the module, and remove the top cover of the control box.
- 6.1.8 Adjustment to the air flow to both dewars, simultaneously, is achieved by changing the setting of the teflon valve located in the center compartment of the control box. This valve setting controls the fraction of the air output from the pump which flows to the two dewars. The remaining fraction of air flows to the electronic compartments of the control box as cooling air.
- 6.1.9 If the fraction of the air flow rate to one dewar has to be changed with respect to the air flow rate to the second dewar of the same assembly, install a small bore tube inside the appropriate TFE air tube at the bottom of the dewar. This additional flow restriction serves to divide the flow between the two dewars coming from a common source upstream. Table 6.1.9-1 summarizes the air flow rates set in all of the IBSE modules at the time of delivery of hardware to NASA.

TABLE 6.1.9-1 MEASURED AIR FLOW RATES TO IBSE DEWARS O OCTOBER 27, 1985

•

	5/N 0001	(J/K)	COLD:	63 ML/M	WARM:	250 ML/M
,		(A/B)	COLD:	100 ML/M	COLD:	94 ML/M
)	S/N 0002		COLD:	50 ML/M	WARM:	240 ML/M
>	S/N 0003	(C/D)			WARM:	263 ML/M
)	S/N 0004	(E/F)	COLD:	75 ML/M		110 NT /N
2	S/N 0005	(G/H)	COLD:	106 ML/M	COLD:	119 ML/M

6.2 Thermal Performance Test

The thermal performance of the dewars can be measured from the exterior of the module without disassembly. There is a YSI model 44006 thermistor mounted on the cold plate of each dewar whose resistance can be read from test jacks on the front electrical subpanel of each module. The thermal performance test procedure using these test jacks is as follows:

- 6.2.1 Mount the IBSE module inside a block of foam similar to the Pyrell foam used in a middeck locker.
- 6.2.2 Electrically attach a resistance measuring meter to each of the two, white, test jacks on the front panel. The left, white, test jack is connected to the left dewar, etc. The third, blue, test jack, located below the two, white, test jacks as viewed from the front of the panel with the primary power connector to the right, is common to both thermistors.
- 6.2.3 Apply 28 v dc power to the module.
- 6.2.4 Record the transient temperature of each cold plate during cooldown.
- 6.2.5 Figure 6.2.5-1 is a plot of the transient cooldown temperature of each of the cold plates in the seven cold dewars (5±1°C steady state) recorded just prior to the time of delivery of the flight hardware to NASA with an ambient air temperature of approximately 70°F.
- 6.2.6 If the transient temperatures of a dewar does not approximately correspond with its previously measured values (see Figure 6.2.5-1), then initiate a detailed inspection and testing of all hardware and electronic subassemblies using engineering and electrical drawings.

7.0 Fault Analysis

The primary parameters to measure for an operating IBSE module are the temperatures of the cold plates inside the dewar (see Figure 6.2.5-1) and the current draw of the locker. With an ambient air temperature of approximately 70°F, the initial current draw of a cold/cold module (at 28 v dc) is approximately 3.5 amperes. The steady state current draw of the same locker for the same air temperature is approximately 1.5 amperes. The initial current draw of a cold/warm module is approximately 3.5 amperes and its steady state value is approximately 1 ampere. The current being drawn by a cold/warm module may show a 5 - 10 second cyclical variation of approximately .1 ampere due to the controller for the warm dewar periodically shutting off. If

60 Q T_{air} = 70°F 80 ELAPSED TIME (MINUTES) ۲ ۲ 70 DEWAR H DEWAR G DEWAR E ပ В 4 DEWAR DEWAR DEWAR 60 1 ı I I ī ł 00 \bigcirc + ⊲ + 50 40 8 0 20 0 \$ 2 С 1s1C-3 0 10 ഹ 20 25

OCTOBER 27, 1985 TRANSIENT COOLDOWN OF FIVE COLD IBSE DEWARS -FIGURE 6.2.5-1

measured current values are significantly different from the values plotted in Figure 6.2.5-1, a detailed inspection of the electrical subsystem should be initiated.

There are three thermostats located inside the control box which will open if local, internal temperatures exceed 118°F. The thermostats are in series and will shut off all power to the module except to the two air cooling fans. There is one additional thermostat on each cooling fan to turn off power to that fan if its temperature exceeds 118°F.

There are three fuses inside the control box which, if blown, will shut off power to one of the two electronic temperature controllers (3a) or to the entire module (7a).

8.0 Operating Instruction

With the control box bolted to the dewar assembly, and a block of foam surrounding the assembly, the only remaining action to be taken by a user to initiate operation of the module is to apply 28 v dc power to the electrical connector on the front of the electrical subpanel. The application of 28 v dc will cause the green LED on the front of the locker to light. A blown fuse inside the locker does not turn the LED off. In a 0-g environment, the air flow flag over the exit air slot or grille will flutter to show that the cooling air fans are operating and, therefore, the module is receiving power.

After operating a cold dewar for an extended period of time, the cold plate/bag support structure should be removed from the dewar assembly and any condensed water removed. A jet of pressurized air can be used to blow out any moisture which has accumulated around the thermoelectric devices which are located in the gap between the hot and cold plates of the dewar lid assembly.

PRECEDING PAGE BLANK NOT FILMED

.

THE CENTER FOR BLOOD RESEARCH

٩.

800 Huntington Avenue, Boston, Massachusetts 02115 (617) 731-6470

November 19, 1985

Mr. William Paton/LSSM CS-SED-4 NASA JFK Space Center Florida 33899

RE: NAS 9-17222/IBSE

Dear Bill:

Enclosed please find a copy of the DD 1149 form with equipment and supply list attached.

For your information we also enclose an itemized bill-of-lading made out by the transportation company, Daley & Wanzer (form #1190).

The equipment will be leaving Boston on Friday, November 22nd, and will be delivered between 8 am and noon on Friday, December 6th.

If there are any changes that are required we are able to reach the driver in transit on a daily basis.

Please be advised that there are four IEC centrifuges that weigh approximately 1100 lbs. each and are palletized for fork-lift unloading at KSC.

Very truly yours,

Richard L. Korn Executive Vice President

RLK:rar Encl.

cc: E. Michel/EX4 VD. Surgenor, CBR F. Lionetti, CBR C. Kelliher, ONR

RECEIVED

NOV 2 0 1985

PRESIDENT'S OFFICE

Ten ts	SHEET TOTAL GRAND TOTAL 20. RECEIVER'S VOUCH
nen ts.	\$ <mark>\$ 8</mark>
tred to	
returne d Resea	
s to be etion of	Antipation of the second of th
r completion	
Equi The afte	
t amon/IEC of Mass. hildren's Hospital enter for Blood Research rthur D. Little signee rn elliher . Ed Michel/EX4)	I: David Almgren)
Equipment Equipment C = Cr C =	ADL (ATTN
ctor Furnished attached sheet ition Distribu No. Origi 2. F 4. h	ON VIA MATS 5.
Contra see Requis Copy	A MANSPORTAT
	<pre>(a)</pre>

Attachment A to DD 1149 19 Nov. 1985 Reg # NAS 9-1722 RQ003 DAMON/IEC DIVISION 15 Boxes-includes: Notor assembly (216) Nylon cushion Conical (315) Shield (320) (326) trunnion ring 354 trunnion nn (356) Shield 21 889) rotors Z tabes (890) Poter Ttubes 9 cap (897) rotor 7 tubes (92) Ĵ head (977) rotor Emoly cotest inst cy to logis Buckets (1024 (1580 metal sealing caps 6000 200/2001 WHZ (2257) eur tubes (2278) DPR auto clear tubes 2360) Centra 7R (115/60 cycle) Centra 7 (115/60 cercle) Case autoclear tube CS autoclear tube 2362)* Centra 7 286) 2997) ≁ (3600) M-25 (3615) 50/60Hz * 1201 Micro-MB Central adapter 7 plac - 6/2 (5707) Central adapter 19 plac red (5719) Central adapter 37 plac yellow (5737) holder afto bucket (5749) adapter (7574 6 ORIGIMAL PAGE IS OF POOR OUALITY

NASA IBSE IMMUNOLOGY: UMASS MED CTR.



LIST OF ITEMS TO BE SHIPPED FROM C.B.R.



achment B to DD 1149

Keg# NAS 9-17222 -RQ-003



Approximate value: \$160.00

BOX	5:	Description	<u>Oty.</u>
POX		20 L plastic carboy Plastic grad. cylinder, 2L' Plastic grad. cylinder, 1h Plastic grad. cylinder, 500 m Plastic grad. cylinder, 100 m Notebook & paper Plastic vol. flask 1L 5% Bovine albumin P.B.S.	1 1 1 1 1 2 500 ml
		3 containers, each with the r	

180g NaCl 22.44g Na₂HPO₄ 6.712g NaH₂PO₄ 20 NaN₃

 $\mathcal{A}_{\mathcal{S}}$ containers, each with 1.5g Na₂EDTA

Attachment C to DD 1149 AND. 1985 Reg# NAS 9-17220RQ-003 Q.Y Dr. Jacobson and Dr. Keny-Children's Hosp. 10 boxes Dr. Ausprunk - Children's Hosp - 3 boxes Dr. Scanlon

Hachment D to DD1149 ANOV. 1985 Reg # NAS 7-17222-RQ-003 p.5 R. Fabian Lionetti (White Cells) Center for Blood Research Dilu. Vials (144) 1 Case Box Sampling Site Complexs 10x 75 mm glass tubes 10 ml disp. pipets (96) 1 Cuse 4 boxes 5 18 G G synnes DONE pipet tips (blue) 300 nel T-packs diff stains 1000 24 H Stains 17×100 ml tabes T Caps bx (125)50Ks NCC. Syrings tubes (crange top) tubes (orange top) cc syninges m 50 nl x 2 : pks 100 0 2-200 ml) hx erlennugr flasks Sidearm Stoppers pipets 53/4" Stenlization filter 500 mc scintillation nalcene Vials うの

Hachment D to DD1149 19 Nov. 1985 Reg # NAS 9-1722-RQ-003 DX 3: Glass beakers (5 ranging from 10,mol -> 500 al) Himen 500 ml glass graduated glinder (PG) dextran (100g) Clean tiverse II Scintillation Cocktail Dimburgered autoclause by on 50 cc. Syrings Kinwipes racks pipetman tex mixer 45 Calculator 5) 12 × 75 mm tubes tube c caps (1000 201 ul pipetman ul pipetman in tabing o Caps (> 130) tr PBS Hite, red, + blue label tape) bx microscope slides) bx cover slips 3 spatulas marker pens

19 Nov 1985 KM # NAS 7-1722= KQ 003 tachment D to DD1149 p.7 BOX # JKKR1 Platelets 41 Countoff Sochium Phosphake abasic 5009 1009 EGTA 4.39 Sodiom uzide SDS 2509 .02m/Ethenol 100ml PMSF 5009 Tris 2 l ACD Sodium meta perpodate 1009 loome triton X-100 Los (00mlNonicet-P-40

tachnut D to D1149 19 NW. 1915 Reg # NAS -9-1722 KB 003 Box # JKKRZ Platelets P8 Folins reagents D CUSOY 18 Boome (3 2% Nazloz LOOOML (3 2% Na/h testak 500ml NaOH Boomd IonlyInal Schety Pacting ٠.

Reg # NAS - 9-1722RQ003 19 Nov. 1985 tach ment D to 1149 DD Box # JKK123 Platelets (P) ncl 5009 5009 Nall 5004 Dextrose SUCTOSE 5009 1 Bottle of Follows Keagent 500ml Call 8 met anticoous 500g # gorlands# ł, ORIGIMAL PAGE IS OF POOR QUALITY

Reg # NAS 9 A202 ROOM Hackment D to DD 1149 M Nor. MK (p.10) Box # JKKR4 Plaklets (p.10) recorder for Dual Channel agregometer • • . . . **.**. . . . -ORIGINAL PAGE IS OF POOR QUALITY

achnest D to 1149 DD 19 Nov 1985 Keg # NMS 9 17322 KQ 863 Box # JKKR5 Platelets (PI) Dual channel aggregometer · -....**.** • • · • . · · · · · . . • • ••• . .

19 Nov 1985 Reg# NAS 9 17222 RQ 03 Hachment D to 1149DD Box # JIKKR 6 Platelets (P12) 2 Bottles Hydrogen Perovice Loonlea in Original Safety Pucking 500g Paraformattehyde Bog magnesion chlorice . . . **. .** Sadium phosphake monobaser 1000g 2509 EDTA boomd glycerin PUE gloves 1 Box large -(Box small latex gloves 1 Box medium small 11 11 large • Solety packing loome. 1Box BME in -----Paralim 1 stir bar trit 10ml PAPetter 8 bays of 25pt • 5ml 11 10 bugs of 25 • Iml " 6 bays & 11 <u>____</u> Ind 11 1 Box of 11 ••• محمد بمنتخذ من من م , · · · and the second -----. . . . **.** . . . **.** . ORIGINAL PAGE IS

OF POOR QUALITY

achnest D to 1149 DD 19 Nov. 1985 Reg # NAS 9 17532 R.Q 003 (p.13) Box # JKKR7 Platelets Hcl Special Safety Packing 500ml Cougolation assay waterbath 1 Box Silrconized Pastuer Pipettes Z Boys Disposable 11 /1 500 Blue pippette tips 1000. Yellow 11 11 3 rachs for 50 cc test tubes 3 100 ml. Graduate Cylinders 1 22 11 11 1 12 les 11 1 500ml 11 11 2 ruchs for 15 ml test tubes. 1 thermometer for Water bath • 6 500ml plastre beaters 2 loooml 11 1 250ml 11 17 2 100 ml 11

• •

•

ORIGINAL PAGE IS OF POOR QUALITY
Hachment D to 1144 DD 19 Nov. 1985 Rey # NAS 9 1732 2 RQ 203

Box # JKKR8 Plaklets (2.14) 100 15 mL Orange TOP test tubes 375 13×100 Blue 11 11 11 - ... · · | | 500 17×75 Blue 11 (1) _____ 100 scintillation Vials **_**___. 1 buy of santillation Vial Caps 12 Reagent Bottles • • • . • •• •• •• ··· - - - -2 ·• · · •• . · · · . ••

.

tachment D to 1149 DD 19 1401955 Reg # NHS9 A222 10 003 Box # JKKR9 Plaklet 200 scintillation viels + tops (7ml) (2.15) 2 usige ractes____ Z_5ml_scrstelf_tybe_(conrul)_ractis____ 2 rolls of green fepe 11 11 White 11_____ 11 11 red 11 11 11 Orange 11 3_test_tube_racks____ 1_Ice_bocket____ to Box of Freher weigh Paper 20 lg_weigh_boats <u>6 Spatulas</u> . 1000 CONTRAL Sarstelt tubes_____ 250 red 11 CLPS 250 Blue 11_____ 250 white 11 Z_Boxes_kim_wifel____

tach ment D to 1149 DD 19 Nov. 1985 Reg # NAS 9 17222 RQ 003 (P.K) anlelt Debbie boxes 3 No etc. Curby Hoists + bikes 2 .

Hackment E to 1149 DD 19Nov. 1985 Reg# NAS 9 1722200 - (- p: 1+ Arthur D. Little includes 2 NASA emptige containers 5 . .

		DALEY 9 1	HOUSEH	OLD 600	DS- DESCRIPTI	VE INVENTORY	PAGE NO	NO, OF PAGES	
W	60105 m	DALET & Y	VANLER, IN	.	ABENT	DHW	CARRIER	3 HEFENENCE NO.	7
JWHER'S	GRADE OR	RATING AND NAME	Cent	R Fr	Blood K	PESEARCH	CONTRAC	T OR GOL. HO.	1
RIGIN L	DADING ADD	ACSS	Sont	Intin	Land AIR	RIT	60VT. SE	RVICE ORDER NO.	-
ESTINAT	ION		- Vai	119		Postal	VAN NUM	DER	-
	Di	SCRIPTIVE SYMB	OLS CHAR	yent -	EXCE	ENTER PTION SYMBOLS		1.004.7104 574801	-
0/W - BLA C - COLOR CP - CARP PBO - PAC CD - CARP	ACK & WHIT I TV HER PACKEL KED BY OW HER DISASS	E TV DBO-DISA PB-PROFE PE-PROFE NER PP-PROFE EMBLED NOTE: THE C	SSEMBLED BY OWNER SSIONAL BOOKS SSIONAL EQUIPMENT SSIONAL PAPERS MISSION OF THES	BE - BENT BR - BROK BU - BURN CH - CHIPF CU - CONTI DITIO E SYMBOLS IN	EN D-DE ED F-FA 50 G-GO ENTS & CON-L-LO N UNKNOWN NDICATES GOOD C	NTED M-MARRED DED MI-MILDEW JGEO MO-MOTHEATEN 35E R-RUBBED RU-RUSTED ON DITION EXCEPT #	SC-SCRATCHED SH-SHORT SO-SOILED T-TORH W-BADLY WORH Z-CRACKED OR NORMAL WEAF	1. AUM 7. REAR 2. BOTTOM 8. RIGHT 3. CORMER 9. SIDE 4. FRONT 10. TOP 5. LEFT 11. VENEEJ 6. LEGS 12. EDGE 2.	
17EM NO.	CR. REF.	AR1	TICLES		COND	TION AT ORIGIN		EXCEPTIONS (IF ANY) AT DESTINATION	_
2		10000	man Red C						-
3		LATCHE U	woven (ente		Caste	7	P. P. BRC	
4	-	SMALL	C	am 1	1 These	JEVEN (A	ens	20 800	
5		Large n	DEN	u 1	to the	-finite	young	T. H. BRC	
6			И	• (Those	items and	ang.	<i>U</i>	
7		MEDIUM	FLAT .	*	Linth	the Ont	top "L		•
8		LARGE	K 1	4	J				
<u> </u>		X-tra La	vige 1	(
1		LARGE	н	+				4	
;		MEDIUM							
- 3		JMTU.	#						
4		<u> </u>	<i>M</i>	<u>;</u>				<u> </u>	
5		MEDINA						4	
6		3-4	<u>//</u> /	1		7	=+		
7		MEDIUM	WHITE #	+		+us car	on -		
8		и.						<u>7</u>	
9		4	(7	
)		Small		d					
		<u> </u>		·					
		LANGE		/				4	
-+		Long Fr	AT 4					<u> </u>	
		<u> </u>	4 4	· <u> </u>				<u> </u>	
		2.1						4	
		Lene				The cash			
			- 7					······	
			4						
		3.1	4		Z	tu conto	$\overline{\mathcal{Q}}$		
		<u> </u>						<u>4</u>	
		LARGE							
			K						
			<u>¥</u>						
		Shapt 4			······			4	
		ARGE	11			······································		001000	
		A TAL	n						. P/
		4 The	<u>ک</u>					- Ur POOR	QU
-L		<u>H h</u>	~ ~					4	-
REM	ARKS/E	CEPTIONS			· · · · · · · · · · · · · · · · · · ·		<u>1</u>	<u> </u>	
A VE CHEC	CRED ALL THE	TEMS LISTED AND NUM	BERED 1 TO	INCLUSIVE A	ACKNOWLEDGE THAT	THIS TAPE LOT NO	TAPF		
RNI	NG-	BEFORE S	SIGNING CHECK SHIP	MENT. COUNT	ITEMS AND DESCRIE	E NOR PRO	COLO	A	
	CONTRACTO	, CARRY ON AUTHO	AIZED APENT (PRIVER)	DATE	OVE.	TRACTOR, CARRIER OR A	THRU	IVER) DATE	
T k	SIGNATURE	Hour	Thaburn	1/19k	DESTI-	SHATURE]			
61N		~~ V1		POATE	NATION 0	ALK OR AUTHORIZED AGE	NT	DATE	

-

ACTOR OR CARRI	IER HOU	ISEHOLD GOOI	DS DESCRIPTIVE II	VENTORY	2	REFERENCE NO.	_
1	DALEY & WANZER	INC.	AGENT Z	HJ		396	
	TING AND NAME	+ -	pho Par	and		RANCE	
R 3 68402 00 100	<u> </u>	we for	CITY	STATE	GOVT. SER	VICE ORDER NO.	
N LOADING ADOR	[35		Art		VAN NUMB	ER	
INATION	Kes	medy)	oace Onles	N EVHENIS		LOCATION SYMBO	
DE - BLACK & WHITE OLOR TV CARRIER PACKED - PACKED BY GWN CARRIER DISASSI	SCRIPTIVE SYMBOLS E TV DBD-DISASSEMBLED BY PB-PROFESSIONAL BOO PE-PROFESSIONAL BOO LER PP-PROFESSIONAL PAP EMBLED	OWNER SE-BEN BE-BEN BE-BEN BU-BUR CH-CHI CH-CHI CH-CHI OIT THESE SYMBOLS	T D-DENTED NED F-FADED PFED 6-600GED ITENTS & CON- L-LOOSE ION UNKNOWN INDICATES GOOD COND	M - MARRED MI - MILDEW MO - MOTHEATEN R - RUBBED RU - RUSTED ITION EXCEPT FO	SC - SCRATCHED SH - SHORT 50 - SOILED T - TORN W - BADLY WORN Z - CRACKED R NORMAL WEAF	1. AEM 7. NLA 2. BOTTOM 8. RIG 3. CORNER 9. SIDI 4. FRONT 10. TOP 5. LEFT 11. VEN 6. LEGS 12. EDG	HT E EER SE
M CR.	NOTE: THE DESIGN C		CON 01710	AT ORIGIN		AT DESTINATION	
o. REF.	No. 1	Cartel				P.B. BA	C
	X-the large	<u>urado</u>					<u> </u>
	<u> </u>	4	VERY FRA	gie		4	
	// //	4	<u>N</u>	n	. <u></u>	4	
	N K	×	k	<u>k</u>			
5	here this					u	
;	H ⁰	<i>u</i>		I Park		4	
9	3.1		D/	<u></u>	q		_
9	4			<u> </u>		4	
r						4	
<u>'</u>	Lange					4	
2						4	
3	<u> </u>						
		4					
<u></u>	" WHITE	4				4	
7	- H - MARCE	и					
8	и	4				+	
9							
10							
1							
2							
3							
4							
3							
8							
9							
70							
1						ORIGINA	
2		<u>-</u>				OF POOR	- FAGE
3						1	QUAL
9							
80							
HEM REMA	RKS/EXCEPTIONS						
		ED 1 TO	WCLUSIVE AND ACKNOWLEDGE	THAT THIS TAPE LOT	NO	TAPE COLOR	
WE HAVE CHECK	OMPLETE LIST OF THE GOODS TENDER	DAND OF THE STATE OF	ITHE GOODS RECEIVED	SCRIBE NOS. FROM	M	THRU	
ARNI	NG LOSS OR DAI	AGE IN SPACE ON TH	AE RIGHT ABOVE	CONTRACTOR, CARR	ER OR AUTHORIZED	AGENT (DRIVER)	DALE
	ONTRACTOR, CARRY OR AUTHORIT	tortown	11/18 DESTI-	(SIGNATURE)	IZED AGENT		DATE
	INNAIUNEI ALLEGO					1	

TABLE I

	PLAN K*	PLAN KE*	PLAN E*
Primary Landing Site	KSC		Edward
Secondary Landing Site		Edward	
Location of Lab & Investigators	KSC	KSC	Boston
Flight Hardwares at the time of launch	KSC	KSC	KSC
Ground samples (Zero time and Control)	KSC	KSC	Boston
Loading of the Dewars	KSC	KSC	KSC
Transportation of ground control after launch	Within KSC	Within KSC	From KSC to Boston in dewars via air-condi- tioned van
Transportation of flight sample after launch	Within KSC	From Edward to KSC via a jet arranged by IBSE team	From Edward to Boston via a jet arranged by IBSE team

SUMMARY OF MAIN IBSE ACTIVITY

Plan K: KSC landing

Plan KE: Edward landing site with KSC as primary landing site Plan E: Edward as the primary landing site

:

Edward as the primary landing site Plan E:

TABLE II

PRE-LAUNCH ACTIVITY

PERSONNEL					KW K			
LOCATION	Boston	ADL Boston	Boston to KSC	KSC	CBR Boston	CBR Boston	Boston to KSC	
TARGET DATE & PERIOD			lits	0 l	t and d to KSC	and	and	
DESCRIPTION	Acceptance of Flight Nardware	Storage of Flight Hardware	Transportation of 5 un of hardware to KSC	Receiving and Testing hardware	Inventory of equipmen labwares to be shippe	Packing of equipment labwares	Shipping of equipment and labwares to KSC	
PLAN* INVOLVEMENT	К, КЕ, Е	K, KE, E	K, KE, E	K, KE, E	K, KE	K, KE	K, KE	
ACTIVITY CODE	I-1	II-2	II-3	H-4	L-1	L-2	L-3	

KSC landing Edward landing with KSC as the primary landing site Edward landing with Edward as the primary landing site * Plan K: Plan KE: Plan E:

• E PERSONNEL -LOCATION Boston KSC Boston KSC Boston KSC KSC KSC KSC KSC TARGET & PERIOD DATE Arrangement for transporta-tion of personnel (Air and Ground) Receiving of equipment and labwares Check-out equipment and Jabwares List of Personnel and Setting up IBSE Lab Room and Board Communication DESCRIPTION clearance INVOLVEMENT PLAN* K, KE ACTIVITY CODE P-3 P-4 P-2 I^{-2} L-6 P-1 L-4

•

PERSONNEL	CFIJB Personnel DV & NJ	FL and Crew DV & MJ CFBB Personnel	DV, NJ, BP	DV, MJ, BP	ſ₩	ſW	BP
LOCATION	CFRB Or lando	CFBB	CFBB	CFBB	Orlando to KSC	Orlando to KSC	Orlando to Boston
TARGET DATE & PERIOD	Jo			ហ្គ	its of	Time	Time
DESCRIPTION	Collection of 36 units c whole blood	Component Preparations	Pooling Blood Samples	Allocating Mlood Sample	Transportation of 56 un Alood samples to KSC	Transportation to Zero Sample	Transportation of Zero Sample
PLAN* INVOLVEMENT	K, KE, E	K, KE, E	K, KE, E	K, KE, E	K, KE, E	K, KE	ш
ACTIVITY CODE	B-1	B-2	l1−.3	l) - 4	B−5	B-6	B-7

1

...

PERSONNEL	EB, DA, WC JY, RI	NJ, FC	DVP, BP	EB, DA, WC JY, RI, NJ, DVP, BF			NASA personnel
IOCATION	KSC	KSC:	KSC	KSC	KSC	KSC	KSC
TARGET DATE & PERIOD	Preloading preparation of the Dewars	Preloading storage of Blood Samples	Coding and distributing of zero time sumples	Loading Blood bags into the dewars	Close off the dewars	Transfer flight lockers with test samples to NASA personnel	Install IRSE dewars/lockers in the space shuttle
INVELVENIENT 11.1.AN★	K, KE, E	К, КЕ, Е	K, KE	K, KE, E	К, КЕ, Е	K, KE, E	K, KE, E
ACTIVITY CODE	[]-1	11-2	1)-3	D-4	1-5	1)-6	1-1

LC.

.

personnel PERSONNEL DVP, BP NASA to Boston LOCATION Buston Boston CBR KSC KSC KSC KSC NSC & PERIOD TARGET **UATE** Transport loaded ground lockers Transport loaded ground lockers to the storage area Connect and check power supply · · Connecting and checking power Coding and distributing zero to the ground/dewnrs lockers with power supply connected Storage of ground lockers Storage of ground lockers supply to the dewars time sumples DESCRIPTION in Boston to Boston I NVOLVEMENT I NVOLVEMENT K, KE, E КE K, KE KE к. К Κ, ω μ تنا ACTIVITY D-10 11-11 E-2 E-3 E-1 CODE 0-0 U∘.8

. _

PERSONNEL								
LOCATION	Boston	Roston	Boston	Buston	Boston	Or Lando KSC Bos ton	Or lando KSC	Or Jando KSC Boston
TARGET DATE & FERIOD								
DESCRIPTION	List of location of laboratories	Request for special and additional equipment	Distribution of equipment	Check out equipment	List of Personnel and clearance	Arrangement for Transportation (Air and Ground)	Rooms and board	Communication
PLAN*	ы	ы	сı	сı	Ŀ	ت د	Ŀ	Ŀ
ACTIVITY CODE	₽-3	5-5	E-6	E-7	E-8	රා - ප	E-10	E-11

OF POOR QUALITY

. .

TABLE IV F LANDING ACTIVITY
rost L

.

•

PERSONNEL	NASA Personnel	NASA Personnel	NASA Personnel	NASA IBSE Personnel		MC
LOCATION	KSC (Plan K) Edward (Plan KE or E)	KSC (Plan K) Edward (Plan KE or E)	KSC (K) Edward (KE or E)	KSC	KSC (Plan K & KE) Boston (Plan E)	KSC (Plan K & KE) Boston (Plan E)
TARGET DATE & PERJOD						
DESCRIPTION	Dewar/locker operational check prior to removal from the Shuttle	Remove the dewars from the Shuttle	Connecting to Battery power supply	Transport Flight lockers to IRSE lab and turn over to IBSE team	Transport ground lockers to JBSE Lab	Locker operational check prior to opening
FLAN INVOLVENENT	K, KE, E	K, KE, E	K, KE, E	Ж	K, KE, E	K, KE, E
ACTIVITY CODE	PL-1	5-11	PL3	PL,-4	ΓL-5	9- 1d

ORIGINAL PAGE IS OF POOR QUALITY

<u>.</u>

,

•

PERSONNEL	<pre></pre>	K & An E)	K& an E)	on IBSE KE Investiga an E) tors	K & an E)	K& an E)	K & an E)
LOCATION	KSC (Plan F KE) Boston (Ple	KSC (Plan KE) Boston (Ple	KSC (Plan KE) Boston (Pla	KSC & Bost (Plan K & Boston (Pl	KSC (Plan KE) Boston (P]	KSC (Plan Boston (Pl	KSC (Plan KE) Boston (Pl
TARGET DATE & PERIOD						шеа	Lion
DESCRIPTION	Opening the dewars lifting the content photography	Remove Blood bags and fill out check list	Coding and distribution of Rlood samples	Testing and Assay	lloist Team on position	Coding and distribution to in position	Investigator team in posi-
F 1.AN I NVOLVENENT	K, KE, E	К, КЕ, Е	K, KE, E	K, KE, E	K, KE, E	K, KE, E	K, KE, E
ACTI VITY CODE	Ն-1վ	F1,-8	6-14	PL-10	FL11	PI,-12	PL-13

ORIGINAL PAGE IS OF POOR QUALITY

:

.....

-

٠

PERSONNEL	NASA Personnel			- - - - - -		
I.OCATION	Edward Air Rase	Edward	Edward	KSC (Plan KE) Boston (Plan E)	KSC (Plan KE) Boston (Plan E)	KSC (Plan KE) Boston (Plan E)
TARGET DATE & PERIOD						i es
DESCRIPTION	Transport flight lockers to MIV Field Airport & turn over to IBSE team	Locker Operational Check	Secure lockers and batteries into the transporting jet	Fly sample locker to INSE lab location	Remove lockers and batteries from the jet	Transport lockers and batter to IBSE lab
LINIISHISHI. LIVN	KE, E	KE, E	KE, E	KE, E	KE, E	KE, E
ACT I VI TY CODE	PL-14	PL-15	11-16	FL-17	PI 18	PL-19

ORIGINAL PAGE IS OF POOR QUALITY

Ξ

٠

Acorn Park Cambridge, Massachusetts 02140-2390 617 864-5770 Telex 921436

Arthur D. Little. Inc.

October 15, 1985 McDonnell Douglas Technical Services, Inc. 16441 Space Center Boulevard Houston, TX 77058 Mr. Stephen Gotch Attention: Mail Code D3 Power Dissipation of the Initial Blood Storage Subject: Experiment ADL Reference 53282-51 Gentlemen: The following steady-state power dissipation data for the Initial Blood Storage Experiment (IBSE) is based on test data recorded during the period of time from October 2, 1985 to October 9, 1985, with a 28V dc power source. Power 1.0 Air pressure: 14.7 psia 1.1 $T_{air} = 65^{\circ}F$ 20W C/C Module (4215-1400) 18W C/W Module (4215-1200) 1.2 $T_{air} = 75^{\circ}F$ 31W C/C Module (4215-1400) 25W C/W Module (4215-1200) 1.3 $T_{air} = 80^{\circ} F$ 50W C/C Module (4215-1400) 52W C/W Module (4215-1200) 1.4 $T_{air} = 90^{\circ}F$ (max.power dissipation capability of Modules @ 28 Vdc) 81W C/C Module (4215-1400) 70W C/W Module (4215-1200)

> ORIGINAL PAGE IS OF POOR QUALITY

Brussels Houston	Madrid Mexico City Milan	San Francisco São Paulo Singapore	Toronto Washington Wiesbaden
London	Milan	Singapore	Wiesbade

Arthur D. Little. Inc.

McDonnell Douglas Technical Services, Inc. Attention: Mr. Stephen Gotch

2.0 Air Pressure: 10.2 psia

2.2
$$T_{air} = -75^{\circ}F$$

C/C Module (4215-1400) 34W
C/W Module (4215-1200) 44W

If you have any questions, please feel free to call me.

Sincerely, David W. Almgren IBSE Engineering Manager

dwa:ln





.

Acorn Park Cambridge. Massachusetts 02140-2390 617 864-5770 Telex 921436

Arthur D. Little, Inc.

October 10, 1985

- To: D. MacN. Surgenor IBSE Principal Investigator
- From: D. Almgren IBSE Engineering Manager

Subject: Quick Look Report for IBSE Thermal Performance Tests

During the one week period from October 2, 1985 to October 9, 1985, a cold/cold IBSE module (S/N 0005) and a cold/warm IBSE module (S/N 0001) were tested over a range of temperatures at two pressures: 14.7 and 10.2 psia. This Quick Look Report represents a summary of the steady state raw data as recorded during the tests.

I.	Steady State Cond	lition No. 1	T air	- 65°F	Pressure = 14 .7 psia
	Date: October	3, 1985 @ 1	200		
	C/C Module:	cold plate cold plate	No. 1: No. 2: T: T _{air} : Power:	4.4°C 4.3°C 1.3°C 6.1°C 20W	(top to bottom) (entering at bottom)
	C/W Module:	cold plate cold plate P	(C): (W): Power:	4.2°C 21.0°C 18W	

OR	ENAL	PAGE	IS
OF	POOR	QUALI	TΥ

id San Francisco co City São Paulo 1 Singapore Tokyo	Toronto Washington Wiesbaden
	rid San Francisco co City São Paulo n Singapore s Tokyo

Arthur D. Little. Inc.

October 10, 1985 Page 2 T_{air} = 75[°]F Pressure = 14.7 psia II. Steady State Condition No. 2 Date: October 4, 1985 @ 0840 4.5°C cold plate No. 1: C/C Module: 4.3°C 1.4°C cold plate No. 2: (top to bottom) T: 6.4°C (entering at bottom) T_{air}: 31W Power: 4.4[°]C 21.3[°]C cold plate (C): C/W Module: cold plate (W): 25W Power: III: Steady State Condition No. 3 T = 90°F Pressure = 14.7 psia Date: October 5, 1985 @ 1030 9.6°C 8.7°C 2.1°C (top to bottom) 11.7°C (entering at bottom) cold plate No. 1: C/C Module: cold plate No. 2: **T**: T_{air}: 81W Power: 10.2°C 24.7°C cold plate (C): C/W Module: cold plate (W): 70W Power: IV. Steady State Condition No. 4 $T_{air} = 80^{\circ}F$ Pressure = 14.7 psia Date: October 7, 1985 @ 1430 4.6°C 4.4°C 1.9°C cold plate No. 1: C/C Module: cold plate No. 2: (top to bottom) Τ: 7.2°C (entering at bottom) T_{air}: 50W Power: 5.6°C cold plate (C): C/W Module: 21.5°C cold plate (W): 52W Power:

Arthur D. Little, Inc.

October 10, 1985 Page 3 V. Steady State Condition No. 5 T = 80°F Pressure = 10.2 psia Date: October 8, 1985 @ 0900 4.6°C 4.4°C 1.7°C (top to bottom) 7°C (entering at bottom) cold plate No. 1: C/C Module: cold plate No. 2: T: T_{air}: 49W Power: 6.4[°]C 21.5[°]C cold plate (C): C/W Module: cold plate (W): 54W Power: VI. Steady State Condition No. 6 $T_{air} = 75^{\circ}F$ Pressure = 10.2 psia Date: October 8, 1985 @ 1330 4.5°C 4.3°C cold plate No. 1: C/C Module: 4.3 C 1.6°C (top to bottom) 6.6°C (entering at bo cold plate No. 2: T: (entering at bottom) T_{air}: 34W Power: 4.5[°]C 21.5[°]C cold plate (C): C/W Module: cold plate (W): 44W Power: VII. Steady State Condition No. 7 $T_{air} = 65^{\circ}F$ **Pressure =** 10.2 psia Date: October 9, 1985 @ 1300 4.4°C 4.2°C cold plate No. 1: C/C Module: cold plate No. 2: 1.4°C (top to bottom) 6.3°C (entering at bottom) T: T_{air}: 22W Power: 4.4°C 21.5°C cold plate (C): C/W Module: cold plate (W): 21W Power:

David W. Almgren

Cc: R. Berthiaume K. Csigi W. Curby

P. Glaser

. .

IBSE Data Analysis

Donald Blevins, Ph.D. Nan Laird, Ph.D.

۰.

January 15,1987



~

JAN 1 6 1987

PRESIDENT'S OFFICE

PRECEDING PAGE BLANK NOT FILMED

Page 2

DATA ANALYSIS

Experimental design:

Each of three different blood components, red cells, white cells and platelets were subjected to microgravity aboard the space shuttle. Each were stored in three different types of storage bags, PVC DEHP, PVC TOTM and Polyolefin and placed in labeled Dewar flasks. Controls matched both for type of storage bag and Dewar position remained on earth.

Upon return to earth, the contents of the bags as well as the contents of the control bags were distributed in randomly numbered samples to individual investigators. The investigators had no knowledge of which were the control samples or of the different storage conditions of the samples. A total of 74 different outcome measurements were reported to us for statistical analysis.

Method of analysis:

The data received from the investigators was listed only by coded sample number. The first step was to match the coded samples with the original bag numbers to determine the information on bag type, orbital status, and storage location. The method of analysis used was analysis of variance, using type of bag and orbital status as the main factors. Most data sets included either 18 or 20 samples, hence allowed investigation of a possible interaction effect as well. The analysis tool used was SAS (Statistical Analysis System, release 82.4, SAS Institute, Cary, NC). For those data sets with fully balanced factors, the SAS ANOVA procedure was used, while for the others the GLM (General Linear Models) procedure was used.

One of the bags (bag number 56) was found to have developed a leak at some time before the samples were taken for distribution to investigators. For experiments involving that bag, the statistical analysis was redone with it omitted.

Results:

The results of the statistical analysis are summarized in the following. The data are listed alphabetically by the name of the principal investigator – often in multiple tables if many measurements were done by the same investigator. The first column of each table is the sample number as given to the investigator, the second number is the number of the bag from which the sample was taken, the third column is the label of the Dewar flask in which the bag was stored. Flasks <u>a</u> to <u>d</u> were on the orbiter while <u>e</u> to <u>h</u> remained on earth. In terms of storage configuration, <u>a</u> and <u>e</u>, <u>b</u> and <u>f</u>, <u>c</u> and <u>g</u>, <u>d</u> and <u>h</u> occupied corresponding positions. The fourth column of the table gives the orbital status; <u>l</u> indicates orbit and <u>0</u> indicates a control. The fifth column of the table indicates the type of bag used, <u>d</u> is PVC DEHP, <u>p</u> is Polyolefin, and <u>t</u> is PVC TOTM. The remaining column(s) contain the reported measurements. Following this is an explanation of the units of the measurement and a table which reports the mean of the measurement and the standard error of the mean (measurement of the accuracy of reported mean) for each combination of bag type and orbital and control samples. In addition, the overall means by bag type and earth/orbit status are reported in the column and row margins of the table. If the analysis of variance indicated statistical significance for one or more factors, the corresponding p-value is reported. If no such notation occurs, the factor was not significant (at the 5% level).

For experiments involving samples from bag number 56 as discussed above, the results of the analysis omitting that sample are reported in the same format, immediately following the original results.

30114c1t 293 7.8 54.3 3 302 6a1d 290 19.6 54.5 2 303 19f0p 468 12.8 57.1 3 304 17b1p 462 10.0 58.2 3 305 13e0t 478 13.6 53.8 3 306 1g0d 432 17.6 58.8 2 307 3c1d 444 3.8 63.5 3 308 8b1t 467 9.4 61.2 2 309 20a1p 158 10.8 60.7 2 310 15e0p 433 7.6 58.7 3311 312 9g0t 457 6.3 60.4 3311 313 16g0p 459 5.7 70.1 2314 4 e0d 388 4.6 69.1 249 314 4e0d 388 4.6 69.1 249	Sample	Bag	Dewar	<u>Orbit=1</u>	Type	<u>V1</u>	<u>2</u>	<u>_V3</u> _	_V4_
315 10 a 1 c 105 5.9 61.5 3 316 5 b 1 d 405 5.9 61.5 3 317 11 f 0 t 341 5.0 56.0 3 318 21 c 1 p 342 3.8 50.3 4	301 302 303 304 305 306 307 308 309 310 311 312 313 314 315 316 317 318	14 6 19 17 13 1 3 8 20 15 2 9 16 4 10 5 11 21	こった ゆきのこ ひっきん ひゅうせん	1 0 1 0 0 1 1 1 0 0 0 0 0 0 0 1 1 1 0	רסטרטסיטטיטטיטיטיט	293 290 468 462 478 432 444 467 158 433 436 457 459 388 269 405 341 342	7.8 19.6 12.8 10.0 13.6 17.6 3.8 9.4 10.8 7.6 8.3 6.3 5.7 4.6 9.7 5.9 5.0 3.8	54.3 54.5 57.1 58.2 53.8 63.5 61.2 60.7 58.7 64.7 60.4 70.1 69.1 56.9 61.5 56.0 50.3	37.9 25.9 30.1 31.8 32.6 23.6 32.7 29.4 28.5 33.7 27.0 33.3 24.2 26.3 33.4 32.6 39.0 45.9

Ausprunk - Transmission Electron Microscopy Data - Red Cells

~

.

V1 = Total Number of Cells Counted

Mean/S.E.M.

0	Type d	ρ	t	
UPDIC \	418 67/51 65	453.33/51.65	425.33/51.65	432.44/29.82
	379 67/51.65	320.67/51.65	343.00/51.65	347.78/29.82
	399.17/36.52	387.00/36.52	384.17/36.52	

V2 = % Normocytes

Mean/S.E.M.

Orbit \ Type d	Ρ	t	2 2 2 1 7 7 2 2
0 10,17/3,066	8.70/3.066	8.30/3.066	9.06/1.770
9 77/3,066	8,20/3.066	<u> </u>	8.98/1.770
9.97/2.168	8.45/2.168	8.63/2.168	

V3 = % Echinocytes type 1 & 11

	Þ	t	170
UFBIT (Type G	08 61.97/2.908	56.73/2.908	60.97/1.6/4
	08 56.40/2.908	57.47/2.908	57.90/1.6/4
62.02/2.0	56 59.18/2.056	57.10/2.056	

V4 = % Echinocytes type III

Mean/S.E.M.

Orbit	Type d	q	t	
	25 63/2 956	29.33/2.956	34.97/2.956	29.98/1.706
	30 40/2 956	35.40/2.956	33.57/2.956	33.12/1.706
	28.02/2.090	32.37/2.090	34.27/2.090	

Page 5 1

<u>Sample</u>	Bag	<u>Dewar</u>	<u>Orbit=1</u>	Type	<u>V1</u>	<u>V2</u>	<u>v3</u>	<u>V4</u>	<u>V5</u>
241	22	а	1	d	0	4	З	2	2.25
241	40	e	0	P	0	2	2	2	1.50
342	24	Ę	Ō	d	2	4	2	2	2.50
343	24	, F	ĩ	t	1	3	3	2	2.25
344	30	0	1	- -	2	2	З	2	2.25
345	39	- C	1	+	1	3	2	1	1.75
346	31	g	0	4	2	2	1	3	2.00
347	25	g	U	Û	1	2	i	2	1.75
348	41	Ь	1	p	1	נ ר	2	2	2 00
349	32	е	0	τ	1	5,	2	1	0 75
350	38	f	0	P	U	1	1	1	2 50
351	23	С	1	d	4	1	1	4	2.00
352	29	а	1	t	3	1	1	2	2.00
353	37	а	1	P	2	3	1	2	2.25
354	28	e	0	d	1	3	1	3	2.00
355	36	q	0	P	2	3	1	2	2.00
356	20	Ť	0	t	2	3	2	2	2.25
350	26	h	1	d	2	3	2	1	2.00
357	20	c C	1	t	2	3	2	3	2.50
300	54	C C	-						

Ausprunk - Transmission Electron Microscopy Data - White Cells

V1 = Degranulated Cells

- 0 = no cells
- 1 = occasional cells
- 2 = half of cells
- 3 = majority of cells
- 4 = all cells

Mean/S.E.M.

	Turne d	D	t	
Urbit 1		0 67/0 638	1,33/0,638	1.22/0.369
0	1.67/0.636	0.87/0.030	2 00/0 638	1.89/0.369
1	2.00/0.638	1.67/0.636	2.00/0.050	
	1.83/0.451	1.17/0.451	1.67/0.451	

V2 = Swollen Cytoplasm - same scoring as above

0	Type d	q	t	
	<u>1996</u> (1	2 00/0.577	3.00/0.577	2.67/0.333
	3.00/0.577	2 67 (0 577	2.33/0.577	2.56/0.333
1	2.67/0.577	2.07/0.07	2 67/0.408	
	2.83/0.408	2.33/0.408	2.0770.400	

V3 = Swollen Nuclear Envelope - same scoring as above

Mean/S.E.M.

	D	t	
<u>Orbit Vivpe d</u>	1 22/0 471	2.00/0.471	1.56/0.272
0 1.33/0.4/1	1.33/0.471	2 00/0 471	1,89/0,272
2.00/0.471	1.67/0.471	2.00/0.4/1	
1.67/0.333	1.50/0.333	2.00/0.333	

V4 = Clumped Chromatin

Mean/S.E.M.

	D	t	
<u>Orbit (Type u</u>	1 67/0 471	1,67/0,471	2.00/0.272
0 2.67/0.4/1	1.6770.471	2 67/0 471	2,44/0,272
2.33/0.471	2.33/0.4/1	2.67/0.471	
2.50/0.333	2.00/0.333	2.17/0.333	

V5 = Average Ranking

- · · ·	D	t	
<u>Drbit Type d</u>	1 42/0 204	2.00/0.204	1.86/0.118
0 2.17/0.204		2 25/0.204	2.19/0.118
1 2.25/0.204	2.08/0.204	2 13/0 144	
2.21/0.144	1.75/0.144	2.15/0.144	

Page 71

Sample	Bag	Dewar	<u>Orbit=1</u>	T∨pe	<u>V1</u>	<u>v2</u>	<u>87</u>	<u>V4</u>	
380	45	h	0	d	0	0	4	4	4.0
301	40 60	h	0	P	1	1	4	3	3.5
301	50	 5	Ô	t	1	1	4	З	3.5
382	<u> </u>	d	1	Ь	1	1	4	3	3.5
383	40 E (d	1	+			-	-	-
384	20	u -	1	n.	2	0	3	2	2.5
385	61	a	1	р 6	1	ñ	4	3	3.5
386	62	đ	1	t T	0	ň	Δ	Ā	4.0
387	47	h	0	0	0	0	4	7	3 5
388	59	h	0	P	1	0	4	1	A 0
389	43	d	1	d	U	0	4	4	4.0 5 E
390	55	h	0	t	1	0	4	3	5.5
391	54	d	1	t	2	0	2	1	1.5
392	63	h	0	P	1	0	4	3	3.5
303	53	Ь	1	t	1	0	4	3	3.5
201	44	h	1	d	0	0	4	4	4.0
394	44	ц Ь	0	d	0	0	4	4	4.0
395	40 E0		1	Ē	1	0	4	3	3.5
396	50	- U -	0	+	1	0	4	3	3.5
397	51	n	0	÷	1	ñ	4	3	3.5
398	98	n	0	ر ب	2	ň	4	2	3.0
399	99	d	1	τ	2	0	-+	<i>L</i>	2.0

Ausprunk - Transmission Electron Microscopy Data - Platelets

V1 = Pseudopods

- 0 = no cells
- 1 = occasional cells
- 2 = half of cells
- 3 = majority of cells
- 4 = all cells

Mean/S.E.M.

O_{2}	a	<u>t</u>	
	1 00/0.225	1.00/0.196	0.70/0.124
	1 33/0 225	1.67/0.226	1.11/0.139
	1.17/0.160	1.29/0.148	
0.17/0.160	1.17/0.160	1.29/0.148	

orbit p = 0.0400type p = 0.0003

V2 = Aggregated - same scoring as above

Orbit \ Type d	ρ	••••	
	0.33/0.231	0.25/0.200	0.20/0.12/
0 0.00/0.2	231 0.00/0.231	0.00/0.231	0.11/0.133
0.17/0.	163 0.17/0.163	0.14/0.151	

Page 8

V3 = Degranulated - same scoring as above

Mean/S.E.M.

0.545	Turne d	Þ	t	
Urbit 1		4 00/0 292	4,00/0,253	4.00/0.160
	4.00/0.292	4.00/0.202	2 23/0 292	3,67/0,169
1	4.00/0.292	3.67/0.292	3.33/0.202	
	4.00/0.207	3.83/0.207	3.71/0.191	

V4 = Swollen or Ruptured Cells - same scoring as above

Mean/S.E.M.

0 1 1 1 1	Turne d	D	t	
Urbit \	1ype u	3 00/0 292	3.00/0.253	3.30/0.160
	4.00/0.292	2 67/0 292	2,00/0,292	2.78/0.169
	3.67/0.292	2.8770.232	2 57/0 191	
	3.83/0.207	2.83/0.207	2.5770.151	

orbit p = 0.0428type p = 0.0014

V5 = Average of V3 and V4

Mean/S.E.M.

0-5:+ \	Type d	Q	t	
UPDIC V	100 (0 277	3 50/0 277	3,50/0,240	3.65/0.152
0	4.00/0.277	3.50/0.277	2 67/0 277	3,22/0,160
1	3.83/0.277	3.17/0.277	2.07/0.277	
<u></u>	3.92/0.196	3.33/0.196	3.14/0.182	

type p = 0.0284
Chao	
------	--

1

Sample	Bag	<u>Dewar</u>	<u>Orbit=1</u>	Type	<u>V1</u>
580	50	h	0	t	1.45
581	48	d	1	d	1.49
582	56	d	1	t	1.39
583	60	h	0	P	1.42
584	45	h	0	d	1.42
585	61	d	1	P	1.37
586	47	h	0	d	1.43
587	55	h	0	t	1.40
588	62	d	1	P	1.32
589	43	d	1	d	1.48
590	59	h	0	P	1.35
591	54	d	1	t	1.36
592	53	d	1	t	1.35
593	63	h	0	P	1.40
594	58	d	1	P	1.41
595	46	h	0	d	1.50
596	44	d	1	d	1.45
597	51	h	0	t	1.43
598	99	d	1	t	1.35
599	98	h	0	t	1.41

V1 = Platelet Count x 10⁹/ml

Mean/S.E.M.

Ombit \	h equi	a	t	
UPDIC (1 39/0 018	1.42/0.016	1.42/0.010
	1.45/0.018	1.37/0.019	1 36/0 016	1.40/0.010
	1.4//0.018	1.37/0.018	1.30/0.011	
	1.46/0.013	1.38/0.013	1.39/0.011	

type p = .0008

Omitting bag number 56: Mean/S.E.M.

Orbit	Type d	P	t	
	1.45/0.018	1.39/0.018	1.42/0.016	1.42/0.010
	1 47/0-018	1,37/0.018	1.35/0.018	1.40/0.010
	1.46/0.013	1.38/0.013	1.39/0.012	

type p = .0011

Chao 2

Sample	Bag	<u>Dewar</u>	<u>Orbit=1</u>	Type	<u>V1</u>	
Sample 580 581 582 583 584 585 586 587 588 589 590 591 592 593 594 595 594 595 596 597	Bad 50 48 56 45 61 47 55 43 54 53 53 53 53 58 44 51 59	Dewar h d d h h d h h d d h d h d h d h d h d	$ \begin{array}{c} 0 \\ 1 \\ 1 \\ 0 \\ 0 \\ 1 \\ 0 \\ 1 \\ 1 \\ 0 \\ 0 \\ 0 \\ 0 \\ 0 \\ 0 \\ 0 \\ 0 \\ 0 \\ 0$		0.000 0.000 0.088 0.000 0.000 0.000 0.000 0.000 0.000 0.000 0.000 0.000 0.000 0.000 0.000 0.000 0.000 0.000 0.000 0.000 0.000 0.000	· .
599	98	h	0	t	0.000	

VI = Platelet response to hypotonic Stress (Δ OD / 2 min)

Insufficient data to analyse

Chao 3

<u>Sample</u>	Bag	<u>Dewar</u>	<u>Orbit=1</u>	Туре	V1
580	50	h đ	0	t	9.1 10.3
581	40 E6	d	1	ť	7.8
582	20	ц р	, n	a	9.5
583	45	h	0	ď	10.6
584	4J 61		1	P	9.1
585	47	h	0	ď	11.3
586	47 55	 b	Õ	t	9.1
507	62	d	1	p	9.2
589	43	đ	1	d	10.1
500	59	h	0	P	9.6
591	54	d	1	t	8.4
592	53	đ	1	t	8.6
593	63	ħ	0	P	9.9
594	58	d	1	Þ	9.5
595	46	h	0	d	10.1
596	44	d	1	d	10.5
597	51	h	0	t	9.2
598	99	d	1	t	8.5
599	98	h	0	t	9.3

V1 = Mean Platelet Volume (µm³)

Mean/S.E.M.

	p	t	
Urbit (IVBE d	9 67/0 182	9,18/0,158	9.77/0.010
0 10.67/0.182	9.67/0.102	9 33/0 158	9,20/0.010
1 10.30/0.182	9.27/0.182	0.33/0.138	
10.48/0.129	9.47/0.129	8.75/0.112	

type p = .0001 orbit p = .0012

Omitting bag number 56: Mean/S.E.M.

O_{ab} is V_{ab} in O_{ab}	Þ	t	
UFBIT TYPE 0	9 67/0 162	9,18/0,141	9.77/0.089
0 10.67/0.162	<u> </u>	9 50/0 162	9.36/0.094
1 10.30/0.162	9.27/0.162	8.30/0.102	
10.48/0.115	9.47/0.115	8.83/0.106	

type p = .0001 orbit p = .0069

Chao

4

Sample	Baq	Dewar	<u>Orbit=1</u>	Type	<u>_v1</u> @
			-		_
580	50	h	0	t	6
581	48	d	1	d	0
582	56	d	1	t	75
583	60	h	0	Þ	0
584	45	h	0	d	0
585	61	đ	1	P	20
586	47	h	0	d	24
587	55	h	0	t	10
588	62	d	1	Þ	36
589	43	d	1	đ	21
590	59	h	0	P	0
591	54	d	1	t	38
592	53	d	1	t	20
593	63	h	0	P	8
594	58	d	1	P	0
595	46	h	0	d	0
596	44	d	1	d	0
597	51	h	0	t	19
598	9 9	d	1	t	21
599	98	h	0	t	0

V1 = Aggregation (by 12.5 µg collagen) % light transmission @ note "trace" interpreted as 0 for purposes of analysis

Mean/S.E.M.

0	Turne d	P	t	
Urbit V		2 67/9 185	8.75/7.954	6.70/5.031
	7 00 (9 185	18 67/9 185	38.50/7.954	23.10/5.031
	7.50/6.495	10.67/6.495	23.63/5.624	

orbit p = .0370

Omitting bag number 56: Mean/S.E.M.

$O=bi+ \lambda$	type d	q	t	
<u>Orbit</u>		2 67/6 731	8.75/5.829	6.70/3.687
	- 8.00/6.731	18 67/6 731	26.33/6.731	17.33/3.886
	7.00/6.731	10 67/4 758	16 29/4 406	
	1.50/4./50	10.0//4./00	10122/41100	

:hao

5

ample	Bag	Dewar	<u>Orbit=1</u>	<u>Type</u>	_v1@_
			_		-
580	50	h	0	t	U
581	48	d	1	d	5
582	56	d	1	t	44
583	60	h	0	P	11
584	45	h	0	d	5
585	61	đ	1	P	0
586	47	h	0	d	10
587	55	h	0	t	5
588	62	d	1	p	11
589	43	d	1	d	0
590	59	h	0	P	8
591	54	d	1	t	45
592	53	d	1	t	9
593	63	ĥ	0	Р	0
594	58	d	1	P	0
595	46	h	0	d	0
596	44	d	1	d	0
597	51	ĥ	Ō	t	0
598	99	d	ī	t	15
590	98	h	0	t	11
333	20		-	•	-

1 = Aggregation (by 100 µM ADP) % light transmission

e note "trace" interpreted as 0 for purposes of analysis

ean/S.E.M.

mbit \	Type d	P	t	
	5 00/5 709	6.33/5.709	4.00/4.944	5.00/3.127
	1 67/5 709	3,67/5,709	28.25/4.944	12.90/3.127
	3.33/4.037	5.00/4.037	16.13/3.496	

mitting bag number 56: ean/S.E.M.

$mhi + \lambda$	Type d	Þ	t	
	5 00/5 159	6.33/5.159	4.00/4.468	5.00/2.826
<u> </u>	1 67/5 159	3 67/5 159	23,00/5,159	9.44/2.979
I	3 33/3 648	5.00/3.648	12.14/3.377	

rbit p = 0.0957ype p = 0.0565nteraction p = 0.0299

Chao 6

<u>Sample</u>	Bag	Dewar	<u>Orbit=1</u>	<u>Type</u>	<u>V1</u>
580 581 582 583 584 585 586 587 586 587 588 589 590 591 592 593 592 593 594 595 596 597 598	50 48 56 60 45 61 47 55 62 43 54 53 63 54 53 63 58 46 44 51 99	,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	0 1 1 0 0 1 0 1 0 1 1 0 1 1 0 1 0 1 0 1 0 1 0 1 0 1 0 1 0 1 0 0 1 1 0 0 1 1 0 0 0 1 1 0 0 1 0 0 1 1 0 0 1 1 0 0 1 1 0 0 1 1 0 0 1 1 0 0 1 1 0 0 1 1 0 0 1 1 0 0 1 1 0 0 1 1 0 0 1 1 0 0 1 1 0 0 1 1 0 0 1 1 0 0 1 1 0 0 1 1 0 1 1 0 1 0 1 0 1 0 1 1 0 1 1 0 1 1 0 1 1 1 1 1 1 1 1 1 1 1 1 1	tdtpdpdtpdpttppddtt	0.00 0.00 2.94 0.00 0.23 0.00 0.29 0.00 0.29 0.00 1.73 0.17 0.00
599	98	h	0	t	0.00

V1 = ATP Release (n mol / 10⁹ platelets)

Mean/S.E.M.

	fune d	q	t	
		0 /0.373	0 /0.323	0 /0.205
<u> </u>		0,17/0,373	1.21/0.323	0.54/0.205
	0 /0.264	0.09/0.264	0.61/0.229	

Dmitting bag number 56: Mean/S.E.M.

<u>ант.</u>	and the second	D	t	
JEDIE V IV		0 (0 219	0 /0.189	0 /0.120
0	0 /0.219	0 /0.210	0 63/0 219	0.27/0.126
1	0 /0.219	0.17/0.219	0.83/0.213	
	0 /0.155	0.09/0.155	0.27/0.143	

Chao 7

Sample	Baq	Dewar	<u>Orbit=1</u>	Type		
580 581 582 583 584 585 586 587 588 587 588 589 590 591 592 593 594 595 594 595 596 597 598	50 48 56 60 45 61 47 55 62 43 59 54 53 63 58 46 44 51 99	ָּרְםָּרָ הַשְּׁרָה שְׁהָ שְׁהַ שְׁ	0 1 1 0 0 1 0 1 1 0 1 1 0 1 1 0 1 0 1 0		> 30.62 > > > > 43.99 > > >	90 90 29.59 90 90 90 90 90 90 90 90 90 90 90 90 90
599	98	n	U	C		-

V1 = Russell's Viperr Venom Time (sec.)

Insufficient data to analyse

<u>Sample</u>	Bag	Dewar	<u>Orbit=1</u>	Type	_ <u>V1</u>
E00	50	ħ	0	t	4.7
500	48	 d	1	d	4.5
581	40 54	d	1	t	57.3
582	20 60	ц Р	0	p	0.0
583	45	н Б	n n	ď	0.0
584	40	d	1	a	1.9
585	47	u b	Î.	ď	0.0
586	4/	н Б	Õ	ť	5.0
587	50	n d	1		0.0
588	62	ц а	1	h	6.1
589	43	ц ь	1	ñ	19.2
590	59	n	1	+	37.9
591	54	a	1	+	n n
592	53	a	1		0.0
593	63	h	U	4	1 9
594	58	d	1	р Г	21
595	46	h	0	a	2.1
596	44	d	1	a	3.2
597	51	h	0	τ	0.0
598	99	d	1	t	11.3
599	98	h	0	t	0.0

V1 = Serotonin Uptake % (10 min. incubation)

Mean/S.E.M.

Chao 8

a tith A Turna d	D	t	
Urbit (Type d	6 40/7 385	2.43/6.396	3.10/4.045
0 0.7077.385	1 27/7 385	26.63/6.396	12.41/4.045
1 4.60/7.385	3 83/5 222	14.53/4.522	
2.65/5.224	J.0J/J.444	,	

Omitting bag number 56: Mean/S.E.M.

	P	t	
Orbit \ Type d	<u> </u>	2.43/4.464	3.10/2.823
0 0.70/5.154	1 27/5 154	16.40/5.154	7.42/2.976
2.65/3.645	3.83/3.645	8.41/3.374	

~

Chao

9

Sample	<u>Bag</u>	<u>Dewar</u>	<u>Orbit=1</u>	Туре	<u></u>
580 581 582 583 584 585 586 587 588 589 590 591 592 593 594	50 48 56 60 45 61 47 55 62 43 59 54 53 63 58	,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	0 1 1 0 0 1 0 1 1 0 1 1 0 1 1 0 1	tdtpdpdtpdpttpp	0.0 0.0 35.9 0.0 0.0 0.0 0.0 0.0 0.0 0.0 29.8 0.0 0.0
595	46	h	0	d	0.0
596	44	d	1	d	0.0
597	. 51	h	U 1	τ +	0.0
598	99	d	1	د +	0.0
599	98	n	U	L	0.0

V1 = Serotonin Release % (4 min. after stimulation)

Insufficient data to analyse.

Chao 10

<u>Sample</u>	Bag	Dewar	<u>Orbit=1</u>	Туре	
580 581 582 583 584 585 586 587 588 589 590 591 592 593 594 593 594 595 596 597 598	50 48 56 60 45 61 47 55 62 43 59 54 53 63 58 46 44 51 99	עמעיעמיעקעמיקי <u>מ</u> ימ	0 1 1 0 0 1 0 1 1 0 1 1 0 1 0 1 0 1 0 1		62.32 53.93 36.96 57.87 65.31 40.78 74.87 59.51 42.98 55.07 43.62 35.04 49.27 52.70 50.46 48.58 82.80 42.65 48.28
599	98	h	0	t	49.32

 $\sqrt{1} = \beta - Thromboglobulin (\mu g / ml)$

Mean/S.E.M.

a a construction of the second s	p	t	
Drbit Vive d	51 40/5 894	53.45/5.104	55.68/3.228
0 62.92/5.094	44 74/5 894	42.39/5.104	49.56/3.228
	48,07/4,168	47.92/3.609	
63.43/4.100	401017 1011		

type p = .0257 Dmitting bag number 56: Mean/S.E.M.

	P	t	
<u>G2 92/6 033</u>	51.40/6.033	53.45/5.225	55.68/3.305
	44.74/6.033	44.20/6.033	50.96/3.483
63.43/4.266	48.07/4.266	49.48/3.950	

type p = .0434

Chao 11

Samole	Bag	Dewar	Orbit=1	Type	<u>vı@</u>
Jampre	<u>5-73</u>	<u></u>			
580	50	h	0	t	5.2
581	48	đ	1	d	4.0
582	56	d	1	t	1.3
583	60	h	0	P	7.4
584	45	h	0	d	16.8
585	61	d	1	P	25.0
586	47	h	0	d	12.0
587	55	h	0	t	5.3
588	62	d	1	p	1.5
589	43	đ	1	d	2.0
590	59	h	0	p	5.2
591	54	d	1	t	1.6
592	53	d	1	t	1.1
593	63	h	0	P	8.0
594	58	d	1	P	1.8
595	46	h	0	d	20.0
596	44	d	1	d	25.0
597	51	h	0	t	20.0
598	99	d	1	t	1.4
599	98	h	0	t	6.8

V1 = Thromboxane B2 Note "> 20" interpreted as 25 for purposes of analysis

Mean/S.E.M.

0	Type d	Þ	t	
UPDIL 1	16 27/4 574	6 87/4 574	9.33/3.961	10.67/2.505
	10.27/4.374	9 43/4 574	1.35/3.961	6.47/2.505
	10.33/4.3/4	9.45/4.5/4	5 34/2 801	
	13.30/3.234	0.15/3.254	5.54/2.001	

Omitting bag number 56: Mean/S.E.M.

$O_{rhit} \setminus T_{VDR} d$	q	t	
01010 16 27/4 746	6 87/4 746	9.33/4.111	10.67/2.600
10.27/4.746	9 43/4 746	1.37/4.746	7.04/2.740
13,30/3,356	8.15/3.356	5.91/3.107	

-

.

Chao 12

Sample	Bag	<u>Dewar</u>	<u>Orbit=1</u>	<u>Type</u>	_V1_
580	50	h	0	t	22.1
581	48	đ	1	đ	22.2
582	56	d	1	t	16.7
583	60	h	0	P	21.7
584	45	h	0	d	23.9
585	61	d	1	P	21.3
586	47	h	0	d	24.0
587	55	h	0	t	23.1
588	62	d	1	P	23.4
589	43	d	1	d	24.4
590	59	h	0	P	23.3
591	54	d	1	t	19.0
592	53	d	1	t	20.8
593	63	h	0	p	24.0
594	58	d	1	P	23.7
595	46	h	0	d	24.8
596	44	d	1	d	22.3
597	51	h	0	t	24.7
598	99	d	1	t	23.2
599	98	h	0	t	24.4
	-				

V1 = Lactate (mM)

Mean/S.E.M.

0-6:+ 1	Type d	D	t	
UPDIC 1	<u>1906</u> <u>0</u> <u>34</u> <u>32</u> (0 <u>937</u>	23 00/0 937	23,58/0.811	23.60/0.513
	22.23/0.937	22 80/0 937	19,93/0.811	21.70/0.513
	23.60/0.663	22.90/0.663	21.75/0.574	

orbit p = 0.0203

Omitting bag number 56: Mean/S.E.M.

$Orbit \land$	Type d	P	tt	
	24 23/0 768	23 00/0.768	23.58/0.665	23.60/0.421
	22.25/0.768	22 80/0 768	21,00/0,768	22.26/0.443
	22.97/0.788	22.00/0.543	22 47/0 503	
	23.60/0.543	22.90/0.545		

orbit p = 0.0465

Chao 13

Sample	<u>Baq</u>	<u>Dewar</u>	<u>Orbit=1</u>	Type	<u></u>
580	50	h	0	t	38
581	48	d	1	d	23
582	56	d	1	t	141
583	60	h	0	p	47
584	45	h	0	d	33
585	61	d	1	p	24
586	47	h	0	d	38
587	55	h	0	t	22
588	62	d	1	P	26
589	43	d	1	d	20
590	59	h	0	P	28
591	54	d	1	t	86
592	53	d	1	t	22
593	63	h	0	P	23
594	58	d	1	P	14
595	46	h	0	d	25
596	44	d	1	d	21
597	51	h	0	t	20
598	99	d	1	t	18
599	98	h	0	t	25

V1 = Glucose (mg / 100 ml)

Mean/S.E.M.

Orbit \	Type	d	P	t	
0,0,0,0	32 0	0/16,149	32,67/16,149	26.25/13.986	29.90/8.845
	21 3	3/16 149	21,33/16,149	66.75/13.986	39.50/8.845
l	26.6	7/11.419	27.00/11.419	46.50/9.889	

•

.

......

Omitting bag number 56: Mean/S.E.M.

Orbit >	Type d	p	t	
	32 00/9 611	32 67/9 611	25.25/8.324	29.90/5.264
	21 33/9 611	21 33/9.611	42.00/9.611	28.22/5.549
<u> </u>	26.67/6.796	27.00/6.796	33.00/6.292	

Curby	i –	Red	Cells	Stress	Test
	•				

Sample	Bag	Dewar	<u>Orbit=1</u>	<u>Type</u>	<u>V1</u>	V2	<u>v3</u>
<u>Sample</u> 401 402 403 404 405 406 407 408 409 410 411 412	Bag 17 1 13 19 14 6 8 2 20 9 3 15	<u>Dewar</u> b g e f c a b f a g c e	<u>Orbit=1</u> 1 0 0 1 1 1 0 1 0 1 0 1 0 0	<u>Type</u> Pdt ptdtd ptd p	V1 27.400 133.000 332.200 248.200 184.600 0.000 149.600 140.400 75.600 148.000 268.600 0.000	V2 14.810 58.620 204.700 123.100 104.000 0.000 91.770 92.630 42.100 108.600 173.300 0.000 60.390	V3 22 23 23 16 22 18 24 30 21 22 23 28 23
413 414 415 416 417 418	11 4 10 16 21 5	f e g c b	0 0 1 0 1 1	t d t p p d	8.199 141.200 256.000 163.200 65.200	5.827 90.080 26.070 58.180 45.020	23 12 20 23 27

V1 = Loss at 20 min., count/m1 x 10^6

Mean/S.E.M.

Drbit \ Typedpt093.87/60.323168.07/60.323197.67/60.323153.20/34.8281111.27/60.32388.73/60.323158.47/60.323119.49/34.828102.57/42.655128.40/42.655178.07/42.655

V2 = Oversized Particle Loss at 20 min., count/ml x 10^4

Mean/S.E.M.

		D	t	
<u>Orbit \</u>		49 72/33,452	124.56/33.452	75.55/19.314
	-52.36/33.452	38.36/33.452	95.28/33.452	68.81/19.314
	62.57/23.654	44.04/23.654	109.92/23.654	

V3 = Peak Channel Shift at 20 min.

a the transferred	D '	t	
Drbit Vivpe 0	21 22/2 5/6	22.67/2.546	23.11/1.470
0 25.33/2.546	21.33/2.546	10 23/2 5/6	21,33/1,470
1 22.67/2.546	22.00/2.546	19.3372.340	
24.00/1.800	21.67/1.800	21.00/1.800	

Sample	Bag	Dewar	<u>Orbit=1</u>	Type	<u>V1</u>	V2	<u> </u>	V4
4.4.1	24	f	0	d	3.688	9.588	90	0.225
441	25	'n	ĩ	t	3.802	10.670	94	1.010
442	20	0	n	÷	2.312	10.100	95	1.041
443	31	g	1	Ă	2 558	9.372	96	0.936
444	22	a	1 0	<u> </u>	5 807	18.630	110	0.635
445	40	e	0	2	3 964	11.080	100	0.712
446	39	С	1	μ +	6 086	15 600	105	1.064
447	29	а	1	ل ر	4 090	10 500	89	0.991
448	25	g	U	0	4.000	20 800	104	0 779
449	32	e	0	τ	7.070	20.800	104	0 514
450	38	f	0	p	3.263	9.100	90	0.314
451	41	Ь	1	p	2.8//	7.132	90	0.372
452	23	С	1	d	2.624	9.736	87	0.000
453	33	f	0	t	4.606	11.500	95	0.724
454	36	g	0	P	3.421	10.240	94	1.040
455	28	e	0	d	7.360	20.440	89	0.785
456	26	b	1	d	2.976	10.350	89	1.107
457	34	-	1	t	4.126	10.070	90	1.060
458	37	a	1	P	6.831	22.390	99	0.726
			7					

Curby 2 - White Cell Distribution

V1 = 1 - 254 count/ml x 10'

Mean/S.E.M.

$Orbit \lambda$	Type d	P	t	
	5 04/1 030	4,16/1,030	4.87/1.030	4.69/0.595
	2 72/1 030	4.56/1.030	4.67/1.030	3.98/0.595
<u>I</u>	3.88/0.728	4.36/0.728	4.77/0.728	

V2 = > 254, count/m1 x 10⁶

Mean/S.E.M.

$Orbit \land$	Type d	P	t	
	12 51/3 067	12 66/3.067	14.13/3.067	13.43/1.771
	13.51/3.00/	13 53/3 067	12.11/3.067	11.82/1.771
<u> </u>	9.82/3.06/	12.10/2.168	13 12/2 168	
	11.66/2.168	13.10/2.100	13.12/20100	

V3 = Peak Channel

Mean/S.E.M.

$Orbit \Lambda$	Type d	q	t	
	89 33/3 413	100.67/3.413	98.00/3.413	96.00/1.970
	90 67/3 413	96.33/3.413	96.33/3.413	94.44/1.970
i	90.00/2.413	98.50/2.413	97.17/2.413	

V4 = Delay Factor - not analysed

Curby 3 - Platelets Cell Distribution

<u>Sample</u>	Bag	<u>Dewar</u>	<u>Orbit=1</u>	Type	<u></u>	<u>V2</u>	<u>v3</u>	
<u>Sample</u> 481 482 483 484 485 486 487 488 489 490 491 492 493 494 495 496 497	<u>Bag</u> 50 48 61 60 45 43 55 47 62 59 54 58 51 46 53 44 63	Dewar h d d h h d h h d d h h d d h .	<u>Orbit=1</u> 0 1 1 0 0 1 0 1 0 1 1 0 0 1 1 1 0 0	<u>Type</u> tdppddtdpptdtdtdp+	V1 1.249 1.632 1.465 1.027 1.306 1.699 2.258 2.406 1.161 1.645 1.496 1.411 1.428 1.663 1.279 1.483 1.280 1.654	V2 3.4 4.6 6.4 5.6 6.8 6.4 3.4 2.8 7.8 6.2 10.2 7.2 5.8 6.2 10.2 5.6 8 6.8 6.2 10.2 5.6 6.8 6.2 10.2 5.6 6.8 6.2 10.2 5.6 6.8 6.2 10.2 5.6 6.8 6.2 10.2 5.6 6.8 6.2 10.2 5.6 6.8 6.2 10.2 5.6 6.8 6.2 10.2 5.6 6.8 6.2 10.2 5.6 6.8 6.2 10.2 5.6 6.8 6.2 10.2 5.6 6.8 6.2 10.2 5.6 6.8 6.2 10.2 5.6 6.8 6.2 10.2 5.6 6.8 6.2 10.2 5.6 6.8 6.2 10.2 5.6 6.8 6.2 10.2 5.6 6.8 6.2 7.2 5.8 6.8 6.2 7.2 5.8 6.8 6.8 6.2 7.2 5.8 6.8 6.8 6.2 7.2 5.8 6.8 6.8 6.2 7.2 5.8 6.8 6.8 6.8 6.2 7.2 5.8 6.8 6.8 6.8 6.2 7.6 6.8 6.8 6.2 7.6 6.8 6.8 6.2 7.2 5.8 6.8 6.8 6.8 6.8 6.8 6.8 6.8 6	$\frac{\sqrt{3}}{12}$ 12 12 14 14 15 15 23 13 14 14 13 14 12 12 12 12	5.215 7.512 6.978 7.764 7.055 7.372 1.475 1.541 6.295 6.338 6.663 2.588 6.590 6.700 7.435 6.211 6.997 6.691
498 499	96 99	đ	1	ť	1.524	9.2	12	5.813

V1 = 1 - 254 count/ml x 10^9 Mean/S.E.M.

Orbit \ Type	h d	P	<u>t</u>	
	79/0 197	1 32/0 197	1.65/0.171	1.59/0.108
	./9/0.197	1 35/0 197	1,43/0,197	1.46/0.144
	.60/0.19/	1 33/0 139	1 56/0.129	
1	./0/0.139	1.33/0.139	1.00/01/20	

V2 = > 254, count/ml x 10^6

Mean/S.E.M.

	Tupo d	Þ	t	
Urbit 1	5 53/1 195	6.33/1.195	6.70/1.035	6.24/0.654
	6 40/1 195	5.13/1.195	7.07/1.195	6.20/0.690
1	5.97/0.845	5.73/0.845	6.86/0.782	

V3 = Peak Channel

Mean/S.E.M.

$Orbit \land$	Type d	p	t	
01010 1	17 00/1.371	13.33/1.371	13.00/1.188	14.30/0.751
	13 00/1.371	13.00/1.371	13.33/1.371	<u>13.11/0.74</u>
<u>i</u>	15.00/0.970	13.17/0.970	13.14/0.898	

V4 = Delay Factor - not analysed

Curby 4 - Red Cells Cell Distribution

Sample	Bag	Dewar	<u>Orbit=1</u>	<u>Type</u>	<u>V1</u>	<u>V2</u>	<u>_V3</u>	_ <u>V4</u>
<u>Samole</u> 401 402 403 404 405 406	<u>Bag</u> 17 13 19 14 6	Dewar b g e f c a	1 0 0 1 1	p d t p t d	3.388 3.871 3.446 4.324 3.435 3.898	7.84 14.02 12.66 67.62 14.26 13.46	110 108 106 140 106 106	1.564 0.491 1.315 1.735 1.331 1.271 1.257
407 408 409 410 411 412 413 414 415 416 417 418	8 20 9 3 15 11 4 10 16 21 5	bf agcefeagcb	1 0 1 0 0 0 1 0 1 1 1	רס סר מ סר מר ם ס מ	3.554 3.109 3.249 3.295 3.737 3.906 3.468 3.361 3.421 3.552 3.839 3.152	14.36 21.06 10.76 16.88 21.64 18.58 12.90 8.74 10.28 10.52 10.80 4.58	112 105 110 112 111 108 108 100 106 108 100	1.633 1.377 1.357 1.552 1.346 1.321 1.374 1.339 1.321 1.539 1.642

V1 = 1 - 254 count/ml x 10^9

Mean/S.E.M.

	h eevT	q	<u>t</u>	
UPDIC 1	<u>- 1996 - 0 177</u>	3 93/0,177	3.40/0.177	3.59/0.102
	3.45/0.177	3 49/0 177	3,46/0,177	3.52/0.102
	3.60/0.177	2.29/0.175	3 43/0 125	
	3.52/0.125	3.71/0.125	5.45/81.20	

V2 = > 254, count/m1 x 10^7

Mean/S.E.M.

0	Turne d	Þ	t	
UPDIC 1		32 24/7 747	14.15/7.747	20.33/4.473
	14.01/7.747	9 80/7 747	13.03/7.747	12.02/4.473
1	13.22/1.141	21 02 (5 478	13 59/5,478	
	13.92/5.4/8	21.02/5.4/0	13:33/3000	

V3 = Peak Channel

Mean/S.E.M.

O-bi+	Type d	p	t	
	109 33/0 177	119,00/0,177	108.00/0.177	112.11/0.102
<u> </u>	105.33/0.177	107.67/0.177	105.33/0.177	106.33/0.102
	107.67/0.125	113.33/0.125	106.67/0.125	

V4 = Delay Factor - not analysed

Jacobson 1 - White Blood Cells

Sample	<u>Baq</u>	Dewar	<u>Orbit=1</u>	Type	<u>V1</u>	<u>V2</u>	<u></u> ¥3
871	22	а	1	d	6.945	42.2	228.1
872	35	ъ	1	t	7.073	25.5	263.7
873	31	g	0	t	7.092	23.9	259.4
874	24	Ē	0	d	6.965	38.7	220.5
875	39	с	1	P	7.051	27.9	259.1
876	40	e	0	P	7.093	24.9	255.3
877	25	a	0	d	6.931	41.5	220.5
878	32	e	0	t	7.101	22.0	263.4
979	29	a	1	t	7.087	24.5	252.6
075	38	f	0	p	7.052	24.5	269.9
000	11	'n	1	p	7.095	24.6	236.5
001	23	c c	1	ď	6.893	44.5	226.6
882	25	د م	1	D	7.031	28.8	261.5
883	27	a 0	n	Ā	6.907	42.2	225.5
884	20	e	1	+	7.079	25.5	259.0
885	34	C E	1	+	7 067	23.8	248.0
886	33	г	0	с Г	7 039	26.5	253.2
887	36	g	U	ч ч	2 030	12 A	240.7
888	26	Ь	1	a	0.930	46.4	240.1

V1 = pH

Mean/S.E.M.

	D	t	
Urbit Type d	7 06/0.015	7.09/0.015	7.03/0.008
	7,06/0,015	7.08/0.015	7.02/0.008
6.93/0.010	7.06/0.010	7.08/0.010	

type p = 0.0001

```
V2 = PCO_2 - mm Hg
```

Mean/S.E.M. .

0-6	 ×	Tuna	

a , · + ,	Turne d	n	t	
Urbit \	10 00 (0 835	25 30/0.835	23.23/0.835	29.77/0.482
	40.00/0.035	27 10/0.835	25.17/0.835	31.77/0.482
	41.92/0.590	26.20/0.590	24.20/0.590	·

type p = 0.0001 orbit p = 0.0129

$$V3 = PO_2 - mm Hg$$

Mean/S.E.M.

O-hit \ Type d	p	t	<u> </u>
UPDIT. 1990 022 17/4 022	259 17/4 923	256.93/4.923	246.1972.842
0 222.17/4.925		259 13/1 923	247.53/2.842
1 231.00/4.923	252.37/4.923	238.43/4.925	
226.98/3.481	255.92/3.481	257.68/3.481	

type p = 0.0001

.

Jacobson 2 - Red Blood Cells

<u>Sample</u>	Bag	Dewar	<u>Orbit=1</u>	Type	<u>V1</u>	<u>_V2</u>	<u>_V3</u>
161	13	е	0	t	6.763	122.1	52.4
162	17	b	1	p	6.707	134.5	49.9
163	19	f	0	P	6.743	133.8	50.1
164	1	a	0	d	6.711	155.2	37.0
165	6	3	1	d	6.738	149.4	36.6
165	14	C C	ī	t	6.730	133.9	48.2
160	, 1	C C	1	d	6.672	169.8	34.1
167	G	C C	0	t	6.747	128.2	48.8
160	0	g S	1	t	6.766	129.0	49.8
169	2	E E	· n	ď	6.710	157.5	36.4
170	15	,	0 0	5	6.740	134.3	45.8
171	10	e	1	г р	6.712	143.4	51.0
172	20	a	1	ч с	6 744	145.8	35.7
173	4	e	0	Ğ	6 758	132.1	51.5
174	16	g	U	р d	6 695	161 2	35.9
175	5	D	1	<u>u</u>	6.075	126 0	47 8
176	11	f	0	τ	6.767	120.0	46.2
177	21	С	1	P	6.728	143.5	40.2
178	10	а	1	t	6.742	134.1	41.6

V1 = pH

Mean/S.E.M.

$Orbit \land$	Type d	p _	<u>t</u>	
	6 72/0 011	6.75/0.011	6.76/0.011	6.74/0.006
	6 70/9 611	6.72/9.611	6.75/9.611	6.72/0.006
1	6.71/0.008	6.73/0.008	6.75/0.008	

type p = 0.0099 orbit p = 0.0336

 $V2 = PCO_2 - mm Hg$

Mean/S.E.M.

Orbit	Type d	p	t	
	152 93/3 243	133 40/3,243	125.43/3.243	137.22/1.872
	102.03/3.243	140 47/3 243	132.33/3.243	144.31/1.872
	160.13/3.245	140.47/3.293	128 88/2 293	
	156.48/2.293	136.93/2.295	120.00/2.293	

type p = 0.0001 orbit p = 0.0201

$$V3 = PO_2 - mm Hg$$

Mean/S.E.M.

Orbit \	Type d	p	<u>t</u>	
010101	36 37/1 164	49.13/1.164	49.67/1.164	45.06/0.6/2
	35 53/1 164	49.03/1.164	48.53/1.164	44.37/0.672
<u>i</u>	35.95/0.823	49.08/0.823	49.10/0.823	

type p = 0.0001

Jacobson 3 - Platelets

.

Sample	Bag	Dewar	<u>Orbit=1</u>	Type	<u>V1</u>	_ <u>V2</u> _	<u>V3</u>
Sample 140 141 142 143 144 145 146 147 148 149 150 151 152 153 154 155 156 157	Bac 61 48 50 56 60 45 55 43 47 62 51 63 84 53 63 84 53 60	<u>ספשר</u> סמדםדעסדמדבדעםעדנ	<u>Orbit=1</u> 1 0 1 0 1 0 1 0 1 0 1 0 1 0 1 1 0 1 1 0 0 1 0 0 1 0 0 1 0 0 1 0 0 1 0 0 1 0 0 0 1 0 0 1 0 0 1 0 0 0 1 0 0 1 0 0 0 1 0 0 0 0 0 1 0 0 0 0 1 0 0 0 0 0 0 0 1 1 0 0 0 0 0 0 0 0 0 0 0 0 0	pdtt pdtt dd ppt ppd t dt	5.811 5.663 5.833 6.544 5.880 5.712 6.237 5.748 5.659 5.734 5.814 5.826 5.732 5.722 5.729 5.662 5.781 5.680 5.763	20.3 10.8 7.3 22.6 11.5 9.2 29.9 7.2 8.7 8.9 18.3 8.0 6.8 6.4 10.6 10.8 12.1 8.6 12.5	170.7 185.2 221.9 186.9 202.1 195.3 145.5 201.4 194.2 207.6 189.1 214.2 211.0 206.8 203.4 184.8 191.6 203.5 189.7
159	98	h	0	t	5.753	6.5	205.0

V1 = pH

Mean/S.E.M.

	an d	D	t	
UPDIL A IS	5 71 (0 104	5 81/0,104	5.77/0.090	5.76/0.057
	<u>5.7170.104</u>	5 78/0.104	6.08/0.090	5.87/0.057
	5.69/0.074	5.80/0.074	5.92/0.064	

omitting bag number 56: Mean/S.E.M.

0-+i+	Type d	p	t	
<u>Orbit</u>	5 71/0 066	5.81/0.066	5.77/0.057	5.76/0.036
	<u> </u>	5,78/0,066	5.93/0.066	5.79/0.038
	5.69/0.047	5.80/0.047	5.84/0.043	

 $V2 = PCO_2 - mm Hg$

Mean/S.E.M.

Orbit \ Type d	P	<u>t </u>	0.01/1.01/2
0 8 90/2 632	8.63/2.632	6.95/2.279	8.04/1.44
10 10/2 632	16,40/2,632	19.28/2.279	15.65/1.442
9.50/1.861	12.52/1.861	13.11/1.612	

orbit p = 0.0022

omitting bag number 56: Mean/S.E.M.

$Orbit \land$	Type d	P	tt	
	8 90/2 661	8,63/2,661	6.95/2.305	8.04/1.458
	10 10/2 661	16 40/2 661	18.17/2.279	14.89/1.536
	0.50/1.882	12 52/1.882	11.76/1.742	
	9.50/1.882	12.52/1.002	11.70/1.74	

orbit p = 0.0065

 $V3 = PO_2 - mm Hg$

Mean/S.E.M.

$Orbit \lambda$	Type d	p	t	
<u>Orbic A</u>	1,02 12/7 620	207 70/7 630	209.83/6.607	206.88/4.179
0	202.13/1.630	201.10/1.030	178 42/6 607	184 11/4,179
1	188.07/7.630	187.73/7.630	178.4278.007	
	195.10/5.395	197.72/5.395	194.13/4.672	

orbit p = 0.0018

omitting	bag	number	56:
Mean/S.E.	Μ.		

<u> </u>	Tuna d		D	t	
UPDIE 1		761	207 70/7 761	209,83/6,721	206.88/4.251
	202.13/1	. /61	207.7077.701	175 60/7 761	183 80/4.481
1	188.07/7	.761	18/./3//./61	1/5.60/7.751	100.007
	195.10/5	.488	197.72/5.488	195.16/5.081	

orbit p = 0.0025

Page 32[°]

.

Jacobson 4

Sample	Bag	<u>Dewar</u>	<u>Orbit=1</u>	<u>Type</u>	<u>V1</u>
81	13	е	0	t	10.0
82	1	a	0	đ	7.0
83	6	a	1	d	7.0
84	19	f	0	P	8.0
85	17	Ь	1	P	7.0
86	14	- C	1	t	10.0
87		b	1	t	10.0
88	15	ē	0	P	8.0
80	9	- a	0	t	10.0
09 00	2	e G	1	d	8.0
90	20	a	1	P	8.0
51	20	f	Ō	ď	10.0
72	5	, h	1	d	10.0
95	10	a	1	t	14.0
94	21	с С	1	Þ	10.0
30	16		Î Î		10.0
96	10	g F	ñ	ť	10.0
97		1	0	Å	7.0
98	4	e	0	u	,

V1 = Plasma Hemoglobin mg/d1

-

Mean/S.E.M.

Orbit \ Type d	p	t	
01010 (19 <u>Pe</u> 00/0 892	8,67/0,892	10.00/0.892	8.89/0.515
0 8.00/0.092	8 33/0 892	11.33/0.892	9.33/0.515
8,17/0,631	8.50/0.631	10.67/0.631	

type p = 0.0324

Sample	Bag	<u>Dewar</u>	<u>Orbit=1</u>	Type	<u>V1</u>	<u>2</u>	<u>_V3</u>		<u>5</u>
<u>Sample</u> 61 62 63 64 65 66 67 68 69 70 71 72 73 74 75	Bag 1 14 6 19 13 17 3 20 8 13 15 16 5 10	<u>Dewar</u> g c af e b c f a b e e g b af	<u>Orbit=1</u> 0 1 1 0 0 1 1 1 0 1 1 0 0 0 1 1 1 0 0 0 1 1 1 0 0	<u>Type</u> dt d pt p d d pt t p p d t t	V1 4.1 3.8 4.2 3.8 3.8 3.8 3.8 3.8 3.8 3.8 3.8 3.8 3.8	V2 3.1 1.8 3.1 3.8 2.5 3.1 4.9 3.8 1.9 2.5 4.4 1.9 2.5 3.8 1.9 2.5 2.5 2.8 2.9 2.5 2.8 2.9 2.5 2.8 2.9 2.5 2.8 2.5 2.8 2.9 2.5 2.8 2.9 2.5 2.8 2.9 2.5 2.8 2.9 2.5 2.8 2.9 2.5 2.8 2.9 2.5 2.8 2.9 2.9 2.8 2.9 2.9 2.8 2.9 2.9 2.9 2.9 2.9 2.9 2.9 2.9	V3 306.7 310.7 323.5 308.0 304.5 313.5 312.7 298.0 304.7 307.0 302.0 302.0 302.0 298.7 305.0 321.3	V4 38.6 38.2 38.5 38.3 38.5 38.4 38.6 38.3 38.6 38.3 38.0 38.5 38.5 38.5 38.5 38.5 38.2 38.2 38.2 38.2	V5 111.0 111.0 112.0
76 77 78	4 21	e	0	d P	4.2	3.1 3.1	337.0 317.0	38.1 38.2	112.0 111.0

Jacobson 5 - Red Blood Cells

V1 = ATP µmoles / gm hgb

Mean/S.E.M.

	Turne	Þ	t	
Urbit 1	<u>1 02 (0, 121</u>	3 70/0.121	3.67/0.121	3.80/0.070
	2.03/0.121	3 93/0.121	3.70/0.121	3.86/0.070
	3.98/0.086	3.82/0.086	3.68/0.086	

V2 = 2,3DPG µmoles / gm hgb

Mean/S.E.M.

	D D	t	
$\frac{\text{Orbit}}{2}$	48 3.13/0.548	2.73/0.548	3.19/0.316
-0 $3.10/0.5$	48 3.13/0.548	2.50/0.548	2.84/0.316
3.30/0.3	87 3.13/0.387	2.62/0.387	

V3 = Glucose mg/dl

0-5:+ \	t eevt	p	t	
UPDIC Y	313 90/6 454	303, 33/6, 454	309.27/6.454	308.83/3.726
<u> </u>	211 62/6 454	311 73/6 454	307.57/6.454	310.31/3.726
	311.03/6.434	307 53/4 564	308.42/4.564	
	312.11/4.004			

V4 = Potassium mlq/liter

Mean/S.E.M.

0-hi+1	Type	Ь	Ρ	t	
	20	43/0 109	38,33/0,109	38.40/0.109	38.39/0.063
		27/0.109	38 30/0 109	38,13/0,109	38.27/0.063
		37/0.10#	29 32/0 077	38 27/0.077	
	38.	40/0.0//	33.32/0.077	50.2770.077	

V5 = Sodium mlq/liter

Orbit	Type (đ	р	t	
0	111.67	/0.272 1	11.67/0.272	111.67/0.272	111.67/0.157
1	112.00	/0.272 1	11.67/0.272	111.00/0.272	111.06/0.13/
	111.83	/0.192 1	11.67/0.192	111.33/0.192	

Jacobson 6 - Red Blood Cells

<u>Sample</u>	Bag	Dewar	<u>Orbit=1</u>	Type	<u></u>	<u>_V2</u>	<u>_V3</u>
81	13	е	0	t	10.0	12.0	167.0
82	1	a	0	d	7.0	11.7	167.0
83	6	a	1	d	7.0	12.7	168.0
84	19	f	0	p	8.0	11.7	167.0
85	17	b	1	P	7.0	12.2	168.0
86	14	C C	1	ť	10.0	11.8	169.0
87	. , A	ĥ	1	t	10.0	12.8	166.5
88	15	6	Ō	р	8.0	11.8	168.0
89	9	ā	Ō	ť	10.0	11.8	169.0
90	- -	E E	1	d	8.0	12.1	168.0
90	20	a	1	P	8.0	12.0	169.0
91	2	f	0	ď	10.0	11.7	169.0
92	5	'n	1	d	10.0	11.9	169.0
93	10	a	1	t	14.0	12.3	168.0
94	21	c C	1	p	10.0	12.6	168.0
95	16	0	0	Þ	10.0	12.5	169.0
90 70	11	9 f	Õ	ť	10.0	12.6	169.0
98	4	ė	õ	d	7.0	12.5	169.0

V1 = Plasma Hemoglobin mgs %

Mean/S.E.M.

$Orbit + \lambda$	Type d	p	t	
	8 00/0 892	8,67/0,892	10.00/0.892	8.89/0.515
	8 33/0 892	8,33/0,892	11.33/0.892	9.33/0.515
	8.17/0.631	8.50/0.631	10.67/0.631	

type p = 0.0324

V2 = Plasma Potassium mlq/liter

Mean/S.E.M.

Orbit >	h equit	p	t	
	11 97/0 246	12.00/0.246	12.13/0.246	12.03/0.142
	12 23/0.246	12.27/0.246	12.30/0.246	12.27/0.142
	12.10/0.174	12.13/0.174	12.22/0.174	

V3 = Plasma Sodium mlq/liter

Orbit \	Type	đ	P	t	
0	168	.33/0.573	168.00/0.573	168.33/0.573	168.22/0.331
	168	33/0.573	168,33/0,573	167.83/0.573	168.17/0.331
1	168	.33/0.405	168.17/0.405	168.08/0.405	

Jacobson 13 Red Cells Phospholipids and Cholesterol

Sample	Baq	Dewar	<u>Orbit=1</u>	<u>Type</u>	<u>V1</u>	<u>2</u>
21	14		1	t	lost	lost
22	17	ь	1	Р	21.84	11.99
23	6	а	1	d	21.20	11.65
24	13	е	0	t	22.01	11.72
25	1	g	0	d	22.25	10.99
26	19	ŕ	0	P	22.46	10.97
27	3	с	1	d	21.43	10.94
28	15	е	0	p	23.42	13.82
29	8	Ь	1	t	23.05	13.16
30	9	g	0	t	22.40	12.65
31	20	a	1	p	21.87	12.38
32	2	f	0	d	24.45	12.83
33	4	е	0	d	23.17	13.32
34	16	q	0	р	26.09	13.32
35	11	Ē	0	t	24.43	11.51
36	10	а	1	t	26.68	12.04
37	21	с	1	P	21.00	12.21
38	5	ь	1	d	21.96	11.02

V1 = Phospholipids $\mu gP/10^8$ cells

Mean/S.E.M.

Orbit \	Type d	P	t	
0	23 29/0 784	23,99/0,784	22.95/0.784	23.41/0.453
	21.53/0.784	21.57/0.784	24.87/0.960	22.38/0.480
	22.41/0.554	22.79/0.554	23.71/0.607	

 $V2 = Cholesterol \mu g / 10^8$ cells

$Orbit \lambda$	Type d	P	t	
	12 28/0 533	12 70/0.533	11.96/0.533	12.35/0.308
	12.30/0.533	12 19/0.533	12,60/0,653	11.92/0.327
I	11.79/0.377	12.45/0.377	12.22/0.413	

Jacobson 14 Plasma Phospholipids and Cholesterol

Sample	Baq	Dewar	Orbit=1	Type	<u></u>	<u>_v2</u>
<u>Jamp r.c</u>	13	e	0	t	1.18	1.48
2	17	Ь	1	р	1.19	1.45
7	19	Ē	0	P	1.27	1.43
4	1	a	0	d	1.22	1.37
5	6	a	1	d	1.24	1.46
6	14	c	1	t	1.20	1.42
7	3	c	1	d	1.20	1.45
8	9	a	0	t	1.20	1.51
9	Ŕ	b	1	t	1.19	1.51
10	2	f	Ō	d	1.25	1.52
10	15	è	Ō	p	1.37	1.34
12	20	a	1	P	1.50	1.36
12	20	- -	0	ď	1.45	1.31
1.4	16	Č	Ő	P	1.33	1.19
14	5	9	1	d	1.27	1.31
15	11	f	0	t	1.36	1.43
10	21		1	ō	1.19	1.23
17	21	ر م	▲ 1	+	1.33	1.31
18	10	a	+	C		

V1 = Phospholipids μ gP/10⁸ cells

Mean/S.E.M.

	vea d	D	t	
UPDITAT		1 32/0 061	1,25/0,061	1.29/0.035
0	1.31/0.061	1.32/0.001	1 24/0 061	1,26/0.035
1	1.24/0.061	1.29/0.081	1.24/0.043	
	1.27/0.043	1.31/0.043	1.24/0.045	

V2 = Cholesterol μ g / 10⁸ cells

Orbit \ Tyr	ne d	P	t	
	1 40/0 056	1,32/0.056	1.47/0.056	1.40/0.033
	1 41/0 056	1.35/0.056	1.41/0.056	1.39/0.033
	1.40/0.040	1.33/0.040	1.44/0.040	

-

Jacobson	15	Plasticizer DEHP				
<u>Sample</u> 4 5 7 10 13 15	<u>Baq</u> 1 6 3 2 4 5	<u>Dewar</u> g a c f e b	<u>Orbit=1</u> 0 1 1 0 0 1	Type d d d d d d	V1 4.65 4.30 4.45 4.25 4.90 4.85	
V1 = Plasticizer mg/100ml plasma Mean (Standard Error of the Mean) Orbit						

•

Samole	Bag	Dewar	<u>Orbit=1</u>	Type	<u></u>	<u>_V2</u>	$\overline{\sqrt{3}}$
500	50	h	0	t	5.8	5.8	5.9
501	60	h	0	p	5.5	5.7	6.0
507	56	Ь	1	t	2.2	3.2	2.7
502	5 0 6 1	Ä	1	p	3.5	2.0	2.9
505	49	d	1	d	10.8	10.1	9.9
504	40	5	1	Ä	10.8	9.8	11.0
505	45	n	0	-	0.5	9 1	94
506	59	n	0	р		101	0 5
507	43	d	1	d	10.3	10.1	9.5
508	54	d	1	t	1.0	1.0	1.5
509	62	đ	1	P	4.5	3.8	3.0
500	47	- -	0	d	7.7	10.0	9.9
510	47		õ	+	7 3	6.7	6.3
511	55	n	0	L L	,	~ ••	

Kenney 1 Surface Proteins - Platelets

 $VI = {}^{125}I - Membrane Proteins$

Mean/S.E.M.

$Orbit \land$	Type d	p	t	
0,0,0,0	9 25/1 419	7.00/1.419	6.55/1.419	7.60/0.819
	10.55/1.419	4.00/1.419	1.60/1.419	5.38/0.819
	9,90/1.004	5.50/1.004	4.07/1.004	

type p = 0.0085

omitting bag number 56: Mean/S.E.M.

$O=bi+ \lambda$	Type d	p	t	
		7 00/1 051	6.55/1.051	7.60/0.607
	9.25/1.051	4.00/1.051	1 00/1.487	6.02/0.665
1	10.55/1.051	4.00/1.031	1.00/1.00	
	9.90/0.743	5.50/0.743	4./0/0.858	

type p = 0.0088

 $V2 = {}^{3}H - Membrane Glycoproteins$

Mean/S.E.M.

$Orbit \wedge$	Type d	p	t	
	9 90/1 263	7.55/1.263	6.25/1.263	7.90/0.729
	10 10/1 263	2 90/1.263	2.10/1.263	5.03/0.729
	10.00/0.893	5.23/0.893	4.17/0.893	

type p = 0.0038orbit p = 0.0240 omitting bag number 56: Mean/S.E.M.

Orbit \ Type d	P	t		
	7 55/0 943	6.25/0.943	7.90/0.544	_
0 9.90/0.945	,	1 00/1 333	5,40/0,596	
1 10.10/0.943	2.90/0.943	1.00/1.335		-
10.00/0.667	5.23/0.667	4.50/0.//0		

type p = 0.0038orbit p = 0.0270

V3 = Cytoskeleton

Mean/S.E.M.

Orbit \	Type d	P	t	
010101	10 45/1 010	7 70/1.010	6.10/1.010	8.08/0.583
	10.43/1.010	0.05/1.010	2 10/1 010	4 92/0.583
1	9.70/1.010	2.95/1.010	2.10/1.010	
	10.08/0.714	5.33/0.714	4.10/0.714	

type p = 0.0008 orbit p = 0.0049

omitting bag number 56: Mean/S.E.M.

$Orbit \lambda$	Type d	P	t	
010101	10 45/0 809	7,70/0,809	6.10/0.809	8.08/0.46/
		2 95/0 809	1.50/1.145	5.36/0.512
	9.7070.803	5 33/0 572	4 57/0.661	
	10.08/0.5/2	5.33/0.372	4.5770.000	

type p = 0.0019 orbit p = 0.0111

Sample	Bag	<u>Dewar</u>	<u>Orbit=1</u>	Type	V1	<u>2</u>	<u></u>
121	19	f	0	P	2.85	7.14	27.14
127	14	Ċ	1	t	2.94	7.35	26.47
122	13	ē	0 0	t	1.40	5.63	22.53
123	17	с Б	1	D	2.63	7.89	28.94
124	1/	0	n	d	2.85	8.57	31.42
125	1 C	9	1	Ā	4.34	7.24	27.53
126	0	a f	0	d	2.94	4.41	35.29
127	2	r -	0	+	A 34	7.24	28.98
128	9	g	0	d d	2 85	7.14	34.28
129	3	C	1		1 11	7 35	32.35
130	8	b	1	ι -	4.41	9.23	36 76
131	15	e	U	p	1.47	0.25	32 35
132	20	а	1	P	1.47	4.41	JZ.JJ
133	10	а	1	t	8.5/	10.00	40.71
134	5	ъ	1	d	6.84	10.95	30.30
135	16	g	0	P	2.89	7.24	36.23
136	4	е	0	d	0.59	5.97	32.83
137	21	С	1	P	1.49	2.98	29.85
138	11	f	0	t	1.21	1.51	28.78

Kevy 1 - Red Cell Osmotic Fragility

V1 = Percent Hemolysis, 0.6% Saline

Mean/S.E.M.

$Orbit \land$	Type d	P	t	
	2 13/1.015	2.40/1.015	2.32/1.015	2.28/0.586
	4 68/1 015	1.86/1.015	5.31/1.015	3.95/0.586
	3.40/0.718	2.13/0.718	3.81/0.718	

V2 = Percent Hemolysis, 0.55% Saline

Mean/S.E.M.

Orbit 1	Type d	P	<u>t</u>	
	6 32/1 224	7,54/1,224	4.79/1.224	<u>6.22/0.70b</u>
	8 44/1 224	5,09/1,224	8.23/1.224	7.26/0.706
	7.38/0.865	6.32/0.865	6.51/0.865	

V3 = Percent Hemolysis, 0.5% Saline

Orbit	Type d	P	t	
0	33, 18/3, 133	33.38/3.133	26.76/3.133	31.11/1.809
	33 39/3,133	30,38/3,133	34.84/3.133	32.87/1.809
	33.28/2.216	31.88/2.216	30.80/2.216	

•

Lionetti 1 - White Cells Recovery and Viability

Sampl	<u>e Bag</u>	Dewar	<u>Orbit=1</u>	Type	<u>V1</u>	<u>_V2</u>	<u>V3</u>
Samp 1 801 802 803 804 805 806 807 808 809 810 811 812 813 814 815 816	e Bag 22 35 31 24 39 40 25 32 29 38 41 23 37 28 33 34	Dewar a b g f c e g e a f b c a e f c	$ \begin{array}{c} 0 \\ 0 \\ 0 \\ 0 \\ 1 \\ 0 \\ 0 \\ 1 \\ 0 \\ 1 \\ 1 \\ 1 \\ 0 \\ 0 \\ 1 \\ 0 \\ 0 \\ 1 \\ 0 \\ 0 \\ 1 \\ 0 \\ 0 \\ 1 \\ 0 \\ 0 \\ 1 \\ 0 \\ 0 \\ 1 \\ 0 \\ 0 \\ 1 \\ 0 \\ 0 \\ 0 \\ 1 \\ 0 \\ 0 \\ 0 \\ 0 \\ 1 \\ 0 \\ 0 \\ 0 \\ 0 \\ 0 \\ 0 \\ 0 \\ 0 \\ 0 \\ 0$	dttdppdttppdpdtt;	9.7 12.1 9.3 8.3 11.8 9.1 13.6 11.0 10.6 11.3 8.0 12.3 10.7 12.4 11.9 11.3	52.9 66.0 50.7 45.3 64.4 49.6 74.2 60.0 57.8 61.6 43.6 67.1 58.4 67.7 64.9 61.6 8	89 83 84 90 93 12 79 85 89 77 88 67 75 88 77 88 77
817 818	36 26	р g	1	p d	9.6	52.4	73
V1 =	count ×	10 ⁶ cell	s/ml	- no	t anal	ysed	

V2 = % recovery

Mean/S.E.M.

Orbit V	Type	Ч	P	t	
	62	40/5 707	60.33/5.707	58.53/5.707	60.42/3.295
	<u> </u>	47/5 707	55 47/5 707	61,80/5,707	58.24/3.295
<u> </u>	<u> </u>	41/5.101	57.00/4.035	60 17/4 035	
	59.	93/4.035	57.90/4.035	80.17/4.000	

V3 = viability % FDA positive

7-b:+ 1	Type d	Þ	t	
	53 67/10 057	86 67/10.057	80.00/10.057	73.44/5.807
	92 33/10 057	82 33/10.057	74.67/10.057	80.11/6.807
	68.50/7.112	84.50/7.112	77.33/7.112	

Sample	Bag	<u>Dewar</u>	<u>Orbit=l</u>	Type	<u>_V1</u>	<u>_V2</u>
801	22	а	1	d	147	141
802	35	ъ	1	t	142	143
803	31	a	0	t	175	195
804	24	f	0	d	171	184
805	39	C	1	P	143	146
806	40	e	0	P	159	178
807	25	a	0	d	147	174
808	32	ē	0	t	157	148
809	29	a	1	t	147	154
810	38	- F	0	р	148	133
811	41	b	1	P	155	178
812	23	- C	1	d	123	130
813	37	a	1	р	141	142
914	28	e	0	d	156	116
014	22	f	0	t	147	154
015	34	ŗ	1	t	163	170
917	36	0	Ō	P	155	162
818	26	b	1	d	149	149
010	<u> </u>					

Lionetti 2 - White Cells Phagocytic Index

V1 = Ingestion of 125 I Staph A, % of N-Ethylmaleimide treated controls, determination 1

Mean/S.E.M.

Orbit \	Type d	P	t	
01010	158,00/6,525	154.00/6.525	159.67/6.525	157.22/3.767
	139 67/6 525	146.33/6.525	150.67/6.525	145.55/3.767
	148.83/4.614	150.17/4.614	155.17/4.614	

orbit p = 0.0490

V2 = Ingestion of 125 I Staph A, % of N-Ethylmaleimide treated controls, determination 2 \cdot

$Orbit \land$	Type	d	P	t	
	158	.00/13.307	157.67/13.307	165.67/13.307	160.44/7.683
	140	.00/13.307	155.33/13.307	155.67/13.307	150.33/7.683
	149	.00/9.408	156.50/9.408	160.67/9.408	

.

Lionetti	3	-	White	Cell	Glucose	Oxidation
	-					

Sample	Bag	Dewar	<u>Orbit=1</u>	<u>Type</u>	<u></u>	<u>_V2</u>
<u>Sample</u> 801 802 803 804 805 806 807 808	Bag 35 31 24 39 40 25 32	<u>Dewar</u> a b g f c e g e	<u>Orbit=1</u> 1 0 0 1 0 0 0	lype dttdppdtt	260 290 320 320 280 410 190 340 290	90 90 110 120 120 120 80 130 130
809 810 811 812 813 814 815 816 817 818	29 38 41 23 37 28 33 34 36 26	a f b c a e f c g b	0 1 1 1 0 0 1 0 1	משסטקרי שס	310 390 200 320 260 360 440 370 280	120 110 60 100 90 120 90 110 90

V1 = ¹⁴C Glucose Oxidation, Stimulated by PMA, 7 over background

Mean/S.E.M.

		n	t	
Orbit \	256.67/32.914	363.33/32.914	340.00/32.914	320.00/19.003
	246.67/32.914 251.67/23.274	<u>330.00/32.914</u> 346.67/23.274	340.00/32.914 340.00/23.274	

type p = 0.0237

V2 = ¹⁴C Glucose Oxidation, Stimulated by FMLP, 7 over background

Mean/S.E.M.

o_{-}	P	t	
<u>OFBIC Type</u> 201	116.67/7.201	120.00/7.201	111.11/4.15/
	103.33/7.201	93.33/7.201	92.22/4.15/
88.33/5.092	110.00/5.092	106.67/5.092	

type p = 0.0230orbit p = 0.0075
Sample	Bag	Dewar	<u>Orbit=1</u>	Type		_ <u>V2</u>	<u>V3</u>	
201 202 203 204 205 206 207 208 209 210 211 212 213 214 215 216 217 218	19 6 17 14 1 13 2 3 9 20 8 15 4 11 21 5 16 10	f a b c g e f c g a b e e f c b g a	0 1 1 1 0 0 0 1 1 1 0 0 1 1 1 0 0 1 1 1 0 0	הם הים הים חים הים הי	14.08 14.20 14.50 14.18 14.48 13.95 14.20 12.28 14.33 14.70 13.83 14.70 16.40 16.25 16.63 16.30 16.60	106.8 72.5 73.8 69.0 73.0 65.2 56.3 74.4 64.6 58.1 63.0 53.9 60.5 71.8 75.2 69.4 62.2 76.3	8.35 8.58 8.58 7.95 6.73 6.73 7.28 7.13 7.85 8.63 9.13 9.25 8.85 9.15 8.68	115.75 116.25 123.50 103.67 110.33 99.50 98.67 96.25 94.50 99.75 99.00 99.25 94.75 94.25 94.25 93.25 96.50 91.00

Szymanski 1 - Red Cells Agglutination

V1 = IgG on intact cells, molecules / cell

Mean/S.E.M.

$\begin{array}{c ccccccccccccccccccccccccccccccccccc$		Type d	P	t	11 70 /0 457
$\begin{array}{c ccccccccccccccccccccccccccccccccccc$	UFDIC 1	14 53/0 791	14,92/0,791	14.89/0.791	14.7870.457
14.37/0.731 .00100		14.35/0.791	15,15/0.791	14.87/0.791	14.80/0.45/
14 45/0 559 15.04/0.559 14.88/0.555			15.04/0.559	14.88/0.559	

V2 = IgG inside cells, molecules/cell

Mean/S.E.M.

	P	t	60 06 / 1 202
$\frac{0}{63} \frac{27}{7} \frac{607}{607}$	74.30/7.607	67.20/7.607	68.26/4.392
$-\frac{0}{72}$ $\frac{0}{10/7}$ $\frac{0}{607}$	69.03/7.607	69.43/7.607	70.1974.392
67.68/5.379	71.67/5.379	68.32/5.379	

V3 = C3c on intact cells, molecules/cell

Mean/S.E.M.

Orbit \ Type d	P	t	0.01/0.309
0 7,55/0,535	8.71/0.535	7.77/0.535	8.01/0.30
1 8.23/0.535	8.56/0.535	8.52/0.535	0.44/0.50
7.89/0.378	8.64/0.378	8.15/0.3/8	

V4 = C3d on intact cells, molecules/cell

Mean/S.E.M.

	T - d	D	t	
<u>Orbit \</u>	Type u	102 75/6 568	100 44/6.568	101.67/3.792
0	100.81/6.568	103.7576.560	104 75/6 568	103.05/3.792
1	102.72/6.568	101.67/6.568	104.75/8.500	100,00976
	101 77/4.644	102.71/4.644	102.60/4.644	
	1011, 1, 1, 1, 1, 1, 1, 1, 1, 1, 1, 1, 1,			
Mean	(Standard Error	of the Mean)		
Type d	: 101.77 (3.20)	p: 102.71 (4.26	5) t:102.60 ((4.97)

Orbit 0:101.67 (2.35) 1:103.05 (4.06)

.

.

<u>Sample</u>	Bag	Dewar	<u>Orbit=1</u>	T∨pe	<u></u>	_V2	<u>_V3</u>	
201	19	f	0	P	21.90	48.85	76.30	54.85
202	6	a	1	d	21.25	47.10	73.10	50.45
202	17	Ь	1	Þ	18.25	44.30	73.30	50.55
203	14	- C	1	t	16.30	40.40	72.95	50.45
204	1	ā	0	d	22.25	57.40	77.25	56.35
205	13	9 e	Õ	t	21.35	48.60	76.00	55.70
200	.2	f	Ō	d	21.60	46.95	74.60	54.70
207	2	, C	1	d	22.20	48.30	75.35	53.60
208	2 Q	9	0	t	14.35	39.00	73.00	49.65
209	20	9	1	p	16.80	42.20	72.80	49.60
210	20	5	1	t	27.85	54.90	78.00	56.15
211	15		î.	Ď	22.10	50.15	77.40	56.80
212	10	e	0	d	26.55	54.85	77.80	56.95
213	11	e F	0	+	23.85	50.75	77.25	56.80
214	21	1	1	n	26.50	52.55	76.95	56.00
215	21		1	ч Б	24 60	51.30	76.55	56.90
216	5	D	1	u 5	26 30	51 65	76.90	56.60
217	16	g	U	ч +	20.30	48 15	76.20	54.40
218	10	а	1	Ľ	20.90	40.15	,	

Szymanski 2 - Red Cells Agglutination

V1 = 7. Agglutination with anti-C3a

Mean/S.E.M.

Orbit \	Type d	P	t	
010101	23 47/2 378	23,43/2,378	19.85/2.378	22.25/1.373
	22 68/2 378	20 52/2 378	21.68/2.378	21.63/1.373
	23.08/1.682	21.98/1.682	20.77/1.682	

V2 = % Agglutination with anti-C3c

Mean/S.E.M.

Orbit \ Type d	P	t	
<u>57 67 C 7995</u> 0 52 07/2 963	50 22/2 963	46.12/2.963	49.80/1.711
1 48 90/2 963	46 35/2 963	47.82/2.963	47.69/1.711
50,98/2,095	48.28/2.095	46.97/2.095	

•

 $\sqrt{3} = \%$ Agglutination with anti-C3d

Mean/S.E.M.

Drhit \	Type d	P	t	
<u> </u>	76 55/1 124	76.87/1.124	75.42/1.124	76.28/0.649
	75 00/1 124	74 35/1,124	75.72/1.124	75.02/0.649
	75.78/0.795	75.61/0.795	75.57/0.795	

V4 = 7 Agglutination with anti-IgG

•

Mean/S.E.M.

Orbit \	Type d	p	t	
0	56 00/1.63	6 56,08/1.636	54.05/1.636	55.38/0.945
1	53 65/1.63	6 52.05/1.636	53.67/1.646	53.12/0.945
ł	54.83/1.15	7 54.07/1.157	53.86/1.157	

.

Szymanski 3 – Platelets

Sample	Bag	Dewar	<u>Orbit=1</u>	Туре	<u></u> V1	_ <u>V2</u> _	<u>_V3</u>
600	60	h	0	P	4018	2085	2049
601	45	h	0	d	6611	2198	2123
602	50	h	0	t	1597	2041	2002
603	61	d	1	p	2840	1574	1840
604	48	d	1	d	6374	2181	2750
605	56	d	1	t	2997	1218	1200
606	54	d	1	t	1844	1148	1207
607	62	d	1	р	2488	1478	1934
608	43	d	1	d	5855	2312	2783
600	47	h	0	d	5425	2302	2560
610	59	h	0	р	5059	2164	2461
611	55	h	0	t	5902	2132	2400
612	46	h	Ō	d	5656	2030	2324
613	53	Ь	1	t	2147	1561	1719
614	63	h	Ō	p	4208	2013	2218
615	51	h	Õ	t	4355	1775	2412
616	58	Ь	1	p	5433	1715	2273
617	44	ď	1	d	6586	2130	2775
017		-					

V1 = Molecules of bound IgG/platelet

Mean/S.E.M.

$0-bi+ \lambda$	h equi	P	t	
	5997/688 9	4428/688.9	3951/688.9	4759/397.8
	6771/688 9	3587/688.9	2329/688.9	4062/397.8
	6084/487.2	4007/487.2	3140/487.2	

type p = 0.0032

omitting bag number 56:

Mean/S.E.M.

3	Tupe d	a	t	
<u>JEDIC A</u>	<u>5897/705 4</u>	4428/705.4	3951/705.4	4759/407.2
	<u> </u>	3587/705.4	1995/863.9	4195/431.9
<u>I</u>	6084/498.8	4007/498.8	3169/546.4	

type p = 0.0052

V2 = Number of C3d molecules/platelet

Mean/S.E.M.

Orbit \	T∨pe d	P	t	
0	2176/85,26	2087/85.26	1982/85.26	2082/49.22
1	2207/85.26	1589/85.26	1309/85.26	1701/49.22
	2192/60.29	1838/60.29	1645/60.29	

type p = 0.0001orbit p = 0.0001interaction p = 0.0037

omitting bag number 56: Mean/S.E.M.

Orbit \	Type d	P	t	
0	2176/86,91	2087/86.91	1982/86.91	2082/50.18
1	2207/86.91	1589/86.91	1354/106.4	1762/53.22
	2192/61.45	1838/61.45	1731/67.32	

type p = 0.0005orbit p = 0.0011interaction p = 0.0085

v3 = Number of C3c molecules/platelet

Mean/S.E.M.

Orbit \	Type d	P	t	
0	2335/126.0	2242/126.0	2271/126.0	2283/72.77
	2769/126.0	2015/126.0	1375/126.0	2053/72.77
	2552/89.12	2129/89.12	1823/89.12	

type p = 0.0003orbit p = 0.0454interaction p = 0.0007

omitting bag number 56: Mean/S.E.M.

Orbit \	Type d	P	t	
0	2335/126.2	2242/126.2	2271/126.2	2293/72.88
1	2769/126.2	2015/126.2	1463/154.6	2160/77.30
······	2552/89.25	2129/89.25	1948/97.77	

type p = 0.0019interaction p = 0.0025

-

Szymanski – Platelets

	D	Davar	0 = bi = 1	Type	V 1	V2
Sample	Bag	Dewar	$\frac{01010-1}{0}$		$37\frac{1}{122}$	33104
600	60	n F	0	ч Б	52517	45906
601	45	n L	0	+	30479	29382
602	50	n,	U ,	- L	20264	25424
603	61	d	1	p	20204	20424
604	48	d	1	d	4/284	40910
605	56	d	1	t	16298	13301
606	54	d	1	t	16234	14390
607	62	d	1	р	26533	24045
608	43	d	1	d	54000	48145
609	47	ĥ	0	d	47809	42384
610	59	h	0	P	37094	32035
611	55	h	0	t	33935	28033
612	46	h	0	d	44638	38982
613	53	d	1	t	24738	22591
614	63	h	0	P	36985	32777
C15	51	 Б	n	t	39813	35458
612	51		1	5	42508	36978
616	29	a	1	4	559/1	19255
617	44	d	1	a	55641	49233

V1 Total Platelet IgG (molecules/platelet)

Mean/S.E.M.

0		Þ	t	
OFBIL	49221/2972	37067/2972	34742/2972	40043/1716
	<u> </u>	32435/2972	19090/2972	34633/1716
1	50348/2101	34751/2101	26916/2101	

type p = 0.0001orbit p = 0.0457interaction p = 0.0200

omitting bag number 56: Mean/S.E.M.

Orbit \	Type d	q	t	
<u> </u>	49321/3047	37067/3047	34742/3047	40043/1759
	52375/3047	32435/3047	18490/3732	36925/1865
1	50348/2154	34751/2154	29039/2360	

type p = 0.0001 interaction p = 0.0342 V2 Total free IgG in Frozen - Thawed platelets (molecules/platelet)

Mean/S.E.M.

		D	t	
UFBIT		22629/2632	30957/2632	35340/1520
0	42424/2632		16761/2632	30560/1520
1	46103/2632	28816/2632	22850/1861	
	44264/1861	30727/1861	23859/1001	

type p = 0.0001orbit p = 0.0461interaction p = 0.0171

omitting bag number 56: Mean/S.E.M.

	h d	P	t	
UPDIC VIVE	42424/2649	32629/2648	30957/2648	35340/1529
	4242472640	28816/2648	18491/3243	32717/1622
	46103/2646	20727/1973	25970/2051	
	44264/18/3	3072771075	2557072001	

type p = 0.0001interaction p = 0.0327

meer deeron

Page 53

Taylor White Cells

Sample	Bag	<u>Dewar</u>	<u>Orbit=1</u>	<u>Type</u>	<u>V1</u>	<u>v2</u>	<u></u>	<u></u>	<u>_V5</u> _	<u></u> V6	<u>7</u>	<u>8</u>
901 904 910 912 913 918	22 40 29 32 37 28	a e a e	1 0 1 0 1 0	d ptt pd	94 98 91 94 98 95	95 95 90 97 97 98	38.2 38.4 52.9 15.6 17.1 19.5	41.5 30.0 75.9 23.4 21.0 26.8	16.9 36.3 47.1 26.0 22.8 60.4	27.0 79.8 98.7 80.3 115.5 79.7	38.6 55.1 88.4 81.0 94.9 74.8	0.78 0.72 2.06 1.15 4.98

V1 = % viability

Mean/S.E.M.

0	h equ	Ð	t	
UPDIC 3	<u>1996</u> 05 00/2 997	98 00/2 887	94.00/2.887	95.67/1.667
<u> </u>	95.00/2.887	99 00/2 887	91,00/2,887	94.33/1.667
	94.00/2.887	98.00/2.00/	02 50/2 041	
	94.50/2.041	98.00/2.041	92.50/2.041	

V2 = % Viability - 72 hours

Mean/S.E.M.

$2-bi+ \sqrt{2}voe$ d	a	t	
<u>99 00/2 261</u>	95,00/2,261	97.00/2.261	96.67/1.599
0 98.00/2.201	97 00/2 261	90.00/2.261	94.33/1.599
96.50/1.985	96.00/1.985	93.50/1.985	

V3 = Protein Synthesis w/o Mitogen

Mean/S.E.M.

3-51+ N	Tupe d	Þ	t	
JEDIC J	19 <u>0</u> 0	29 40/15 37	15,60/15.37	24.50/8.874
	19.50/15.37	30.40/13.37	E2 90/15 37	36,07/8,874
1	38.20/15.37	17.10/15.37	52.30/13.37	
	28.85/10.87	27.75/10.87	34.25/10.8/	

/4 = Protein Synthesis in response to phytohemagglutinin

1ean/S.E.M.

3-5 i + 3	Type d	Ð	t	
<u>JIDIC (</u>	26 90/19 75	30 00/19.75	23.40/19.75	26.73/11.41
	26.80/19.75	21.00/10.75	75 90/19 75	46,13/11,41
1	41.50/19.75	21.00/19.75	10.55/13.07	
	34.15/13.97	25.50/13.9/	49.65/13.9/	

Continued on next page)

Page 54

/5 = Blastogenesis w/o Mitogen

1ean/S.E.M.

Indit 1	Type d	P	t	
<u>,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,</u>	60.40/16.85	36.30/16.85	26.00/16.85	40.90/9.729
	16,90/16,85	22.80/16.85	47.10/16.85	28.93/9.729
	38.65/11.92	29.55/11.92	53.86/11.92	

/6 = Blastogenesis in response to PHA

fean/S.E.M.

hebit \	Type d	P	t	
<u>, , , , , , , , , , , , , , , , , , , </u>	79 70/34 23	79.80/34.23	80.30/34.23	79.93/19.19
	27 00/34.23	115.5/34.23	98.70/34.23	80.40/19.19
	53.35/23.50	97.65/23.50	89.50/23.50	

/7 = Blastogenesis in response to Pokeweed mitogen

lean/S.E.M.

Indit \ Type d	P	t	
0 74 80/23 78	55,10/23,78	81.00/23.78	70.30/13.73
	94,90/23,78	88.40/23.78	73.97/13.73
56.70/16.82	75.00/16.82	84.70/16.82	

'8 = Flow cytometry - phycoerythrin controls - % autoflourescence

lean/S.E.M.

Indit \	Type d	Ρ	t	
<u></u>	4.98/1.766	0.78/1.766	2.06/1.766	2.61/1.020
		1.15/1.766	0.72/1.766	0.94/1.249
	4.98/1.766	0.97/1.249	1.39/1.249	



Addendum to Statistical Report

It was suggested by Dr. Surgenor that certain of the measured variables fell into natural clusters, and that it might be helpful to attempt to analyse the influence of the experimental factors on these groups as a whole. The clusters suggested by Dr. Surgenor were analysed using multivariate analysis of variance methods. The SAS procedure GLM was again used, utilizing the MANOVA option.

The first group analysed was a white cell group consisting of the following variables:

Variable	<u>Table in original report</u>
pH	Jacobson 1
pCO_	"
p0 ₂	"
viability, % FDA positive	Lionetti l
Phagocytic index 1	Lionetti 2
Phagocytic index 2	"
Glucose oxidation, PMA	Lionetti 3
Glucose oxidation, FMLP	"

For this cluster, the orbit factor was not significant at the 5% level (p = .09), the bag type factor was significant (p = .0002) and there was no significant interaction effect.

The second group was a red cell metabolism group:

Variable	Table in original report
рН рСО ₂	Jacobson 2
ATP 2,3DPG Glucose	Jacobson 5 "

For this cluster, the orbit factor was not significant, the bag type factor was significant (p = .002) and there was no significant interaction effect.

Page 2 Blevins 5/11/87

The third cluster was a red cell structure group:

Variable

Table in original report

Cell Potassium	Jacobson 5 "
Cell Sodium Plasma Potassium Plasma Sodium	Jacobson 6
Plasma hemoglobin Osmotic fragility	**
.5%, .55%, and .6% saline Phospholipids	Kevy 1 Jacobson 13 "
Choresceror	

None of the factors were significant for this cluster.

The fourth cluster was a red cell membrane binding group:

<u>Variable</u>

Table in original report

tec on intact cells	Szymanski	1
	**	
IgG inside cells	**	
C3c on intact cells	**	
C3c inside cells	Savmaneki	2
Agglutination with anti-C3a	SZYINALISKI	-
Agglutination with anti-C3C		
Agglutination with anti-C3d		
Agglutination with anti-IgG		

None of the factors were significant for this cluster.

Page 3 Blevins 5/12/87

The fifth cluster was a platelet metabolism group:

_ ... ____

Variable	Table in original report
Glucose pH pCO ₂	Chao 13 Jacobson 3 "
PO ₂	**
Lactate Platelet count Mean platelet volume	Chao 12 Chao 1 Chao 3
The orbit factor was signific (p = .0002).	ant (p = .03) as was bag type
The sixth cluster was a p	platelet structure group:
	This is estained report

<u>Variable</u>

Table in original report

* Membrane protein	Kenney l
* Membrane glycoproteins	89
* Cytoskeleton	
Bound IgG	Szymanski s
Bound C3d	17
Bound C3c	
Average transmission	Ausprunk
electron microscopy score	Auspi unk

* dropped from the analysis due to missing data.

.

The orbit factor was significant (p = .033) as was the bag type (p = .048).

C.

Page 4 Blevins 5/12/87

The final cluster was a platelet function group:

Variable

Table in original report

The above together with: Aggregation by collagen Aggregation by ADP ATP release Serotonin uptake β - Thromboglobulin Thromboxane B₂

Chao 4 Chao 5 Chao 6 Chao 8 Chao 10 Chao 11

None of the factors were significant for this group.

EPISTEMOS, INC. 28 OLD MILL ROAD QUAKER HILL. CT. 06375 203 442-7764

RECEIVED

JUN 2 9 1987

PRESIDENT'S OFFICE

June 27,1987

Dr. Douglas Surgenor The Center for Blood Research 800 Huntington Avenue Boston, MA 02115

Dear Doug,

Here is a report on the analysis of the data on position of storage of the samples. As you can see, some of the signicance of the orbit or bag type effects change, and there is some effect of position on many but not most of the measured outcomes. I don't know if these effects are consistent or random, but they definitely exist. Taking them into account in the analysis lends somewhat more weight to our conclusions, and it seems obvious that it would be well worth some effort in any future experiments to attempt to dampen these effects, either by designing the storage system to some way reduce the effects and/or by carefully matching on storage

I am going to be out of the country during the entire month of July. I hope Nan can clear up any questions you may have until then, and I will be in touch with you when I return.

Sincerely,

Donald K. Blevins, Ph.D.

PRECEDING PAGE BLANK NOT FILMED

Statistical Report 3 TBSE

We will discuss in this report, an additional analysis which we have performed on the data from the Initial Blood Storage Experiment. This analysis was an attempt to evaluate the effect if any of the physical location of storage of the blood bags in the storage module.

There were identical storage modules for the orbiter and on Earth. In each of these, the red and white cell bags were stored in three different Dewars, labeled a, b, and c in for the ground unit and correspondingly labeled e, f, and g in the flight unit. Within each Dewar, the bags were stored in one of three positions - front, middle or back (see diagram 1). The platelet bags were all stored in one Dewar, labeled d and h respectively, in each storage unit. Five platelet bags were stored in front and five in back. We therefor introduced for red and white cells two new factors, Dewar and position, and for platelets one new factor, position (front or back). We included these factors together with the previous factors of bag type and orbital status in an analysis of variance as before. Due to the larger number of experimental factors, there was now insufficient data to include any interaction effects, so we used only a main effects model. In the majority of cases, the inclusion of these factors caused no significant change in the conclusions of the previous analyses, but as will be described in the following, in a number of cases one or both of the following effects were noted. First, for some of the experimental variables, the Dewar and/or position in which the sample was stored had a statistically significant effect on the outcome. Second, the inclusion these new factors in the model in some cases changed the significance levels of the original factors. This can occur if the changes in the experimental variable due to the new factors was large enough to partially mask the change caused by the others.

Following is a summary of the new results. Any variables not noted were not significantly changed in this analyses as compared to the previous analysis. (page numbers refer to the previously distributed complete report)

Ausprunk - Transmission Electron Microscopy Data - Red Cells V1 Total Number of Cells Counted - orbit p = 0.0459 Ausprunk - Transmission Electron Microscopy Data - White Cells (page 6) V5 Average Ranking - orbit p = 0.0565 Chao 4 (page 12) V1 Aggregation (by collagen) - orbit p = 0.0202 position p = 0.0195 Chao 5 (page 13) V1 Aggregation (by ADP) - type p = 0.0634 position p = 0.0206

> ORIGINAL PAGE IS OF POOR QUALITY

IBSE Report 3

```
Chao 6 (page 14)
          V1 ATP Release - orbit p = 0.0692
                            position p = 0.0464
Chao 8 (page 16)
          V1 Serotonin Uptake - position p = 0.0506
     13 (page 21)
Chao
          V1 Glucose - position p = 0.0344
         - White Cell Distribution
Curby
      2
          V1 - Dewar p = 0.0032
V2 - Dewar p = 0.0006
         V2
          V3 Peak Channel - type p = 0.0049
                            Dewar p = 0.0049
         1 White Blood Cells (page 26)
Jacobson
          V1 \text{ pH} - type p = 0.0001
                    position p = 0.0001
                    Dewar p = 0.0577
          V2 pCO_2 - orbit p = 0.0001
                   type p = 0.0001
                     position p = 0.0016
                     Dewar p = 0.0044
Jacobson 2 - Red Blood Cells (page 28)
         V1 pH - type p = 0.0001
                     orbit p = 0.0005
                     position p = 0.0004
                     Dewar p = 0.0428
          V2 pCO_{2} - orbit p = 0.0060
                     type p = 0.0001
Jacobson
         3 - Platelets (page 30)
          V1 pH - type p = 0.0589
                     position p = 0.0277
Jacobson 5 - Red Blood Cells (page 33)
          V3 Glucose - position p = 0.0359
Jacobson 6 - Red Blood Cells (page 35)
          V2 Plasma Potassium - orbit p = 0.0006
                                  position p = 0.0001
          V3 Plasma Sodium - position p = 0.0548
Lionetti 1 - White Cells Recovery and Viability (page 42)
          V1 Count - position p = 0.0010
                       Dewar p = 0.0250
          V2 % recovery - position p = 0.0010
                            Dewar p = 0.0249
```

```
page 2
```

~

Lionetti 2 - White Cells Phagocytic Index (page 43) orbit p = 0.0104V1 position p = 0.0079position p = 0.0090V2 Lionetti 3 - White Cells Glucose Oxidation (page 44) ¹⁴C Glucose Oxidation, Stimulated by PMA -V1 type p = 0.0068position p = 0.0370v2 ¹⁴C Glucose Oxidation, Stimulated by FMLP type p = 0.0018orbit p = 0.0005 position p = 0.0764Dewar p = 0.0114Szymanski 1 - Red Cells Agglutination (page 45) V1 IgG on intact cells - position p = 0.0099 Szymanski 2 - Red Cells Agglutination (page 47) V1 % agglutination with anti-C3a - position = 0.0171 V2 % agglutination with anti-C3c - position = 0.0595

orbit							
	red	white	red	whit	e red	wh	ite
A	6	22	10	29	20	37	
B	8	41	17	26	5	35	
С	21	32	3	31	14	2]	
D		56 62 43 54 44			48 53 58 19 51		

earth							
r	rd	white	red	white	ved	nh.	Ŧe
F	4	40	13	32	15	z <i>8</i>	
F	n	24	19	38	2	33	
G	16	3]	1	25	9	36	
17		60 47 51 59			50 98 63		
	•	~ I 45	.	4	55		

Diagram 1

EPISTEMOS, INC.

28 OLD MILL ROAD QUAKER HILL, CT. 06375

203 442-7764 203 444-073]

September 8, 1987

Dr. Douglas Surgenor The Center for Blood Research 800 Huntington Avenue Boston, MA 02115

Dear Doug;

Here is the analysis of the TOTM platelet bags that Nan said you were interested in. I think it is self explanatory, but as usual feel free to cll me if you have questions.

Sincerely Donald K. Blevins, Ph.D.



SEP 1 0 1987

Statistical Report 4 IBSE

We have done an analysis of the effect of microgravilty versus controls considering only the platelet bags of type PVC TOTM. Since the only variable factor was the two level earth/orbit factor, a simple t-test suffices for the analysis. The SAS TTEST procedure was used. A basic assumption of the ttest is that the underlying population of each of the two groups are normally distributed with equal variances. If the equal variance assumption is not satisfied an approximate t-test can be used. The SAS procedure TTEST reports both results together with a test statistic for the equal variance assumption. For these data, the assumption of equal variance was statistically acceptable (p > .10) except for a few of the outcome variables. Therefor in the following summary, we will report only the results of the standard t-test except for those few items, for which we will report both the standard statistic and the approximate statistic. Unless otherwise noted the degree of freedom for the t-test is 5. The degree of freedom for the approximate t-test depends on the calculated sample variances.

Results

Ausprunk	data:						
	Psuedopods	Τ=	-2.39,	P =	0.06		
	Aggregated	T =	0.85,	p =	0.44		
	Degranulated	T =	1.20,	p =	0.29		
	Swollen/Ruptured	T =	2.07,	p =	0.09		
	Average Score	T =	1.66,	p =	0.16		
Chao data	3:						
	Platelet count	T =	5.16,	P =	0.0036		
	Hypotonic stress	T =	-1.41,	P =	0.22		
	Mean volume	T =	9.07,	P =	0.0003		
	Aggregation (collagen)						
		T =	-2.59,	P =	0.0489		
	Aggregation (ADP)	T =	-1.94.	p =	0.11		
	(Unequal var. p =	0.0	626, ap	prox	T = -1.66	$\mathbf{p} = 0$.23)
	ATP release	Τ=	-1.38,	P ≃	0.23		
	Russell's Viper venom time						
		T ≖	1.20,	P =	0.29		
	Serotonin uptake	T =	-1.46,	P =	0.20		-
	(Unequal var. p =	0.0	105, ap	prox	T = -1.23	, p = 0	.34)

```
Serotonin release T = -1.20, p = 0.29
                             T = 1.40, p = 0.22
         8-Thromboglobin
                              T = 1.88, p = 0.12
         Thromboxane B<sub>2</sub>
           (Unequal var. p = 0.0025, approx. T = 2.22, p = 0.11)
                              T = 2.07, p = 0.09
         Lactate
                              T = -0.83, p = 0.45
         Glucose
           (Unequal var. p = 0.0318, approx. T = -0.70, p = 0.55)
Curby data:
         1-254/m1 \times 10^9
                              T = 0.80, p = 0.46
         > 254/ml x 10<sup>9</sup>
                              T = -0.19, p = 0.86
                              T = -0.33, p = 0.75
         Peak channel
                              T = -1.10, p = 0.32
         Delay factor
Jacobson data:
                              T = -1.21, p = 0.28
         pH
           (Unequal var. p = 0.0165, apprx. T = -1.02, p = 0.41)
                              T = -2.28, p = 0.07
         pC0,
           (Unequal var. p = 0.0001, approx. T = -1.91, p = 0.20)
                              T = 2.50, p = 0.0543
         P<sup>0</sup>2
Kenney data:
         1251 - membrane protein
                              T = 4.27 df = 1, p = 0.14
         <sup>3</sup>H - membrane glycoprotein
                              T = 6.74 df = 1, p = 0.09
                              T = 13.28 df = 1, p = 0.0479
         Cytoskeleton
Szymanski data:
                              T = 1.20 df = 3, p = 0.32
         Bound IgG
         C3d molecules / platelet
                                   3.03 df = 3, p = 0.0560
                              T =
         C3c molecules / platelet
                              T = 3.13 df = 3, p = 0.0520
                              T = 3.01 df = 3, p = 0.0572
         Total IgG
         Total free IgG in frozen - thawed platelets
                              T = 2.94 df = 3, p = 0.06
```