Light Microscopy Module Light Microscopy Module International Space Station Premier Automated Microscope





SPACE





AND

PHYSICAL SCIENCES










Fluorescently tagged root of Arabidopsis thaliana (thale cress) from Characterizing Arabidopsis Root Attractions (CARA, ISS SpaceX-3), provided by Anna-Lisa Paul and Rob Ferl, University of Florida

- Fluorescently labeled actin (red) and focal adhesion (green in rabbit synovial fibroblasts, provided by Eduardo Almeida, NASA Ames Research Center.
- 3. ACE-T1 studies the self-assembly of Janus particles and how they can be used to spontaneously form precisely organized 3D structures. This work has the promise of ultimately providing efficient and affordable manufacturing processes for functional devices and materials with novel or enhanced properties.
- from Lawrence DeLucas, Year 2 Progress Report, NASA grant
- osite of CVB data, Cover-Transport Phenomena Fundamentals, 3rd Edition, Joe Plawsky, NASA grant number NNX09AL98G
- 6. Supercubes, Cover-Soft Matter, July 2011, Rossi, Sacanna, Irvine,
- Three-dimensional model of colloidal gelation, Peter Lu, Harvard, Peter Lu, org.

Characteristics and Features



1 IMPERX GEV-B2020 cameras

- 7 LMM control box
- 8 Equipment transfer module



For more information about the Light Microscopy Module please visit ://spaceflightsystems.grc.nasa.gov/SOPO/ICHO/IRP/FCF/Investigations/LM

Points of Contact

Ron Sicker—LMM Payloads—NASA Glenn Research Center 21000 Brookpark Road MS 77–7 | Cleveland, OH 44135 216–433–6498 Ronald.J.Sicker@nasa.gov Kelly Bailey—7IN Technologies. Inc.

Kelly Bailey—ZIN Technologies, Inc. 745 Engle Road | Cleveland, OH 44130 619-322-5686 baileyk@zin-tech.com

illiam Foster—LMM Facility—NASA Glenn Research Center 21000 Brookpark Road MS 77–7 | Cleveland, OH 44135 216–433–2368 William.M.Foster@nasa.gov

Andrew Suttles—NASA Glenn Research Center 000 Brookpark Road MS 77-7 ↓ Cleveland, OH 44135 216-433-8328 Andrew.C.Suttles@nasa.gov

Francic D. Chierana I. D. rancis P. Chiaramonte, Ph.D. ace Life and Physical Science NASA Headquarters 300 E Street SW Vashington, DC 20546–0001 202–358–0693 nncis.P.Chiaramonte@nasa.g

An Advanced Automated Microscope

Innovative Microscopy Research Capability from Space

Innovation and imagination are all that are required to use the Light Microscopy Module (LMM) as a laboratory microscope to perform research aboard the International Space Station (ISS). The LMM is a remotely controllable, automated microscope that gives scientists the ability to study—in real time— the effects of the space environment on physics and biology. Specimens can be studied without the need to return the samples to Earth.

Microscope Modified for Space Research

The LMM flight unit features a modified commercial laboratory Leica RXA microscope configured to operate in an automated mode with interaction from the ground support staff. Its core capabilities include a level of containment, white light imaging (available now), fluorescence, confocal microscopy (available in 2016 to 2017), and an imaging capability from a Q-Imaging Retiga 1300 camera.

LMM Supported in the Fluids Integrated Rack

The LMM operates in the Fluids Integrated Rack (FIR), which is located in the U.S. Destiny Laboratory of the ISS. The FIR provides the LMM with the laboratory infrastructure common to most investigations, including an optics bench, temperature control, power control, illumination, imaging and frame capture, data processing, and other resources. The FIR also provides isolation from vibrations on the station to allow for a more stable environment to obtain high-resolution images. The LMM, in conjunction with the FIR, will help fulfill the vision of a true laboratory in space, which is ideal for low-cost payload development

Critical Research Enabled by LMM

How matter is organized and moves on the microscopic level profoundly affects the macroscopic world. Understanding these processes will help scientists and engineers build more efficient machines and consumer products both on Earth and space applications. A suite of experiments is enabled by the LMM to allow for a detailed characterization of fluids colloids two-phase media and biological samples. In the future, the LMM could be used to assist in maintenance of station crew health to advance knowledge of the effects of space on biology, and to contribute to long-term mission space exploration.

Modified Microscope

Engineers at NASA Glenn Research Center modified a Leica RXA laboratory-grade microscope by adding 23 micromotors to permit remote control by scientists on the ground and to meet the demands of space flight and crew-tended operations. As such, it contains all of the necessary optical components for use as a fully functional microscope. The microscope can house many different lenses corresponding to magnifications of 2.5 \times , 4 \times , 10 \times , 20 \times , 40 \times , 50 \times , 63 \times (air), 63 \times and 100 \times oil-coupled objectives. Present capabilities include brightfield and epi-illumination microscopy. Future planned capabilities include high-resolution color video microscopy, condenser assembly, confocal microscopy, and possibly laser tweezers.






Sample Modules

Biological samples for the LMM launched on the Space Shuttle Discovery's STS-133 mission on February 24, 2011, included fixed slides containing yeast, bacteria, a leaf, a fly (Drosophila), a butterfly wing, tissue sections and blood, and six containers of live C. elegans worms. The wing was from "Butterflies in Space" (a previous study that involved students from around the country) that was flown into space in 2009 on STS-129. In addition, some of the worms were descendants of those that survived the Space Shuttle Columbia STS-107 accident. These experiments were operated using OptiCells™ and sample slides.

Physical Science Sample Cells

On-going experiments in thermophysics are conducted in the Constrained Vapor Bubble (CVB) sample modules. In development are "specialty cells" that use physical properties (temperature, E-field) to implement specific colloidal science test matrices. The colloidal science specialty cells will have a geometry and complexity similar to that of the CVB sample module. Prototype temperature, electric field, and biology specialty cells are in





Pre-Advanced Colloids Experiment (PACE) LED Base

Present operations use a sample module mounted on a PACE light-emitting diode (LED) base. The ACE-M sample module has 15 small wells (1.2 µl) to observe colloid experiments in epi-illumination. The current experiment with Procter & Gamble will test product stability. A planned experiment with the University of Pennsylvania will test temperature sensitive particles. The PACE module will be modified to include a heater for this experiment. ACE-T uses the same mechanical interface as the PACE LED base. When ACE-T launched in 2016 it became the new sample module. The petri base design can be used with Velcro® attachment for plant biology and macromolecular biophysics experiments.





Characteristics and Features

Engineers at NASA Glenn Research Center modified a Leica RXA laboratory-grade microscope by adding 23 micromotors to permit remote control by scientists on the ground and to meet the demands of space flight and crew-tended operations. As such, it contains all of the necessary optical components for use as a fully functional microscope. The microscope can house many different lenses corresponding to magnifications of $2.5 \times .4 \times .10 \times .20 \times .40 \times .50 \times .63 \times$ (air), $63 \times$ and $100 \times$ oil-coupled objectives. Present capabilities include brightfield and epi-illumination microscopy. Future planned capabilities include high-resolution color video microscopy, condenser assembly, confocal microscopy, and possibly laser tweezers

LMM Control Box (LCB)

The electrical design of the LMM uses parts of the existing electronics of the Leica microscope and supplementary internal and external electronics that support enhanced automation and imaging capabilities. Motors and linear actuators have been added to motorize the manual functions of the Leica microscope. The LCB provides 16 axes of control for stepper motors and 4 axes of control for servo motors



microscope with powerful laser-diagnostic

hardware and interfaces.





Two cameras can be mounted on the headpiece of the microscope; one coaxially with the viewing axis of the





microscope and one mounted at an angle on the confocal tube assembly. The two present cameras employed are identical Q-Imaging Retiga 1300 units. In addition to these two cameras, there is a small surveillance camera that can be mounted inside the AFC (shown on the next page). The surveillance camera has a fixed window size of 640imes480 pixels/frame. The Q-Imaging 1300 cooled monochrome camera (6.7imes6.7 μ m) has a maximum window size of 1280×1024 pixels/frame. A camera upgrade was initiated in 2014 using the IMPERX GEV–B2020. The resolution of this camera is 2048 imes 2048 pixels/frame

Auxiliary Fluids Container (AFC)

The LMM provides an enclosed work area called the AFC, which is the main work area for sample cell processing and containment for fluids and shatterable materials. The AFC consists of two sealed glove ports, gloves, and an attachment port for the equipment transfer module (ETM) used for transporting experiment samples from stowage to the LMM. The AFC is fastened to the microscope body and sealed to provide a clean working space and one level of containment. Glove ports allow access to the sample area for cleaning before opening the box and experiment sample changeout or reconfiguration. The ETM can be configured to support various experiment modules and is located below the AFC, which has a pass-through for the samples. Materials are thus transferred without the risk of contamination release. The ETM is loaded with experiment modules on the ground and provides contained storage until the samples are utilized in the experiment.

During all testing in the AFC, the experiment sample is always moved relative to the "fixed" objective lens of the microscope. Movement is accomplished by first mounting the sample on the translation stage assembly, which automates the movements in the plane normal to the main axis of the microscope (i.e., in the X-Y plane). Movement along the main axis of the microscope (Z direction) is accomplished through the stage mount, which is connected to a Z-drive mechanism in the body of the microscope. Translation in the X-Y plane is done using motors on the translation stage itself. These motors receive drive signals through the single electrical fitting located at the end of the X-axis housing.


The LMM has the ability to image custom made microscopic slides for biological tests and sample wells containing colloids, along with commercial OptiCells™, which are a unique cell culture format for growing, monitoring, and transporting biological cells. OptiCells™ contain two parallel, gas-permeable, cell culture treated, polystyrene membranes attached to a standard microliter plate-sized frame. Each side has a growth area of 50 cm², total 100 cm², 75-µm thick membranes, 2 mm apart. A standard petri dish can be used for plant biology investigations.


Current Data and Future Capabilities

In situ mixing allows for the observation of samples shortly after mixing and for samples to be reinitialized at any time without intervention by the ISS crew. In situ mixing is a feature of the ACE-T specialty cell.

ocal Microscopy (2016—2020)

Confocal microscopy will be implemented using a 532-nm frequency- doubled Nd:YAG laser, a Nipkow disk confocal scanner, and an 8- to 12-bit digital camera. The scanner will allow up to 24 frames per second of confocal images to be collected by the camera. The crystal's threedimensional structure will be reconstructed by assembling the image slices with an image analysis program from which colloidal growth, structure, and dynamics can be determined. The confocal module will be attached and aligned to the side of the LMM and will access the sample through a camera port on the Leica RXA. The microscope's reflected light turret will contain a reflecting mirror to direct the light to and from the sample. The confocal upgrade includes a custom image processor, two cameras, and a confocal scanner as shown on the cover.





1.8 µm particles.



Epi-illuminated image of live C. elegans using a 10× objective. The C. elegans were cultured in an OptiCell™ on

Fluorescent image of Arabidopsis (root) taken during LMM ground testing for CARA/Petri Plant demonstration planned for SpaceX-3 mission.