

38P

N64-18119
Code 1...

FINAL PROGRESS REPORT
AND MANUAL

RUGGEDIZED INFRARED CARBON DIOXIDE ANALYZER

(NASA CONTRACT NAS2-711)

Date Submitted: March 15, 1963

Submitted To: National Aeronautics and Space Administration
Ames Research Center
Life Science Research Laboratory
Moffett Field, California

Submitted By: ⁰³⁶⁸²⁸² Beckman Instruments, Inc.
Special Projects Division
2400 Harbor Blvd.
Fullerton, California

(NASA CR-53359) OTS:
\$3.60 p/h,
\$134
m/f

OTS PRICE

XEROX \$ 3.60
MICROFILM \$ 1.34

Prepared By: Robert R. Pinter
Robert R. Pinter, Project Engineer

15 Mar. 1963 38 p

Approved By: Raymond L. Madsen
Raymond L. Madsen, Chief Engineer

Approved By: Max D. Liston
Max D. Liston, Manager, Special Projects Division

AMES IR CO₂ ANALYZER
FINAL REPORT AND MANUAL
CONTRACT NAS2-711

CONTENTS

	<u>PAGE NO.</u>
1.0 INTRODUCTION	1
1.1 References	1
1.2 Principles Of Operation Of The Pneumatic MONOBEAM Infrared Detector	1
1.2.1 CO ₂ Measurement System Description	5
2.0 INSTALLATION	8
2.1 Mounting of Detector and Electronics	8
2.2 Electrical Connections	11
2.2.1 Input Power	11
2.2.2 Recorder Output	11
2.2.3 Detector	11
2.3 Sampling System	12
3.0 OPERATION	12
3.1 Controls and Switches	12
3.1.1 Power On-Off Switch	14
3.1.2 Fine Zero	14
3.1.3 Span	14
3.1.4 Coarse Zero	14
3.1.5 Calibrate Adjust	15
3.1.6 Calibrate Test Button	15

CONTENTS (CONT'D.)

	<u>PAGE NO.</u>
3.1.7 Loop Gain Test Button	15
3.2 Warm Up	16
3.3 Sample Gas	16
3.4 Calibration	16
4.0 MAINTENANCE	18
4.1 Sample Cell Cleaning	18
4.1.1 Removal of Source and Sample Cell Block	18
4.1.2 Cleaning the Detector Windows	18
4.1.3 Cleaning of the Sample Compartment in the Source Block	20
4.1.4 Reassembly	20
4.2 Quadrature Adjustment	21
4.3 Phase Adjustment	23
4.4 Loop Gain Adjust	27
4.5 Detector Compensation	27
4.6 Resistance-Voltage Checks	30

LIST OF FIGURES

<u>FIGURE</u>	<u>TITLE</u>
1	Beckman CO ₂ Analyzer
2	MONOBEAM Pneumatic Detector Diagram
3	CO ₂ Analyzer Block Diagram
4	Suggested Installation For CO ₂ Analyzer Electronics
5	Suggested Detector Installation
6	CO ₂ Analyzer Sampling System Diagram
7	Source Removal Diagram
8	Magnetic Valve Adjustment Diagram
9	Electronics, Top View
10	Phase Adjustment Waveforms
11	Compensation Circuit Schematic
12	Schematic

AMES IR CO₂ ANALYZER

1.0 INTRODUCTION

This is the final report and instruction manual for the Ruggedized Infrared Carbon Dioxide Analyzer, NAS2-711. This report will be concerned with, 1) description of the instrument, including the theory of operation; 2) installation; 3) operating instructions and; 4) maintenance procedures.) A photograph of the completed instrument is shown in Figure 1.

1.1 References:- Beckman Proposal CS-61-131, June, 1961

Progress Report #1, January 31, 1962

Progress Report #2, March 1, 1962

Progress Report #3, April 4, 1962

Progress Report #4, May 10, 1962


Progress Report #5, June 5, 1962

Progress Report #6, July 18, 1962


Progress Report #7, August 20, 1962

Progress Report #8, September 18, 1962

1.2 Principles Of Operation Of The Pneumatic MONOBEAM Infrared Detector:- The CO₂ sensor for the instrument is a pneumatic detector shown in diagram form in Figure 2. This detector is unique in that it is composed of a sample detector and a reference detector which are optically in tandem. Modulated infrared energy from the source is transmitted through the gas sample being measured in cell #1, or the sample detector. Cell #1 is charged with CO₂. The infrared energy in the CO₂ absorption band is virtually eliminated in this first cell. The remaining energy is transmitted



BECKMAN CO₂ ANALYZER
FIGURE 1



MONOBEAM PNEUMATIC DETECTOR

FUNCTIONAL DIAGRAM

FIGURE 2

into cell #2, or the reference detector. Cell #2 is charged with COS. In cell #2, the energy at the edge of the CO₂ absorption band is absorbed by the COS. } The result is a CO₂ sensor with a sample detector highly selective to infrared energy changes in the center of the CO₂ absorption band and a reference detector, highly selective to energy changes on the CO₂ band edge. When CO₂ is present in the gas sample being measured, cell #1 experiences a greater loss in received energy than cell #2. This relative change of absorbed energy is ultimately detected and used to indicate the concentration of the CO₂ in the gas sample being measured. The energy changes are detected by measuring the effect of a volume change produced in each of the cells. A thin diaphragm, indicated by D1 in cell #1, and D2 in cell #2, moves in response to the volume change and acts as the plate of a varying capacitor formed by the diaphragms and the fixed buttons indicated by B1 and B2 respectively. Each of the capacitors is biased with a d.c. voltage and produces an electrical signal in response to a capacitance change. In order to prevent damage to the thin diaphragms, from static pressure differences between the front and the back of the diaphragms, a by-pass hole is provided connecting the volumes on each side of the diaphragms. Thus, each cell will respond only to time varying changes in the incoming infrared energy. It is important to note at this point the selection of a time-varying infrared source is not for this purpose alone but primarily to enable synchronous detection methods to be used for optimum signal to noise performance. Since the input energy is time varying, an a.c. signal is produced in each detector cell. By biasing each of the diaphragm capacitors in opposition, the a.c. signals can be

made to cancel under the reference condition of no CO_2 in the gas sample being measured. When CO_2 is present the a.c. signal from the first cell drops due to the reduced infrared energy. The second cell remains essentially unaffected, and an unbalance is produced causing a signal to appear. This signal is then amplified and synchronously rectified to indicate concentration of CO_2 . Since the detector functions with a single infrared source and utilizes a common optical path for both the sample detector and the reference detector, the name, "Pneumatic MONOBEAM Infrared Detector", is applied. Since the MONOBEAM Pneumatic Detector contains both a sample and reference detector in tandem, an inherent zero stability is obtained. This is due to the fact that both the sample detector and reference detector respond in the same manner to flat spectral changes in infrared energy due to source deterioration and intensity variation or dirt on the detector and source windows. Since the output of the MONOBEAM detector is the difference between signals obtained from the reference detector and the sample detector, these above factors will not influence the zero signal or balanced detector conditions.

1.2.1 CO_2 Measurement System Description:- A block diagram of the Ruggedized Infrared CO_2 Analyzer is shown in Figure 3. The pneumatic detector receives modulated infrared energy from the source after having been passed through the gas sample cell. The detector is biased by a voltage (E), of which a voltage E_1 is impressed on cell #1 and a voltage E_2 on cell #2. These voltages are varied to obtain a zero or balanced electrical output from the detector. The electrical output of the detector is directed into a high impedance preamplifier containing two Nuistor electrometer amplifiers and a third transistor amplifier stage. The preamplifier




FIGURE 3

CO₂ ANALYZER BLOCK DIAGRAM


is connected to a main amplifier having an adjustable gain potentiometer which is used to set the loop gain of the electronic feedback system. The output of the amplifier is connected to a synchronous detector or rectifier which is switched under the control of the clock. The clock also turns the infrared source on and off, consequently establishing synchronism between the modulated infrared energy and the synchronous detector. The synchronous detector allows the detection or rectification of the electrical signals having a frequency equal to that of the source variation. This enables it to discriminate against other sources of electrical signals which are not of the same frequency or in synchronism with the switching of the synchronous detector. This method of detection is quite commonly used in order to detect a signal of known frequency and phase in the presence of incoherent noise or other sources of electrical interference. The output of the synchronous detector is a d.c. voltage proportional to the magnitude of the time varying electrical signal output of the pneumatic detector. This d.c. is directed to the output terminals through a potentiometer which is used to control the range of the output voltage or the span of the instrument. This d.c. is also fed back to the zero adjust potentiometer to complete a feedback loop from which a ratio system is achieved. By means of this feedback loop the output voltage of the synchronous detector becomes an error signal proportional to the unbalance of the pneumatic detector. This signal is characterized by the following features:

- a. It is insensitive to variations in intensity of infrared energy due to deterioration of the source and loss in


MOUNTING PROCEDURES

1. REMOVE 4 FASTENING SCREWS - (A).
2. REMOVE AMPLIFIER FROM ITS CASE.
3. INSTALL CASE AS SUGGESTED.
4. SLIDE CHASSIS BACK IN CASE AND MAKE SURE TO ENGAGE PROPERLY OVER STOPPINS - (B).
5. REINSTALL THE 4 PANEL FASTENING SCREWS.

NOTE: DO NOT REMOVE DOUBLER - (C) and (D).



SUGGESTED INSTALLATION FOR CO₂ ANALYZER AMPLIFIER



SUGGESTED DETECTOR INSTALLATION

FIGURE 5

2.2 Electrical Connections:- There are three receptacles provided on the front panel of the control unit. One is for power input, one for output and one for connection to the detector.

2.2.1 Input Power:- This instrument has been designed to operate on 380 to 420 c.p.s., 110 to 120 VRMS only. The instrument has an isolation transformer so that a grounded power supply may be used. The power receptacle connections are tabulated below.

<u>Pin</u>	<u>Connection</u>
A	Chassis Ground
B	Unused
C	400 ± 20 c.p.s., 115 ± 5 VRMS
D	400 ± 20 c.p.s., 115 ± 5 VRMS
E	Unused

The chassis ground, pin A, may be used to bond the chassis to ground if required. The instrument is protected by a 1/2 amp. 3 AG fuse. Accidental connection of the instrument to 60 c.p.s., 115 V power will cause the fuse to blow but will not cause other damage.

2.2.2 Recorder Output:- The instrument output is available at pins C and D of the recorder receptacle. Pin C is the high side and pin D is the low side. The remaining pins are unused. The output is in parallel with the voltmeter on the panel. Full scale on the voltmeter is approximately 2 volts d.c. The impedance of the recorder output is about 20,000 ohms.

2.2.3 Detector:- The detector connector has an interlock so that the input power is interrupted if the detector is disconnected while the instrument is on. This protects the preamplifier filament power supply Zener which will become overheated if the instrument is operated with

the filaments disconnected. The electronics may be operated without the detector connected, by shorting pins K and L and connecting a 60 ohm, 1 watt resistor across pin H and J of the twelve pin connector located in the lower center of the front panel.

- 2.3 Sampling System:- A sampling system is required to drive the sample gas through the detector sample cell. A typical sampling system is shown in Figure 6. It consists of a suitable vacuum pump, a restrictor valve for flow control, a flow gauge and a catheter sampling tube. The catheter tube should have a 1/16" I.D. and may be up to six feet in length. Larger inside diameters will result in an apparent loss in response due to mixing in the tube. Smaller diameters will result in excessive pressure drops causing a possible instrument span error.

IMPORTANT NOTE: ALL RESTRICTIONS IN THE SAMPLING SYSTEM FOR THE PURPOSE OF FLOW CONTROL SHOULD BE LOCATED DOWNSTREAM OF THE DETECTOR. ANY RESTRICTION OF THE SAMPLE GAS UPSTREAM OF THE DETECTOR MAY CAUSE A REDUCTION OF PRESSURE IN THE SAMPLE CELL, RESULTING IN SERIOUS SPAN ERRORS.

The sample vacuum pump should be capable of drawing a continuous sample of 1000 ml per minute. A sample rate of 500 ml per minute is recommended for normal operation.

3.0 OPERATION

- 3.1 Controls and Switches:- The instrument is equipped with the following controls:

Power On-Off Switch

Fine Zero




FIGURE 6

Span

Coarse Zero

Calibrate Test Button

Loop Gain Test Button

- 3.1.1 Power On-Off Switch:- The power switch is a double pole single throw. When it is off, both sides of the line are disconnected from the instrument. Power can only be applied to the instrument when the detector plug is connected because of an electrical interlock included in the detector connector. The input power circuit is completed through pins K and L of the detector connector. Application of power is indicated by the red indicator lamp on the front panel.
- 3.1.2 Fine Zero:- The fine zero is used to zero the instrument under the presence of 0% CO₂ in the sample cell. The panel meter provides the visual read-out for making the zero adjustment. Final zero adjustments should be made after allowing the instrument to warm up for at least one hour.
- 3.1.3 Span:- The span control is used to adjust the instrument for a particular output voltage or meter reading for a given percentage of CO₂ in the sample cell. The instrument is limited by the sample cell to about 10% CO₂. The span control is capable of providing a full scale reading for about 4% CO₂ in the sample cell. The span control was factory adjusted to yield a full scale deflection with 8% CO₂ in the sample cell.

NOTE: ALL SPAN CONTROL ADJUSTMENTS SHOULD BE MADE WITH THE RECORDER CONNECTED TO THE RECORDER CONNECTOR BECAUSE OF THE LOADING EFFECTS OF THE RECORDER.

- 3.1.4 Coarse Zero:- The coarse zero is a screw driver adjustment available through a hole in the front panel. This control provides essentially the same action as the fine zero but with less resolution. The coarse zero

was factory adjusted for instrument zero with the fine zero in its mid position. The coarse zero should require little attention since most of the instrument drift will be well within the range of the fine zero adjustment.

3.1.5 Calibrate Adjust:- The calibrate adjust is a screw driver adjustment available through a hole in the front panel. The range of this control will vary between instruments but should yield about 1/2 scale deflection when the calibrate test button is pressed. The calibrate control is adjusted with the calibrate test button pressed to provide a given output indication with 0% CO₂ in the sample cell. This adjustment is usually made after the span adjust has been set, using no CO₂ gas mixtures in the sample cell. After the span is set, the calibrate control is adjusted to provide a given meter or output indication within its range with 0% CO₂ in the sample cell.

3.1.6 Calibrate Test Button:- The calibrate test when pressed produces an electrical signal equivalent to CO₂ in the sample cell. The exact amount of CO₂ that it represents is determined by the span and calibrate adjust controls. The purpose of this test button is to determine the status of the instrument or make span adjustments in the absence of known CO₂ sample gas mixtures.

3.1.7 Loop Gain Test Button:- This control indicates the instrument loop gain. This test should only be done with 0% CO₂ in the sample cell. The output indication when this test button is pressed is an inverse function of the loop gain. For this instrument, loop gain indications should be between 2% and 4% of full scale for satisfactory operation.

3.2 Warm Up:- Maximum accuracy and stability will be attained after the instrument has reached a stable operating temperature. This will take approximately one hour if the instrument has been installed in accordance with paragraph 2.1. If the detector is not mounted to the suggested heat sink, the detector will experience a greater temperature rise (30°C) and a longer period of time will be required to reach a stable operating point (4 to 6 hours). Operation of the detector in this manner will not cause damage. The electronics portion of the instrument does not require a heat sink under normal operating conditions.

3.3 Sample Gas:- The sample gas is drawn through the detector sample cell with a small vacuum pump (not supplied with the instrument). A pump that can continuously draw a sample of 1000 milliliters per minute with a suitable length of 1/16" I.D. catheter is recommended. The response time of 0.2 seconds to 90% of final reading is obtainable when the sample gas flow is at least 500 milliliters per minute. Lower sample gas flow rates will result in slower responses.

IMPORTANT NOTE: ADJUSTMENT OF THE FLOW RATE SHOULD BE ACCOMPLISHED BY A THROTTLING VALVE INSERTED BETWEEN THE DETECTOR AND PUMP. ADJUSTMENT OF THE FLOW RATE BY RESTRICTING THE SAMPLE GAS INLET SIDE OF THE DETECTOR MAY RESULT IN A REDUCTION OF PRESSURE IN THE SAMPLE CELL WHICH COULD CAUSE A CONSIDERABLE SPAN ERROR.

3.4 Calibration:- This instrument was designed for operation into a recorder. The meter on the front panel was provided to aid in zero

adjustment and initial setup of the instrument only. All calibration should be accomplished with the intended recorder connected.

Calibration of the instrument is accomplished by the measurement of known concentrations of CO_2 . Calibration should be done after the instrument is installed as outlined in Section 2.0 and allowed to reach stable operation as described in paragraph 3.2. For best results, the calibration should be performed by duplicating, as closely as possible, the conditions under which the instrument is to be used, such as, mounting, ambient temperature, flow rate, and catheter lengths. The calibrating gas should be allowed to enter the sample cell through the catheter. Care should be taken in connecting the catheter to the calibrating gas supply so that no restriction occurs upstream of the detector. This can be accomplished by placing the catheter loosely in a larger tube coming from the calibrating gas supply. The calibration gas flow must be greater than that drawn by the catheter. Another method would be to connect the calibrating gas supply to the catheter through a plastic bag to prevent pressure changes in the sample cell due to unequal flow rates. The CO_2 sensor is nonlinear so more than one calibrate gas mixture will be required to determine the transfer characteristic. The number of mixtures required will be determined by the expected accuracy. The diluent gas will affect the calibration of the instrument because of an effect referred to as "collision broadening", where the absorption characteristic is affected by the diluent gas. However, the maximum expected error in going from a diluent of 100% O_2 to 100% N_2 is less than 1%.

The instrument zero can be accomplished by evacuation of the sample cell or by using 100% N_2 as a sample gas. The use of normal room temperature may also be employed as the CO_2 content is usually

only about .05%. The use of the span and calibrate adjust and push button is described in paragraph 3.1.

4.0 MAINTENANCE

The following procedures, with the exception of paragraph 4.1, describe critical factory adjustments which should not be attempted unless absolutely necessary and only when the proper equipment is available.

4.1 Sample Cell Cleaning:- Continued use of the instrument will result in an accumulation of dirt in the sample cell. This reduces the magnitude of the infrared energy reaching the detector, resulting in a reduction of loop gain. The calibration of the instrument will not be seriously affected, however, the response of the instrument will be reduced by a proportional amount. The sample cell may be easily cleaned with the aid of the following procedure:


4.1.1 Removal of Source and Sample Cell Block:-

- a. Disconnect sampling tubing to Luer Lock connectors.
- b. Unscrew one (1) cap screw (6-40 thread) at the cable end of the cover.
- c. Unscrew both screws (10-32 thread) from cable end of detector. NOTE: NEVER REMOVE SCREWS ON OTHER END.
- d. Pull out the source and sample cell block approximately 1/16", then tilt it and remove as shown in Figure 7.

4.1.2 Cleaning the Detector Windows:- The window on the detector is now exposed for cleaning. Wipe it with a tooth brush and detergent solution. In case of hard crust deposit, such as dry mucus, a razor blade can be used to remove the deposit. Due to the hardness of the window material (synthetic sapphire) no damage will occur. Rinse with clean

FIGURE 7

SOURCE REMOVAL DIAGRAM



AW
12/12/62

water and dry with a clean cloth. CAUTION: PREVENT SOLUTION FROM RUNNING INTO PREAMPLIFIER SECTION UNDER THE COVER.

- 4.1.3 Cleaning of the Sample Compartment in the Source Block:- Remove the "O" ring and dispose of it in case of damage or deterioration. The complete block can now be cleaned by immersion into a detergent solution. The window, in case of hard deposit, can be scraped with a pen knife or similar tool. Final clean and dry the same as the detector window. Replace "O" ring. CAUTION: CARE SHOULD BE TAKEN SO AS NOT TO ACCIDENTALLY BEND THE SOURCE CONTACT PINS.


The bores, .055 diameter, on inlet and outlet of the sample cell can be cleaned with a wire or suitable size drill. The Luer Lock fittings can be unscrewed. When reinstalling them, turn until seated against the nylon washer then tighten from 1/4 to 1/2 turn. Replacements for the Luer Lock insert and its nylon seal washer can be obtained from Beckman Instruments, Inc.

- 4.1.4 Reassembly:- Follow the disassembly sequence in reverse. In case of difficulties, loosen the other cover screw by one to two turns. Do not remove the cover unless absolutely necessary. Make sure the spring contacts to the source pins are properly engaged.

The soft metal gasket (indium) between the cover and detector body is not an absolute necessity but it will help to keep dirt out of the preamplifier section. The material can be obtained from Beckman Instruments, Inc., Special Projects Division. Insert narrow strip, approximately 1/8" wide, between the cover and body. Tighten the cover down evenly with both screws and cut off the excess with a

pocket knife. CAUTION: BE CAREFUL NOT TO GET STRIP IN CONTACT WITH EXPOSED WIRES INSIDE.

4.2 Quadrature Adjustment:- A quadrature voltage will exist when the detector is balanced if the phase angle between the detector and reference cells is not exactly 180° . The quadrature voltage is undesirable. It serves no useful purpose and causes an increase of ripple at the output of the instrument. The quadrature voltage could also result in amplifier saturation if it becomes excess. The quadrature voltage is minimized by adjustment of one of two magnetically operated pneumatic valves located in each cell of the detector. The quadrature voltage is excessive if the 40 c.p.s. ripple at the output of the instrument exceeds .050 volts peak to peak. It is possible to reduce the ripple to less than .025 volts peak to peak by adjustment of one of the magnetic valves with a strong horseshoe permanent magnet, having an air gap of approximately $1/4$ ". The magnetic valves are located under nonmagnetic stainless steel covers located on the back and side of the detector. The valve stem has a magnetic bar located directly under the round stainless steel cover. The valve may be rotated with the magnet as shown in Figure 8. Each of the two valves is rotated while the 40 c.p.s. ripple on the instrument output is observed. Each valve is rotated independently until the ripple is minimized. There is an infinite combination of the two relative valve positions that will result in minimum ripple. The most desirable combination is that which results when one of the valves is fully closed as this results in maximum detector signal to noise ratio. It may be necessary to readjust the phase and loop gain controls after the quadrature voltage has been minimized.



MAGNETIC VALVE ADJUSTMENT
DIAGRAM

FIGURE 8

An alternate method for minimizing the quadrature and one that may have to be used if the quadrature signal is excessive is given below.

- a. Prepare the instrument for normal operation, except that the electronics should be removed from its enclosure.
- b. Remove the demodulator board completely from its socket.
- c. Connect a cathode ray oscilloscope (CRO) between the base of Q9 on the amplifier and the instrument common.
- d. Apply power to the instrument and allow it to warm up for several minutes.
- e. Close fully both magnetic valves, using the magnet as shown in Figure 8.
- f. Adjust the fine zero, and coarse zero if necessary for minimum 20 c.p.s. signal on the CRO.
- g. Rotate one of the magnetic valves until there is a noticeable change in the signal appearing on the CRO.
- h. Readjust the fine zero for minimum signal. If this value is less than the value observed in "f", repeat "g" and "h" until a minimum minimum is achieved. If this value is larger than that of step "f", return the valve to the fully closed position and go back to step "g", rotating the other valve.

This procedure should result in a quadrature signal of 1/2 volt peak to peak or less. This adjustment is somewhat difficult because of the harmonic content of the waveform.

After the quadrature has been minimized, using the alternate procedure, the phase and loop gain adjustments should be checked.

ture component. The phasing is adjusted with the aid of the following procedure:

- a. Prepare the instrument for normal operation, except that the electronics should be removed from its enclosure.
- b. Remove the demodulator board completely from its socket, Figure 9.
- c. Connect a dual trace CRO as follows: Channel A between the emitter end of R52 and instrument common and channel B between clock board pin A and instrument common. A single trace CRO may be used if the square wave, ± 12 volts, 20 c.p.s., available at pin A of the clock board and instrument common, is used to trigger the CRO horizontal sweep.
- d. Adjust the fine zero and coarse zero, if necessary, to obtain 3 volts peak to peak on channel A, as shown in Figure 10. Two positions of the zero pot, one on either side of balance, will yield the required signals, which will be about 180° apart. CW rotation of the zero control should yield curve (1). The phase pot should be rotated until curves (1) and (2) are symmetrically located with respect to the square wave on channel B, Figure 10.

CAUTION: EXTREME ROTATION OF THE PHASE POT MAY RESULT IN A COMPLETE PHASE REVERSAL, THAT IS, CW ROTATION OF THE ZERO POT FROM BALANCE YIELDING CURVE (2) FIGURE 10, INSTEAD OF CURVE (1). THE INSTRUMENT WOULD BE INOPERATIVE UNDER THESE CONDITIONS AS IT WOULD BE IMPOSSIBLE TO OBTAIN

Demodulator Board




FIGURE 9




FIGURE 10

Proper phasing of the instrument will result in maximum obtainable loop gain. After phasing it may be necessary to retouch the loop gain control as described in paragraph 4.4.

- 4.4 Loop Gain Adjust: - The instrument loop gain is controlled by the Trimpot located on the amplifier board as shown in Figure 9. The loop gain should be increased until the response of the instrument is critically damped. This can be accomplished by rotating the loop gain pot CW until the output "rings" slightly when the catheter is quickly placed in a stream of 8% CO₂. The flow rate through the sample cell should be 500 ml per minute and the output displayed on a CRO, or other suitable recorder. (Do not attempt to use the panel meter to determine the response). The gain should then be reduced until the ringing in response to stepped changes in CO₂ just disappears. Do not attempt loop gain adjustments unless the phasing, paragraph 4.3, and quadrature, paragraph 4.2, are known to be correct.

- 4.5 Detector Compensation: - The compensation network consists of C16, C17, C36, and R15. The purpose of this network is to minimize the response of the preamplifier to changes in the cell bias due to feedback. A simplified schematic of the compensation network and associated circuitry is shown in Figure 11. D₁ and D₂ are the detector and reference cell capacities respectively. R13, R15, and C36, are located in the preamplifier while C16 and C17 are located with the electronics. C_g represents the grid to ground capacities associated with V₁, and any strays that may exist. When the network is properly adjusted, that is, characteristics of the network made up of C16, C17, C36, and R15, are identical to those of the network made up of D₁, D₂,




FIGURE 11

C_8 , and R_{13} , the voltages appearing at the grid and cathode of V_1 , as a result of a voltage applied at F.B., are equal and do not appear on the plate of V_1 . The changes in detector output appearing at the grid only of V_1 , as a result of changes in the voltage appearing at point F.B., are not affected by the compensation network. Misadjustment of the compensation network will result in erratic operation of the instrument and possible oscillation. The compensation may be adjusted with the aid of the following procedure.

- a. Prepare the instrument for normal operation except that the electronics should be removed from its enclosure.
- b. Remove the demodulator and clock boards completely from their sockets.
- c. Connect a CRO between the base of Q_9 and instrument common.
- d. Connect an audio oscillator between the center tap of the coarse zero and instrument common.
- e. Remove completely C_{16} and C_{17} , located on the power supply board, and temporarily connect in their places, two capacitor decade boxes variable from .001 to 1 mfd.
- f. Adjust the decade capacitor in place of C_{16} to .25 mfd and the decade capacitor in place of C_{17} to .15 mfd.
- g. Adjust the output of the audio oscillator until a signal of about 15 volts peak to peak appears on the CRO. (A slight clipping may occur.)
- h. Adjust C_{17} for minimum signal.
- i. Increase C_{16} by .01 mfd and readjust C_{17} for minimum

signal.

- j. If the signal resulting from step "i" is less than that resulting from step "h", repeat step "i" until a minimum minimum signal is achieved. If the comparison shows that it is larger, repeat step "i", decreasing C16 by .01 mfd. Increase the output of the signal generator as the correct adjustment is approached, until the output of the signal generator is 3 volts peak to peak. Correct adjustment of the compensation network is achieved when voltage on the base of Q9 is between 3 and 6 volts peak to peak and the output of the signal generator is 3 volts peak to peak. It may be necessary to use smaller increments (.001 mfd) in step "i" to achieve this.
- k. Replace the decade capacitors with fixed values made up of 2 or 3 paper or mylar dielectric capacitors as required, and recheck to see if the compensation is still within the limits outlined in step "j". (A difference may occur when the fixed values are installed, due to the dissipation factors of the decade and fixed capacitors.)

4.6

Resistance-Voltage Checks:- The following resistance-voltage check list may be useful in isolating a failure that may occur. Table I shows some of the power supply voltages that should exist at the printed circuit board sockets. These measurements may be taken with all boards removed. The clock output should be measured with all boards in place. The clock output should be a ± 12 volt squarewave.

board should be equal and 180° out of phase.

Table II shows the resistance and continuity checks that may be made at the connector of the detector assembly.

Board Pin	Amp	Demodulator	Clock
A	- 12	Clock	Clock
B	E _O	Clock*	Blank
C	E _O *	Blank	- 25
D	+ 12	Common	+ 3
E	Blank	Key	Common
F	Blank	E _{out}	Key
H	Key	+ 12	Source Drive
J	Common	- 12	+ 12
K	Blank	E _{in} *	Clock*
L	E _{in}	E _{in}	- 12

PRINTED CIRCUIT BOARD CONNECTOR PIN IDENTIFICATION

TABLE I

* Inverted

Connector Pin	Resistance To Case	Remarks
A	100K	
B	>100 Meg	<.1 ohm to pin F
C	>100 Meg	
D	22 to 26K	
E	50K to 10 Megs	This reading depends upon type of ohmmeter used and polarity of connection.
F	>100 Meg	
H	10 to 20 ohms	Preamplifier filaments, cold.
J	<.1 ohm	
K	>100 Meg	<.1 ohm to pin L
L	>100 Meg	
M	>100 Meg	10 to 11 ohms to pin N
N	>100 Meg	

DETECTOR ASSEMBLY RESISTANCE AND CONTINUITY CHART

TABLE II

