

N66 26657

FACILITY FORM 602

(ACCESSION NUMBER)	(THRU)
68	1
(PAGES)	(CODE)
CR 15723	04
(NASA CR OR TMX OR AD NUMBER)	

JET PROPULSION LABORATORY
CALIFORNIA INSTITUTE OF TECHNOLOGY
PASADENA, CALIFORNIA

GPO PRICE \$ _____

CFSTI PRICE(S) \$ _____

Hard copy (HC) 3.00

Microfiche (MF) .75

PROGRESS REPORT

Microorganism Study. JPL Contract No. 950783

Prof. W. B. Bollen, Microbiologist

Oregon State University, Corvallis, Oregon

April 9, 1965

This work was performed for the Jet Propulsion Laboratory,
California Institute of Technology, sponsored by the
National Aeronautics and Space Administration under
Contract NAS7-100.

10560

PROGRESS REPORT

Microorganism Study. CIT Contract No. 950783

W. B. Bollen, Microbiologist

Oregon State University, Corvallis, Oregon

April 9, 1965

Work on 19 Bacillus cultures has been completed, except for precise growth temperature range determinations which await modification and adaptation of a borrowed temperature gradient apparatus. Studies on the other cultures submitted in the first two lots are nearly complete and work has been started on the third collection delivered February 25. A graduate assistant has been assigned to the project and is now giving full time to the systematic studies. Progress will now be more rapid.

A descriptive chart for each Bacillus is presented herewith. The species have been identified according to Berkey's Manual, and Smith, Gordon, and Clark's (Aerobic Sporeforming Bacteria) characters considered most significant. Important variations are indicated.

Each descriptive chart is supplemented by an attached sheet giving a more complete characterization of the colony characteristics.

The morphology descriptions are supplemented with photomicrographs of nigrosin, Gram, and flagella stains, taken at 1000X. Each division on the micrometer scale has a value of 0.5 micron. Dimensions given are from measurements of 10 or more cells on the nigrosin slide. Size range as well as the mean is shown under Cell Morphology. Where colors are named and preceded by a number and letters--e.g., 2ba Pearl-- this indicates the identification given in the Color Harmony Manual, published by the Container Corporation of America. I have found this system to be accurate, rapid, and adaptable to microbial cultures and it is the most comprehensive set of color standards available.

Other details on the descriptive charts are self-explanatory. Ages, temperatures and kind of media are indicated. Unless otherwise stated, trypticase soy media were used.

A full discussion of these isolants, their characteristics, and the methods will be presented at a later date. Several of the cultures have been very frustrating to work with because of apparent contamination. Now, however, we believe the differently appearing colonies are evidence of dissociation. Absolute proof of this would require single cell isolations, but we do not have the necessary micromanipulators.

W.B. Bollen




4-9-65 -

Name of organism *Bacillus megaterium* Studied by Dr. W.B. Bollen Culture No 10-2-B

Source.

Habitat Soil

Date 4-8-65

Descriptions (Underline required terms.)		Sketches X1000
CELL MORPHOLOGY Medium: Trypticase soy agar Temp. 28 °C. Vegetative cells: Age: 24 hr. Form and arrangement: streptococci, diplococci, micrococci, sarcines, rods, commas, spirals, branched rods, filaments. Motility in broth: Flagella: Peritrichous Size: 1.9 x 3.6 μ. Irregular forms: 1.5 - 2.5 x 1.5 - 6.5 μ. Sporangia: none, rods, spindles, elliptical, cleavate, drumstick. Age: Endospores: Shape: spherical, ellipsoid, cylindrical. Position: central to eccentric, terminal, subterminal.	N: grossin	 $\delta \text{ dia.} = 0.5 \mu$
STAINING CHARACTERISTICS Gram: + Age: 24 hr. Method: Kopeloff's (modified) Special stains:	Gram stain	
AGAR STROKE Age: 5 day	Temp. 28°C.	 24 hr
AGAR COLONIES Age: 3 day	Temp. 28 °C.	 48 hr
NUTRIENT BROTH Age: 24 hour	Temp. 28°C.	
Surface growth: none, ring, pellicle, flocculent, membranous. Subsurface growth: none, turbid, granular. Amount of growth: scanty, moderate, abundant. Sediment: none, granular, flocculent, viscid, flaky.		
GELATIN STAB Age: 24 hour	Temp. 25°C.	
Liquefaction: none, crateriform, infundibuliform, napiform, saccate, stratiform. Rate: slow, moderate, rapid.		
OTHER MEDIA Age: Temp. °C.		
Potato slant: 2db Ivory		
Soybean agar slant:		Growth abundant

FERMENTATION		Temp. 25 °C.				
Medium Nutrient broth	Carbohydrate: 1%	Glucose	Lactose	Sucrose	Xylose	Mannitol
Indicator: BCP						
Acid in 3 days	+					
Acid in days						
Gas in 3 days	-					
Gas in days						

ACTION ON MILK		Temp. 28 °C.		
Indicator:	Days			
Litmus	4	30		
Reaction	=	=		
Acid curd	=	=		
Rennet curd	=	=		
Peptonization	+	+		
Reduction (before coagulation)	=	=		

ACTION ON NITRATES

Medium: 1% KNO₃ broth Temp. 28 °C.
 Nitrite: + 2 d. ; d. ; d.
 Gas (N): - 3 d. ; d. ; d.

INDOLE PRODUCTION

Medium: Tryptophan broth Age: 10 days
 Method: Kovac's Temp. 28 °C.
 Indole: present, absent.

HYDROGEN SULFIDE PRODUCTION

Medium: Lead Acetate &
Thio-iron
 H₂S: present, absent.

Age 10 days Temp. 28 °C.

RELATION TO FREE OXYGEN

TEMPERATURE RELATIONS

Growth in refrigerator (°C.): present, absent.
 Growth at room temperature (°C.): present, absent.
 Growth at 37° C.: present, absent.
 Growth at 50° C.: present, absent.

Medium: Dextrose Trypt-soy Age: 10 days
 Method: Shake Tubes Temp. 28 °C.

Aerobic growth: absent, present, better than anaerobic growth, poorer
 than anaerobic growth.

Anaerobic growth: present, absent.

Pasteurization survival, 80°C. 10 minutes: Positive

ADDITIONAL TESTS

Casein Hydrolysis: Positive
 Fat Hydrolysis: Negative
 Gelatin Hydrolysis: Positive
 Starch Hydrolysis: Negative (+ in Peptone)
 Urea Hydrolysis: Negative

NH₄ from Peptone: Positive
 Metabolism: Oxidative - Glucose
 Acetyl methyl carbinol: Negative

NH₄ as sole Nitrogen source: Positive

Sole Carbon sources: Citrate - Positive
 Glucose - Positive
 Sucrose - Positive
 Xylose - Positive

Methylene blue reduction: Negative

Crystalline Dextrins: Not observed

Salt tolerances: 2% - Positive
 7% - Positive
 10% - Positive

Cultural Characteristics of Bacterial Colonies

Culture No. 10.B

Sypticase - soy agar
8 day

I. Surface Colonies

a. Macroscopic appearance.

1. Size, mm 4mm

2. Shape: Outline round, oval, irregular, filamentous.

Elevation- effuse, flat, raised, convex, rugose, papillate, umbonate, pulvinate.

Topography- smooth, rough, wrinkled, contoured, striated Habit- compact, spreading.

3. Optical properties:

(a) Color: Color Harmony Manual No. 26a. (Pence)

(b) Appearance by reflected light- dull, opalescent, iridescent, glistening, fluorescent.

(c) Appearance by transmitted light- transparent, translucent, opaque.

b. Microscopic appearance (X100).

1. Margin- entire, granular, cleft, lobed, undulate, crenate, close, ciliate, rhizoid.

2. Internal structure- amorphous, dense, granules (fine, coarse), filamentous, striated, interlaced.

c. Consistency- moist, slimy, soft, butyrous, waxy, tough, adherent, brittle.

II. Deep colonies

a. Size, mm


b. Shape

c. Color

Name of organism Bacillus pumilus Studied by Dr. W. P. Boller Culture No. 14-2-B

Source Habitat Soil Date 4-3-65

Descriptions (Underline required terms.)		Sketches
<p>CELL MORPHOLOGY Medium: Trypticase soy agar Temp. 28 °C.</p> <p>Vegetative cells: Age: 24 hr.</p> <p>Form and arrangement: streptococci, diplococci, micrococci, sarcinae, rods, commas, spirals, branched rods, filaments.</p> <p>Motility in broth: + Flagella: Peritrichous Size: 0.65 x 2.1μ Irregular forms:</p> <p>Sporangia: none, rods, spindles, elliptical, clavate, drumstick. Age:</p> <p>Endospores: Shape: spherical, ellipsoid, cylindrical.</p> <p>Position: central to eccentric, terminal, subterminal.</p>		

<p>STAINING CHARACTERISTICS</p> <p>Gram: + Age: 24 hr. Method: Kopeloff's (modified)</p>		
<p>AGAR STROKE Age: 5 day.</p> <p>Amount of growth: scanty, moderate, abundant.</p> <p>Form: filiform, echinulate, beaded, spreading, rhizoid.</p> <p>Consistency: buoyant, viscid, membranous, brittle.</p> <p>Chromogenesis: fluorescent, iridescent, phosphorescent.</p> <p>Color: Brick red</p>	Temp. 28 °C. Gram stain	24 hr 48 hr

<p>AGAR COLONIES Age: 3 day</p> <p>Form: punctiform, circular, filamentous, rhizoid, irregular.</p> <p>Elevation: effuse, flat, raised, convex.</p> <p>Surface: smooth, contoured, radiate, concentric, rugose.</p> <p>Margin: entire, undulate, erose, filamentous, curled.</p> <p>Density: opaque, translucent.</p>	Temp. 28 °C.
<p>NUTRIENT BROTH Age: 1 day</p> <p>Surface growth: none, ring, pellicle, flocculent, membranous.</p> <p>Subsurface growth: none, turbid, granular.</p> <p>Amount of growth: scanty, moderate, abundant.</p> <p>Sediment: none, granular, flocculent, viscid, flaky.</p>	Temp. 28 °C.

<p>GELATIN STAB Age: 1 day</p> <p>Liquefaction: none, crateriform, infundibuliform, napiform, saccate, stratiform.</p> <p>Rate: slow, moderate, rapid.</p>	Temp. 25 °C.
---	--------------

<p>OTHER MEDIA</p> <p>Potato slant: 6pe Dk. Lacquer Red</p> <p>Soybean agar slant: No growth</p>	Age: Temp. °C.
---	----------------

FERMENTATION		Temp. 25 °C.		
Medium Nutrient broth	Carbohydrate: %	Glucose	Lactose	Sucrose
Indicator: BCP	1%	=	=	=
Acid in 3 days		=	=	=
Acid in days				
Gas in 3 days		=	=	=
Gas in days				

ACTION ON MILK		Temp. 28 °C.				
Indicator:	Days					
Litmus	14					
Reaction	=					
Acid curd	=					
Rennet curd	=					
Peptonization	=					
Reduction (before coagulation)	=					

ACTION ON NITRATES

Medium: 1% KNO_3 broth Temp 28 °C.
 Nitrite: -10 d. ; d.
 Gas (N): -10 d. ; d.

INDOLE PRODUCTION

Medium: Tryptophan broth Age: 10 da.
 Method: Kovac's Temp 28 °C.
 Indole: present, absent.

HYDROGEN SULFIDE PRODUCTION

Medium: Lead acetate &
 H_2S : present, absent. Thio iron

Age: 10 da. Temp 28 °C.

RELATION TO FREE OXYGEN - Catalase: Positive
 Medium: Dextrose Trypt. Soy Age: 10 day

Method: Kovac's Temp 28 °C.

Aerobic growth: absent, present, better than anaerobic growth, poorer than anaerobic growth.

Anaerobic growth: present, absent.

TEMPERATURE RELATIONS

Growth in refrigerator (°C.): present, absent.
 Growth at room temperature (28 °C.): present, absent.
 Growth at 37° C.: present, absent.
 Growth at 50° C.: present, absent.

Pasteurization survival, 80°C. 10 minutes: Positive
 ADDITIONAL TESTS

Casein Hydrolysis: Positive
 Fat Hydrolysis: Positive
 Gelatin Hydrolysis: Positive
 Starch Hydrolysis: Negative
 Urea Hydrolysis: Negative

NH_4 from Peptone: Positive

Metabolism: Oxidative

Acetyl methyl carbinol: Positive

NH_4 as sole Nitrogen source: Positive

Sole Carbon sources: Citrate - Positive
 Glucose - Positive
 Sucrose - Positive
 Xylose - Positive

Methylene blue reduction: Positive

Crystalline Dextrins: Not observed

Salt tolerances: 2% - Positive
 7% - Positive
 10% - Positive

Cultural Characteristics of Bacterial Colonies

Culture No. 146

Sypticase soy agar
3 day

I. Surface Colonies

a. Macroscopic appearance.

1. Size, mm 2 mm
2. Shape: Outline- round, oval, irregular, filamentous.
Elevation- effuse, flat, raised, convex, rugose, papillate, umbonate, pulvinate.
Topography- smooth, rough, wrinkled, contoured, striated. Habit- compact, spreading.
3. Optical properties:
 - (a) Color: Color Harmony Manual No. Spe Serra Cotta
 - (b) Appearance by reflected light- dull, opalescent, iridescent, glistening, fluorescent.
 - (c) Appearance by transmitted light- transparent, translucent, opaque.

b. Microscopic appearance (X100).


1. Margin- entire, granular, cleft, lobed, undulate, crenate, erose, ciliate, rhizoid.
2. Internal structure- amorphous, dense, granular (fine, coarse), filamentous, striated, interlaced.
- c. Consistency- moist, slimy, soft, butyrous, waxy, tough, adherent, brittle.

II. Deep colonies

- a. Size, mm
- b. Shape
- c. Color

Name of organism Bacillus pumilus Studied by Dr. N. B. Holloman Culture No. 14-2-D

Source Habitat Soil Date 4-9-65

Descriptions (Underline required terms.)	Sketches
<p>CELL MORPHOLOGY Medium: Trypticase soy agar Temp. 28°C.</p> <p>Vegetative cells: Age: 24 hr. Form and arrangement: streptococci, diplococci, micrococci, sarcinae, rods, commas, spirals, branched rods, filaments.</p> <p>Motility in broth: + Flagella: Peritrichous - Size: $0.5 \times 1.0 \approx 2.5 \mu$ Irregular forms:</p> <p>Sporangia: none, rods, spindles, elliptical, clavate, drumstick. Age:</p> <p>Endospores: Shape: spherical, ellipsoid, cylindrical. Position: central to eccentric, terminal, subterminal.</p>	
<p>STAINING CHARACTERISTICS</p> <p>Gram: + Age: 24 hr Method: Kopeloff's Special stains: (modified)</p>	 
<p>AGAR STROKE Age: 5 day</p> <p>Amount of growth: scanty, moderate, abundant. Form: filiform, echinulate, beaded, spreading, rhizoid. Consistency: butyrous, viscid, membranous, brittle. Chromogenesis: Red brown fluorescent, iridescent, photogenic.</p>	Temp. 28°C.  
<p>AGAR COLONIES Age: 3 day</p> <p>Form: punctiform, circular, filamentous, rhizoid, irregular. Elevation: effuse, flat, raised, convex. Surface: smooth, contoured, radiate, concentric, rugose. Margin: entire, undulate, erose, filamentous, curled. Density: opaque, translucent.</p>	Temp. 28°C.
<p>NUTRIENT BROTH Age: 24 hr.</p> <p>Surface growth: none, ring, pellicle, flocculent, membranous. Subsurface growth: none, turbid, granular. Amount of growth: scanty, moderate, abundant. Sediment: none, granular, flocculent, viscid, sticky.</p>	Temp. 28°C.
<p>GELATIN STAB Age: 24 hr.</p> <p>Liquefaction: none, crateriform, infundibuliform, napiform, saccate, stratiform. Rate: slow, moderate, rapid.</p>	Temp. 25°C.
<p>OTHER MEDIA</p> <p>Potato slant: 6 pe Dark Lacquer Red Soybean agar slant: No growth</p>	Temp. °C.

FERMENTATION		Temp. 25 °C.		
Medium: Nutrient broth		Glucose	Lactose	Sucrose
Carbohydrate: 1 %		=	=	=
Indicator: BCP				
Acid in 3 days		=	=	=
Acid in days				
Gas in 3 days		=	=	=
Gas in days				

ACTION ON MILK		Temp. 28°C.				
Indicator:	Days					
Litmus	14					
Reaction	=					
Acid curd	=					
Rennet curd	=					
Peptonization	=					
Reduction (before coagulation)	=					

ACTION ON NITRATES

Medium: $\frac{1}{2}$ KNO_3 broth Temp. 28°C.
Nitrite: d. ; d. ; -10 d.
Gas (N): d. ; d. ; -10 d.

INDOLE PRODUCTION

Medium: Tryptophan broth Age: 10 da.
Method: Kovac's Temp. 28°C.
Indole: present, absent.

HYDROGEN SULFIDE PRODUCTION

Medium: Lead acetate & H₂S: present, absent. Age: 10 da.
thio iron Temp. 28°C.

TEMPERATURE RELATIONS

Growth in refrigerator (°C.): present, absent.
Growth at room temperature (28°C.): present, absent.
Growth at 37° C.: present, absent.
Growth at 50° C.: present, absent.

RELATION TO FREE OXYGEN - Catalase: Positive

Medium: Dextrose Trypt. Soy Age: 10 da.
Method: Shake tubes Temp 28 °C.

Aerobic growth: absent, present, better than anaerobic growth, poorer than anaerobic growth.

Anaerobic growth: present, absent.

Pasteurization survival, 80°C. 10 minutes: Positive
ADDITIONAL TESTS

Casein Hydrolysis: Positive
Fat Hydrolysis: Positive
Gelatin Hydrolysis: Positive
Starch Hydrolysis: Negative
Urea Hydrolysis: Negative

NH₄ from Peptone: Positive

Metabolism: Oxidative

Acetyl methyl carbinol: Positive

NH₄ as sole Nitrogen source: Positive

Sole Carbon sources: Citrate - Positive
Glucose - Positive
Sucrose - Positive
Xylose - Positive

Methylene blue reduction: Positive

Crystalline Dextrins: Not observed

Salt tolerances: 2% - Positive
7% - Positive
10% - Positive

Cultural Characteristics of Bacterial ColoniesCulture No. 14dDrypticase soy agar
3 dayI. Surface Colonies

a. Macroscopic appearance.

1. Size, mm 2 mm
2. Shape: Outline- round, oval, irregular, filamentous.
Elevation- effuse, flat, raised, convex, rugose, papillate, umbonate, pulvinate.
Topography- smooth, rough, wrinkled, contoured, striated. Habit- compact, spreading.
3. Optical properties:

- (a) Color: Color Harmony Manual No. Spe Serra Cotta
- (b) Appearance by reflected light- dull, opalescent, iridescent, glistening, fluorescent.
- (c) Appearance by transmitted light- transparent, translucent, opaque.

b. Microscopic appearance (X100).

1. Margin- entire, granular, cleft, lobed, undulate, crenate, erose, ciliate, rhizoid.
2. Internal structure- amorphous, dense, granular (fine, coarse), filamentous, striated, interlaced.
- c. Consistency- moist, slimy, soft, butyrous, waxy, tough, adherent, brittle.


II. Deep colonies

- a. Size, mm
- b. Shape
- c. Color

Name of organism Bacillus pumilus Studied by Dr. W. B. Bollen Culture No. 15-2-C

Source Habitat Soil Date

Descriptions (Underline required terms.)		Sketches
CELL MORPHOLOGY Medium: Trypticase soy agar Vegetative cells: Age: 24 hr. Form and arrangement: streptococci, diplococci, micrococci, sarcinae, rods, commas, spirals, branched rods, filaments. Motility in broth: + Flagella: Peritrichous Size: 0.65 x 2.34 μ regular forms: 0.9 - 0.75 x 1.5 - 3.5 μ irregular forms: Sporangia: None, rods, spindles, elliptical, clavate, drumstick. Age: Endospores: Shape: spherical, ellipsoid, cylindrical. Position: central to eccentric, terminal, subterminal.	Temp. 28°C.	

STAINING CHARACTERISTICS Gram: + Age: 24 hr. Method: Kopeloff's (modified) Special stains:	Temp. 25°C.		
---	-------------	---	---

AGAR STROKE Age: 5 day	Temp. 25°C.		24 hr
------------------------	-------------	---	-------

AGAR COLONIES Age: 3 day	Temp. 28°C.		48 hr
--------------------------	-------------	---	-------

NUTRIENT BROTH Age: 24 hour	Temp. 28°C.	
-----------------------------	-------------	---

GELATIN STAB Age: 24 hr.	Temp. 25°C.	
--------------------------	-------------	---

OTHER MEDIA Age:	Temp. °C.	
Potato slant: 6pe Dark	Lacquer Red	
Soybean agar slant: No growth		

FERMENTATION		Temp. 25 °C.		
Medium: Nutrient broth		Glucose	Lactose	Sucrose
Carbohydrate: 1 %		-	-	-
Indicator: B C P				
Acid in 3 days		-	-	-
Acid in days				
Gas in 3 days		-	-	-
Gas in days				

ACTION ON MILK Temp. 28 °C.	
Indicator:	Days
Litmus	14
Reaction	=
Acid curd	=
Rennet curd	=
Peptonization	=
Reduction (before coagulation)	=

ACTION ON NITRATES

Medium: 1% KNO₃ broth Temp 28 °C.
 Nitrite: d. ; d. ; -10d.
 Gas (N): d. ; d. ; -10d.

INDOLE PRODUCTION

Medium: Tryptophan broth Age: 10 day
 Method: Kovac's Temp: 28 °C.
 Indole: present, absent.

HYDROGEN SULFIDE PRODUCTION

Medium: Lead Acetate & H₂S: present, absent. Age: 10 da. Temp: 28 °C.
 Thio iron

TEMPERATURE RELATIONS

Growth in refrigerator (10°C): present, absent. (+ -)
 Growth at room temperature (28°C): present, absent.
 Growth at 37° C.: present, absent.
 Growth at 50° C.: present, absent. (+ -)

RELATION TO FREE OXYGEN - Catalase: Positive

Medium: Dextrose Trypt. soy Age: 10 day
 Method: Shake tubes Temp: 28 °C.

Aerobic growth: absent, present, better than anaerobic growth, poorer than anaerobic growth.

Anaerobic growth: present, absent.

Pasteurization survival, 80°C. 10 minutes: Positive

ADDITIONAL TESTS

Casein Hydrolysis: Positive
 Fat Hydrolysis: Positive
 Gelatin Hydrolysis: Positive
 Starch Hydrolysis: Negative
 Urea Hydrolysis: Negative

NH₄ from Peptone: Positive

Metabolism: Oxidative

Acetyl methyl carbinol: Positive

NH₄ as sole Nitrogen source: Positive

Sole Carbon sources: Citrate - Positive
 Glucose - Negative
 Sucrose - Negative
 Xylose - Positive

Methylene blue reduction: Positive

Crystalline Dextrins: Not observed

Salt tolerances: 2% - Positive
 7% - Positive
 10% - Positive

Cultural Characteristics of Bacterial Colonies

Culture No.

15c

I. Surface Colonies

a. Macroscopic appearance.

1. Size, mm 2 mm2. Shape: Outline- round, oval, irregular, filamentous.Elevation- effuse, flat, raised, convex, rugose, papillate, umbonate, pulvinate.Topography- smooth, rough, wrinkled, contoured, striated. Habit- compact, spreading.3. Optical properties:(a) Color: Color Harmony Manual No. ~~5~~ ^{Sp} Sereia Cotta(b) Appearance by reflected light: dull, opalescent, iridescent, glistening, fluorescent.(c) Appearance by transmitted light- transparent, translucent, opaque.


b. Microscopic appearance (X100).

1. Margin- entire, granular, cleft, lobed, undulate, ericate, erose, ciliate, rhizoid.2. Internal structure- amorphous, dense, granular (fine, coarse), filamentous, striated, interlaced.c. Consistency- moist, slimy, soft, butyrous, waxy, tough, adherent, brittle.II. Deep coloniesa. Size, mmb. Shapec. ColorSoy agar
3 day

Name of organism Bacillus subtilis var. niger Studied by Dr. W.B. Bollen Culture No. 16C
 Source White Mountain Habitat Soil Date _____

Descriptions (Underline required terms.)		Sketches
CELL MORPHOLOGY Medium: <u>Trypticase soy</u> Vegetative cells: Age: <u>24 hrs</u> Form and arrangement: <u>streptococci, diplococci, micrococci, sarcines, rods, commas, spirals, branched rods, filaments.</u> Motility in broth: + Flagella: <u>Peritrichous</u> Size: <u>0.75 x 2.8 μ</u> Irregular forms: <u>oval, round, rod, bacillus, elliptical, coccoid, drumstick.</u> Age: Sporangia: <u>none, rod, bacillus, elliptical, coccoid, drumstick.</u> Age: Endospores: Shape: <u>spherical, ellipsoid, cylindrical.</u> Position: <u>central to eccentric, terminal, subterminal.</u>		Temp. 28 °C. 

STAINING CHARACTERISTICS		
Gram: + - Age: 24 hr Method: Modified Kopeloff Special stains:		MIGROSIN 

AGAR STROKE		
Age: 24 hrs.		Temp. 25 °C. 24 hrs. GRAM STRAIN 

AGAR COLONIES		
Age: _____		Temp. °C.

NUTRIENT BROTH		
Age: 2 days		Temp. 25 °C.

GELATIN STAB		
Age: 24 hours		Temp. 25 °C.

OTHER MEDIA		
Potato slant:	Age: 2ge	Covert Tan
Soybean agar slant:	3lg	Adobe Brown - abundant
Glucose nitrate slant:		- scant growth
Tyrosine agar slant:		Black pigment
Fat agar:	3nl	Dark Brown

FERMENTATION		Temp. 25 °C.				ACTION ON MILK		Temp. 28 °C.	
Medium: Nutrient broth	Carbohydrate: 1% Indicator: BCP	Glucose	Lactose	Sucrose	Maltose				
Acid in 2 days		-	-	+	-				
Acid in 8 days		-	-	+	-				
Gas in 2 days		-	-	-	-				
Gas in 8 days		-	-	-	-				

Indicator:	Days	
Litmus	2	4
Reaction	=	ALKALINE
Acid curd		
Rennet		
Peptonization	+	
Reduction (before coagulation)	+	

ACTION ON NITRATES

Medium: 1% KNO₃ broth Temp. 28°C.
 Nitrite: + 3 d. : d. ; d.
 Gas (N): - 3 d. : d. ; d.

INDOLE PRODUCTION

Medium: Tryptophan broth Age: 10 days
 Method: Kovac's Temp. 28 °C.
 Indole: present, absent.

HYDROGEN SULFIDE PRODUCTION

Medium: Lead Acetate & Thio iron
 H₂S: present, absent.

Age 10 days
 Temp. 28°C.

RELATION TO FREE OXYGEN - Catalase: Positive
 Medium: Dextrose Nutrient a Age: 10 days
 Method: Shake Tubes Temp. 28°C.

TEMPERATURE RELATIONS

Growth in refrigerator (10°C.): present, absent.
 Growth at room temperature (28°C.): present, absent.
 Growth at 37° C.: present, absent.
 Growth at 50° C.: present, absent. (sl.)

Aerobic growth: absent, present, better than anaerobic growth, poorer than anaerobic growth.
 Anaerobic growth: present, absent.

Pasteruization survival, 80°C. 10 minutes: Positive

ADDITIONAL TESTS

Casein Hydrolysis: Positive
 Fat Hydrolysis: Negative
 Gelatin Hydrolysis: Positive
 Starch Hydrolysis: Positive
 Urea Hydrolysis: Negative

NH₄ from Peptone: Positive

Metabolism: Fermenter (glucose and sucrose) anaerogenic
 Nonoxidizer-nonfermenter (lactose and xylose)

Acetyl methyl carbinol: Positive

NH₄ as sole Nitrogen source: Positive

Sole Carbon sources: Citrate - Positive
 Glucose - Positive
 Sucrose - Positive
 Xylose - Positive

Methylene blue reduction: Positive

Crystalline Dextrins: Not observed

Salt tolerances: 2% - Positive
 7% - Positive
 10% - Positive

Cultural Characteristics of Bacterial ColoniesCulture No. 16cSugarcane - soy agar
6 dayI. Surface Colonies



a. Macroscopic appearance.

1. Size, mm 2mm2. Shape: Outline- round, oval, irregular, filamentous.Elevation- effuse, flat, raised, convex, rugose, papillate, umbonate, pulvinate.Topography- smooth, (rough) wrinkled, contoured, striated. Habit- compact, spreading.3. Optical properties:(a) Color: Color Harmony Manual No. 3pa (Chocolate Brown)(b) Appearance by reflected light- dull, opalescent, iridescent, glistening, fluorescent.(c) Appearance by transmitted light- transparent, translucent, opaque.

b. Microscopic appearance (X100).

1. Margin- entire, granular, cleft, lobed, undulate, eremate, erose, ciliate, rhizoid.2. Internal structure- amorphous, dense, granular (fine, coarse), filamentous, striated, interlaced.c. Consistency- moist, slimy, soft, butterous, waxy, tough, adherent, brittle.II. Deep coloniesa. Size, mmb. Shapec. Color

Name of organism Bacillus subtilis Studied by Dr. W. B. Bollen Culture No. 16C1
 var. niger
 Source White Mountain Habitat Soil Date 4-8-65

Descriptions (Underline required terms.)		Sketches	
CELL MORPHOLOGY	Medium: <u>Trypticase soy</u> Vegetative cells: Age: 24 hrs Form and arrangement: streptococci, diplococci, micrococci, sarcines, rods, commas, spirals, branched rods, filamentous. Motility in broth: + Flagella: Peritrichous Size: $0.96 \times 2.8 \mu$ Irregular forms: $0.75 - 1.0 \times 3.0 - 4.5 \mu$ Sporangia: none rods, spindles, elliptical, clavate, drumstick. Age: Endospores: Shape: spherical, ellipsoid, cylindrical. Position: central to eccentric, terminal, subterminal, nor swollen	Temp. 28 °C.	 Nigroso
STAINING CHARACTERISTICS	Gram: + - Age: 24 hrs Method: Kopeloff's (modified) Special stains:		
AGAR STROKE	Age: 24 hrs. Amount of growth: scanty, moderate, abundant. Form: filiform, echinulate, beaded, spreading, rhizoid. Consistency: butyrous, viscid, membranous, brittle. Chromogenesis: white; fluorescent, iridescent, photogenic.	Temp. 25°C.	 24 hr 48 hr
AGAR COLONIES	Age: _____	Temp. °C.	
NUTRIENT BROTH	Age: 2 days Surface growth: none, ring, pellicle, flocculent, membranous. (islands) Subsurface growth: none, turbid, granular. Amount of growth: scanty, moderate, abundant. Sediment: none, granular, flocculent, viscid, flaky.	Temp. 25 °C.	 Flagella
GELATIN STAB	Age: 24 hrs Liquefaction: none, crateriform, infundibuliform, napiform, saccate, stratiform. Rate: slow, moderate, rapid.	Temp. 25 °C.	(7)
OTHER MEDIA	Age: _____	Temp. °C.	
Potato slant:	31e	Biscue	
Soybean agar slant:	31g	Adobe Brown	
Glucose Nitrate slant:			scant growth
Tyrosine agar slant:		Black	
Fat agar:	31g	to white (mainly white)	

FERMENTATION		Temp. 25 °C.			
Medium: Nutrient broth		Glucose	Lactose	Sucrose	Xylose
Carbohydrate: 1 %		+	=	+	=
Indicator: BCP					
Acid in 2 days					
Acid in days					
Gas in 2 days		-	-	-	-
Gas in days					

ACTION ON MILK		Temp. 28 °C.
Indicator:	Days	
Litmus	2	
Reaction	-	
Acid curd		
Rennet curd		
Peptonization	+	
Reduction (before coagulation)	+	

- (aeroline)

ACTION ON NITRATES

Medium: 1% KNO₃ broth
 Nitrite: 3 d. ; d.
 Gas (N): 3 d. ; d.

Temp. 28°C.
 ; d.
 ; d.

INDOLE PRODUCTION

Medium: Tryptophan broth
 Method: Kovac's
 Indole: present, absent.

Age: 10 days
 Temp. 28°C.

HYDROGEN SULFIDE PRODUCTION

Medium: Lead Acetate &
Thio iron
 H₂S: present, absent.

Age: 10 days
 Temp. 28°C.

RELATION TO FREE OXYGEN - Catalase: Positive
 Medium: Dextrose Nutrient a Age: 10 days
 Method: Shake tubes Temp. 28°C.

TEMPERATURE RELATIONS

Growth in refrigerator (10°C.): present, absent.
 Growth at room temperature (28°C.): present, absent.
 Growth at 37° C.: present, absent.
 Growth at 50° C.: present, absent.

Aerobic growth: absent, present, better than anaerobic growth, poorer than anaerobic growth.
 Anaerobic growth: present, absent.

Pasteurization survival, 80°C 10 minutes: Positive
 ADDITIONAL TESTS

Casein Hydrolysis: Positive
 Fat Hydrolysis: Negative
 Gelatin Hydrolysis: Positive
 Starch Hydrolysis: Positive
 Urea Hydrolysis: Negative

NH₄ from Peptone: Positive

Metabolism: Fermenter - anaerogenic (Glucose and Sucrose)
 Nonoxidizer - nonfermenter (Lactose and Xylose)

Acetyl methyl carbinol: Positive

NH₄ as sole Nitrogen source: Positive

Sole Carbon sources: Citrate - Positive
 Glucose - Positive
 Sucrose - Positive
 Xylose - Positive

Methylene blue reduction: Positive

Crystalline Dextrans: Not observed

Salt tolerances: 2% - Positive
 7% - Positive
 10% - Positive

Cultural Characteristics of Bacterial Colonies

Culture No. 16C1

Soybean soy agar
5 day

I. Surface Colonies

a. Macroscopic appearance.

1. Size, mm 1 mm

2. Shape: Outline- round, oval, irregular, filamentous.

Elevation- effuse, flat, raised, convex, rugose, papillate, umbonate, pulvinate.

Topography- smooth, rough, wrinkled, contoured, striated. Habit- compact, spreading.

3. Optical properties:

(a) Color: Color Harmony Manual No. 3gc St. San

(b) Appearance by reflected light- dull, opalescent, iridescent, glistening, fluorescent.

(c) Appearance by transmitted light- transparent, translucent, opaque.

b. Microscopic appearance (X100).

1. Margin- entire, granular, cleft, lobed, undulate, cremate, erose, ciliate, rhizoid.

2. Internal structure- amorphous, dense, granular (fine, coarse), filamentous, striated, interlaced.

c. Consistency- moist, slimy, soft, butyrous, waxy, tough, adherent, brittle.

II. Deep colonies




a. Size, mm

b. Shape

c. Color

Name of organism Bacillus pumilus Studied by Dr. W.B. Bollen Culture No. 18

Source White Mountain Habitat Soil Date 4-8-65

Descriptions (<u>Underscore required terms.</u>)	Sketches
CELL MORPHOLOGY Medium: Trypticase soy agar Temp. 30 °C. Vegetative cells: Age: 24 hrs. Form and arrangement: streptococci, diplococci, micrococci, sarcines, rods, commas, spirals, branched rods, filaments. Motility in broth: + Flagella: Peritrichous Size: $0.85 \times 2.15 \mu$ Irregular forms: $0.75 - 1.0 \times 1.5 - 3.5 \mu$ Sporangia: none, rods, spindles, elliptical, clavate, drumstick. Age: Endospores: Shape: spherical, ellipsoid, cylindrical. Position: central to eccentric, terminal, subterminal.	
STAINING CHARACTERISTICS Gram: + - Age: 24 hr Method: Kopeloff's (modified) Special stains:	
AGAR STROKE Age: 24 hour Amount of growth: scanty, moderate, abundant. Form: filiform, echinulate, beaded, spreading, rhizoid. Consistency: butyrous, viscid, membranous, brittle. Chromogenesis: ; fluorescent, iridescent, phosphorescent.	Temp. 25 °C. 
AGAR COLONIES Age: 2 day Form: punctiform, circular, filamentous, rhizoid, irregular. Elevation: effuse, flat, raised, convex. Surface: smooth, contoured, radiate, concentric, rugose. Margin: entire, undulate, erose, filamentous, curled. Density: opaque, translucent.	Temp. 28 °C. 
NUTRIENT BROTH Age: 12 days Surface growth: none, ring, pellicle, flocculent, membranous. islands Subsurface growth: none, turbid, granular. Amount of growth: scanty, moderate, abundant. Sediment: none, granular, flocculent, viscid, fleshy.	Temp. 25 °C. 
GELATIN STAB Age: Liquefaction: none, crateriform, infundibuliform, napiform, saccate, striform. Rate: slow, moderate, rapid.	Temp. 25 °C.
OTHER MEDIA Potato slant: 21e Mustard Soybean agar slant: 3ca Pearl Pink Fat agar: Clear to lt. yellow	Temp. °C.

FERMENTATION		Temp. 25 °C.			
Medium: Nutrient broth		Glucose	Lactose	Sucrose	Xylose
Carbohydrate: 1 %					
Indicator: BCP					
Acid in 2 days	+	=	+	=	
Acid in days					
Gas in 2 days	-	-	-	-	
Gas in days					

Indicator:	Days		
	3	5	7
Litmus			
Reaction	-	=	-
Acid curd	=	=	=
Rennet curd	=	+	+
Peptonization	+	+	+
Reduction (before coagulation)	+	+	+

ACTION ON NITRATES

Medium: 1% KNO₃ broth Temp 28 °C.
 Nitrite: d. ; d. ; -3 d.
 Gas (N): d. ; d. ; -3 d.

INDOLE PRODUCTION

Medium: Tryptophan
 Method: Kovac's
 Indole: present, absent.

Age: 10 day
 Temp. 28°C.

HYDROGEN SULFIDE PRODUCTION

Medium: Lead acetate & H₂S: present, absent. Age: 10 day
 Thio iron Temp. 28°C.

TEMPERATURE RELATIONS

Growth in refrigerator (10°C.): present, absent.
 Growth at room temperature (28°C.): present, absent.
 Growth at 37° C.: present, absent.
 Growth at 50° C.: present, absent.

RELATION TO FREE OXYGEN - Catalase: Positive

Medium: Dextrose Nutrient a Age: 10 day
 Method: Shake tubes Temp. 28°C.

Aerobic growth: absent, present, better than anaerobic growth, poorer than anaerobic growth.
 Anaerobic growth: present, absent.

Pasteurization survival, 80°C. 10 minutes: Positive

ADDITIONAL TESTS

Casein Hydrolysis: Positive
 Fat Hydrolysis: Positive
 Gelatin Hydrolysis: Positive
 Starch Hydrolysis: Negative
 Urea Hydrolysis: Negative

NH₄ from Peptone: Positive

Metabolism: Fermenter - anaerogenic - Glucose, Sucrose, and Xylose
 Nonoxidizer-nonfermenter - Lactose

Acetyl methyl carbinol: Positive

NH₄ as sole Nitrogen source: Positive

Sole Carbon sources: Citrate - Positive
 Glucose - Positive
 Sucrose - Positive
 Xylose - Positive

Methylene blue reduction: Positive

Crystalline Dextrins: Not observed

Salt tolerances: 2% - Positive
 7% - Positive
 10% - Positive

Cultural Characteristics of Bacterial Colonies

Culture No. 18

I. Surface ColoniesSoybean-casein agar
2 days


a. Macroscopic appearance.

1. Size, mm 2mm2. Shape: Outline- round, oval, irregular, filamentous.Elevation- effuse, flat, raised, convex, rugose, papillate, umbonate, pulvinate.Topography- smooth, rough, wrinkled, contoured, striated. Habit- compact, spreading.3. Optical properties:(a) Color: Color Harmony Manual No. 26b (Bamboo)(b) Appearance by reflected light- dull, opalescent, iridescent, glistening, fluorescent.(c) Appearance by transmitted light- transparent, translucent, opaque.

b. Microscopic appearance (X100).

1. Margin- entire, granular, cleft, lobed, undulate, crenate, erose, ciliate, rhizoid.2. Internal structure- amorphous, dense, granular (fine, coarse), filamentous, striated, interlaced.c. Consistency- moist, slimy, soft, butyrous, waxy, tough, adherent, brittle.II. Deep coloniesa. Size, mmb. Shapec. Color

Name of organism Bacillus cereus Studied by Dr. W. B. Bollen Culture No. 18-2-C
 Source White Mountain Habitat Soil Date 4-8-65

Descriptions (Underline required terms.)		Sketches
CELL MORPHOLOGY	Medium: Trypticase soy agar Temp. 30 °C. Vegetative cells: Age: 24 hrs Form and arrangement: streptococci, diplococci, micrococci, sarcinae, rods, commas, spirals, branched rods, filaments. Motility in broth: + Flagella: Peritrichous Size: $2 \times 4.7\mu$ Irregular forms: $2 \times 3.0 - 7.0\mu$ Sporangia: none, rods, spindles, elliptical, elevate, drumstick. Age: Endospores: Shape: spherical, ellipsoid, cylindrical. Position: central to excentric, terminal, subterminal.	N: gross: 
STAINING CHARACTERISTICS	Gram: + Age: 24 hr Method: Kopeloff's (modified) Special stains:	 
AGAR STROKE	Age: 5 day	Temp. 28 °C. Gram stain 24 hr 48 hr
AGAR COLONIES	Age: 6 day	Temp. 28 °C.
NUTRIENT BROTH	Age: 24 hour	Temp. 28 °C.
GELATIN STAB	Age: 24 hour	Temp. 25 °C.
OTHER MEDIA	Age:	Temp. °C.
Potato slant:	2ec	Biscuit
Soybean agar slant:		Abundant growth

FERMENTATION		Temp. 25 °C.		
Medium: Nutrient broth		Glucose	Lactose	Sucrose
Carbohydrate: 1%				
Indicator: BCP				
Acid in days	+			
Acid in days				
Gas in days	-			
Gas in days				

ACTION ON MILK		Temp. 25 °C.		
Indicator:	Litmus	Days		
Reaction	4	10	14	
Acid curd	=	=	=	
Rennet curd	=	=	=	
Peptonization	=	+	+	
Reduction (before coagulation)	+	+	+	

ACTION ON NITRATES

Medium: 1% KNO₃ broth Temp. 28°C.
 Nitrite: d. ; d. +2 d.
 Gas (N): d. ; d. -2 d.

INDOLE PRODUCTION

Medium: Tryptophan broth Age: 10 day
 Method: Kovac's Temp. 28°C.
 Indole: present, absent.

HYDROGEN SULFIDE PRODUCTION

Medium: Lead acetate &
 H₂S: present, absent. Thio iron Age: 10 da. Temp. 28°C.

RELATION TO FREE OXYGEN - Catalase: Negative
 Medium Dextrose Trypt. sov aAge: 10 da.
 Method: Shake tubes Temp. 28°C.

TEMPERATURE RELATIONS

Growth in refrigerator (°C.): present, absent.
 Growth at room temperature (28°C.): present, absent.
 Growth at 37° C.: present, absent.
 Growth at 50° C.: present, absent.

Aerobic growth: absent, present, better than anaerobic growth, poorer than anaerobic growth.
 Anaerobic growth: present, absent.

Pasteurization survival, 80°C. 10 minutes: Positive

ADDITIONAL TESTS

Casein Hydrolysis: Positive
 Fat Hydrolysis: Negative
 Gelatin Hydrolysis: Positive
 Starch Hydrolysis: Positive
 Urea Hydrolysis: Negative

NH₄ from Peptone: Positive

Metabolism: Nonoxidizer-nonfermenter - Glucose

Acetyl methyl carbinol: Positive

NH₄ as sole Nitrogen source: Positive

Sole Carbon sources: Citrate - Positive
 Glucose - Positive
 Sucrose - Positive
 Xylose - Positive

Methylene blue reduction: Negative

Crystalline Dextrins: Not observed

Salt tolerances: 2% - Positive
 7% - Positive
 10% - Positive

Cultural Characteristics of Bacterial ColoniesCulture No. 18cSypticase soy agar
6 days.I. Surface Colonies

a. Macroscopic appearance.




1. Size, mm 5 mm2. Shape: Outline- round, oval, irregular, filamentous.
Elevation- effuse, flat, raised, convex, rugose,
papillate, umbonate, puvinate.Topography- smooth, rough, wrinkled, contoured,
striated. Habit- compact, spreading.3. Optical properties:(a) Color: Color Harmony Manual No. 1/2ca (Cream)(b) Appearance by reflected light- dull, opalescent,
iridescent, glistening, fluorescent.(c) Appearance by transmitted light- transparent,
translucent, opaque.

b. Microscopic appearance (X100).

1. Margin- entire, granular, cleft, lobed, undulate,
cremate, erose, ciliate, rhizoid.2. Internal structure- amorphous, dense, granular (fine,
coarse), filamentous, striated, interlaced.c. Consistency- moist, slimy, soft, butyrous, waxy, tough,
adherent, brittle.II. Deep coloniesa. Size, mmb. Shapec. Color

Name of organism Bacillus pumilus Studied by Dr. W.B. Bollen Culture No. 19

Source White Mountain Habitat Soil Date 4-8-65

Descriptions (Underline required terms.)		Sketches
CELL MORPHOLOGY Medium: Trypticase soy agar temp. 30 °C. Vegetative cells: Age: 24 hr Form and arrangement: streptococci, diplococci, micrococci, sarcines, rods, commas, spirals, branched rods, filaments. Motility in broth: + Flagella: Peritrichous Size: 0.74 x 2.5 μ Irregular forms: Sporangia: none, rods, spindles, elliptical, cleistospores, drumstick. Age: Endospores: Shape: spherical, ellipsoid, cylindrical. Position: central to eccentric, terminal, subterminal.	Nigrosporin	
STAINING CHARACTERISTICS Gram: + - Age: 24 hr Method: Kopeloff's (modified) Special stains:	Gram stain	
AGAR STROKE Age: 24 hr.	Temp. 25 °C	 24 hr
Amount of growth: scanty, moderate, abundant. Form: filiform, echinulate, beaded, spreading, rhizoid. Consistency: glutinous, viscid, membranous, brittle. Chromogenesis: white; fluorescent, iridescent, photogenic.		48 hr
AGAR COLONIES Age: 5 day	Temp. 28 °C.	
Form: punctiform, circular, filamentous, rhizoid, irregular. Elevation: effuse, flat, raised, convex. Surface: smooth, contoured, radiate, concentric, rugose. Margin: entire, undulate, erose, filamentous, curled. Density: opaque, translucent.		
NUTRIENT BROTH Age: 2 day	Temp. 25 °C.	
Surface growth: none, ring, pellicle, flocculent, membranous. Subsurface growth: none, turbid, granular. Amount of growth: scanty, moderate, abundant. Sediment: none, granular, flocculent, viscid, sticky.		Flagella
GELATIN STAB Age: 5 day	Temp. 25 °C.	
Liquefaction: none, crateriform, infundibuliform, napiform, saccate, stratiform. Rate: slow, moderate, rapid.		
OTHER MEDIA Age: Temp. °C.		
Potato slant: 3ie Camel Soybean agar slant: 3ca Pearl pink Fat agar: clear with yellow-orange tint		scant to moderate

FERMENTATION		Temp. 25 °C.				
Medium Nutrient broth	%	Glucose	Lactose	Sucrose	Xylose	
Indicator: BCP 1						
Acid in 12 days	+	=	+	=		
Acid in days						
Gas in 12 days	-	-	-	-		
Gas in days						

ACTION ON MILK		Temp. 22 °C.	
Indicator:	Days		
Litmus	3 14		
Reaction	= -		
Acid curd	= =		
Rennet curd	= =		
Peptonization	+ +		
Reduction (before coagulation)	+ +		

ACTION ON NITRATES

Medium: 1% KNO_3 broth Temp. 28°C.
 Nitrite: d. ; d. ; -3 d.
 Gas (N): d. ; d. ; -3 d.

INDOLE PRODUCTION

Medium: Tryptophan broth Age: 10 day
 Method: Kovac's
 Indole: present, absent.

HYDROGEN SULFIDE PRODUCTION

Medium: Lead acetate & Thiocyanate Age: 16 day
 H_2S : present, absent. Temp. 28°C.

TEMPERATURE RELATIONS

Growth in refrigerator (10°C.): present, absent.
 Growth at room temperature (28°C.): present, absent.
 Growth at 37°C.: present, absent.
 Growth at 50°C.: present, absent. (+-)

RELATION TO FREE OXYGEN - Catalase: Positive

Medium: Dextrose Nutrient A Age: 16 day
 Method: Shake tubes Temp. 28°C.

Aerobic growth: absent, present, better than anaerobic growth, poorer than anaerobic growth.

Anaerobic growth: present, absent.

Pasteurization survival, 80°C 10 minutes: Positive
 ADDITIONAL TESTS

Casein Hydrolysis: Positive
 Fat Hydrolysis: Positive
 Gelatin Hydrolysis: Positive
 Starch Hydrolysis: Negative
 Urea Hydrolysis: Negative

NH_4 from Peptone: Positive

Metabolism: Fermentative - anaerogenic - Glucose and Sucrose
 Nonoxidative-nonfermentative - Lactose and Xylose

Acetyl methyl carbinol: Positive

NH_4 as sole Nitrogen source: Positive

Sole Carbon sources: Citrate - Positive
 Glucose - Positive
 Sucrose - Positive
 Xylose - Positive

Methylene blue reduction: Positive

Crystalline Dextrans: Not observed

Salt tolerances: 2% - Positive
 7% - Positive
 10% - Positive

Cultural Characteristics of Bacterial Colonies

Culture No. 19

Soy agar
5 day.I. Surface Coloniesa. Macroscopic appearance.1. Size, mm 3mm.2. Shape: Outline- round, oval, irregular, filamentous.Elevation- effuse, flat, raised, convex, rugose, papillate, umbonate, pulvinate.Topography- smooth, rough, wrinkled, contoured, striated. Habit- compact, spreading.3. Optical properties:(a) Color: Color Harmony Manual No. 41c, Pastel Orange ^{some colonies darker and mottled}(b) Appearance by reflected light- dull, opalescent, iridescent, glistening, fluorescent.(c) Appearance by transmitted light- transparent, translucent, opaque.b. Microscopic appearance (X100).1. Margin- entire, granular, cleft, lobed, undulate, crenate, erose, ciliate, rhizoid.2. Internal structure- amorphous, dense, granular (fine, coarse), filamentous, striated, interlaced.c. Consistency- moist, slimy, soft, butyrous, waxy, tough, adherent, brittle.II. Deep coloniesa. Size, mmb. Shapec. Color

Name of organism Bacillus pumilis Studied by Dr. W.B. Bollen Culture No. 19A1

Source White Mountain Habitat Soil Date 4-8-65

Descriptions (Underline required terms.)		Sketches
CELL MORPHOLOGY Medium: Trypticase soy Vegetative cells: Age: 24 hr Form and arrangement: streptococci, diplococci, micrococci, sarcinae, rods, commas, spirals, branched rods, filaments. Motility in broth: + Flagella: Peritrichous Size: 1.0 x 2.3 μ Irregular forms: 1.0 x 2.0-3.0 μ Sporangia: none, rods; spindles, elliptical, clavate, drumstick. Age: 60 Hr.	Temp. 28 °C.	
ENDOSPORES: Shape: spherical, ellipsoid, cylindrical. Position: central to eccentric, terminal, subterminal.	Nigrosin	
STAINING CHARACTERISTICS Gram: + Age: 24 Method: Kopeloff's (modified)		
AGAR STROKE Age: 24 hr. Amount of growth: scanty, moderate, abundant. Form: filiform, echinulate, beaded, spreading, rhizoid. Consistency: butyrous, viscid, membranous, brittle. Chromogenesis: white ; fluorescent, iridescent, photogenic.	Temp. 25 °C.	
AGAR COLONIES Age: 5 da.	Temp. °C.	
NUTRIENT BROTH Age: 2 day Surface growth: none, ring, pellicle, flocculent, membranous. Subsurface growth: none, turbid, granular. Amount of growth: scanty, moderate, abundant. Sediment: none, granular, flocculent, viscid, fleshy.	Temp. 25 °C.	
GELATIN STAB Age: 48 hr Liquefaction: none, crateriform, infundibuliform, napiform, saccate, stratiform. Rate: slow, moderate, rapid.	Temp. 25 °C.	
OTHER MEDIA Potato slant: Black pigment Soybean agar slant: 3ie Camel Glucose nitrate slant: Clear to white shiny (moderate growth) Tyrosine agar slant: Black water soluble pigment Fat agar: 2ca Lt. Ivory	Temp. °C.	

FERMENTATION		Temp. 25 °C.				
Medium: Nutrient broth		Glucose	Lactose	Sucrose	Maltose	
Carbohydrate: 1 %		=	=	+	=	
Indicator: BCP						
Acid in 2 days		=	=	+	=	
Acid in days						
Gas in 2 days		=	=	=	=	
Gas in days						

Indicator: Litmus	ACTION ON MILK Temp. 28 °C.		
	3	5	7
Reaction		-	
Acid curd			
Rennet curd		+	
Peptonization	+		
Reduction (before coagulation)	+		

ACTION ON NITRATES

Medium: $\frac{1}{5}$ KNO₃ broth Temp. 28°C.
 Nitrite: -3 d. ; d.
 Gas (N): -3 d. ; d.

INDOLE PRODUCTION

Medium: Tryptophan broth Age: 10 day
 Method: Kovac's
 Indole: present, absent.

HYDROGEN SULFIDE PRODUCTION

Medium: Lead Acetate & Thio iron
 H₂S: present, absent.

Age: 10 day Temp. 28°C.

TEMPERATURE RELATIONS

Growth in refrigerator (10°C.): present, absent. ()
 Growth at room temperature (20°C.): present, absent.
 Growth at 37° C.: present, absent.
 Growth at 50° C.: present, absent.

RELATION TO FREE OXYGEN - Catalase: Negative
 Medium: Dextrose Nutrient a Age: 10 day

Method: Shake tubes Temp. 28°C.

Aerobic growth: absent, present, better than anaerobic growth, poorer than anaerobic growth.

Anaerobic growth: present, absent.

Pasteurization survival , 80°C 10 minutes: Positive

ADDITIONAL TESTS

Casein Hydrolysis: Positive
 Fat Hydrolysis: Negative
 Gelatin Hydrolysis: Positive
 Starch Hydrolysis: Negative
 Urea Hydrolysis: Negative

NH₄ from Peptone: Positive

Metabolism: Fermenter - anaerogenic (Glucose and Sucrose)
 Nonoxidizer-nonfermenter (Lactose and Xylose)

Acetyl methyl carbinol: Positive

NH₄ as sole Nitrogen source: Positive

Sole Carbon sources: Citrate - Positive
 Glucose - Positive
 Sucrose - Positive
 Xylose - Positive

Methylene blue reduction: Positive

Crystalline Detrins: Not observed

Salt tolerances: 2% - Positive
 7% - Positive
 10% - Positive

Cultural Characteristics of Bacterial Colonies

Culture No. 19A1

Drypticase soy agar
5 days

I. Surface Colonies

a. Macroscopic appearance.

1. Size, mm 4 mm

2. Shape: Outline- round, oval, irregular, filamentous.

Elevation- effuse, flat, raised, convex, rugose, papillate, umbonate, pulvinate.

Topography- smooth, rough, wrinkled, contoured, striated. Habit- compact, spreading.

3. Optical properties:

(a) Color: Color Harmony Manual No. 200 St. Quarry

(b) Appearance by reflected light- dull, opalescent, iridescent, glistening, fluorescent.

(c) Appearance by transmitted light- transparent, translucent, opaque.

b. Microscopic appearance (X100).

1. Margin- entire, granular, cleft, lobed, undulate, crenate, close, ciliate, rhizoid.

2. Internal structure- amorphous, dense, granular (fine, coarse), filamentous, striated, interlaced.

c. Consistency- moist, slimy, soft, butyrous, waxy, tough, adherent, brittle.

II. Deep colonies

a. Size, mm

b. Shape

c. Color

Name of organism Bacillus pumilus Studied by Dr. W.B. Bollen Culture No. 19A2

Source Habitat Soil Date 4-3-65

Descriptions (Underline required terms.)		Sketches.
CELL MORPHOLOGY Medium: Trypticase soy agar Temp. 28 °C. Vegetative cells: Age: 24 hr. Form and arrangement: streptococci, diplococci, micrococci, sarcinae, rods, commas, spirals, branched rods, filaments. Motility in broth: + Flagella: Peritrichous Size: 0.65×2.45 Irregular forms: $0.3 \times 1.2 \times 1.9$ Sporangia: none, rods, spindles, elliptical, clavate, drumstick. Age: Endospores: Shape: spherical, ellipsoid, cylindrical. Position: central to eccentric, terminal, subterminal.	Nigrosin	
STAINING CHARACTERISTICS Gram: + Age: 24 hr Method Kopeloff's (modified) Special stains:		
AGAR STROKE Age: 24 hr. Amount of growth: scanty, moderate, abundant. Form: gliform, echinulate, beaded, spreading, rhizoid. Consistency: butyrous, viscid, membranous, brittle. Chromogenesis: fluorescent, iridescent, photogenic.	Temp. 25 °C. Gram Stain	24 hr 48 hr
AGAR COLONIES Age: 7 day Form: punctiform, circular, filamentous, rhizoid, irregular. Elevation: effuse, flat, raised, convex. Surface: smooth, contoured, radiate, concentric, rugose. Margin: entire, undulate, erose, filamentous, curled. Density: opaque, translucent.	Temp. 25 °C.	
NUTRIENT BROTH Age: 2 day Surface growth: none, ring, pellicle, flocculent, membranous. Subsurface growth: none, turbid, granular. Amount of growth: scanty, moderate, abundant. Sediment: none, granular, flocculent, viscid, sticky.	Temp. 25 °C.	Flagella
GELATIN STAB Age: 2 day Liquefaction: none, crateriform, infundibuliform, napiform, saccate, striform. Rate: slow, moderate, rapid.	Temp. 25 °C.	
OTHER MEDIA Age: Temp. °C.		
Potato slant: 3le Cinnamon Soybean agar slant: 2db Ivory Fat agar:		(abundant) scant

FERMENTATION		Temp. 25 °C.			
Medium: Nutrient broth	Carbohydrate: 1%	Glucose	Lactose	Sucrose	
Indicator: BCP					
Acid in 1 days	+	=	+	-	
Acid in 12 days	+	+	+	-	
Gas in 1 days	=	=	=	=	
Gas in 12 days	=	=	=	=	

ACTION ON MILK		Temp 28 °C.	
Indicator:	Days		
Litmus	8 35		
Reaction	= +		
Acid curd	= =		
Rennet curd	= =		
Peptonization	= =		
Reduction (before coagulation)	+	+	

ACTION ON NITRATES

Medium: 1% KNO_3 broth Temp. 28°C.
 Nitrite: -3 d. ; d.
 Gas (N): -3 d. ; d.

INDOLE PRODUCTION

Medium: Tryptophan broth Age 10 day
 Method: Kovac's
 Indole: present, absent.

HYDROGEN SULFIDE PRODUCTION

Medium: Lead acetate & H₂S: present, absent. Thio iron
 Age: 10 day Temp. 28°C.

TEMPERATURE RELATIONS

Growth in refrigerator (10°C.): present, absent. (+-)
 Growth at room temperature (28°C.): present, absent.
 Growth at 37°C.: present, absent.
 Growth at 50°C.: present, absent. (-+>)

RELATION TO FREE OXYGEN - Catalase: Positive

Medium: Dextrose Nutrient Age: 10 day
 Method: Shake tubes Temp. 28°C.

Aerobic growth: absent, present, better than anaerobic growth, poorer than anaerobic growth.

Anaerobic growth: present, absent.

Pasteurization survival, 80°C. 10 minutes: Positive
 ADDITIONAL TESTS

Casein Hydrolysis: Positive
 Fat Hydrolysis: Positive
 Gelatin Hydrolysis: Positive
 Starch Hydrolysis: Negative
 Urea Hydrolysis: Negative

NH_4 from Peptone: Positive

Metabolism: Fermentative - anaerogenic - Glucose and Sucrose
 Nonoxidative-nonfermentative - Lactose and Xylose

Acetyl methyl carbinol: Positive

NH_4 as sole Nitrogen source: Positive

Sole Carbon sources: Citrate-- Positive
 Glucose - Positive
 Sucrose - Negative
 Xylose - Positive

Methylene blue reduction: Positive

Crystalline Dextrins: Not observed

Salt tolerances: 2% - Positive
 7% - Positive
 10% - Positive

Cultural Characteristics of Bacterial Colonies

Culture No. 19A2

Soypticas - soy agar
2 dayI. Surface Colonies

a. Macroscopic appearance.



1. Size, mm 2mm2. Shape: Outline- round, oval, irregular, filamentous.Elevation- effuse, flat, raised, convex, rugose, papillate, umbonate, pulvinate.Topography- smooth, rough, wrinkled, contoured, striated. Habit- compact, spreading.3. Optical properties:(a) Color: Color Harmony Manual No. 2hb (Mauve)(b) Appearance by reflected light- dull, opalescent, iridescent, glistening, fluorescent.(c) Appearance by transmitted light- transparent, translucent, opaque.

b. Microscopic appearance (X100).

1. Margin- entire, granular, cleft, lobed, undulate, cremate, erose, ciliate, rhizoid.2. Internal structure- amorphous, dense, granular (fine, coarse), filamentous, striated, interlaced.c. Consistency- moist, slimy, soft, butyrous, waxy, tough, adherent, brittle.II. Deep coloniesa. Size, mmb. Shapec. Color

Name of organism Bacillus pumilus Studied by Dr. W.B. Bollen Culture No. 19DL

Source Habitat Soil Date Apr. 26, 1965

Descriptions (Underline required terms.)		Sketches
CELL MORPHOLOGY	Medium: Trypticase soy agar Temp. 28 °C. Vegetative cells: Age: 24 hr. Form and arrangement: streptococci, diplococci, micrococci, sarcines, rods, commas, spirals, branched rods, filaments. Motility in broth: + Flagella: Peritrichous Size: 0.65×2.5 Irregular forms: $0.5 - 0.75 \times 1.5 - 3.5$ Sporangia: none, rods, spindles, elliptical, clavate, drumstick. Age: Endospores: Shape: spherical, ellipsoid, cylindrical. Position: central to eccentric, terminal, subterminal.	
STAINING CHARACTERISTICS	Gram: Age: 24 hr. Method: Kopeloff's (modified) Special stains:	 
AGAR STROKE	Age: 24 hour Amount of growth: scanty, moderate, abundant. Form: filiform, echinulate, beaded, spreading, rhizoid. Consistency: buoyant, viscid, membranous, brittle. Chromogenesis: White ; fluorescent, iridescent, photogenic.	Temp. 25 °C. Gram stain 24 hr. 48 hr.
AGAR COLONIES	Age: 5 day Form: punctiform, circular, filamentous, rhizoid, irregular. Elevation: effuse, flat, raised, convex. Surface: smooth, contoured, radiate, concentric, rugose. Margin: entire, undulate, erose, filamentous, curled. Density: opaque, translucent.	Temp. 28 °C.
NUTRIENT BROTH	Age: 2 day Surface growth: none, ring, pellicle, flocculent, membranous. islands Subsurface growth: none, turbid, granular. Amount of growth: scanty, moderate, abundant. Sediment: none, granular, flocculent, viscid, sticky.	Temp. 25 °C.
GELATIN STAB	Age: 5 day Liquefaction: none, crateriform, infundibuliform, napiform, saccate, stratiform. Rate: slow, moderate, rapid.	Temp. 25 °C. Flagella
OTHER MEDIA	Age: Temp. °C.	
Potato slant:	3ng Yellow	maple
Soybean agar slant:	3ca	
Fat agar:	2gc	Clear to

FERMENTATION		Temp. 25 °C.			
Medium Nutrient broth:	Indicator: 1% RCP	Glucose	Lactose	Sucrose	Maltose
Carbohydrate:		=	=	+	=
Indicator:					
Acid in 2 days					
Acid in days					
Gas in 2 days		-	-	-	-
Gas in days					

ACTION ON MILK		Temp. 28 °C.			
Indicator: Litmus	Days	8	21		
Reaction		=	=		
Acid curd		=	=		
Rennet curd		=	=		
Peptonization		=	+		
Reduction (before coagulation)		+	+		

ACTION ON NITRATES

Medium: 1% KNO₃ broth Temp. 28°C.
 Nitrite: d. : d. : -3 d.
 Gas (N): d. : d. : -3 d.

INDOLE PRODUCTION

Medium: Tryptophan broth Age: 10 da.
 Method: Kovac's
 Indole: present, absent

HYDROGEN SULFIDE PRODUCTION

Medium: Lead Acetate & Age: 10 da.
 H₂S: present, absent. Thio iron Temp. 28°C.

RELATION TO FREE OXYGEN - Catalase: Positive
 Medium: Dextrose Nutrient a Age: 10 da.
 Method: Shake tubes Temp. 28°C.

TEMPERATURE RELATIONS

Growth in refrigerator (10°C.): present, absent. ()
 Growth at room temperature (25°C.): present, absent.
 Growth at 37° C.: present, absent.
 Growth at 50° C.: present, absent. ()

Aerobic growth: absent, present, better than anaerobic growth, poorer
 than anaerobic growth.
 Anaerobic growth: present, absent.

Pasteurization survival, 80°C. 10 minutes: Positive

ADDITIONAL TESTS

Casein Hydrolysis: Positive
 Fat Hydrolysis: Positive
 Gelatin Hydrolysis: Positive
 Starch Hydrolysis: Negative
 Urea Hydrolysis: Negative

NH₄ from Peptone: Positive

Metabolism: Fermentative - anaerogenic - Glucose and Sucrose
 Nonoxidative-nonfermentative - Lactose and Xylose

Acetyl methyl carbinol: Positive

NH₄ as sole Nitrogen source: Positive

Sole Carbon sources: Citrate - Positive
 Glucose - Positive
 Sucrose - Positive
 Xylose - Positive

Methylene blue reduction: Positive

Crystalline Dextrans: Not observed

Salt tolerances: 2% - Positive
 7% - Positive
 10% - Positive

Cultural Characteristics of Bacterial Colonies

Culture No. 19d1

Soyticasce soy agar
5 dayI. Surface Colonies

a. Macroscopic appearance.

1. Size, mm 3mm

2. Shape: Outline- round, oval, irregular, filamentous.
Elevation- effuse, flat, raised, convex, rugose, papillate, umbonate, pulvinate.

Topography- smooth, rough, wrinkled; contoured, striated. Habit- compact, spreading.

3. Optical properties:

(a) Color: Color Harmony Manual No. ~~265~~ Bamboo

(b) Appearance by reflected light- dull, opalescent, iridescent, glistening, fluorescent.

(c) Appearance by transmitted light- transparent, translucent, opaque.

b. Microscopic appearance (X100).

1. Margin- entire, granular, cleft, lobed, undulate, crenate, erose, ciliate, rhizoid.

2. Internal structure- amorphous, dense, granular (fine, coarse), filamentous, striated, interlaced.

c. Consistency- moist, slimy, soft, butyrous, waxy, tough, adherent, brittle.



II. Deep colonies

a. Size, mm

b. Shape

c. Color

Name of organism Bacillus subtilis
var. niger Studied by Dr. W.B. Bollen Culture No. 19G2-2
Source..... Habitat Soil Date 4-8-65

Descriptions (Underline required terms.)	Sketches
CELL MORPHOLOGY Medium: Trypticase soy agar Temp. 28 °C. Vegetative cells: Age: 24 hrs. Form and arrangement: streptocci, diplococci, micrococci, sarcines, rods, communi, spirals, branched rods, filaments. Motility in broth: + Flagella: Peritrichous Size: $0.75 \times 2.45\mu$ Irregular forms: $0.75 \times 2.0 - 3.0\mu$ Sporangia: none, rods, spindles, elliptical, dense, dramatic. Age: Endospores: Shape: spherical, ellipsoid, cylindrical. Position: central to eccentric, terminal, subterminal.	
STAINING CHARACTERISTICS Gram: + Age: 18 hr. Method Kopeloff's (modified)	 
AGAR STROKE Age: 18 hr. Amount of growth: scanty, moderate, abundant. Form: filamentous, schizomycetous, spreading, rhizoid. Consistency: butyrous, viscid, membranous, brittle. Chromogenesis: fluorescent, iridescent, photogenic. White becoming black in 10 days. Black water soluble pigment.	Temp. 25°C. Gram stain 24 hr 48 hr
AGAR COLONIES Age: 24 hr. Form: punctiform, circular, filamentous, rhizoid, irregular. Elevation: effuse, flat, raised, convex. Surface: smooth, contoured, radiculate, concentric, rugose. Margin: entire, undulate, erose, filamentous, curled. Density: opaque, translucent.	Temp. 28 °C.
NUTRIENT BROTH Age: 2 day ⁸ Surface growth: none, ring, pellicle, flocculent, membranous. islands Subsurface growth: none, turbid, granular. Amount of growth: scanty, moderate, abundant. Sediment: none, granular, flocculent, viscid, fleshy.	Temp. 25 °C.  Flagella
GELATIN STAB Age: 2 day Liquefaction: none, crateriform, infundibuliform, napiform, saccate, stratiform. Rate: slow, moderate, rapid.	Temp. 25 °C.
OTHER MEDIA Potato slant: 31e Camel Soybean agar slant: 31i Beaver abundant Glucose nitrate slant: scant Tyrosine agar slant: Black pigment Fat agar: 3nl Dark Brown	Temp. °C.

FERMENTATION		Temp. 25 °C.				
Medium Nutrient broth		Glucose	Lactose	Sucrose	Xylose	
Carbohydrate: 1%		+	=	+	-	
Indicator: BCP						
Acid in 2 days	+	=	+	-		
Acid in days						
Gas in 2 days	-	-	-	-		
Gas in days						

Indicator:	ACTION ON MILK						Temp. 28 °C.
	Days						
Litmus	4	7					
Reaction	-						
Acid curd							
Rennet curd							
Peptonization	+						
Reduction (before coagulation)		+					

ACTION ON NITRATES

Medium: 1% KNO_3 broth Temp. 28c.
 Nitrite: 3 d. ; d.
 Gas (N): 3 d. ; d.

INDOLE PRODUCTION

Medium: Tryptophan broth Age: 16 day
 Method: Kovac's
 Indole: present, absent

HYDROGEN SULFIDE PRODUCTION

Medium: Lead acetate & Age: 10 da. H₂S: present, absent
 Thio iron Temp. 28c.

TEMPERATURE RELATIONS

Growth in refrigerator (10°C.): present, absent.
 Growth at room temperature (28°C.): present, absent.
 Growth at 37° C.: present, absent.
 Growth at 50° C.: present, absent.

RELATION TO FREE OXYGEN - Catalase: Positive

Medium: Dextrose Nutrient a Age: 10 day
 Method: Shake tubes Temp. 28c.

Aerobic growth: absent, present, better than anaerobic growth, poorer than anaerobic growth.

Anaerobic growth: present, absent.

Pasteurization survival, 80°C. 10 minutes: Positive

ADDITIONAL TESTS

Casein Hydrolysis: Positive
 Fat Hydrolysis: Negative
 Gelatin Hydrolysis: Positive
 Starch Hydrolysis: Positive
 Urea Hydrolysis: Negative

NH₄ from Peptone: Positive

Metabolism: Fermenter (anaerogenic) Glucose and Sucrose
 Nonoxidizer-nonfermenter Lactose and Xylose

Acetyl methyl carbinol: Positive

NH₄ as sole Nitrogen source: Positive

Sole Carbon sources: Citrate - Positive
 Glucose - Positive
 Sucrose - Positive
 Xylose - Positive

Methylene blue reduction: Positive

Crystalline Dextrins: Not observed

Salt tolerances: 2% - Positive
 7% - Positive
 10% - Positive

Cultural Characteristics of Bacterial ColoniesCulture No. 1962-2Sypticase soy agar
24 hourI. Surface Colonies

a. Macroscopic appearance.



1. Size, mm 3mm2. Shape: Outline- round, oval, irregular, filamentous.Elevation- effuse, flat, raised, convex, rugose, papillate, umbonate, pulvinate.Topography- smooth, rough, wrinkled, contoured, striated. Habit- compact, spreading.3. Optical properties:(a) Color: Color Harmony Manual No. 39c St. Sow(b) Appearance by reflected light- dull, opalescent, iridescent, glistening, fluorescent.(c) Appearance by transmitted light- transparent, translucent, opaque.

b. Microscopic appearance (X100).

1. Margin- entire, granular, cleft, lobed, undulate, ericate, erose, ciliate, rhizoid.2. Internal structure- amorphous, dense, granular (fine, coarse), filamentous, striated, interlaced.c. Consistency- moist, slimy, soft, buryrous, waxy, tough, adherent, brittle.II. Deep coloniesa. Size, mmb. Shapec. Color

Name of organism Bacillus subtilis Studied by Dr. W. B. Bollen Culture No. 19H1

Source Habitat Soil Date 4-9-65

Descriptions (<u>Underscore required terms.</u>)	Sketches
CELL MORPHOLOGY Medium: Trypticas soy agar Temp. 28 °C. Vegetative cells: Age: 24 hrs. Form and arrangement: streptococci, diplococci, micrococci, sarcines, rods, commas, spirals, branched rods, filaments. Motility in broth: + Flagella: Peritrichous Size: $0.6 \times 2.2\mu$ Irregular forms: $0.5 - 0.75 \times 1.5 - 2.0\mu$ Sporangia: none, rods, spindles, elliptical, deviate, drumstick. Age: Endospores: Shape: spherical, ellipsoid, cylindrical. Position: central to eccentric, terminal, subterminal.	
STAINING CHARACTERISTICS Gram: + -- Age: 24 hr. Method Kopeloff's (modified) Special stains:	 
AGAR STROKE Age: 24 hrs. Amount of growth: scanty, moderate, abundant. Form: filiform, echinulate, beaded, spreading, rhizoid. Consistency: butyrous, viscid, membranous, brittle. Chromogenesis: ; fluorescent, iridescent, phosphorescent. white	Temp. 25 °C. Gram stain 24 hr 48 hr
AGAR COLONIES Age: 7 days Form: punctiform, circular, filamentous, rhizoid, irregular. Elevation: effuse, flat, raised, convex, center. Surface: smooth, contoured, radiate, concentric, rugose. Margin: entire, undulate, erose, filamentous, curled. Density: opaque, translucent.	Temp. 25 °C.
NUTRIENT BROTH Age: 12 days Surface growth: none, ring, pellicle, flocculent, membranous. Subsurface growth: none, turbid, granular. Amount of growth: scanty, moderate, abundant. Sediment: none, granular, flocculent, viscid, sticky.	Temp. 25 °C.
GELATIN STAB Age: 6 days Liquefaction: none, crateriform, infundibuliform, napiform, saccate, straiform. Rate: slow, moderate, rapid.	Temp. 25 °C.
OTHER MEDIA Age: Temp. °C.	
Potato slant: 3ne Topaz Soybean agar slant: 2ca Lt. Ivory Fat agar: 2ea Lt. wheat	

FERMENTATION		Temp. 25 °C.			
Medium Nutrient broth	1%	Glucose	Lactose	Sucrose	<u>Lysozyme</u>
Carbohydrate:	%	+	+	+	-
Indicator: BCD	1				
Acid in 2 days		+	+	+	-
Acid in days					
Gas in 2 days		-	-	-	-
Gas in days					

Indicator:	Days		
	4	5	7
Litmus	-	-	-
Reaction	=	=	=
Acid curd	=	=	=
Rennet curd	=	=	=
Peptonization	=	=	+
Reduction (before coagulation)	=	+	+

ACTION ON NITRATES

Medium 1% KNO_3 broth Temp 28 °C.
 Nitrite: -3 d. ; d.
 Gas (N): -3 d. ; d.

INDOLE PRODUCTION

Medium Tryptophan broth Age 10 day
 Method: Kovac's
 Indole: present, absent.

HYDROGEN SULFIDE PRODUCTION

Medium Lead acetate & Thiocolorin Age: 10 day
 H₂S: present, absent Temp. 28°C.

TEMPERATURE RELATIONS

Growth in refrigerator (10 °C.): present, absent.
 Growth at room temperature (28°C.): present, absent.
 Growth at 37° C.: present, absent.
 Growth at 50° C.: present, absent.

RELATION TO FREE OXYGEN - Catalase: Negative

Medium Dextrose Nutrient a Age: 10 day
 Method: Shake tubes Temp. 28°C.

Aerobic growth: absent, present, better than anaerobic growth, poorer
 than anaerobic growth.

Anaerobic growth: present, absent.

Pasteurization survival, 80°C. 10 minutes: Positive

ADDITIONAL TESTS

Casein Hydrolysis: Positive
 Fat Hydrolysis: Positive
 Gelatin Hydrolysis: Positive
 Starch Hydrolysis: Positive
 Urea Hydrolysis: Negative

NH₄ from Peptone: Positive

Metabolism: Fermenter - anaerogenic - Glucose and Sucrose
 Nonoxidizer- nonfermenter - Lactose and Xylose

Acetyl methyl carbinol: Positive

NH₄ as sole Nitrogen source: Positive

Sole Carbon sources: Citrate - Positive
 Glucose - Negative
 Sucrose - Positive
 Xylose - Negative

Methylene blue reduction: Positive

Crystalline Dextrins: Not observed

Salt tolerances: 2% - Positive
 7% - Positive
 10% - Positive

Cultural Characteristics of Bacterial Colonies

Culture No. 19H1

I. Surface Colonies

a. Macroscopic appearance.

1. Size, mm 2mm2. Shape: Outline- round, oval, irregular, filamentous.Elevation- effuse, flat, raised, convex, rugose, papillate, umbonate, pulvinate.Topography- smooth, rough, wrinkled, contoured, striated. Habit- compact, spreading.3. Optical properties:(a) Color: Color Harmony Manual No. 2 age (St. Sennon Yellow)(b) Appearance by reflected light- dull, opalescent, iridescent, glistening, fluorescent.(c) Appearance by transmitted light- transparent, translucent, opaque.

b. Microscopic appearance (X100).

1. Margin- entire, granular, cleft, lobed, undulate, eremate, erose, ciliate, rhizoid.2. Internal structure- amorphous, dense, granular (fine, coarse), filamentous, striated, interlaced.c. Consistency- moist, slimy, soft, butyrous, waxy, tough, adherent, brittle.II. Deep coloniesa. Size, mmb. Shapec. ColorSnyptcase soy agar
2 days

Name of organism Bacillus subtilis Studied by Dr. W. B. Bollen Culture No. 19H2
 var. niger Source Habitat Soil Date

Descriptions (Underline required terms.)		Sketches
CELL MORPHOLOGY	Medium: Trypticase soy agar Temp. 28 °C. Vegetative cells: Age: 24 hr. Form and arrangement: streptococci, diplococci, micrococci, sarcines, rods, commas, spirals, branched rods, filaments. Motility in broth: + Flagella: Peritrichous Size: 0.75 x 2.65 μ Irregular forms: E. I. ments Sporangia: none, rods, spores, elliptical, ovoid, dense, drumstick. Age: Endospores: Shape: spherical, ellipsoid, cylindrical. Position: central to eccentric, terminal, subterminal.	Nigrosin
STAINING CHARACTERISTICS	Gram: + - Age: 18 Method: Kopeloff's (modified) Special stains:	
AGAR STROKE 73	Age: 18 hrs.	Temp 25 °C. Gram stain
		24 hrs 48 hrs
AGAR COLONIES 73	Age: 5 days	Temp. °C.
		Flagella
NUTRIENT BROTH	Age: 2 days	Temp. 25 °C.
GELATIN STAB	Age: 2 days	Temp. 25 °C.
OTHER MEDIA	Age:	Temp. °C.
Potato slant:	3ie	Camel
Soybean agar slant:	3lg	Adobe Brown
Glucose nitrate slant:		White Abundant growth
Tyrosine agar slant:		Not blackened - slightly darkened.
Fat agar:	3nl	Dark Brown

FERMENTATION		Temp. 25 °C.				
Medium: Nutrient broth		Glucose	Lactose	Sucrose	Xylose	
Carbohydrate: 1 %		+	-	+	-	
Indicator: BCP						
Acid in 2 days		+	-	+	-	
Acid in days						
Gas in 2 days		-	-	-	-	
Gas in days						

ACTION ON MILK		Temp. 28 °C.				
Indicator:		Days				
Litmus		3	4			
Reaction	=					
Acid curd						
Rennet curd						
Peptonization			+			
Reduction (before coagulation)			+			

ACTION ON NITRATES

Medium: 1% KNO₃ broth
 Nitrite:-3d. ;d.
 Gas (N):-3d. ;d.

Temp. 28°C.
 ;d.
 ;d.

INDOLE PRODUCTION

Medium: Tryptophan broth
 Method: Kovac's
 Indole: present, absent

Age: 10 days
 Temp. 28°C.

HYDROGEN SULFIDE PRODUCTION

Medium: Lead Acetate &
 Thio Iron
 H₂S: present, absent

Age: 10 day
 Temp. 28°C.

RELATION TO FREE OXYGEN - Catalase: Positive-
 Medium: Dextrose Nutrient Age: 10 da. Negative
 Method: Shake tubes Temp. 28°C.
 Aerobic growth: absent, present, better than anaerobic growth, poorer
 than anaerobic growth.
 Anaerobic growth: present, absent.

TEMPERATURE RELATIONS

Growth in refrigerator (10°C.): present, absent. (+ -)
 Growth at room temperature (28°C.): present, absent.
 Growth at 37°C.: present, absent.
 Growth at 50°C.: present, absent.

Pasteurization survival, 80°C. 10 minutes: Positive

ADDITIONAL TESTS

Casein Hydrolysis:	Positive
Fat Hydrolysis:	Negative
Gelatin Hydrolysis:	Positive
Starch Hydrolysis:	Positive
Urea Hydrolysis:	Negative

NH₄ from Peptone: Positive

Metabolism: Fermenter - anaerogenic (Glucose and Sucrose)
 Nonoxidizer-nonfermenter (Lactose & Xylose)

Acetyl methyl carbinol: Positive

NH₄ as sole Nitrogen source: Positive

Sole Carbon sources: Citrate - Positive
 Glucose - Positive
 Sucrose - Positive
 Xylose - Positive

Methylene blue reduction: Positive

Crystalline Dextrins: Not observed

Salt tolerances: 2% - Positive
 7% - Positive
 10% - Positive

Cultural Characteristics of Bacterial ColoniesCulture No. 19H2Soybean casein agar
5 dayI. Surface Colonies

a. Macroscopic appearance.

1. Size, mm 6 mm.
2. Shape: Outline- round, oval, irregular, filamentous.
Elevation- effuse, flat, raised, convex, rugose, papillate, umbonate, pulvinate.
Topography- smooth, rough, wrinkled, contoured, striated. Habit- compact, spreading.
3. Optical properties:
 - (a) Color: Color Harmony Manual No. 3gc St. Ivory
 - (b) Appearance by reflected light- dull, opalescent, iridescent, glistening, fluorescent.
 - (c) Appearance by transmitted light- transparent, translucent, opaque.

b. Microscopic appearance (X100).



1. Margin- entire, granular, cleft, lobed, undulate, ericate, erose, ciliate, rhizoid.
2. Internal structure- amorphous, dense, granular (fine, coarse), filamentous, striated, interlaced.
3. Consistency- moist, slimy, soft, butyrous, waxy, tough, adherent, brittle.

II. Deep colonies

- a. Size, mm
- b. Shape
- c. Color

Name of organism Bacillus polymyxa Studied by Dr. W.B. Boller Culture No. 191

Source White Mountain Habitat Soil Date

Descriptions (Underline required terms.)	Sketches
<p>CELL MORPHOLOGY Medium: Trypticase soy agar temp. 28 °C.</p> <p>Vegetative cells: Age: 24 hrs</p> <p>Form and arrangement: streptococci, diplococci, micrococci, sarcines, rods, commens., spirals, branched rods, filaments.</p> <p>Motility in broth: + Flagella: Peritrichous</p> <p>Size: $0.75 - 1 \times 1.5 - 4.0 \mu$ Irregular forms:</p> <p>Sporangia: none, rods; spores, elliptical, clavate, drumstick. Age:</p> <p>Endospores:</p> <p>Shape: spherical, ellipsoid, cylindrical.</p> <p>Position: central to eccentric, terminal, subterminal.</p>	 <p>N: gros: n</p>
<p>STAINING CHARACTERISTICS</p> <p>Gram: + - Age: 24 hr Method: Kopeloff (modified)</p>	 
<p>AGAR STROKE Age: 24 hrs.</p> <p>Amount of growth: scanty, moderate, abundant.</p> <p>Form: filiform, echinulate, beaded, spreading, rhizoid.</p> <p>Consistency: butyrous, viscid, membranous, brittle.</p> <p>Chromogenesis: Ivory fluorescent, iridescent, photogenic.</p>	<p>Temp. 28 °C.</p> <p>GRAM STAIN</p> <p>24 hr. 48 hr.</p>
<p>AGAR COLONIES Age: 48 hours</p> <p>Form: punctiform, circular, filamentous, rhizoid, irregular.</p> <p>Elevation: effuse, flat, raised, convex.</p> <p>Surface: smooth, contoured, radiate, concentric, rugose.</p> <p>Margin: entire, undulate, erose, filamentous, curled.</p> <p>Density: opaque, translucent.</p>	<p>Temp. 28 °C.</p>
<p>NUTRIENT BROTH Age: 24 hrs.</p> <p>Surface growth: none, ring, pellicle, flocculent, membranous.</p> <p>Subsurface growth: none, turbid, granular.</p> <p>Amount of growth: scanty, moderate, abundant.</p> <p>Sediment: none, granular, flocculent, viscid, sticky.</p>	<p>Temp. 28 °C.</p> <p>Flagella</p>
<p>GELATIN STAB Age: 5 days</p> <p>Liquefaction: none, crateriform, infundibuliform, napiform, saccate, stratiform.</p> <p>Rate: slow, moderate, rapid.</p>	<p>Temp. 25 °C.</p>
<p>OTHER MEDIA Age: 10 days</p> <p>Potato slant: 2db Ivory</p> <p>Soybean agar slant: 3ec Bisque (abundant)</p> <p>Glucose nitrate slant: White (abundant)</p> <p>Fat agar: White</p> <p>Tyrosine agar slant: No pigment.</p>	<p>Temp. 28 °C.</p>

FERMENTATION Temp. 28 °C.				
Medium: Nutrient broth	Glucose	Lactose	Sucrose	Xylose
Carbohydrate: 1%	+	+	+	+
Indicator: BCP				
Acid in 2 days	+	+	+	+
Acid in days				
Gas in 2 days	+	+-	+	+
Gas in days				

ACTION ON MILK Temp. 28 °C.	
Indicator:	Days
Litmus	3
Reaction	
Acid curd	+ (coag. at 28°)
Rennet curd	
Peptonization	+
Reduction (before coagulation)	+

ACTION ON NITRATES

Medium: 1% KNO₃ broth Temp. 28 °C.
 Nitrite: -3 d. ; d.
 Gas (N): -3 d. ; d.

INDOLE PRODUCTION

Medium: Tryptophan
 Method: Kovac's
 Indole: present, absent.

Age: 10 days
 Temp. 30 °C.

HYDROGEN SULFIDE PRODUCTION

Medium: Lead Acetate &
 Thio-iron
 H₂S: present, absent.

Age: 10 days RELATION TO FREE OXYGEN - Catalase: Positive
 Temp. 30 °C. Medium: Dextrose Nutrient A Age: 10 days
 Method: Shake Tubes Temp. 30°C.

TEMPERATURE RELATIONS

Growth in refrigerator (10°C.): present, absent.
 Growth at room temperature (28°C.): present, absent.
 Growth at 37°C.: present, absent.
 Growth at 50°C.: present, absent.

Aerobic growth: absent, present, better than anaerobic growth, poorer
 than anaerobic growth.
 Anaerobic growth: present, absent.

Pasteurization survival, 80°C. 10 minutes: Positive

ADDITIONAL TESTS

Casein Hydrolysis: Positive
 Fat Hydrolysis: Negative
 Gelatin Hydrolysis: Positive
 Starch Hydrolysis: Positive
 Urea Hydrolysis: Negative

NH₄ from Peptone: Negative

Metabolism: Fermenter - aerogenic

Acetyl methyl carbinol: Positive

NH₄ as sole Nitrogen source: Positive

Sole Carbon sources: Citrate - Positive
 Glucose - Positive
 Sucrose - Positive
 Xylose - Positive

Methylene blue reduction: Positive

Crystalline Dextrins: Not observed

Salt tolerances: 2% - Positive
 7% - Negative
 10% - Negative

Reduces Brom Thymol Blue.

Cultural Characteristics of Bacterial Colonies

Culture No. 191

Sypticase soy agar
2 days

I. Surface Colonies

a. Macroscopic appearance.

1. Size, mm 3

2. Shape: Outline- round, oval, irregular, filamentous.

Elevation- effuse, flat, raised, convex, rugose, papillate, umbonate, pulvinate.

Topography- smooth, rough, wrinkled, contoured, striated. Habit- compact, spreading.

3. Optical properties:

(a) Color: Color Harmony Manual No. 100 (pale yellow)

(b) Appearance by reflected light- dull, opalescent, iridescent, glistening, fluorescent.

(c) Appearance by transmitted light- transparent, translucent, opaque.

b. Microscopic appearance (X100).

1. Margin- entire, granular, cleft, lobed, undulate, cremate, erose, ciliate, rhizoid.

2. Internal structure- amorphous, dense, granular (fine, coarse), filamentous, striated, interlaced.

c. Consistency- moist, slimy, soft, butyrous, waxy, tough, adherent, brittle.

II. Deep colonies



a. Size, mm

b. Shape

c. Color

Name of organism Bacillus polymyxa* Studied by Dr. W.B. Bollen Culture No. 118D

Source Habitat Soil Date 4-8-65

Descriptions (Underline required terms.)	Sketches
CELL MORPHOLOGY Medium: Trypticase soy agar Temp. 28 °C. Vegetative cells: Age: 24 Hrs. Form and arrangement: streptococci, diplococci, micrococci, sarcines, rods, commas, spirals, branched rods, filaments. Motility in broth: + Flagella: Peritrichous Size: 1.1 X 6.0 μ Irregular forms: 1.0 - 1.5 X 9.0 - 10.0 μ Sporangia: none, rods, spindles, elliptical, dense, drumstick. Age: Endospores: Shape: spherical, ellipsoid, cylindrical. Position: central to eccentric, terminal, subterminal.	 <i>Nigrosin</i>
STAINING CHARACTERISTICS Gram: + - - Age: 18 hr Method: Kopeloff's Special stains: (modified)	 <i>Gram stain</i>
AGAR STROKE Age: 18 hr. Amount of growth: scanty, moderate, abundant. Form: filiform, echinulate, beaded, spreading, rhizoid. Consistency: butyrous, viscid, membranous, brittle. Chromogenesis: white fluorescent, iridescent, photogenic.	Temp. 25 °C. 
AGAR COLONIES Age: 10 day Form: punctiform, circular, filamentous, rhizoid, irregular. Elevation: effuse, flat, raised, convex. Surface: smooth, contoured, radiate, concentric, rugose. Margin: entire, undulate, erose, filamentous, curled. Density: opaque, translucent.	Temp. 28 °C. 24 hr 48 hr 
NUTRIENT BROTH Age: 8 day Surface growth: none, ring, pellicle, flocculent, membranous. Subsurface growth: none, turbid, granular. Amount of growth: scanty, moderate, abundant. Sediment: none, granular, flocculent, viscid, fleshy.	Temp. 25 °C.  <i>Flagella</i>
GELATIN STAB Age: 10 day Liquefaction: none, crateriform, infundibuliform, napiform, saccate, stratiform. Rate: slow, moderate, rapid.	Temp. 25 °C.
OTHER MEDIA Potato slant: 2db Ivory Soybean agar slant: White abundant Glucose nitrate slant: scant Tyrosine agar slant: Black pigment Fat agar: 2ec Biscuit	Temp. °C.

FERMENTATION		Temp.	25 °C.
Medium: Nutrient broth		Glucose	
Carbohydrate: 1%		Lactose	
Indicator: BCP		Sucrose	Xylose
Acid in 2 days	+	=	+
Acid in days			-
Gas in 2 days	+		+
Gas in days			

ACTION ON MILK		Temp. 28 °C.
Indicator: Litmus	Days	
Reaction	+	5 35
Acid curd	±	(whey spot)
Rennet curd	=	
Peptonization	=	
Reduction (before coagulation)	+	

+ = acid

* Varies from typical B. polymyxa in casein and gelatin hydrolysis (negative), and salt tolerances (2% negative).

ACTION ON NITRATES

Medium: 1%KNO₃ broth Temp. 28°C.
 Nitrite: - 3 d. ; d.
 Gas (N): - 3 d. ; d.

INDOLE PRODUCTION

Medium: Tryptophan broth Age: 10 day
 Method: Kovac's Temp. 28°C.
 Indole: present absent.

HYDROGEN SULFIDE PRODUCTION

Medium: Lead acetate & H₂S: present, absent. Age: 10 da. Temp. 28°C.
 Thio iron

RELATION TO FREE OXYGEN - Catalase: Positive

Medium:Dextrose Nutrient a. Age: 10 day
 Method:Shake tubes Temp. 28°C.

Aerobic growth: absent, present, better than anaerobic growth, poorer than anaerobic growth.

Anaerobic growth: present, absent.

TEMPERATURE RELATIONS

Growth in refrigerator (°C.): present, absent.
 Growth at room temperature (°C.): present, absent.
 Growth at 37° C.: present, absent.
 Growth at 50° C.: present, absent.

Pasteurization survival, 80°C. 10 minutes: Positive

ADDITIONAL TESTS

Casein Hydrolysis: Negative
 Fat Hydrolysis: Negative
 Gelatin Hydrolysis: Negative
 Starch Hydrolysis: Negative
 Urea Hydrolysis: Negative

NH₄ from Peptone: Negative

Metabolism: Fermentative, aerogenic

Acetyl methyl carbinol: Positive

NH₄ as sole Nitrogen source: Positive

Sole Carbon sources: Citrate - Negative
 Glucose - Positive
 Sucrose - Positive
 Xylose - Positive

Methylene blue reduction: Positive

Crystalline Dextrins: Not observed

Salt tolerances: 2% - Negative
 7% - Negative
 10% - Negative

Reduces Brom Cresol Purple and Brom Thymol Blue.

Cultural Characteristics of Bacterial ColoniesCulture No. 118dSyptreas soy agar
10 days.I. Surface Colonies

a. Macroscopic appearance.



1. Size, mm 1mm.2. Shape: Outline- round, oval, irregular, filamentous.Elevation- effuse, flat, raised, convex, rugose, papillate, umbonate, pulvinate.Topography- smooth, rough, wrinkled, contoured, striated. Habit- compact, spreading.3. Optical properties:(a) Color: Color Harmony Manual No. 125a Cegam(b) Appearance by reflected light- dull, opalescent, iridescent, glistening, fluorescent.(c) Appearance by transmitted light- transparent, translucent, opaque.

b. Microscopic appearance (X100).

1. Margin- entire, granular, cleft, lobed, undulate, crenate, erose, ciliate, rhizoid.2. Internal structure- amorphous, dense, granular (fine, coarse), filamentous, striated, interlaced.c. Consistency- moist, slimy, soft, butyrous, waxy, tough, adherent, brittle.II. Deep coloniesa. Size, mmb. Shapec. Color

Name of organism Bacillus circulans # Studied by Dr. W. B. Bollen Culture No. 122A

Source Habitat Soil Date 4-8-65

Descriptions (Underline required terms.)	Sketches
<p>CELL MORPHOLOGY Medium: Trypticase soy agar Temp. 28 °C. Vegetative cells: Age: 24 hr. Form and arrangement: streptococci, diplococci, micrococci, sarcines, rods, commas, spirals, branched rods, filaments. Motility in broth: + Flagella: Peritrichous Size: 0.65 x 2.0μ Irregular forms: Sporangia: none, rods, spindles, elliptical, clavate, drumstick. Age: Endospores: Shape: spherical, ellipsoid, cylindrical. Position: central to eccentric, terminal, subterminal.</p>	 <i>Nigrospore</i>
<p>STAINING CHARACTERISTICS Gram: + -- Age: 18 hr. Method: Kopeloff's (modified) Special stains:</p>	 <i>Gram stain</i>
<p>AGAR STROKE Age: 18 hr. Amount of growth: scanty, moderate, abundant. Form: filiform, echinulate, beaded, spreading, rhizoid. Consistency: buoyant, viscid, membranous, brittle. Chromogenesis: White; fluorescent, iridescent, photogenic.</p>	 <i>24 hr</i>
<p>AGAR COLONIES Age: 2 day Form: punctiform, circular, filamentous, rhizoid, irregular. Elevation: effuse, flat, raised, convex. Surface: smooth, contoured, radiate, concentric, rugose. Margin: entire, undulate, erose, filamentous, curled. Density: opaque, translucent.</p>	 <i>48 hr</i>
<p>NUTRIENT BROTH Age: 8 day Surface growth: none, ring, pellicle, flocculent, membranous. Subsurface growth: none, turbid, granular. Amount of growth: scanty, moderate, abundant. Sediment: none, granular, flocculent, viscid, flaky.</p>	 <i>Temp. 25°C.</i>
<p>GELATIN STAB Age: 2 day Liquefaction: none, crateriform, infundibuliform, napiform, saccate, stratiform. Rate: slow, moderate, rapid.</p>	 <i>Temp. 25°C.</i>
<p>OTHER MEDIA Age: Potato slant: 3ba Pearl Soybean agar slant: White Glucose nitrate slant: No pigment Tyrosine agar slant: Covert Tan Fat agar: 2gc</p>	 <i>Temp. 25°C.</i>

FERMENTATION		Temp. 25 °C.			
Medium: Nutrient broth	1%	Glucose	Lactose	Sucrose	Xylose
Carbohydrate: 1%		+	+	+	+
Indicator: BCP					
Acid in 2 days		+	+	+	+
Acid in days					
Gas in 2 days		-	-	-	-
Gas in days					

ACTION ON MILK		Temp. 28 °C.	
Indicator:	Days		
Litmus	2 13		
Reaction	+		
Acid curd		+ (coag & whey sep.)	
Rennet curd			
Peptonization			
Reduction (before coagulation)	+		

* small variant.

ACTION ON NITRATES

Medium: 1% KNO_3 broth Temp. 28°C.
 Nitrite: - 3 d. ; d. ; d.
 Gas (N): - 3 d. ; d. ; d.

INDOLE PRODUCTION

Medium: Tryptophan broth Age: 10 day
 Method: Kovac's
 Indole: present, absent.

HYDROGEN SULFIDE PRODUCTION

Medium: Lead acetate & H₂S: Age: 10 day
 Thio iron Temp. 28°C.

RELATION TO FREE OXYGEN - Catalase: Positive
 Medium: Dextrose Nutrient a Age: 10 day
 Method: Shake tubes Temp. 28°C.

TEMPERATURE RELATIONS

Growth in refrigerator (°C): present, absent.
 Growth at room temperature (°C): present, absent.
 Growth at 37° C.: present, absent.
 Growth at 50° C.: present, absent.

Aerobic growth: absent, present, better than anaerobic growth, poorer than anaerobic growth.
 Anaerobic growth: present, absent.

Pasteurization survival, 80°C. 10 minutes: Positive

ADDITIONAL TESTS

Casein Hydrolysis: Positive
 Fat Hydrolysis: Positive
 Gelatin Hydrolysis: Positive
 Starch Hydrolysis: Positive
 Urea Hydrolysis: Negative

NH₄ from Peptone: Negative

Metabolism: Fermenter - anaerogenic

Acetyl methyl carbinol: Negative

NH₄ as sole Nitrogen source: Positive

Sole Carbon sources: Citrate - Positive
 Glucose - Negative
 Sucrose - Positive
 Xylose - Positive

Methylene blue reduction: Negative

Crystalline Dextrins: Observed

Salt tolerances: 2% - Positive
 7% - Negative
 10% - Negative

Cultural Characteristics of Bacterial Colonies

Culture No. 122A

Soybean soy agar
2 days.I. Surface Colonies

a. Macroscopic appearance.

1. Size, mm 1mm.
2. Shape: Outline- round, oval, irregular, filamentous.
Elevation- effuse, flat, raised, convex, rugose, papillate, umbonate, pulvinate.
Topography- smooth, rough, wrinkled, contoured, striated. Habit- compact, spreading.
3. Optical properties:
 - (a) Color: Color Harmony Manual No. 1½ ca (cream)
 - (b) Appearance by reflected light- dull, opalescent, iridescent, glistening, fluorescent.
 - (c) Appearance by transmitted light- transparent, translucent, opaque.

b. Microscopic appearance (X100).

1. Margin- entire, granular, cleft, lobed, undulate, eremate, erose, ciliate, rhizoid.
2. Internal structure- amorphous, dense, granular (fine, coarse), filamentous, striated, interlaced.
3. Consistency- moist, slimy, soft, butyrous, waxy, tough, adherent, brittle.



II. Deep colonies

- a. Size, mm
- b. Shape
- c. Color

Name of organism *Bacillus circulans* Studied by Dr. W.B. Bollen Culture No. 126

Source Habitat Soil Date

Descriptions (Underline required terms.)		Sketches
CELL MORPHOLOGY Medium: <u>Trypticase soy</u> Vegetative cells: Age: 24 hr. Form and arrangement: <u>streptococci, diplococci, micrococci, sarcines, rods, commas, spirals, branched rods, filaments.</u> Motility in broth: + Flagella: <u>Peritrichous</u> Size: <u>1.5 - 2.0 μ</u> Irregular forms: Sporangia: none, rods spindles, elliptical, elevate, drumstick. Age: Endospores: Shape: <u>spherical, ellipsoid, cylindrical.</u> Position: <u>central to eccentric, terminal, subterminal.</u>	Temp. 28 °C.	

STAINING CHARACTERISTICS Gram: + - - Age: 18 hr. Method: Kopeloff (modified) Special stains:	Temp. 25 °C.	
AGAR STROKE Age: 18 hrs. Amount of growth: <u>scanty, moderate, abundant.</u> Form: <u>filiform, echinulate, beaded, spreading, rhizoid.</u> Consistency: <u>butyrous, viscid, membranous, brittle.</u> Chromogenes: White : <u>fluorescent, iridescent, photogenic.</u>	Temp. 25 °C.	

AGAR COLONIES Age:	Temp. °C.	
Form: <u>punctiform, circular, filamentous, rhizoid, irregular.</u> Elevation: <u>effuse, flat, raised, convex.</u> Surface: <u>smooth, contoured, radiate, concentric, rugose.</u> Margin: <u>entire, undulate, erose, filamentous, curled.</u> Density: <u>opaque, translucent.</u>		

NUTRIENT BROTH Age: 8 days	Temp. 25°C.	
Surface growth: <u>none, ring, pellicle, flocculent, membranous.</u> Subsurface growth: <u>none, turbid, granular.</u> Amount of growth: <u>scanty, moderate, abundant.</u> Sediment: <u>none, granular, flocculent, viscid, flaky.</u>		

GELATIN STAB Age:	Temp. °C.	
Liquefaction: <u>none, crateriform, infundibuliform, napiform, saccate, striform.</u> Rate: <u>slow, moderate, rapid.</u>		

OTHER MEDIA Age:	Temp. °C.	
Potato slant: 31g	Nude tan	
Soybean agar slant:	Grey-white	(scant)
Glucose nitrate slant:	Scant	
Fat agar: 3nl	Dark Brown	(abundant) medium darkened
Tyrosine agar slant:	No pigment	

Glucose nutrient agar slant: growth as good or better than nutrient agar

FERMENTATION		Temp. 28 °C.			
Medium: Nutrient broth	Glucose	Lactose	Sucrose	Galactose	
Carbohydrate: 1 %					
Indicator: BCP					
Acid in 2 days	+	+	+	+	
Acid in days					
Gas in 2 days	-	-	-	-	
Gas in days					

Indicator:	ACTION ON MILK		Temp. 28 °C.	agar
	Days	Days		
Litmus	1	20		
Reaction	ACID			
Acid curd		+ (coag. & whey sep.)		
Rennet curd				
Peptonization				
Reduction (before coagulation)	+			

ACTION ON NITRATES

Medium: 1% KNO₃ broth Temp. 28°C.
 Nitrite: -3 d. ; d.
 Gas (N): -3 d. ; d.

INDOLE PRODUCTION

Medium: Tryptophan broth Age: 10 days
 Method: Kovac's
 Indole: present, absent.

HYDROGEN SULFIDE PRODUCTION

Medium: Lead Acetate &
 Thio-iron
 H.S.: present, absent.

Age: 10 days
 Temp. 28 °C.

RELATION TO FREE OXYGEN - Catalase: Positive - Neg
 Medium: Dextrose-Nutrient A. Age: 10 days
 Method: Shake Tubes Temp. 28°C.

TEMPERATURE RELATIONS

Growth in refrigerator (10°C.): present, absent. (+-)
 Growth at room temperature (28°C.): present, absent.
 Growth at 37°C.: present, absent.
 Growth at 50°C.: present, absent.

Aerobic growth: absent, present, better than anaerobic growth, poorer
 than anaerobic growth.
 Anaerobic growth: present, absent.

Pasteurization survival, 80°C. 10 minutes: Positive

ADDITIONAL TESTS

Casein Hydrolysis: Positive
 Fat Hydrolysis: Negative
 Gelatin Hydrolysis: Positive
 Starch Hydrolysis: Positive
 Urea Hydrolysis: Negative

NH₄ from Peptone: Negative

Metabolism: Fermenter - anaerogenic

Acetyl methyl carbinol: Negative

NH₄ as sole Nitrogen source: Positive

Sole Carbon sources: Citrate - Negative
 Glucose - Positive
 Sucrose - Positive
 Xylose - Positive

Methylene blue reduction: Positive 1 day: Negative 5 days.

Crystalline Dextrins: Not observed

Salt tolerances: 2% - Positive
 7% - Positive
 10% - Negative

Cultural Characteristics of Bacterial Colonies

Culture No. 126

*Syntexis soy agar
2 days.*

I. Surface Colonies

a. Macroscopic appearance.

1. Size, mm 1 mm

2. Shape: Outline- round, oval, irregular, filamentous.

Elevation- effuse, flat, raised, convex, rugose, papillate, umbonate, pulvinate.

Topography- smooth, rough, wrinkled, contoured, striated. Habit- compact, spreading.

3. Optical properties:

(a) Color: Color Harmony Manual No. 1ea (pale yellow)

(b) Appearance by reflected light- dull, opalescent, iridescent, glistening, fluorescent.

(c) Appearance by transmitted light- transparent, translucent, opaque.

b. Microscopic appearance (X100).

1. Margin- entire, granular, cleft, lobed, undulate, cremate, erose, ciliate, rhizoid.

2. Internal structure- amorphous, dense, granular (fine, coarse), filamentous, striated, interlaced.

c. Consistency- moist, slimy, soft, butyrous, waxy, tough, adherent, brittle.

II. Deep colonies

a. Size, mm

b. Shape

c. Color