

NASA CONTRACTOR  
REPORT

NASA CR-61240

April 29, 1968

NASA CR-61240

GPO PRICE \$ \_\_\_\_\_

CFSTI PRICE(S) \$ \_\_\_\_\_

Hard copy (HC) 3.00

Microfiche (MF) .65

ff 653 July 65

EVALUATION OF ANALYTICAL STANDARDS BY DIFFERENTIAL  
THERMAL ANALYSIS AND DIFFERENTIAL  
SCANNING CALORIMETRY

Prepared under Contract No. NAS 8-20073 by  
J. P. Evans and K. G. Scroggins

BROWN ENGINEERING COMPANY

N 68-35618

FACILITY FORM 602

(ACCESSION NUMBER)  
68  
(PAGES)  
CR-61240  
(NASA CR OR TMX OR AD NUMBER)

(THRU)  
1  
(CODE)  
06  
(CATEGORY)



For

NASA-GEORGE C. MARSHALL SPACE FLIGHT CENTER  
Huntsville, Alabama

April 29, 1968

NASA CR-61240

EVALUATION OF ANALYTICAL STANDARDS BY  
DIFFERENTIAL THERMAL ANALYSIS AND  
DIFFERENTIAL SCANNING CALORIMETRY

By

J. P. Evans and K. G. Scrogham

Prepared under Contract No. NAS 8-20073 by

BROWN ENGINEERING COMPANY  
Huntsville, Alabama

For

Propulsion and Vehicle Engineering Laboratory

Distribution of this report is provided in the interest of  
information exchange. Responsibility for the contents  
resides in the author or organization that prepared it.

NASA-GEORGE C. MARSHALL SPACE FLIGHT CENTER

TECHNICAL REPORT BSVD-M-68-TR-001

EVALUATION OF ANALYTICAL STANDARDS BY  
DIFFERENTIAL THERMAL ANALYSIS AND  
DIFFERENTIAL SCANNING CALORIMETRY

BY

J. P. Evans and K. G. Scrogham  
Thermal Measurements Section  
Chemistry Branch  
Materials Department  
Space Vehicle Division  
Brown Engineering Company, Inc.

ABSTRACT

The effects of heating rate and furnace atmosphere on four analytical standards were investigated by DTA and DSC. The results are presented in both graphical and tabular form.

## TABLE OF CONTENTS

	Page
SUMMARY.....	1
INTRODUCTION.....	1
EXPERIMENTAL.....	2
Differential Thermal Analysis.....	2
Differential Scanning Calorimetry .....	4
Heats of Fusion.....	6
RESULTS AND DISCUSSION .....	8
DTA.....	8
DSC .....	13
Heats of Fusion.....	13
CONCLUSIONS.....	19
APPENDIX A .....	20
List of Thermograms (DTA) .....	20
APPENDIX B .....	45
List of Thermograms (DSC).....	45

## LIST OF TABLES

Table	Title	Page
I	DTA of Copper Sulphate Pentahydrate	9
II	DTA of Calcium Oxalate Monohydrate	10
III	DTA of Potassium Nitrate	11
IV	DTA of Silver Nitrate	12
V	DSC of Copper Sulphate Pentahydrate	14
VI	DSC of Potassium Nitrate	15
VII	DSC of Silver Nitrate	16
VIII	Heats of Inversion and Fusion	17

# EVALUATION OF ANALYTICAL STANDARDS BY DIFFERENTIAL THERMAL ANALYSIS AND DIFFERENTIAL SCANNING CALORIMETRY

BY

J. P. Evans and K. G. Scrogham  
Brown Engineering Company, Inc.

## SUMMARY

Differential thermal analyses have been run on copper sulfate pentahydrate, calcium oxalate monohydrate, potassium nitrate and silver nitrate under the conditions requested. Thermograms are included for each of the analyses.

DSC runs were also made on the copper, potassium, and silver salts. Heats of inversion and fusion were determined for potassium nitrate and silver nitrate and comparison made between experimental and literature values.

## INTRODUCTION

Several factors are known to affect differential thermal analysis results. These factors are dependent on two types of variables, instrumental and sample characteristics. The effects of two of the instrumental factors, heating rate and furnace atmosphere, on four inorganic salts were investigated. Copper sulphate pentahydrate, calcium oxalate monohydrate, potassium nitrate, and silver nitrate were analyzed on the DuPont 900 Differential Thermal Analyzer using standard procedures.

These same salts with the exception of calcium oxalate monohydrate, were also analyzed on the Perkin-Elmer Differential Scanning Calorimeter. DTA results were compared with the DSC results.

The heats of fusion and inversion for potassium nitrate and the heat of fusion of silver nitrate were determined at a heating rate of 5°C./min. using the Perkin-Elmer Differential Scanning Calorimeter. These data were compared with available literature values.

The standard analytical compounds used in this study were supplied by Sadtler Research Laboratories, Philadelphia, Pennsylvania.

## EXPERIMENTAL

### Differential Thermal Analysis

The DuPont 900 Differential Thermal Analyzer was used for all DTA experiments. The instrument was equipped with the Standard DTA Cell containing the heating block for micro (2 mm) sample tubes.

In all analyses the following experimental procedure was followed:

- 1) The sample tube was tared and the sample introduced.
- 2) The sample and sample tube were weighed ( $\pm 0.05$  mg).
- 3) The sample tube was then tapped until the sample occupied all visible voids.
- 4) The thermocouple was inserted into the center of the sample and the sample tube was placed in the heating block of the Cell.
- 5) The procedure for preparation of the reference material ( $\text{Al}_2\text{O}_3$ ) was identical to that above.
- 6) For those runs requiring a dynamic atmosphere, the Cell was purged for at least fifteen minutes with dry nitrogen gas prior to the start of the analysis.

The specifics for the analysis of each material are outlined below.

1) Copper sulfate pentahydrate

a) Run at 5, 10, and 15°C. /minute rate of rise from ambient to 500°C. in static atmosphere (air).

b) Run at 5, 10, and 15°C. /minute rate of rise from ambient to 500°C. in dynamic atmosphere (N<sub>2</sub> at 1 SCFH for envelope).

c) All sample weights were 1.15 mg  $\pm$  0.05 mg.

2) Calcium oxalate monohydrate

a) Run at 5, 10, and 15°C. /minute rate of rise from ambient to 500°C. in dynamic atmosphere (N<sub>2</sub> at 1 SCFH for envelope).

b) Run at 5, 10, and 15°C. /minute rate of rise from ambient to 500°C. in dynamic atmosphere (air at 1 SCFH for envelope).

c) All sample weights were 1.60 mg  $\pm$  0.15 mg.

3) Potassium nitrate

a) Run at 5°C. /minute rate of rise from ambient to 350°C. in static atmosphere (air). Sample was allowed to cool to ambient temperature with thermocouple in place and rerun under identical experimental conditions.



b) Procedure a) was repeated with a new sample for heating rates of 10 and 15°C. /minute.

c) All sample weights were 1.65 mg.

4) Silver nitrate

a) Run at 5°C. /minute rate of rise from ambient to 250°C. in static atmosphere (air). Sample was allowed to cool to ambient temperature with thermocouple in place and rerun under identical experimental conditions.

b) Procedure a) was repeated with a new sample for heating rates of 10 and 15°C. /minute.

c) All sample weights were 3.40 mg  $\pm$  0.05 mg.

Differential Scanning Calorimetry

The Perkin-Elmer DSC-1 Differential Scanning Calorimeter was used for all DSC experiments. The instrument was calibrated at each heating rate using indium, tin, and lead. The melting points of these materials are 429°K., 505°K., and 600°K. respectively. Calibration was carried out to within 1° of these temperatures.

In all analyses the following experimental procedure was followed:

- 1) The sample pan was tared and the sample introduced.
- 2) The sample and sample pan were weighed ( $\pm 0.05$  mg).
- 3) The pan lid was crimped into place as per the manufacturer's instructions.
- 4) The reference side of the cell contained an empty pan and lid for all runs.
- 5) For those runs requiring a dynamic atmosphere, the cell was purged for at least fifteen minutes with dry nitrogen gas prior to the start of the analysis.

The specifics for the analysis of each material are outlined below.

1) Copper sulfate pentahydrate

- a) Run at 5, 10, and 20°C. /minute rate of rise from ambient to 500°C. in static atmosphere (air).
- b) Run at 5, 10, and 20°C. /minute rate of rise from ambient to 500°C. in dynamic atmosphere (N<sub>2</sub> at 30 cc/minute for envelope).
- c) All sample weights were 3.5 mg  $\pm$  0.1 mg.

2) Potassium nitrate

- a) Run at 5°C. /minute rate of rise from ambient to 350°C. in static atmosphere (air). Sample was allowed to cool to ambient temperature in the cell and rerun under identical experimental conditions.

b) Procedure a) was repeated with a new sample for heating rates of 10 and 20°C./minute.

c) All sample weights were 3.5 mg  $\pm$  0.1 mg.

3) Silver nitrate

a) Run at 5°C./minute rate of rise from ambient to 250°C. in static atmosphere (air). Sample was allowed to cool to ambient temperature in the cell and rerun under identical experimental conditions.

b) Procedure a) was repeated with a new sample for heating rates of 10 and 20°C./minute.

c) All sample weights were 3.5 mg  $\pm$  0.1 mg.

Heats of Fusion

The heats of inversion and fusion for potassium nitrate were determined along with the heat of fusion of silver nitrate. These determinations were made in duplicate but only at a heating rate of 5°C./minute. Lack of "standards" prevented further investigation at higher rates. Hermetically sealed sample pans were used in these experiments to prevent volatilization. Samples were weighed before and after analysis to assure constant weight. Range setting (4x) and chart speed (15 or 30 cm/hr.) were selected to give peak areas in the range 25 to 85 cm<sup>2</sup>. All areas were determined using a K & E Compensating Polar

Planimeter, each area being measured at least twice.  
The arithmetic average of all measurements on each peak was taken as the peak area.

Calibration of the power input to the sample was performed by measurements of the heat of fusion of indium ( $\Delta H_f = 6.80 \text{ cal/gm}$ ).

## RESULTS AND DISCUSSION

### DTA

The results of all DTA experiments are tabulated in Tables I through IV and shown on pages 21 through 44 in Appendix A. Analysis of the data on copper sulfate pentahydrate shows the dehydration to be insensitive to either atmosphere or heating rate. This is not the case, however, for the decomposition reaction. An increase in heating rate causes an elevation in peak temperature in both air and nitrogen atmospheres.

Analysis of the calcium oxalate monohydrate data shows the peak temperatures to be relatively insensitive to either atmosphere or heating rate.

Potassium nitrate was run in air and then rerun at each heating rate. The inversion temperature and the melting temperature are insensitive to heating rate and reproducible from run to run.

Data analysis shows this is not the case with silver nitrate. The inversion is permanent and non-reversible and therefore is not observed upon rerunning the salt. The fusion temperature, however, is found to be constant.

TABLE I

DTA of Copper Sulfate Pentahydrate

Run No.	Wt. (mg.)	Rate (°C. /Min.)	Atmosphere		Temperature °C. <sup>a</sup>	
			Air	N <sub>2</sub>	Onset	Peak
SRL-1	1.20	5	X		96 101 112 239	98 104 123 256
SRL-2	1.15	10	X		94 100 116 244	97 104 126 261
SRL-3	1.10	15	X		96 102 120 250	98 105 126 266
SRL-4	1.20	5		X	96 98 109 220	98 103 123 240
SRL-5	1.20	10		X	96 100 113 226	98 104 124 247
SRL-6	1.15	15		X	96 103 116 233	99 108 129 253

<sup>a</sup>  $\Delta T$  Sensitivity, 0.5°C. /in.

TABLE II

DTA of Calcium Oxalate Monohydrate

Run No.	Wt. (mg. )	Rate (°C. /Min. )	Atmosphere		Temperature °C. <sup>a</sup>	
			Air	N <sub>2</sub>	Onset	Peak
SRL-10	1. 70	5	X		224	239
SRL-11	1. 55	10	X		226	240
SRL-12	1. 60	15	X		234	244
SRL-7	1. 45	5		X	216	236
SRL-8	1. 45	10		X	232	243
SRL-9	1. 45	15		X	233	244

<sup>a</sup>  $\Delta T$  Sensitivity, 0.5°C. /min.

TABLE III

DTA of Potassium Nitrate<sup>a</sup>

Run No.	Wt. (mg. )	Rate (°C. /Min. )	First Run	Rerun	Temperature °C. <sup>b</sup>	
					Onset	Peak
SRL-13	1.65	5	X		128 331	131 333
SRL-14	1.65	5		X	128 331	129 333
SRL-15	1.65	10	X		129 331	133 334
SRL-16	1.65	10		X	129 331	130 334
SRL-17	1.65	15	X		129 331	133 334
SRL-18	1.65	15		X	129 331	130 334

<sup>a</sup> Atmosphere, Static Air

<sup>b</sup>  $\Delta T$  Sensitivity, 0.5°C. /in.



TABLE IV

DTA of Silver Nitrate<sup>a</sup>

Run No.	Wt. (mg. )	Rate (°C. /Min. )	First Run	Rerun	Temperature °C. <sup>b</sup>	
					Onset	Peak
SRL-19	3. 50	5	X		164 207	166 209
SRL-20	3. 40	5		X	--- 207	--- 209
SRL-21	3. 40	10	X		164 207	167 209
SRL-22	3. 40	10		X	--- 206	--- 208
SRL-23	3. 45	15	X		165 208	167 210
SRL-24	3. 45	15		X	--- 207	--- 209

<sup>a</sup> Atmosphere, Static Air

<sup>b</sup>  $\Delta T$  Sensitivity, 0.5°C. /min.

## DSC

The results of all DSC experiments are summarized in Tables V through VII and shown on pages 46 through 63 in Appendix B. Analysis of the data on copper sulfate pentahydrate shows a variation in the temperatures at which the dehydration is observed. According to the workers at Perkin-Elmer<sup>a</sup> this is to be expected since the peaks may occur at any temperature within the range of stability of the respective hydrates.

The analyses of potassium nitrate are very reproducible as shown in the table. No significant deviation is observed in either heating rate or repeated analysis.

The excellent reproducibility is also observed in the silver nitrate data. Again it is shown that the inversion is permanent and non-reversible.

## Heats of Fusion

The heat of fusion was determined for both  $\text{KNO}_3$  and  $\text{AgNO}_3$  and is shown in Table VIII. The heat associated with the inversion of both salts was also determined.

<sup>a</sup> Thermal Analysis Newsletter No. 7, 1967,  
Perkin-Elmer Corporation, Norwalk, Connecticut

TABLE V

DSC of Copper Sulfate Pentahydrate

Run No.	Wt. (mg.)	Rate (°C/Min.)	$\Delta T$ Sensitivity (Mcal./Sec.)	Atmosphere		Peak Temperature, °C.
				Air	N <sub>2</sub>	
259	3.5	5	4	X		79 105 248
266	3.5	10	8	X		77 107 247
270	3.6	20	8	X		92 119 256
258	3.5	5	4		X	75 100 221
265	3.5	10	8		X	78 107 238
269	3.5	20	8		X	92 118 243

TABLE VI  
DSC of Potassium Nitrate<sup>a</sup>

Run No.	Wt. (mg.)	Rate (°C. /Min.)	$\Delta T$ Sensitivity (Mcal. /Sec.)	First Run	Rerun	Peak Temperature, °C.
260A	3.5	5	8	X		133 334
260B	3.5	5	8		X	128 334
267A	3.4	10	16	X		130 336
267B	3.4	10	16		X	129 335
271A	3.6	20	16	X		133 335
271B	3.6	20	16		X	130 336

<sup>a</sup> Atmosphere, Static Air

TABLE VII  
DSC of Silver Nitrate<sup>a</sup>

Run No.	Wt. (mg.)	Rate (°C./Min.)	$\Delta T$ Sensitivity (Mcal./Sec.)	First Run	Rerun	Peak Temperature, °C.
262A	3.4	5	4	X		171 211
262B	3.4	5	4		X	--- 212
268A	3.6	10	8	X		169 211
268B	3.6	10	8		X	--- 212
272A	3.5	20	8	X		172 211
272B	3.5	20	8		X	--- 212

<sup>a</sup> Atmosphere, Static Air

TABLE VIII  
Heats of Inversion and Fusion

Run No.	Sample	$\Delta H_f$ (Cal/gm)	$T_f$ °C.	Literature $\Delta H_f$ T	
261	$\text{KNO}_3^a$	10.5	130	---	---
		11.0	129	---	---
	$\text{KNO}_3$	19.9	336	25.4 <sup>b</sup>	308
		19.6	335	27.7 <sup>c</sup>	337
263	$\text{AgNO}_3^a$	2.3	168	---	---
264		2.8	169	---	---
263	$\text{AgNO}_3$	15.1	211	16.7 <sup>d</sup>	212
		14.5	210	17.7 <sup>b</sup>	208
				16.2 <sup>c</sup>	210

<sup>a</sup> Inversion

<sup>b</sup> Handbook of Chemistry and Physics, 40th. Edition, 1959

<sup>c</sup> Handbook of Chemistry, Lange, 9th. Edition

<sup>d</sup> Handbook of Differential Thermal Analysis, Smothers and Chang, 1966

Analysis of the data on fusion shows a reproducibility of 4% or better and a variation from available literature sources on the order of 20%. No data was available for comparison of the heats of inversion.

## CONCLUSIONS

Differential thermal analysis has shown that the dehydration of copper sulphate pentahydrate is not affected by either heating rate or atmosphere. However, the decomposition reaction was affected by these two variables. The dehydration of calcium oxalate monohydrate was also shown to be relatively insensitive to either atmosphere or heating rate. The inversion and the melting of potassium nitrate were insensitive to both variables as was the melting of silver nitrate. However, the silver nitrate inversion was found to be permanent and non-reversible.

Differential scanning calorimetry results displayed a variation in the dehydration temperatures of copper sulphate pentahydrate as should be expected. The analyses of potassium nitrate and silver nitrate by DSC were very reproducible in the range of rates investigated.


The heats of fusion determined by DSC for potassium nitrate and silver nitrate had a reproducibility of 4% and a variation of 20% from available literature sources.




## Appendix A

### List of Thermograms (DTA)


	<u>Page</u>
Copper Sulphate Pentahydrate in Air at 5°C. /min.	21
Copper Sulphate Pentahydrate in Air at 10°C. /min.	22
Copper Sulphate Pentahydrate in Air at 15°C. /min.	23
Copper Sulphate Pentahydrate in N <sub>2</sub> at 5°C. /min.	24
Copper Sulphate Pentahydrate in N <sub>2</sub> at 10°C. /min.	25
Copper Sulphate Pentahydrate in N <sub>2</sub> at 15°C. /min.	26
Calcium Oxalate Monohydrate in N <sub>2</sub> at 5°C. /min.	27
Calcium Oxalate Monohydrate in N <sub>2</sub> at 10°C. /min.	28
Calcium Oxalate Monohydrate in N <sub>2</sub> at 15°C. /min.	29
Calcium Oxalate Monohydrate in Air at 5°C. /min.	30
Calcium Oxalate Monohydrate in Air at 10°C. /min.	31
Calcium Oxalate Monohydrate in Air at 15°C. /min.	32
Potassium Nitrate at 5°C. /min.	33
Potassium Nitrate Rerun at 5°C. /min.	34
Potassium Nitrate at 10°C. /min.	35
Potassium Nitrate Rerun at 10°C. /min.	36
Potassium Nitrate at 15°C. /min.	37
Potassium Nitrate Rerun at 15°C. /min.	38
Silver Nitrate at 5°C. /min.	39
Silver Nitrate Rerun at 5°C. /min.	40
Silver Nitrate at 10°C. /min.	41
Silver Nitrate Rerun at 10°C. /min.	42
Silver Nitrate at 15°C. /min.	43
Silver Nitrate Rerun at 15°C. /min.	44




<b>SAMPLE:</b> $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$  <b>ORIGIN:</b> Sadtler Research Laboratories	<b>SIZE</b> 1.15 mg in microtube <b>REF.</b> $\text{Al}_2\text{O}_3$ <b>PROGRAM MODE</b> Heat <b>RATE</b> 10 $\frac{^\circ\text{C}}{\text{min}}$ <b>START</b> 26 $^\circ\text{C}$	<b>ATM.</b> Air, Static <b>T</b> 50 $\frac{^\circ\text{C}}{\text{min}}$ <b>SCALE SETTING</b> 0.5 $\frac{^\circ\text{C}}{\text{mm}}$ Baseline Slope -0.2	<b>RUN NO.</b> SRL-2 <b>DATE</b> 2/1/68 <b>OPERATOR</b> J. Evans




<b>SAMPLE:</b> $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$  <b>ORIGIN:</b> Sadtler Research Laboratories	<b>SIZE</b> 1.10 mg in microtube <b>REF.</b> $\text{Al}_2\text{O}_3$ <b>PROGRAM MODE</b> Heat <b>RATE</b> 15 $\frac{\%}{\text{min}}$ , <b>START</b> 24 $^{\circ}\text{C}$	<b>ATM.</b> Air, Static <b>T</b> $^{\circ}\text{C}$ <b>SCALE</b> 50 $\frac{\%}{\text{min}}$ <b>SETTING</b> Baseline Slope -0.2	<b>RUN NO.</b> SRL-3 <b>DATE</b> 2/1/68 <b>OPERATOR</b> J. Evans




<b>SAMPLE:</b> $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$  <b>ORIGIN:</b> Sadtler Research Laboratories	<b>SIZE</b> 1.20 mg. in microtube <b>REF.</b> $\text{Al}_2\text{O}_3$ <b>PROGRAM MODE</b> Heat <b>RATE</b> 5 $\frac{^\circ\text{C}}{\text{min}}$ <b>START</b> 24 $^\circ\text{C}$	<b>ATM.</b> $\text{N}_2$ , 1 SCFH <table border="1"> <tr> <td><b>T</b></td> <td><b><math>\Delta T</math></b></td> </tr> <tr> <td>50 <math>\frac{^\circ\text{C}}{\text{min}}</math></td> <td>0.5 <math>\frac{^\circ\text{C}}{\text{min}}</math></td> </tr> </table>	<b>T</b>	<b><math>\Delta T</math></b>	50 $\frac{^\circ\text{C}}{\text{min}}$	0.5 $\frac{^\circ\text{C}}{\text{min}}$	<b>RUN NO.</b> SRL-4 <b>DATE</b> 2/23/68 <b>OPERATOR</b> J. Evans
	<b>T</b>	<b><math>\Delta T</math></b>					
	50 $\frac{^\circ\text{C}}{\text{min}}$	0.5 $\frac{^\circ\text{C}}{\text{min}}$					
	<b>SCALE SETTING</b> Baseline Slope -0.2						




T. °C (CORRECTED FOR CHROMEL ALUMEL THERMOCOUPLES)







<b>SAMPLE:</b> $\text{CaC}_2\text{O}_4$ $\text{H}_2\text{O}$  <b>ORIGIN:</b> Sadtler Research Laboratories	<b>SIZE</b> 1.45 mg in microtube <b>REF.</b> 1.85 mg. $\text{AL}_2\text{O}_3$ <b>PROGRAM MODE</b> Heat <b>RATE</b> 5 $\frac{\%}{\text{min}}$ , <b>START</b> 24 $^{\circ}\text{C}$	<b>ATM.</b> $\text{N}_2$ , 1 SCFH <b>T</b> $\Delta T$ <b>SCALE SETTING</b> 50 $\frac{\%}{\text{min}}$ 0.5 $\frac{\%}{\text{min}}$ Baseline Slope -0.2	<b>RUN NO.</b> SRL-7 <b>DATE</b> 3/15/68 <b>OPERATOR</b> J. Evans









<b>SAMPLE:</b> $\text{CaC}_2\text{O}_4$ $\text{H}_2\text{O}$  <b>ORIGIN:</b> Sadtler Research Laboratories	<b>SIZE</b> <u>1.45 mg. in microtube</u> <b>REF.</b> <u>1.85 mg. <math>\text{Al}_2\text{O}_3</math></u>	<b>ATM.</b> <u><math>\text{N}_2</math></u> <u>1 SCFH</u> <b>T</b> <u>50</u> <b><math>\Delta T</math></b> <u>0.5</u>	<b>RUN NO.</b> <u>SRL-9</u> <b>DATE</b> <u>3/15/68</u>
	<b>PROGRAM MODE</b> <u>Heat</u> <b>RATE</b> <u>15</u> <b><math>\frac{\text{°C}}{\text{min}}</math></b> <b>START</b> <u>25</u> <b><math>\text{°C}</math></b>	<b>SCALE SETTING</b> <u>0.5</u> Baseline Slope <u>-0.2</u>	<b>OPERATOR</b> <u>J. Evans</u>




<b>SAMPLE:</b> $\text{CaC}_2\text{O}_4 \cdot \text{H}_2\text{O}$  <b>ORIGIN:</b> Sadtler Research Laboratories	<b>SIZE</b> 1.70 mg. in microtube <b>REF.</b> 1.85 mg. $\text{AL}_2\text{O}_3$	<b>ATM.</b> Air, 1 SCFH <div style="display: flex; justify-content: space-around;"> <span>T</span> <span><math>\Delta T</math></span> </div>	<b>RUN NO.</b> SRL-10 <b>DATE</b> 3/19/68 <b>OPERATOR</b> J. Evans	
	<b>PROGRAM MODE</b> Heat <b>RATE</b> 5 $\frac{^\circ\text{C}}{\text{min}}$	<b>SCALE</b> 50 $\frac{^\circ\text{C}}{\text{mm}}$ <b>SETTING</b> 0.5 $\frac{^\circ\text{C}}{\text{mm}}$		
	<b>START</b> 25 $^\circ\text{C}$ Baseline Slope -0.2			




<b>SAMPLE:</b> $\text{CaC}_2\text{O}_4$ $\text{H}_2\text{O}$	<b>SIZE</b> 1.55 mg. in microtube	<b>ATM.</b> Air, 1 SCFH		<b>RUN NO.</b> SRL-11
	<b>REF.</b> 1.85 mg. $\text{AL}_2\text{O}_3$	<b>T</b>	<b><math>\Delta T</math></b>	<b>DATE</b> 3/18/68
	<b>PROGRAM MODE</b> Heat	<b>SCALE</b> 50 $\frac{\%}{\text{min}}$	<b>0.5 <math>\frac{\%}{\text{min}}</math></b>	<b>OPERATOR</b> J. Evans
	<b>RATE</b> 10 $\frac{\%}{\text{min}}$ , START 25 $^{\circ}\text{C}$	<b>SETTING</b> Baseline Slope -0.2		
<b>ORIGIN:</b> Sadtler Research Laboratories				




T,  $^{\circ}\text{C}$  (CORRECTED FOR CHROMEL ALUMEL THERMOCOUPLES)




<b>SAMPLE:</b> $\text{KNO}_3$	<b>SIZE</b> 1.65 mg	<b>ATM.</b> Air, Static		<b>RUN NO.</b> SRL-13
	<b>REF.</b> 1.85 mg. $\text{AL}_2\text{O}_3$	<b>T</b>	$\Delta T$	<b>DATE</b> 3-21-68
	<b>PROGRAM MODE</b> Heat	<b>SCALE</b> 50 $\frac{\mu\text{V}}{\text{mm}}$	0.5 $\frac{\mu\text{V}}{\text{mm}}$	<b>OPERATOR</b> J. Evans
	<b>RATE</b> 5 $\frac{\%}{\text{min}}$ , <b>START</b> 25 $^{\circ}\text{C}$	<b>SETTING</b> Baseline Slope -0.2		
<b>ORIGIN:</b> Sadtler Research Laboratories				



T,  $^{\circ}\text{C}$  (CORRECTED FOR CHROMEL ALUMEL THERMOCOUPLES)




<b>SAMPLE:</b> $\text{KNO}_3$	<b>SIZE</b> 1.65 mg.	<b>ATM.</b> Air, Static	<b>RUN NO.</b> SRL-15
	<b>REF.</b> 1.85 mg. $\text{Al}_2\text{O}_3$	<b>T</b> $\Delta T$	<b>DATE</b> 3-22-68
	<b>PROGRAM MODE</b> Heat	<b>SCALE</b> 50 $\frac{\text{mV}}{\text{m}}$ 0.5 $\frac{\text{mV}}{\text{m}}$	<b>OPERATOR</b> J. Evans
	<b>ORIGIN:</b> Sadtler Research Laboratories	<b>SETTING</b> Baseline Slope -0.2	
	<b>RATE</b> 10 $\frac{\text{mV}}{\text{min}}$ <b>START</b> 23 $^{\circ}\text{C}$		




T,  $^{\circ}\text{C}$  (CORRECTED FOR CHROMEL ALUMEL THERMOCOUPLES)




<b>SAMPLE:</b> $\text{KNO}_3$ Rerun	<b>SIZE</b> 1.65 mg. in microtube	<b>ATM.</b> Air, Static		<b>RUN NO.</b> SRL-16
	<b>REF.</b> 1.85 mg. $\text{AL}_2\text{O}_3$	<b>T</b>	<b><math>\Delta T</math></b>	<b>DATE</b> 3-22-68
	<b>PROGRAM MODE</b> Heat	<b>SCALE</b> 50 $\frac{\%}{\text{min}}$ 0.5 $\frac{\%}{\text{min}}$		<b>OPERATOR</b> J. Evans
	<b>ORIGIN:</b> Sadtler Research Laboratories	<b>RATE</b> 10 $\frac{\%}{\text{min}}$ , <b>START</b> 30 $^\circ\text{C}$		
<b>SETTING</b> Baseline Slope -0.2				




<b>SAMPLE:</b> $\text{KNO}_3$  <b>ORIGIN:</b> Sadtler Research Laboratories	<b>SIZE</b> 1.65 mg. <b>REF.</b> 1.85 mg. $\text{Al}_2\text{O}_3$ <b>PROGRAM MODE</b> Heat <b>RATE</b> 15 $\frac{^\circ\text{C}}{\text{min}}$ , <b>START</b> 23 $^\circ\text{C}$	<b>ATM.</b> Air, Static <b>T</b> $^\circ\text{C}$ $\Delta T$ <b>SCALE</b> 50 $\frac{^\circ\text{C}}{\text{mm}}$ 0.5 $\frac{^\circ\text{C}}{\text{mm}}$ <b>SETTING</b> Baseline Slope -0.2	<b>RUN NO.</b> SRL-17 <b>DATE</b> 3-22-68 <b>OPERATOR</b> J. Evans




T,  $^\circ\text{C}$  (CORRECTED FOR CHROMEL ALUMEL THERMOCOUPLES)




<b>SAMPLE:</b> $\text{AgNO}_3$  <b>ORIGIN:</b> Sadtler Research Laboratories	<b>SIZE</b> 3.40 mg. in microtube <b>REF.</b> 1.85 mg. $\text{AL}_2\text{O}_3$	<b>ATM.</b> Air, Static <div style="display: flex; justify-content: space-around;"> <span>T</span> <span><math>\Delta T</math></span> </div>	<b>RUN NO.</b> SRL-19 <b>DATE</b> 4/5/68 <b>OPERATOR</b> J. Evans
	<b>PROGRAM MODE</b> Heat <b>RATE</b> 5 $\frac{^\circ\text{C}}{\text{min}}$ , <b>START</b> 24 $^\circ\text{C}$	<b>SCALE</b> 50 $\frac{^\circ\text{C}}{\text{in}}$ <b>SETTING</b> 0.5 $\frac{^\circ\text{C}}{\text{in}}$ Baseline Slope -0.2	




<b>SAMPLE:</b> $\text{AgNO}_3$ Rerun  <b>ORIGIN:</b> Sadtler Research Laboratories	<b>SIZE</b> 3.40 mg. in microtube <b>REF.</b> 1.85 mg. $\text{Al}_2\text{O}_3$ <b>PROGRAM MODE</b> Heat <b>RATE</b> 5 $\frac{^\circ\text{C}}{\text{min}}$ , START 30 $^\circ\text{C}$	<b>ATM.</b> Air, Static <b>T</b> <b>SCALE</b> 50 $\frac{^\circ\text{C}}{\text{mm}}$ <b>SETTING</b> 0.5 $\frac{^\circ\text{C}}{\text{mm}}$ Baseline Slope -0.2	<b>RUN NO.</b> SRL-20 <b>DATE</b> 4/5/68 <b>OPERATOR</b> J. Evans
--	--	---	---






<b>SAMPLE:</b> AgNO <sub>3</sub> Rerun  <b>ORIGIN:</b> Sadtler Research Laboratories	<b>SIZE</b> 3.40 mg. in microtube <b>REF.</b> 1.85 mg. AL <sub>2</sub> O <sub>3</sub> <b>PROGRAM MODE</b> Heat <b>RATE</b> 10 $\frac{^\circ\text{C}}{\text{min}}$ , <b>START</b> 25 $^\circ\text{C}$	<b>ATM.</b> Air, Static <table border="1"> <tr> <td>T</td> <td><math>\Delta T</math></td> </tr> <tr> <td>50 <math>\frac{^\circ\text{C}}{\text{min}}</math></td> <td>0.5 <math>\frac{^\circ\text{C}}{\text{min}}</math></td> </tr> </table>	T	$\Delta T$	50 $\frac{^\circ\text{C}}{\text{min}}$	0.5 $\frac{^\circ\text{C}}{\text{min}}$	<b>RUN NO.</b> SRL-22 <b>DATE</b> 4/4/68 <b>OPERATOR</b> J. Evans
	T	$\Delta T$					
	50 $\frac{^\circ\text{C}}{\text{min}}$	0.5 $\frac{^\circ\text{C}}{\text{min}}$					
	<b>SCALE SETTING</b> Baseline Slope -0.2						




T, °C (CORRECTED FOR CHROMEL ALUMEL THERMOCOUPLES)

SAMPLE: $\text{AgNO}_3$	SIZE 3.45 mg.	ATM. Air, Static	RUN NO. SRL-23
	REF. 1.85 mg. $\text{Al}_2\text{O}_3$	T $\Delta T$	DATE 4/3/68
ORIGIN: Sadtler Research Laboratories	PROGRAM MODE Heat	SCALE 50 $\frac{\%}{\text{min}}$	OPERATOR J. Evans
	RATE 15 $\frac{\%}{\text{min}}$ , START 30 $^{\circ}\text{C}$	SETTING Baseline Slope -0.2	





<b>SAMPLE:</b> AgNO <sub>3</sub> Rerun  <b>ORIGIN:</b> Sadtler Research Laboratories	<b>SIZE</b> 3.45 mg. <b>REF.</b> 1.85. mg. <b>PROGRAM MODE</b> Heat <b>RATE</b> 15 $\frac{^\circ\text{C}}{\text{min}}$ , <b>START</b> 33 $^\circ\text{C}$	<b>ATM.</b> Air, Static <table border="1"> <tr> <td><b>T</b></td> <td><b><math>\Delta T</math></b></td> </tr> <tr> <td>50 <math>\frac{^\circ\text{C}}{\text{min}}</math></td> <td>0.5 <math>\frac{^\circ\text{C}}{\text{min}}</math></td> </tr> </table>	<b>T</b>	<b><math>\Delta T</math></b>	50 $\frac{^\circ\text{C}}{\text{min}}$	0.5 $\frac{^\circ\text{C}}{\text{min}}$	<b>RUN NO.</b> SRL-24 <b>DATE</b> 4/3/68 <b>OPERATOR</b> J. Evans
	<b>T</b>	<b><math>\Delta T</math></b>					
	50 $\frac{^\circ\text{C}}{\text{min}}$	0.5 $\frac{^\circ\text{C}}{\text{min}}$					
	<b>SCALE SETTING</b> Baseline Slope -0.2						





T, °C (CORRECTED FOR CHROMEL ALUMEL THERMOCOUPLES)


## Appendix B

### List of Thermograms (DSC)

	<u>Page</u>
Copper Sulphate Pentahydrate in Air at 5°C. /min.	46
Copper Sulphate Pentahydrate in Air at 10°C. /min.	47
Copper Sulphate Pentahydrate in Air at 20°C. /min.	48
Copper Sulphate Pentahydrate in N <sub>2</sub> at 5°C. /min.	49
Copper Sulphate Pentahydrate in N <sub>2</sub> at 10°C. /min.	50
Copper Sulphate Pentahydrate in N <sub>2</sub> at 20°C. /min.	51
Potassium Nitrate at 5°C. /min.	52
Potassium Nitrate Rerun at 5°C. /min.	53
Potassium Nitrate at 10°C. /min.	54
Potassium Nitrate Rerun at 10°C. /min.	55
Potassium Nitrate at 20°C. /min.	56
Potassium Nitrate Rerun at 20°C. /min.	57
Silver Nitrate at 5°C. /min.	58
Silver Nitrate Rerun at 5°C. /min.	59
Silver Nitrate at 10°C. /min.	60
Silver Nitrate Rerun at 10°C. /min.	61
Silver Nitrate at 20°C. /min.	62
Silver Nitrate Rerun at 20°C. /min.	63









DSC

Sample .....  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$  .....

Date ..... 1/31/68 ..... Run # ..... 258 .....


Sample Wt. .... 3.5 ..... mg Ref. Wt. .... mg


Ref. .... PANALIP .....

Scan Rate ..... 5 .....  $^{\circ}\text{C}/\text{min}$  Temp. ....  $^{\circ}\text{C}$  .....

Atmos. ....  $\text{N}_2$  ..... Flow ..... 30 ..... ml/min

Operator ..... WALLACE ..... Log Book ..... p. 59 .....





DSC

Sample .....  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$  .....

Date ..... 3/4/68 ..... Run # ..... 265 .....


Sample Wt. .... 3.5 ..... mg Ref. Wt. .... mg

Ref. .... PAN. LID .....


Scan Rate .... 10 ..... °C/min Temp. IN °C .....


Atmos. ....  $\text{N}_2$  ..... Flow .... 30 ..... ml/min


Operator WALLACE ..... Log Book ..... p. 60 .....














DSC

Sample ..... KNO<sub>3</sub> .....  
Date ..... 3/5/68 ..... Run # ..... 267 B .....  
Sample Wt. ..... 3.4 ..... mg Ref. Wt. .... mg  
Ref. .... PAN, LID .....  
Scan Rate ..... 10 ..... °C/min Temp. .... °C .....  
Atmos. .... AIR ..... Flow .... STATIC ..... ml/min  
Operator WALLACE ..... Log Book ..... p. 60 .....

CENTIGRADE

50

100

150


200


250


300

350











DSC

Sample .....  $\text{AgNO}_3$  .....  
Date ..... 2/27/68 ..... Run # ..... 262 A .....  
Sample Wt. .... 3.4 ..... mg Ref. Wt. .... mg .....  
Ref. .... PAN, LID .....  
Scan Rate ..... 5 ..... °C/min Temp. in °C .....  
Atmos. .... AIR ..... Flow STATIC ..... ml/min .....  
Operator WALLACE ..... Log Book p. 60 .....









DSC


Sample .....  $\text{AgNO}_3$  .....  
Date ..... 3 / 6 / 68 ..... Run # ..... 268 A .....  
Sample Wt. .... 3.6 ..... mg Ref. Wt. .... mg .....  
Ref. .... PAN. LID .....  
Scan Rate ..... 10 ..... °C/min Temp. IN ..... °C .....  
Atmos. .... Air ..... Flow STATIC ..... ml/min .....  
Operator WALLACE ..... Log Book ..... p. 60 .....




DSC


Sample ..... AgNO<sub>3</sub> .....  
 Date ..... 3/6/68 ..... Run # ..... 268 B .....  
 Sample Wt. ..... 3.6 ..... mg Ref. Wt. .... mg  
 Ref. .... PAN, LID .....  
 Scan Rate ..... 10 ..... °C/min Temp. .... °C  
 Atmos. .... Air ..... Flow ..... STATIC ..... ml/min  
 Operator ..... WALLACE ..... Log Book ..... p. 60 .....





DSC  
 Sample .....  $\text{AgNO}_3$  .....  
 Date ..... 3/12/68 ..... Run # ..... 272 A .....  
 Sample Wt. .... 3.5 ..... mg Ref. Wt. .... mg  
 Ref. .... PAN., LID .....  
 Scan Rate ..... 20 ..... °C/min Temp. .... °C .....  
 Atmos. .... AIR ..... Flow ..... STATIC ..... ml/min  
 Operator ..... WALLACE ..... Log Book ..... p. 60 .....





DSC

Sample .....  $\text{AgNO}_3$  .....  
Date ..... 3/12/68 ..... Run # ..... 272 B .....  
Sample Wt. .... 3.5 .....mg Ref. Wt. ....mg .....  
Ref. .... PAN, LID .....  
Scan Rate .... 20 .....°C/min Temp. .... °C .....  
Atmos. .... Air ..... Flow STATIC ..... ml/min .....  
Operator .... WALLACE ..... Log Book .... p. 60 .....

CHARTWORK

APPROVAL


EVALUATION OF ANALYTICAL STANDARDS BY DIFFERENTIAL THERMAL  
ANALYSIS AND DIFFERENTIAL SCANNING CALORIMETRY


By


J. P. Evans and K. G. Scrogam

The information in this report has been reviewed for security classification. Review of any information concerning Department of Defense or Atomic Energy Commission programs has been made by the MSFC Security Classification Officer. This report, in its entirety, has been determined to be unclassified.

This document has also been reviewed and approved for technical accuracy.

  
\_\_\_\_\_  
J. R. Nunnelley  
Chief, Chemistry Branch

  
\_\_\_\_\_  
J. E. Kingsbury  
Chief, Materials Division

  
\_\_\_\_\_  
W. R. Lucas  
Director, Propulsion and Vehicle Engineering Laboratory