

Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

HMS Critical Readiness Review:

- AED (Automated External Defibrillator)
- ISS Medical Checklist
- Med Kit Redesign
- ALSP (Advanced Life Support Pack)
- AMP (Ambulatory Medical Pack)
- IMAK (ISS Medical Accessory Pack)



Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

HMS Critical Readiness Review:

- HASP (HMS Ancillary Support Pack)
- PHS (Periodic Health Status)
- PanOptic Digital Ophthalmoscope
- IFEP



Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

Critical Readiness Review:

AED (Automated External Defibrillator)

- Moved to LAB1O5_D1, extra battery velcroed to locker door
- Resupply all on 39P (8/31/10)

Med Kit Redesign

- Ongoing project
- Scheduled for delivery 43P (April, 2011)





Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

ISS Medical Checklist

- Will be updated with Med Kit Redesign
- Consolidate hard copies on orbit to one book
- 2 copies on orbit
- Possibly go to xml for most procedures





Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

ALSP (Advanced Life Support Pack)

- Moved to LAB1O5_D1
- Resupplied (piece part) on 19A (ALSP Emergency Surg Pack),
- Resupplied IV fluids Saline, Dextrose (32P)
- Drug Pack (ULF3)
- EpiPens (3) stowed on inside D1 locker door.





Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

AMP (Ambulatory Medical Pack)

- Moved to LAB1O5_D2
- EC8+ Cartridges resupplied (36P, ULF4, 39P)
- Surgical Pack piece parts on 19A
- Dental Pack piece parts on 19A
- Pallet 4 (36P)
- Pallet 1 (ULF6)





Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

IMAK (ISS Medical Accessory Pack)

- Started sending new Teflon disposable IMAKs (2JA)
- Nomex IMAKs still on orbit
- IMAK 1014 contains extra and resupplied medication (in LAB1O5_D2), HITT has inventory





Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

HASP (HMS Ancillary Support Pack)

Moved to LAB1O5_G2

PHS (Periodic Health Status)

- PCBA Control Solution Info
 - Received comments from crew
 - BUN consistently low on red (low) controls
 - pH frequently out on red (low) and blue (high) controls
 - Updated procedure to eliminate repeat of control if out of range



Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

PanOptic Digital Ophthalmoscope

- Welch Allyn hardware with camera attachment
- Flew on 17A
- Ability to downlink digital photos and video
- Requirements being added to MRID for regular exam









Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

·IFEP

- No current issues
- Forward Work: CDCA (Common Data Collection Application) software
 - in development for FFQ and IFEP
 - will replace current software on T61p



Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

Back-Up Slides



Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

Automated External Defibrillator (AED)

Current HW Status

- AED successfully activated and checked out in June, 2008
- Requires monthly inspection (5 minute task)
- Resupply every 2 years or soonest opportunity if used
- Procedures updated to reflect AED use
- Stowed in LAB1O4_D1





Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

Automated External Defibrillator (AED)

Functionality:

- ✓ Nominally functions in automated mode— Crewmember simply powers AED on and follows audio/visual prompts
- ✓ Display screen shows cardiac rhythms, shock count, and elapsed time
- ✓ Powered by non-rechargeable battery
- ✓One AED unit on Station
- ✓Unit has one prime and one backup battery, 2 sets of electrodes






Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10



Maintenance:

Performs daily and monthly self checks and indicates on status display if service required

Data:

No real time data downlink. Data stored in AED – downloaded to a computer upon return. NO MORE COMMANDING!

Provides ECG rhythm, defibrillator incrementor, CPR countdown clock

AED



Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10



ISS Medical Checklist Planned/Recent Changes

- Incorporated:
 - AED incorporation into emergency procedures.
 - New ACLS flip book sent up on 1J, incorporated with full book revision
- Book Status
 - Hybrid Book
 - Current version flew 1J
 - PCN published after ULF2



Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

ISS Medical Checklist

Planned/Recent Changes (continued)

- Recent Changes:
 - Update with medication change (Mucinex-D substituted for Entex PSE) (30P)
 - CHeCS rack reference updates
 - Ultrasound procedure updates
 - Misc Op Nom changes

Future Changes

- Med Kit Redesign will require procedure changes
- Possible elimination of hard copy, change in ACLS Flip Book



Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

Med Kit Redesign

- CR written and approved for redesign
 - Driven by changes in medication/configuration/availability, carpuject issues
 - System Requirement Review in December, 2009
 - Tentative delivery ULF5
 - Will affect Medical Checklist and Training



Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

ALSP

Current HW Status

 ALSP Drug Pack, Airway Pack, and Emergency Surgery Pack flown on 1E

Manifest of new HW

- Packs are certified to fly on their own to facilitate resupply.
- Dextrose (0.5 L), NS (0.5 L, 1L) (32P)
- ALSP IV Admin Pack (32P)
- ALSP Airway Piece Parts (15A)
- ALSP Emergency Surgery Piece Parts- lidocaine (18S)
- ALSP Drug Pack Nitro stat (18S)



Misc Change

- None

Timeline Changes/IDRD Req

No updates



Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

AMP

- Current HW Status
 - Full AMP flown on 13A.1
- Manifest of new HW
 - AMP Pallet 1, AMP Pallet 2, and AMP Injectable Pack resupplied (30P) due to expiring medications
 - AMP Pocket 2 resupplied (ULF2)
 - AMP Pallet 4 resupplied (31P)
 - AMP Dental Piece Parts resupplied (15A)
 - AMP Surgical Piece Parts lidocaine (18S)
- Timeline Changes/IDRD Req.
 - No updates





Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

IMAK

Current HW Status

- Flown on all Russian and US flights for crew attributable medical items
- Used as catch-all for resupply

New HW or HW Changes

 CR submitted to substitute ziplock bags for IMAKs for efficiency, ease of use

Manifest of new HW

On each flight, crew request items & some resupply

Timeline Changes/IDRD Req.

 Added IMAK unpack to 'As Needed' section of Annex 4





Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

IMAK

- New HW Design
 - Size 12" x 12"
 - Teflon coated for fire retardancy
 - Disposable
 - Opaque for privacy
 - Est STS127/2JA deployment





Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

HASP

- Current HW Status
 - HASP flown on 9A
- Manifest of new HW
 - Resupplied NS on 31P (6-1 L, 3-0.5L bags) due to expiration dates
- Timeline Changes/IDRD Req
 - No updates





Medical Operations

Presenter

Date

SD/Kathy Pierpoline

02/04/10

PHS

Current HW Status

- Required hardware included in the AMP Blood Analyzer Pack, AMP Physical Exam Pack, AMP Pallet 3, and ALSP
- EC8 Cartridges resupplied inside IMAK or other container











Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

Manifest of new HW

- Many changes due to decrease in upmass
- Control Solutions Kit flying on ULF3
- New EC-8 cartridges flying on 33P, 34P
- EC-8 cartridges are sent up in Ziplock, qty = 25, no more use of Nomex container







Medical Operations

Presenter


SD/Kathy Pierpoline

Date

02/04/10

Current Issue – Carpuject Syringe Cap Failure:

- Issue
 - » Carpuject cap failed during E16
 - » Needle cap crumbled when carpuject used
- Action
 - » Procedure written providing instructions for taping needle cap for crew protection for non-emergent events (on IPV Uplink Shelf 12-0520)
 - » All carpujects are resupplied except Rocephin kit
 - » Caps are taped on ground with clear tape





Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

- Current Issue Carpuject Syringe Cap Failure:
 - » Long Term:
 - Med Kit Redesign ongoing
 - Provide most medications in manufacturer supplied carpujects
 - Supplied without needles –
 eliminates cap crumbling issues
 - Maintains manufacturer's expiration date
 - Medications not available in carpujects will be supplied in vials
 - Evaluating medication handling system





Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10



HRF Ultrasound

Current HW Status

- Power problems in Increment16
- Changed drawer location and is now operational
- HRF available for contingency and Payload use
- Moved to COL module
- Working on CMRS certification in COL module



Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

Validation of On-Orbit Methodology for the Assessment of Cardiac Function and Changes in the Circulating Volume Using Ultrasound and Braslet-M Occlusion Cuffs (SDTO 17011)





Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

Overview



- Collaborative effort between NASA and Russian FSA
 - » PI J. Michael Duncan MD (NASA), V. Bogomolov MD (IBMP)
 - » Co-Investigators I. Alferova (IBMP), S. Dulchavsky MD (Henry Ford Health System, Detroit), D. Ebert PhD (Wyle), D. Hamilton MD, PhD (Wyle), V. Matveev MD (GCTC), A. Sargsyan MD (Wyle)



Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

Overview (continued)

- Utilizes HRF Ultrasound and current Russian countermeasures (Braslet-M cuffs)
- Zero Upmass experiment
- Applicable Increments
 - » Increment 16 (5 scans)
 - » Increment 18 (5 scans)
 - » Increment 19 (requesting 2 paired scans)





Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

- Overview (continued)
 - Increment 19/20 Braslet Study
 - » Want to conduct 2 pair of scans in Incr 19/20
 - » M. Barratt, T. Kopra, K. Wakata, G. Padalka











Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

Objectives

- Develop and validate methodology to assess circulating volume changes and cardiac function in spaceflight conditions, through short term hemodynamic modification using Braslet-M occlusion cuffs
- Tests the performance of occlusion cuffs in modifying fluid shift, including physiological measurements of Braslet cuff use





Medical Operations

Presenter

SD/Kathy Pierpoline

Date


02/04/10

Benefit/Rationale

- Augments the capability of the onboard ultrasound system, including
 - » A newly validated operating mode
 - » Adds methodologies available
 - » Hones telemedicine experience
- Increased understanding of the Braslet countermeasures









Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

Resources

- Crew Time (nominal)
 - » Set-up 55 minutes
 - » Prep 5 minutes
 - » Scan 60 minutes total
 - Operator and Subject
 - Scan 25 minutes
 - Apply Braslet (Subject) wait at least 10 minutes
 - Scan 25 minutes
 - Cardiac, femoral vein, jugular scans
 - Valsalva and Mueller maneuvers
 - » Stow 20 minutes





Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

Resources (continued)

- Ground Time
 - » Before Scan
 - HRF Rack Power-up commanding (< 90 minutes before scan)
 - Marshall Activation commanding (< 90 minutes before scan)
 - PI team Preparation
 - MPSR/BME Console support
 - » After scan
 - HRF Downlink commanding (≤ 2 hours)
 - TSC provide videos to PI team
 - Marshall Deactivation command (< 2 hours)
 - BME Still photo download, MDRF for video if necessary





Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

Resources (continued)

- Hardware
 - » Ultrasound system (part of HRF)
 - » Braslet-M cuffs
 - » Ultrasound gel
 - » Straw
 - » Cameras PD 100/G1 and DCS 760/D2X (video and still image)





Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

Scheduling Constraints

- No exercise within 2 hrs of scan session.
- No Braslet use within 24 hrs of scan session
- No caffeine 12 hrs prior to scan session
- No heavy meals 4 hrs prior to scan session
- No food at all 2 hrs prior to scan session
- Acoustic Constraint 4 hours









Medical Operations

Presenter

SD/Kathy Pierpoline

Date

02/04/10

Questions?