

Paper 2038 1
Session 2.06
IEEE Aerospace Conference
Big Sky, Montana
March, 2018



## **EVA Swab Tool**

to Support Planetary Protection and Astrobiology Evaluations

NASA/Johnson Space Center/ Michelle Rucker, Drew Hood, Mary Walker
NASA-CalTech Jet Propulsion Laboratory/ Dr. Kasthuri Venkateswaran
University of Florida/ Dr. Andy Schuerger



### Issue: We have knowledge gaps!

 Weather/how microbes are released from crewed pressure systems Strategic Knowledge
Gap B5
Forward Contamination
for Mars

### Why do we care?

- Informs Mars operational concepts
- Informs architecture decisions
- Informs landing site selection decisions

- How to protect the science
- Open vs. closed life support systems?
- How close can we operate to where life may be present?



### First, We Need A Toolbox

# Microbial sampling on Earth is simple

Sterile, wetted swab

to pressure cabin




# Sampling in space is more complicated

Difficult to handle swab with bulky gloves
Wetted swab won't work in a vacuum
Sample container has to maintain

integrity when transiting from vacuum



Extravehicular Activity (EVA) Gloves





### 8-Sample Swab Kit

0.22 micron pore Microbial filter stackin bottom of each sample canister allows canisters to vent but remain sterile



6 canisters on top, 2 on bottom (not shown)

# NASA

# Tool Form, Fit, Function Tests Suit differential pressure



Lab Environment (Not Vacuum Chamber)

#### Form, Fit, Function Testing

Suited subject in a Mark III suit at differential suit pressure

### Swab Tip Environmental Testing

- Macrofoam + two other types of swab
- Bend and pull tests at temperatures from -73.3°C (-100°F) to 37.8°C (100°F)
- Testing indicated swabs would maintain integrity under expected loading, even at temperature extremes

#### **Swab Tip Effectiveness**

 Tests at JPL and JSC demonstrated dry macrofoam swab was at least as effective as standard wetted swabs

Tool Handle

We need to characterize the

suits before we send suited

crew to sample anything else

Mini Work Station



# Extravehicular Mobility Unit (EMU) Suit differential pressure



Saved cost by piggy-backing onto EMU-suited ISS crew training runs

IEEE Aerospace 2018, Session 2.06 Paper 2038 1

## NASA

# EMU Wrist Joint Swabbing 4.3 psi differential suit pressure



Initial tests: did not modify any suit cleaning or handling protocols



## Suit joints & vents are the most likely microbial escape paths

- Microbes only need 0.5 to 1.0 µm gap
- Vents can be filtered, but joints can't

## **Evaluations with EMU Suits**




Click here to play video



### External Vacuum Swab Test

- Piggy-backed onto a 6-test Orion project series
  - Two different types of suits



Orion Crew Survival System




Modified Advanced Crew Escaped System



### **Design Improvements**

- End Effector-to-Canister Seal Redesign
  - Seal retaining ring, alternate seal profile, or alternate seal type
- End Effector Locking Device
  - To ensure end effectors remain inside canisters during handling or shipping



- Swab Re-Use Prevention
  - Ratcheting device to prevent inadvertent re-use of a swab
- Additional Caddy Handle
  - For better operator control

Will be incorporated into flight design



### Next Steps Analysis, flight certification, and ISS swab

#### Culture and DNA analysis in progress





- Will feed results back to suit designers and publish test data
- Working with CASIS to identify potential commercial partners
  - Companies interested in looking for extremophiles outside ISS



ISS is seeking ideas for additional research that could be conducted after USOS increases to 4 crew (after September, 2017)



**From:** InnoCentive [mailto:no-reply@innocentive.com]

**Sent:** Thursday, August 24, 2017 4:44 PM

**Subject:** Award Announcement for NASA@work Challenge: Submit Your Research Idea to be

Conducted on ISS!: Congratulations!

Congratulations, your submission to NASA@work Challenge: 2270 - has been awarded!

Payload integration agreement in place for 2018



### Conclusions



- Ground testing validated tool effectiveness and interface with pressurized EVA suits
  - Test subjects found tool straightforward to operate
  - Design improvements identified
  - Minimized development costs by repurposing retired hardware and piggy-backing onto flight crew training
  - Additional potential uses identified: astrobiology research and micrometeoroid/orbital debris evaluation
- Baseline suit microbial data collected
  - Supports development of human mission planetary protection protocols



### Acknowledgements

The authors thank Douglas Terrier and the Science and Technology Mission Directorate, and a crossorganizational test and analysis team, including Mary Sue Bell, Alex Horvath, Justin Connolly, Bekki Bruce, Christian Castro, Dr. Aaron Regberg, Dr. Ganesh Babu Malli Mohan, and Dr. Camilla Urbaniak

## Questions?

NASA Johnson Space Center XM/Michelle Rucker Michelle.a.rucker@nasa.gov

14