




NASA's BioSentinel mission: lessons learned & what's next

Sergio R. Santa Maria ARC Space Biosciences

sergio.santamaria@nasa.gov





Recent publications

ASTROBIOLOGY

Volume 20, Number 8, 2020 © Mary Ann Liebert, Inc. DOI: 10.1089/ast.2019.2068

Introduction to BioSentinel

BioSentinel:

A Biological CubeSat for Deep Space Exploration

Sofia Massaro Tieze, 1.2 Lauren C. Liddell, 2.3 Sergio R. Santa Maria, 2.4 and Sharmila Bhattacharya 2

ASTROBIOLOGY Volume 20, Number 8, 2020 Mary Ann Liebert, Inc. DOI: 10.1089/ast.2019.2073

Research Article

BioSentinel:

Long-Term Saccharomyces cerevisiae Preservation for a Deep Space Biosensor Mission

Sergio R. Santa Maria, 1.2 Diana B. Marina, 2.3 Sofia Massaro Tieze, 2.4 Lauren C. Liddell, 2.5 and Sharmila Bhattacharya 2

ASTROBIOLOGY


Volume 21, Number 5, 2021 Mary Ann Liebert, Inc. DOI: 10.1089/ast.2020.2305

Research Article

BioSentinel:

A Biofluidic Nanosatellite Monitoring Microbial Growth and Activity in Deep Space


Michael R. Padgen,¹ Lauren C. Liddell,^{1,2} Shilpa R. Bhardwaj,^{1,3} Diana Gentry,¹ Diana Marina,^{1,4} Macarena Parra,¹ Travis Boone,^{1,5} Ming Tan,^{1,6} Lance Ellingson,^{1,5} Abraham Rademacher,^{1,5} Joshua Benton,^{1,7} Aaron Schooley,^{1,5} Aliyeh Mousavi,⁸ Charles Friedericks,¹ Robert P. Hanel,¹ Antonio J. Ricco,¹ Sharmila Bhattacharya,¹ and Sergio R. Santa Maria^{1,9}



DOI. No. 10.1109/MAES.2019.2953760

CubeSats for microbiology and astrobiology research

Luis Zea^a, Sergio R. Santa Maria^b, and Antonio J. Ricco^b
^aBioServe Space Technologies, University of Colorado, Boulder, CO, United States,
^bNASA Ames Research Center, Moffett Field, CA, United States







Proceedings

Developing Technologies for Biological Experiments in Deep Space †

Elizabeth M. Hawkins 1,2,3, Ada Kanapskyte 1,4 and Sergio R. Santa Maria 5,6,*

Feature Article:

BioSentinel: A 6U Nanosatellite for Deep-Space Biological Science

Antonio J. Ricco, Sergio R. Santa Maria, Robert P. Hanel, Sharmila Bhattacharya, BioSentinel Team, NASA Ames Research Center Radworks Group, Johnson Space Center







erspective

Space Biology Research and Biosensor Technologies: Past, Present, and Future †


Ada Kanapskyte ^{1,2}, Elizabeth M. Hawkins ^{1,3,4}, Lauren C. Liddell ^{5,6}, Shilpa R. Bhardwaj ^{5,7}, Diana Gentry ⁵ and Sergio R. Santa Maria ^{5,8,*}



What are we going to encounter beyond LEO?



What are we going to encounter beyond LEO?




- Both interplanetary & modulated by the solar cycle
- GCRs: high-energy protons and highly charged, energetic heavy particles
- GCRs not effectively shielded & can fragment

<u>Challenges</u>: SPEs – unpredictable; large doses in short time


GCRs – biology effects poorly understood (but most hazardous)

Cell nuclei showing DNA damage






The limits of life in space - as we know it - is 12.5 days on a lunar round trip or 1.2 years in LEO. As we send people further into space, we need to understand the biological risks and how they can be addressed





NASA Ames biological CubeSat missions







NASA Artemis I mission












Artemis 1: secondary payloads (6U CubeSats)














BioSentinel mission

Objectives:

- Develop a deep space tool with autonomous life support technologies
- o Study the space radiation environment beyond LEO and its effect on biology









BioSentinel: key components

Subsystem	Description	Vendor	Subsystem	Description	Vendor
4U BioSensor			LET spectrometer		
Fluidics	18 PC cards & manifolds with integrated valves, pumps, bubble traps & desiccant chambers	NASA ARC	Radiation sensor	Timepix-based sensor (TID and LET particle characterization)	NASA JSC
Optics	3 LED lights & TAOS sensor per fluidic well (16 wells per card)		Bus ADCS		
Thermal	Dedicated thermal control per fluidic card		Reaction Wheels / Star Tracker / IMU	XACT 3-axis Attitude-Control System & Star Tracker	Blue Canyon Technologies
Sensors	Temperature, RH, pressure		Sun Sensors	5 sun sensors	SolarMEMS
Bus C&DH			Bus Power		
Processor	UT700 LEON 3	Space Dynamics	EPS	TI-MSP430 FRAM	NASA ARC
Data Storage	MRAM & Flash		Batteries	18650B Li-ion	Panasonic
Bus Comm			Solar Array	4 HaWK solar arrays	MMA Design
Radio	Iris v2 transponder	JPL	Propulsion		
Antenna	Medium & low gain		Thruster	7-nozzle 3D printed system	Lightsey Space Research




BioSentinel: science instrumentation
























BioSentinel: microfluidics

























BioSentinel: microfluidics



















BioSentinel: optical detection system

Dedicated 3-color optical system at each well to track growth *via* optical density and cell metabolic activity *via* dye color changes


LEDs:


570 nm (green, measures pink)


630 nm (red, measures blue)

850 nm (infrared, measures growth)










alamarBlue turns pink when cells are metabolically active



HR repair defective cells show sensitivity to ionizing radiation



BioSentinel: optical detection system

Dedicated 3-color optical system at each well to track growth *via* optical density and cell metabolic activity *via* dye color changes


LEDs:


570 nm (green, measures pink)


630 nm (red, measures blue)

850 nm (infrared, measures growth)










alamarBlue turns pink when cells are metabolically active




HR repair defective cells show sensitivity to GCRsim ionizing radiation



BioSentinel: ISS experiment

- ISS instrument launched on SpaceX CRS-24 on Dec 21, 2021, and installed into SABL incubator
- First set of cards activated in late-Jan 2022; last one completed in mid-June
- ISS payload returned on SpX-25 in late-August
- Data used as LEO control and to validate instrumentation.

Card 8 optical readings after ~72 hours







Bio activity ~48 hours after activation (and continued for ~7 days)

Reagents (growth medium and alamarBlue dye) working properly

ISS LET spectrometer data from April 14






Ground control unit at NASA ARC




BioSentinel: ISS experiment





BioSentinel: how did we get here?





- Model organism & strain selection
- How to fly the biology? Dry/wet?
- Long-term viability & biocompatibility
- Materials, sterilization method, reagents, fluidic components...
- Multiple fluidic cards; 288 wells
- LEO control (and ground)
- Data telemetry & processing





BioSentinel: radiation environment characterization

- LET spectrometer mounted outside the BioSensor enclosure
- Reports cumulative dose, dose rate & particle distribution based on LET
- Performing nominally; unaffected by BioSensor issues

Science objectives:

- Characterize radiation environment
- Compare ISS, deep space, and ground data
- Characterize contribution from solar and GCR pa
- Characterize solar particle events (SPEs)

≺	p	ar	U	CI

rticles	
	(tra


based on LET	(mGy/day)	(6-mo mission)	
Ground	~ 0.0003	0.055 mGy	
ISS (Jan 31 – Jun 17, 2022)	0.217	-	
Free-flyer	0.306	55.78 mGy	
Free-flyer (transport modeling)	$\textbf{0.25} \pm \textbf{0.04}$	-	

Absorbed dose rate


Total absorbed dose

LET spectrometer -


Ground LET spectrometer - Dec 16



ISS LET spectrometer - April 14



Free-flyer LET spectrometer - Dec 16






BioSentinel: radiation environment characterization

- LET spectrometer mounted outside the BioSensor enclosure
- Reports cumulative dose, dose rate & particle distribution based on LET
- Performing nominally; <u>unaffected by BioSensor issues</u>

Science objectives:

- Characterize radiation environment
- Compare ISS, deep space, and ground data
- Characterize contribution from solar and GCR particles
- Characterize solar particle events (SPEs)




ased on LET	(mGy/day)	(6-mo mission)	
Ground	~ 0.0003	0.055 mGy	
ISS (Jan 31 – Jun 17, 2022)	0.217	-	
Free-flyer	0.306	55.78 mGy	
Free-flyer (transport modeling)	$\textbf{0.25} \pm \textbf{0.04}$	-	

Absorbed dose rate

Total absorbed dose


LET spectrum (in Si; 20M events); first 49 days







BioSentinel: radiation environment characterization

- CME event on March 13
- >100 MeV proton intensity enhancement
- Steady return to background levels
- BioS LET spectrometer detected event









BioSentinel: top risks & lessons learned

	Risk Description	Mitigation Approach	
1	Long-term pre-launch integration time for biology	Turn spacecraft over as late as possible	
2	DSN antenna oversubscription	Use non-DSN antennas (e.g., ESA)	
3 Flight software update after deployment		SIM tests; revert to previous version if issues arise	



BioSentinel: top risks & lessons learned

	Risk Description	Mitigation Approach	
1	Long-term pre-launch integration time for biology	Turn spacecraft over as late as possible	
2	DSN antenna oversubscription	Use non-DSN antennas (e.g., ESA)	
3 Flight software update after deployment		SIM tests; revert to previous version if issues arise	

Top Issue: multiple launch delays

How the project managed the issue?

Long-term testing (reagents, biology, etc.)

Environmental monitoring requests at KSC

Constant communication with partners

Results: Loss of personnel

Constant replanning (e.g., moving trajectories)

Additional battery recharging

Effects of long-term storage inside BioSensor enclosure

Recommendations: Constant communication between project, NASA, and SLS

Plan for such complex mission (e.g., personnel, battery charging)

Improve environmental monitoring if biology involved



BioSentinel: top risks & lessons learned

	Risk Description	Mitigation Approach	
1	Long-term pre-launch integration time for biology	Turn spacecraft over as late as possible	
DSN antenna oversubscriptionFlight software update after deployment		Use non-DSN antennas (e.g., ESA)	
		SIM tests; revert to previous version if issues arise	

Issues, lessons learned & mitigations

	What Happened?	What did we learn from it?	Proposed Mitigation Strategies	Recommended Actions
1	Lower yeast viability	Long duration effects inside sealed enclosure (个 RH)	Late integration and/or payload swap	Modify design for future missions to include late load & integration; ↑ desiccant inside enclosure
2	Unknown launch date	Loss of personnel & constant changes in trajectory	Analyze risks associated with unique mission conditions	Document for future personnel on why decisions were made
3	Late payload EVT not ideal	EVT provided many lessons learned that could not be implemented for this mission	Work around issues without design changes	Conduct early/scaled-down EVT to feed potential design changes




How are we improving LEIA based on lessons learned from BioSentinel?




LEIA lunar surface mission

LEIA: Lunar Explorer Instrument for space biology Applications










- 1. Determine cellular sensitivity to the lunar environment
- 2. Test production of antioxidant carotenoids for use in crew dietary supplementation
- 3. Engineer genetic countermeasures to test for enhanced tolerance to the lunar surface
- 4. Measure biologically relevant radiation at the south pole landing site






LEIA improvements

LEIA prototype (open configuration)








LEIA improvements

How are we preventing high humidity inside the BioSensor enclosure?




Increased desiccant

Improved seals & protected connectors



LEIA improvements



Redesigned bag tray & re-routed plumbing

Late load installation



BLEO instrumentation: future & ongoing

A flexible design that will be used in different space environments

