

AEC-NASA Tech Briefs announce new technology derived from the research and development program of the U.S. AEC or from AEC-NASA interagency efforts. They are issued to encourage commercial application. Tech Briefs are available on a subscription basis from the National Technical Information Service, Springfield, Virginia 22151. Requests for individual copies or questions relating to the Tech Brief program may be directed to the Technology Utilization Office, NASA, Code KT, Washington, D.C. 20546.

## Improved Intensifying Screen Reduces X-Ray Exposure



Spectral Emission of CaWO<sub>4</sub> and Gd <sub>2</sub>O<sub>2</sub>S:Tb

A terbium-activated gadolinium oxysulfide  $(Gd_2O_2S:Tb)$  intensifying screen increases the sensitivity of x-ray film by a factor of 2 when compared to a standard calcium tungstate  $(CaWO_4)$  screen. Optimizing the film to this new x-ray phosphor can increase film sensitivity by a factor of 20. Consequently, x-ray exposure can be reduced comparably without changing the quality of present images.

Terbium-activated gadolinium oxysulfide has over five times the intrinsic luminescense efficiency of  $CaWO_4$ , with its emissions primarily in the green spectral range near 5400 Å. Calcium tungstate, on the other hand, emits a broadband which peaks in the blue spectral range near 4200 Å. The figure compares the two spectra; the CaWO<sub>4</sub> intensity scale is expanded for clarity.

## (continued overleaf)

This document was prepared under the sponsorship of the Atomic Energy Commission and/or the National Aeronautics and Space Administration. Neither the United States Government nor any person acting on behalf of the United States Government assumes any liability resulting from the use of the information contained in this document, or warrants that the use of any information, apparatus, method, or process disclosed in this document may not infringe privately owned rights.

The spectral sensitivity of medical x-ray film (see figure) has developed over the years to match the blue emission from  $CaWO_4$ . The sensitivity range of this blue-sensitive film does not include the principle emission from  $Gd_2O_2S$ :Tb. There does exist, however, a green-sensitive film with the spectral sensitivity range shown. Using the green-sensitive film with a  $Gd_2O_2S$ :Tb screen, 1/18 of the x-ray exposure is needed to get film darkening equivalent to that obtained using a  $CaWO_4$  screen. Using the blue-sensitive film with a  $Gd_2O_2S$ :Tb screen, 1/2 the exposure is needed for equivalent darkening. An even larger exposure reduction (about 1/20 of that presently required) will be possible when a new photographic film with maximum sensitivity near 5400 Å is developed for the new screens.

The new x-ray intensifying screen may make possible new radiographic procedures where detection speed and x-ray tube power have previously been the limiting factors. It will also significantly reduce the total population exposure to harmful radiation in the United States. Note:

Requests for further information may be directed

to:

Mr. Glenn K. Ellis Technology Utilization Officer Office of Information Services U.S. Atomic Energy Commission Washington, D.C. 20545 Reference: TSP72-10232

## Patent status:

Inquiries concerning rights for commercial use of this information may be made to:

Mr. George H. Lee, Chief Chicago Patent Group U.S. Atomic Energy Commission Chicago Operations Office 9800 South Cass Avenue Argonne, Illinois 60439

> Source: R. A. Buchanan Lockheed Aircraft Corporation under contract to Atomic Energy Commission (AEC-10090)