


# NASA TECH BRIEF

## *Marshall Space Flight Center*



NASA Tech Briefs announce new technology derived from the U.S. space program. They are issued to encourage commercial application. Tech Briefs are available on a subscription basis from the National Technical Information Service, Springfield, Virginia 22151. Requests for individual copies or questions relating to the Tech Brief program may be directed to the Technology Utilization Office, NASA, Code KT, Washington, D.C. 20546.

### Remote Weighing Device



#### The problem:

Present electrobalances can determine a fraction of a milligram change in specimens weighing several hundred grams. In general, these balances will handle only specimens up to approximately 200-g and of limited physical size.

#### The solution:

A remote-weighing device has been designed and developed that is capable of measuring a weight change of 3-mg of a 500-g specimen. Because of its sensitivity, the instrument has been used for a moisture diffusion test. In general, the instrument can be used for any testing where small weight changes in large samples are required; both positive and negative changes can be measured.

#### How it's done:

The device, as shown in the figure, is constructed with a stainless steel flexure that supports a beam on one end and to which the specimen is attached on the opposite end. The flexure has four strain gages that measure bending stress and cancel tension and compression stresses. The beam is attached to the flexure so that the distance from the point of attachment to the specimen is six times the distance from the attach point to the balancing weight. Any change in the specimen weight is then multiplied by a factor of six and creates a bending stress in the flexure that is measured by the strain gage. This ratio can be varied to obtain greater resolutions or to measure heavier specimens. As the specimen holder is not part of the device, it can be custom designed to satisfy corrosion, contamination, and support features of the specimen.

(continued overleaf)

The device is compact and can accept specimens of greater physical size or total weight than electrobalances presently on the market. Possible fields of study in which the instrument might be used include corrosion studies, biological weight gains or losses, nuclear experiments, evaporation rates, and weight changes in hazardous environments.

**Patent status:**

No patent action is contemplated by NASA.

Source: J. P. Valinsky of  
McDonnell Douglas Corp.  
under contract to  
Marshall Space Flight Center  
(MFS-21556)

**Notes:**

1. Information concerning this innovation may be of interest to research or testing laboratories and instrument manufacturers.
2. Requests for further information may be directed to:  
Technology Utilization Officer  
Marshall Space Flight Center  
Code A&TS-TU  
Huntsville, Alabama 35812  
Reference: B72-10325